## Annex : Calculation of Mean and Standard Deviation

• A cholesterol control is run 20 times over 25 days yielding the following results in mg/dL:

192, 188, 190, 190, 189, 191, 188, 193, 188, 190, 191, 194, 194, 188, 192, 190, 189, 189, 191, 192.

- Using the cholesterol control results, follow the steps described below to establish QC ranges. An example is shown on the next page.
- 1. Make a table with 3 columns, labeled A, B, C.
- 2. Insert the data points on the left (column A).
- 3. Add Data in column A.
- 4. Calculate the mean: Add the measurements (sum) and divide by the number of measurements (n).

$Mean = \sum \underline{x_1 + x_2 + x_3 + \dots + x_n}$	<u>3809</u> = <b>190.5 mg/dL</b>
Ν	20

- 5. Calculate the variance and standard deviation: (see formulas below)
  - a. Subtract each data point from the mean and write in column B.
  - b. Square each value in column B and write in column C.
  - c. Add column C. Result is 71 mg/dL.
  - d. Now calculate the variance: Divide the sum in column C by n-1 which is 19. Result is **4 mg/dL**.
  - e. The variance has little value in the laboratory because the units are squared.
  - f. Now calculate the SD by taking the square root of the variance.
  - g. The result is **2 mg/dL**.

A	В	С
Data points.	—	$(-)^{2}$
X <sub>1</sub> -Xn	$\mathbf{x}_{i} - \mathbf{x}$	$(\mathbf{x}_i - \mathbf{x})$
192 mg/dL	1.5	$2.25 \text{ mg}^2/\text{dL}^2$
188 mg/dL	-2.5	$6.25 \text{ mg}^2/\text{dL}^2$
190 mg/dL	-0.5	$0.25 \text{ mg}^2/\text{dL}^2$
190 mg/dL	-0.5	$0.25 \text{ mg}^2/\text{dL}^2$
189 mg/dL	-1.5	$2.25 \text{ mg}^2/\text{dL}^2$
191 mg/dL	0.5	$0.25 \text{ mg}^2/\text{dL}^2$
188 mg/dL	-2.5	$6.25 \text{ mg}^2/\text{dL}^2$
193 mg/dL	2.5	$6.25 \text{ mg}^2/\text{dL}^2$
188 mg/dL	-2.5	$6.25 \text{ mg}^2/\text{dL}^2$
190 mg/dL	-0.5	$0.25 \text{ mg}^2/\text{dL}^2$
191 mg/dL	0.5	$0.25 \text{ mg}^2/\text{dL}^2$
194 mg/dL	3.5	$12.25 \text{ mg}^2/\text{dL}^2$
194 mg/dL	3.5	$12.25 \text{ mg}^2/\text{dL}^2$
188 mg/dL	-2.5	$6.25 \text{ mg}^2/\text{dL}^2$
192 mg/dL	1.5	$2.25 \text{ mg}^2/\text{dL}^2$
190 mg/dL	-0.5	$0.25 \text{ mg}^2/\text{dL}^2$
189 mg/dL	-1.5	$2.25 \text{ mg}^2/\text{dL}^2$
189 mg/dL	-1.5	$2.25 \text{ mg}^2/\text{dL}^2$
191 mg/dL	0.5	$0.25 \text{ mg}^2/\text{dL}^2$
192 mg/dL	1.5	$2.25 \text{ mg}^2/\text{dL}^2$

 $\sum x=3809$   $\sum = -1$   $\sum (x_i - \overline{x})^2$  Sum of Col C is 71 mg<sup>2</sup>/dL<sup>2</sup>  $SD = \sqrt{S^2} = \sqrt{\frac{\sum (X_i - \overline{X})^2}{n-1}} mg/dL$   $SD = \sqrt{S^2} = \sqrt{71}/19 = 2mg/dL$ 

## The square root returns the result to the original units.

The sum of the squared differences of each value from the mean (column C) is 71.

Notes:

a) In the calculations for variance, n-1 is used rather than n. This has been shown to reduce bias and provide a more true measure of variation. Therefore, for 20 data points, n-1 = 19. b)  $S^2$  is the variance, SD is the square root.

## **Calculate the Ranges**

The mean of these data is 190.5, and the SD is 2.

To calculate the acceptable ranges for use in quality control decisions:

1. Range for 1 SD: Subtract the SD from the mean (190.5 - 2 = 188.5)

Add the SD to the mean (190.5 + 2 = 192.5)

→ Range for 1 SD is 188.5 - 192.5.

2. Range for 2 SD: Multiply the SD by  $2 (2 \times 2 = 4)$ Add and subtract 4 from the mean (190.5)

→ Range for 2 SD is 186.5 - 194.5.

3. Range for 3 SD: Multiply the SD by 3  $(2 \times 3 = 6)$ Add and subtract 6 from the mean (190.5)

→ Range for 3 SD is 184.5 – 196.5.

Next make Levey-Jennings charts by plotting the mean and SD. See content sheets 7-4 and 7-5 for details.