Mpox disease Emergency Use Listing Procedure (EUL) for IVDs
Product: VIASURE Monkeypox virus Real Time PCR Detection Kit
EUL Number: MPXV-12649-243-00
Outcome: Accepted.

The EUL process is intended to expedite the availability of in vitro diagnostics needed in public health emergency situations and to assist interested UN procurement agencies and Member States in determining the acceptability of using specific products in the context of a Public Health Emergency of International Concern (PHEIC), based on an essential set of available quality, safety, and performance data. The EUL procedure includes the following:

- Quality Management Systems Review and Plan for Post-Market Surveillance: a desktop review of the manufacturer's Quality Management System documentation and specific manufacturing documents.
- Product Dossier Review: assessment of the documentary evidence of safety and performance. This evaluation of limited scope is to verify critical analytical and performance characteristics.

The VIASURE Monkeypox virus Real Time PCR Detection Kit, with product codes VS-MPX101L, VS-MPX101H, VS-MPX106L, VS-MPX106H, VS-MPX112L, VS-MPX112H, VS-MPX113L, VS-MPX113H, VS-MPX101, VS-MPX136, VS-MPX172, and VS-MPX196T, CE-mark regulatory version manufactured by CerTest Biotec, S.L, located at Pol. Industrial Rio Gallego II. Calle J, N°1 50840 San Mateo de Gallego, Zaragoza, Spain, was listed as eligible for WHO procurement on 9 September 2025.

Intended use:

According to the claim of intended use from CerTest Biotec, S.L., *"VIASURE Monkeypox virus Real Time PCR Detection Kit is a real-time PCR test designed for the qualitative identification of DNA from Monkeypox virus in lesion swabs from individuals suspected of Monkeypox virus infection by their healthcare professional (HCP). This test is intended to be used as an aid in the diagnosis of Monkeypox virus infection in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, amplified using real-time PCR and detected using fluorescent reporter dye probes specific for Monkeypox virus.*

The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real time PCR and in vitro diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and nucleic acid extraction system)."

Validated specimen type:

Human skin lesion swabs collected using sterile synthetic fibre swabs with plastic placed immediately into a sterile transport tube containing Universal transport medium (UTM Universal Transport Medium (Copan) or in Viral Transport Medium (-VTM- Virus Transport and Preservation Medium (Non-inactivated) (biocomma)).

Test kit contents:

- VS-MPX101L: 1 x 8-well strips, low profile
- VS-MPX101H: 1 x 8-well strips, high profile
- VS-MPX106L: 6 x 8-well strips, low profile
- VS-MPX106H: 6 x 8-well strips, high profile
- VS-MPX112L: 12 x 8-well strips, low profile
- VS-MPX112H: 12 x 8-well strips, high profile
- VS-MPX113L: 96-well plate, low profile
- VS-MPX113H: 96-well plate, high profile
- VS-MPX101: 2 x 4-well strips, Rotor-Gene
- VS-MPX136: 9 x 4-well strips, Rotor-Gene
- VS-MPX172: 18 x 4-well strips, Rotor-Gene
- VS-MPX196T: 4 tubes x 24 reactions

Items required but not provided:

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e., individual tubes, well-strips and/or microplates). Only for Tubes format.
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30° C to - 10° C and/or $\leq - 70^{\circ}$ C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 μ L, 20-200 μ L, 100-1000 μ L).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene instruments). Only for Rotor Gene format.
- Blank sample.

Validated instruments for the different workflows

1. MagDEA Dx SV Kit / magLEAD 12gC instrument (Precision System Science Co.) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad).
2. MagDEA Dx SV Kit / magLEAD 12gC instrument (Precision System Science Co.) + DTPrime real-time PCR instrument (DNA-Technology).
3. MagDEA Dx SV Kit / magLEAD 12gC instrument (Precision System Science Co.) + Rotor-Gene Q (QIAGEN).
4. MagDEA Dx SV Kit / magLEAD 12gC instrument (Precision System Science Co.) + 7500 Fast Real Time PCR System (Applied Biosystems).
5. QIAamp DNA Mini Kit (QIAGEN) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad).

6. QIAamp DNA Mini Kit (QIAGEN) + DTPPrime real-time PCR instrument (DNA-Technology).
7. QIAamp DNA Mini Kit (QIAGEN) + Rotor-Gene Q (QIAGEN).
8. QIAamp DNA Mini Kit (QIAGEN) + 7500 Fast Real Time PCR System (Applied Biosystems).
9. MagCore Viral Nucleic Acid Extraction Kit, using the MagCore automated nucleic acid extractors (RBC Bioscience) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad).

Validated thermocycler software:

- CFX Maestro 2.3 software, versions 5.0.021.0616 and 5.3.022.1030, in Bio-Rad instruments.
- Realtime_PCR software, v7.9, specifically version 7.9.9.25 (last edition), in DNA-Technology thermocycler.
- Rotor-Gene Q Series Software 2.3.5 (Build 1), used in Rotor-Gene Q (QIAGEN) instrument.
- 7500 Software v2.0.6, applied in the Applied Biosystems thermocycler.

Validated automated extraction systems software:

- 1.0.3 firmware version 1.2A, belonging to the magLEAD 12gC instrument (Precision System Science Co.).
- L0-0205-D008, included in the MagCore automated nucleic acid extractors (RBC Bioscience).

Storage:

The test kit must be stored at 2-40 °C.

Shelf-life upon manufacture:

The shelf-life is currently assigned 24 months dating.

Warnings/limitations:

Please refer to the Instructions for Use attached to this public report.

Product dossier assessment

CerTest Biotec, S.L. submitted the product dossier for the VIASURE Monkeypox virus Real Time PCR Detection Kit in alignment with the Instructions and requirements for Emergency Use Listing (EUL) Submission: In vitro diagnostics detecting Monkeypox virus nucleic acid (PQDx_457). The WHO reviewed the information provided in the dossier.

The risk-benefit assessment conclusion was acceptable.

Quality Management Systems Review

To establish eligibility for WHO procurement, CerTest Biotec, S.L. was asked to provide up-to-date information about the status of its quality management system.

Based on the WHO's review of the submitted quality management system documentation, CerTest Biotec, S.L. provided sufficient information to fulfil the requirements described in the Instructions and requirements for EUL Submission: In vitro diagnostics detecting Monkeypox virus nucleic acid (PQDx_457).

The conclusion of the quality management system assessment was acceptable.

Plan for Post-Market Surveillance

Post-market surveillance, including monitoring all customer feedback, detecting and acting on adverse events, product problems, non-conforming goods and processes is a critical component of minimising the potential harm of an IVD listed for emergency use.

The following post-EUL activities are required to maintain the EUL status:

1. Notification to WHO of any planned changes to a prequalified product, in accordance with *“Reportable changes to WHO prequalified and emergency use listed in vitro diagnostics”*¹; and
2. Post-market surveillance activities, in accordance with *“WHO guidance on post-market surveillance of in vitro diagnostics” (ISBN 978 92 4 150921 3)*².

CerTest Biotec, S.L. is also required to submit an annual report summarising sales data and all complaints. Certain complaints and changes to the product must be notified immediately to WHO, as per the above-mentioned documents. The sales data will serve as denominator data to guide the frequency of re-inspection.

The manufacturer has committed to ensuring that post-emergency use listing safety, quality, and performance monitoring activities are in place, which are in accordance with WHO guidance on post-market surveillance of in vitro diagnostics.

Scope and duration of procurement eligibility

The VIASURE Monkeypox virus Real Time PCR Detection Kit, with product codes VS-MPX101L, VS-MPX101H, VS-MPX106L, VS-MPX106H, VS-MPX112L, VS-MPX112H, VS-MPX113L, VS-MPX113H, VS-MPX101, VS-MPX136, VS-MPX172, and VS-MPX196T, manufactured by CerTest Biotec, S.L., is eligible for WHO procurement for 12 months from the day of listing. The assay detects the monkeypox virus nucleic acid, including Clade I and Clade II. This listing does not infer that the product meets WHO prequalification requirements and does not mean that the product is listed as WHO-prequalified. As part of the ongoing requirements for listing

¹ <https://iris.who.int/handle/10665/381373>

² <https://iris.who.int/handle/10665/337551>

as eligible for WHO procurement, CerTest Biotec, S.L. must engage in post-market surveillance activities to ensure that the product continues to meet safety, quality, and performance requirements. CerTest Biotec, S.L. is required to notify WHO of any complaints, including adverse events related to the use of the product, within 10 days.

WHO reserves the right to rescind eligibility for WHO procurement if additional information on the safety, quality, and performance during post-market surveillance activities and if new data becomes available to WHO that changes the risk-benefit balance.

Labelling


- 1. Labels**
- 2. Instructions for use**

1.0 Labels

1.1 Low Profile

1.1.1 VS-MPX101L Well Strip (low profile)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , Low Profile

certest.es/
viasure/labeling

REF VS-MPX101L

LOT MPX101L-XXX

yy **yy** **mm**

1 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
1 x Tear-off 8-cap strips

(01)08435440240556
(17)991231
(10)MPX101L000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit


CERTEST BIOTEC S.L Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX101L v.00

1.1.2 VS-MPX106L Well Strip (low profile)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , Low Profile

certest.es/
viasure/labeling

REF VS-MPX106L

LOT MPX106L-XXX

yy **yy** **mm**

6 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
6 x Tear-off 8-cap strips


(01)08435440240471
(17)991231
(10)MPX106L000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX106L v.00

1.1.3 VS-MPX112L Well Strip (low profile)**VIASURE *Monkeypox virus***
Real Time PCR Detection Kit , Low Profile**REF** VS-MPX112L**LOT** MPX112L-XXX

yyyy-mm


12 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
12 x Tear-off 8-cap strips

**UDI**

(01)08435440240495
(17)991231
(10)MPX112L000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit **CERTEST BIOTEC S.L.** Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX112L v.00

1.1.4 VS-MPX1SL Well Strip (low profile)**VIASURE *Monkeypox virus***
Real Time PCR Detection Kit , Low Profile**REF** VS-MPX112L**LOT** MPX112L-XXX

yyyy-mm

12 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
12 x Tear-off 8-cap strips


**UDI**

(01)08435440240495
(17)991231
(10)MPX112L000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit **CERTEST BIOTEC S.L.** Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX112L v.00

1.1.5 VS-MPX1SL Well Strip (low profile) in aluminium pouch


1.1.6 VS-MPX113L Well Plate (low profile)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , Low Profile

certest.es/viasure/labeling

REF VS-MPX113L **LOT** MPX113L-XXX  yyyy-mm

1 x MPX 96-well plate
 1 x Rehydration Buffer (1,8 mL)
 1 x Positive Control
 1 x Negative Control (1 mL)
 1 x Water RNase/DNase free (1 mL)
 12 x Tear-off 8-cap strips


 (01)08435440240518
 (17)991231
 (10)MPX113L000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit


CERTEST BIOTEC S.L Pol. Industrial Río Gállego II, Calle J, Nº 1,
 50840 San Mateo de Gállego, Zaragoza (SPAIN)


VS-MPX113L v.00

1.1.7 VS-MPX1PL Plate (low profile) in aluminium pouch


1.1.8 Vial labels

1.1.9 VS-MPX101H Well Strip (high profile)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , High Profile

1 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
1 x Tear-off 8-cap strips

UDI

(01)08435440240563
(17)991231
(10)MPX101H000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX101H v.00


1.1.10 VS-MPX106H Well Strip (high profile)

VIASURE *Monkeypox virus*

Real Time PCR Detection Kit , High Profile

48

[certest.es/
viasure/labeling](https://certest.es/viasure/labeling)**REF** VS-MPX106H**LOT** MPX106H-XXX

yyyy-mm

6 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
6 x Tear-off 8-cap strips

(01)08435440240488
(17)991231
(10)MPX106H000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX106H v.00


1.1.11 VS-MPX112H Well Strip (high profile)

VIASURE *Monkeypox virus*

Real Time PCR Detection Kit , High Profile

96

[certest.es/
viasure/labeling](https://certest.es/viasure/labeling)**REF** VS-MPX112H**LOT** MPX112H-XXX

yyyy-mm

12 x MPX 8-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
12 x Tear-off 8-cap strips

(01)08435440240501
(17)991231
(10)MPX112H000


VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX112H v.00


1.1.12 VS-MPX1SH Well Strip (high profile)

1.1.12 VS-MPX113H Well Strip (high profile)

VIASURE *Monkeypox virus*

Real Time PCR Detection Kit , High Profile

[certest.es/
viasure/labeling](http://certest.es/viasure/labeling)

REF VS-MPX113H

LOT MPX113H-XXX

yyyy-mm

1 x MPX 96-well plate
 1 x Rehydration Buffer (1,8 mL)
 1 x Positive Control
 1 x Negative Control (1 mL)
 1 x Water RNase/DNase free (1 mL)
 12 x Tear-off 8-cap strips

(01)08435440240525
 (17)991231
 (10)MPX113H000


VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
 50840 San Mateo de Gállego, Zaragoza (SPAIN)

VS-MPX113H v.00

1.1.13 VS-MPX1PH Well Strip (high profile)

Vial labels

**Rehydration
Buffer**

REF VS-RB02
LOT RB02-xxx
VOL 1,8 mL
CERTEST BIOTEC S.L.

yyyy-mm
2°C 40°C
VS-RB02 v.03

**Water
RNAse/DNAse free**

REF VS-H2O
LOT H2O-xxx
VOL 1 mL
CERTEST BIOTEC S.L.

yyyy-mm
2°C 40°C
VS-H2O v.03

Negative Control


REF VS-NC1
LOT NC1-xxx
VOL 1 mL
CERTEST BIOTEC S.L.

yyyy-mm
2°C 40°C
VS-NC1 v.03

**MPX
Positive Control**

REF VS-MPX1C
LOT MPX1C-xxx
yyyy-mm

20
2°C 40°C
VS-MPX1C v.00

1.1.14 VS-MPX101 Well Strip (Rotor Gene)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , Rotor-Gene

CE IVD 8 2°C to 40°C [certest.es/
viasure/labeling](http://certest.es/viasure/labeling)

REF VS-MPX101 LOT MPX101-XXX yyyy-mm

2 x MPX 4-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNAse free (1 mL)
2 x Tear-off 4-cap strips


UDI
(01)08435440240570
(17)991231
(10)MPX101000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L Pol. Industrial Río Gállego II, Calle J, N° 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)


VS-MPX101 v.00

1.1.15 Well Strip (Rotor-Gene) in aluminium pouch

1.1.16 VS-MPX136 Well Strip (Rotor Gene)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , Rotor-Gene

[certest.es/
viasure/labeling](http://certest.es/viasure/labeling)

REF VS-MPX172

LOT MPX172-XXX

Exp yyyy-mm

18 x MPX 4-well strips
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)
18 x Tear-off 4-cap strips


(01)08435440240549
(17)991231
(10)MPX172000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit


 CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)


VS-MPX172 v.00

1.1.17 VS-MPX1RG Well Strip (Rotor-Gene) in aluminium pouch

1.1.18 Vial labels

1.1.19 VS-MPX196T Tube (Reaction-Mix Tube)

VIASURE *Monkeypox virus* Real Time PCR Detection Kit , Tube Format

CE IVD certest.es/viasure/labeling

REF VS-MPX196T LOT MPX196T-XXX yyyy-mm

4 x MPX Reaction-Mix tubes
1 x Rehydration Buffer (1,8 mL)
1 x Positive Control
1 x Negative Control (1 mL)
1 x Water RNase/DNase free (1 mL)

UDI
(01)08435440240587
(17)991231
(10)MPX196T000

VIASURE *Monkeypox virus* Real Time PCR Detection Kit

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)



VS-MPX196T v.00

1.1.20 Vial labels

**MPX
Reaction-Mix tube**

REF VS-MPX1TB
LOT MPX1TB-xxx

yyy-mm

**Rehydration
Buffer**

REF VS-RB02 yyy-mm

LOT RB02-xxx

VOL 1,8 mL

CERTEST BIOTEC S.L.


**Water
RNase/DNase free**

REF VS-H2O yyy-mm

LOT H2O-xxx

VOL 1 mL

CERTEST BIOTEC S.L.


**Negative Control**

REF VS-NC1 yyy-mm

LOT NC1-xxx

VOL 1 mL



CERTEST BIOTEC S.L.

**MPX
Positive Control**

REF VS-MPX1C

LOT MPX1C-xxx

yyy-mm

2.0 Instructions for Use³

³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MPX101L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MPX101H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MPX106L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MPX106H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MPX112L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MPX112H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MPX113L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MPX113H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MPX101
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MPX136
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MPX172

Table A 1. References for Open and Rotor-Gene format with internal control products. / *Referencias para productos Open y Rotor-Gene Format con control interno.*

TUBE FORMAT WITH INTERNAL CONTROL (SEE ANNEX 2) / FORMATO TUBO CON CONTROL INTERNO (VER ANEXO 2)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MPX196T

Table A 2. References for Tube format with internal control products. / *Referencias para productos formato Tubo con control interno.*

Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1. Intended use	6
2. Summary and Explanation.....	6
3. Principle of the procedure	7
4. Reagents provided	8
5. Reagents and equipment to be supplied by the user	8
6. Transport, storage and use conditions	10
7. Precautions for users	11
8. Test procedure	13
8.1. Specimen collection, transport and storage	13
8.2. DNA extraction	15
9. Result interpretation	16
9.1. References with internal control (references in Annexes 1 and 2)	16
10. Limitations of the test.....	18
11. Quality control	20
12. Analytical performance characteristics	20
12.1. Analytical linearity	20
12.2. Analytical sensitivity. Limit of Detection (LoD).....	21
12.3. Measuring range.....	22
12.4. Accuracy	22
12.4.1. Trueness (Veracity)	22
12.4.2. Precision.....	23
12.4.3. Accuracy	25
12.5. Analytical specificity and reactivity	26
12.5.1. Analytical specificity.....	26
12.5.2. Analytical reactivity	39
12.6. Metrological traceability	40
13. Clinical performance characteristics	40
ANNEX 1	43
A1.1 Principle of the procedure	43

A1.2 Reagents provided	43
A1.3 Test procedure	45
A1.3.1 Lyophilised positive control	45
A1.3.2 PCR protocol	45
ANNEX 2	47
A2.1 Principle of the procedure	47
A2.2 Reagents provided	47
A2.3 Test procedure	48
A2.3.1 Lyophilised positive control	48
A2.3.2 Lyophilised Reaction-Mix tube	48
A2.3.3 PCR protocol	48

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ENGLISH

1. Intended use

VIASURE *Monkeypox virus* Real Time PCR Detection Kit is a real-time PCR test designed for the qualitative identification of DNA from Monkeypox virus in lesion swabs from individuals suspected of Monkeypox virus infection by their healthcare professional (HCP). This test is intended to be used as an aid in the diagnosis of Monkeypox virus infection in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, amplified using real-time PCR and detected using fluorescent reporter dye probes specific for Monkeypox virus.

The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and nucleic acid extraction system).

2. Summary and Explanation

Monkeypox virus (MPXV) is a linear, double-stranded DNA virus which belongs to the genus *orthopoxvirus* in the family of *Poxviridae*, known for being the responsible of causing the monkeypox disease (mpox) (Chauhan et al., 2023; Lu et al., 2023; Qudus et al., 2023; WHO, 2024). This disease is a zoonotic viral condition with symptoms very similar to the ones caused by other *orthopoxvirus* (such as variola virus or cowpox virus), which includes eruptive skin lesions, rash, fever, fatigue and lymphadenopathy (Chauhan et al., 2023; Lu et al., 2023; Qudus et al., 2023). Mpox is described as self-limiting and less severe than the close condition smallpox, but may become complicated depending on the viral strain and the individual immune status; for example, some population risk groups are children, pregnant women or immunocompromised individuals (Lu et al., 2023; Qudus et al., 2023; WHO, 2024). Two clades have been described so far, the West African (currently known as clade II) and the Congo Basin (Central African, currently known as clade I) (Lu et al., 2023; Qudus et al., 2023; WHO, 2024). These are very different in relation to their epidemiologic, clinical characteristics and geographical localisation (Qudus et al., 2023).

Person-to-person MPXV transmission may occur directly by contact with skin lesions, body fluids, or respiratory secretions from infected individuals, or indirectly by contact with infected objects or linens (Lu et al., 2023; Qudus et al., 2023; WHO, 2024). Then, this virus enters the healthy individual's body through mucous membranes or broken skin, and subsequently, it spreads via tissue-resident immune cells and draining lymph nodes (Lu et al., 2023; Qudus et al., 2023), giving rise to an incubation period of 6-13 days usually, although possible to range from 5-21 days (WHO, 2024).

Since the first mpox cases in humans described in 1970s in African countries, infections have been steadily increased throughout these countries until its spreading to non-endemic regions, including USA and European countries (Chauhan et al., 2023; Qudus et al., 2023). In 2022, multi-country outbreaks were

described for the first time, due to human-to-human transmission through sexual contact (Chauhan et al., 2023). According to the WHO and as reproduce the ECDC, since the beginning of mpox surveillance in 2022 and up to July 2024, more than 100,000 confirmed cases due to clades I and II, including more than 200 deaths, were reported in more than 120 countries globally (ECDC, 2024; WHO, 2024a, 2024c). In 14th August of 2024, the WHO declared again the mpox as Public Health Emergency of International Concern (PHEIC) based on the detection of a new subclade (clade Ib) in Easter Democratic Republic of the Congo, which it is spreading rapidly and being detected in the neighbouring countries that has not previously reported mpox, and has the potential to extend outside Africa and beyond (WHO, 2024c).

The similarities in clinical symptoms between diseases caused by *orthopoxvirus* and other viral infections such as varicella-zoster virus make the diagnosis challenging (Chauhan et al., 2023). Real-time PCR assays are required for controlling outbreaks and reducing human-to-human transmission (Chauhan et al., 2023).

NOTE: to ease the comprehension, we will use the following abbreviations throughout the entire document:

- Mpox: Monkeypox disease (former nomenclature).
- MPXV: Monkeypox virus.

3. Principle of the procedure

VIASURE *Monkeypox virus* Real Time PCR Detection Kit is a non-automatized real-time PCR test designed for the qualitative detection of DNA from Monkeypox virus (MPXV), including Clades I and II, in lesion swabs. After DNA isolation, the identification of MPXV is performed by the amplification of a conserved region of the *G2R* and *F3L* genes, using specific primers and a fluorescent-labelled probes. This test may be performed using several validated workflows (see section 5), which allow for greater versatility and adaptability of use.

VIASURE *Monkeypox virus* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *Monkeypox virus* Real Time PCR Detection Kit contains all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilised¹ format. The assay uses a human housekeeping gene as an **Endogenous Internal Control (EIC)** (haemoglobin subunit beta gene, *HBB* gene) to follow-up the integrity of the sample, to monitor the extraction process and/or discard the inhibition of the polymerase activity. Human housekeeping genes are involved in basic cell maintenance

¹ Please note that both "stabilised", and "lyophilised" terms are used indistinguishable and as synonyms throughout the entire document.

and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels. If there is an improper collection, storage, and/or processing of the specimens, this control may become negative.

Moreover, this kit contains two additional controls: Positive Control (PC) and Negative Control (NC). The PC is a non-infectious synthetic DNA positive for the target under study to be used for the monitoring of the PCR amplification in each execution. On the other hand, the NC is used to detect unspecific amplifications during the amplification step. Please consult section 11. Quality control for further information.

4. Reagents provided

VIASURE *Monkeypox virus* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for Open and Rotor-gene format products with Endogenous Internal Control (EIC) and in Annex 2 for Tube format products with Endogenous Internal Control (EIC).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e., individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or $\leq -70^{\circ}\text{C}$.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 μL , 20-200 μL , 100-1000 μL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1).
- Blank sample (See section 8.2 DNA extraction).

VIASURE *Monkeypox virus* Real Time PCR Detection Kit was tested with different workflows, although nine main workflows were selected to carry out the analytical and/or clinical evaluation and validation of the device. These instruments are considered as validated for the use of the VIASURE *Monkeypox virus* Real Time PCR Detection Kit:

	Instruments
1	MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad) ²
2	MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + DTPRime real-time PCR instrument (DNA-Technology)
3	MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + Rotor-Gene® Q (QIAGEN)
4	MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + 7500 Fast Real Time PCR System (Applied Biosystems™)
5	QIAamp DNA Mini Kit (QIAGEN) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad) ²
6	QIAamp DNA Mini Kit (QIAGEN) + DTPRime real-time PCR instrument (DNA-Technology)
7	QIAamp DNA Mini Kit (QIAGEN) + Rotor-Gene® Q (QIAGEN)
8	QIAamp DNA Mini Kit (QIAGEN) + 7500 Fast Real Time PCR System (Applied Biosystems™)
9	MagCore® Viral Nucleic Acid Extraction Kit, using the MagCore® automated nucleic acid extractors (RBC Bioscience) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad) ²

Table 1. Extraction-amplification methods (workflows) validated using VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

The software and/or firmware versions used during the performance validation of the VIASURE kit were the following ones:

- Thermocyclers:
 - CFX Maestro 2.3 software, versions 5.0.021.0616 and 5.3.022.1030, in Bio-Rad instruments.
 - Realtime_PCR software, v7.9, specifically version 7.9.9.25 (last edition), in DNA-Technology thermocycler.
 - Rotor-Gene Q Series Software 2.3.5 (Build 1), used in Rotor-Gene® Q (QIAGEN) instrument.
 - 7500 Software v2.0.6, applied in the Applied Biosystems™ thermocycler.
- Automated extraction systems:
 - 1.0.3 firmware version 1.2A, belonging to the magLEAD® 12gC instrument (Precision System Science Co.).
 - L0-0205-D008, included in the MagCore® automated nucleic acid extractors (RBC Bioscience).

To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

² Both, the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) and the CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad) are equivalent and have the same operating characteristics. The CFX96™ Real-Time PCR Detection System (Bio-Rad) thermocycler also has similar characteristics and can be considered equivalent.

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Samples selected for the evaluation of the device were lesion swabs collected with sterile synthetic fibre swabs with plastic placed immediately into a sterile transport tube containing Universal transport medium (UTM® Universal Transport Medium™ (Copan) (workflows 1-8) or in Viral Transport Medium (-VTM- Virus Transport and Preservation Medium (Non-inactivated) (biocomma ®)) (workflow 9).

6. Transport, storage and use conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the kit label.
- Avoid vibrations during transport to prevent liquid leakage.
- Once the positive control has been re-suspended, store it at -20°C until the expiration date stated on the kit label. It is recommended to separate it into aliquots (for instance, 4 aliquots of 25 µL each) to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- For Tube format kits: Once the *Monkeypox virus* Reaction-Mix tube has been reconstituted, it may be kept at 25±5°C or 4±2°C for up to 4 hours. For a longer period, it is recommended store at -20°C until the expiration date stated on the kit label. Reaction-Mix in tube format has been validated as still being stable after 6 freeze-thaw cycles. Once the tube has been opened/or rehydrated, it must be stored to keep the vial away from light.
- For Open format kits (strips and plates): For proper use of the kit, it is recommended not to rehydrate additional unused wells. If all the provided wells in an 8-well strip or 96-well plate are not used, the required number of reaction wells can be easily cut from 8-well strips or 96-well plates. If all wells are not used, it is recommended that unused wells are stored in the provided foil bags with the silica gel inside to prevent contact with light until use. The open bag with not rehydrated wells at 2-40°C is stable until the shelf life stated on the kit label.
- For strips and plates: Once the required wells are rehydrated, the addition of nucleic acid and performance of the Real-Time PCR must be carried out immediately.
- After nucleic acid extraction from the clinical sample, the eluted could be stored for up to 24 hours at 4-25°C or six months at -20°C prior to addition to the Real-Time PCR assay.

The following table summarises the conditions for transport, storage and use of the overall kit and each component:

Component	Transport Conditions	Storage Conditions	In-use conditions
Entire VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit.	2-40°C during the shelf life stated in the kit label.	Before use: 2-40°C during the shelf life stated in the kit label.	* See in-use conditions of each component.
<i>Monkeypox virus</i> well-strip or well plate (See Annex 1).		Before use: 2-40°C during the shelf life stated in the kit label. Non-rehydrated wells in an open pouch with the silica gel: 2-40°C during the shelf life stated in the kit label.	Room temperature (18-31°C)*.
<i>Monkeypox virus</i> Reaction-Mix tube (see Annex 2).		Before use: 2-40°C during the shelf life stated in the kit label. Rehydrated tube: 25±5°C or 4±2°C for up to 4 hours. For a longer period, store directly at -20°C during the shelf life stated in the kit label.	Immediate use at Room temperature (18-31°C)*.
<i>Monkeypox virus</i> Positive Control vial		Before use: 2-40°C during the shelf life stated in the kit label. Rehydrated control: -20°C during the shelf life stated in the kit label.	Immediate use at Room temperature (18-31°C)*.
Rehydration Buffer		Before use: 2-40°C during the shelf life stated in the kit label. Once opened: 2-40°C during the shelf life stated in the kit label.	Room temperature (18-31°C)*.
Negative Control vial		Before use: 2-40°C during the shelf life stated in the kit label. Once opened: 2-40°C during the shelf life stated in the kit label.	Room temperature (18-31°C)*.
Water RNase/DNase free vial		Before use: 2-40°C during the shelf life stated in the kit label. Once opened: 2-40°C during the shelf life stated in the kit label.	Room temperature (18-31°C)*.

Table 2. Summary of the conditions for transport, storage and use of the VIASURE *Monkeypox virus* Real Time PCR Detection Kit and each component.

* Room temperature = 18°C - 31°C. This range is established after performing tests where the ambient temperature conditions (at 30°C±2) and/or the exposure time of the reaction mixture to these temperatures (immediate use vs. 15 min) were modified before and after rehydration, both cases with the protective seal open. PCR amplification was performed at the ambient temperature established by the thermocycler manufacturers, corresponding to a range of 18°C-30°C.

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of Real-Time PCR and *in vitro* diagnostic procedures (including training on the Real-Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Instructions for use must be read carefully before using the VIASURE *Monkeypox virus* Real Time PCR Detection Kit. Do not perform the PCR assay until the information about procedures, safety precautions and limitations described in the Instructions for use have been understood.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.

- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Once the protective aluminium seal has been removed from the plates or strips, or the tubes have been opened, the end user must use the plates, strips or tubes in the shortest time possible to protect the master mix from sunlight.
- Please note: In case the number of reactions required/needed is less than supplied on one strip or plate, before removing the protective seal, carefully cut out the required wells and store the rest inside the bag with the desiccant.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-MPX113L, VS-MPX113H, VS-MPX136, and VS-MPX172). Remove any excess air in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- To avoid label deterioration, do not use the product near solvents.
- Do not mix reagents from different pouches and/or kits and/or lots and/or another supplier.
- Protect reagents against humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-MPX101, VS-MPX136, and VS-MPX172 (compatible for use with Qiagen/Corbett Rotor-Gene® instruments), use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the Reaction-Mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Once the pouch of Positive Control or Reaction-Mixture in Tube format has been opened, rehydrate the vials immediately, use them in the PCR reaction or store them properly. It is not recommended to leave the vial open without resuspending, as this may cause deterioration of the reagent.
- Once the protective aluminium seal has been removed from the plate or strips, use them immediately for the Real-Time PCR (rehydrate the wells, add the samples, and place the strips/plate into the thermocycler). It is not recommended to leave the strips or plate open without resuspending. Likewise, once the required wells are rehydrated, the addition of nucleic acid and performance of the Real-Time PCR must be carried out immediately.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Monkeypox virus* Real Time PCR Detection Kit and any additional reagents or equipment required for testing are not contaminated. Always avoid microbial and

ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.

- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment, because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration. A statement declaring no requirement of Material Safety Data Sheet could be requested to Certest Biotec.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions, and procedures.
- Make sure that the configuration and conditions of the settings on the thermocycler is done following the instructions in the section 'PCR protocol' of Annexes 1 and 2 (DNA/RNA protocol, number of cycles, temperatures etc.).
- The certificate of analysis is not included with the device; however, it could be requested to Certest Biotec in case of need.

8. Test procedure

Please see Annex 1 for Open and Rotor-gene format products with Endogenous Internal Control Test Procedure and Annex 2 for Tube format products with Endogenous Internal Control Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *Monkeypox virus* Real Time PCR Detection Kit was tested using lesion swabs, including vesicles, pustules or scabs, collected with sterile synthetic fibre swabs with plastic placed immediately into

a sterile transport tube containing Universal transport medium (UTM® Universal Transport Medium™ (Copan)) or Viral Transport Medium (-VTM- Virus Transport and Preservation Medium (Non-inactivated) (biocomma®)), depending on the sample. Other types of samples must be validated by the user.

The equivalence between the use of Copan's UTM and biocomma's VTM for Monkeypox virus sample collection lies primarily in their functional roles as viral transport media. Both are designed to stabilise viral samples for diagnostic testing, particularly during transport from the collection site to the laboratory. UTM/VTM media generally consists of a balanced salt solution and antimicrobial agents to maintain viral integrity while preventing bacterial or fungal growth. The Copan's UTM has been widely validated and used for a range of viral pathogens, including MPXV, as it preserves the viability of viruses for molecular testing. Similarly, biocomma's VTM performs a comparable function, ensuring the preservation of viral nucleic acid, which is crucial for PCR-based diagnostic tests. Both media are designed to meet WHO and CDC guidelines for viral sample collection and transport.

Collection, storage, and transport of specimens should be maintained per the conditions validated by the user. Overall, clinical samples should be collected as soon as possible, during the acute phase of illness or at any time during the clinical course, and before a therapy begin, if possible. Specimens should be obtained by avoiding biological secretions that are not of interest, and subsequently collected in sterile, clean containers with or without transport media (depending on sample type). All specimens must be labelled with the name and identification number of the person from whom the specimen was collected, the source of the specimen, and the date and time it was collected. After collection, specimens should be placed in a biohazard bag and should be transported and processed as soon as possible to guarantee the quality of the test. Specimens submitted for molecular testing must be stored in controlled conditions so that nucleic acids do not degrade during storage. It is recommended to use fresh specimens for the test. For both transport (if necessary) and storage of clinical samples until processing, it should be considered the temperature at which they have to be stored and for how long (also following local and national regulations for transport of pathogenic material). In the case of lesion swab specimens, they should be maintained at 2 to 8°C for up to 7 days³, or frozen at -20°C or lower, ideally at -80°C, for conservation if time until processing exceeds 7 days³ (with a maximum of 2 freeze/thaw cycles).

The clinical specimens must be collected, transported and stored according to appropriate laboratory guidelines and/or laboratory policy manuals. As example, refer to the *Guidelines for Collecting and Handling Specimens for Mpox Testing*. (<https://www.cdc.gov/mpox/hcp/diagnosis-testing/collecting-specimens.html>).

Please note: Specimen collection, transport, and storage conditions indicated above are suggested based on the recommendations for lesion swabs (collected in transport medium (UTM or VTM)) intended to be used for nucleic acid detection, which appear in the referenced CDC guidelines. Nevertheless, we

³ *Guidelines for Collecting and Handling Specimens for Mpox Testing*. (<https://www.cdc.gov/mpox/hcp/diagnosis-testing/collecting-specimens.html>).

recommend following laboratory guidelines, and/or microbiology laboratory policy manual for proper transport and preservation of samples^{3,4}.

An internal specimen stability study was conducted with VIASURE *Monkeypox virus* Real Time PCR Detection Kit, using contrived lesion swabs (negative skin lesion swabs in UTM® Universal Transport Medium™ (Copan) spiked with MXPV) at two different concentrations (2xLoD and 5xLoD). Samples were examined at 24, 48 and 168 hours after being stored at 4°C±2°C, and at 24 hours after being maintained at 25°C±2°C. Results showed a good performance of samples stored at all conditions tested, meeting the acceptance criteria initially stated and complying with the sample storage criteria recommended by the CDC guidelines.

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

Please note: Positive and Negative control must not be extracted. However, to monitor or control the extraction process and discard possible contaminations, a blank sample (not provided in the kit) could be extracted following the same procedure used for clinical sample. Blank sample should consist of the same matrix (real or simulated) as clinical samples, previously characterized as negative for all targets.

VIASURE *Monkeypox virus* Real Time PCR Detection Kit was validated with lesion swabs collected in UTM® Universal Transport Medium™ (Copan) extracted with MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) or QIAamp DNA Mini Kit (QIAGEN), and analysed on CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad), DTPPrime real-time PCR instrument (DNA-Technology), Rotor®-Gene Q (QIAGEN) or 7500 Fast Real Time PCR System (Applied Biosystems™). In addition to these workflows, VIASURE *Monkeypox virus* Real Time PCR Detection Kit was validated with lesion swabs collected in VTM Virus Transport and Preservation Medium (Non-inactivated) (biocomma ®), extracted with MagCore® Viral Nucleic Acid Extraction Kit / MagCore® automated nucleic acid extractors (RBC Bioscience) and analysed on CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad).

In case of using a different workflow than those mentioned above, the validation of the new workflow shall be performed by the end-user following the manufacturer's instructions for use.

⁴ [Surveillance, case investigation and contact tracing for Monkeypox: Interim guidance \(who.int\)](https://www.who.int/publications-detail/surveillance-case-investigation-and-contact-tracing-for-monkeypox-interim-guidance)

9. Result interpretation

9.1. References with internal control (references in Annexes 1 and 2)

The result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Endogenous Internal Control (EIC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. The threshold must be established in each run. Please, if your equipment sets the threshold automatically, check and verify that it adjusts to the positive control or adjust it manually. For each channel, select the well corresponding the positive control, and fix the threshold within the exponential phase of the fluorescence curve and above any background signal (below the baseline). Once the threshold has been established, the rest of the samples of the same run can be interpreted.

The threshold value of the different channels could be different due to the chemical nature of the different fluorophores. The threshold value for different instruments may vary due to different signal intensities.

It is also recommended to include a blank, consisting of a confirmed negative sample (the same matrix as tested samples) for the targets detected in order to set the baseline of the assay.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for MPXV in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	MPXV (FAM) ¹	Endogenous Internal Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	>40 or no signal	Valid

Table 3. Expected Performance of Controls. Ct values. no signal = no amplification curves.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for the target channel), all results are reported as 'Invalid', and retesting is required.

2 The positive template control includes human housekeeping *HBB* gene target; therefore, amplification signals are observed in all target channels, including the Endogenous Internal Control (EIC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, **select in the thermocycler only the channels where targets are detected**. Then, use the following table, read, and analyse the results:

MPXV (FAM)	Endogenous Internal Control (HEX)	Interpretation for patient's individual samples	
≤40	≤40 ¹	Valid	MPXV DNA Detected
>40 or no signal	≤40 ²	Valid	Target not Detected ²
>40 or no signal	>40 or no signal ²	Invalid	Test Failure – Repeat Testing ²

Table 4. Interpretation of individual patient sample results. Ct values no signal = no amplification curves.


1 The Endogenous Internal Control (EIC) must show an amplification signal (Ct ≤40). The endogenous internal control is a human housekeeping gene that is presented in all nucleated cells. In case the target shows an amplification signal and the Endogenous Internal Control does not show any, the reaction should be considered 'Invalid' and it is recommended to repeat the qPCR diluting the DNA sample 1:10 following the laboratory guidelines and/or microbiology laboratory policy manuals. Ct values between 35-40 may indicate suboptimal sample collection, storage, and/or processing, so test repetition with recollected and/or reprocessed samples may be considered based on clinical judgement.

2 In the case of MPXV target genes negative, EIC must show an amplification signal with Ct ≤ 40. The Ct value could be very variable due to the Endogenous Internal Control is a human housekeeping gene that is present in all human nucleated cells in the original samples. If there is an absence of signal or Ct value > 40 of the Endogenous Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues, following the laboratory guidelines and/or microbiology laboratory policy manuals. Possible causes of a Ct value > 40 of the Endogenous Internal Control might be improper collection, storage, and/or processing of the specimens, failure of one or more of the reagents or the PCR thermal cycling process, as well as the presence of inhibitors.

In case of a continued ambiguous result, it is recommended to review the instructions for use; the extraction process used by the user; to verify the correct performance of each qPCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- a) repeat qPCR with the same isolated DNA sample, or
- b) re-extract and retest another aliquot of the same specimen or,
- c) obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad).

10. Limitations of the test

- The results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests.
- The Endogenous Internal Control (EIC) must show an amplification signal ($Ct \leq 40$). Ct values between 35-40 may indicate suboptimal sample collection, storage, and/or processing, so test repetition with recollected and/or reprocessed samples may be considered based on clinical judgement.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from lesion swabs.
- Possible crosstalk might be observed in empty channels of the thermocycler if there is no target to be detected, so it is required to select only the channels where these are amplified when interpretation of results is performed. Please contact to viasuresupport@certest.es for any queries.
- When using Qiagen/Corbett Rotor-Gene® instruments, only one VIASURE product must be used at a time to avoid potential crosstalk between channels. Please contact to viasuresupport@certest.es for any queries.
- The operational use of this product has been validated only within the recommended temperature range of 18–31 °C, which aligns with the operational specifications of commonly used real-time PCR thermocyclers (consult Table 2). In case of using a different temperature condition, the validation of the new operational temperature shall be performed by the end-user.
- The quality of the test depends on the quality of the sample: nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present. It is not possible to correlate the Ct values obtained by PCR with the concentration of the sample as they depend on the thermocycler used and the run itself.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.

- Please note the intended measurement range of the assay, as samples with concentrations above or below this range may give erroneous results.
- There is a possibility of false positive results due to cross-contamination by MPXV, either by *Monkeypox virus* Positive Control, by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown MPXV clades.
 - A viral load in the specimen below the limit of detection for the assay.
 - The presence of Real-Time PCR inhibitors or other types of interfering substances. The impacts of vaccines, some antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection have not been evaluated.
 - The effect of interfering substances has only been evaluated for those indicated in section 12.5.1 of this Instructions for Use (Interferences and inhibitors of PCR section). Please, see this section to check the most common endogenous and exogenous substances that induce a total or partial interference of Real-Time PCR reaction. Other substances not indicated in this part could lead to erroneous results.
 - Failure to follow instructions for use and the assay procedure.
- Potential microbial interference from the orthopoxviruses Camelpox and Variola could be obtained, since $\geq 80\%$ homology with one of the two targets detected with VIASURE assay is observed. The results of the test should be evaluated by a healthcare professional in the context of virus incidence, clinical symptoms, and other diagnostic tests.
- *In silico* analysis showed $\geq 80\%$ sequence homology with certain *Lactobacillus species*, *Staphylococcus epidermidis*, *Streptococcus pyogenes*, and *Pseudomonas aeruginosa*. Nevertheless, based on analytical, biological, and epidemiological plausibility assessments, the likelihood of non-specific amplification leading to interference is considered very low and acceptable. Consult section 12.5.1 Analytical specificity for more information.
- A positive test result does not necessarily indicate the presence of viable viruses and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targeted viral sequences.
- Negative results do not preclude MPXV infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral

levels during infections caused by MPXV have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the virus.

- If diagnostic tests for other poxvirus-related illnesses are negative and the patient's clinical presentation and epidemiological information suggest that MPXV infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment (even being the same model), extraction system, type of sample, previous treatment of the sample, etc... among others.
- Positive control and Negative control must not be extracted, but a blank sample is recommended to be included throughout the extraction and qPCR workflow to discard possible contaminations. Potential false positive results might remain undetected if blank sample is omitted.
- Some samples may fail to exhibit *HBB* amplification curves due to low human cell number in the original clinical sample. A negative EIC signal does not preclude the presence of MPXV in clinical specimens.

11. Quality control


VIASURE *Monkeypox virus* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. The Positive Control (PC) is a non-infectious synthetic DNA positive for the target under study to be used for the monitoring of the PCR amplification in each execution. The use of the PC is solely intended for amplification control purposes and has no quantitative objective. PC might not amplify due to the PCR inhibition by the presence of an artifact/inhibitor, due to a failure of one or more of the reagents or at the well-rehydration step, or because of a failure at the PCR thermal cycling process. Thus, if the PC does not amplify, the run is invalidated. On the other hand, the Negative Control (NC) is used to detect unspecific amplifications during the amplification step. If amplifications are observed, the run is invalidated. Finally, the Endogenous Internal Control (EIC) in each well confirms the correct performance of the technique and the extraction process, and its presence approves the quality and the integrity of the sample processing. A negative EIC indicates that no human DNA is present in the sample owing to insufficient sample collection or due to an improper sample collection, storage, and/or processing.

12. Analytical performance characteristics

12.1. Analytical linearity

The linearity of the assay was determined and confirmed by testing a series of ten-fold dilutions containing a known concentration (ranging from 10^7 to 10^1 copies per reaction, that is, from 2×10^9 to 2×10^3 copies/mL) of specific and synthetic DNA belonging to MPXV. Example of the amplification plot resulting from an assay run on the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) is included below:

Figure 2. Dilution series of MPXV (10^7 - 10^1 copies/rxn) template run on the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).

The linearity was also tested on DTPPrime real-time PCR instrument (DNA-Technology), Rotor-Gene® Q (QIAGEN) and 7500 Fast Real Time PCR System (Applied Biosystems™) with series of ten-fold dilutions containing a known concentration (ranging from 10^7 to 10^1 copies per reaction (2×10^9 to 2×10^3 copies/mL)) of specific and synthetic DNA belonging to MPXV. Results lead to confirm that the linearity, linear regression, efficiency and level of detection observed for the target detected with VIASURE *Monkeypox virus* Real Time PCR Detection Kit, using the four thermocyclers mentioned, meet the acceptance criteria.

12.2. Analytical sensitivity. Limit of Detection (LoD)

Analytical sensitivity or limit of detection (LoD) of VIASURE *Monkeypox virus* Real Time PCR Detection Kit was analysed using skin lesion swabs in Copan®'s UTM (for workflows 1-8), and skin lesion swabs in biocomma®'s VTM (workflow 9), spiked with reference material of MPXV. Samples were extracted and analysed following nine workflows:

	Extraction	Amplification
1	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.)	CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
2	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.)	DTPPrime real-time PCR instrument (DNA-Technology)
3	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.)	Rotor-Gene® Q (QIAGEN)
4	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.)	7500 Fast Real Time PCR System (Applied Biosystems™)
5	QIAamp DNA Mini Kit (QIAGEN)	CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
6	QIAamp DNA Mini Kit (QIAGEN)	DTPPrime real-time PCR instrument (DNA-Technology)
7	QIAamp DNA Mini Kit (QIAGEN)	Rotor-Gene® Q (QIAGEN)
8	QIAamp DNA Mini Kit (QIAGEN)	7500 Fast Real Time PCR System (Applied Biosystems™)
9	MagCore® Viral Nucleic Acid Extraction Kit / MagCore® automated nucleic acid extractors (RBC Bioscience)	CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)

Table 5. Extraction and amplification methods (workflows) validated using VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

VIASURE *Monkeypox virus* Real Time PCR Detection Kit showed an experimental detection limit indicated in the following table, with a positive rate of $\geq 95\%$:

skin lesion swabs in Copan®'s UTM								skin lesion swabs in biocomma®'s VTM
Workflow 1	Workflow 2	Workflow 3	Workflow 4	Workflow 5	Workflow 6	Workflow 7	Workflow 8	Workflow 9
MPXV	MPXV	MPXV	MPXV	MPXV	MPXV	MPXV	MPXV	MPXV
25 copies/mL	50 copies/mL	50 copies/mL	25 copies/mL	25 copies/mL	50 copies/mL	25 copies/mL	50 copies/mL	15 copies/mL

Table 6. Experimental LoD results of VIASURE *Monkeypox virus* Real Time PCR Detection Kit using nine workflows and skin lesion swabs in Copan®'s UTM or biocomma®'s VTM spiked with reference material for MPXV.

After using the Probit Model to analyse the results, VIASURE *Monkeypox virus* Real Time PCR Detection Kit showed the detection limit indicated in the following table, with a positive rate of $\geq 95\%$:

skin lesion swabs in Copan®'s UTM								skin lesion swabs in biocomma®'s VTM
Workflow 1	Workflow 2	Workflow 3	Workflow 4	Workflow 5	Workflow 6	Workflow 7	Workflow 8	Workflow 9
MPXV	MPXV	MPXV	MPXV	MPXV	MPXV	MPXV	MPXV	MPXV
14 copies/mL	28 copies/mL	48 copies/mL	23 copies/mL	24 copies/mL	32 copies/mL	23 copies/mL	28 copies/mL	14 copies/mL

Table 7. LoD results of VIASURE *Monkeypox virus* Real Time PCR Detection Kit obtained after Probit analysis, using nine workflows and skin lesion swabs in Copan®'s UTM or biocomma®'s VTM spiked with reference material for MPXV.

The obtained LoDs per workflow comply with the acceptance criteria initially stated. Nevertheless, experimental LoDs were considered as reference for the rest of performance analyses.

12.3. Measuring range

The measuring range of the assay was determined by testing a series of ten-fold dilutions containing a known concentration of specific and synthetic DNA belonging to MPXV, and further refined with the LoD assay. Results allowed to confirm the correct detection of the target from 2×10^9 copies/mL 25 copies/mL, complying with the acceptance validation criteria initially stated.

12.4. Accuracy

12.4.1. Trueness (Veracity)

The veracity of VIASURE *Monkeypox virus* Real Time PCR Detection Kit was evaluated by testing reference material listed below. Reference strain were extracted with MagDEA Dx SV Kit, using the magLEAD® 12gC instrument (Precision System Science Co.), and analysed with the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) / CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad) and the DTPrime real-time PCR instrument (DNA-Technology). In the case of reference DNA/RNA, they were directly analysed with the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad) or the DTPrime real-time PCR instrument (DNA-Technology).

1. Synthetic DNA fragments:

- Name of synthetic DNA fragment for MPXV (*GR2* gene): MPXXPC, FAM channel.
- Name of synthetic DNA fragment for MPXV (*F3L* gene): MPXXPC, FAM channel.

2. ZeptoMetrix Corporation:

- NATtrol™ Monkeypox Virus External Run Control, isolate USA/MA001/2022* (ZeptoMetrix Ref: NATMPXV-ERC).

* This reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: Monkeypox Virus, hMPXV/USA/MA001/2022 (Lineage B.1, Clade IIb), NR-58622.

3. Control materials:

External Reference	Pathogen Name	Variety
MBC146-R	AMPLIRUN® MONKEYPOX VIRUS DNA CONTROL	West African clade

Table 8. Monkeypox virus reference material from Vircell S.L.

4. EQA Programmes

Procedence	Pathogen Name	Samples correctly detected
QCMD 2022 and 2023 Poxviruses EQA Pilot Study	Monkeypox Clade I (Congo Basin clade /Central African isolate)	POX22S2-01, POX22S2-06, POX22S2-07, POX22S2-08, POX22S2-10, POX23C1A-01, POX23C1A-04, POX23C1A-05, POX23C1A-07, POX23C1B-07 and POX23C1B-10
QCMD 2022 and 2023 Poxviruses EQA Pilot Study	Monkeypox Clade IIb (West African clade/2022 outbreak isolate)	POX22S2-03, POX22S2-04, POX23C1A-03, POX23C1A-09, POX23C1A-10, POX23C1B-02, POX23C1B-04, POX23C1B-05, POX23C1B-08 and POX23C1B-09

Table 9. Reference material from EQA programmes.

12.4.2. Precision

To determine the precision of VIASURE *Monkeypox virus* Real Time PCR Detection Kit, intra-assay (repeatability), inter-assay, inter-batch, and inter-thermocycler assays (reproducibility) were performed. The panel of samples were analysed following the workflow MagDEA Dx SV Kit, using the magLEAD® 12gC instrument (Precision System Science Co.) / DTPPrime real-time PCR instrument (DNA-Technology).

Intra-assay

The intra-assay was tested by analysing replicates of all samples in the same run using VIASURE *Monkeypox virus* Real Time PCR Detection Kit and the selected workflow. The following table shows the results obtained in this assay.

Sample	Target	Viasure channel	$\bar{x}$ (Ct)	σ	CV %
Positive 1 (6xLoD)	MPXV	FAM	32.14	0.32	1.01
Positive 2 (3xLoD)	MPXV	FAM	33.48	0.50	1.49
Negative sample	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	26.56	0.27	1.00
Positive Control	MPXV	FAM	24.30	0.09	0.38
Negative Control	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	Neg	n.a.	n.a.

Table 10. Intra-assay results of VIASURE *Monkeypox virus* Real Time PCR Detection Kit. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-assay

The inter-assay was tested by analysing the different samples on three different days by three different operators with the VIASURE *Monkeypox virus* Real Time PCR Detection Kit run on DTPPrime real-time PCR instrument (DNA-Technology). A summary of results is shown in the table below.

Sample	Target	Viasure channel	$\bar{x}$ (Ct)	σ	CV %
Positive 1 (6xLoD)	MPXV	FAM	32.59	0.52	1.61
Positive 2 (3xLoD)	MPXV	FAM	33.33	0.30	0.89
Negative sample	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	26.66	0.71	2.67
Positive Control	MPXV	FAM	24.40	0.20	0.81
Negative Control	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	Neg	n.a.	n.a.

Table 11. Inter-assay results of VIASURE *Monkeypox virus* Real Time PCR Detection Kit. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-batch

The inter-batch values were determined with three replicates of the different samples by using three batches of VIASURE *Monkeypox virus* Real Time PCR Detection Kit, run on DTPPrime real-time PCR instrument (DNA-Technology). The results are described below in the following table.

Sample	Target	Viasure channel	$\bar{x}$ (Ct)	σ	CV %
Positive 1 (6xLoD)	MPXV	FAM	32.30	0.58	1.79
Positive 2 (3xLoD)	MPXV	FAM	33.42	0.48	1.43
Negative sample	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	26.86	0.45	1.67
Positive Control	MPXV	FAM	24.46	0.40	1.65
Negative Control	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	Neg	n.a.	n.a.

Table 12. Inter-batch results of VIASURE *Monkeypox virus* Real Time PCR Detection Kit. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-thermocycler

The inter-thermocycler values were determined with three replicates of the same samples used for intra-assay, inter-assay and inter-batch assay, using the VIASURE *Monkeypox virus* Real Time PCR Detection

Kit. The DNA extracted from samples and positive and negative controls were tested all together in the same run, on CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad), 7500 Fast Real Time PCR System (Applied Biosystems™) and DTPrime real-time PCR instrument (DNA-Technology). The results are described below in the following table.

Sample	Target	Viasure channel	$\bar{x}$ (Ct)	σ	CV %
Positive 1 (6xLoD)	MPXV	FAM	32.19	0.26	0.82
Positive 2 (3xLoD)	MPXV	FAM	33.13	0.31	0.95
Negative sample	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	26.51	0.69	2.61
Positive Control	MPXV	FAM	24.41	0.26	1.05
Negative Control	MPXV	FAM	Neg	n.a.	n.a.
	EIC	HEX	Neg	n.a.	n.a.

Table 13. Inter-thermocycler results of VIASURE *Monkeypox virus* Real Time PCR Detection Kit. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

12.4.3. Accuracy

The accuracy of VIASURE *Monkeypox virus* Real Time PCR Detection Kit was determined considering the combination of the trueness and precision properties. Taking the obtained data of the precision experiment for the 3xLoD concentration point as reference of the expected Ct values (Trueness), 3 different batches were used to test samples at the same concentration used in precision assays. A total of 180 data points (Precision) were analysed: three runs were performed in three different thermocyclers, and in each run, a total of 60 replicates (20 replicates amplified per batch) were included. Samples were extracted with MagDEA DX SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) and analysed with the CFX Opus 96™ Real-Time PCR Detection System (Bio-Rad), 7500 Fast Real Time PCR System (Applied Biosystems™) and DTPrime real-time PCR instrument (DNA-Technology).

Sample	Target	Viasure channel	Accuracy reference values (Trueness)			Accuracy results (Precision)		
			$\bar{x}$ (Ct)	σ	CV%	$\bar{x}$ (Ct)	σ	CV%
Sample 1 (3xLoD)	MPXV	FAM	32.40	0.50	1.55	32.85	0.35	1.05

Table 14. Accuracy results in comparison to accuracy reference values (precision assay results → Trueness) of VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

The results obtained showed that the mean Ct value obtained from the 180 replicates was within the corresponding range of the mean Ct value obtained in the precision experiments $\pm$ the standard deviation for the target in study. In addition to this, the deviation and the variation coefficient meet the acceptance requirements.

12.5. Analytical specificity and reactivity

The analytical specificity and analytical reactivity were evaluated for the VIASURE *Monkeypox virus* Real Time PCR Detection Kit *in silico* and experimentally, using different starting material such as certified reference strains, certified reference RNA/DNAs and material from the EQAs programmes.

12.5.1. Analytical specificity

Analytical specificity is the assay's ability to detect the intended target. There are two components to be considered for analytical specificity: Cross-reactivity and Interference. Cross-reactivity may occur when genetically related genes/sequences are present in a patient specimen, while interference may happen if the presence of specific substances potentially present in the sample matrixes affects the performance of the qPCR.

Cross-reactivity *in silico* analysis

The cross-reactivity was assessed by using sequences of other microorganisms related, but different from the detected target. Sequences downloaded from the NCBI GenBank (<https://www.ncbi.nlm.nih.gov/genbank/>) were evaluated, while alignment of primers and probes were performed with BLAST (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>) and an in-house bioinformatic analysis software. Results shown below only indicate sequences with a homology percentage higher than 80%. Aligned sequences with a percentage of alignment less than 80% of homology were considered unlikely to be detected. The most relevant results are shown in the following table:

Microorganism	ID Reference in NCBI or RefSeq number	% Homology
Borealpox virus ^a	MN240300	85%
Camelpox virus ^b	KC841322, KC841321, KC841320, KC841319, U87839, U87837, U87838, U87840	91.25%
Cowpox virus ^c	OM460002	97.26%
	LT896721, LR812035	95.89%
	KY549148, KY549144, HQ420900	94.52%
	KC813494, KC813493, FJ769342, FJ769341, FJ769340, FJ769339, FJ769337, AY623479, U90235, U90234, U90233, U90231, U90228, KY569019, KC813505, KC813496, FJ769336, U90227, U90226, KY549147, KY549146, KU861504, KC813510, KC813503	91.25%
	KC813509, KC813495, KC813492, KC813491, HQ328780, U90232, U90230, U90229, U90225, KC813506, KC813498, KC813497, FJ769338, KY549151	90%
Taterapox virus ^d	U86874	93.33%
Vaccinia virus ^e	MN369532, BK013340	90%
Variola major virus ^f	L22579	93.51%
Variola minor virus ^g	Y16780	93.51%
Variola virus ^h	LR800245, LR800244	94.81%
	BK010317, LT706529, LT706528, KY358055, NC_001611, OL468962, OL468961, DQ437592, DQ437591, DQ437590, DQ437589, DQ437588, DQ437587, DQ437586, DQ437585, DQ437584, DQ437583, DQ437582,	93.51%

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Table 15. Percentage of homology of Monkeypox virus primers and probe set in the cross-reactivity analysis, comparing reference sequences of Monkeypox virus-like pathogens from the NCBI database.

^a In the case of the Borealex virus (formerly Alaskapox), the detected sequence is the only Borealex virus sequence available in the NCBI Nucleotide database, and aligns with a very low homology against the monkeypox detection primers and probe set, being its detection very unlikely. Moreover, only a few cases of borealex have been confirmed, with no evidence of human-to-human infection⁵. However, the possible detection of the pathogen has been analysed experimentally by testing a series of ten-fold dilutions containing a known concentration (ranging from 10⁷ to 10¹ copies per reaction) of specific and synthetic DNAs belonging to the mismatches observed in the Borealex virus sequences, and results obtained showed that VIASURE assay barely detects Borealex virus, discarding possible cross-reactivity issues.

^b Relating to Camelpox disease, the virus responsible of this condition is usually restricted to camels ^{6,7}, affecting humans only in punctual cases, with no reported human-to-human transmission.

^{c,d} Sporadic cases of Cowpox⁸ and Taterapox⁹ have been reported in Europe, mostly linked to handling infected animals, usually rodents and cats. There is not known evidence of human-to-human transmission. In addition, in the case of the Cowpox virus it was verified the non-detection in a wet analysis conducted with the human-like specimen containing the Cowpox, OPV 88-Lunge (indicated in section “Cross-reactivity wet testing”).

^e Vaccinia virus is the pathogen used for the smallpox vaccine, being responsible of its eradication. Therefore, human cases are extremely rare, with no human-to-human infections notified ^{10,11}. Again, the non-detection of this virus was verified in a wet analysis conducted with a human-like specimen containing the Vaccinia virus, Elstree (indicated in section “Cross-reactivity wet testing”).

⁵ Borealex. (n.d.). Retrieved April 18, 2024, from <https://health.alaska.gov/dph/Epi/id/Pages/Borealex.aspx>

⁶ Balamurugan V, Venkatesan G, Bhanuprakash V, Singh RK. Camelpox, an emerging orthopox viral disease. Indian J Virol. 2013 Dec;24(3):295-305. doi: 10.1007/s13337-013-0145-0. Epub 2013 Jul 16. PMID: 24426291; PMCID: PMC3832703.

⁷ Dahiya, S. S., Kumar, S., Mehta, S. C., Narnaware, S. D., Singh, R., & Tuteja, F. C. (2016). Camelpox: A brief review on its epidemiology, current status and challenges. Acta tropica, 158, 32–38. <https://doi.org/10.1016/j.actatropica.2016.02.014>

⁸ Krankowska, D. C., Woźniak, P. A., Cybula, A., Izdebska, J., Suchacz, M., Samelska, K., Wiercińska-Drapała, A., & Szaflik, J. P. (2021). Cowpox: How dangerous could it be for humans? Case report. International journal of infectious diseases : IJID : official publication of the International Society for Infectious Diseases, 104, 239–241. <https://doi.org/10.1016/j.ijid.2020.12.061>

⁹ Parker, S., Crump, R., Hartzler, H., & Buller, R. M. (2017). Evaluation of Taterapox Virus in Small Animals. Viruses, 9(8), 203. <https://doi.org/10.3390/v9080203>

¹⁰ Xiang, Y., & Lane, R. K. (2021). Vaccinia Virus (Poxviridae). Encyclopedia of Virology: Volume 1-5, Fourth Edition, 1–5, 854–859. <https://doi.org/10.1016/B978-0-12-814515-9.00047-3>

¹¹ Henderson D. A. (2011). The eradication of smallpox--an overview of the past, present, and future. Vaccine, 29 Suppl 4,D7–D9. <https://doi.org/10.1016/j.vaccine.2011.06.080>

^{f,g,h} Finally, Smallpox (Variola virus) was declared eradicated (eliminated) by the World Health Assembly in 1980, and no cases of naturally occurring smallpox have happened since ¹².

A more detailed information about sequences homologies $\geq 80\%$ with the forward, reverse and probe oligonucleotides included in VIASURE reaction mix can be found in the following tables. Table 16 includes the organisms that align with the forward, reverse and probe oligonucleotides for the G2R target detection with a percentage of alignment higher than 80%, but with a homology percentage lower than 80% with the forward, reverse and probe oligonucleotides for the F3L target detection (data not shown). Table 17 shows the list of organisms that align with both sets of primers and probe for the *G2R* and *F3L* target detection with a percentage of alignment higher than 80%. Finally, Table 18 includes the analysed microorganisms and its percentage of alignment (homology) with G2R forward, G2R reverse, G2R probe, F3L forward, F3L reverse and F3L probe with a percentage greater than 80%. Please note that PCR amplification is expected only if the three oligonucleotides belonging to one set (F3L set or G2R set) achieve a percentage of homology greater than 80%, and/or if the amplicon is lower than 500 nucleotides. If not, cross-reactivity can be discarded.

FASTA Name	Homology <i>G2R</i> target- Forward (%)	Homology <i>G2R</i> target- probe (%)	Homology <i>G2R</i> target- reverse (%)	Total Homology (%)
AY623479.1 Cowpox virus strain CoPV/DEN/98 tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	96.15	83.33	95.83	91.25
BK010317.1 TPA_asm: Variola virus isolate VD21, complete genome	88.46	86.67	100	91.25
BK013340.1 TPA_asm: Vaccinia virus isolate VK02, complete genome	92.31	83.33	95.83	90
DQ437580.1 Variola virus strain Afghanistan 1970 Variolator 4, complete genome	88.46	86.67	100	91.25
DQ437581.1 Variola virus strain Bangladesh 1975 v75-550 Banu, complete genome	88.46	86.67	100	91.25
DQ437582.1 Variola virus strain China Horn 1948, complete genome	88.46	86.67	100	91.25
DQ437583.1 Variola virus strain Congo 1970, complete genome	88.46	86.67	100	91.25
DQ437584.1 Variola virus strain Germany 1958 Heidelberg, complete genome	88.46	86.67	100	91.25
DQ437585.1 Variola virus strain India 1964 7124 Vellore, complete genome	88.46	86.67	100	91.25
DQ437586.1 Variola virus strain India 1964 7125 Vellore, complete genome	88.46	86.67	100	91.25
DQ437587.1 Variola virus strain Iran 1972 2602 Tabriz, complete genome	88.46	86.67	100	91.25
DQ437588.1 Variola virus strain Nepal 1973, complete genome	88.46	86.67	100	91.25
DQ437589.1 Variola virus strain Pakistan 1969 (Rafiq Lahore), complete sequence	88.46	86.67	100	91.25
DQ437590.1 Variola virus strain Somalia 1977, complete genome	88.46	86.67	100	91.25
DQ437591.1 Variola virus strain Sumatra 1970 V70-222, complete genome	88.46	86.67	100	91.25
DQ437592.1 Variola virus strain Syria 1972 V72-199, complete genome	88.46	86.67	100	91.25

¹² <https://www.cdc.gov/smallpox/index.html>

DQ441416.1 Variola virus strain Benin, Dahomey 1968 (v68-59), complete genome	88.46	86.67	100	91.25
DQ441417.1 Variola virus strain Botswana 1972 (v72-143), complete genome	88.46	86.67	100	91.25
DQ441418.1 Variola virus strain Botswana 1973 (v73-225), complete genome	88.46	86.67	100	91.25
DQ441419.1 Variola virus strain Brazil 1966 (v66-39 Sao Paulo), complete genome	88.46	86.67	100	91.25
DQ441420.1 Variola virus strain Bangladesh 1974 (nur islam), complete genome	88.46	86.67	100	91.25
DQ441421.1 Variola virus strain Bangladesh 1974 (Shahzaman), complete genome	88.46	86.67	100	91.25
DQ441423.1 Variola virus strain Congo 9 1970 (v74-227 Gispen), complete genome	88.46	86.67	100	91.25
DQ441424.1 Variola virus strain Ethiopia 1972 (Eth16 R14-1X-72 Addis), complete genome	88.46	86.67	100	91.25
DQ441425.1 Variola virus strain Ethiopia 1972 (Eth17 R14-1X-72 Addis), complete genome	88.46	86.67	100	91.25
DQ441426.1 Variola virus strain Guinea 1969 (005), complete genome	88.46	86.67	100	91.25
DQ441427.1 Variola virus strain India 1953 (Kali-Muthu-M50 Madras), complete genome	88.46	86.67	100	91.25
DQ441428.1 Variola virus strain India 1953 (New Delhi), complete genome	88.46	86.67	100	91.25
DQ441429.1 Variola virus strain Japan 1946 (Yamada MS-2(A) Tokyo), complete genome	88.46	86.67	100	91.25
DQ441430.1 Variola virus strain Japan 1951 (Harper, Masterseed), complete genome	88.46	86.67	100	91.25
DQ441431.1 Variola virus strain Japan 1951 (Stillwell, Masterseed), complete genome	88.46	86.67	100	91.25
DQ441432.1 Variola virus strain Korea 1947 (Lee, Masterseed), complete genome	88.46	86.67	100	91.25
DQ441433.1 Variola virus strain Kuwait 1967 (K1629), complete genome	88.46	86.67	100	91.25
DQ441434.1 Variola virus strain Niger 1969 (001, importation from Nigeria), complete genome	88.46	86.67	100	91.25
DQ441435.1 Variola virus strain South Africa 1965 (102 Natal, Ingwavuma), complete genome	88.46	86.67	100	91.25
DQ441436.1 Variola virus strain South Africa 1965 (103 T'vaal, Nelspruit), complete genome	88.46	86.67	100	91.25
DQ441437.1 Variola virus strain Sierra Leone 1969 (V68-258), complete genome	88.46	86.67	100	91.25
DQ441438.1 Variola virus strain Somalia 1977 (V77-1252), complete genome	88.46	86.67	100	91.25
DQ441439.1 Variola virus strain Somalia 1977 (V77-1605), complete genome	88.46	86.67	100	91.25
DQ441440.1 Variola virus strain Sudan 1947 (Juba), complete genome	88.46	86.67	100	91.25
DQ441441.1 Variola virus strain Sudan 1947 (Rumbec), complete genome	88.46	86.67	100	91.25
DQ441442.1 Variola virus strain Sumatra 1970 V70-228, complete genome	88.46	86.67	100	91.25
DQ441443.1 Variola virus strain Tanzania 1965 kembula, complete genome	88.46	86.67	100	91.25
DQ441444.1 Variola virus strain United Kingdom 1946 Harvey, complete genome	88.46	86.67	100	91.25
DQ441445.1 Variola virus strain United Kingdom 1946 Hinden (Middlesex), complete genome	88.46	86.67	100	91.25
DQ441446.1 Variola virus strain United Kingdom 1947 Higgins (Staffordshire), complete genome	88.46	86.67	100	91.25

DQ441447.1 Variola virus strain United Kingdom 1952 Butler, complete genome	88.46	86.67	100	91.25
DQ441448.1 Variola virus strain Yugoslavia 1972 V72-164, complete genome	88.46	86.67	100	91.25
FJ769336.1 Cowpox virus isolate FIN/T2000 cytokine response modifier B (CrmB) gene, complete cds	92.31	83.33	95.83	90
FJ769337.1 Cowpox virus isolate No-H2 cytokine response modifier B (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
FJ769338.1 Cowpox virus isolate No-F1 cytokine response modifier B (CrmB) gene, complete cds	92.31	83.33	95.83	90
FJ769339.1 Cowpox virus isolate No-F2 cytokine response modifier B (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
FJ769340.1 Cowpox virus isolate No-H1 cytokine response modifier B (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
FJ769341.1 Cowpox virus isolate Swe-H1 cytokine response modifier B (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
FJ769342.1 Cowpox virus isolate Swe-H2 cytokine response modifier B (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
HQ328780.1 Cowpox virus isolate IT-Llama09 tumor necrosis factor receptor II-like protein (crmB) gene, partial cds	92.31	83.33	95.83	90
KC813491.1 Cowpox virus strain BeaBer04/1, complete genome	92.31	83.33	95.83	90
KC813492.1 Cowpox virus strain HumLan08/1, complete genome	92.31	83.33	95.83	90
KC813493.1 Cowpox virus strain HumLit08/1, complete genome	96.15	83.33	95.83	91.25
KC813494.1 Cowpox virus strain HumLue09/1, complete genome	96.15	83.33	95.83	91.25
KC813495.1 Cowpox virus strain HumMag07/1, complete genome	92.31	83.33	95.83	90
KC813496.1 Cowpox virus strain HumPad07/1, complete genome	92.31	83.33	95.83	90
KC813497.1 Cowpox virus strain JagKre08/1, complete genome	92.31	83.33	95.83	90
KC813498.1 Cowpox virus strain JagKre08/2, complete genome	92.31	83.33	95.83	90
KC813503.1 Cowpox virus strain RatGer09/1, complete genome	92.31	83.33	95.83	90
KC813505.1 Cowpox virus strain RatKre08/2, complete genome	92.31	83.33	95.83	90
KC813506.1 Cowpox virus strain CatPot07/1, complete genome	92.31	83.33	95.83	90
KC813509.1 Cowpox virus strain HumBer07/1, complete genome	92.31	83.33	95.83	90
KC813510.1 Cowpox virus strain HumGra07/1, complete genome	92.31	83.33	95.83	90
KC841319.1 Camelpox virus isolate CMLV-DELHI tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92.31	86.67	95.83	91.25
KC841320.1 Camelpox virus isolate CMLV-BIKANER tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92.31	86.67	95.83	91.25
KC841321.1 Camelpox virus isolate CMLV-BARMER tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92.31	86.67	95.83	91.25
KC841322.1 Camelpox virus isolate CMLV-JAISALMER tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92.31	86.67	95.83	91.25
KU861504.1 Cowpox virus isolate Kostroma/2015 cytokine response modifier B gene, complete cds	92.31	86.67	95.83	91.25
KY358055.1 Variola virus strain VD21, 17th century, complete genome	88.46	86.67	100	91.25
KY549146.1 Cowpox virus isolate CPXV_K779, complete genome	92.31	83.33	95.83	90
KY549147.1 Cowpox virus isolate CPXV_K780, complete genome	92.31	83.33	95.83	90
KY549151.1 Cowpox virus isolate Norwayfeline, complete genome	92.31	83.33	95.83	90
KY569019.1 Cowpox virus strain CPXV CheNova_DK_2014, complete genome	96.15	83.33	95.83	91.25
L22579.1 Variola major virus (strain Bangladesh-1975) complete genome	88.46	86.67	100	91.25

LR800244.1 Variola virus genome assembly, complete genome: monopartite	92.31	86.67	100	92.5
LR800245.1 Variola virus genome assembly, complete genome: monopartite	92.31	86.67	100	92.5
LR800246.1 Variola virus genome assembly, complete genome: monopartite	88.46	86.67	95.83	90
LR800247.1 Variola virus genome assembly, complete genome: monopartite	88.46	86.67	95.83	90
LT706528.1 Variola virus isolate VARV_V563 genome assembly, chromosome: I	88.46	86.67	100	91.25
LT706529.1 Variola virus isolate VARV_V1588 genome assembly, chromosome: I	88.46	86.67	100	91.25
MN240300.1 Boreapox virus, complete genome	88.46	80	87.5	85
MN369532.1 Vaccinia virus isolate VK01, complete genome	92.31	83.33	95.83	90
NC_001611.1 Variola virus, complete genome	88.46	86.67	100	91.25
OL468961.1 Variola virus strain GHA68, complete genome	88.46	86.67	100	91.25
OL468962.1 Variola virus strain Congo 1970 v70-60, complete genome	88.46	86.67	100	91.25
U18339.1 Variola virus Garcia-1966 right near-terminal region	88.46	86.67	100	91.25
U18341.1 Variola virus Somalia-1977 right near-terminal region	88.46	86.67	100	91.25
U86874.1 Taterapox virus tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92.31	86.67	91.67	90
U87837.1 Camelpox virus strain Somalia-1978 tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92.31	86.67	95.83	91.25
U87838.1 Camelpox virus CP1 tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92.31	86.67	95.83	91.25
U87839.1 Camelpox virus strain Saudi-M3, tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92.31	86.67	95.83	91.25
U87840.1 Camelpox virus CP5 tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92.31	86.67	95.83	91.25
U88145.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U88146.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U88147.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U88148.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U88150.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U88151.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U88152.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88.46	86.67	100	91.25
U90225.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92.31	83.33	95.83	90
U90226.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92.31	83.33	95.83	90
U90227.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92.31	83.33	95.83	90
U90228.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
U90229.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92.31	83.33	95.83	90
U90230.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92.31	83.33	95.83	90

U90231.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
U90232.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92.31	83.33	95.83	90
U90233.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
U90234.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
U90235.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96.15	83.33	95.83	91.25
X67117.1 Variola virus (XhoI-F,O,H,P,Q genome fragment) genes	88.46	86.67	100	91.25
X69198.1 Variola virus DNA complete genome	88.46	86.67	100	91.25
X70841.1 Variola virus genes for ORF1R, ORF2R and ORF3R	88.46	86.67	100	91.25
Y16780.1 variola minor virus complete genome	88.46	86.67	100	91.25

Table 16. Percentage homology of forward, reverse and probe for the *G2R* target detection (VIASURE Monkeypox Virus Real Time PCR Detection Kit) of the aligned sequences (>80%).

FASTA Name	Homology <i>G2R</i> target-forward (%)	Homology <i>G2R</i> target - probe (%)	Homology <i>G2R</i> target - reverse (%)	Homology <i>F3L</i> target - forward (%)	Homology <i>F3L</i> target - probe (%)	Homology <i>F3L</i> target-reverse (%)
HQ420900.1 Cowpox virus strain UK2000_K2984, complete genome	92.31	86.67	95.83	100.00	92.00	91.30
KY549144.1 Cowpox virus isolate CPXV_Catpox5wv1, complete genome	92.31	86.67	95.83	100.00	92.00	91.30
KY549148.1 Cowpox virus isolate CPXV_1639, complete genome	92.31	86.67	95.83	100.00	92.00	91.30
LR812035.1 Cowpox virus isolate GerMygEK 938/17 genome assembly, complete genome: monopartite	96.15	83.33	95.83	96.00	96.00	95.65
LT896721.1 Cowpox virus isolate Ger 2010 MKY genome assembly, complete genome: monopartite	96.15	83.33	95.83	96.00	96.00	95.65
OM460002.1 Cowpox virus isolate No-H2, complete genome	96.15	83.33	95.83	100.00	96.00	95.65

Table 17. Percentage homology of forward, reverse and probe for the *G2R* and *F3L* target detection (VIASURE Monkeypox Virus Real Time PCR Detection Kit) of the aligned sequences (>80%).

Pathogens/microorganisms	Taxonomy ID	Accession ID	<i>F3L</i> gene-FWD/ *	<i>F3L</i> gene-RV/ *	<i>F3L</i> gene-Probe/ *	<i>G2R</i> gene-FWD/ *	<i>G2R</i> gene-RV/ *	<i>G2R</i> gene-Probe/ *
Variola virus (smallpox)	10255	NC_001611.1	NA	86.96 %	80.00 %	88.46 %	100%	86.67 %
Molluscum contagiosum virus	10280	NC_001731.1	NA	NA	NA	NA	NA	NA
Parapoxvirus (e.g., Orf virus) Parapoxvirus (txid:10257)	129726	NC_013804.1	NA	NA	NA	NA	NA	NA
	1579460	NC_025963.1	NA	NA	NA	NA	NA	NA
	10258	NC_005336.1	NA	NA	NA	NA	NA	NA
	129727	NC_005337.1	NA	NA	NA	NA	NA	NA
	187984	NC_035188.1	NA	NA	NA	NA	NA	NA
Vaccinia virus Orthopoxvirus vaccinia (txid:10245)	10245	NC_006998.1	NA	86.96 %	80.00 %	92.31 %	95.83 %	NA
Cowpox virus	10243	NC_003663.2	NA	91.30 %	80.00 %	92.31 %	95.83 %	86.67 %
Ectromelia (mousepox) virus	12643	NC_004105.1	NA	86.96 %	80.00 %	96.15 %	NA	NA
Camelpox virus	28873	NC_003391.1	NA	86.96	NA	92.31	95.83	86.67

				%		%	%	%
Buffalopox	28873	MG599038.1	NA	86.96 %	80.00 %	92.31 %	95.83 %	NA
Herpes simplex virus (HSV-1 and HSV-2)	10298	NC_001806.2	NA	NA	NA	NA	NA	NA
	10310	NC_001798.2	NA	NA	NA	NA	NA	NA
Varicella-zoster virus (Chickenpox)	10335	NC_001348.1	NA	NA	NA	NA	NA	NA
<i>Streptococcus mitis</i>	28037	CP012646.1	NA	NA	NA	NA	NA	NA
<i>Staphylococcus aureus</i>	93061	NC_007795.1	NA	NA	NA	NA	NA	NA
<i>Staphylococcus epidermidis</i>	1282	CP035288.1	80.00 %	NA	NA	NA	NA	NA
		CP035289.1	NA	NA	NA	NA	NA	NA
		CP035290.1	NA	NA	NA	NA	NA	NA
<i>Streptococcus pyogenes</i>	1314	LS483338.1	80.00 %	NA	NA	NA	NA	NA
<i>Streptococcus agalactiae</i>	1311	CP012480.1	NA	NA	NA	NA	NA	NA
		CP012742.1	NA	NA	NA	NA	NA	NA
<i>Pseudomonas aeruginosa</i>	208964	AE004091.2	NA	82.61 %	NA	NA	NA	NA
<i>Trichophyton rubrum</i>	5551	GG700682.1	NA	NA	NA	NA	NA	NA
		GG700681.1	NA	NA	NA	NA	NA	NA
		GG700680.1	NA	NA	NA	NA	NA	NA
		GG700679.1	NA	NA	NA	NA	NA	NA
		GG700678.1	NA	NA	NA	NA	NA	NA
		GG700677.1	NA	NA	NA	NA	NA	NA
		GG700676.1	NA	NA	NA	NA	NA	NA
		GG700675.1	NA	NA	NA	NA	NA	NA
		GG700674.1	NA	NA	NA	NA	NA	NA
		GG700673.1	NA	NA	NA	NA	NA	NA
		GG700672.1	NA	NA	NA	NA	NA	NA
		GG700671.1	NA	NA	NA	NA	NA	NA
		GG700670.1	NA	NA	NA	NA	NA	NA
		GG700669.1	NA	NA	NA	NA	NA	NA
		GG700668.1	NA	NA	NA	NA	NA	NA
		GG700667.1	NA	NA	NA	NA	NA	NA
		GG700666.1	NA	NA	NA	NA	NA	NA
		GG700665.1	NA	NA	NA	NA	NA	NA
		GG700664.1	NA	NA	NA	NA	NA	NA
		GG700663.1	NA	NA	NA	NA	NA	NA
		GG700662.1	NA	NA	NA	NA	NA	NA
		GG700661.1	NA	NA	NA	NA	NA	NA
		GG700660.1	NA	NA	NA	NA	NA	NA
		GG700659.1	NA	NA	NA	NA	NA	NA
		GG700658.1	NA	NA	NA	NA	NA	NA
		GG700657.1	NA	NA	NA	NA	NA	NA
		GG700656.1	NA	NA	NA	NA	NA	NA
		GG700655.1	NA	NA	NA	NA	NA	NA
		GG700654.1	NA	NA	NA	NA	NA	NA
		GG700653.1	NA	NA	NA	NA	NA	NA
		GG700652.1	NA	NA	NA	NA	NA	NA
		GG700651.1	NA	NA	NA	NA	NA	NA

			GG700650.1	NA	NA	NA	NA	NA	NA
			GG700649.1	NA	NA	NA	NA	NA	NA
			GG700648.1	NA	NA	NA	NA	NA	NA
			GG730037.1	NA	NA	NA	NA	NA	NA
<i>Corynebacterium jeikeium</i>		38289	CP063195.1	NA	NA	NA	NA	NA	NA
<i>Candida albicans</i>		237561	CP017623.1	NA	NA	NA	NA	NA	NA
			CP017624.1	NA	NA	NA	NA	NA	NA
			CP017625.1	NA	NA	NA	NA	NA	NA
			CP017626.1	NA	NA	NA	NA	NA	NA
			CP017627.1	NA	NA	NA	NA	NA	NA
			CP017628.1	NA	NA	NA	NA	NA	NA
			CP017629.1	NA	NA	NA	NA	NA	NA
			CP017630.1	NA	NA	NA	NA	NA	NA
Human Genomic DNA		9606	GRCh38.p14	NA	NA	NA	NA	NA	NA
<i>Lactobacillus species</i>	<i>Lactobacillus acetotolerans</i>	1600	CP044496.1	80.00 %	NA	NA	NA	NA	NA
	<i>Lactobacillus acidophilus</i>	1579	CP139575.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus amylolyticus</i>	83683	CP031835.1	NA	NA	NA	NA	NA	NA
			CP031836.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus apis</i>	303541	CP029476.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus casei</i>	1423732	AP012544.1	NA	NA	NA	NA	NA	NA
			AP012545.1	NA	NA	NA	NA	NA	NA
			AP012546.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus crispatus</i>	47770	CP039266.1	NA	NA	NA	NA	NA	NA
			CP039267.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus delbrueckii</i>	1217420	CP018218.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus gasseri</i>	324831	CP000413.1	80.00 %	NA	NA	NA	NA	NA
	<i>Lactobacillus kullabergensis</i>	1218493	CP029477.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus kefiranofaciens</i>	267818	CP061341.1	80.00 %	NA	NA	NA	NA	NA
			CP061343.1	NA	NA	NA	NA	NA	NA
			CP061342.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus johnsonii</i>	33959	CP059055.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus jensenii</i>	109790	CP018809.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus isalae</i>	2993455	NZ_OX443569.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus iners</i>	147802	CP049230.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus helveticus</i>	1587	CP151471.1	NA	NA	NA	NA	NA	NA

	<i>Lactobacillus helsingborgensis</i>	1218494	CP029544.1	NA	NA	NA	NA	NA	NA
			CP029545.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus xylocopicola</i>	2976676	NZ_AP026803.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus terrae</i>	2269374	CP024610.1	80.00 %					
	<i>Lactobacillus taiwanensis</i>	508451	CP059278.1	NA	NA	NA	NA	NA	NA
			CP059277.1	NA	NA	NA	NA	NA	NA
			CP059276.1	NA	NA	NA	NA	NA	NA
	<i>Lactobacillus paragasseri</i>	2107999	AP018549.1	80.00 %	NA	NA	NA	NA	NA
	<i>Lactobacillus panisapium</i>	2012495	CP048268.1	NA	NA	NA	NA	NA	NA
<i>Lactobacillus mulieris</i>	2508708	CP160088.1	NA	NA	NA	NA	NA	NA	
<i>Escherichia coli</i>		386585	BA000007.3	NA	NA	NA	NA	NA	NA
		471	AB011548.2	NA	NA	NA	NA	NA	NA
<i>Acinetobacter calcoaceticus</i>				AB011549.2	NA	NA	NA	NA	NA
			NZ_CP020000.1	NA	NA	NA	NA	NA	NA
<i>Bacteroides fragilis</i>		272559	CR626927.1	NA	NA	NA	NA	NA	NA
			CR626928.1	NA	NA	NA	NA	NA	NA
<i>Enterococcus faecalis</i>		1351	NZ_AP026714.1	NA	NA	NA	NA	NA	NA
<i>Streptococcus Group C</i>	<i>Streptococcus equi subsp. zooepidemicus</i>	40041	CP065054.1	NA	NA	NA	NA	NA	NA
			CP065055.1	NA	NA	NA	NA	NA	NA
<i>Streptococcus Group G</i>		1320	GQ184567.1	NA	NA	NA	NA	NA	NA
<i>Corynebacterium diphtheriae</i>		1717	NZ_LN831026.1	NA	NA	NA	NA	NA	NA
<i>Neisseria gonorrhoeae</i>		485	AP023069.1	NA	NA	NA	NA	NA	NA
			AP023070.1	NA	NA	NA	NA	NA	NA
<i>Chlamydia trachomatis</i>		272561	NC_000117.1	NA	NA	NA	NA	NA	NA
<i>Mycoplasma pneumoniae</i>		2104	NZ_LR214945.1	NA	NA	NA	NA	NA	NA
			NZ_LR214946.1	NA	NA	NA	NA	NA	NA
<i>Mycoplasma genitalium</i>		2097	CP159789.1	NA	NA	NA	NA	NA	NA
Human papillomavirus (HPV)		10566	NC_027779.1	NA	NA	NA	NA	NA	NA
<i>Trichomonas vaginalis</i>		412133	CM044617.1	NA	NA	NA	NA	NA	NA
			CM044618.1	NA	NA	NA	NA	NA	NA
			CM044619.1	NA	NA	NA	NA	NA	NA
			CM044620.1	NA	NA	NA	NA	NA	NA

		CM044621.1	NA	NA	NA	NA	NA	NA
		CM044622.1	NA	NA	NA	NA	NA	NA
<i>Treponema pallidum</i>	243276	NC_021490.2	NA	NA	NA	NA	NA	NA

Table 18. Percentage homology of G2R forward (FWD), G2R reverse (RV), G2R probe, F3L forward (FWD), F3L reverse (RV) and F3L probe (VIASURE *Monkeypox Virus* Real Time PCR Detection Kit) of the aligned sequences (>80%). NA = not applicable: homology percentage below 80%. * Results at $\geq 80\%$ homology.

Overall, as already seen in Table 15, *in-silico* cross-reactivity was observed for Variola, Cowpox and Camelpox viruses. As explained before, the non-detection of Cowpox virus was verified in a wet analysis conducted with the human-like specimen containing the Cowpox, OPV 88-Lunge (indicated in section “Cross-reactivity wet testing”), and in the case of Variola and Camelpox, low-risk cross-reactivity is expected, since the eradication declaration of Variola by the World Health Assembly in 1980, and the punctual cases of human infection of Camelpox virus, with no reported human-to-human transmission.

Regarding the results presented in Table 18, for those microorganisms that showed $\geq 80\%$ homology with one of the targets design included in the device (e.g., *S. epidermidis*, *S. pyogenes*, *P. aeruginosa*, *L. acetotolerans*, *L. gasseri*, *L. kefirifaciens*, *L. terrae* and *L. paragasseri*), a microbial interference risk assessment was performed to evaluate the potential cross-reactivity of the assay’s primers and probes with non-target organisms. According to the established evaluation criteria, if only one or two oligonucleotides of a primer/probe set (F3L or G2R) show $\geq 80\%$ homology, cross-reactivity can be excluded since PCR amplification is not expected. Likewise, when the potential amplicon length exceeds 500 nucleotides, cross-reactivity risk can also be discarded. Taking these criteria into account, together with the analytical, biological, and epidemiological plausibility assessments, the overall risk of microbial interference was determined to be very low and acceptable.

In conclusion, the VIASURE Monkeypox virus target designs should not cause false positives in detecting this virus when other organisms are present.

Analytical specificity wet testing

Cross-reactivity wet testing

The cross-reactivity of VIASURE *Monkeypox virus* Real Time PCR Detection Kit was confirmed by testing certified reference strains and/or nucleic acid from reference strains associated to skin lesions and/or transmission by close contact to body fluids, negative to the target under study. Additionally, a microbial interference study was conducted using certified reference strains and/or nucleic acid from reference strains associated to skin lesions and/or transmission by close contact to body fluids, negative to the target under study, and in presence of the target under study (reference strain of Monkeypox virus). When possible and when the concentration data were available, the microorganisms were evaluated at medically relevant levels (10^6 CFU/mL for bacteria and 10^5 PFU/mL for viruses).

No cross-reactivity was detected between any of the following microorganisms tested, excepted for the targeted microorganism when present. Likewise, no microbial interference was observed, as no significant modification of the Ct obtained with respect to the control sample was observed.

Cross-reactivity testing					
Control sample (MPXV) *	-/+	Varicella-Zoster Virus strain rubrum	-	<i>Pseudomonas aeruginosa</i> strain ATCC 27853	-
<i>Candida albicans</i> strain NCYC 597	-	Varicella-Zoster Virus strain Ellen	-	<i>Treponema pallidum</i>	-
<i>Enterococcus faecalis</i> strain NCIMB 775	-	Varicella-Zoster Virus strain Ellen / MPXV *	-/+	<i>Trichomonas vaginalis</i>	-
<i>Escherichia coli</i> 0.1285;O18:H7:K1	-	<i>Chlamydia trachomatis</i> [§]	-	Cowpox OPV 88-Lunge	-
Herpes simplex virus 1 (HSV-1) Strain MacIntyre	-	<i>Chlamydia trachomatis</i> [§] / MPXV *	-/+	Vaccinia Virus Elstree	-
Herpes simplex virus 1 (HSV-1) Strain MacIntyre / MPXV *	-/+	<i>Bacteroides fragilis</i>	-	<i>Staphylococcus epidermidis</i> PCI 1200	-
Herpes simplex virus 2 (HSV-2) strain MS	-	Human papillomavirus 16	-	<i>Lactobacillus plantarum</i> strain Lp 39	-
Herpes simplex virus 2 (HSV-2) strain MS / MPXV *	-/+	Human papillomavirus 18	-	<i>Mycoplasma pneumoniae</i> strain PI 1428	-
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> strain ATCC 25923	-	<i>Mycoplasma genitalium</i> M30	-	<i>Streptococcus pyogenes</i> strain MGAS 8232	-
<i>Streptococcus agalactiae</i> strain Z019	-	<i>Neisseria gonorrhoeae</i> strain NCTC 8375 [B 5025]	-	Human Genomic DNA	-
Monkeypox Clade I (Congo Basin clade/Central African isolate)	-/+	Monkeypox Clade IIb (West African clade/2022 Outbreak isolate)	-/+		

Table 19. Reference microorganisms used in this study.

* This reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: Monkeypox Virus, hMPXV/USA/MA001/2022 (Lineage B.1, Clade IIb), NR-58622.

[§] A panel comprising 12 clinical isolates representing serovars associated with trachoma and urogenital infections. *Chlamydia trachomatis*, Serovar A (ATCC VR-571B). *Chlamydia trachomatis*, Serovar B (ATCC VR573). *Chlamydia trachomatis*, Serovar Ba (ATCC VR-347). *Chlamydia trachomatis*, Serovar C (ATCC VR-1477). *Chlamydia trachomatis*, Serovar D (ATCC VR-885). *Chlamydia trachomatis*, Serovar E (ATCC VR-348B). *Chlamydia trachomatis*, Serovar F (ATCC VR-346). *Chlamydia trachomatis*, Serovar G (ATCC VR-878). *Chlamydia trachomatis*, Serovar H (ATCC VR-879). *Chlamydia trachomatis*, Serovar I (ATCC VR-880). *Chlamydia trachomatis*, Serovar J (ATCC VR-886). *Chlamydia trachomatis*, Serovar K (ATCC VR-887).

Interferences and inhibitors of PCR

A study of interfering substances was performed to test the possible interfering and cross-reactive effect of endogenous and exogenous substances on VIASURE *Monkeypox virus* Real Time PCR Detection Kit. A total of 19 potentially interfering substances were evaluated. The samples were extracted with MagDEA Dx SV Kit, using the magLEAD® 12gC instrument, and run on DTprime real-time PCR instrument (DNA-Technology). The rehydrated master mix without any interfering substance was used as a control for the test, to which the nucleic acid extracted from the simulated sample was added (Mixture control or Positive Matrix Control, PMC). In addition to this, the negative matrix without potentially interfering substances nor added microorganism/ reference material was tested as Negative Matrix Control (NMC).

Evaluation of the interference in qPCR assay alone

The following results were obtained when potentially interfering substance was added directly to the rehydrated master mix, which also contains the nucleic acid previously extracted from the simulated samples:

Substance Name	Concentration tested	Result
Acyclovir	6.60x10 ⁻² mg/mL	N.I
Diphenhydramine (Neosayomol cream/Benadryl cream)	0.625 mg/mL	N.I
EDTA whole blood	0.625%	N.I
Petrolatum containing lip/skin care	0.625 mg/mL	N.I
Casein	7 mg/mL	N.I
Cornstarch	2.5 mg/mL	N.I
Faeces	0.22 %	N.I
Female urine	1%	N.I
Benzocaine containing local anaesthetic	0.625 mg/mL	N.I
Human Serum Albumin	10 mg/mL	N.I
Hydrocortisone cream	0.625 mg/mL	N.I
Male urine	1 %	N.I
Mucin	2.5 mg/mL	N.I
Bacitracin	0.625 mg/mL	N.I
Seminal fluid	7 %	N.I
Lidocaine cream	0.65 mg/mL	N.I
Douche	0.625 mg/mL	N.I
Vaginal Lubricant Liquid – water based	2.5 mg/mL	N.I
Zinc Oxide ointment	0.625 mg/mL	N.I

Table 20. Potential interference substances. N.I: No reportable interference.

Evaluation of the interference in the whole process (extraction and the qPCR)

The following results were obtained when potentially interfering substances were added to the simulated sample that was going to be extracted, from which the obtained eluate was added to the qPCR assay:

Substance Name	Concentration tested	Result
Acyclovir	6.60x10 ⁻² mg/mL	N.I
Diphenhydramine (Neosayomol cream/Benadryl cream)	2.5 mg/mL	N.I
EDTA whole blood	5 %	N.I
Petrolatum containing lip/skin care	2.5 mg/mL	N.I
Casein	7 mg/mL	N.I
Cornstarch	2.5 mg/mL	N.I
Faeces	0.22 %	N.I
Female urine	1%	N.I
Benzocaine containing local anaesthetic	2.5 mg/mL	N.I
Human Serum Albumin	10 mg/mL	N.I
Hydrocortisone cream	2.5 mg/mL	N.I
Male urine	1 %	N.I
Mucin	2.5 mg/mL	N.I
Bacitracin	2.5 mg/mL	N.I
Seminal fluid	7 %	N.I
Lidocaine cream	2.6 mg/mL	N.I
Douche	2.5 mg/mL	N.I
Vaginal Lubricant Liquid – water based	2.5 mg/mL	N.I
Zinc Oxide ointment	2.5 mg/mL	N.I

Table 21. Potential interference substances. N.I: No reportable interference.

In conclusion, different potentially interfering substances, both endogenous and exogenous, were tested on VIASURE *Monkeypox virus* Real Time PCR Detection Kit. The results obtained lead to conclude that, at the concentrations tested, no interference of any the evaluated substances is observed at the concentrations indicated for each condition.

12.5.2. Analytical reactivity

Analytical reactivity can be defined as the percentage of target microbial strains or DNA/RNA samples that give the correct positive result. Analytical reactivity was studied *in silico* and wet analyses.

Analytical reactivity *in silico* analysis

The Analytical Reactivity was assessed by using sequences downloaded from NCBI GenBank (<https://www.ncbi.nlm.nih.gov/genbank/>), GISAID EpiPox¹³ database, and an in-house bioinformatic analysis software, in order to demonstrate that MPXV can be correctly detected by the device under study.

The following sequences of the target pathogen were obtained:

- 12,779 sequences of Monkeypox virus, submitted up to December 05th, 2024. These sequences include the following:
 - 620 sequences of Monkeypox virus clade Ia (formerly clade I).
 - 195 sequences of Monkeypox virus clade Ib.
 - 48 sequences of Monkeypox virus clade IIa.
 - 11,499 sequences of Monkeypox virus clade IIb.

The results obtained after the analysis of sequences are shown in the following table:

Target	% sequences experimentally detected without mismatches	N° aligned sequences
MPXV clade Ia	100	540
MPXV clade Ib	100	176
MPXV clade IIa	50	4
MPXV clade IIb	99.32	11,460

Table 22. Analytical Reactivity *in silico* assay.

In conclusion, the analytical reactivity analysis showed correct detection of the MPXV with the VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

¹³ Khare, S., et al (2021) GISAID's Role in Pandemic Response. China CDC Weekly, 3(49): 1049-1051. doi: 10.46234/ccdcw2021.255
PMCID: 8668406

Analytical reactivity wet testing

The analytical reactivity of VIASURE *Monkeypox virus* Real Time PCR Detection Kit was evaluated against the following certified reference strains and material, including from the EQAs programmes, with a positive result:

- NATtrol™ Monkeypox Virus External Run Control, isolate USA/MA001/2022* (ZeptoMetrix Ref: NATMPXV-ERC).

* This reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: Monkeypox Virus, hMPXV/USA/MA001/2022 (Lineage B.1, Clade IIb), NR-58622.

- AMPLIRUN® MONKEYPOX VIRUS DNA CONTROL, West African clade (Vircell, Ref: MBC146-R).
- Monkeypox Clade I (Congo Basin clade/Central African isolate).
- Monkeypox Clade IIb (West African clade/2022 Outbreak isolate).

In conclusion, all the material used in the analytical reactivity study was correctly detected, meeting the established acceptance criteria.

12.6. Metrological traceability

This assay is not designed for measuring purposes.

13. Clinical performance characteristics

The clinical performance of VIASURE *Monkeypox virus* Real Time PCR Detection Kit was tested using clinical samples (mainly lesion swabs) from patients with suspected MPXV infection. In order to determine the clinical diagnostic accuracy, different multicenter evaluations were conducted through collaboration with national and international entities. A summary of the sites, sample type and workflow are included in the following table.

	Site	Sample type	Workflow	Target
1	Hospital Universitario Príncipe de Asturias (HUPA), Madrid, Spain	Lesion swabs in UTM	MagCore HF16 automated extraction system (RBC bioscience) + CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)	MPXV
2	North-West London Pathology (NWLP), UK	Lesion/vesicle swabs in UTM® Universal Transport Medium™ (Copan) *	EZ1 instrument (QIAGEN) + 7500 Fast real time PCR instrument (Applied Biosystems™)	MPXV
3.A	South-West London Pathology (SWLP), UK	Throat and lesion swabs in UTM® Universal Transport Medium™ (Copan) §	QIAGEN EZ1&2 Virus Mini Kit 2.0 (QIAGEN) / EZ1 instrument (QIAGEN) + Rotor-Gene® Q (QIAGEN)	MPXV
3.B		Throat swabs in UTM® Universal Transport Medium™ (Copan) §		
3.C		Lesion swabs in UTM® Universal Transport Medium™ (Copan) §		
4	Viapath Analytics laboratory at St Thomas' Hospital, London, UK	Lesion swabs in UTM® Universal Transport Medium™ (Copan) #	QIAsymphony DSP Virus/Pathogen Mini Kit (QIAGEN) / QIAsymphony SP instrument (QIAGEN) + 7500 Fast real time PCR instrument (Applied Biosystems™)	MPXV

5	Hospital Universitario Miguel Servet (HUMS), Zaragoza, Spain	Ulcerative lesion (rash or growth) swabs in VTM (Non-inactivated) (biocomma®) ~	MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad) and MagCore Viral Nucleic Acid Extraction Kit (Low PCR Inhibition) / MagCore® automated nucleic acid extractors (RBC Bioscience) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	MPXV
6	Certest Biotec S.L, Zaragoza, Spain	Ulcerative lesion swabs in UTM® Universal Transport Medium™ (Copan)	MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	MPXV
7	ONELABT S.A, Santa Elena, Ecuador	Skin lesion/vesicle swabs in UTM® Universal Transport Medium™ (Copan)	QIAamp DNA Mini Kit (QIAGEN) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	MPXV

Table 23. Site, sample type, workflow, and target.

* Previously inactivated with MagNA Pure 96 External Lysis Buffer (Roche Life Science) at 1:1.

\$ Previously mixed with NucliSENS® easyMAG® Lysis Buffer (LB) (bioMérieux) for lysing at 1:1.

Previously lysed by heat treatment with Lysis Buffer (QIAGEN) (200 µL of sample + 430µL of buffer).

~ Previously inactivated with Lysis Buffer Universal LB (Seegene) at 1:1.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values and the likelihood ratios (LR) for VIASURE *Monkeypox virus* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV	LR+	LR-
1	Reference laboratory method *	MPXV	3	6	0	0	1 (0.29-1)	1 (0.54-1)	1 (0.43-1)	1 (0.61-1)	12.25 (0.82-181.65)	0.135 (0.010-1.814)
2	CDC MPX specific assay Li <i>et al</i> 2010 ¹⁴ + Sanger sequencing (only positive samples)	MPXV	39	21	1	0	1 (0.91-1)	0.95 (0.77-0.99)	0.97 (0.86-0.99)	1 (0.83-1)	15.1 (3.2-71.17)	0.013 (0.001-0.21)
3.A	CDC MPX specific assay Li <i>et al</i> 2010 ¹⁴ + Sanger sequencing (only positive samples)	MPXV	20	34	0	0	1 (0.83-1)	1 (0.89-1)	1 (0.79-1)	1 (0.87-1)	68.3 (4.35-1,071)	0.024 (0.002-0.317)
3.B			1	3	0	0	1 (0.02-1)	1 (0.29-1)	1 (0.05-1)	1 (0.3-1)	6 (0.39-90.4)	0.28 (0.01-0.377)
3.C			19	31	0	0	1 (0.82-1)	1 (0.88-1)	1 (0.79-1)	1 (0.86-1)	62.4 (3.98-977)	0.025 (0.002-0.39)
4	CDC MPX specific assay Li <i>et al</i> 2010 ¹⁴ + Sanger sequencing (only positive samples)	MPXV \$	42	13	2	0	1 (0.91-1)	0.86 (0.59-0.98)	0.95 (0.84-0.99)	1 (0.75-1)	6.32 (2.02-19.76)	0.014 (0.001-0.21)
5 ~	RealStar Orthopoxvirus PCR kit 1.0 (Altona Diagnostics)	MPXV	69	96	0	0	1 (0.97-1)	1 (0.98-1)	1 (0.93-1)	1 (0.95-1)	192.61 (12.133-3057.7)	0.007 (0.000-0.114)
6	RealStar Orthopoxvirus PCR kit 1.0 (Altona Diagnostics)	MPXV	38	348	1	0	1(0.9-1)	0.99(0.98-1)	0.97(0.86-0.99)	1(0.98-1)	230(46.6-1137)	0.013(0.001-0.2)
7	WHO MPX	MPXV	102	160	0	1	0.99	1	1	0.99	317.3	0.01 (0.003-

¹⁴ Real-time PCR assays for the specific detection of monkeypox virus West African and Congo Basin strain DNA. Yu Li, Hui Zhao, Kimberly Wilkins, Christine Hughes, Inger K. Damon. Journal of Virological Methods. Volume 169, Issue 1, October 2010, Pages 223-227.

	molecular assay (Li et al 2010) ¹⁴					(0.94-1)	(0.97-1)	(0.96-1)	(0.96-1)	(19.9-5052)	0.071)
	Monkeypox Virus (MPXV) Fast Real Time PCR Kit (BioPerfectus)	102	161	0	0	1 (0.96-1)	1 (0.97-1)	1 (0.97-1)	1 (0.97-1)	322.4 (20.28-5133)	0.005 (0-0.077)
	CDC MPX specific assay Li et al 2010 ¹⁴	102	161	0	0	1 (0.96-1)	1 (0.97-1)	1 (0.97-1)	1 (0.97-1)	322.4 (20.28-5133)	0.005 (0-0.077)

Table 24. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Positive Predictive Values (PPV), Negative Predictive Values (NPV) and the likelihood ratios (LR) for VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

* Two orthopoxvirus-generic PCRs and another MPXV-specific PCR, based on three publications^{15,16,17} + Sanger sequencing (for confirmation of MPXV positive samples) + Allplex™ Genital Ulcer Assay (Seegene) (only with MPXV negative samples).

\$ Kappa statistic value of 0.9055.

~ The results were identical between the two workflows, verifying that the device test has excellent performance as an open format product.

Results show high agreement to detect MPXV using VIASURE *Monkeypox virus* Real Time PCR Detection Kit. In addition, the use of other types of samples in some clinical evaluations helped to observe that in cases where vesicle swab samples were unavailable, certain alternative sample types showed comparable positive results when the viral load was sufficient.

¹⁵ Loparev, V. N., Massung, R. F., Esposito, J. J., & Meyer, H. (2001). Detection and differentiation of old world orthopoxviruses: restriction fragment length polymorphism of the crmB gene region. *Journal of clinical microbiology*, 39(1), 94–100. <https://doi.org/10.1128/JCM.39.1.94-100.2001>

¹⁶ Li, Y., Olson, V. A., Laue, T., Laker, M. T., & Damon, I. K. (2006). Detection of monkeypox virus with real-time PCR assays. *Journal of clinical virology: the official publication of the Pan American Society for Clinical Virology*, 36(3), 194–203. <https://doi.org/10.1016/j.jcv.2006.03.012>

¹⁷ Kulesh, D. A., Loveless, B. M., Norwood, D., Garrison, J., Whitehouse, C. A., Hartmann, C., Mucker, E., Miller, D., Wasieleski, L. P., Jr, Huggins, J., Huhn, G., Miser, L. L., Imig, C., Martinez, M., Larsen, T., Rossi, C. A., & Ludwig, G. V. (2004). Monkeypox virus detection in rodents using real-time 3'-minor groove binder TaqMan assays on the Roche LightCycler. *Laboratory investigation; a journal of technical methods and pathology*, 84(9), 1200–1208. <https://doi.org/10.1038/labinvest.3700143>

ANNEX 1

OPEN FORMAT AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MPX101L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MPX101H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MPX106L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MPX106H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MPX112L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MPX112H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MPX113L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MPX113H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MPX101
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MPX136
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MPX172

Table A1 1. References

A1.1 Principle of the procedure

VIASURE *Monkeypox virus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilised¹ format. The assay uses a human housekeeping gene (haemoglobin subunit beta gene, *HBB* gene) as an **Endogenous Internal Control (EIC)** to follow-up the integrity of the sample, to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
Monkeypox virus (MPXV)	FAM	<i>G2R</i> and <i>F3L</i> genes
Endogenous Internal Control (EIC)	HEX*	<i>HBB</i> gene

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel; to check most common detection channels consult the website www.certtest.es.

A1.2 Reagents provided

VIASURE *Monkeypox virus* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilised PCR Reaction-Mix could be placed inside different wells and could be

¹ Please note that both “stabilised”, and “lyophilised” terms are used indistinguishable and as synonyms throughout the entire document.

marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible instruments. Table A1.4 includes materials and reagents to be used with 96-well plate compatible instruments. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on Certest's website www.certest.es).

Reagent/Material	Description	Concentration Range	Colour	Amount
<i>Monkeypox virus</i> 8-well strips	Lyoprotectors and Stabilizers	±6 g/100 mL*	White (Opaque)	1/6/12 x 8-well strip
	Nucleotide triphosphate (dNTPs)	±1 mM*		
	Primers and Probes	0.2-1 nMol/μL*		
	Enzymes	10-100 U/rxn*		
Rehydration Buffer	Saline Solution Mixture	±13 mM	Blue	1 vial x 1.8 mL
	Buffer (TRIS, pH)	±67 mM		
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilised ¹ DNA	1.9x10 ⁴ copies/μL*	Red	1 vial
Negative control	RNAse/DNAse free water	1 g/mL	Violet	1 vial x 1 mL
Water RNAse/DNAse free	RNAse/DNAse free water	1 g/mL	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	-	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref. VS-MPX101L, VS-MPX101H, VS-MPX106L, VS-MPX106H, VS-MPX112L and VS-MPX112H.

* For component in stabilised format, the concentration range means after rehydration.

Reagent/Material	Description	Concentration Range	Colour	Amount
<i>Monkeypox virus</i> 96-well plate	Lyoprotectors and Stabilizers	±6 g/100 mL*	White (Opaque)	1 plate
	Nucleotide triphosphate (dNTPs)	±1 mM*		
	Primers and Probes	0.2-1 nMol/μL*		
	Enzymes	10-100 U/rxn*		
Rehydration Buffer	Saline Solution Mixture	±13 mM	Blue	1 vial x 1.8 mL
	Buffer (TRIS, pH)	±67 mM		
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilised ¹ DNA	1.9x10 ⁴ copies/μL*	Red	1 vial
Negative control	RNAse/DNAse free water	1 g/mL	Violet	1 vial x 1 mL
Water RNAse/DNAse free	RNAse/DNAse free water	1 g/mL	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	-	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref VS-MPX113L and VS-MPX113H.

* For component in stabilised format, the concentration range means after rehydration.

Reagent/Material	Description	Concentration Range	Colour	Amount
<i>Monkeypox virus</i> 4-well strips	Lyoprotectors and Stabilizers	±6 g/100 mL*	Transparent	2/9/18 x 4-well strip
	Nucleotide triphosphate (dNTPs)	±1 mM*		
	Primers and Probes	0.2-1 nMol/μL*		
	Enzymes	10-100 U/rxn*		
Rehydration Buffer	Saline Solution Mixture	±13 mM	Blue	1 vial x 1.8 mL
	Buffer (TRIS, pH)	±67 mM		
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilised ¹ DNA	1.9x10 ⁴ copies/μL*	Red	1 vial
Negative control	RNAse/DNAse free water	1 g/mL	Violet	1 vial x 1 mL
Water RNAse/DNAse free	RNAse/DNAse free water	1 g/mL	White	1 vial x 1 mL
4-cap strip	Optical caps for sealing plate during thermal cycling	-	Transparent	2/9/18 x 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref. VS-MPX101, VS-MPX136 and VS-MPX172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

* For component in stabilised format, the concentration range means after rehydration.

A1.3 Test procedure

A1.3.1 Lyophilised positive control

Monkeypox virus Positive Control contains high copies of the template; the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilised *Monkeypox virus* Positive Control (red vial) by adding 100 μL of the supplied Water RNAse/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it into aliquots to minimize freeze and thaw cycles. For instance, Positive Control may be divided into 4 aliquots of 25 μL each one, but end user should choose the number of aliquots to prepare considering the maximum freeze/thaw cycles verified. Positive control has been validated as still being stable after 6 freeze-thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay.

Please note: In case the number of reactions required/needed is less than supplied on one strip or plate, before removing the protective seal, carefully cut out the required wells and store the rest inside the bag with the desiccant. Peel off protective aluminium seal from plates or strips you need for the assay.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Monkeypox virus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

3) Set up the thermocycler (consult thermocycler compatibility on Certest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (MPXV) and HEX (Endogenous Internal Control (EIC)). Depending on the equipment used, select in the thermocycler only the proper channels where targets are detected (to check most common detection channels consult website www.certest.es).

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MPX196T

Table A2. 1.References.

A2.1 Principle of the procedure

VIASURE *Monkeypox virus* Real Time PCR Detection Kit contains in Reaction-Mix tube well all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay uses a human housekeeping gene (haemoglobin subunit beta gene, *HBB* gene) as an **Endogenous Internal Control (EIC)** to follow-up the integrity of the sample, to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
Monkeypox virus (MPXV)	FAM	<i>G2R</i> and <i>F3L</i> genes
Endogenous Internal Control (EIC)	HEX*	<i>HBB</i> gene

Table A2. 2.Target, channel and genes.

*Depending on the equipment used select the proper detection channel; to check most common detection channels consult website www.certtest.es.

A2.2 Reagents provided

VIASURE *Monkeypox virus* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Concentration Range	Colour	Amount
<i>Monkeypox virus</i> Reaction-Mix tube	Lyoprotectors and Stabilizers	±6 g/100 mL*	White (Opaque)	4 vials
	Nucleotide triphosphate (dNTPs)	±1 mM*		
	Primers and Probes	0.2-1 nMol/μL*		
	Enzymes	10-100 U/rxn*		
Rehydration Buffer	Saline Solution Mixture	±13 mM	Blue	1 vial x 1.8 mL
	Buffer (TRIS, pH)	±67 mM		
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilised ¹ DNA	1.9x10 ⁴ copies/μL*	Red	1 vial
Negative control	RNAse/DNAse free water	1 g/mL	Violet	1 vial x 1 mL
Water RNAse/DNAse free	RNAse/DNAse free water	1 g/mL	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref. VS-MPX196T.

* For component in stabilised format, the concentration range means after rehydration.

A2.3 Test procedure

A2.3.1 Lyophilised positive control

Monkeypox virus Positive Control contains high copies of the template; the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilised *Monkeypox virus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it into aliquots to minimize freeze and thaw cycles. For instance, Positive Control may be divided into 4 aliquots of 25µL each one, but end user should choose the number of aliquots to prepare considering the maximum freeze/thaw cycles verified. Positive control has been validated as still being stable after 6 freeze-thaw cycles.

A2.3.2 Lyophilised Reaction-Mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Monkeypox virus* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilised Reaction-Mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 4°C±2°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Monkeypox virus* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Monkeypox virus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on Certest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (MPXV) and HEX (Endogenous Internal Control (EIC)). Depending on the equipment used, select in the thermocycler only the proper channels where targets are detected (to check most common detection channels consult website www.certest.es).

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HQ328780.1 Cowpox virus isolate IT-Llama09 tumor necrosis factor receptor II-like protein (crmB) gene, partial cds	92,31	83,33	95,83	90
KC813491.1 Cowpox virus strain BeaBer04/1, complete genome	92,31	83,33	95,83	90
KC813492.1 Cowpox virus strain HumLan08/1, complete genome	92,31	83,33	95,83	90
KC813493.1 Cowpox virus strain HumLit08/1, complete genome	96,15	83,33	95,83	91,25
KC813494.1 Cowpox virus strain HumLue09/1, complete genome	96,15	83,33	95,83	91,25
KC813495.1 Cowpox virus strain HumMag07/1, complete genome	92,31	83,33	95,83	90
KC813496.1 Cowpox virus strain HumPad07/1, complete genome	92,31	83,33	95,83	90
KC813497.1 Cowpox virus strain JagKre08/1, complete genome	92,31	83,33	95,83	90
KC813498.1 Cowpox virus strain JagKre08/2, complete genome	92,31	83,33	95,83	90
KC813503.1 Cowpox virus strain RatGer09/1, complete genome	92,31	83,33	95,83	90
KC813505.1 Cowpox virus strain RatKre08/2, complete genome	92,31	83,33	95,83	90
KC813506.1 Cowpox virus strain CatPot07/1, complete genome	92,31	83,33	95,83	90
KC813509.1 Cowpox virus strain HumBer07/1, complete genome	92,31	83,33	95,83	90
KC813510.1 Cowpox virus strain HumGra07/1, complete genome	92,31	83,33	95,83	90
KC841319.1 Camelpox virus isolate CMLV-DELHI tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92,31	86,67	95,83	91,25
KC841320.1 Camelpox virus isolate CMLV-BIKANER tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92,31	86,67	95,83	91,25
KC841321.1 Camelpox virus isolate CMLV-BARMER tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92,31	86,67	95,83	91,25
KC841322.1 Camelpox virus isolate CMLV-JAISALMER tumor necrosis factor receptor II-like protein (crmB) gene, complete cds	92,31	86,67	95,83	91,25
KU861504.1 Cowpox virus isolate Kostroma/2015 cytokine response modifier B gene, complete cds	92,31	86,67	95,83	91,25
KY358055.1 Variola virus strain VD21, 17th century, complete genome	88,46	86,67	100	91,25
KY549146.1 Cowpox virus isolate CPXV_K779, complete genome	92,31	83,33	95,83	90
KY549147.1 Cowpox virus isolate CPXV_K780, complete genome	92,31	83,33	95,83	90
KY549151.1 Cowpox virus isolate Norwayfeline, complete genome	92,31	83,33	95,83	90
KY569019.1 Cowpox virus strain CPXV CheNova_DK_2014, complete genome	96,15	83,33	95,83	91,25
L22579.1 Variola major virus (strain Bangladesh-1975) complete genome	88,46	86,67	100	91,25
LR800244.1 Variola virus genome assembly, complete genome: monopartite	92,31	86,67	100	92,5
LR800245.1 Variola virus genome assembly, complete genome: monopartite	92,31	86,67	100	92,5
LR800246.1 Variola virus genome assembly, complete genome: monopartite	88,46	86,67	95,83	90
LR800247.1 Variola virus genome assembly, complete genome: monopartite	88,46	86,67	95,83	90
LT706528.1 Variola virus isolate VARV_V563 genome assembly, chromosome: I	88,46	86,67	100	91,25
LT706529.1 Variola virus isolate VARV_V1588 genome assembly, chromosome: I	88,46	86,67	100	91,25
MN240300.1 Borealpox virus, complete genome	88,46	80	87,5	85
MN369532.1 Vaccinia virus isolate VK01, complete genome	92,31	83,33	95,83	90
NC_001611.1 Variola virus, complete genome	88,46	86,67	100	91,25
OL468961.1 Variola virus strain GHA68, complete genome	88,46	86,67	100	91,25
OL468962.1 Variola virus strain Congo 1970 v70-60, complete genome	88,46	86,67	100	91,25

U18339.1 Variola virus Garcia-1966 right near-terminal region	88,46	86,67	100	91,25
U18341.1 Variola virus Somalia-1977 right near-terminal region	88,46	86,67	100	91,25
U86874.1 Taterapox virus tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92,31	86,67	91,67	90
U87837.1 Camelpox virus strain Somalia-1978 tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92,31	86,67	95,83	91,25
U87838.1 Camelpox virus CP1 tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92,31	86,67	95,83	91,25
U87839.1 Camelpox virus strain Saudi-M3, tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92,31	86,67	95,83	91,25
U87840.1 Camelpox virus CP5 tumor necrosis factor receptor II homolog (crmB) gene, complete cds	92,31	86,67	95,83	91,25
U88145.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U88146.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U88147.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U88148.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U88150.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U88151.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U88152.1 Variola virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	88,46	86,67	100	91,25
U90225.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92,31	83,33	95,83	90
U90226.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92,31	83,33	95,83	90
U90227.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92,31	83,33	95,83	90
U90228.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96,15	83,33	95,83	91,25
U90229.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92,31	83,33	95,83	90
U90230.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92,31	83,33	95,83	90
U90231.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96,15	83,33	95,83	91,25
U90232.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	92,31	83,33	95,83	90
U90233.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96,15	83,33	95,83	91,25
U90234.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96,15	83,33	95,83	91,25
U90235.1 Cowpox virus tumor necrosis factor receptor II homolog (CrmB) gene, complete cds	96,15	83,33	95,83	91,25
X67117.1 Variola virus (Xhol-F,O,H,P,Q genome fragment) genes	88,46	86,67	100	91,25
X69198.1 Variola virus DNA complete genome	88,46	86,67	100	91,25
X70841.1 Variola virus genes for ORF1R, ORF2R and ORF3R	88,46	86,67	100	91,25
Y16780.1 variola minor virus complete genome	88,46	86,67	100	91,25

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PRODUCTO	REFERENCIA
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MPX101L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MPX101H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MPX106L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MPX106H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MPX112L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MPX112H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MPX113L
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MPX113H
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MPX101
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MPX136
VIASURE <i>Monkeypox virus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MPX172

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Bibliography/Bibliografía

- Bawage, S. S., Tiwari, P. M., Pillai, S., Dennis, V., & Singh, S. R. (2013). Recent advances in diagnosis, prevention, and treatment of human respiratory syncytial virus. In *Advances in Virology* (Vol. 2013). <https://doi.org/10.1155/2013/595768>
- Chauhan, R. P., Fogel, R., & Limson, J. (2023). Overview of Diagnostic Methods, Disease Prevalence and Transmission of Mpox (Formerly Monkeypox) in Humans and Animal Reservoirs. In *Microorganisms* (Vol. 11, Issue 5). MDPI. <https://doi.org/10.3390/microorganisms11051186>
- French, C. E., McKenzie, B. C., Coope, C., Rajanaidu, S., Paranthaman, K., Pebody, R., Nguyen-Van-Tam, J. S., Higgins, J. P. T., & Beck, C. R. (2016). Risk of nosocomial respiratory syncytial virus infection and effectiveness of control measures to prevent transmission events: a systematic review. *Influenza and Other Respiratory Viruses*, 10(4), 268–290. <https://doi.org/10.1111/irv.12379>
- Kuo, R. L., Yang, S. L., Liu, Y. C., Chen, L. T., Mok, C. K., Kuo, S. M., Shih, S. R., & Tsao, K. C. (2014). Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. *Journal of Virological Methods*, 208, 41–46. <https://doi.org/10.1016/j.jviromet.2014.07.020>
- Lu, J., Xing, H., Wang, C., Tang, M., Wu, C., Ye, F., Yin, L., Yang, Y., Tan, W., & Shen, L. (2023). Mpox (formerly monkeypox): pathogenesis, prevention, and treatment. In *Signal Transduction and Targeted Therapy* (Vol. 8, Issue 1). Springer Nature. <https://doi.org/10.1038/s41392-023-01675-2>
- Neumann, G., & Kawaoka, Y. (2015). Transmission of influenza A viruses. *Virology*, 479–480, 234–246. <https://doi.org/10.1016/j.virol.2015.03.009>
- Paul Glezen, W., Schmier, J. K., Kuehn, C. M., Ryan, K. J., & Oxford, J. (2013). The Burden of Influenza B: A Structured Literature Review. *American Journal of Public Health*, 103(3), e43–e51. <http://www.ajph.org>
- Qudus, M. S., Cui, X., Tian, M., Afaq, U., Sajid, M., Qureshi, S., Liu, S., Ma, J., Wang, G., Faraz, M., Sadia, H., Wu, K., & Zhu, C. (2023). The prospective outcome of the monkeypox outbreak in 2022 and characterization of monkeypox disease immunobiology. *Frontiers in Cellular and Infection Microbiology*, 13. <https://doi.org/10.3389/fcimb.2023.1196699>
- WHO. (2024). *Monkeypox virus*. <https://www.who.int/emergencies/disease-outbreak-news/item/2022-DON385>

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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certest

F-566 rev.03

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.
The products, services and data set out in this document may suffer changes
and/or variations on the texts and pictures shown.