## WHO Emergency Use Assessment SARS-CoV-2 IVDs PUBLIC REPORT

**Product: STANDARD Q COVID-19 Ag Test**

**Manufacturer: SD Biosensor, Inc.**

**EUL Number: EUL 0563-117-00**

**Outcome: Accepted**

The EUL process is intended to expedite the availability of in vitro diagnostics needed in public health emergency situations and to assist interested UN procurement agencies and the Member States in determining the acceptability of using specific products in the context of a Public Health Emergency of International Concern (PHEIC), based on an essential set of available quality, safety and performance data. The EUL procedure includes the following:

- Quality Management Systems Review and Plan for Post-Market Surveillance: desk-top review of the manufacturer’s Quality Management System documentation and specific manufacturing documents;
- Product Dossier Review: assessment of the documentary evidence of safety and performance.

STANDARD Q COVID-19 Ag Test, product codes 09COV30D, 09COV31D, 09COV32D, 09COV33D and 10COVC11 (COVID-19 Ag Control swab), CE-mark regulatory version, manufactured by SD Biosensor, Inc., C 4th and 5th, 16 Deogyong-daero, 1556 beon-gil Suwon-si, Geonggi-do, 16690, Republic of Korea, was listed on 22 September 2020.

### Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the EUL product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the table below.

Version	Summary of amendment	Date of report amendment
2.0	Addition of a warning that STANDARD COVID-19 Ag Control (product code:10COVC10) was not assessed with STANDARD Q COVID-19 Ag Test and is not part of EUL.	6 November 2020
3.0	Correction of the warning added in version to read as, <i>“STANDARD COVID-19 Ag Control (product code:10COVC10) was assessed and found not acceptable.”</i>	18 November 2020
4.0	Change the outside packaging box in the public report to align with the approved IFU.	15 January 2021
5.0	1. Changed the shape of the plastic container cover for the nozzle cap (Blister - flat).	29 June 2022

	<ol style="list-style-type: none"> <li>2. Addition of a new buffer tube rack.</li> <li>3. Addition of a new manufacturing site.</li> <li>4. Addition of a new warehousing site.</li> <li>5. Addition of new suppliers for components.</li> <li>6. Change the Site plan of BioNote.</li> <li>7. Introduction of new quality control swabs.</li> <li>8. Addition of a specimen type (Nasal swab).</li> <li>9. Change of the product IFU and addition of new IFUs.</li> <li>10. Change of the product package and addition of new packages.</li> </ol>	
--	---	--

### Intended use

According to the claim of intended use from SD Biosensor, Inc., “*STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens of SARS-CoV-2 present in human nasopharyngeal specimens. This product is intended for healthcare professionals at the clinical setup and point of care sites, as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms of SARS-CoV-2 infection. It provides only an initial screening test result. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.*”

### Specimen types that were validated

Nasopharyngeal and nasal swab specimens (please pay attention to the product codes).

### Test kit contents:

Component	25 tests (product code 09COV30D)	25 tests (product code 09COV31D)	25 tests (product code 09COV32D)	25 tests (product code 09COV33D)	10 tests/kit COVID-19 Ag Control swab (product code 10COVC11)
Test device (individually in a foil pouch with desiccant)	25	25	25	25	N/A
Extraction buffer tube	25	25	25	25	N/A
Nozzle cap	25	25	25	25	N/A
Sterile swab	25 (Nasopharyngeal)	25 (Nasal)	25 (Nasopharyngeal)	25 (Nasal)	N/A
Negative and Positive control swabs	N/A	N/A	1 of each	1 of each	10 of each
Instructions for use	1	1	1	1	1

### Items required but not provided

- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
- Timer
- Biohazard container
- STANDARD F Analyzer

### Storage

2-30°C.

### Shelf-life upon manufacture

24 months for STANDARD Q COVID-19 Ag Test (real-time stability studies are ongoing).

30 months for STANDARD COVID-19 Ag Control Swab.

### Warnings/limitations:

STANDARD COVID-19 Ag Control (product code:10COVC10) was assessed and found not acceptable.

## Product dossier assessment

SD Biosensor, Inc. submitted a product dossier for the STANDARD Q COVID-19 Ag Test for detecting severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) as per the *“Instructions for Submission Requirements: In vitro diagnostics (IVDs) Detecting SARS-CoV-2 Nucleic Acid and rapid diagnostics tests detecting SARS-CoV-2 antigens (PQDx\_0347 version 4)”*. The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

### Post listing commitment for EUL

As a commitment to listing, the manufacturer is required to:

1. Provide the real-time stability studies report by 31 July 2022. The commitment is under review.
2. Submit the final transport and real-time stability studies report of the new quality control swab as part of the next EUL renewal application.
3. Submit nasal swab storage stability studies at -20 °C for 4 months and -70 °C for 6 months as part of the next EUL renewal application.
4. Submit the results of the accelerated stability testing up to week 20 as part of the next EUL renewal application.
5. Amend the IFU’s font size to make it more user-friendly in the following IFU review. Submit the amended version of the IFU as part of the next EUL renewal application.

Risk-benefit assessment is acceptable.

## Quality Management Systems Review

To establish eligibility for WHO procurement, SD Biosensor, Inc. was asked to provide up-to-date information about the status of its quality management system.

Based on the review of the submitted quality management system documentation by WHO staff, it was established that sufficient information was provided by SD Biosensor, Inc. to fulfil the requirements described in the *“Instructions for Submission Requirements: In vitro diagnostics (IVDs) Detecting SARS-CoV-2 Nucleic Acid and rapid diagnostics tests detecting SARS-CoV-2 antigens (PQDx\_347 version 4)”*.

Quality management documentation assessment is acceptable.

## Plan for Post-Market Surveillance

Post-market surveillance, including monitoring all customer feedback, detecting and acting on adverse events, product problems, non-conforming goods and processes is a critical component of minimizing potential harm of an IVD listed for emergency use.

The following post-EUL activities are required to maintain the EUL listing status:

1. Notification to WHO of any planned changes to a EUL product, in accordance with *“WHO procedure for changes to a WHO prequalified in vitro diagnostic”* (document number PQDx\_121); and
2. Post-market surveillance activities, in accordance with *“Guidance for post-market surveillance and market surveillance of medical devices, including in vitro diagnostics”* (ISBN 978-92-4-001531-9).

SD Biosensor, Inc. is also required to submit an annual report that details sales data and all categories of complaints in a summarized form. There are certain categories of complaints and changes to the product that must be notified immediately to WHO, as per the above-mentioned documents.

The manufacturer has committed to ensuring that post-emergency use listing safety, quality and performance monitoring activities are in place in accordance with WHO guidance *“WHO guidance on post-market surveillance of in vitro diagnostics”*.<sup>1</sup>

## Scope and duration of procurement eligibility

STANDARD Q COVID-19 Ag Test, product codes 09COV30D, 09COV31D, 09COV32D, 09COV33D and 10COVC11 (COVID-19 Ag Control swab), manufactured by SD Biosensor, Inc. is considered to be eligible for WHO procurement until further notice. The assay may be used

---

<sup>1</sup> Available on the web page

<https://www.who.int/publications/i/item/guidance-for-post-market-surveillance-and-market-surveillance-of-medical-devices-including-in-vitro-diagnostics>

to detect the Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) antigens. This listing does not infer that the product meets WHO prequalification requirements and does not mean that the product is listed as WHO prequalified.

As part of the ongoing requirements for listing as eligible for WHO procurement, SD Biosensor, Inc. must engage in post-market surveillance activities to ensure that the product continues to meet safety, quality and performance requirements. SD Biosensor, Inc. is required to notify WHO of any complaints, including adverse events related to the use of the product, within 7 days and any changes to the product.

WHO reserves the right to rescind eligibility for WHO procurement if additional information on the safety, quality, and performance during post-market surveillance activities and if new data becomes available to WHO that changes the risk-benefit balance.

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


**STANDARD**

**1.1.2 Product code 09COV31D**

- 1) Test device (individually in a foil pouch with desiccant) x 25
- 2) Desiccant buffer table x 25
- 3) Nasal swab x 25
- 4) Storage media x 25
- 5) Instructions for use x 1

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
SD Biosensor, Inc. is a subsidiary of SD Biosensor GmbH, a company registered in Germany.

SD Biosensor GmbH  
Promed Consulting GmbH  
Friedrichstraße 100, 10117 Berlin, Germany  
Tel: +49 30 2500 1000  
www.sdbiosensor.com

**1.1.3 Product code 09COV32D**

7) Instructions for use x 1  
8) STANDARD COVID-19 Ag Negative Control swab x 1  
9) STANDARD COVID-19 Ag Positive Control swab x 1  
10) Device swab x 25  
11) Instructions for use x 25  
12) Standard buffer x 25  
13) Test device (individually in a foil pouch with desiccant) x 25

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25 TESTS  
K I T

## STANDARD Q COVID-19 Ag Test

Manufactured by  
**SD Biosensor, Inc.**  
10000 E. Harvard Ave., Suite 100  
Denver, CO 80231, USA  
www.sdbiosensor.com  
Manufacturing site:  
10000 E. Harvard Ave., Suite 100  
Denver, CO 80231, USA  
Authorized Representative:  
**WT Promed Consulting GmbH**  
Wittmannstraße 10, 10585 Berlin, Germany  
The kit is used only for the control and not for the final  
diagnosis of COVID-19. For more information, please  
contact us at [info@sdbiosensor.com](mailto:info@sdbiosensor.com)

SDCOV19 Ag KIT  
Issue date: 2022.09

#### **1.1.4 Product code 09COV33D**

www.sdbiosensor.com

SD BIOSENSOR

CE, RoHS, REACH, ISO 13485, ISO 9001

7) Instructions for use x 1  
 8) STANDARD COVID-19 Ag Negative Control Swab x 1  
 9) STANDARD COVID-19 Ag Positive Control Swab x 1  
 10) Nasal Swab x 25  
 11) Media (Individually in a pouch with desiccant) x 25  
 12) Expansion buffer (x 25)

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www.sdbiosensor.com

RESOURCES: WHO  
 Issue Date: 2020.05

Authorized by SD Biosensor, Inc.  
 SD Biosensor, Inc. is a U.S. company with a U.S. FDA  
 registration number: 1412227. SD Biosensor, Inc.  
 is a U.S. company with a U.S. FDA registration  
 number: 1412227. SD Biosensor, Inc. is a U.S.  
 company with a U.S. FDA registration number:  
 1412227. SD Biosensor, Inc. is a U.S. company  
 with a U.S. FDA registration number: 1412227.


Authorized by SD Biosensor, Inc.  
 SD Biosensor, Inc. is a U.S. company with a U.S. FDA  
 registration number: 1412227. SD Biosensor, Inc.  
 is a U.S. company with a U.S. FDA registration  
 number: 1412227. SD Biosensor, Inc. is a U.S.  
 company with a U.S. FDA registration number:  
 1412227. SD Biosensor, Inc. is a U.S. company  
 with a U.S. FDA registration number: 1412227.

**1.1.5 Product code 10COVC11 (COVID-19 Ag Control swab)**

SD BIOSENSOR

www.sdbiosensor.com

CONTENTS  
1) STANDARD COVID-19 Ag Positive Control swab x 10  
2) STANDARD COVID-19 Ag Negative Control swab x 10  
3) Instructions for Use x 1

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## STANDARD COVID-19 Ag Control swab

Manufactured by  
**SD Biosensor, Inc.**

Head office  
C-4thB3B, 16, Deogyong-dong 153Bonsil, Yeongdeung-gu,  
Seoul 0, Ganggu-do 16220, REPUBLIC OF KOREA


Manufacturing site  
74, Donggansongmyeong 4-ro, Ojeong-eup, Heungdeok-gu,  
Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Authorized Representative  
**MT Promed Consulting GmbH**  
Altenhofstrasse 93 66386 St. Ingbert Germany  
Phone : +49 6894 581020, Fax : +49 6894 581021

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Classification : Tool for biopsy

Product-License No : 14-2811

This product is a class II medical device


FDA Registered and CE Approved

Peel off the plastic film

Remove the protective tube

Please check the breakpoint

STERILITY GUARANTEED UNLESS PACKAGE HAS BEEN OPENED OR DAMAGED.  
REMOVED THE PROTECTIVE TUBE OF THE SWABS AND USE (NFS-1, NFS-2, NFS-3, NFS-4, NFS-5, NFS-6)  
DO NOT EXCESSIVE FORCE THE PRODUCT AT THE TIME OF SAMPLE COLLECTION (ESPECIALLY PRODUCT  
WITH BREAK POINT)  
THIS PRODUCT IS A CLASS II MEDICAL DEVICE LICENSED IN THE REPUBLIC OF KOREA  
MDD CE 2292 CE-Marking according to directive 93/42/EEC.

R. = REF L. = LOT E. = ⓧ

STERILE R

Registered to ISO 13485 & GMP  
MDD CE 2292

EC REP S.B. PHARMA GMBH

Bunsenstr. 14 53121 Bonn, Germany  
TEL: +49 (0) 228 52266916 FAX: +49 (0) 228 52266918

주식회사 노블바이오  
**Noble Biosciences, Inc.**  
Clinical Diagnostic Products

경기도 화성시 정남면 신백길 13-50  
13-50, Sinbaek-gil, Jeongnam-myeon, Hwaseong-si,  
Gyeonggi-do, 18521, Republic of Korea  
www.noblebio.co.kr E-mail : info@noblebio.co.kr  
TEL : (82)-31-291-0044 FAX : (82)-31-291-0046

**Noble Bio**  
Clinical Diagnostic Products

**PEEL HERE**


**NFS-SWAB APPLICATOR™**  
STERILE NYLON (FLOCKED) SWABS

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
Classification : Tool for biopsy  
Product-License No : 14-2811  
This product is a class II medical device  
FDA Registered and CE Approved

Peel off the plastic film

Remove the protective tube

Please check the breakpoint

MOLDED BREAKPOINT

STERILITY GUARANTEED UNLESS PACKAGE HAS BEEN OPENED OR DAMAGED  
REMOVE THE PROTECTIVE TUBE OF THE SWABS AND USE (NFS-1, NFS-2, NFS-3, NFS-4, NFS-5, NFS-6)  
DO NOT EXCESSIVE FORCE THE PRODUCT AT THE TIME OF SAMPLE COLLECTION (ESPECIALLY PRODUCT WITH BREAK POINT)  
THIS PRODUCT IS A CLASS II MEDICAL DEVICE LICENSED IN THE REPUBLIC OF KOREA

MDD CE 2252 CE-Marking according to directive 93/42/EEC.

R.= [REF] L.= [LOT] E.= [E]

R. NFS-1

L. NFS-XXXXX-XX

E. YYYY/MM

STERILE R

Registered to ISO 13485 & GMP  
MDD CE 2292

EC REP S.B. PHARMA GMBH

Bunsenstr. 14 53121 Bonn, Germany  
TEL: +49 (0) 228 52266916 FAX: +49 (0) 228 52266918



주식회사 노블바이오  
Noble Biosciences, Inc.  
Clinical Diagnostic Products

경기도 화성시 정남면 신백길 13-50  
13-50, Siribaek-gil, Jeongnam-myeon, Hwaseong-si,  
Gyeonggi-do, 18521, Republic of Korea  
www.noblebio.co.kr E-mail: info@noblebio.co.kr  
TEL: (82)-31-291-0044 FAX: (82)-31-291-0046

VP-007-2000(SWAB) Rev.04 Date:2020.12

10.3

189mm

159.8

18.9

PEEL HERE

**EZ-FINDER SWABS™**  
FLOCKED SWABS

STERILE EO

Sterile only if peel-pouch is intact.  
This product is a disposable medical device

CE FDA

Product Code: IZ-01  
MFG: YYYY.MM.DD  
EXP: YYYY.MM.DD  
LOT: XXXXXXXX  
Made in China

30

200.00 mm

Peel here

### Miraclean Technology Co.,Ltd.


Room 301, Building A, No.18, Rongshuxia Industrial Zone, Tongxin Community, Baolong Street, Longgang District, Shenzhen, 518116, Guangdong, P.R. China

Tel: 86-755-89616772  
Web: www.mantacc.com


P/N 96000

### Disposable Sampling Swab 1pc

1. Peel apart the package

2. Remove the sterilized swab to collect samples

⚠ Use with caution if allergic to nylon fibre and ABS material.

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


Type:A-04

FDA Registered

- 1.Please read the IFU carefully before use.
- 2.The product is sterile.
- 3.Don't use if package is damaged.

 **JIANGSU HANHENG MEDICAL TECHNOLOGY CO., LTD.**  
16-B4,#1 North Qingyang Road, Tianning District, 213017  
Changzhou, Jiangsu, China

EXECUTOR: **EXTREME MACHINERY MANUFACTURE IMP&EXP CO., LIMITED**  
Flat A, Floor 15, Manly Commercial Building, 16 Soy Street,  
Mong Kok, Kowloon, Hong Kong

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STERILE EO

FDA Registered

- 1.Please read the IFU carefully before use.
- 2.The product is sterile.
- 3.Don't use if package is damaged.

JIANGSU HANHENG MEDICAL TECHNOLOGY CO., LTD  
16-B4,#1 North Qingyang Road, Tianning District, 213017  
Changzhou, Jiangsu, China

EXECUTOR: EXTREME MACHINERY MANUFACTURE IMP&EXP CO., LIMITED  
Flat A, Floor 15, Manly Commercial Building, 16 Soy Street,  
Mong Kok, Kowloon, Hong Kong

XXXXXXXX  
YYYYMMDD  
YYYYMMDD

#### **1.4 Nasal swab labels**

① ⑤

PEEL DOWN

## FA SPONGE SWAB

Specimen Collection Swab  
REF FASPS02

CONTENTS  
SINGLE TIP  
APPLICATORS

- STERILITY ONLY IF PEEL-POUCH IS INTACT
- DISPOSABLE USE.
- DESTROYED AFTER USE
- PLEASE DO NOT USE IT IF EXPIRED OR PACKING DAMAGE

CE 1639 **STERILE EO**

USA FDA Registration  
Number : 3011555966

L = LOT

E = EXP.DATE(YYYY-MM)

MFR. : FA INC.  
10-5, Myeonghaksan-  
danseo-ro, Yeondong-  
myeon, Sejong-si,  
30068 Rep. of KOREA  
[www.facompany.co.kr](http://www.facompany.co.kr)

EC REP EU Representative :  
MT Promedt Consulting  
GmbH  
Aletnhofstrasse 80, 66386  
St, Ingbert, Germany

① ⑤

PEEL DOWN

## FA SPONGE SWAB

Specimen Collection Swab  
REF: FASPS02

CONTENTS  
SINGLE TIP  
APPLICATORS

- STERILITY ONLY IF  
PEEL-POUCH IS INTACT
- DISPOSABLE USE.
- DESTROYED AFTER USE
- PLEASE DO NOT USE IT  
IF EXPIRED OR PACKING  
DAMAGE

CE 1639

STERILE EO

USA FDA Registration  
Number : 3011555966

L =

E =

L.XXXXXXXXXXX

E.YYYY-MM

MFR. : FA INC.  
10-5, Myeonghaksar-  
danseo-ro, Yeondong-  
myeon, Sejong-si,  
30068 Rep. of KOREA  
[www.facompany.co.kr](http://www.facompany.co.kr)

EU REP: EU Representative :  
MT Promedt Consulting  
GmbH  
Aletnhofstrasse 80, 66386  
St, Ingbert, Germany

200.00 mm

Peel here

### Miraclean Technology Co.,Ltd.

No.18, Rongshuxia Industrial Zone, Tongle Community, Longgang District, Shenzhen, China, 518116

Tel: 86-755-89616772  
Web: www.mantacc.com


P/N 93050

### Disposable Sampling Swab 1pc

1. Peel apart the package

2. Remove the sterilized swab to collect samples

⚠ Use with caution if allergic to sponge and POM material.

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CE 0413 FDA Registered  
Medico Technology Co., Ltd

Address: Room 201 of Building 14th and Building 17th, Hengyi Lane, Yuanhu Road, Zhangbei Industrial Park, Longcheng Street, Longgang district, Shenzhen, Guangdong, China  
www.medicoswab.com

Wellkang Ltd(www.CE-marking.eu)  
Enterprise Hub,NW Business Complex,1  
BeraghmoreRd,Derry,BT488SE,N.Ireland.UK

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CE 0413 FDA Registered  
Medico Technology Co., Ltd

Address: Room 201 of Building 14th and Building 17th, Hengyi Lane, Yuanhu Road, Zhangbei Industrial Park, Longcheng Street, Longgang district, Shenzhen, Guangdong, China  
www.medicoswab.com

EC REP

Wellkang Ltd(www.CE-marking.eu)  
Enterprise Hub,NW Business Complex,1  
BeraghmoreRd,Derry,BT488SE,N.Ireland.UK

Product code: MFS-9700KQ

YYYY.MM.DD

YYYY.MM.DD

LOT XXXXXX

190mm

# Disposable sampling swab

EC REP Luxus Lebenswelt GmbH  
Kochstr. 1, 47877, Willich, Germany

Type: G-014

STERILE EO

FDA Registered


1. Please read the IFU carefully before use.
2. The product is sterile.
3. Don't use if package is damaged.

Jiangsu HanHeng Medical Technology Co., Ltd.  
16-B4, #1 North Qingyang Road, Tianning District, Changzhou,  
215017 Jiangsu P.R. China

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STERILE EO

FDA Registered


- 1. Please read the IFU carefully before use.
- 2. The product is sterile.
- 3. Don't use if package is damaged.

Jiangsu HanHeng Medical Technology Co., Ltd.  
16-B4, #1 North Qingyang Road, Tianning District, Changzhou,  
215017 Jiangsu P.R. China

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
20mm

### **1.5 Controls pouch label**

# STANDARD COVID-19 Ag Negative Control swab

195mm x 50


 Printing part  
(black)

# STANDARD COVID-19 Ag Positive Control swab

195mm x 50

 Printing part  
(black)

## **1.6 Extraction tube label**

# Extraction Buffer Tube

 SD BIOSENSOR

2°C  
(36°F)

30°C  
(86°F)

ML27RT1ENR2  
Issue date: 2021.01

# Extraction Buffer Tube

STANDARD Q COVID-19 Ag Test

LOT No. : XXXXXXXXXXXX  
MFG DATE : YYYY.MM.DD.  
EXP DATE : YYYY.MM.DD.  
Quantity : 25pcs

 SD BIOSENSOR

2°C  
(36°F)

30°C  
(86°F)

ML27RT1ENR2  
Issue date: 2021.01

## **1.7 Nozzle label**

**Nozzle cap**  
**80 x 12**

80mm x 12mm

**Nozzle cap**

Code

LOT

Qty 25 pcs

 SD BIOSENSOR

L4PRTZENR2  
Issue date : 2021.01

**Nozzle cap**  
**80 x 12**

80mm x 12mm

**Nozzle cap**

Code

MSH10-2

Qty 25 pcs

Lot

XXXXXX-XXXXX

YYYY.MM.DD

SD BIOSENSOR

L4PRTZENR2  
Issue date : 2021.01

## **2.0 Instructions for use<sup>2</sup>**

---

<sup>2</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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
# STANDARD Q COVID-19 Ag Test

STANDARD™ Q COVID-19 Ag Test

READ AND UNDERSTAND CAREFULLY BEFORE PERFORMING THE TEST

SD BIOSENSOR

## KIT CONTENTS

STANDARD COVID-19 Ag Positive Control swab  
STANDARD COVID-19 Ag Negative Control swab  
Instructions for use


## SPECIMEN COLLECTION AND PREPARATION

### Specimen preparation (Nasal swab)

1. Tilt patient's head back slightly.
2. While rotating the swab, insert swab less than one inch (about 2cm) into nostril until resistance is met at turbinate.
3. Rotate the swab 4 times against nasal wall.
4. Repeat in other nostril using the same swab.

Only one swab is used per test. Rub both nostrils with the same swab. Using proper pressure, rub the swab against all of the inside nasal walls of your nostril. Make at least 4 big circles along the nasal wall. Do not just spin the swab.

1. Apply 4 drops of extraction buffer to the specimen well.
2. Read the test result in 15-30 minutes.


Place the test device on a flat surface. Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles. Do not read test results after 30 minutes. It may give false results.

1. Insert the swab into an extraction buffer tube. While squeezing the buffer, rotate the swab more than 5 times.
2. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
3. Press the nozzle cap tightly onto the tube.
4. Specimen should be tested as soon as possible after collection.
5. Specimens may be stored at room temperature (15 – 25°C or 2-8°C/ 36-64°F for up to 4 hours prior to testing.

Without the tube squeezing process, improper results may occur due to the large amount of buffer absorbed by the swab. If the specimen storage condition is out of instructions as below, do not use. The Nasal swab is stored in extraction buffer for more than 4 hours at 5±3°C or 20±5°C. Freezing and thawing of Nasal swab is more than 1 cycle.

## PREPARATION AND TEST PROCEDURE

### Preparation

1. Carefully read instructions for using STANDARD Q COVID-19 Ag Test.
2. Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
3. Check the test device and the desiccant pack in the foil pouch.

### Test Procedure

1. Apply 4 drops of extracted specimen to the specimen well of the test device.
2. Read the test result in 15-30 minutes.

Place the test device on a flat surface. Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles. Do not read test results after 30 minutes. It may give false results.

## PERFORMANCE CHARACTERISTICS

### Clinical evaluation

Clinical performance of the STANDARD Q COVID-19 Ag Test was evaluated using nasal swab samples from 696 subjects in a prospective study at a clinical center in Germany. The study cohort included adults at high risk for SARS-CoV-2 infection according to clinical suspicion, 311 subjects underwent nasal sampling performed by healthcare professionals and 385 subjects followed instructions to obtain a nasal swab sample by themselves. Self-collection was performed under the supervision of healthcare workers without interference or assistance. Test procedures and result reading were always performed by healthcare professionals. RT-PCR tests (Roche cobas® SARS-CoV-2 and TibMolbiol SARS-CoV-2 E-gene assay) using combined nasopharyngeal/oropharyngeal swab samples were used as the comparator methods. Nasal sampling always preceded the combined NP/OP sampling.

### Test sensitivity & specificity

The following tables summarize the patient and performance characteristics of the STANDARD Q COVID-19 Ag Test. The relative sensitivity was 89.6% (CI: 70.3% to 95.9% CI: 79.7% - 95.7%) for professionally collected samples, and 89.1% (CI value ≤ 30; 95% CI: 78.8% - 95.5%) for self-collected samples. For patients for whom data post symptom onset was known, and was 0-5 days, the relative sensitivity in comparison to RT-PCR was 86.7% (95% CI: 75.4% - 94.1%) for professionally collected nasal samples and 88.9% (95% CI: 77.4% - 95.8%) for self-collected nasal samples. The relative specificity in comparison to RT-PCR was 99.1% (95% CI: 96.9% - 99.9%) for professionally collected nasal samples and 99.0% (95% CI: 97.2% - 99.8%) for self-collected nasal samples. In total, nasal swab samples from 150 PCR-positive and 546 PCR-negative individuals were evaluated using the STANDARD Q COVID-19 Ag Test. The relative sensitivity and relative specificity were 82.7% (95% CI: 75.6% - 88.4%) and 99.1% (95% CI: 97.3% - 99.7%), respectively.

## INTERPRETATION OF TEST RESULTS

Test result	Example	Description
Negative		1. A purple colored band will appear in the upper section of the result window. This band is control line (C). 2. A purple colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T). 3. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.
Positive		1. A purple colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T). 2. A purple colored band will appear in the upper section of the result window. This band is control line (C). 3. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.
Invalid		1. The presence of any line no matter how faint the result is considered positive. Positive results should be considered in conjunction with the clinical history and other data available.

The presence of any line no matter how faint the result is considered positive. Positive results should be considered in conjunction with the clinical history and other data available.

## CONTROL PREPARATION AND TEST PROCEDURE

### Positive/Negative control

#### Preparation

1. Put the positive or negative control swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
2. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
3. Press the nozzle cap tightly onto the tube.

#### Test Procedure

1. Apply 4 drops of extraction buffer to the specimen well.
2. Read the test result in 15-30 minutes.

Place the test device on a flat surface. Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles. Do not read test results after 30 minutes. It may give false results.

1. Put the positive or negative control swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
2. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
3. Press the nozzle cap tightly onto the tube.

#### Test Procedure

1. Apply 4 drops of extraction buffer to the specimen well.
2. Read the test result in 15-30 minutes.

Place the test device on a flat surface. Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles. Do not read test results after 30 minutes. It may give false results.

## INTERPRETATION OF CONTROL TEST RESULT

STANDARD COVID-19 Ag Control Positive swab: Positive	Result	Interpretation	Follow up
Test (T) Line Positive	PASS	-	-
Test (T) Line Negative	FAIL	Retest	-
No Control (C) Line	Invalid	-	Retest

STANDARD COVID-19 Ag Negative Control swab: Negative	Result	Interpretation	Follow up
Test (T) Line Negative	PASS	-	-
Test (T) Line Positive	FAIL	Retest	-
No Control (C) Line	Invalid	-	Retest

## EXPLANATION AND SUMMARY

### Introduction

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or SARS-CoV-2 (COVID-19), was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as Severe Acute Respiratory Syndrome (SARS). These kits are helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

### Intended use

STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific nucleocapsid protein antigen from SARS-CoV-2 present in human nasal or nasopharyngeal specimens. This product is intended for healthcare professionals at the clinical setup and point of care sites, as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms of SARS-CoV-2 infection. It provides only an initial screening test result. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

### Test principle

STANDARD Q COVID-19 Ag test device has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-SARS-CoV-2 antibody conjugated with color particles are used as detectors for SARS-CoV-2 antigenic during. During the test, SARS-CoV-2 antigen in the specimen interact with the monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making antigen-antibody color particle complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the mouse monoclonal anti-SARS-CoV-2 antibody. A colored test line would be visible in the result window if SARS-CoV-2 antigens are present in the specimen. The intensity of colored test line will vary

depending upon the amount of SARS-CoV-2 antigen present in the specimen. If SARS-CoV-2 antigens are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

### Kit contents (09COV33D)

- 1. Test device (individually in a foil pouch with desiccant) x 25
- 2. Extraction buffer tube x 25
- 3. Nozzle cap x 25
- 4. Sterile swab x 25
- 5. STANDARD COVID-19 Ag Positive Control swab x 1
- 6. STANDARD COVID-19 Ag Negative Control swab x 1
- 7. Instructions for use x 1

## MATERIALS REQUIRED BUT NOT PROVIDED

1. Personal protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
2. Timer
3. Biohazard container

## KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

## WARNINGS AND PRECAUTIONS

1. Bring the kit contents and the specimens to room temperature before testing.
2. Do not re-use the test kit.
3. Do not use the test kit if the pouch is damaged or the seal is broken.
4. Do not use the extraction buffer tube of another lot.
5. Do not smoke, drink or eat while handling specimen.
6. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
7. Clean up spills quickly using an appropriate disinfectant.
8. Handle all specimens as if they contain infectious agents.
9. Observe established precautions against microbiological hazards throughout testing procedures.
10. Dispose of all specimens and materials used to perform the test as biohazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
12. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:  
Warning: H317 May cause an allergic skin reaction.  
H412 Harmful to aquatic life with long lasting effects.  
H319 Causes serious eye irritation.  
Prevention: P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P273 Avoid release to the environment.  
P280 Wear eye protection/face protection.  
Response: P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.  
P337 + P313 If eye irritation persists: Get medical advice/attention.  
P362 + P364 Take off contaminated clothing and wash it before reuse.  
For customers in the European Economic Area: Contains SVHC: octyl/nonylphenol ethoxylates.  
For use as part of an IVD method and under controlled conditions only – acc. to Art. 56.3 and 3.23 REACH Regulation.

## LIMITATION OF TESTS

1. The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
2. The test should be used for the detection of SARS-CoV-2 antigen in human nasal swab specimens only, other specimen types have not been validated.
3. This test can not be used for quantifying SARS-CoV-2 antigen concentration. Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
4. The test result must always be evaluated with other data available to the physician.
5. A negative result may occur if the concentration of antigen in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or molecular assay.
7. Positive test results do not rule out co-infections with other pathogens.

## ES REF Q-NCOV-016 Cat No. 09COV33D


# STANDARD Q COVID-19 Ag Test

STANDARD™ Q COVID-19 Ag Test

READ AND UNDERSTAND CAREFULLY BEFORE PERFORMING THE TEST

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
**COVID-19 Ag Control swab**

STANDARD™ COVID-19 Ag Control swab

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU  
PERFORM THE TEST**EXPLANATION AND SUMMARY****■ Intended use**


STANDARD COVID-19 Ag Control swab is intended for use as an external quality control material to monitor the performance of STANDARD COVID-19 Ag products. It is important to perform control tests with positive and negative control materials to assure STANDARD COVID-19 Ag testing system is working properly.

**■ Principle**

This product was designed for use with STANDARD COVID-19 Ag products for purpose of monitoring assay performance and maintaining quality assurance.

**KIT CONTENTS**

STANDARD COVID-19 Ag Control swab is contained 10 Positive control swabs and 10 Negative control swabs.

STANDARD COVID-19  
Ag Positive Control  
swabSTANDARD COVID-19  
Ag Negative Control  
swabInstructions for  
use**MATERIALS REQUIRED BUT NOT PROVIDED**

- STANDARD COVID-19 Ag Products


**PREPARATION AND CONTROL TEST PROCEDURE**

This product should be treated the same as patient specimens and run in accordance with instruction accompanying the instrument and kit being used. It is recommended that positive and negative controls be run:


- Once for each new lot,
- Once for each untrained operator,
- Once for each new shipments of test kits,
- As required by test procedures in this instructions and in accordance with local, state and federal regulations of accreditation requirements.

**■ Preparation**

1. Check the expiry date at the label of a STANDARD COVID-19 Ag Control swab package and Front of the foil pouch of a STANDARD COVID-19 Ag Control swab. Use another lot, if it has expired.


2. Insert the positive or negative control swab into an extraction buffer tube which is in the STANDARD COVID-19 Ag products. Stir the swab at least five times.

3. Remove the swab while squeezing the sides of the tube to extract the liquid from swab. Dispose of the used swab in accordance with your biohazard waste disposal protocol.

4. Press the nozzle cap tightly onto the tube.
5. Testing and interpreting in accordance with the instruction accompanying the STANDARD COVID-19 Ag products.

**■ Test procedure**

1. Test should be run following the instructions for use provided by the test kit for unknown specimens.
2. Results should be determined in the same manner as used for unknown specimens when tested using the test kits.
3. Results may vary among methodologies, among products, and among different lots of the same kit.

**INTERPRETATION OF TEST RESULTS**

Test should be interpreted following the instructions for use provided by the test kit for unknown specimens.

**STORAGE AND STABILITY**

Store the STANDARD COVID-19 Ag Control swab at 2-30°C/36-86°F. Kit materials are stable until expiration date printed on the outer box.

**WARNINGS AND PRECAUTIONS**

1. If there is evidence of microbial contamination in the reconstituted control, discard the control.
2. Wear protective clothing and gloves when handling specimens or reagents.
3. Clean any spillage by immediately and thoroughly wiping up with a suitable disinfectant such as 1% sodium hypochlorite solution.
4. Handle and dispose of all specimens, controls and materials used in testing as though they contain infectious agents.
5. Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local distributor of SD Biosensor, Inc.

## LIMITATIONS

1. This product is provided for quality assurance purposes and must not be used for calibration or as primary reference preparations in any test procedure.
2. Adverse storage conditions or use of outdated reagents may produce erroneous results.
3. This product should not be used past the expiration date.
4. Alterations in physical appearance may indicate instability or deterioration of this product. If there is evidence of microbial contamination in this product, discard it.

**Manufactured by SD Biosensor, Inc.**

**Head office :** C-4th&5th, 16, Deogyong-daero 1556beongil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA

**Manufacturing site :** 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

**Authorized Representative**

**MT Promedt Consulting GmbH**

Altenhofstrasse 80 66386 St. Ingbert Germany

Phone : +49 6894 581020, Fax : +49 6894 581021

Please contact us for any complaints/inquiries/suggestions via email ([sales@sdbiosensor.com](mailto:sales@sdbiosensor.com)), website ([www.sdbiosensor.com](http://www.sdbiosensor.com)).

L24COVC11ENR2  
Issue date: 2022.05

Reference number

In vitro Diagnostics

Consult Instructions for Use

Contains Sufficient for <math>n</math>-Tests

To indicate the temperature limitations in which the transport package has to be kept and handled.

Do not re-use

Use by

Batch code

Manufacturer

Date of manufacture

Indicate that you should keep the product dry

Keep away from sunlight

Fulfill the requirements of Directive 98/79/EC on in vitro diagnostic medical devices

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© SD BIOSENSORS