## WHO Emergency Use Assessment SARS-CoV-2 IVDs PUBLIC REPORT

**Product: Panbio COVID-19 Ag Rapid Test Device (NASAL)**

**Manufacturer: Abbott Rapid Diagnostics Jena GmbH**

**EUL Number: EUL 0587-032-00**

**Outcome: Accepted**

The EUL process is intended to expedite the availability of in vitro diagnostics needed in public health emergency situations and to assist interested UN procurement agencies and the Member States in determining the acceptability of using specific products in the context of a Public Health Emergency of International Concern (PHEIC), based on an essential set of available quality, safety, and performance data. The EUL procedure includes the following:

- Quality Management Systems Review and Plan for Post-Market Surveillance: a desktop review of the manufacturer’s Quality Management System documentation and specific manufacturing documents.
- Product Dossier Review: assessment of the documentary evidence of safety and performance.

Panbio COVID-19 Ag Rapid Test Device (NASAL), product codes 41FK11 and 41FK21, CE marked regulatory version, manufactured by Abbott Rapid Diagnostics Jena GmbH, Orlaweg 1, 07743 Jena, Germany, was listed on 19 November 2020. The product was assessed for EUL renewal and met EUL renewal requirements on 6 June 2022. Therefore, the product is eligible for WHO procurement until further notice.

### Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the accepted product under EUL, for which WHO has been notified and has undertaken a review. Amendments to the report are summarized, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Addition of a new kit with reference number 41FK21 that contains test devices with a 2D barcode printed on the test device, which encodes traceability information for the product, and changes to the labelling.	20 January 2021
3.0	1. Validation and implementation of additional nasal swabs from a different supplier to Panbio COVID-19 Ag Rapid Test Device (Nasal).	4 May 2021

	2. Amendments to the IFUs to include data from a study on asymptomatic individuals. 3. Amendments to the nasal product IFU to add the results of a study of supervised self-collection of nasal swab specimens <sup>1</sup> .	
4.0	Extension of shelf life to 24 months based on accelerated stability data. Real-time studies are ongoing.	9 November 2021
5.0	Updating of IFU to version A6.	13 January 2022
6.0	Extension of product's eligibility for procurement after EUL renewal assessment until further notice.	13 June 2022
7.0	Change in kit box size and kit components, including pre-filled buffer tubes, new Panbio COVID-19 Ag Rapid Test Device Positive Control Kit.	19 October 2022

**Intended use:**

According to the claim of intended use from Abbott Rapid Diagnostics Jena GmbH, "*Panbio COVID-19 Ag Rapid Test Device is an in vitro diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag) in human nasal swab specimens from individuals who meet COVID-19 clinical and / or epidemiological criteria. Panbio COVID-19 Ag Rapid Test Device is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation. The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 infection and they cannot be used as the sole basis for treatment or other management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information. The test is not intended to be used as a donor screening test for SARS-CoV-2.*"

**Validated specimen type:**

Nasal swab specimens.

---

<sup>1</sup> The WHO EUL assessment was abridged based on reliance on the product authorization under the USFDA Emergency Use Authorization. WHO does not recommend the use of health care provider instructed/observed self-collected nasal, nasopharyngeal or oropharyngeal swab specimens. Please refer to the WHO guidance on Diagnostic testing for SARS-CoV-2 at the following link.

<https://apps.who.int/iris/discover?query=Diagnostic+testing+for+SARS-CoV-2%3A+Interim+guidance>

**Test kit contents:**

Component	25 tests (product code 41FK11)	25 tests (product code 41FK21)
Test device (individually in a foil pouch with desiccant)	25	25 with a 2D barcode printed on the test device.
Buffer tubes	25	25
Sterilized nasal swabs for sample collection	25	25
Tube rack	1	1
Instructions for use with quick reference guide	1	1

**Items required but not provided**

- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
- Timer
- Biohazard container
- Panbio COVID-19 Ag Rapid Test Device Positive Control Kit (41FK11CON)

**Storage:**

2-30°C.

**Shelf-life upon manufacture:**

24 months (real-time stability studies are ongoing).

**Warnings/limitations:**

Refer to the instructions for use (IFU).

**Product dossier assessment**

Abbott Rapid Diagnostics Jena GmbH submitted a change notification for their listed product, Panbio COVID-19 Ag Rapid Test Device (Nasopharyngeal), to introduce nasal swab specimens as an additional claimed specimen type. The new configuration was adapted from the corresponding accepted EUL product (Panbio COVID-19 Ag Rapid Test Device (Nasopharyngeal)), for which a WHO EUL assessment had already been conducted and accepted. Additional data was generated to meet particular requirements for EUL as set out in the “*Instructions for Submission Requirements: In vitro diagnostics (IVDs) Detecting SARS-CoV-2 Nucleic Acid and rapid diagnostics tests detecting SARS-CoV-2 nucleic acid or antigens (PQDx\_0347)*”. The information (data and documentation) submitted in the change notification was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

### Post listing Commitments for EUL:

As commitments to listing, Abbott Rapid Diagnostics Jena GmbH committed to,

1. Provide additional microbial interference study report by 31 December 2020. The commitment was fulfilled. Issue closed.
2. Provide accelerated stability data by 31 January 2021, 31 March 2021, and 31 July 2021. The commitment was fulfilled. Issue closed.
3. Provide interim and final real-time stability reports by 31 July 2022 and 31 October 2022, respectively.
4. Provide additional clinical study information by 31 March 2021. The commitment was fulfilled. Issue closed.

The risk-benefit assessment is acceptable.

### Quality Management Systems Review

To establish eligibility for WHO procurement, Abbott Rapid Diagnostics Jena GmbH was asked to provide up-to-date information about the status of its quality management system.

Based on the review of the submitted quality management system documentation by WHO staff, it was established that Abbott Rapid Diagnostics Jena GmbH provided sufficient information to fulfil the requirements described in the *“Instructions for Submission Requirements: In vitro diagnostics (IVDs) Detecting SARS-CoV-2 Nucleic Acid and rapid diagnostics tests detecting SARS-CoV-2 nucleic acid or antigens (PQDx\_347)”*.

Quality management documentation assessment is acceptable.

### Plan for Post-Market Surveillance

Post-market surveillance, including monitoring all customer feedback, detecting and acting on adverse events, product problems, non-conforming goods and processes is critical component of minimizing potential harm of an IVD listed for emergency use.

The following post-EUL activities are required to maintain the EUL listing status:

1. Notification to WHO of any planned changes to a EUL product, in accordance with *“WHO procedure for changes to a WHO prequalified in vitro diagnostic”* (document number PQDx\_121); and
2. Post-market surveillance activities, in accordance with *“Guidance for post-market surveillance and market surveillance of medical devices, including in vitro diagnostics”* (ISBN 978-92-4-001531-9).

Abbott Rapid Diagnostics Jena GmbH is also required to submit all categories of complaints in a summarized form. Certain categories of complaints and changes to the product must be notified immediately to WHO, as per the above-mentioned documents.

The manufacturer has committed to ensuring that post-emergency use listing safety, quality and performance monitoring activities are in accordance with WHO guidance "*Guidance for post-market surveillance and market surveillance of medical devices, including in vitro diagnostics*" (ISBN 978-92-4-001531-9).<sup>2</sup>

### **Scope and duration of procurement eligibility**

Panbio COVID-19 Ag Rapid Test Device (NASAL), product codes 41FK11 and 41FK21, manufactured by Abbott Rapid Diagnostics Jena GmbH, is considered eligible for WHO procurement until further notice. The assay may detect the Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) antigen. This listing does not infer that the product meets WHO prequalification requirements and does not mean that the product is listed as WHO-prequalified.

As part of the ongoing requirements for listing as eligible for WHO procurement, Abbott Rapid Diagnostics Jena GmbH must engage in post-market surveillance activities to ensure that the product continues to meet safety, quality, and performance requirements. Abbott Rapid Diagnostics Jena GmbH is required to notify WHO of any complaints, including adverse events related to the use of the product, within 7 days and any changes to the product.

WHO reserves the right to rescind eligibility for WHO procurement if additional information on the safety, quality, and performance during post-market surveillance activities and if new data becomes available to WHO that changes the risk-benefit balance.

---

<sup>2</sup> Available on the web page

<https://www.who.int/publications/i/item/guidance-for-post-market-surveillance-and-market-surveillance-of-medical-devices-including-in-vitro-diagnostics>.

## **Labelling**


### **1.0 Labels**

### **2.0 Instructions for Use (IFU)**

## **1.0 Product labels**

## **1.1 Outside box label**

### **1.1.1 Product code 41FK11**

REF 41FK11 IVD CE 20°C 30°C 36°F 86°F

# Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

## Abbott Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

**EN**  
In vitro diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag)

- 25 Test devices
- 25 Buffer tubes
- 25 Sterilized nasal swabs for sample collection
- 1 Tube rack
- 1 Instructions for Use with Quick Reference Guide (Nasal)

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Made in Korea: Abbott Diagnostics Korea Inc.  
65, Borahgal-ro, Gheung-gu, Yongin-si,  
Gyeonggi-do, 17099, Republic of Korea

## Abbott Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

Abbott Rapid Diagnostics Jena GmbH  
Orlweg 1, D-07743 Jena, Germany  
www.globalpointofcare.abbott


LOT XXXXXXXX Sub : X  
buffer : XXXXXXXX  
Mfg : YYYY.MM.DD YYYY.MM.DD  
(YY)MMDD (YY)MMDD  
(10)XXXXXXXX (240)XXXXXXXX

## Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

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
PRODUCT - Panbio™ COVID-19 Ag  
Rapid Test Device (Nasal), 25T  
SIZE - 230 X 125 X 80MM

PMS 303 C  
PMS 2269 C

PN - 41FK11P-70-A0  
REV. DATE - 2022. 03. 10  
CREATOR - THINKBOX

**1.1.2 Product code 41FK21**

REF 41FK21 IVD CE 20°C 36°F 25

# Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

## Abbott Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

**EN**  
In vitro diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag)

- 25 Test devices
- 25 Buffer tubes
- 25 Sterilized nasal swabs for sample collection
- 1 Tube rack
- 1 Instructions for Use with Quick Reference Guide (Nasal)

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Made in Korea: Abbott Diagnostics Korea Inc. 65, Borahagal-ro, Gheung-gu, Yongin-si, Gyeonggi-do, 17099, Republic of Korea

## Abbott Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)


Abbott Rapid Diagnostics Jena GmbH  
Orlweg 1, D-07743 Jena, Germany  
www.globalpointofcare.abbott

LOT XXXXXXXX Sub : X  
buffer : XXXXXXXX  
Mfg : YYYY.MM.DD (YYYY.MM.DD)

(01) 04051892006751  
(17) YYMMDD  
(10) XXXXXXXX  
(240) XXXXXXXX

## Abbott Panbio COVID-19 Ag RAPID TEST DEVICE (NASAL)

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© 2022 Abbott. All rights reserved.  
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

PRODUCT - Panbio™ COVID-19 Ag  
Rapid Test Device (Nasal), 25T


SIZE - 230 X 125 X 80MM

PMS 303 C  
PMS 2269 C

PN - 41FK21P-70-A0  
REV. DATE - 2022. 03. 22  
CREATOR - THINKBOX

### **1.1.3 COVID19 Ag Rapid Test Device Positive Control Kit (Product code 41FK11CON)**

2°C 30°C 36°F 86°F 10 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79 80 81 82 83 84 85 86 87 88 89 90 91 92 93 94 95 96 97 98 99 100 REF 41FK11CON IVD

# COVID-19 AG RAPID TEST DEVICE POSITIVE CONTROL KIT

**EN**  
Panbio™ COVID-19 Ag Rapid Test Device Positive Control Kit contains quality control swabs for use with Panbio™ COVID-19 Ag Rapid Test Device only.

- 10 Positive Control Swabs
- 1 Instructions for Use with Quick Reference Guide

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
Mfg : YYYY.MM.DD    YYY.MM.DD (YYYY.MM.DD)

(01) 04051892006782  
(17) YYMMDD  
(10) XXXXXXXX  
(240) AF41FK11CON

## COVID-19 AG RAPID TEST DEVICE POSITIVE CONTROL KIT

2°C 30°C 36°F 86°F 10 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79 80 81 82 83 84 85 86 87 88 89 90 91 92 93 94 95 96 97 98 99 100 REF 41FK11CON IVD

**PT**  
O Panbio™ COVID-19 Ag Rapid Test Device Positive Control Kit contém swabs de controle de qualidade para uso apenas com o Panbio™ COVID-19 Ag Rapid Test Device.

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© 2022 Abbott. All rights reserved.  
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Made in Korea: Abbott Diagnostics Korea Inc.  
65, Borahagal-ro, Giheung-gu, Yongin-si,  
Gyeonggi-do, 17099, Republic of Korea



PRODUCT - Panbio™ COVID-19 Ag Rapid Test Device Positive Control Kit, 10T

PMS 303 C  
PMS 2269 C



PN. - 41FK11CON-70-A0  
REV. DATE - 2022. 03. 22  
CREATOR - THINKBOX

SIZE - 190 X 68 X 40 MM

## 1.2 Test device labels (front and back side) Big Pouch


### 1.3 Test device labels (front and back side) Small Pouch


### 1.4. Testing device printing


#### 1.4.1 Product code 41FK11




#### 1.4.2 Product code 41FK21


### 1.5 Extraction tubes (bag label)

1.6 Sterilized swabs

Buffer tube label printed on the tube sealing

### 1.7 Positive Control swabs

## **2.0 Instructions for use<sup>3</sup>**

---

<sup>3</sup> English version of the IFU was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

**REF** 41FK11/41FK21

Panbio™  
**COVID-19 Ag Rapid  
Test Device**  
(NASAL)

**Abbott Rapid Diagnostics Jena GmbH**  
Orlaweg 1, D-07743 Jena, Germany  
[www.globalpointofcare.abbott](http://www.globalpointofcare.abbott)

Date Issued : 2022.03  
41FK11/41FK21-CP-07-A0

© 2022 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

## About the Test

### Introduction

The novel coronaviruses belong to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

### Test Principle

Panbio™ COVID-19 Ag Rapid Test Device contains a membrane strip, which is pre-coated with immobilized anti-SARS-CoV-2 antibody on the test line and mouse monoclonal anti-chicken IgY on the control line. Two types of conjugates (human IgG specific to SARS-CoV-2 Ag gold conjugate (binds to the nucleocapsid protein) and chicken IgY gold conjugate) move upward on the membrane chromatographically and react with anti-SARS-CoV-2 antibody and pre-coated mouse monoclonal anti-chicken IgY respectively. For a positive result, human IgG specific to SARS-CoV-2 Ag gold conjugate and anti-SARS-CoV-2 antibody will form a test line in the result window. Neither the test line nor the control line are visible in the result window prior to applying the patient specimen. A visible control line is required to indicate a test result is valid.

### Intended Use

Panbio™ COVID-19 Ag Rapid Test Device is an *in vitro* diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag) in human nasal swab specimens from individuals who meet COVID-19 clinical and / or epidemiological criteria. Panbio™ COVID-19 Ag Rapid Test Device is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation. The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 infection and they cannot be used as the sole basis for treatment or other management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information. The test is not intended to be used as a donor screening test for SARS-CoV-2.

### Kit Variants

- **41FK11** No 2D barcode printed on the contained test devices
- **41FK21** Contains test devices with a 2D barcode printed on the test device, which encodes traceability information for the product

## Materials Provided

- 25 Test devices with desiccant in individual foil pouch
- 25 Buffer tubes
- 25 Sterilized nasal swabs for sample collection
- 1 Tube rack
- 1 Instructions for Use with Quick Reference Guide (Nasal)

## Materials Required but not Provided

- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves), Timer, Biohazard container

## Active Ingredients of Main Components

- **1 Test device** Gold conjugate: Human IgG specific to SARS-CoV-2 Ag gold colloid and Chicken IgY - gold colloid, Test line: Mouse monoclonal anti-SARS-CoV-2, Control line: Mouse monoclonal anti-Chicken IgY
- **Buffer** Tricine, Sodium Chloride, Tween 20, Sodium Azide (<0.1%), Proclin 300

## Storage and Stability

1. The test kit should be stored at a temperature between 2-30 °C. Do not freeze the kit or its components.

**Note:** When stored in a refrigerator, all kit components must be brought to room temperature (15-30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.

2. Perform the test immediately after removing the test device from the foil pouch.
3. Do not use the test kit beyond its expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.
6. Direct swab specimens should be tested immediately after collection.

## Warnings

1. For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
2. These instructions must be strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instruction prior to performing a test.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and buffer.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially

contaminated materials (i.e. swab, buffer tube, test device) in a biohazard container as if they were infectious waste and dispose according to applicable local regulations.

8. Do not mix or interchange different specimens.
9. Do not mix reagent of different lots or those for other products.
10. Do not store the test kit in direct sunlight.
11. To avoid contamination, do not touch the head of provided swab when opening the swab pouch.
12. The sterilized swabs should be used only for nasal specimen collection.
13. To avoid cross-contamination, do not reuse the sterilized swabs for specimen collection.
14. Do not dilute the collected swab with any solution except for the provided buffer.
15. The buffer contains <0.1% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with a large volume of water.

## **Test Procedure (Refer to Figure)**

### **Nasal swab Specimens**

**Note:** Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.

### **Test Preparation**

1. Allow all kit components to reach a temperature between 15-30 °C prior to testing for 30 minutes.
2. Remove the test device from the foil pouch prior to use. Place on a flat, horizontal and clean surface.
3. Keep buffer tube upright and remove foil. Place the buffer tube in the tube rack.

### **Nasal Mid-Turbinate (NMT) Specimen Collection & Extraction**

Specimens are collected by the professional user as described below. Alternatively, nasal specimen collection steps 1-3 can be completed by the patient according to oral instructions and under supervision of the professional user. For supervised patient self-collection, the swab is handed to the patient by the professional user and after sampling, the patient hands the swab back to the professional user to complete the remaining steps of the procedure.

1. Tilt the patient's head back 70 degrees. While gently rotating the swab, insert swab less than one inch (about 2 cm) into nostril (until resistance is met at the turbinates).  
⚠ **Caution:** Ensure that the patient's head is kept still during nasal specimen collection, as sudden movements may cause swab stick breakage.
2. Rotate the swab five times against the nasal wall then slowly remove from the nostril.
3. Using the same swab repeat the collection procedure with the second nostril.  
**Note:** Ensure a minimum waiting period of 24 hours before performing a new nasal sampling from both nostrils (e.g. for a repeat test).  
⚠ **Caution:** If the nasal swab stick breaks prior to obtaining a nasal specimen, repeat specimen collection with a new swab. If the nasal swab breaks during sampling, consultation with a medical healthcare professional is recommended to determine and initiate necessary treatment and monitoring.
4. Swirl the swab tip in the buffer fluid inside the buffer tube, pushing into the wall of the buffer tube at least five times and then squeeze out the swab by squeezing the buffer tube with your fingers.
5. Break the swab at the breakpoint and close the buffer tube with the attached nozzle.

### Reaction with Test Device

1. Dispense 5 drops of extracted specimens vertically into the specimen well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.  
⚠ **Caution:** Bubbles that occur in the Buffer tube can lead to inaccurate results. If you are unable to create sufficient drops, this may be caused by clogging in the dispensing nozzle. Shake the tube gently to release the blockage until you observe free drop formation.
2. Dispose of the buffer tube containing the used swab according to your local regulations and biohazard waste disposal protocol.
3. Start timer. Read result at 15 minutes. Do not read results after 20 minutes.
4. Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

### Test Interpretation (Refer to Figure)

1. **Negative result:** The presence of only the control line (C) and no test line (T) within the result window indicates a negative result.
2. **Positive result:** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a positive result.  
⚠ **Caution:** The presence of any test line (T), no matter how faint, indicates a positive result.
3. **Invalid result:** If the control line (C) is not visible within the result window after performing the test, the result is considered invalid.

## Test Limitations

1. The contents of this kit are to be used for the professional and qualitative detection of SARS-CoV-2 antigen from nasal swab. Other specimen types may lead to incorrect results and must not be used.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. A negative test result may occur if the specimen was collected, extracted or transported improperly. A negative test result does not eliminate the possibility of SARS-CoV-2 infection and should be confirmed by viral culture or a molecular assay.
4. Positive test results do not rule out co-infections with other pathogens.
5. Test results must be evaluated in conjunction with other clinical data available to the physician.
6. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
7. Panbio™ COVID-19 Ag Rapid Test Device is not intended to detect from defective (non-infectious) virus during the later stages of viral shedding that might be detected by PCR molecular tests.
8. Positive results may occur in cases of infection with SARS-CoV.

## Quality Control

### 1. Internal Quality Control:

The test device has a test line (T) and a control line (C) on the surface of the test device. Neither the test line nor the control line are visible in the result window before applying a specimen. The control line is used for procedural control and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

### 2. External Quality Control:

Good laboratory practices recommends the use of external control materials. Users should follow federal, state and local guidelines concerning the frequency of external control testing. Dedicated external control material can be ordered by the following catalog number: 41FK11CON

## Performance Characteristics

### 1. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Symptomatic)

Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing 104 positive nasal swab specimens and 404 negative specimens for SARS-CoV-2 antigen (Ag) to have a sensitivity of 98.1% (95% CI: 93.2-99.8%) and a specificity of 99.8% (95% CI: 98.6-100.0%). Clinical specimens were determined to be positive or negative using an FDA EUA RT-PCR reference method. The individuals on which the reported sensitivity and specificity are based also had a nasopharyngeal swab taken, which was tested in the FDA EUA approved RT-PCR.

## Panbio™ COVID-19 Ag Rapid Test Device Results

		Nasal PCR Test Result		
		Positive	Negative	Total
Panbio™ COVID-19 Ag Rapid Test Device Result (nasal swab specimens)	Positive	102	1	103
	Negative	2	403	405
	Total	104	404	508
		<b>Sensitivity</b>	<b>Specificity</b>	<b>Overall Percent Agreement</b>
		<b>98.1%</b> [93.2%; 99.8%]	<b>99.8%</b> [98.6%; 100.0%]	<b>99.4%</b> [98.3%; 99.9%]

- Performance data was calculated from a study of individuals suspected of exposure to COVID-19 or who have presented with symptoms in the last 7 days.
- Stratification of the positive specimens post onset of symptoms or suspected exposure between 0-3 days has a sensitivity of 100.0% (95% CI: 92.3-100.0%; n=46) and 4-7 days has a sensitivity of 96.6% (95% CI: 88.1-99.6%; n=58).
- Positive agreement of the Panbio™ COVID-19 Ag Rapid Test Device is higher with samples of Ct values  $\leq 30$  with a sensitivity of 100.0% (95% CI: 96.0-100.0%) and Ct values  $\leq 33$  with a sensitivity of 99.0% (95% CI: 94.5-100.0%). Patients with Ct value  $> 30$  are no longer contagious.
- The clinical performance data was also calculated vs nasopharyngeal swab specimens using an FDA EUA RT-PCR reference and has a sensitivity of 91.1% (95% CI: 84.2-95.6%) and specificity of 99.7% (95% CI: 98.6-100.0%).

### 2. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Asymptomatic)

Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing 483 asymptomatic subjects for SARS-CoV-2 antigen (Ag). Clinical specimens were determined to be positive or negative using an FDA EUA RT-PCR reference method.

The positive results (n=50) were stratified by the comparator method cycle threshold (Ct) counts and assessed to better understand the correlation of product performance, as a surrogate for the amount of virus present in the clinical sample. A lower Ct value corresponds to a higher virus concentration. As presented in the table below, the positive agreement increases with lower Ct values.

The specificity (n=433) was 100% with 95% CI [99.2%; 100.0%].

The results for sensitivity are summarized in the following table:

	All Nasal PCR Positive Samples (n=50)	Ct values ≤ 33 (n=40)	Ct values ≤ 30 (n=32)
Sensitivity [CI 95%]	66.0% [51.2%; 78.8%]	80.0% [64.4%; 90.9%]	93.8% [79.2%; 99.2%]

Patients with Ct value >30 are no longer contagious.

### 3. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Self-Collected Swab)

The clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was assessed in 287 symptomatic subjects (≥16 years of age) who collected their swab specimen (self swabbing) under the direction and supervision of a trained professional. The swab was then handed to the trained professional who executed the remaining steps of the procedure. The trained professional also collected a nasopharyngeal swab from each subject to be used as a reference specimen. The reference specimen was tested on the Panbio™ COVID-19 Ag Rapid Test Device.

The results are summarized in the following table:

		Panbio™ COVID-19 Ag Rapid Test Device (Nasopharyngeal)		
		Positive	Negative	Total
Panbio™ COVID-19 Ag Rapid Test Device (Nasal) – Self-collected Swab	Positive	110	0	110
	Negative	2	175	177
	Total	112	175	287
		Positive Agreement	Negative Agreement	Overall Percent Agreement
		98.2% [93.7%; 99.8%]	100.0% [97.9%; 100.0%]	99.3% [97.5%; 99.9%]

### 4. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Pediatric)

Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing a total of 93 positive nasal swab specimens and 318 negative specimen for SARS-CoV-2 antigen (Ag) from pediatric symptomatic and asymptomatic subjects between 0 and 15 years who were suspected of exposure to COVID-19 or who have presented with symptoms in the last 7 days, to have a sensitivity of 82.8% (95% CI: 73.6-89.8%) and a specificity of 100% (95% CI: 98.8-100%). Clinical specimens were determined to be positive or negative using a nasal swab specimen with an FDA EUA RT-PCR reference method. Lower Ct value corresponds to a higher virus concentration.

The specificity (n=318) was 100% with 95%CI [98.8-100%].

Overall results for sensitivity are summarized in the following table according to age:

		All Positive Samples	Ct values ≤ 33	Ct values ≤ 30
Sensitivity [CI 95%]	Total	82.8% [73.6%; 89.8%] (n=93)	87.2% [73.5%; 89.8%] (n=86)	93.2% [84.7%; 97.7%] (n=73)
	Age 0-5	79.4% [62.1%; 91.3%] (n=34)	87.1% [70.2%; 96.4%] (n=31)	91.3% [72.0%; 99.0%] (n=23)
	Age 6-10	84.8% [68.1%; 94.9%] (n=33)	87.0% [71.0%; 96.5%] (n=32)	93.3% [77.9%; 99.2%] (n=30)
	Age 11-15	84.6% [65.1%; 95.6%] (n=26)	87.0% [66.4%; 97.2%] (n=23)	95.0% [75.1%; 99.8%] (n=20)

The following table presents results for RT-PCR positive symptomatic and asymptomatic pediatric cohorts:

		All Nasal PCR Positive Samples	Ct values ≤ 33	Ct values ≤ 30
Sensitivity [CI 95%]	Symptomatic	87.0% [77.4%; 93.6%] (n=77)	91.5% [82.5%; 96.8%] (n=71)	95.1% [86.3%; 99.0%] (n=61)
	Asymptomatic	62.5% [35.4%; 84.8%] (n=16)	66.7% [38.4%; 88.2%] (n=15)	83.3% [51.6%; 97.9%] (n=12)

#### 5. Detection Limit

Panbio™ COVID-19 Ag Rapid Test Device was confirmed to detect  $2.5 \times 10^{1.8}$  TCID<sub>50</sub>/ml of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in Korea.

#### 6. Hook Effect

There is no hook effect at  $1.0 \times 10^{5.8}$  TCID<sub>50</sub>/ml of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in Korea.

#### 7. Cross Reactivity

Cross-reactivity of Panbio™ COVID-19 Ag Rapid Test Device was evaluated by testing 46 viruses and 21 other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Table below. The following viruses and other microorganisms except the Human SARS-coronavirus Nucleoprotein have no effect on the test results of Panbio™ COVID-19 Ag Rapid Test Device.

Panbio™ COVID-19 Ag Rapid Test Device has cross-reactivity with Human-SARS-coronavirus Nucleoprotein at a concentration of 25 ng/ml or more because SARS-CoV has high homology (79.6%) to the SARS-CoV-2.

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
1	Virus	Adenovirus Type 1	1.54 X 10 <sup>7</sup> PFU/ml	No cross reaction
2		Adenovirus Type 5	4.0 X 10 <sup>8</sup> PFU/ml	No cross reaction
3		Adenovirus Type 7	2.0 X 10 <sup>9</sup> PFU/ml	No cross reaction
4		Enterovirus (EV68)	2.0 X 10 <sup>7</sup> PFU/ml	No cross reaction
5		Echovirus2	7.0 X 10 <sup>5.5</sup> PFU/ml	No cross reaction
6		Echovirus11	3.5 X 10 <sup>6.25</sup> PFU/ml	No cross reaction
7		Enterovirus D68	2.0 X 10 <sup>7</sup> PFU/ml	No cross reaction
8		Human herpesvirus (HSV) 1	3.5 X 10 <sup>7.5</sup> PFU/ml	No cross reaction
9		Human herpesvirus (HSV) 2	3.5 X 10 <sup>5.75</sup> PFU/ml	No cross reaction
10		Mumps Virus Ag	1.1 X 10 <sup>5</sup> PFU/ml	No cross reaction
11		Influenza virus A (H1N1) Strain (A/Virginia/ATCC1/2009)	2.6 X 10 <sup>5</sup> PFU/ml	No cross reaction
12		Influenza virus A (H1N1) Strain (A/WS/33)	3.5 X 10 <sup>7.25</sup> PFU/ml	No cross reaction
13		Influenza virus A(H1N1) Strain (A/California/08/2009/pdm09)	1.1 X 10 <sup>8</sup> PFU/ml	No cross reaction
14		Influenza virus B Strain (B/Lee/40)	3.5 X 10 <sup>6.25</sup> PFU/ml	No cross reaction
15		Parainfluenza Type 1	2.1 X 10 <sup>8</sup> PFU/ml	No cross reaction

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
16	Virus	Parainfluenza Type 2	$3.5 \times 10^5$ PFU/ml	No cross reaction
17		Parainfluenza Type 3	$4.6 \times 10^7$ PFU/ml	No cross reaction
18		Parainfluenza Type 4A	$2.0 \times 10^7$ PFU/ml	No cross reaction
19		Respiratory syncytial virus (RSV) type A	$3.0 \times 10^5$ PFU/ml	No cross reaction
20		Respiratory syncytial virus (RSV) type B	$3.9 \times 10^5$ PFU/ml	No cross reaction
21		Rhinovirus A16	$8.8 \times 10^5$ PFU/ml	No cross reaction
22		HCoV-HKU1	1.5mg/ml	No cross reaction
23		HCoV-NL63	$1.2 \times 10^5$ PFU/ml	No cross reaction
24		HCoV-OC43	$6.2 \times 10^5$ PFU/ml	No cross reaction
25		HCoV-229E	$1.1 \times 10^6$ PFU/ml	No cross reaction
26		Human SARS-coronavirus Nucleoprotein	25 ng/ml	Cross Reaction
27		MERS-CoV Nucleoprotein	0.25 mg/ml	No cross reaction
28		Human Metapneumovirus (hMPV) 16 Type A1	$1.1 \times 10^6$ PFU/ml	No cross reaction
29		Adenovirus Type 2	$1.96 \times 10^7$ PFU/ml	No cross reaction
30		Adenovirus Type 3	$1.4 \times 10^{6.5}$ PFU/ml	No cross reaction
31	Adenovirus Type 4	$3.5 \times 10^{6.5}$ PFU/ml	No cross reaction	
32	Enterovirus C	$6.0 \times 10^7$ PFU/ml	No cross reaction	
33	Influenza virus A(H3N2) Strain (A/Hong Kong/8/68)	$3.5 \times 10^{5.5}$ PFU/ml	No cross reaction	

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
34	Virus	Influenza virus A(H5N1)	1.5 mg/ml	No cross reaction
35		Influenza virus B Strain (Victoria)	$5.46 \times 10^6$ PFU/ml	No cross reaction
36		Rhinovirus 14	$1.6 \times 10^8$ PFU/ml	No cross reaction
37		Human cytomegalovirus	$7.0 \times 10^5$ PFU/ml	No cross reaction
38		Norovirus	$7.14 \times 10^7$ PFU/ml	No cross reaction
39		Varicella-zoster virus	$1.96 \times 10^4$ PFU/ml	No cross reaction
40		Measles virus	$6.1 \times 10^5$ PFU/ml	No cross reaction
41		EB virus	$5.6 \times 10^8$ copies/ml	No cross reaction
42		Influenza virus(H7N9)	1.5mg/ml	No cross reaction
43		Influenza virus B Strain (Yamagata)	$2.73 \times 10^{10}$ PFU/ml	No cross reaction
44		Rhinovirus 54	$3.5 \times 10^{5.67}$ PFU/ml	No cross reaction
45		Rotavirus	$1.12 \times 10^7$ PFU/ml	No cross reaction
46		Adenovirus type 11	$3.0 \times 10^6$ PFU/ml	No cross reaction

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
1	Other Microorganism	<i>Staphylococcus saprophyticus</i>	7.9 X 10 <sup>7</sup> CFU/ml	No cross reaction
2		<i>Neisseria sp. (Neisseria lactamica)</i>	6.8 X 10 <sup>8</sup> CFU/ml	No cross reaction
3		<i>Staphylococcus haemolyticus</i>	1.4 X 10 <sup>10</sup> CFU/ml	No cross reaction
4		<i>Streptococcus salivarius</i>	7.84 X 10 <sup>7</sup> CFU/ml	No cross reaction
5		<i>Hemophilus parahaemolyticus</i>	8.8 X 10 <sup>8</sup> CFU/ml	No cross reaction
6		<i>Proteus vulgaris</i>	2.9 X 10 <sup>7</sup> CFU/ml	No cross reaction
7		<i>Moraxella catarrhalis</i>	1.9 X 10 <sup>8</sup> CFU/ml	No cross reaction
8		<i>Klebsiella pneumoniae</i>	2.0 X 10 <sup>7</sup> CFU/ml	No cross reaction
9		<i>Fusobacterium necrophorum</i>	7.0 X 10 <sup>8</sup> CFU/ml	No cross reaction
10		<i>Mycobacterium tuberculosis</i>	10mg/ml	No cross reaction
11		Pooled human nasal wash	N/A*	No cross reaction
12		<i>Streptococcus pyogenes</i>	3.6 X 10 <sup>7</sup> CFU/ml	No cross reaction
13		<i>Mycoplasma pneumoniae</i>	4.0 X 10 <sup>8</sup> CFU/ml	No cross reaction

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
14	Other Microorganism	<i>Staphylococcus aureus</i>	1.3 X 10 <sup>8</sup> CFU/ml	No cross reaction
15		<i>Escherichia coli</i>	6.8 X 10 <sup>6</sup> CFU/ml	No cross reaction
16		<i>Chlamydia pneumoniae</i>	9.1 X 10 <sup>7</sup> IFU/ml	No cross reaction
17		<i>Haemophilus influenzae</i>	3.4 X 10 <sup>8</sup> CFU/ml	No cross reaction
18		<i>Legionella pneumophila</i>	1.2 X 10 <sup>6</sup> CFU/ml	No cross reaction
19		<i>Streptococcus pneumoniae</i>	1.3 X 10 <sup>6</sup> CFU/ml	No cross reaction
20		<i>Bordetella pertussis</i>	4.4 X 10 <sup>9</sup> CFU/ml	No cross reaction
21		<i>Pneumocystis jirovecci</i> (PJP)	1.0 X 10 <sup>8</sup> nuclei/ml	No cross reaction

\* No concentration provided by supplier. Undiluted stock solution was tested.

## 8. Interfering Substances

The following 43 potentially interfering substances have no impact on Panbio™ COVID-19 Ag Rapid Test Device. The final test concentrations of the interfering substances are documented in the Table below.

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
1	Endogenous Substance	Mucin	0.5%	No Interference
2		Hemoglobin	100 mg/L	No Interference
3		Triglycerides	1.5 mg/L	No Interference
4		Icteric (Bilirubin)	40 mg/dL	No Interference
5		Rheumatoid factor	200 IU/ml	No Interference
6		Anti-nuclear antibody	>1:40	No Interference
7		Pregnant	10-fold dilution	No Interference
8	Exogenous Substance	Guaiacol glyceryl ether	1 µg/ml	No Interference
9		Albuterol	0.005 mg/dL	No Interference
10		Ephedrine	0.1 mg/ml	No Interference
11		Chlorpheniramine	0.08 mg/dL	No Interference
12		Diphenhydramine	0.08 mg/dL	No Interference
13		Ribavirin	26.7 µg/ml	No Interference
14		Oseltamivir	0.04 mg/dL	No Interference

No.	Types of Specimen	Tested substance	Final Test Concentration	Test Result
15	Exogenous Substance	Zanamivir	17.3 µg/ml	No Interference
16		Phenylephrine hydrochloride	15% v/v	No Interference
17		Oxymetazolin hydrochloride	15% v/v	No Interference
18		Amoxicillin	5.4 mg/dL	No Interference
19		Acetylsalicylic acid	3 mg/dL	No Interference
20		Ibuprofen	21.9 mg/dL	No Interference
21		Chlorothiazide	2.7 mg/dL	No Interference
22		Indapamide	140 ng/ml	No Interference
23		Glimepiride (Sulfonylureas)	0.164 mg/dL	No Interference
24		Acarbose	0.03 mg/dL	No Interference
25		Ivermectin	4.4 mg/L	No Interference
26		Lopinavir	16.4 µg/L	No Interference
27		Ritonavir	16.4 µg/L	No Interference
28		Chloroquine phosphate	0.99 mg/L	No Interference
29		Sodium chloride with preservatives	4.44 mg/ml	No Interference
30		Beclomethasone	4.79 ng/ml	No Interference
31		Dexamethasone	0.6 µg/ml	No Interference
32		Flunisolide	0.61 µg/ml	No Interference
33		Triamcinolone	1.18 ng/ml	No Interference
34		Budesonide	2.76 ng/ml	No Interference
35		Mometasone	1.28 ng/ml	No Interference
36		Fluticasone	2.31 ng/ml	No Interference
37		Sulfur	9.23 µg/ml	No Interference
38		Benzocaine	0.13 mg/ml	No Interference
39		Menthol	0.15 mg/ml	No Interference
40		Mupirocin	10 µg/ml	No Interference
41		Tobramycin	24.03 µg/ml	No Interference
42		Biotin	1.2 µg/ml	No Interference
43	HAMA	63.0 ng/ml	No Interference	


## 9. **Repeatability & Reproducibility**

Repeatability & Reproducibility of Panbio™ COVID-19 Ag Rapid Test Device was established using in-house reference panels containing negative specimens and a range of positive specimens. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.


## PREPARATION

- 1 Allow all kit components to reach a temperature between 15-30 °C prior to testing for 30 minutes.  
**Note:** Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.
- 2 **Open the package and look for the following:**
  1. Test device with desiccant in individual foil pouch
  2. Buffer tube
  3. Sterilized nasal swabs for sample collection
  4. Tube rack
  5. Instructions for Use with Quick Reference Guide (Nasal)
- 3 Carefully read these instructions prior to using Panbio™ COVID-19 Ag Rapid Test Device kit.
- 4 Look at the expiration date of the kit box. If the expiration date has passed, use another kit.
- 5 **Open the foil pouch and look for the following:**
  1. Result window
  2. Specimen wellThen, label the device with the patient identifier.

◆◆◆◆◆: 5 drops of the extracted specimen


- 6 Keep buffer tube upright and remove foil. Place the buffer tube in the tube rack before proceeding to the next step.  
**Note:** Be careful not to spill the tube contents.


## TEST PROCEDURE

- 1 Tilt the patient's head back 70 degrees. While gently rotating the swab, insert swab less than one inch (about 2 cm) into nostril (until resistance is met at the turbinates). Rotate the swab five times against the nasal wall. Using the same swab repeat the collection procedure with the second nostril. Slowly remove swab from the nostril.

**⚠ Caution:** If the swab stick breaks during specimen collection, repeat specimen collection with a new swab.


- 2 Insert the swab specimen in the buffer tube. Swirl the swab tip in the buffer fluid inside the buffer tube, pushing into the wall of the buffer tube at least five times and then squeeze out the swab by squeezing the buffer tube with your fingers.


## TEST PROCEDURE

- 3 Break the swab at the breakpoint and close the buffer tube with the attached nozzle.

- 4 Dispense 5 drops of extracted specimens vertically into the specimen well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.

**⚠ Caution:** Bubbles that occur in the buffer tube can lead to inaccurate results. If you are unable to create sufficient drops, this may be caused by clogging in the dispensing nozzle. Shake the tube gently to release the blockage until you observe free drop formation.

## TEST PROCEDURE

- 5** Dispose of the buffer tube containing the used swab according to your local regulations and biohazard waste disposal protocol.

- 6** Start timer. Read result at 15 minutes. Do not read results after 20 minutes.


- 7** Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

## TEST INTERPRETATION

### NEGATIVE


The presence of only the control line (C) and no test line (T) within the result window indicates a negative result.

### POSITIVE


The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a positive result.

**⚠ Caution:** The presence of any test line (T), no matter how faint, indicates a positive result.

### INVALID

If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly. It is recommended to read the IFU again before re-testing the specimen with a new test device.

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
- MEMO -

- MEMO -

Panbio™  
**COVID-19 Ag RAPID TEST DEVICE** (NASAL)  
 (NASAL)

**Abbott**

**REF** 41FK11CON

**Panbio™  
COVID-19 Ag  
Rapid Test Device  
Positive  
Control Kit**

## Introduction

The Panbio™ COVID-19 Ag Rapid Test Device Positive Control Kit are specifically formulated and manufactured to ensure performance of the Panbio™ COVID-19 Ag Rapid Test Device and used to verify the user's ability to properly perform the test and interpret the result. The Positive Control will produce positive results and has been manufactured to produce visible test line (T) and control line (C).

## Intended Use

Panbio™ COVID-19 Ag Rapid Test Device Positive Control Kit contains positive quality control swabs for use with Panbio™ COVID-19 Ag Rapid Test Device only.

It is the responsibility of each laboratory using Panbio™ COVID-19 Ag Rapid Test Device to establish an adequate quality assurance program to ensure the performance of the device under its specific locations and conditions of use.

Perform the Panbio™ COVID-19 Ag Rapid Test Device Control Kit at periodic intervals as indicated by user facility, such as:

- When training a new operator prior to them performing tests on patient specimens
- When opening a new test lot
- Whenever a new shipment of test kits is received

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- Panbio™ COVID-19 Ag Rapid Test Device kit
- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
- Timer
- Biohazard container

## Storage and Stability

1. The control kit should be stored at a temperature between 2-30 °C. Do not freeze the kit.  
**Note:** When stored in a refrigerator, the control kit must be brought to room temperature (15-30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
2. Perform the test immediately after removing the positive control swab from the foil pouch.
3. Do not use the control kit beyond its expiration date.

4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the control kit if the pouch is damaged or the seal is broken.

## Warnings

1. For *in vitro* diagnostic use only.
2. Do not use the positive control swab for specimen collection.
3. These instructions must be strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instruction prior to performing a test.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Do not use Positive Controls for other products.
7. Do not store the test kit in direct sunlight.
8. To avoid contamination, do not touch the head of provided swab when opening the swab pouch.
9. Do not dilute the positive control swab with any solution except for the provided buffer.

## Limitations

Panbio™ COVID-19 Ag Rapid Test Device Positive Control Kit contains quality control swabs for use with Panbio™ COVID-19 Ag Rapid Test Device only.

## Test Procedure (Refer to Figures on Other Side)

### Test Preparation

1. Allow all test device and control kit components to reach a temperature between 15-30 °C prior to testing for 30 minutes.
2. Remove the test device from the foil pouch prior to use. Place on a flat, horizontal and clean surface.

### Reaction with Test Device

1. Keep the buffer tube upright and remove the foil. Place the buffer tube in the tube rack.  
**Note:** Be careful not to spill the tube contents.
2. Remove the positive control swab from the foil pouch and insert it into buffer tube.
3. Swirl the swab tip in the buffer fluid inside the buffer tube, pushing into the wall of the buffer tube at least five times and then squeeze out the swab by squeezing the buffer tube with your fingers.
4. Dispose of the used control swab in accordance with your local regulations and biohazard waste disposal protocol.

5. Press the nozzle cap tightly onto the buffer tube.
6. Dispense 5 drops of extracted specimens vertically into the specimen well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.
7. Dispose of the buffer tube according to your local regulations and biohazard waste disposal protocol.
8. Start timer. Read result at 15 minutes. Do not read result after 20 minutes.
9. Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

## Test Interpretation (Refer to Figures on Other Side)

1. **Positive result:** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a positive result.  
**Note:** All lines should appear.
2. **Invalid result:** If the control line (C) is not visible within the result window after performing the test and/or the presence of only the control line (C) and no test line (T) within the result window are considered invalid. Instructions may not have been followed correctly. It is recommended to read the IFU again before re-testing the positive control swab with a new test device.

## Quality Control

The Positive Control is specifically formulated and manufactured to ensure performance of the Panbio™ COVID-19 Ag Rapid Test Device and are used to verify the user's ability to properly perform the test and interpret the result. The Positive Control will procedure positive test results and has been manufactured to produce visible test line (T) and control line (C).

Good laboratory practice suggests the use of Controls to ensure that:

- Test reagents are working,
- The test is correctly performed

The Controls can be run under any of the following circumstances:

- When training a new operator prior to them performing tests on patient specimens
- When opening a new test lot
- Whenever a new shipment of test kits is received

# PREPARATION

- 1 Allow all test and control kit components to reach a temperature between 15-30 °C prior to testing for 30 minutes.

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- 2 **Open the package and look for the following:**


1. Positive Control Swab
2. Instructions for Use with Quick Reference Guide

- 3 Carefully read these instructions prior to using Panbio™ COVID-19 Ag Rapid Test Device Control Kit.

- 4 Look at the expiration date of the kit box. If the expiration date has passed, use another kit.

- 5 **Open the foil pouch and look for the following:**

1. Result window
2. Specimen well (S)


●●●●● : 5 drops of the extracted specimen

## TEST PROCEDURE

- 1 Keep the buffer tube upright and remove the foil. Place the buffer tube in the tube rack.


**Note:** Be careful not to spill the tube contents.

- 2 Remove the positive control swab from the foil pouch and insert it into buffer tube.


Swirl the swab tip in the buffer fluid inside the buffer tube, pushing into the wall of buffer tube at least five times and then squeeze out the swab by squeezing the buffer tube with your fingers.

**Note:** Dispose of the used control swab in accordance with your local regulations and biohazard waste disposal protocol.


# 3

Press the nozzle cap tightly onto the buffer tube.

# 4

Dispense 5 drops of extracted specimens vertically into the specimen well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.

# 5

Dispose of the buffer tube according to your local regulations and biohazard waste disposal protocol.

# 6

Start timer. Read result at 15 minutes. Do not read result after 20 minutes.

**WAIT 15 MINUTES**

**READ AT  
15-20 MINUTES**

# 7


Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

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**Note:** All lines should appear.

### INVALID

If the control line (C) is not visible within the result window after performing the test and/or the presence of only the control line (C) and no test line (T) within the result window are considered invalid. Instructions may not have been followed correctly. It is recommended to read the IFU again before re-testing the positive control swab with a new test device.

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