

**WHO Emergency Use Assessment SARS-CoV-2 IVDs  
PUBLIC REPORT**

**Product: Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)**

**Manufacturer: Guangzhou Wondfo Biotech Co., Ltd**

**EUL Number: EUL 0676-004-00**

**Outcome: Accepted**

The EUL process is intended to expedite the availability of in vitro diagnostics needed in public health emergency situations and to assist interested UN procurement agencies and Member States in determining the acceptability of using specific products in the context of a Public Health Emergency of International Concern (PHEIC), based on an essential set of available quality, safety and performance data. The EUL procedure includes the following:

- Quality Management Systems Review and Plan for Post-Market Surveillance: desk-top review of the manufacturer's Quality Management System documentation and specific manufacturing documents.
- Product Dossier Review: assessment of the documentary evidence of safety and performance. This evaluation of limited scope is to verify critical analytical and performance characteristics.

Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) with product codes W634P0013, W634P0014, W634P0015, rest of world regulatory version manufactured by Guangzhou Wondfo Biotech Co., Ltd, No. 8 Lizhishan Road, Science City, Luogang District, 510663 Guangzhou, China, was listed as eligible for WHO procurement on 18 July 2022.

**Intended use**

According to the claim of intended use from Guangzhou Wondfo Biotech Co., Ltd, "*Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an immunochromatographic assay for rapid, qualitative detection of novel coronavirus antigen extracted from the nasal swab specimen collected from the individuals who are suspected of novel coronavirus infection, with or without symptoms. The test is to be used as an aid in the diagnosis of coronavirus disease (COVID-19), which is caused by novel coronavirus. Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is for professional use only. The test may be used in any laboratory and non-laboratory settings that meet the requirements specified in the instructions*

*for use and local regulation. The test provides preliminary test results. Negative results cannot exclude novel coronavirus infection and they cannot be used as the sole basis for treatment or other management decision. For in vitro diagnostic use only. ”*

**Specimen type that was validated:** Nasal swab specimens.

**Assay description**

According to the claim of assay description from Guangzhou Wondfo Biotech Co., Ltd, *“Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is based on the principle of Immunochromatography, using colloidal gold as the label for detection of novel coronavirus antigen extracted from the nasal swab specimen. When the extracted specimen is added into the test device, the specimen is absorbed into the device by capillary action, reacts with the novel coronavirus antibody-dye conjugate and flows across the precoated membrane.*

*When the novel coronavirus antigen level in the specimen is at or above the target cutoff (the detection limit of the test), it will react with the antibody labeled with colloidal gold to form a complex. Through capillary action, the reaction complex moves forward along the nitrocellulose membrane to the detection region (T), where it reacts with the precoated antibodies on the nitrocellulose membrane to form an antibody-antigen-antibody complex, and this produces a colored test band that indicates a positive result. When the novel coronavirus antigen level in the specimen is zero or below the target cutoff, there is no a visible colored band in the detection region (T) of the device. This indicates a negative result.*

*To serve as a procedure control, a colored band will appear at the control region (C) if the test has been performed properly.”*

**Test kit contents**

<b>Components</b> \ <b>Catalog No.</b>	<b>W634P0013</b>	<b>W634P0014</b>	<b>W634P0015</b>
Sealed Pouches* (pcs)	20	100	20
Extraction Buffer (6mL/vial)	2	10	2
Sample Extraction Tube (pcs)	20	100	20
Drippers (pcs)	20	100	20
Nasal Swab (pcs)	20	100	20
Positive Control Swab (pcs)	\	1	1
Negative Control Swab (pcs)	\	1	1
IFU (pcs)	1	1	1

**Note:** \*Each sealed pouch contains: 1 Test Cassette and 1 Desiccant Pouch

### Items required but not provided

1. Timer
2. Test tube rack
3. Personal protective equipment, such as protective gloves, medical mask, goggles and lab coat.
4. Appropriate biohazard waste container and disinfectants.

### Storage

The test kit should be stored at 2-30 °C.

### Shelf-life upon manufacture

Nine months (real-time stability studies are ongoing)

### Warnings/limitations

Please refer to the attached instructions of use (IFU).

## Product dossier assessment

Guangzhou Wondfo Biotech Co., Ltd submitted a product dossier for the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) as per the “*Instructions for Submission Requirements: In vitro diagnostics (IVDs) Detecting SARS-CoV-2 nucleic acid or antigen (PQDx\_0347)*”. The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and an external assessor appointed by WHO.

### Post listing Commitments for EUL

As a requirement for listing, the manufacturer is required to:

1. When available, assess the traceability of the materials used to validate the product with the WHO International Standard for SARS-CoV-2 antigen and provide the report to WHO.
2. Estimate the limit of detection (LoD) of the product with the WHO International Standard for SARS-CoV-2 antigen when available and provide the report to WHO within one month of completion of the study.
3. Partake in an independent performance evaluation conducted by a laboratory commissioned by WHO. Any such performance evaluation testing will be performed using the protocol and technical criteria established by WHO.
4. Provide evidence supporting the equivalence of inactivated virus samples with samples containing the live virus, per PQDx\_347 requirements.

5. Provide a study protocol and report characterizing the inactivated virus stocks by PCR to be consistent with the PQDx\_347 requirements, and the stock concentration will be given in copies/mL.
6. Provide clinical evidence supporting the required proportion of samples collected at different times post-onset symptoms, consistent with PQDx\_347 requirements. Additionally, Wondfo commits to providing the clinical protocol/reports of the clinical studies and the claimed performance in IFU to be adjusted accordingly.
7. Provide a revised IFU including the in-use (after opened) stability claim for the extraction buffer to include the claim supported by the conclusion presented in the respective stability study report.
8. Provide a revised IFU, which will include information about enrolled subjects in the clinical performance to match the information provided in the trial design outlined in section 4 of the clinical studies.

Risk-benefit assessment conclusion is acceptable.

### Quality Management Systems Review

To establish the eligibility for WHO procurement, Guangzhou Wondfo Biotech Co., Ltd was asked to provide up-to-date information about the status of their quality management system.

Based on the review of the submitted quality management system documentation by WHO staff and external technical experts (assessors), it was established that sufficient information was provided by Guangzhou Wondfo Biotech Co., Ltd to fulfil the requirements described in the *“Instructions for Submission Requirements: In vitro diagnostics (IVDs) Detecting SARS-CoV-2 nucleic acid or antigen( PQDx\_ 347)”*.

The quality management documentation assessment conclusion is acceptable.

### Plan for Post-Market Surveillance

Post-market surveillance, including monitoring all customer feedback, detecting and acting on adverse events, product problems, non-conforming goods, and processes is a critical component of minimizing potential harm of an IVD listed for emergency use.

The following post-prequalification activities are required to maintain the prequalification status:

1. Notification to WHO of any planned changes to a prequalified product, in accordance with “*WHO procedure for changes to a WHO prequalified in vitro diagnostic*” (document number PQDx\_121); and
2. Post-market surveillance activities, in accordance with “*Guidance for post-market surveillance and market surveillance of medical devices, including in vitro diagnostics*” (ISBN 978-92-4-001531-9).

Guangzhou Wondfo Biotech Co., Ltd is also required to report complaints related to the product. There are certain categories of complaints and changes to the product that must be notified immediately to WHO, as per the above-mentioned documents.

The manufacturer has committed to ensuring that post-emergency use listing safety, quality, and performance monitoring activities are in place, which are in accordance with WHO guidance “*Guidance for post-market surveillance and market surveillance of medical devices, including in vitro diagnostics.*”<sup>1</sup>

### **Scope and duration of procurement eligibility**

The Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) with product codes W634P0013, W634P0014, and W634P0015, manufactured by Guangzhou Wondfo Biotech Co., Ltd, is considered eligible for WHO procurement for 12 months from the day of listing. The assay may detect the Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) antigen. This listing does not infer that the product meets WHO prequalification requirements and does not mean that the product is listed as WHO prequalified.

As part of the ongoing requirements for listing as eligible for WHO procurement, Guangzhou Wondfo Biotech Co., Ltd must engage in post-market surveillance activities to ensure that the product meets safety, quality, and performance requirements, Guangzhou Wondfo Biotech Co., Ltd must notify WHO of any complaints, including adverse events related to the use of the product, within 7 days and any changes made to the product.

WHO reserves the right to rescind eligibility for WHO procurement if additional information on the safety, quality, and performance during post-market surveillance activities and if new data becomes available to WHO that changes the risk-benefit balance.

---

<sup>1</sup> <https://www.who.int/publications/i/item/9789240015319>

## **Labelling**


- 1. Labels**
- 2. Instructions for use**

## **1.0 Labels**

### **1.1 Product code W634P0013 (20 T/kit)**

#### **1.1.1. Kit box design**

**Wondfo**  
Leading POCT Manufacturer



**Wondfo**  
Leading POCT Manufacturer

**FOR IN VITRO DIAGNOSTIC USE ONLY**  
**FOR PROFESSIONAL USE ONLY**



**Wondfo**  
Leading POCT Manufacturer



Rev. A1 Rel.:2021/06/28



**Wondfo**  
Leading POCT Manufacturer

REF W634P0013  
LOT XXXXXXXXX  
XXXX-XX-XX  
XXXX-XX-XX



**Wondfo**  
Leading POCT Manufacturer

We Are Working For Your Health


IVL

Guangzhou Wondfo Biotech Co., Ltd.  
No.8 Lishihe Road, Shimen Town, Dongying District, 510663, Guangzhou, P.R.China  
Tel: +86-20-33222919 Fax: +86-20-33222923  
E-mail: gk@wondfo.com.cn Website: www.wondfo.com.cn

Guangzhou Wondfo Biotech Co., Ltd.  
No.8 Lishihe Road, Shimen Town, Dongying District, 510663, Guangzhou, P.R.China  
Tel: +86-20-33222919 Fax: +86-20-33222923  
E-mail: gk@wondfo.com.cn Website: www.wondfo.com.cn



### 1.1.2 kit box label



### 1.1.3 Test cassette pouch

18mm

**Wondfo**  
Leading POCT Manufacturer



## Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)

**LOT** XXXXXXXXX

 XXXX-XX-XX

**REF** W634P0013

 XXXX-XX-XX

Rev. A2 Rel.: 2021.03.03

**FOR IN VITRO DIAGNOSTIC USE ONLY**

**Wondfo**




 **Guangzhou Wondfo Biotech Co., Ltd.**  
No.8 Lizhishan Road, Science City, Luogang  
District, 510663, Guangzhou, P.R.China  
Tel: +86-20-32053962 Fax: +86-20-32296063  
E-mail: [global@wondfo.com.cn](mailto:global@wondfo.com.cn) Website: [www.wondfo.com.cn](http://www.wondfo.com.cn)

### 1.1.4 Round sticker




### 1.1.5 Extraction buffer label




### 1.1.6 Sample extraction tube + dripper label

Front



Back



### **1.1.7 Nasal swab label**

38.0 mm

25.0 mm



PEEL HERE

Disposable swabs

CF 150-P 3 B

STERILE

EO



- 1, Single use only.
- 2, There is a risk of cross- infection if re-use
- 3, Do not use if package is damaged

EC REP

CF101 Rev.00 Date201701

Lins Service & Consulting GmbH  
Obere Seegasse 34/2,69124  
Heidelberg, Germany  
Email: info@lins-service.com

6.5

CE 0197

LOT 20210206JZ



2024/02/05



Jiangsu Changfeng Medical  
Industry Co., Ltd  
Touqiao Town, Guangling  
District, Yangzhou,  
Jiangsu 225109 China


10

205.0 mm

185.0 mm

## **1.2 Product code W634P0014 (100 T/kit)**

### **1.2.1 Kit box design**



**Wondfo**  
Leading POCT Manufacturer



Guangzhou Wondfo Biotech Co., Ltd.  
No. 8 Lishiwan Road, Science City, Luogang District, 510663, Guangzhou, P.R.China  
Tel: +86-20-32032662 Fax: +86-20-32036663  
E-mail: global@wondfo.com.cn Website: www.wondfo.com.cn

**Wondfo**  
Leading POCT Manufacturer

**FOR *IN VITRO* DIAGNOSTIC USE ONLY**  
**FOR PROFESSIONAL USE ONLY**

**Wondfo**  
Leading POCT Manufacturer



Guangzhou Wondfo Biotech Co., Ltd.  
No. 8 Lishiwan Road, Science City, Luogang District, 510663, Guangzhou, P.R.China  
Tel: +86-20-32032662 Fax: +86-20-32036663  
E-mail: global@wondfo.com.cn Website: www.wondfo.com.cn

Rev. A1 Rel:2021/04/15

**Wondfo**  
Leading POCT Manufacturer



WB34P0014  
XXXXXXXXXX  
XXXX-XX-XX  
XXXX-XX-XX

**FOR *IN VITRO* DIAGNOSTIC USE ONLY**  
**FOR PROFESSIONAL USE ONLY**


Health

IVD

**Wondfo**  
Leading POCT Manufacturer



### 1.2.2 Kit box label



### 1.2.3 Test cassette pouch

18mm

**Wondfo**  
Leading POCT Manufacturer



## Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)

**LOT** XXXXXXXXXX

 XXXX-XX-XX

**REF** W634P0014

 XXXX-XX-XX

Rev. A2 Rel.: 2021.03.03

**FOR IN VITRO DIAGNOSTIC USE ONLY**

**Wondfo**



**Guangzhou Wondfo Biotech Co., Ltd.**

No.8 Lizhishan Road, Science City, Luogang

District, 510663, Guangzhou, P.R.China


Tel: +86-20-32053962 Fax: +86-20-32296063

E-mail: [global@wondfo.com.cn](mailto:global@wondfo.com.cn) Website: [www.wondfo.com.cn](http://www.wondfo.com.cn)

### 1.2.4 Round sticker




### 1.2.5 Extraction buffer label




### 1.2.6 Sample extraction tube + dripper label

Front



Back



### **1.2.7 Nasal swab label**

38.0 mm

25.0 mm



PEEL HERE

Disposable swabs

CF 150-P 3 B

STERILE

EO



- 1, Single use only.
- 2, There is a risk of cross- infection if re-use
- 3, Do not use if package is damaged

CF101 Rev.00 Date201701

EC REP

Lins Service & Consulting GmbH  
Obere Seegasse 34/2,69124  
Heidelberg, Germany  
Email: info@lins-service.com

6.5

CE 0197

LOT 20210206JZ



2024/02/05




Jiangsu Changfeng Medical  
Industry Co., Ltd  
Touqiao Town, Guangling  
District, Yangzhou,  
Jiangsu 225109 China

10


205.0 mm

185.0 mm

### 1.2.8 Positive control swab label



### 1.2.9 Negative control swab label



### **1.3 Product code W634P0015 (20 T/kit)**

#### **1.3.1 Kit box design**



**Wondfo**  
Leading POCT Manufacturer



**Wondfo**  
Leading POCT Manufacturer




**Wondfo**  
Leading POCT Manufacturer



**Wondfo**  
Leading POCT Manufacturer

REF W634P0015  
LOT XXXXXXXXX  
XXXX-XX-XX  
XXXX-XX-XX



**FOR IN VITRO DIAGNOSTIC USE ONLY  
FOR PROFESSIONAL USE ONLY**

Guangzhou Wondfo Biotech Co., Ltd.  
No.8 Lishihe Road, Science City, Tieding District, 510663, Guangzhou, P.R.China  
Tel: +86-20-33222919 Fax: +86-20-33222923  
E-mail: gk@wondfo.com.cn Website: www.wondfo.com.cn

Guangzhou Wondfo Biotech Co., Ltd.  
No.8 Lishihe Road, Science City, Luogang District, 510663, Guangzhou, P.R.China  
Tel: +86-20-33222919 Fax: +86-20-33222923  
E-mail: gk@wondfo.com.cn Website: www.wondfo.com.cn


Rev. A1 Rel.:2021/06/28

We Are Working For Your Health

IVL

**Wondfo**  
Leading POCT Manufacturer

### 1.3.2 kit box label



### 1.3.3 Test cassette pouch

18mm

**Wondfo**  
Leading POCT Manufacturer



## Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)

**LOT** XXXXXXXXXX

 XXXX-XX-XX

**REF** W634P0015

 XXXX-XX-XX

Rev. A2 Rel.: 2021.03.03

**FOR IN VITRO DIAGNOSTIC USE ONLY**

**Wondfo**



**Guangzhou Wondfo Biotech Co., Ltd.**

No.8 Lizhishan Road, Science City, Luogang

District, 510663, Guangzhou, P.R.China


Tel: +86-20-32053962 Fax: +86-20-32296063

E-mail: [global@wondfo.com.cn](mailto:global@wondfo.com.cn) Website: [www.wondfo.com.cn](http://www.wondfo.com.cn)

### 1.3.4 Round sticker




### 1.3.5 Extraction buffer label




### 1.3.6 Sample extraction tube + dripper label

Front



Back



### **1.3.7 Nasal swab label**

38.0 mm

25.0 mm



PEEL HERE

Disposable swabs

CF 150-P 3 B

STERILE

EO



- 1, Single use only.
- 2, There is a risk of cross- infection if re-use
- 3, Do not use if package is damaged

EC REP

CF101 Rev.00 Date201701

Lins Service & Consulting GmbH  
Obere Seegasse 34/2,69124  
Heidelberg, Germany  
Email: info@lins-service.com

6.5

CE 0197

LOT 20210206JZ



2024/02/05




Jiangsu Changfeng Medical  
Industry Co., Ltd  
Touqiao Town, Guangling  
District, Yangzhou,  
Jiangsu 225109 China

10


205.0 mm

185.0 mm

### 1.3.8 Positive control swab



### 1.3.9 Negative control swab





## **2.0 Instructions for use (IFU)<sup>2</sup>**

---

<sup>2</sup>

English version of the IFU was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

### CATALOG NO.

W634P0013 W634P0014 W634P0015

### INTENDED USE

Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an immunochromatographic assay for rapid, qualitative detection of novel coronavirus antigen extracted from the nasal swab specimen collected from the individuals who are suspected of novel coronavirus infection, with or without symptoms. The test is to be used as an aid in the diagnosis of coronavirus disease (COVID-19), which is caused by novel coronavirus.

Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is for professional use only. The test may be used in any laboratory and non-laboratory settings that meet the requirements specified in the instructions for use and local regulation. The test provides preliminary test results. Negative results cannot exclude novel coronavirus infection and they cannot be used as the sole basis for treatment or other management decision.

For *in vitro* diagnostic use only.

### SUMMARY

COVID-19 is an acute respiratory infectious disease. Novel coronavirus transmission occurs with high efficacy and infectivity mainly through the respiratory route, such as oral fluid, sneezing, physical contact, and other air droplets [1-3]. Novel coronavirus can be transmitted by symptomatic or asymptomatic infected individuals [4]. The current epidemiological investigation suggests a mean incubation period of 5 to 7 days and a median incubation period of 3 days (range: 0–24 days) [5]. In symptomatic patients, the clinical manifestations of the disease usually start after less than a week, consisting of fever, cough, nasal congestion, fatigue, and other signs of upper respiratory tract infections [6]. The availability of a cost-effective, rapid point of care diagnostic test is critical to enable healthcare professionals to aid in the diagnosis of patients and prevent further spread of the virus. As a rapid test with a 15-20 min testing time, Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an useful diagnosis tool for coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

### PRINCIPLE

Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is based on the principle of Immunochromatography, using colloidal gold as the label for detection of novel coronavirus antigen extracted from the nasal swab specimen. When the extracted specimen is added into the test device, the specimen is absorbed into the device by capillary action, reacts with the novel coronavirus antibody-dye conjugate and flows across the pre-coated membrane.

When the novel coronavirus antigen level in the specimen is at or above the target cutoff (the detection limit of the test), it will react with the antibody labeled with colloidal gold to form a complex. Through capillary action, the reaction complex moves forward along the nitrocellulose membrane to the detection region (T), where it reacts with the pre-coated antibodies on the nitrocellulose membrane to form an antibody-antigen-antibody complex, and this produces a colored test band that indicates a positive result. When the novel coronavirus antigen level in the specimen is zero or below the target cutoff, there is no a visible colored band in the detection region (T) of the device. This indicates a negative result.

To serve as a procedure control, a colored band will appear at the control region (C), if the test has been performed properly.

### PRECAUTION

1. This kit is for *in vitro* diagnostic use only.
2. All specimens should be treated as contagious. Use appropriate precautions in the collection, handling, storage and disposal of patient specimens and used kit contents.
3. Wear appropriate personal protective equipment (e.g. protective gloves, medical mask, goggles and lab coat) when handling the contents of this kit.
4. Proper specimen collection, storage, and transport are critical to the performance of this test.
5. Discard after first use. Do not reuse the test device and kit components.
6. Avoid excessively high or low temperature in the experiment environment. When stored in a refrigerator, all kit components must be brought to room temperature for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
7. **Do not** touch the reaction area of the test strip.
8. **Do not** use if the test kit is beyond expiration date.
9. **Do not** use the kit if the pouch is punctured or not well sealed.
10. **Do not** dilute the collected swab with any solution except with the provided extraction buffer.
11. Components from different lots must not be mixed or used together.

12. The test result should be interpreted by the professional along with clinical findings and other laboratory test results.
13. **DISPOSAL OF THE USED KIT:** All specimens and the used kit have the infectious risk. The process of disposing used kit must follow the local infectious disposal law or laboratory regulation.

### MATERIALS

#### Materials Provided

Components	REF	W634P0013	W634P0014	W634P0015
Sealed Pouches* (pcs)		20	100	20
Extraction Buffer (6mL/vial)		2	10	2
Sample Extraction Tube (pcs)		20	100	20
Drippers (pcs)		20	100	20
Nasal Swab (pcs)		20	100	20
Positive Control Swab (pcs)	/		1	1
Negative Control Swab (pcs)	/		1	1
IFU (pcs)		1	1	1

*Note: \*Each sealed pouch contains: 1 Test Cassette and 1 Desiccant Pouch*

#### Materials Required but Not Provided

1. Timer
2. Test tube rack
3. Personal protective equipment, such as protective gloves, medical mask, goggles and lab coat.
4. Appropriate biohazard waste container and disinfectants.


### STORAGE AND STABILITY

1. The test kit should be stored at a temperature between 2-30 °C. The validity period of this product is 24 months. **Do not** freeze.
2. The test cassette should be used within 1 hour after taking out from the sealed pouch. Bottled buffer solution should be re-capped after use. Store the opened buffer at 2-30 °C for no more than 10 weeks.
3. Keep away from sunlight, moisture, and heat.
4. Kit contents are stable until the expiration date printed on the outer package.
5. The manufacturing date is printed on the outer package.

### SPECIMEN COLLECTION AND PREPARATION

Bring the test components to room temperature. Take out the nasal swab from the notch of the package.

1. Tilt the head of the patient backwards (about 70 degrees).
2. Gently twist the swab, insert the entire absorbent tip of the nasal swab into a nostril to about 1.5 cm deep.
3. Perform the first sampling by rubbing the nasal wall firmly with the nasal swab, turning it five times against the nasal walls so that the absorbent surface of the nasal swab is wetted all round.  
**Note:**  
a) *Be careful not to hurt the patient.*  
b) *This process may take about 15 seconds.*
4. Slowly remove the nasal swab from the first nostril.
5. Repeat the collection process with the same nasal swab in the other nostril.



Nasal Swab Collection

**Note:** Simply twirling the swab against one part of the inside of the nose or leaving the swab in the nose for 15 seconds is **not** a proper technique and may result in an insufficient specimen.


### TEST PROCEDURE

Please read the instructions for use carefully before performing the test.

#### Nasal Swab Testing: For W634P0013/W634P0014/ W634P0015

1. Transfer 10 drops (about 400  $\mu$ L) extraction buffer to the sample extraction tube, holding buffer vial vertically.
2. Insert the moist swab tip with nasal secretions into the sample extraction tube; place the nasal swab tip against the bottom of the extraction tube and rotate gently for 10 times to release the

3. secretion from the tip of the nasal swab. Place the sample extraction tube on the test tube rack (if applicable) and leave the swab in the extraction buffer for 1 minute.
4. While squeezing the middle of the extraction tube, slowly pull out the nasal swab to extract as much liquid as possible from the tip of the nasal swab. Discard the used swab in accordance with the biohazard waste disposal protocol.
5. Tightly close the extraction tube with dropper.  
**Note:** Swab specimens should be tested immediately after collection. If it is not possible to test immediately, the swab specimen can be kept in an extraction tube filled with extraction buffer at room temperature (10–40 °C) for up to 2 hours prior to testing. And used within 3 hours when stored at 2 °C ~8 °C.
6. Remove a test cassette from the sealed pouch by tearing at the notch and place it on a level surface.
7. Turn the sample extraction tube upside down, hold it vertically and add 4 drops (about 80  $\mu$ L) processed specimen into the small, round, white sample well. The result window (C/T) now turns purple.
8. Start the timer. Read the results at 15–20 minutes. **Do not read results after 20 minutes.**



#### Positive/Negative Control Swab Testing: For W634P0014/ W634P0015

**Caution:** The positive and negative control swab are for control use only. **Do not** use the control swabs for specimen collection.

Transfer about 10 drops (about 400  $\mu$ L) extraction buffer to the sample extraction tube vertically. Take out the positive or negative swab and insert the swab into the extraction buffer, rotate the swab tip 10 times against the bottom and sides of the extraction tube to release the specimen from the swab tip. Return the sample extraction tube to the test tube rack (if applicable) and leave the swab in the extraction buffer for 1 minute. Then follow the above test procedure c-g of [Nasal Swab Testing]

#### Note:

1. Please refer to the **RESULT INTERPRETATION** section of this instructions for use to interpret the results.
2. Please refer to the **External Quality Control** section of this instructions for use for the frequency of testing control swabs.

### RESULT INTERPRETATION

#### Positive Result


Colored bands appear at both test line (T) and control line (C) indicating a positive result for the novel coronavirus antigen in the specimen.

#### Negative Result

Colored band appears at control line (C) only indicating that the concentration of the novel coronavirus antigen is zero or below the detection limit of the test.

#### Invalid Result

No visible colored band appears at control line after performing the test. The directions may have not been followed correctly or the test kit may have deteriorated. It is recommended to re-sample and re-test.



## QUALITY CONTROL

### 1. Internal Quality Control

A procedural control is included in the test. A colored band appearing in the control region (C) is considered an internal procedural control. If the test procedure is performed properly and the test reagents of the control line are working, then the control line will always appear and the test result is considered valid.

### 2. External Quality Control

Good laboratory practice recommends the use of the control swabs to ensure that the test reagent is working and the test is properly performed.

**It is recommended that the control swabs are performed at any of the following scenarios:**

- (1) New users run the control swabs prior to testing patient specimens;
- (2) When receiving a new test reagent lot;
- (3) When receiving a new test shipment;
- (4) Users should follow the appropriate federal state, and local guidelines concerning the frequency of assaying external quality control swabs.

## LIMITATIONS OF PROCEDURE

1. This reagent is designed to detect novel coronavirus antigen N protein in human nasal swab specimen.
2. The specimen collection process will affect the accuracy of the test, such as improper specimen collection, improper specimen storage, or repeated freezing and thawing of the specimen etc.
3. This reagent is a qualitative assay. It is not designed to determine the quantitative concentration of novel coronavirus antigen. If you need to test the quantitative concentration, please refer to the instructions of relevant professional instruments or reagents for quantitative assay.
4. The test results of this reagent are for clinical reference only and should not be used as the sole basis for clinical diagnosis and treatment. The clinical management of patients should be comprehensively considered based on their symptoms / signs, medical history, other laboratory examinations and treatment response.
5. Limited by the method of antigen test reagents, it is recommended to use nucleic acid detection or virus culture identification methods for review and confirmation for negative test results.
6. Positive test results do not rule out co-infections with other pathogens. A negative result of this reagent can be caused by improper specimen collection, improper specimen transfer or handling.

## PERFORMANCE CHARACTERISTICS

### A. Clinical Performance

#### 1. Clinical Trial of Wondfo 2019-nCoV Antigen Test (Symptomatic)

The clinical performance of the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) was evaluated by testing nasal swab specimens collected from symptomatic subjects who were suspected of COVID-19. A total of 793 clinical specimens were tested using the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) and were compared to the results of oropharyngeal swab/Nasopharyngeal swab tested with molecular (RT-PCR) assay.

Clinical Performance of Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) vs RT-PCR Assay.

Wondfo 2019-nCoV Antigen Test (Nasal swab)	RT-PCR Test		Total
	Positive	Negative	
Positive	191	5	196
Negative	14	583	597
Total	205	588	793

Sensitivity: 93.17%(95%CI:88.86-95.89%)  
Specificity: 99.15%(95%CI:98.03-99.64%)  
Overall Percent Agreement: 97.60%(95%CI:96.29-98.46%)

Stratification of the positive specimens post onset of symptoms or suspected exposure between 0-3 days has a sensitivity of 91.78% (n=73) and 4-7 days has a sensitivity of 92.31% (n=91). Positive agreement of the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is higher with specimens of Ct values ≤30 with a sensitivity of 97.04%.

#### 2. Clinical Trial of Wondfo 2019-nCoV Antigen Test (Asymptomatic)

The clinical performance of the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) was evaluated by testing nasal swab specimens collected from asymptomatic subjects. A total of 509 clinical specimens were tested using the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) and were compared to the results of oropharyngeal swab tested with an Emergency Use Authorized molecular (RT-PCR) assay.

A total of 78 specimens with PCR positive results are stratified based on the Ct value of the target gene amplification as presented in the table below:

The specificity (n=431) was 100% with 95% CI [99.12%-100%].

Clinical Performance of Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) vs FDA EUA RT-PCR Assay.

	Total (n)	Sensitivity (%)	95% CI
CT≤25	29	96.55	82.82-99.39
CT≤30	63	92.06	82.73-96.56
CT≤33	72	83.33	73.09-90.20
RT-PCR positive	78	78.21	67.84-85.92

### B. Limit of Detection

The Limit of Detection (LoD) of the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) was determined using limiting dilutions of UV-inactivated novel coronavirus (ZeptoMetrix 0810587UV). The ZeptoMetrix material is a preparation of SARS-Related Coronavirus 2, isolate USA-WA1/2020. The material was supplied frozen at a concentration of 4.57 x10<sup>6</sup> TCID<sub>50</sub>/mL.

The Limit of Detection (LoD) of the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) was confirmed as 1.5×10<sup>2</sup> TCID<sub>50</sub>/mL.

### C. Analytical Reactivity/Inclusivity

The analytical reactivity of the novel coronavirus in the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) were evaluated with currently available novel coronavirus strains (see table below).

Virus Strain/WHO Label	Lineage#	Isolate	Concentration
UV-Inactivated SARS-CoV-2 virus Culture Fluid	/	USA-WA1/2020	1.5×10 <sup>6</sup> TCID <sub>50</sub> /mL
Alpha	B.1.1.7	IQTC34925	4.0×10 <sup>5</sup> PFU/mL
Beta	B.1.351	IQTC22239	4.0×10 <sup>5</sup> PFU/mL
Delta	B.1.617.2	GDPC2	4.0×10 <sup>5</sup> PFU/mL

### D. Interference

The test result of Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is not influenced by the presence of following substances/conditions:

No.	Types of Specimen	Interfering Substances	Concentration	No.	Types of Specimen	Interfering Substances	Concentration
1	Nasal gel	Naso GEL (NeilMed)	5% v/v	17	Anti-viral drugs	Ribavirin	10.0 µg/mL
2	Nasal sprays	Oxymetazoline hydrochloride	15% v/v	18		Oseltamivir	0.5 mg/mL
3	or drops	CVS Nasal Spray (Cromolyn)	15% v/v	19		Zanamivir	10.0 µg/mL
4		CVS Nasal Drops (Phenylephrine)	15% v/v	20		Ritonavir	2.0 mg/mL
5		Budesonide	15% v/v	21		Loprenavir	0.5 mg/mL
6	Nasal corticosteroids	Beclomethasone	15% v/v	22		Interferon alpha	5% v/v
7		Fluticasone Propionate	5% v/v	23		Palamivir	0.15 mg/mL
8	Throat churgs	Chloraseptic (Menthol/Benzocaine)	1.5 mg/mL	24	Allergy medication	Abidor	1.0 mg/mL
9		Sore Throat Phenol Spray	15% v/v	25		Diphenhydramine	10.0 µg/mL
10		Tobramycin	24.0 µg/mL	26	Anti-inflammatory medication	Chlorpheniramine	10.0 µg/mL
11		Levofloxacin	0.1 mg/mL	27		Acetylsalicylic acid	1.0 mg/mL
12	Antibiotic and Antibacterial	Azithromycin	1.0 mg/mL	28	Others	Ibuprofen	1.0 mg/mL
13		Ceftriaxone	0.8 mg/mL	29		Mucin	0.5% v/v
14		Meropenem	0.5 mg/mL	30		Whole Blood	4% v/v
15		Mupirocin	10.0 mg/mL	31		Biotin	0.1 mg/mL
16		Amoxicillin	0.25 mg/mL	32		HAMA	0.5% v/v

### E. Cross-reactivity and Microbial Interference

No cross-reactivity or interference was observed with the following microorganisms when tested at the concentrations presented in the table below.

No.	Types of Specimen	Cross Reaction Substance	Concentration
1	Other high priority pathogens from the same virus family	Human coronavirus 229E	4.68×10 <sup>4</sup> TCID <sub>50</sub> /mL
2		Human coronavirus OC43	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
3		Human coronavirus NL63	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
4		MERS-coronavirus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
5		Human coronavirus HKU1	0.25 mg/mL
6		SARS-coronavirus	0.25 mg/mL
7		Adenovirus-1	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
8		Adenovirus-2	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
9		Adenovirus-3	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
10		Adenovirus-4	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
11		Adenovirus-5	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
12	Adenovirus-7	1×10 <sup>5</sup> TCID <sub>50</sub> /mL	
13	Adenovirus-55	1×10 <sup>5</sup> TCID <sub>50</sub> /mL	
14	Human Metapneumovirus (hMPV-3)	1×10 <sup>5</sup> TCID <sub>50</sub> /mL	
15	Other high priority organisms	Parainfluenza virus 1	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
16		Parainfluenza virus 2	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
17		Parainfluenza virus 3	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
18		Parainfluenza virus 4	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
19		Epstein-Barr virus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
20		Measles virus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
21		Human Cytomegalovirus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
22		Rotavirus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
23		Norovirus	0.1 mg/mL
24		Mumps virus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
25		Varicella-zoster virus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
26	Influenza A H1N1	1×10 <sup>5</sup> TCID <sub>50</sub> /mL	
27	Influenza A H3N2	1×10 <sup>5</sup> TCID <sub>50</sub> /mL	

No.	Types of Specimen	Cross Reaction Substance	Concentration
28	Other high priority organisms	Influenza B Victoria	4.68×10 <sup>4</sup> TCID <sub>50</sub> /mL
29		Influenza B Yamagata	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
30		Enterovirus Group A (CA16)	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
31		Enterovirus Group B (CB1)	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
32		Enterovirus Group C (CA24)	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
33		Enterovirus Group D (EV68)	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
34		Respiratory syncytial virus (RSV) type A	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
35		Respiratory syncytial virus (RSV) type B	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
36		Rhinovirus-A	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
37		Rhinovirus-B	1×10 <sup>5</sup> TCID <sub>50</sub> /mL
38		Chlamydia pneumoniae	1×10 <sup>5</sup> IFU/mL
39		Haemophilus influenzae	1×10 <sup>5</sup> CFU/mL
40		Legionella pneumophila	1×10 <sup>5</sup> CFU/mL
41		Mycobacterium tuberculosis	1.0 mg/mL
42		Streptococcus pneumoniae	1×10 <sup>5</sup> CFU/mL
43		Streptococcus pyogenes	1×10 <sup>5</sup> CFU/mL
44		Bordetella pertussis	1×10 <sup>5</sup> CFU/mL
45		Mycoplasma pneumoniae	1×10 <sup>5</sup> CFU/mL
46		Pneumocystis jirovecii (PJP)	N/A
47		Pooled human nasal wash	N/A

\*Pneumocystis jirovecii (PJP) have not been tested. There is no significant similarities in Amino Acid between the nucleocapsid phosphoprotein of SARS-CoV-2 and the PJP.

### F. Hook Effect

Within the titer range of the following novel coronavirus culture, there is no hook effect in the test results of this product.

Virus Strain/WHO Label	Lineage#	Isolate	Concentration
UV-Inactivated SARS-CoV-2 virus Culture Fluid	/	USA-WA1/2020	4.57×10 <sup>6</sup> TCID <sub>50</sub> /mL
Alpha	B.1.1.7	IQTC34925	8.6×10 <sup>5</sup> PFU/mL
Beta	B.1.351	IQTC22239	2.0×10 <sup>5</sup> PFU/mL
Delta	B.1.617.2	GDPC2	2.0×10 <sup>5</sup> PFU/mL

### G. Precision

Repeatable & Reproducibility of Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) was established using in-house reference panels containing negative specimens and a range of positive specimens. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.

## BIBLIOGRAPHY

- [1] World Health Organization. WHO 2020, Statement on the first meeting of the International Health Regulations (2005) Emergency Committee regarding the outbreak of novel coronavirus (2019-nCoV). (Available from: [https://www.who.int/news/item/23-01-2020-statement-on-the-meeting-of-the-international-health-regulations-\(2005\)-emergency-committee-regarding-the-outbreak-of-novel-coronavirus-\(2019-ncov\)](https://www.who.int/news/item/23-01-2020-statement-on-the-meeting-of-the-international-health-regulations-(2005)-emergency-committee-regarding-the-outbreak-of-novel-coronavirus-(2019-ncov)))
- [2] Leung N H L, Chu D K W, Shiu E Y C, et al. Respiratory virus shedding in exhaled breath and efficacy of face masks[J]. Nature medicine, 2020, 26(5): 676-680.
- [3] Han Q, Lin Q, Ni Z, et al. Uncertainties about the transmission routes of 2019 novel coronavirus[J]. Influenza and Other Respiratory Viruses, 2020, 14(4): 470-471.
- [4] Johansson M A, Quandelacy T M, Kada S, et al. SARS-CoV-2 transmission from people without COVID-19 symptoms[J]. JAMA network open, 2021, 4(1): e2035057-e2035057.
- [5] Ciotti M, Ciccozzi M, Terrinoni A, et al. The COVID-19 pandemic[J]. Critical reviews in clinical laboratory sciences, 2020, 57(6): 365-388.
- [6] World Health Organization. WHO 2020, Coronavirus. (Available from: [https://www.who.int/health-topics/coronavirus/coronavirus#tab=tab\\_1](https://www.who.int/health-topics/coronavirus/coronavirus#tab=tab_1))

## INDEX OF SYMBOL

IVD In Vitro Diagnostic Use	See Instruction for Use	Expiry Date	Tests per Kit
Manufacturing Date	Keep Dry	LOT Batch Number	REF Catalog number
Keep away from Sunlight	Manufacturer	Do not reuse	Store between 2-30°C
Positive Control	Negative Control	CONTROL USE ONLY DO NOT USE for specimen collection	

**Guangzhou Wondfo Biotech Co., Ltd.**  
No. 8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R.China  
Tel: +86-20-32053962 400-888-5268(Toll free)  
Fax: +86-20-32296063 E-mail: [global@wondfo.com.cn](mailto:global@wondfo.com.cn)  
Website: [www.en.wondfo.com.cn](http://www.en.wondfo.com.cn)