WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Abbott RealTime HCV
WHO reference number: PQDx 0450-027-00

Abbott RealTime HCV with product codes **4J86-90, 4J86-80, 4J86-70, and 04J86-10¹**, manufactured by **Abbott Molecular Inc, CE-Marked regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 10 December 2019.

Summary of WHO prequalification assessment for the Abbott RealTime HCV

	Date	Outcome
Prequalification listing	10 December 2019	listed
Dossier review	N/A	N/A
Site inspection(s) of quality management system	28 October 2024	MR
Product performance evaluation	Quarter (Q) 3 and Q4 of 2019	MR

MR: Meet Requirements

N/A: Not Applicable

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarised in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	The manufacturer has submitted a change request to extend the claim of the Abbott RealTime HCV assay (4J86) to include quantification of hepatitis C viral RNA (HCV RNA) in whole blood spotted on cards as dried blood spots (DBS) (i.e. obtained via venipuncture or capillary blood) from HCV-infected individuals.	16 September 2020.
3.0	Updated Abbott's European Authorized Representative (EC Rep) legal entity name from Abbott GmbH & Co. KG to Abbott GmbH. Labelling	1 November 2021.

¹ product code 04J86-10 was added to allow for the use of dried blood spot (DBS) specimens in addition to plasma specimens.

	changes to comply with the labelling requirements for products registered under IVDR.	
4.0	The addition of the DBS list number (04J86-10) was omitted in version 3 of the public report.	9 August 2024.
5.0	Closure of IFU-related post-prequalification commitment, which was included in the product dossier assessment section of this public report.	16 September 2024.
6.0	Administrative changes to the labelling to meet the requirements of the new IVDR 2017/746/EU.	27 October 2025

Intended use:

According to the claim of intended use from Abbott Molecular Inc, *“Abbott RealTime HCV assay is an in vitro reverse transcription polymerase chain reaction (RT-PCR) assay for the quantitation of hepatitis C viral ribonucleic acid (HCV RNA) in human serum, plasma and whole blood spotted on cards as dried blood spots (DBS) (ie, obtained via venipuncture or capillary blood) from HCV-infected individuals. The Abbott RealTime HCV assay is intended for use as an aid in the management of HCV-infected patients undergoing antiviral therapy. The Abbott RealTime HCV assay is not for screening blood, plasma, serum or tissue donors for HCV, or to be used as a diagnostic test to confirm the presence of HCV infection.*

The intended users for the Abbott RealTime HIV-1 assay are laboratory and healthcare professionals.”

Assay description:

According to the claim of assay description from Abbott Molecular Inc, *“Abbott RealTime HCV assay consists of 3 reagent kits:*

- *Abbott RealTime HCV Amplification Reagent Kit*
- *Abbott RealTime HCV Control Kit*
- *Abbott RealTime HCV Calibrator Kit*

The Abbott RealTime HCV assay uses RT-PCR to generate amplified product from the RNA genome of HCV in clinical specimens. An RNA sequence that is unrelated to the HCV target sequence is introduced into each specimen at the beginning of sample preparation. This unrelated RNA sequence is simultaneously amplified by RT-PCR and serves as an internal control (IC) to demonstrate that the process has proceeded correctly for each sample. The amount of HCV target sequence that is present at each amplification cycle is measured through the use of fluorescent-labeled oligonucleotide probes on the Abbott m2000rt instrument. The probes do not generate signal unless they are specifically bound to the amplified product. The amplification cycle at which fluorescent signal is detected by the Abbott m2000rt is proportional to the log of the HCV RNA concentration present in the original sample”.

Test kit contents:

Component	Number of tests and product code
Abbott RealTime HCV Amplification Reagent Kit	96 tests (product code 4J86-90)
Abbott RealTime HCV Internal Control	4 vials x 1.2mL (4J86Y)
Abbott RealTime HCV Amplification Reagent Pack	4 packs X 24 tests/pack (4J86)
Abbott RealTime HCV Control Kit	product code 4J86-80
Abbott RealTime HCV Negative Control	8 vials x 1.8mL/ vial (4J86Z)
Abbott RealTime HCV Low Positive Control	8 vials x 1.3ml/vial (4J86W)
Abbott RealTime HCV High Positive Control	8 vials x 1.3mL/vial (4J86X)
Abbott RealTime HCV Calibrator Kit	Product code 4J86-70
Abbott RealTime HCV Calibrator A	12 vials x 1.3 mL/vial (4J86A)
Abbott RealTime HCV Calibrator B	12 vials x 1.3 mL/vial (4J86B)
Abbott RealTime HCV Amplification kit for DBS	Product code 4J86-10
Abbott RealTime HCV Supplement for Dried Blood Spot (DBS)	Supplement

Items required but not provided:

Component	Product code and description
Sample preparation area	
Abbott m2000sp instrument	9K14-02
Abbott m1000sp instrument	4J72-01
Abbott m24sp instrument	3N06-01
Abbott RealTime HCV m2000 ROW System Combined Application CD-ROM	1L69-07 or higher
Abbott m2000sp software with version 6.0 or higher	50-148400 or higher
5 mL Reaction Vessels	General Lab material
Abbott mSample Preparation System	04J70-24
Calibrated precision pipettes	20 µL to 1 000 µL
Aerosol barrier pipette tips	20 µL to 1 000µL
13 to 16 mm sample tubes	13 to 16 mm sample tubes

200 µL disposable tips	04J71-17
1000 µL disposable tips	04J71-17
Vortex mixer	General lab equipment
Abbott Optical Adhesive Covers	04J71-75
Abbott Adhesive Cover Applicator	9K32-01
Abbott Splash-Free Support Base	09K31-01
200 mL Reagent Vessels	4J71-60
Master Mix Tube	04J71-80
Abbott 96-Deep-Well Plate	04J71-30
Abbott 96-Well Optical Reaction Plate	04J71-70
Centrifuge capable of 2000g	General lab equipment
1.4 mL Micro Vial 15 mm Caps (optional)	3N20-01
Sample preparation area for DBS	
Abbott mSample Preparation System DBS Buffer Kit	09N02-001
m2000 System 13mm DBS-PoS Set	09N03-001
Amplification area	
Abbott m2000rt	9K15-01
Abbott RealTime HCV m2000 ROW System Combined Application CD-ROM	1L69-07 or higher
Abbott m2000rt Optical Calibration Kit	List No. 4J71-93
Amplification area for DBS	
Abbott RealTime HCV m2000 ROW System Combined Application CD-ROM	01L69-040 or higher
Other Materials	
Biological safety cabinet approved for working with infectious materials	General lab equipment
Sealable plastic bags	General lab material
RNase-free water (Eppendorf or equivalent) †	General lab material
1.7 mL molecular biology grade microcentrifuge tubes (Dot Scientific, Inc. or equivalent) †	
Cotton Tip Applicators (Puritan or equivalent) †	General lab material
Other Materials for DBS	
Sealable bag for DBS card	General lab material
Desiccant pack for DBS card	General lab material
Ahlstrom-Munktel TFN paper card and/or Whatman 903	General lab material

†NOTE: These 3 items are used in the procedure for Monitoring the Laboratory for the Presence of Contamination. Refer to the QUALITY CONTROL PROCEDURES section of this package insert.

Storage:

- The Abbott RealTime HCV Amplification Reagent Pack and Internal Control vials(4J86-90) must be stored at – 10°C or colder when not in use. Care must be taken to separate the Abbott RealTime HCV Amplification Reagent Pack that is in use from direct contact with samples, calibrators, and controls.
- The Abbott RealTime HCV Negative and Positive Controls (4J86-80) must be stored at – 10°C or colder.
- The Abbott RealTime HCV Calibrator A and Calibrator B (4J86-70) must be stored at – 10°C or colder.

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18 months.

Warnings/limitations:

Refer to the current instructions for use.

Prioritization for prequalification

Based on the established eligibility criteria, Abbott RealTime HCV was given priority for WHO prequalification assessment.

Product dossier assessment

In accordance with the WHO procedure for abridged prequalification assessment, Abbott Molecular Inc. was not required to submit a product dossier for the Abbott RealTime HCV as per the *“Instructions for compilation of a product dossier”* (PQDx_018 version 3). Notwithstanding, certain aspects of the product dossier previously submitted for stringent regulatory review were reviewed by an assessor during the site inspection. Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for the Abbott RealTime HCV for prequalification.

Manufacturing site inspection

At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality. Routine inspections of the Manufacturing site will be conducted with copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current assessment performed at a manufacturing site for in vitro diagnostic products and summarises the assessment findings.

<https://extranet.who.int/pqweb/vitro-diagnostics/who-public-inspection-reports>

All published WHOPIRs are with the agreement of the manufacturer. Based on the site inspection and corrective action plan review, the quality management system for the Abbott RealTime HCV meets WHO prequalification requirements.

Product performance evaluation

The Abbott RealTime HCV assay is an in vitro reverse transcription polymerase chain reaction (RT-PCR) assay for the quantitation of hepatitis C viral ribonucleic acid (HCV RNA) in human serum and EDTA/ACD plasma from HCV-infected individuals. The Abbott RealTime HCV assay is intended for use as an aid in the management of HCV-infected patients undergoing antiviral therapy. A volume of 0.2 mL or 0.5 mL of specimen is used to perform the assay (however, at least 0.7 – 1.3 mL of specimen should be available for use with the *m1000* and *m2000sp* systems). This type of assay requires laboratory equipment and cannot be performed in laboratories with limited facilities.

The Abbott RealTime HCV assay was evaluated at the National Serology Reference Laboratory (NRL), Melbourne, Australia, on behalf of WHO in Q3 and Q4 2019. From this evaluation, we drew the following conclusions.

Analytical evaluation

The precision (intra- and inter-assay) of measurement was found to be acceptable. At 10^2 IU/mL, CV% were <8% and at 10^4 IU/mL, CV% were < 3%.

The limit of detection of the assay was estimated at 31 IU/mL (95% CI: 11-52 IU/mL), which is consistent with the claim (30 IU/mL with the 0.2 mL sample preparation procedure).

There was no cross-contamination when high positive and negative specimens were ran together.

All 6 HCV genotypes included in the 4th HCV RNA Genotype Panel for Nucleic Acid Amplification (NIBSC reference: 14/209) and in the NRL HCV Mixed Genotype Panel were detected by the assay.

Clinical evaluation

In this limited performance evaluation on a panel of 197 ACD plasma specimens, using a sample volume of 0.2 mL, we found a bias of $-0.05 \log_{10} \text{ IU/mL}$ (limits of agreement: -0.42 to $0.32 \log_{10} \text{ IU/mL}$) compared to the reference results using the cobas HCV quantitative nucleic acid test with cobas 6800 System (Roche Diagnostics GmbH). The sensitivity (N=99) was 100% (95% CI: 96.3% – 100%) and specificity (N=98) was 100% (95% CI: 96.2% - 100%) compared to the reference results.

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(en) For *In Vitro* Diagnostic Use. The Abbott RealTime HCV assay is an *in vitro* reverse transcription-polymerase chain reaction (RT-PCR) assay for the quantitation of hepatitis C viral ribonucleic acid (HCV RNA) in human serum and plasma from HCV-infected individuals. The Abbott RealTime HCV assay is intended for use as an aid in the management of HCV-infected patients undergoing antiviral therapy. The Abbott RealTime HCV assay is not for screening blood, plasma, serum or tissue donors for HCV, or to be used as a diagnostic test to confirm the presence of HCV infection.

Contents:

1. **[INTERNAL CONTROL]** Abbott RealTime HCV Internal Control (4 vials, 1.2 mL per vial). <0.01% noninfectious Armored RNA® with internal control sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HbsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin® 300 and 0.15% ProClin 950.
2. **[AMPLIFICATION REAGENT PACK]** Abbott RealTime HCV Amplification Reagent Pack (4 packs, 24 tests/pack). Each Reagent Pack contains:
 - 1 Bottle (0.141 mL) ThermoStable rTth Polymerase Enzyme (2.9 to 3.5 Units/µL) in buffered solution.
 - 1 Bottle (1.10 mL) HCV Oligonucleotide Reagent. <0.1% synthetic oligonucleotides (4 primers and 2 probes), and <0.3% dNTPs in a buffered solution with a reference dye. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
 - 1 Bottle (0.40 mL) Activation Reagent. 30 mM manganese chloride solution. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.

ProClin is a registered trademark of Rohm and Haas.

Armored RNA is a registered trademark of Ambion.

Abbott RealTime is a trademark of Abbott.

(en) **CAUTION:** Handle human sourced materials as potentially infectious. Consult instructions for use. / (de) **ACHTUNG:** Humanmaterial gilt als potentiell infektiös und muss mit der entsprechenden Vorsicht gehandhabt werden. Gebrauchsanweisung beachten. / (fr) **ATTENTION :** Manipuler les produits d'origine humaine comme s'ils étaient potentiellement infectieux. Consulter les instructions d'utilisation. / (es) **ATENCIÓN:** maneje los productos de origen humano como potencialmente infecciosos. Consulte las instrucciones de uso. / (it) **ATTENZIONE:** Trattare i materiali di origine umana come potenzialmente infettivi. Consultare le istruzioni per l'uso. / (pt) **ATENÇÃO:** manusear os materiais de origem humana como potencialmente infecciosos. Consultar as instruções de utilização.

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




1.3.2 Abbott RealTime HCV CAL B vial label (4J86B)

2. Instructions for use²

² English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

HCV Calibrators

NOTE: Changes Highlighted

Key to Symbols Used	
REF	Reference Number
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
	Use by
CAL A	Calibrator A
CAL B	Calibrator B
	Upper limit of temperature
	Consult instructions for use
	Warning
	CAUTION
EC REP	Authorized representative in the European Community
	Manufacturer

Notice to User

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the Customer service section or Technical assistance section of these instructions.

Intended Use

The Abbott RealTime HCV Calibrators are for calibration of the Abbott RealTime HCV assay when used for the quantitative determination of Hepatitis C Virus (HCV) RNA in human serum and plasma from HCV infected individuals.

Intended User

The intended users for the Abbott RealTime HCV Calibrators are laboratory professionals.

Contents

- CAL A Abbott RealTime HCV Calibrator A (List No. 4J86A)** (12 vials, 1.3 mL per vial). Noninfectious Armored RNA® with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HIV RNA, HCV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin® 300 and 0.15% ProClin 950.
 - CAL B Abbott RealTime HCV Calibrator B (List No. 4J86B)** (12 vials, 1.3 mL per vial). Noninfectious Armored RNA with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HIV RNA, HCV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
- Calibrator concentrations are specified in each Abbott RealTime HCV Calibrator Kit Card.
 - The Abbott RealTime HCV Calibrator Kit must only be used with the Abbott RealTime HCV assay (List No. 4J86-90).

Standardization

Abbott manufactures internal reference standards for the Abbott RealTime HCV assay. These internal standards are referenced to the Second WHO International Standard for Hepatitis C Virus RNA (NIBSC Code 96/798)¹ at each concentration level. The Abbott RealTime HCV calibrators are manufactured against these internal standards.

Precautions

- IVD** In Vitro Diagnostic Medical Device
- For In Vitro Diagnostic Use Only**
- Do not use beyond expiration date.

CAUTION: This preparation contains human sourced and/or potentially infectious components. Components sourced from human blood have been tested and found to be nonreactive by FDA-licensed tests for antibody to HCV, antibody to HIV-1, antibody to HIV-2, and HBsAg. The material is also tested and found to be negative by FDA-licensed PCR methods for HIV-1 RNA and HCV RNA. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. These reagents and human specimens should be handled as if infectious using safe laboratory procedures, such as those outlined in Biosafety in Microbiological and Biomedical Laboratories,² OSHA Standards on Bloodborne Pathogens,³ CLSI Document M29-A4,⁴ and other appropriate biosafety practices.⁵ Therefore all human sourced materials should be considered infectious.

These precautions include, but are not limited to, the following:

- Wear gloves when handling specimens or reagents.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% sodium hypochlorite or other suitable disinfectant.²
- Decontaminate and dispose of all potentially infectious materials in accordance with local, state and federal regulations.⁵

Components of the Abbott RealTime HCV Calibrator Kit (List No. 4J86-70) contain the following components:

- 2-Methyl-2H-isothiazol-3-one
- Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one (EC no. 247-500-7) and 2-methyl-2H-isothiazol-3-one (EC no. 220-239-6)(3:1)
- Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one (EC no. 247-500-7) and 2-methyl-4-isothiazolin-3-one (EC no. 220-239-6)(3:1)

The following warnings apply:

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing mist / vapours / spray.
P280	Wear protective gloves / protective clothing / eye protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P501	Dispose of contents / container in accordance with local regulations.

Consult instructions for use

Store at $\leq -10^{\circ}\text{C}$

Shipping Conditions

Ship on dry ice.

BIBLIOGRAPHY

1. Saldanha J, Lelie N, Heath A. Establishment of the first international standard for nucleic acid amplification technology (NAT) assays for HCV RNA. *Vox Sang.* 1999;76:149-58.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009. [Also available online. Type> www.cdc.gov, search>BMBL5>look up sections III and IV.]
3. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
4. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
5. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva, Switzerland: World Health Organization; 2004.

Technical Assistance

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

SUMMARY OF SAFETY AND PERFORMANCE STATEMENT

A summary of safety and performance (SSP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSP location after the launch of European Database on Medical Devices. Search for device using UDI-DI provided on the outer packaging of the device.

Armored RNA is a registered trademark of Ambion.

ProClin is a registered trademark of Rohm and Haas.

Abbott RealTime is a trademark of Abbott.

The Abbott RealTime HCV Calibrators Kit is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

©2005, 2024 Abbott. All Rights Reserved.

February 2024

HCV Controls

NOTE: Changes Highlighted

Key to Symbols Used	
REF	Reference Number
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
CONTROL -	Negative Control
CONTROL L	Low Positive Control
CONTROL H	High Positive Control
 -10°C	Upper limit of temperature
	Warning
	Consult instructions for use
	CAUTION
EC REP	Authorized representative in the European Community
	Manufacturer

Notice to User

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the Customer service section or Technical assistance section of these instructions.

Intended Use

The Abbott RealTime HCV Controls are used to establish run validity of the Abbott RealTime HCV assay when used for the quantitation of Hepatitis C Virus (HCV) RNA in human serum and plasma from HCV infected individuals.

Intended User

The intended users for the Abbott RealTime HCV Controls are laboratory professionals.

Contents

- CONTROL - Abbott RealTime HCV Negative Control** (List No. 4J86Z) (8 vials, 1.8 mL per vial). Negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HIV RNA, HCV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin® 300 and 0.15% ProClin 950.
- CONTROL L Abbott RealTime HCV Low Positive Control** (List No. 4J86W) (8 vials, 1.3 mL per vial). Noninfectious Armored RNA® with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HIV RNA, HCV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
- CONTROL H Abbott RealTime HCV High Positive Control** (List No. 4J86X) (8 vials, 1.3 mL per vial). Noninfectious Armored RNA with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HIV RNA, HCV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.

- Control concentrations are specified in each Abbott RealTime HCV Control Kit Card.
- The Abbott RealTime HCV Control Kit must only be used with the Abbott RealTime HCV assay (List No. 4J86-90).

Precautions

- IVD In Vitro Diagnostic Medical Device**
- For In Vitro Diagnostic Use Only**
- Do not use beyond expiration date.

CAUTION: This preparation contains human sourced and/or potentially infectious components. Components sourced from human blood have been tested and found to be nonreactive by FDA-licensed tests for antibody to HCV, antibody to HIV-1, antibody to HIV-2, and HBsAg. The material is also tested and found to be negative by FDA-licensed PCR methods for HIV-1 RNA and HCV RNA. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. These reagents and human specimens should be handled as if infectious using safe laboratory procedures, such as those outlined in Biosafety in Microbiological and Biomedical Laboratories,¹ OSHA Standards on Bloodborne Pathogens,² CLSI Document M29-A4,³ and other appropriate biosafety practices.⁴ Therefore all human sourced materials should be considered infectious.

These precautions include, but are not limited to, the following:

- Wear gloves when handling specimens or reagents.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% sodium hypochlorite or other suitable disinfectant.¹
- Decontaminate and dispose of all potentially infectious materials in accordance with local, state and federal regulations.⁴

Components of the Abbott RealTime HCV Control Kit (List No. 4J86-80) contain the following components:

- 2-Methyl-2H-isothiazol-3-one
- Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one (EC no. 247-500-7) and 2-methyl-2H-isothiazol-3-one (EC no. 220-239-6)(3:1)
- Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one (EC no. 247-500-7) and 2-methyl-4-isothiazolin-3-one (EC no. 220-239-6)(3:1)

The following warnings apply:

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing mist / vapours / spray.
P280	Wear protective gloves / protective clothing / eye protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P501	Dispose of contents / container in accordance with local regulations.

Consult instructions for use

Store at $\leq -10^{\circ}\text{C}$

Shipping Conditions

Ship on dry ice.

BIBLIOGRAPHY

1. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
[Also available online. Type> www.cdc.gov, search>BMBL5>look up sections III and IV.]
2. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
3. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
4. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva, Switzerland: World Health Organization; 2004.

Technical Assistance

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

SUMMARY OF SAFETY AND PERFORMANCE STATEMENT

A summary of safety and performance (SSP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSP location after the launch of European Database on Medical Devices. Search for device using UDI-DI provided on the outer packaging of the device.

Armored RNA is a registered trademark of Ambion.

ProClin is a registered trademark of Rohm and Haas.

Abbott RealTime is a trademark of Abbott.

The Abbott RealTime HCV Controls Kit is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

© 2005, 2024 Abbott. All Rights Reserved.
February 2024

REF 4J86

REF 4J86

51-602124/R13

51-602124/R13

NOTE: Changes Highlighted

Key to Symbols Used	
REF	Reference Number
LOT	Lot Number
IVD	In Vitro Diagnostic Medical Device
	Use By
CONTROL -	Negative Control
CONTROL L	Low Positive Control
CONTROL H	High Positive Control
CAL A	Calibrator A
CAL B	Calibrator B
INTERNAL CONTROL	Internal Control
AMPLIFICATION REAGENT PACK	Amplification Reagent Pack
	Upper Limit of Temperature
	Consult instructions for use
	Caution
	Warning
	Manufacturer
EC REP	Authorized Representative in the European Community

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

NOTICE TO USER

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the Customer service section or Technical assistance section of these instructions.

CUSTOMER SERVICE

INTERNATIONAL: CALL YOUR ABBOTT REPRESENTATIVE

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

Abbott RealTime HCV

INTENDED USE

The Abbott RealTime HCV assay is an in vitro reverse transcription-polymerase chain reaction (RT-PCR) assay for the quantitation of hepatitis C viral ribonucleic acid (HCV RNA) in human serum and plasma from HCV-infected individuals. The Abbott RealTime HCV assay is intended for use as an aid in the management of HCV-infected patients undergoing antiviral therapy. The Abbott RealTime HCV assay is not for screening blood, plasma, serum or tissue donors for HCV, or to be used as a diagnostic test to confirm the presence of HCV infection.

INTENDED USER

The intended users for the Abbott RealTime HCV assay are laboratory professionals.

SUMMARY AND EXPLANATION OF THE TEST

HCV is a single-stranded RNA virus, with a genome of 9,500 nucleotides.¹ HCV has been identified as the major etiological agent for post-transfusion non-A, non-B hepatitis worldwide.^{2,3} HCV has been transmitted primarily through intravenous drug use and through blood products. Sensitive serological tests for HCV antibodies have greatly reduced the incidence of new infections from donated blood. About 85% of HCV-infected individuals develop chronic hepatitis, with up to 20% of chronically infected individuals developing cirrhosis. In patients with cirrhosis, the incidence of hepatocellular carcinoma is 1-4% per year.^{4,5} Quantitation of HCV RNA has been instrumental in understanding the effectiveness of antiviral response to interferon monotherapy, interferon plus ribavirin combination therapy, and peginterferon plus ribavirin combination therapy.⁶⁻¹⁰ Current guidelines for the management and treatment of HCV recommend quantitative testing for HCV RNA before the start of antiviral therapy, during therapy, and after the conclusion of treatment. The objective of treatment is a sustained virological response (SVR), defined as the absence of HCV RNA detectable by a sensitive test after the end of treatment. SVR is almost always preceded by an early virological response (EVR), defined as a 2-log or greater decrease in HCV viral load after 12 weeks of therapy.^{4,5,11}

HCV RNA in serum or plasma can be quantitated using nucleic acid amplification or signal amplification technologies.¹²⁻¹⁴ The Abbott RealTime HCV assay uses RT-PCR technology combined with homogeneous real time fluorescent detection for the quantitation of HCV RNA. The selection of a conserved region of the HCV genome provides for the detection of genotypes 1 through 6. The assay is standardized against the Second WHO International Standard for Hepatitis C Virus RNA (NIBSC Code 96/798)¹⁵ and results are reported in international units/mL (IU/mL).

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The Abbott RealTime HCV assay consists of 3 reagent kits:

- Abbott RealTime HCV Amplification Reagent Kit
- Abbott RealTime HCV Control Kit
- Abbott RealTime HCV Calibrator Kit

The Abbott RealTime HCV assay uses RT-PCR¹⁶ to generate amplified product from the RNA genome of HCV in clinical specimens. An RNA sequence that is unrelated to the HCV target sequence is introduced into each specimen at the beginning of sample preparation.

This unrelated RNA sequence is simultaneously amplified by RT-PCR and serves as an internal control (IC) to demonstrate that the process has proceeded correctly for each sample. The amount of HCV target sequence that is present at each amplification cycle is measured through the use of fluorescent-labeled oligonucleotide probes on the Abbott *m2000rt* instrument. The probes do not generate signal unless they are specifically bound to the amplified product. The amplification cycle at which fluorescent signal is detected by the Abbott *m2000rt* is proportional to the log of the HCV RNA concentration present in the original sample.

Sample Preparation

The purpose of sample preparation is to extract and concentrate the target RNA molecules, to make the target accessible for amplification, and to remove potential inhibitors of amplification from the extract.

The Abbott *mSample* Preparation System (4 × 24 Preps) uses magnetic particle technology to capture nucleic acids and washes the particles to remove unbound sample components. The bound nucleic acids are eluted and transferred to output tubes or a 96 deep-well plate. The nucleic acids are then ready for amplification. The IC is taken through the entire sample preparation procedure along with the calibrators, controls, and specimens.

The Abbott *m2000sp* System can be used to prepare samples for the Abbott RealTime HCV assay. The Abbott *m2000sp* provides automated sample eluate transfer and reaction assembly in the Abbott 96-Well Optical Reaction Plate. Alternatively, samples can be prepared manually using the Abbott *mSample* Preparation System followed by manual reaction assembly.

Reagent Preparation and Reaction Plate Assembly

The Abbott *m2000sp* combines the Abbott RealTime HCV amplification reagent components (HCV Oligonucleotide Reagent, Thermostable rTth Polymerase Enzyme, and Activation Reagent). The Abbott *m2000sp* dispenses the resulting master mix to the Abbott 96-Well Optical Reaction Plate along with aliquots of the nucleic acid samples prepared by the Abbott *m2000sp*. The plate is ready, after manual application of the optical seal, for transfer to the Abbott *m2000rt*.

Manual sample preparation method users manually combine the Abbott RealTime HCV amplification reagent components to create the amplification master mix and transfer aliquots of the master mix and sample eluates to the reaction plate. The plate is ready, after manual application of the optical seal and centrifugation, for transfer to the Abbott *m2000rt*.

Amplification

During the amplification reaction on the Abbott *m2000rt*, the target RNA is converted to cDNA by the reverse transcriptase activity of the thermostable rTth DNA polymerase. First, the HCV and IC reverse primers anneal to their respective targets and are extended during a prolonged incubation period. After a denaturation step, in which the temperature of the reaction is raised above the melting point of the double-stranded cDNA:RNA product, a second primer anneals to the cDNA strand and is extended by the DNA polymerase activity of the rTth enzyme to create a double-stranded DNA product.

During each round of thermal cycling, amplification products dissociate to single strands at high temperature, allowing primer annealing and extension as the temperature is lowered. Exponential amplification of the product is achieved through repeated cycling between high and low temperatures, resulting in a billion-fold or greater amplification of target sequences. Amplification of both targets (HCV and IC) takes place simultaneously in the same reaction.

The target sequence for the Abbott RealTime HCV assay is in the 5'UTR region of the HCV genome. This region is specific for HCV and is highly conserved.¹⁷ The primers are designed to hybridize to the 5'UTR region with the fewest possible mismatches among HCV genotypes 1 through 6. The IC target sequence is derived from the hydroxypyruvate reductase gene from the pumpkin plant, *Cucurbita pepo*, and is delivered in an Armored RNA[®] particle that has been diluted in negative human plasma.

Detection

During the read cycles of amplification on the Abbott *m2000rt*, the temperature is lowered further to allow fluorescent detection of amplification products as the HCV and IC probes anneal to their targets (real-time fluorescence detection). The HCV and IC probes are single-stranded DNA oligonucleotides consisting of a probe sequence with a fluorescent moiety that is covalently linked to the 5' end of the probe and a quenching moiety that is covalently linked to the 3' end of the probe. In the absence of the HCV or IC target sequences, probe fluorescence is quenched. In the presence of HCV or IC target sequences, probe

hybridization to complementary sequences separates the fluorophore and the quencher and allows fluorescent emission and detection. The HCV and IC probes are each labeled with a different fluorophore, thus allowing for simultaneous detection of both amplified products at each cycle. The amplification cycle at which fluorescent signal is detected by the Abbott *m2000rt* is proportional to the log of the HCV RNA concentration present in the original sample.

PREVENTION OF NUCLEIC ACID CONTAMINATION

The possibility of nucleic acid contamination is minimized because:

- Reverse transcription, PCR amplification, and oligonucleotide hybridization occur in a sealed Abbott 96-Well Optical Reaction Plate.
- Detection is carried out automatically without the need to open the Abbott 96-Well Optical Reaction Plate.
- Pipettes with aerosol barrier tips or disposable transfer pipettes are used for all pipetting. The disposable pipettes or pipette tips are discarded after use.
- Separate dedicated areas are used to perform the Abbott RealTime HCV assay. Refer to the **SPECIAL PRECAUTIONS** section of this package insert.

REAGENTS

Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90)

1. **INTERNAL CONTROL** Abbott RealTime HCV Internal Control (List No. 4J86Y) (4 vials, 1.2 mL per vial)
< 0.01% noninfectious armored RNA with internal control sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin[®] 300 and 0.15% ProClin 950.
2. **AMPLIFICATION REAGENT PACK** Abbott RealTime HCV Amplification Reagent Pack (List No. 4J86) (4 packs, 24 tests/pack)
 - 1 bottle (0.141 mL). Thermostable rTth Polymerase Enzyme (2.9 to 3.5 units/μL) in buffered solution.
 - 1 bottle (1.10 mL). HCV Oligonucleotide Reagent. < 0.1% synthetic oligonucleotides (4 primers and 2 probes) and < 0.3% dNTPs in a buffered solution with a reference dye. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
 - 1 bottle (0.40 mL). Activation Reagent. 30 mM manganese chloride solution. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.

Abbott RealTime HCV Control Kit (List No. 4J86-80)

1. **CONTROL -** Abbott RealTime HCV Negative Control (List No. 4J86Z) (8 vials, 1.8 mL per vial)
Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
2. **CONTROL L** Abbott RealTime HCV Low Positive Control (List No. 4J86W) (8 vials, 1.3 mL per vial)
Noninfectious armored RNA with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
3. **CONTROL H** Abbott RealTime HCV High Positive Control (List No. 4J86X) (8 vials, 1.3 mL per vial)
Noninfectious armored RNA with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.

Abbott RealTime HCV Calibrator Kit (List No. 4J86-70)

1. **CAL A** Abbott RealTime HCV Calibrator A (List No. 4J86A) (12 vials, 1.3 mL per vial)
Noninfectious armored RNA with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.
2. **CAL B** Abbott RealTime HCV Calibrator B (List No. 4J86B) (12 vials, 1.3 mL per vial)
Noninfectious armored RNA with HCV sequences in negative human plasma. Negative human plasma tested and found to be nonreactive for HBsAg, HBV DNA, HCV RNA, HIV RNA, anti-HIV-1/HIV-2, and anti-HCV. Preservatives: 0.1% ProClin 300 and 0.15% ProClin 950.

WARNINGS AND PRECAUTIONS

IVD

For In Vitro Diagnostic Use

The Abbott RealTime HCV assay is not for screening blood, plasma, serum or tissue donors for HCV, or to be used as a diagnostic test to confirm the presence of HCV infection.

Safety Precautions

Refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals, Hazard Section, or the Manual Sample Preparation for Abbott RealTime RNA Assays, Handling Precaution section, for instructions on safety precautions.

CAUTION: This preparation contains human sourced and/or potentially infectious components. Components sourced from human blood have been tested and found to be nonreactive by FDA-licensed tests for antibody to HCV, antibody to HIV-1, antibody to HIV-2, and HBsAg. The material is also tested and found to be negative by FDA-licensed PCR methods for HIV-1 RNA and HCV RNA. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. These reagents and human specimens should be handled as if infectious using safe laboratory procedures, such as those outlined in Biosafety in Microbiological and Biomedical Laboratories,¹⁸ OSHA Standards on Bloodborne Pathogens,¹⁹ CLSI Document M29-A4,²⁰ and other appropriate biosafety practices.²¹ Therefore all human sourced materials should be considered infectious.

These precautions include, but are not limited to, the following:

- Wear gloves when handling specimens or reagents.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% sodium hypochlorite or other suitable disinfectant.¹⁸
- Decontaminate and dispose of all potentially infectious materials in accordance with local, state, and federal regulations.²¹

Components of the Abbott RealTime HCV Calibrator Kit (List No. 4J86-70), the Abbott RealTime HCV Control Kit (List No. 4J86-80), and the Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90) contain the following components:

- 2-Methyl-2H-isothiazol-3-one
- Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one (EC no. 247-500-7) and 2-methyl-2H-isothiazol-3-one (EC no. 220-239-6)(3:1)
- Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one (EC no. 247-500-7) and 2-methyl-4-isothiazolin-3-one (EC no. 220-239-6)(3:1)

The following warnings apply:

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing mist/vapours/spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves/protective clothing/eye protection.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash before reuse.
P501	Dispose of contents / container in accordance with local regulations.

SPECIAL PRECAUTIONS

Handling Precautions

The Abbott RealTime HCV assay is only for use with human serum and plasma specimens that have been handled and stored in capped tubes as described in the **SPECIMEN COLLECTION, STORAGE, AND TRANSPORT TO THE TEST SITE** section.

During preparation of samples, compliance with good laboratory practices is essential to minimize the risk of cross-contamination between samples, and the inadvertent introduction of ribonucleases (RNases) into samples during and after the extraction procedure.

Proper aseptic technique should always be used when working with RNA.

Amplification technologies such as PCR are sensitive to accidental introduction of product from previous amplification reactions. Incorrect results could occur if either the clinical specimen or the RealTime reagents used in the amplification step become contaminated by accidental introduction of even a few molecules of amplification product. Measures to reduce the risk of contamination in the laboratory include physically separating the activities involved in performing PCR in compliance with good laboratory practices.

Work Areas

Use 3 dedicated areas within the laboratory for performing the Abbott RealTime HCV assay with manual sample preparation using the Abbott mSample Preparation System and the Abbott m2000rt:

- The **Reagent Preparation Area** is dedicated to combining the Abbott RealTime HCV amplification reagent components to create the amplification master mix and transferring aliquots of the master mix to the reaction plate. Laboratory coats, pipettes, pipette tips, and vortexers used in the **Reagent Preparation Area** must remain in this area and not be moved to either the **Sample Preparation Area** or the **Amplification Area**.
- The **Sample Preparation Area** is dedicated to processing samples (specimens, Abbott RealTime HCV Controls, and Calibrators), and to adding processed samples, controls, and calibrators to the Abbott 96-Well Optical Reaction Plate. **All reagents used in the Sample Preparation Area should remain in this dedicated area at all times. Laboratory coats, pipettes, pipette tips, and vortexers used in the Sample Preparation Area must remain in this area and not be moved to either the Reagent Preparation Area or the Amplification Area. Do not bring amplification product into the Sample Preparation Area.**
- The **Amplification Area** is dedicated to the amplification and detection of amplified product. Laboratory coats and equipment used in the **Amplification Area** must remain in this area and not be moved to either the **Reagent Preparation Area** or the **Sample Preparation Area**.

Only 2 dedicated areas, Sample Preparation Area and Amplification Area, are recommended when the Abbott m2000sp and Abbott m2000rt are used.

Components contained within a kit are intended to be used together. Do not mix components from different kit lots. For example, do not use the negative control from control kit lot X with the positive controls from control kit lot Y.

Do not use kits or reagents after the dates shown on kit labels.

Work areas and instrument platforms must be considered potential sources of contamination. Change gloves after contact with potential contaminants (specimens, eluates, and/or amplified product) before handling unopened reagents, negative control, positive controls, calibrators, or specimens. Refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals for instrument cleaning procedures.

If the Abbott *m2000sp* instrument run is aborted, dispose of all commodities and reagents according to the Abbott *m2000sp* Operations Manual. If the Abbott *m2000sp* master mix addition protocol is aborted, seal the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose according to the Abbott *m2000sp* Operations Manual, Hazards section, along with the gloves used to handle the plate.

If the Abbott *m2000rt* instrument run is interrupted or aborted, seal the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.

Decontaminate and dispose of all potentially biohazardous materials in accordance with local, state, and federal regulations.²¹ All materials should be handled in a manner that minimizes the chance of potential contamination of the work area.

NOTE: Autoclaving the sealed Reaction Plate will not degrade the amplified product and may contribute to the release of the amplified product by opening the sealed plate. The laboratory area can become contaminated with amplified product if the waste materials are not carefully handled and contained before and after processing.

Aerosol Containment

To reduce the risk of nucleic acid contamination due to aerosols formed during manual pipetting, aerosol barrier pipette tips must be used for all manual pipetting. The pipette tips must be used only 1 time. Clean and disinfect spills of specimens and reagents as stated in the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals.

Contamination and Inhibition

The following precautions should be observed to minimize the risks of RNase contamination, cross-contamination between samples, and inhibition:

- Wear appropriate personal protective equipment at all times.
- Use powder-free gloves.
- Change gloves after having contact with potential contaminants (specimens, eluates, and/or amplified product).
- To reduce the risk of nucleic acid contamination due to aerosols formed during pipetting, pipettes with aerosol barrier tips must be used for all pipetting. The length of the tip should be sufficient to prevent contamination of the pipette barrel. While pipetting, care should be taken to avoid touching the pipette barrel to the inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.
- Change aerosol barrier pipette tips between ALL manual liquid transfers.
- The Abbott *m* Sample Preparation System (4 × 24 Preps) reagents are single use only. Use new reagent troughs or vessels, reaction vessels, and newly opened reagents for every new Abbott RealTime HCV assay run. At the end of each run, discard all remaining reagents from the worktable as stated in the Abbott *m2000sp* Operations Manual and the Abbott *m*Sample Preparation System (4 × 24 Preps) product information sheet.

STORAGE INSTRUCTIONS

Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90).

 -10°C The Abbott RealTime HCV Amplification Reagent Pack and Internal Control vials must be stored at -10°C or colder when not in use. Care must be taken to separate the Abbott RealTime HCV Amplification Reagent Pack that is in use from direct contact with samples, calibrators, and controls.

Abbott RealTime HCV Control Kit (List No. 4J86-80).

 -10°C The Abbott RealTime HCV Negative and Positive Controls must be stored at -10°C or colder.

Abbott RealTime HCV Calibrator Kit (List No. 4J86-70).

 -10°C The Abbott RealTime HCV Calibrator A and Calibrator B must be stored at -10°C or colder.

SHIPPING CONDITIONS

- Abbott RealTime HCV Amplification Reagent Kit: Ship on dry ice.
- Abbott RealTime HCV Control Kit: Ship on dry ice.
- Abbott RealTime HCV Calibrator Kit: Ship on dry ice.

INDICATION OF INSTABILITY OR DETERIORATION OF REAGENTS

When a positive or negative control value is out of the expected range, it may indicate deterioration of the reagents. Associated test results are invalid and samples must be retested. Assay recalibration may be necessary.

INSTRUMENT PROCEDURE

The Abbott RealTime HCV application files must be installed on the Abbott *m2000sp* and Abbott *m2000rt* systems from the Abbott RealTime HCV *m2000* ROW System Combined Application CD-ROM prior to performing the assay. For detailed information on application file installation, refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals, Operating Instructions section.

SPECIMEN COLLECTION, STORAGE, AND TRANSPORT TO THE TEST SITE

Specimen Collection and Storage

Human serum and plasma (EDTA and ACD-A) specimens may be used with the Abbott RealTime HCV assay. Follow the manufacturer's instructions for processing collection tubes.

Freshly drawn specimens (whole blood) may be held at 2°C to 30°C for up to 6 hours prior to centrifugation.

Separate serum or plasma from cells by centrifugation.

After centrifugation, serum or plasma may be removed from cells. Serum or plasma specimens may be stored:

- at 15°C to 30°C for up to 24 hours
- at 2°C to 8°C for up to 3 days
- at -70°C or colder for longer term²²⁻²⁵

Multiple freeze/thaw cycles should be avoided. If frozen, thaw specimens at 15°C to 30°C or at 2°C to 8°C. Once thawed, if specimens are not being processed immediately, they can be stored at 2°C to 8°C for up to 6 hours.

NOTE: Serum and plasma specimens should not be frozen in non-gel blood collection tubes.

Specimen Transport

Ship specimens according to recommended storage temperature and time listed in the **Specimen Collection and Storage** section above. For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical, diagnostic, or biological specimens.

ABBOTT REALTIME HCV ASSAY PROCEDURE

This Abbott RealTime HCV package insert contains 2 assay protocols:

- Samples prepared for amplification using the manual sample preparation method, follow **ASSAY PROTOCOL I**.
- Samples prepared for amplification using the Abbott *m2000sp* instrument follow **ASSAY PROTOCOL II**.
The Abbott RealTime HCV assay provides 2 sample volume options (0.2 mL and 0.5 mL). (See assay protocol step 6 and **INTERPRETATION OF RESULTS** section.)

Materials Provided

- Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90)

Materials Required But Not Provided

- Abbott RealTime HCV Calibrator Kit (List No. 4J86-70)
- Abbott RealTime HCV Control Kit (List No. 4J86-80)

For manual sample preparation, refer to the Materials and Equipment Required Section of the Manual Sample Preparation for Abbott RealTime RNA Assays procedure (List No. 06L73).

For Abbott *m2000sp* Instrument Sample Preparation Area

- Abbott *m2000sp* instrument
- 5 mL Reaction Vessels
- Abbott *m*Sample Preparation System (4 × 24 Preps) (List No. 04J70-24)
- Calibrated precision pipettes capable of delivering 20 to 1000 µL
- Aerosol barrier pipette tips for 20 to 1000 µL pipettes
- 13 to 16 mm sample tubes
- 200 µL and 1000 µL disposable tips
- Vortex mixer
- Abbott Optical Adhesive Covers (List No. 04J71-75)
- Abbott Adhesive Cover Applicators
- Abbott Splash-Free Support Base (List No. 9K31-01)
- Master Mix Vial
- 200 mL reagent vessels
- Abbott 96-Deep-Well Plate (List No. 04J71-30)
- Abbott RealTime HCV *m2000* ROW System Combined Application CD-ROM (List No. 1L69)
- Abbott 96-Well Optical Reaction Plate (List No. 04J71-70)
- Centrifuge capable of 2,000g

For Manual Sample Preparation Reagent Preparation Area

- PCR cooler, either StrataCooler® 96 Benchtop Cooler or Eppendorf PCR-Cooler
- Abbott 96-Well Optical Reaction Plate (List No. 04J71-70)
- Calibrated precision pipettes capable of delivering 20 to 1000 µL
- Aerosol barrier pipette tips for 20 to 1000 µL pipettes
- Single-use RNase/DNase-free tube or container
- Vortex mixer
- Abbott Optical Adhesive Covers (List No. 04J71-75)

For Abbott *m2000rt* Instrument Amplification Area

- Abbott *m2000rt* instrument
- Abbott RealTime HCV *m2000* ROW System Combined Application CD-ROM (List No. 1L69)
- Abbott *m2000rt* Optical Calibration Kit (List No. 04J71-93)

- Abbott Adhesive Cover Applicators
- Abbott Splash-Free Support Base (List No. 9K31-01)

Other Materials

- Biological safety cabinet approved for working with infectious materials
- Sealable plastic bags
- RNase-free water (Eppendorf or equivalent)*
- 1.7 mL molecular biology grade microcentrifuge tubes (Dot Scientific, Inc. or equivalent)*
- Cotton tip applicators (Puritan or equivalent)*

*NOTE: These 3 items are used in the procedure for Monitoring the Laboratory for the Presence of Contamination. Refer to the QUALITY CONTROL PROCEDURES section of this package insert.

Procedural Precautions

Read the instructions in this package insert carefully before processing samples.

The Abbott RealTime HCV Calibrators, Internal Control, Negative Control, and Low and High Positive Control vials are intended for single-use only and should be discarded after use.

Use aerosol barrier pipette tips or disposable pipettes only 1 time when pipetting specimens, IC, or amplification reagents. To prevent contamination to the pipette barrel while pipetting, care should be taken to avoid touching the pipette barrel to the inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.

Monitoring procedures for the presence of amplification product can be found in the QUALITY CONTROL PROCEDURES section in this package insert.

To reduce the risk of nucleic acid contamination, clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% sodium hypochlorite or other suitable disinfectant.

The Abbott RealTime HCV Calibrators and Controls must be prepared in conjunction with the specimens to be tested. The use of the Abbott RealTime HCV Controls and Calibrators is integral to the performance of the Abbott RealTime HCV assay. Refer to the QUALITY CONTROL PROCEDURES section of this package insert for details.

ASSAY PROTOCOL I: MANUAL SAMPLE PREPARATION METHOD AND ABBOTT *m2000rt* INSTRUMENT

For a detailed description of how to perform an Abbott *m2000rt* protocol refer to the Abbott *m2000rt* Operations Manual, Operating Instructions section.

Laboratory personnel must be trained to operate the Abbott *m2000rt* instrument. The operator must have a thorough knowledge of the software applications and must follow good laboratory practices.

1. Thaw assay controls and IC at 15°C to 30°C or at 2°C to 8°C. Thaw calibrators at 15°C to 30°C or at 2°C to 8°C only if performing a calibration run; see QUALITY CONTROL PROCEDURES section of this package insert.
 - Once thawed, assay controls, IC, and calibrators can be stored at 2°C to 8°C for up to 24 hours before use.
 - Vortex each assay calibrator and each control 3 times for 2 to 3 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
2. Thaw amplification reagents at 15°C to 30°C or at 2°C to 8°C and store at 2°C to 8°C until required for the amplification master mix procedure.
 - Once thawed, the amplification reagents can be stored at 2°C to 8°C for up to 24 hours if not used immediately.

Sample Preparation Area

For the manual sample preparation method refer to Extraction Protocol Section of the Manual Sample Preparation for Abbott RealTime RNA Assays procedure.

3. Gently invert the Abbott *mSample* Preparation bottles to ensure a homogeneous solution. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.
4. Vortex each IC 3 times for 2 to 3 seconds before use.

5. Use a calibrated precision **PIPETTE DEDICATED FOR INTERNAL CONTROL USE ONLY** to add 500 µL of IC to each bottle of *mLysis* Buffer. Mix by gently inverting the container 5 to 10 times to minimize foaming.
6. A negative control, a low positive control, and a high positive control are included in each run.
 - The Abbott RealTime HCV assay provides 2 sample volume options for manual sample preparation (0.2 mL and 0.5 mL).
 - If frozen, thaw specimens at 15°C to 30°C or at 2°C to 8°C. Once thawed, specimens can be stored at 2°C to 8°C for up to 6 hours if not processed immediately.
 - Vortex each specimen 3 times for 2 to 3 seconds. **Specimens showing particulate matter or turbidity should be clarified by centrifugation at 2,000g for 5 minutes prior to testing.** Aliquot each specimen into clean tubes or vials if necessary. Avoid touching the inside of the cap when opening tubes.
 - **The assembly of the amplification master mix and sample eluates into the Abbott 96-Well Optical Reaction Plate (step 13) must be initiated within 1 hour after completion of Sample Preparation.**

Amplification Area

7. Switch on and initialize the Abbott *m2000rt* instrument.

NOTE: The Abbott *m2000rt* instrument requires 15 minutes to warm up.

8. Create the Abbott *m2000rt* test order. Refer to the Operating Instructions section of the Abbott *m2000rt* Operations Manual. From the Protocol screen, select the appropriate application file corresponding to the sample volume being tested.
 - Enter calibrator (needed if a calibration curve has not been stored on the Abbott *m2000rt*) and control lot specific values in the test order for accurate calibration and control evaluation. Lot-specific values are specified in each Abbott RealTime HCV Calibrator and Control Kit Card.

Reagent Preparation Area

All reagent preparation must take place in the dedicated Reagent Preparation Area. Refer to the Handling Precautions section of this package insert before preparing reagents.

NOTE: Change gloves before handling the amplification reagents.

9. Prepare the amplification master mix.
 - Each amplification reagent pack supports up to 24 reactions.
 - Prior to opening the amplification reagents, ensure that the contents of the vials are at the bottom by tapping the vials in an upright position on the bench to bring the liquid to the bottom of the vials.
 - Prepare the master mix by using a **PIPETTE DEDICATED FOR REAGENT USE ONLY** to add 271 µL of the HCV Activation Reagent (Reagent 1) and 949 µL of the HCV Oligonucleotide Reagent (Reagent 2) together in the Thermostable rTth DNA Polymerase Enzyme bottle (Reagent 3).
 - If performing 25 to 48 reactions, prepare a second amplification master mix with a second amplification reagent pack.
 - **The Abbott *m2000rt* protocol (step 16) must be initiated within 40 minutes of the addition of Activation Reagent into the first rTth Enzyme Reagent bottle (step 9).**
10. Pipette the contents of the master mix from the enzyme bottle(s) into a single-use RNase/DNase-free tube and vortex to mix.
11. Place an Abbott 96-Well Optical Reaction Plate in a PCR cooler, stored as indicated in the PCR cooler instruction manual. Using a **DEDICATED PIPETTE**, dispense 50 µL aliquots of the amplification master mix into the Abbott 96-Well Optical Reaction Plate. A calibrated repeat pipettor may be used. Visually verify that 50 µL has been dispensed into each well.
12. Transfer the Abbott 96-Well Optical Reaction Plate on the PCR cooler to the Sample Preparation Area.

Sample Preparation Area

13. In the Sample Preparation Area, transfer 50 µL of sample eluate to the Abbott 96-Well Optical Reaction Plate on the PCR cooler. **Use a separate pipette tip for each sample eluate transfer.** During the transfer of each sample, mix the reaction by pipetting up and down 3 to 5 times. Visually verify that 100 µL has been dispensed into each well.

- Seal the Abbott 96-Well Optical Reaction Plate according to the instructions in the Abbott *m2000rt* Operations Manual.
- Remove the Abbott 96-Well Optical Reaction Plate from the PCR cooler to the Abbott Splash-Free Support Base. Centrifuge the Abbott 96-Well Optical Reaction Plate in the Abbott Splash-Free Support Base at 500*g* for 5 minutes. Transfer to the Amplification Area.

NOTE: Do not transfer the PCR cooler to the Amplification Area.

Amplification Area

- Place the Abbott 96-Well Optical Reaction Plate in the Abbott *m2000rt* instrument. From the Protocol screen, select the appropriate application file corresponding to the sample volume being tested. Initiate the Abbott RealTime HCV protocol, as described in the Abbott *m2000rt* Operations Manual, Operating Instructions Section.

POST PROCESSING PROCEDURES

- Clean the PCR cooler as described in the PCR cooler instruction manual and return to the Reagent Preparation Area.
- Place the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.
- Clean the Splash-Free Support Base before next use, according to the Abbott *m2000rt* Operations Manual.
- For the manual sample preparation method users, refer to the Clean Up section of the Manual Sample Preparation for Abbott RealTime RNA Assays Procedure (List No. 06L73).

ASSAY PROTOCOL II: ABBOTT *m2000sp* INSTRUMENT AND ABBOTT *m2000rt* INSTRUMENT

For a detailed description of how to perform an Abbott *m2000sp* instrument and Abbott *m2000rt* instrument protocol refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals, Operating Instructions sections. The 96-sample capability requires Abbott *m2000sp* Software Version 4.0 or higher. Please follow Abbott *m2000sp* Operations Manual (List 09K20-04 or higher) and addendum or addenda.

Laboratory personnel must be trained to operate the Abbott *m2000sp* and Abbott *m2000rt* instruments. The operator must have a thorough knowledge of the applications run on the instruments and must follow good laboratory practices.

- Thaw assay controls and IC at 15°C to 30°C or at 2°C to 8°C. Thaw calibrators at 15°C to 30°C or at 2°C to 8°C only if performing a calibration run; see **QUALITY CONTROL PROCEDURES** section of this package insert.
 - Once thawed, assay controls, IC, and calibrators can be stored at 2°C to 8°C for up to 24 hours before use.
 - Vortex each assay calibrator and each control 3 times for 2 to 3 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
- Thaw amplification reagents at 15°C to 30°C or at 2°C to 8°C and store at 2°C to 8°C until required for the amplification master mix procedure.
 - Once thawed the amplification reagents can be stored at 2°C to 8°C for up to 24 hours if not used immediately.

NOTE: Use 1 bottle of *m*Lysis Buffer, 1 vial of IC, and 1 Abbott RealTime HCV Amplification Reagent Kit to support up to 24 reactions. Use a second set of reagents to support 25 to 48 reactions, a third set of reagents to support 49 to 72 reactions and a fourth set of reagents to support 73 to 96 reactions WITH THE EXCEPTION OF *m*MICROPARTICLES. USE ONLY TWO BOTTLES OF *m*MICROPARTICLES WHEN PROCESSING 25 TO 96 REACTIONS.

- Gently invert the Abbott *m*Sample Preparation bottles to ensure a homogeneous solution. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.
- Vortex each IC 3 times for 2 to 3 seconds before use.
- Use a calibrated precision **PIPETTE DEDICATED FOR INTERNAL CONTROL USE ONLY** to add 500 µL of IC to each bottle of *m*Lysis Buffer. Mix by gently inverting the container 5 to 10 times to minimize foaming.

- A total of 96 samples can be processed in each run.** A negative control, a low positive control, and a high positive control are included in each run, therefore allowing a maximum of 93 specimens to be processed per run.

- The Abbott RealTime HCV assay minimum sample volume and associated rack requirements on the Abbott *m2000sp* are:

		Abbott RealTime HCV Minimum Sample Volume Assay Application	
Rack	Tube Diameter ^a	0.2 mL	0.5 mL
13 mm	11.6 mm - 14.0 mm	0.7 mL	1.0 mL
16 mm	15.0 mm - 16.0 mm	1.0 mL	1.3 mL

^a Refers to sample tube outer diameter

- If frozen, thaw specimens at 15°C to 30°C or at 2°C to 8°C. Once thawed, specimens can be stored at 2°C to 8°C for up to 6 hours, if not processed immediately.
 - Vortex each specimen 3 times for 2 to 3 seconds before loading on the Abbott *m2000sp* worktable. **Specimens showing particulate matter or turbidity should be clarified by centrifugation at 2,000*g* for 5 minutes prior to testing.** Aliquot each specimen into clean tubes or vials if necessary. Refer to the Abbott *m2000sp* Operations Manual for tube sizes. Avoid touching the inside of the cap when opening tubes.
- Place the low and high positive controls, the negative control, the calibrators, if applicable, and the patient specimens into the Abbott *m2000sp* sample rack.
 - Place the 5 mL Reaction Vessels into the Abbott *m2000sp* 1 mL subsystem carrier.
 - Load the Abbott *m*Sample Preparation System reagents and the Abbott 96 Deep-Well Plate on the Abbott *m2000sp* worktable as described in the Abbott *m2000sp* Operations Manual, Operating Instructions.
 - From the Protocol screen, select the appropriate application file corresponding to the sample volume being tested. Initiate the sample extraction protocol as described in the Abbott *m2000sp* Operations Manual, Operating Instruction.
 - Enter calibrator (needed if a calibration curve has not been stored on the Abbott *m2000rt*) and control lot specific values in the **Sample Extraction: Worktable Setup, Calibrator and Control** fields. Lot-specific values are specified in each Abbott RealTime HCV Calibrator and Control Kit Card.
 - The Abbott *m2000sp* Master Mix Addition protocol (step 12) must be initiated within 1 hour after completion of Sample Preparation.**

NOTE: Change gloves before handling the amplification reagents.

- Load the amplification reagents and the master mix vial on the Abbott *m2000sp* worktable after sample preparation is completed.
 - Each amplification reagent pack supports up to 24 reactions.
 - Prior to opening the amplification reagents, ensure that the contents are at the bottom of the vials by tapping the vials in an upright position on the bench.
 - Remove and discard the amplification vial caps.
 - A second amplification reagent pack is required if performing 25 to 48 reactions.
 - A third amplification reagent pack is required if performing 49 to 72 reactions.
 - A fourth amplification reagent pack is required if performing 73 to 96 reactions.
- Select the appropriate deep well plate that matches the corresponding sample preparation extraction. Initiate the Abbott *m2000sp* Master Mix Addition protocol. Follow the instructions as described in the Abbott *m2000sp* Operations Manual, Operating Instructions section.

NOTE: The operator should not manually fill any empty/unfilled wells in the Abbott 96-Well Optical Reaction Plate.

- After sample extraction is complete, the Abbott *m2000sp* automatically fills any empty wells in the Abbott 96-Well Optical Reaction Plate when there are greater than 48 samples processed within a run. Plate fill is not performed for runs containing 48 samples or fewer.

- If prompted by the instrument, Reagent Carrier 2 should remain in place, minimally containing the reagent vessel for *mElution* Buffer (Reagent Carrier 2, location 6). If this reagent vessel has been unloaded, place a new reagent vessel with the *mElution* Buffer label into Reagent Carrier 2, location 6. System fluid will be added to the reagent vessel and used to fill empty wells. Once this process is complete, the system will continue with the master mix addition.

NOTE: System instructions for use of the automated plate-filling feature are found in the Abbott *m2000sp* Operations Manual (List No. 9K20-04 or higher), section 5, Operating Instructions, Sample Extraction—Closed Mode.

- The Abbott *m2000rt* protocol (step 16) must be started within 50 minutes of the initiation of the Master Mix Addition protocol (step 12).

NOTE: If the run is aborted for any reason subsequent to step 12, a new 96-well PCR plate must be used if the Abbott *m2000sp* Master Mix Addition Protocol (Step 12) will be repeated.

13. Switch on and initialize the Abbott *m2000rt* instrument in the amplification area.

NOTE: The Abbott *m2000rt* requires 15 minutes to warm-up.

NOTE: Remove gloves before returning to the sample preparation area.

14. Seal the Abbott 96-Well Optical Reaction Plate after the Abbott *m2000sp* instrument has completed addition of samples and master mix according to the Abbott *m2000sp* Operations Manual, Operating Instructions section.
15. Place the sealed optical reaction plate into the Splash-Free Support Base for transfer to the Abbott *m2000rt* instrument.
16. Place the Abbott 96-Well Optical Reaction Plate in the Abbott *m2000rt* instrument. From the Protocol screen, select the appropriate application file corresponding to the sample volume being tested. Initiate the Abbott RealTime HCV protocol, as described in the Abbott *m2000rt* Operations Manual, Operating Instructions section.

NOTE: If creating the Abbott *m2000rt* test order manually, enter sample IDs in the corresponding PCR tray locations according to the “Wells for Selected Plate” grid, found on the detail screen of the “PCR Plate Results” on the Abbott *m2000sp*. See Section 5 of the Abbott *m2000sp* Operations Manual.

POST PROCESSING PROCEDURES

1. Remove the Abbott 96 Deep-Well Plate from the worktable and dispose of according to the Abbott *m2000sp* Operations Manual.
2. Place the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose according to the Abbott *m2000rt* Operations Manual, along with the gloves used to handle the plate.
3. Clean the Splash-Free Support Base before next use, according to the Abbott *m2000rt* Operations Manual.

QUALITY CONTROL PROCEDURES

Abbott *m2000rt* Optical Calibration

Refer to the Calibration Procedures section in the Abbott *m2000rt* Operations Manual for a detailed description of how to perform an Abbott *m2000rt* Optical Calibration.

Optical calibration of the Abbott *m2000rt* instrument is required for the accurate measurement and discrimination of dye fluorescence during the Abbott RealTime HCV assay.

The following Abbott *m2000rt* Optical Calibration Plates are used to calibrate the Abbott *m2000rt* instrument for the Abbott RealTime HCV assay:

- FAM™ Plate (Carboxyfluorescein)
- ROX™ Plate (Carboxy-X-rhodamine)
- VIC® Plate (Proprietary dye)

Assay Calibration

For a detailed description of how to perform an assay calibration refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals, Operating Instructions sections.

A calibration curve is required to quantitate the HCV RNA concentration of specimens and controls. Two assay calibrators are run in replicates of 3 to generate a calibration curve (HCV concentration versus the threshold cycle [C_t] at which a reactive level of fluorescent signal is detected). The calibration curve slope and intercept are calculated and

stored on the instrument. The concentration of HCV RNA in a sample is calculated from the stored calibration curve. Results are automatically reported on the Abbott *m2000rt* workstation.

Follow the procedure for sample extraction, master mix addition, amplification and detection protocols as stated in the Abbott *m2000sp* Operations Manual and the Abbott *m2000rt* Operations Manual.

Once an Abbott RealTime HCV calibration is accepted and stored, it may be used for 6 months. During this time, all subsequent samples may be tested without further calibration unless:

- An Abbott RealTime HCV Amplification Reagent Kit with a new lot number is used.
- An Abbott *mSample* Preparation System (4 × 24 Preps) with a new lot number is used.
- An Abbott RealTime HCV application file for a different sample volume is used.

Detection of Inhibition

An IC threshold cycle (C_t) assay validity parameter is established during a calibration run.

A defined, consistent quantity of IC is introduced into each specimen, calibrator, and control at the beginning of sample preparation and measured on the Abbott *m2000rt* instrument to demonstrate proper specimen processing and assay validity. The IC is composed of an RNA sequence unrelated to the HCV target sequence.

The median amplification cycle at which the IC target sequence fluorescent signal is detected in calibration samples establishes an IC C_t validity range to be met by all subsequent processed specimens.

An error control flag is displayed when a specimen or control fails to meet this specification. Refer to the Abbott *m2000rt* Operations Manual for an explanation of the corrective actions for the error control flag. Specimens whose IC C_t value exceeds the established range must be retested starting with sample preparation.

Negative and Positive Controls

A negative control, a low positive control, and a high positive control are included in each test order to evaluate run validity.

The lot specific values for the low positive control and high positive control are specified on each Abbott RealTime HCV Control Kit Card and must be entered into the assay test order when a run is performed.

An error control flag is displayed when a control result is out of range. Refer to the Abbott *m2000rt* Operations Manual for an explanation of the corrective actions for the error control flag. If negative or positive controls are out of range, all of the specimens and controls from that run must be reprocessed, beginning with sample preparation.

The presence of HCV must not be detected in the negative control. HCV detected in the negative control is indicative of contamination by other samples or by amplified product introduced during sample preparation or during preparation of the Abbott 96-Well Optical Reaction Plate. To avoid contamination, clean the Abbott *m2000sp* instrument and the Abbott *m2000rt* instrument and repeat sample processing for controls and specimens following the **Procedural Precautions**. If negative controls are persistently reactive, contact your Abbott representative.

Monitoring the Laboratory for the Presence of Contamination

It is recommended that this test be done at least once a month to monitor laboratory surfaces and equipment for contamination by amplification product. It is very important to test all areas that may have been exposed to processed specimens, controls, and calibrators, and/or amplification product. This includes routinely handled objects such as pipettes, the Abbott *m2000sp* and Abbott *m2000rt* function keys, laboratory bench surfaces, microcentrifuges, and centrifuge adaptors.

1. Add 0.8 mL RNase-free water to a 1.7 mL molecular biology grade microcentrifuge tube.
2. Saturate the cotton tip of an applicator (Puritan or equivalent) in the RNase-free water from the microcentrifuge tube.
3. Using the saturated cotton tip of the applicator, wipe the area to be monitored using a sweeping motion. Place the applicator into the microcentrifuge tube.
4. Swirl the cotton tip in RNase-free water 10 times, and then press the applicator along the inside of the tube so that the liquid drains back into the solution at the bottom of the microcentrifuge tube. Discard the applicator.
5. Pipette 0.5 mL of *mWash* 1 buffer to a clean tube using the pipette dedicated for Internal Control use.
6. Add 20 µL of the *mWash* 1 buffer to each microcentrifuge tube.
7. Cap the microcentrifuge tube.

8. Test this sample according to the assay procedure section of this package insert.
 - Transfer liquid from the microcentrifuge tube to a 5 mL Reaction Vessel.
 - Bring the volume to 1.5 mL with RNase-free water.
9. The presence of contamination is indicated by the detection of HCV nucleic acid in the swab samples.
10. If HCV nucleic acid is detected on equipment, follow the cleaning and decontaminating guidelines given in that equipment's operations manual. If HCV nucleic acid is detected on surfaces, clean the contaminated areas with 1.0% (v/v) sodium hypochlorite solution, followed by 70% ethanol or water.

NOTE: Chlorine solutions may pit equipment and metal. Use sufficient amounts or repeated applications of 70% ethanol or water until chlorine residue is no longer visible.

11. Repeat testing of the contaminated area by following steps 1 through 10.

RESULTS

Calculation

The concentration of viral HCV RNA in a sample or control is calculated from the stored calibration curve. The Abbott *m2000rt* instrument automatically reports the results on the Abbott *m2000rt* workstation. Assay results can be reported in IU/mL or Log IU/mL.

INTERPRETATION OF RESULTS

Sample Volume	Result	Interpretation
0.5 mL	Not Detected	Target not detected
	< 1.08 Log IU/mL ^a	Detected
	1.08 to 8.00 Log IU/mL	
	> 8.00 Log IU/mL	> ULQ ^c
0.2 mL	Not Detected	Target not detected
	< 1.48 Log IU/mL ^b	Detected
	1.48 to 8.00 Log IU/mL	
	> 8.00 Log IU/mL	> ULQ

^a 12 IU/mL

^b 30 IU/mL

^c ULQ = upper limit of quantitation

LIMITATIONS OF THE PROCEDURE

- **FOR IN VITRO DIAGNOSTIC USE**
- Optimal performance of this test requires appropriate specimen collection, handling, preparation, and storage (refer to the **SPECIMEN COLLECTION, STORAGE, AND TRANSPORT TO THE TEST SITE** section of this package insert).
- Human serum and plasma specimens (ACD-A or EDTA) may be used with the Abbott RealTime HCV assay. The use of other anticoagulants has not been validated with the Abbott RealTime HCV assay.
- Use of the Abbott RealTime HCV assay is limited to personnel who have been trained in the procedures of a molecular diagnostic assay and the Abbott *m2000sp* and Abbott *m2000rt* instruments.
- The instruments and assay procedures reduce the risk of contamination by amplification product. However, nucleic acid contamination from the calibrators, positive controls, or specimens must be controlled by good laboratory practice and careful adherence to the procedures specified in this package insert.
- As with any diagnostic test, results from the Abbott RealTime HCV assay should be interpreted in conjunction with other clinical and laboratory findings. A specimen with the result of "Not Detected" cannot be presumed to be negative for HCV RNA.

SPECIFIC PERFORMANCE CHARACTERISTICS

The following performance characteristics were determined using the Abbott RealTime HCV assay with the Abbott *m2000sp* and the 0.5 mL sample preparation procedure, unless otherwise specified.

Limit of Detection (LoD)

The LoD is defined as the HCV RNA concentration detected with a probability of 95% or greater.

Limit of Detection, 0.5 mL Sample Volume

The LoD of the Abbott RealTime HCV assay is 12 IU/mL with the 0.5 mL sample preparation procedure.

The LoD was determined by testing dilutions of the Second WHO International Standard for Hepatitis C Virus RNA (NIBSC 96/798) prepared in HCV negative human plasma. Testing was performed with 3 lots of amplification reagents on 3 instrument systems. The results, representative of the analytical sensitivity of the Abbott RealTime HCV assay, are summarized in Table 1.

Table 1. Detection Rates for 0.5 mL Sample Volume (LoD)

(IU/mL)	Number Tested	Number Detected	Percent Detected
25.0	57	57	100
20.0	57	57	100
15.0	57	55	96
12.5	57	53	93
10.0	57	56	98
7.5	57	51	89
5.0	57	46	81
2.5	57	33	58

Probit analysis of the data determined that the concentration of HCV RNA detected with 95% probability was 10.5 IU/mL (95% CI 8.6-14.0 IU/mL).

Limit of Detection, 0.2 mL Sample Volume

The LoD of the Abbott RealTime HCV assay is 30 IU/mL with the 0.2 mL sample preparation procedure and was determined as described for the 0.5 mL sample preparation procedure. The results, representative of the analytical sensitivity performance of the Abbott RealTime HCV assay, are summarized in Table 2.

Table 2. Detection Rates for 0.2 mL Sample Volume (LoD)

(IU/mL)	Number Tested	Number Detected	Percent Detected
50.0	60	59	98
35.0	60	60	100
25.0	60	59	98
20.0	59 ^a	48	81
15.0	60	53	88
12.5	60	54	90
10.0	60	39	65
7.5	60	27	45


^a One sample was inhibited and was excluded from the data analysis.

Probit analysis of the data determined that the concentration of HCV RNA detected with 95% probability was 23.8 IU/mL (95% CI 17.4-59.7 IU/mL).


Linear Range

The upper limit of quantitation (ULQ) for the Abbott RealTime HCV assay is 100 million IU/mL and the lower limit of quantitation is equivalent to LoD (12 IU/mL for the 0.5 mL sample preparation procedure, 30 IU/mL for the 0.2 mL sample preparation procedure).

A 9-member panel prepared by diluting HCV armored RNA from 8.21 log IU/mL to 0.91 log IU/mL in HCV negative human plasma was tested. The results, representative of the Abbott RealTime HCV assay linearity, are shown in Figure 1.

Figure 1. Linearity in Plasma

A 9-member panel prepared by diluting HCV armored RNA from 8.33 log IU/mL to 1.03 log IU/mL in HCV negative human serum was tested. The results, representative of the Abbott RealTime HCV assay linearity, are shown in Figure 2.

Figure 2. Linearity in Serum

Linearity analysis was performed following the CLSI EP06-A²⁶ guideline. The Abbott RealTime HCV assay was shown to be linear in serum and plasma across the range of HCV RNA concentrations tested.

Precision

The precision of the Abbott RealTime HCV assay was evaluated for the 0.5 mL sample preparation procedure using the Abbott *m2000sp* instrument and the manual sample preparation method. An 8-member HCV RNA panel was prepared. Panel members 1, 3, and 5 were dilutions of an HCV viral stock in negative human plasma. Panel members 2 and 4 were dilutions of an HCV viral stock in negative human serum. Panel members 6 through 8 were prepared by diluting HCV armored RNA in negative human plasma. One lot of amplification reagents was assigned to each of 3 Abbott *m2000* instruments. For the precision study using the Abbott *m2000sp*, the panel members were tested in replicates of 4 in the first run on each instrument pair and replicates of 5 in each subsequent run, for a total of 15 runs. For the precision study using the manual sample preparation method, new dilutions of panel members 1, 2, 3, 5, 7, and 8 were prepared. These panel members were tested in replicates of 2 or 3 in 5 runs on each Abbott *m2000rt* instrument. Analysis was performed according to CLSI EP10-A2.²⁷ Within-run, between-run, and inter-assay (within-run and between-run) standard deviations were determined. The Abbott RealTime HCV assay was designed to achieve an inter-assay standard deviation (SD) of less than or equal to 0.25 log IU/mL of HCV RNA for samples containing HCV concentrations from 100 to 100 million IU/mL. The results, representative of the precision of the Abbott RealTime HCV assay, are summarized in Tables 3 and 4.

Table 3. Precision with the Abbott *m2000sp* and Abbott *m2000rt* Instruments

Panel Member	n	Conc. Mean (IU/mL)	Conc. Mean (Log IU/mL)	Within-Run Component SD ^a	Between-Run Component SD ^a	Inter-Assay SD ^{a,b}
1	65 ^c	8	0.88	0.25	0.00	0.25
2	72	91	1.96	0.08	0.04	0.09
3	72	564	2.75	0.05	0.03	0.06
4	72	8,853	3.95	0.03	0.02	0.04
5	71 ^d	61,320	4.79	0.04	0.02	0.04
6	71 ^d	925,118	5.97	0.04	0.02	0.04
7	71 ^d	10,045,506	7.00	0.07	0.01	0.07
8	71 ^d	109,257,360	8.04	0.04	0.01	0.04

^a Standard deviations (SD) are in Log IU/mL.

^b Inter-assay contains within-run and between-run components.

^c HCV RNA was not detected in 7 replicates.

^d One replicate was excluded from the data analysis due to an instrument error.

Table 4. Precision with Manual Sample Preparation Method

Panel Member	n	Conc. Mean (IU/mL)	Conc. Mean (Log IU/mL)	Within-Run Component SD ^a	Between-Run Component SD ^a	Inter-Assay SD ^{a,b}
1	41 ^c	8	0.93	0.20	0.00	0.20
2	42	79	1.90	0.09	0.03	0.10
3	42	651	2.81	0.06	0.08	0.10
5	41 ^d	51,114	4.71	0.09	0.08	0.12
7	42	9,628,784	6.98	0.07	0.01	0.07
8	42	95,761,398	7.98	0.05	0.04	0.06

^a Standard deviations (SD) are in log IU/mL.

^b Inter-assay contains within-run and between-run components.

^c HCV RNA was not detected in 1 replicate.

^d One replicate was excluded from the data analysis due to an instrument error.

Potentially Interfering Substance

The susceptibility of the Abbott RealTime HCV assay to interference by elevated levels of endogenous substances was evaluated. HCV negative samples and samples containing 1,000 IU/mL of HCV RNA were tested according to Table 5.

Table 5. Potential Interfering Substances

Potential Interfering Substance	No. Negative Samples Tested	No. Positive Samples Tested
Bilirubin	10	10
Bilirubin Diluent ^a	10	10
Hemoglobin	10	10
Protein	10	10
Liposyn	10	10
Liposyn Diluent ^a	10	10

^a Control condition

No interference in the performance of the Abbott RealTime HCV assay was observed in the presence of the following substances for all HCV positive and negative samples tested:

- Hemoglobin 500 mg/dL
- Liposyn (Triglycerides) 3000 mg/dL
- Bilirubin 20 mg/dL
- Protein 9 g/dL

Antivirals and antibiotics at concentrations in excess of peak plasma or serum levels were tested in 5 pools. No interference in the performance of the Abbott RealTime HCV assay was observed in the presence of the drug pools listed in Table 6 for all HCV positive and negative samples tested.

Table 6. Potentially Interfering Therapeutic Drugs by Drug Pool

Drug Pool	Drugs Tested
1	Zidovudine, Saquinavir, Ritonavir, Clarithromycin, Interferon 2a, Interferon 2b
2	Abacavir sulfate, Amprenavir, Peginterferon 2a, Peginterferon 2b, Ribavirin
3	Tenofovir disoproxil fumarate, Lamivudine, Indinavir sulfate, Ganciclovir, Valganciclovir hydrochloride, Acyclovir
4	Stavudine, Efavirenz, Lopinavir, Enfuvirtide, Ciprofloxacin
5	Zalcitabine, Nevirapine, Nelfinavir, Azithromycin, Valacyclovir

The assay was also evaluated by testing 70 specimens that had been either obtained from individuals diagnosed or screened for an autoimmune disorder (systemic lupus erythematosus [SLE], anti-nuclear antibodies [ANA], and rheumatoid factor [RF]) or serologically characterized as positive for the following viral markers: HBsAg, anti-HTLV-I/II, anti-HIV-1, anti-HIV-2. An additional 10 specimens that were PCR positive for a flavivirus (West Nile virus [n=4] and GB virus C [n=6]) were tested.

HCV RNA was not detected in 77 of the 80 specimens tested. HCV RNA was detected, but less than LoD, in 3 specimens (1 RF, 1 SLE and 1 anti-HIV-1). Insufficient volume did not allow for confirmation.

Specificity

The target specificity of the Abbott RealTime HCV assay is greater than or equal to 99.5% after resolution. In a representative study, specificity was evaluated by testing 56 HCV seronegative serum and 56 HCV seronegative plasma specimens. The specimens were tested on 3 Abbott m2000 instrument systems with 3 lots of amplification reagents. HCV RNA was not detected, resulting in 100% (112/112) specificity (95% CI 96.76-100.00%).

Cross-Reactivity

The viruses and microorganisms listed in Table 7 were evaluated for potential cross-reactivity in the Abbott RealTime HCV assay. Purified nucleic acid or viral lysate from each microorganism or virus was added to HCV RNA negative samples and samples that contained 1,000 IU/mL HCV RNA.

Table 7. Potentially Cross-Reactive Microorganisms/Viruses


Microorganism / Virus	Microorganism / Virus
Human immunodeficiency virus 1	Vaccinia virus
Human immunodeficiency virus 2	BK human polyomavirus
Human T-lymphotropic virus 1	Human papilloma virus 16
Hepatitis B virus	Human papilloma virus 18
Epstein-Barr virus	<i>Neisseria gonorrhoeae</i>
Herpes simplex virus 1	<i>Chlamydia trachomatis</i>
Herpes simplex virus 2	<i>Candida albicans</i>
Cytomegalovirus	<i>Staphylococcus aureus</i>
Human herpesvirus 6B	<i>Staphylococcus epidermidis</i>
Human herpesvirus 8	<i>Mycobacterium gordonae</i>
Varicella-zoster virus	<i>Mycobacterium smegmatis</i>
Dengue virus 4	

No interference in the performance of the Abbott RealTime HCV assay was observed in the presence of the potential cross-reactants for all positive and negative samples tested.

Detection and Quantitation of HCV Genotypes

The ability of the Abbott RealTime HCV assay to detect and quantitate HCV genotypes was evaluated in 3 studies. In the first study, dilution linearity was demonstrated by diluting 6 specimens, one of each genotype 1 through 6, to target concentrations of 5.0 log IU/mL, 3.5 log IU/mL, 2.0 log IU/mL, and 1.0 log IU/mL. Four replicates were tested at each concentration of each genotype. These representative data are shown in Figure 3. The correlation coefficients ranged from 0.994 to 0.998.

Figure 3. Linearity Across HCV Genotypes


In the second study, a total of 75 specimens representative of each genotype were tested with the Abbott RealTime HCV assay. The results are summarized in Table 8. The correlation between the Abbott RealTime HCV assay results and those of a comparator assay is shown in Figure 4. Method comparison analysis was performed according to CLSI Document EP09-A2.²⁶

Table 8. HCV Genotype Detection

HCV Genotypes	Number of Specimens Tested	Number of Specimens Detected	Percent Detected
1	14	14	100
1a	3	3	100
1b	12	12	100
2	9	9	100
2a/2c	1	1	100
2b	2	2	100
3a	12	12	100
3b	1	1	100
3c	1	1	100
4	9	9	100
4a	3	3	100
4c/4d	1	1	100
5a	4	4	100
6	3	3	100

Figure 4. HCV Genotype Detection - Correlation between Abbott RealTime HCV and Comparator

In the third study, purified RNA transcripts from genotypes 1 through 6 were diluted to a target concentration of 1,000 IU/mL and tested. The results are summarized in **Table 9**. The results demonstrated equivalent quantitation of the genotypes tested.


Table 9. HCV Genotype Transcript Quantitation

HCV Genotype	Mean Concentration (log IU/mL)
1	3.00
2	3.04
3	2.98
4	2.80
5	2.82
6	3.14

Correlation

Specimens from 131 HCV infected patients were tested with the Abbott RealTime HCV assay and a comparator assay. Method comparison analysis was performed following **CLSI EP09-A2**.²⁸ The correlation plot is shown in **Figure 5**. One specimen was identified as an outlier and upon re-analysis excluding the outlier, the correlation coefficient was 0.950.

Figure 5. Assay Correlation between Abbott RealTime HCV and Comparator

BIBLIOGRAPHY

- Clarke B. Molecular virology of hepatitis C virus. *J Gen Virol.* 1997;78 (Pt 10):2397-2410. doi:10.1099/0022-1317-78-10-2397.
- Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW, Houghton M. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science.* 1989;244(4902):359-362. doi:10.1126/science.2523562.
- Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science.* 1989;244(4902):362-364. doi:10.1126/science.2496467.
- EASL International Consensus Conference on Hepatitis C. Paris, 26-28, February 1999, Consensus Statement. European Association for the Study of the Liver. *J Hepatol.* 1999;30(5):956-961.
- NIH Consensus Statement on Management of Hepatitis C: 2002. *NIH Consens State Sci Statements.* 2002;19(3):1-46.
- McHutchison JG, Gordon SC, Schiff ER, et al. Interferon alpha-2b alone or in combination with ribavirin as initial treatment for chronic hepatitis C. Hepatitis Interventional Therapy Group. *N Engl J Med.* 1998;339(21):1485-1492. doi:10.1056/NEJM199811193392101.
- Davis GL, Esteban-Mur R, Rustgi V, et al. Interferon alpha-2b alone or in combination with ribavirin for the treatment of relapse of chronic hepatitis C. International Hepatitis Interventional Therapy Group. *N Engl J Med.* 1998;339(21):1493-1499. doi:10.1056/NEJM199811193392102.
- Manns MP, McHutchison JG, Gordon SC, et al. Peginterferon alpha-2b plus ribavirin compared with interferon alpha-2b plus ribavirin for initial treatment of chronic hepatitis C: a randomised trial. *Lancet.* 2001;358(9286):958-965. doi:10.1016/s0140-6736(01)06102-5.
- Fried MW, Shiffman ML, Reddy KR, et al. Peginterferon alpha-2a plus ribavirin for chronic hepatitis C virus infection. *N Engl J Med.* 2002;347(13):975-982. doi:10.1056/NEJMoa020047.
- Hadziyannis SJ, Sette H Jr, Morgan TR, et al. Peginterferon-alpha2a and ribavirin combination therapy in chronic hepatitis C: a randomized study of treatment duration and ribavirin dose. *Ann Intern Med.* 2004;140(5):346-355. doi:10.7326/0003-4819-140-5-200403020-00010.
- Strader DB, Wright T, Thomas DL, Seeff LB; American Association for the Study of Liver Diseases. Diagnosis, management, and treatment of hepatitis C. *Hepatology.* 2004;39(4):1147-1171. doi:10.1002/hep.20119.
- Gerken G, Rothaar T, Rumi MG, et al. Performance of the COBAS AMPLICOR HCV MONITOR test, version 2.0, an automated reverse transcription-PCR quantitative system for hepatitis C virus load determination. *J Clin Microbiol.* 2000;38(6):2210-2214. doi:10.1128/JCM.38.6.2210-2214.2000.
- Germer JJ, Heimgartner PJ, Ilstrup DM, Harmsen WS, Jenkins GD, Patel R. Comparative evaluation of the VERSANT HCV RNA 3.0, QUANTIPLEX HCV RNA 2.0, and COBAS AMPLICOR HCV MONITOR version 2.0 Assays for quantification of hepatitis C virus RNA in serum. *J Clin Microbiol.* 2002;40(2):495-500. doi:10.1128/JCM.40.2.495-500.2002.
- Leckie G, Schneider G, Abravaya K, et al. Performance attributes of the LCx HCV RNA quantitative assay. *J Virol Methods.* 2004;115(2):207-215. doi:10.1016/j.jviromet.2003.10.001.
- Saldanha J, Lelie N, Heath A. Establishment of the first international standard for nucleic acid amplification technology (NAT) assays for HCV RNA. WHO Collaborative Study Group. *Vox Sang.* 1999;76(3):149-158. doi:10.1159/000031040.
- Myers TW, Gelfand DH. Reverse transcription and DNA amplification by a *Thermus thermophilus* DNA polymerase. *Biochemistry.* 1991;30(31):7661-7666. doi:10.1021/bi00245a001.
- Smith DB, Mellor J, Jarvis LM, et al. Variation of the hepatitis C virus 5' non-coding region: implications for secondary structure, virus detection and typing. The International HCV Collaborative Study Group. *J Gen Virol.* 1995;76 (Pt 7):1749-1761. doi:10.1099/0022-1317-76-7-1749.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories, Fifth edition. US Government Printing Office; December 2009.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Occupational Exposure to Bloodborne Pathogens.

20. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Clinical and Laboratory Standards Institute; 2014.
21. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. World Health Organization; 2004.
22. Davis GL, Lau JY, Urdea MS, et al. Quantitative detection of hepatitis C virus RNA with a solid-phase signal amplification method: definition of optimal conditions for specimen collection and clinical application in interferon-treated patients. *Hepatology*. 1994;19(6):1337-1341.
23. Cuypers HT, Bresters D, Winkel IN, et al. Storage conditions of blood samples and primer selection affect the yield of cDNA polymerase chain reaction products of hepatitis C virus. *J Clin Microbiol*. 1992;30(12):3220-3224. doi:10.1128/jcm.30.12.3220-3224.1992.
24. Ginocchio CC, Wang XP, Kaplan MH, et al. Effects of specimen collection, processing, and storage conditions on stability of human immunodeficiency virus type 1 RNA levels in plasma. *J Clin Microbiol*. 1997;35(11):2886-2893. doi:10.1128/jcm.35.11.2886-2893.1997.
25. Sebire K, McGavin K, Land S, Middleton T, Birch C. Stability of human immunodeficiency virus RNA in blood specimens as measured by a commercial PCR-based assay. *J Clin Microbiol*. 1998;36(2):493-498. doi:10.1128/JCM.36.2.493-498.1998.
26. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Clinical and Laboratory Standards Institute (formerly NCCLS): Wayne, PA; 2003.
27. Clinical Laboratory Standards Institute. *Preliminary Evaluation of Quantitative Clinical Laboratory Methods; Approved Guideline – Second Edition*. CLSI Document EP10 A2. CLSI: Wayne, PA, 2006.
28. Clinical Laboratory Standards Institute. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition*. CLSI Document EP09-A2. CLSI: Wayne, PA, 2002.

TECHNICAL ASSISTANCE:

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

SUMMARY OF SAFETY AND PERFORMANCE STATEMENT

A summary of safety and performance (SSP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSP location after the launch of European Database on Medical Devices. Search for device using UDI-DI provided on the outer packaging of the device.

This product is licensed from Chiron Corporation and is covered by one or more of the following US patents 5,714,596, 5,712,088, 5,863,719, 5,372,928, 5,851,759, 6,074,816 and their foreign equivalents.

THE PURCHASE OF THIS PRODUCT ALLOWS THE PURCHASER TO USE IT FOR AMPLIFICATION OF NUCLEIC ACID SEQUENCES AND FOR DETECTION OF NUCLEIC ACID SEQUENCES FOR HUMAN *IN VITRO* DIAGNOSTICS. NO GENERAL PATENT OR OTHER LICENSE OF ANY KIND OTHER THAN THIS SPECIFIC RIGHT OF USE FROM PURCHASE IS GRANTED HEREBY. THIS PROVISION DOES NOT PROHIBIT THE RESALE OF THIS PRODUCT.

Armored RNA® is a patented technology jointly developed by Ambion, Inc. and Cenetron Diagnostics, LLC. US patents #5,677,124, #5,919,625, #5,939,262 and patents pending.

Armored RNA is a registered trademark of Ambion, Inc.

ProClin is a registered trademark of Rohm and Haas.

StrataCooler is a registered trademark of Stratagene.

FAM and ROX are trademarks of Life Technologies Corporation or its subsidiaries in the US and/or certain other countries.

VIC is a registered trademark of Life Technologies Corporation or its subsidiaries in the US and/or certain other countries.

Abbott *m*, *m2000*, *RealTime*, *m2000rt*, and *m2000sp* are trademarks of Abbott.

Abbott Molecular Inc. is the legal manufacturer of the

Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90)

Abbott RealTime HCV Control Kit (List No. 4J86-80)

Abbott RealTime HCV Calibrator Kit (List No. 4J86-70)

The Abbott RealTime HCV Amplification Reagent Kit is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

© 2005, 2024 Abbott. All Rights Reserved.

February 2024

51-602124/R13

Abbott RealTime HCV

Supplement for Dried Blood Spot (DBS)

en

GTIN

REF

04J86-010

IVD

REF 04J86-010

51-608440/R3

51-608440/R3

NOTE: Changes Highlighted

Key to Symbols Used

GTIN

Global Trade Item Number

REF

Reference Number

IVD

In Vitro Diagnostic Medical Device

Maximum Time Allowed

Upper Limit of Temperature

Temperature Limit

EC REP

Authorized Representative in the
European Community/European Union

Manufacturer

CUSTOMER SERVICE INTERNATIONAL: CALL YOUR ABBOTT REPRESENTATIVE

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

This Instructions For Use (IFU) is a supplement to the Abbott RealTime HCV assay (List No. 4J86) package insert, specifically for Dried Blood Spot (DBS) sample type.

NOTICE TO USER

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the customer service or technical assistance section of these instructions.

NAME

Abbott RealTime HCV

INTENDED USE FOR DBS

The Abbott RealTime HCV assay is an in vitro reverse transcription polymerase chain reaction (RT-PCR) assay for the quantitation of hepatitis C viral RNA (HCV RNA) in whole blood spotted on cards as dried blood spots (DBS) (ie, obtained via venipuncture or capillary blood) from HCV-infected individuals. The Abbott RealTime HCV assay is intended for use as an aid in the management of HCV-infected patients undergoing antiviral therapy. The Abbott RealTime HCV assay is not for screening blood, plasma, serum or tissue donors for HCV, or to be used as a diagnostic test to confirm the presence of HCV infection.

INTENDED USER

The intended users for the Abbott RealTime HCV assay are laboratory professionals.

SUMMARY AND EXPLANATION OF THE TEST

Hepatitis C virus (HCV) is one of the major causes of liver disease. It was estimated that about 1.1% of the population are chronically infected with hepatitis C¹ and approximately 700,000 people die each year from hepatitis C-related liver disease.² HCV is a positive-strand RNA virus and transmitted primarily through intravenous drug use and through blood products. About 55 to 85% of HCV-infected individuals develop chronic hepatitis, with up to 30% of chronically infected individuals developing

cirrhosis.¹ In patients with cirrhosis, the incidence of hepatocellular carcinoma is 2 to 4% per year.¹

Quantitative measurement of HCV RNA level in peripheral blood has been shown to be an essential parameter in management of various anti-HCV therapies.⁶⁻¹³ With the advent of direct acting antiviral regimens, the treatments for HCV infection are advanced dramatically. The objective of anti-HCV therapies is to achieve a sustained virologic response (SVR), presented as continued absence of detectable HCV RNA.²⁻⁵ Measurements of HCV viral load have been used to individualize treatment duration and response-guided therapy or to determine and predict sustained or non-sustained virologic response of antiviral therapy.²⁻⁵

The majority of global HCV burden is found in low- and middle-income countries (LMICs).¹⁴ Identification of patients with active HCV infection is very limited in these resource-constrained areas. Whole blood collected as Dried Blood Spot (DBS) can be used as an alternative to serum or plasma for HCV RNA testing. DBS represents a feasible option that eliminates many logistical and technical limitations associated with serum or plasma collection, processing, storage, and shipment. World Health Organization (WHO) Guidelines on Hepatitis C testing include recommendations for using DBS in HCV virological testing in resource-limited areas.¹⁴

The Abbott RealTime HCV assay for DBS uses reverse transcription polymerase chain reaction (RT-PCR) technology combined with homogeneous real time fluorescent detection for the quantitation of HCV RNA in DBS. The selection of a conserved region of the HCV genome provides for the detection of genotypes 1 through 6. The assay is standardized against the 2nd WHO International Standard for Hepatitis C Virus RNA (NIBSC Code 96/798)¹⁵ and results are reported in international units/mL (IU/mL) or its log base 10 equivalent.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The Abbott RealTime HCV assay for DBS requires 3 reagent kits:

- Abbott RealTime HCV Amplification Reagent Kit
- Abbott RealTime HCV Control Kit
- Abbott RealTime HCV Calibrator Kit

The Abbott RealTime HCV assay for DBS uses RT-PCR to generate amplified product from the HCV RNA genome in human DBS specimens. The steps of the Abbott RealTime HCV assay for DBS consist of pre-analytical DBS treatment, sample preparation, RT-PCR assembly, amplification, detection, and result reporting. Pre-analytical DBS treatment is performed using the Abbott mSample Preparation System DBS Buffer Kit, prior to specimen loading on the Abbott m2000sp system. Following the pre-analytical DBS treatment, all subsequent steps are executed automatically by the Abbott m2000sp and Abbott m2000rt systems.

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert for more details on sample preparation, RT-PCR assembly, amplification, and detection.

AMPLIFICATION REAGENT EXTENDED USE

The optional amplification reagent extended use feature allows amplification reagent packs and IC to be used a total of 2 times. Amplification reagent packs containing prepared master mix may be stored at -10°C or colder, capped and protected from light, for up to 14 days before the second use. The internal control (IC) may be stored capped at -10°C or colder for up to 14 days before the second use. The second use should not exceed the expiration date of the reagent.

An overview of this feature is provided in Appendix 1 of this IFU.

Throughout this IFU, amplification reagent packs and IC that have not yet been used will be referred to as new amplification reagent packs and IC (ie, initial use). Amplification reagent packs that have been used once and contain prepared master mix will be referred to as partial amplification reagent packs. IC vials that have been used once will be referred to as partial IC vials.

PREVENTION OF NUCLEIC ACID CONTAMINATION

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

REAGENTS

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

WARNINGS AND PRECAUTIONS

IVD

For In Vitro Diagnostic Use

The Abbott RealTime HCV assay for DBS is not for screening blood, plasma, serum, or tissue donors for HCV, or to be used as a diagnostic test to confirm the presence of HCV infection.

SAFETY PRECAUTIONS

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Abbott RealTime HCV assay for DBS only uses Abbott *m2000sp* and Abbott *m2000rt* instruments.

SPECIAL PRECAUTIONS

Handling Precautions for DBS Specimens

The Abbott RealTime HCV assay for DBS is only for use with whole blood or DBS specimens that have been handled and stored as described in the **SPECIMEN COLLECTION, STORAGE, AND TRANSPORT TO THE TEST SITE** section.

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Handling precautions specific for DBS:

- Ensure that DBS cards are fully dried prior to storage in bag with desiccants. Inappropriate spotting, drying, storage and handling of DBS may lead to inaccurate test results.
- To avoid carryover contamination, use a new pipette tip for each DBS sample when pushing the spot out from the perforated DBS card. If a non-perforated DBS card is used, ensure that the tool used to excise the spot is decontaminated prior to and after contact with the DBS sample.

Work Areas

Only two dedicated areas, Sample Preparation Area and Amplification Area, are recommended when the Abbott *m2000sp* and Abbott *m2000rt* are used.

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Aerosol Contamination

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Contamination and Inhibition


Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Contamination From External dU-Containing Amplified Product

Laboratories that use or have used HCV amplification assays that include post-PCR processing of the amplified product may be contaminated by dU-containing amplified product. Such contamination may cause inaccurate results in the Abbott RealTime HCV assay for DBS. Refer to the **Monitoring the Laboratory for the Presence of Contamination** section of the Abbott RealTime HCV assay (List No. 4J86) package insert. When negative controls are persistently reactive or where contamination with dU-containing HCV amplified product is likely to have occurred, it is recommended that the laboratory use uracil-N-glycosylase (UNG) contamination control procedure if decontamination of the laboratory is unsuccessful.

STORAGE INSTRUCTIONS

Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90)

- **New** and **Partial** Abbott RealTime HCV Amplification Reagent Packs and Internal Control (IC) vials must be stored at -10°C or colder when not in use. Care must be taken to separate the Abbott RealTime HCV Amplification Reagent Pack that is in use from direct contact with samples, calibrators, and controls.
- **Partial** amplification reagent packs and IC must be stored at -10°C or colder, capped, upright, and protected from light, following initial use. If stored this way, **partial** amplification reagent packs with prepared master mix may be used a second time within 14 days of initial use. IC may also be used a second time within 14 days of being thawed, if stored capped at -10°C or colder. After 2 uses, partial amplification reagent packs and IC should be discarded. The second use should not exceed the expiration date of the reagent.

Abbott RealTime HCV Control Kit (List No. 4J86-80)

- The Abbott RealTime HCV Negative and Positive Controls must be stored at -10°C or colder.

Abbott RealTime HCV Calibrator Kit (List No. 4J86-70)

- The Abbott RealTime HCV Calibrator A and Calibrator B must be stored at -10°C or colder.

SHIPPING CONDITIONS

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

INDICATION OF INSTABILITY OR DETERIORATION OF REAGENTS

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

INSTRUMENT PROCEDURE

The Abbott RealTime HCV for DBS application file must be installed on the Abbott *m2000sp* and Abbott *m2000rt* systems from the Abbott RealTime HCV *m2000* ROW System Combined Application CD-ROM (List No. 01L69-040 or higher) prior to performing the assay. For detailed information on application file installation, refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals, Operating Instructions section.

SPECIMEN COLLECTION, STORAGE, AND TRANSPORT TO THE TEST SITE

Specimen Collection and Storage

- To prepare DBS, use finger prick or EDTA whole blood. EDTA whole blood should be maintained under controlled temperature conditions during storage and shipment (2°C to 8°C for up to 48 hours, or ambient temperature not exceeding 30°C for up to 24 hours).
- EDTA whole blood collected in a tube should be mixed thoroughly by gentle inversions prior to spotting. Spot the blood onto the one-half-inch (12 millimeter) circles on Ahlstrom-Munktel TFN paper card or Whatman 903, ensuring that the entire circle is covered. It is recommended that at least 70 μL of blood (approximately 3 to 5 blood drops) be used for each circle to ensure full coverage.
- Air dry the card at ambient temperature for a minimum of 3 hours. Package each card in a sealable bag with desiccant packs (at least 2 to 3 desiccant packs).
- DBS specimen cards can be stored under the conditions as specified in the table below prior to assay processing.

DBS Specimen Card Storage

Condition	Maximum Storage Time
2°C to 30°C at ambient humidity	18 weeks
$\leq 37^{\circ}\text{C}$ at ambient humidity	12 weeks
$\leq 45^{\circ}\text{C}$ at ambient humidity	8 weeks
$\leq 38^{\circ}\text{C}$ at high humidity ($\leq 85\%$)	4 weeks
-25°C to -15°C at ambient humidity	18 weeks
-70°C or colder at ambient humidity	18 weeks

Specimen Transport

Ship specimens according to the recommended storage temperature and time listed in the **Specimen Collection and Storage** section above. For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical, diagnostic, or biological specimens.

ABBOTT REALTIME HCV ASSAY PROCEDURE FOR DBS

Materials Required But Not Provided

- Abbott RealTime HCV Amplification Reagent Kit (List No. 4J86-90)
- Abbott RealTime HCV Calibrator Kit (List No. 4J86-70)
- Abbott RealTime HCV Control Kit (List No. 4J86-80)

For Abbott *m2000sp* Instrument Sample Preparation Area

- Abbott *mSample* Preparation System DBS Buffer Kit (List No. 09N02-001)
- Abbott *m2000sp* instrument (List No. 09K14-01, 09K14-02)
- Abbott *mSample* Preparation System (4 $\times$ 24 Preps) (List No. 04J70-24)
- Calibrated precision pipettes capable of delivering 20 to 1000 μL (Calibrated precision pipettes capable of delivering < 20 μL may be required if using UNG.)

- Aerosol barrier pipette tips for 20 to 1000 µL pipettes (Aerosol barrier pipette tips capable of delivering < 20 µL may be required if using UNG.)
- Vortex Mixer
- Master Mix Tube (List No. 04J71-80)
- Abbott 96-Deep-Well Plate (List No. 04J71-30)
- 5 mL Reaction Vessels (List No. 04J71-20)
- 200 mL reagent vessels (List No. 04J71-60)
- Abbott 96-Well Optical Reaction Plate (List No. 04J71-70)
- Abbott Optical Adhesive Covers (List No. 04J71-75)
- Abbott Adhesive Cover Applicators (List 9K32-01)
- Abbott Splash-Free Support Base (List No. 9K31-01)
- Abbott RealTime HCV *m2000* ROW System Combined Application CD-ROM (List No. 01L69-040 or higher)
- *m2000* System 13mm DBS-PoS Set (List No. 09N03-001)
- 1.4 mL Micro Vial 15mm Caps (List No. 03N20-01) optional
- Uracil N-Glycosylase (UNG)^a (List No. 6L87-02)

^aIf needed, per the **Contamination From External dU-Containing Amplified Product** section of this IFU.

For Abbott *m2000rt* Instrument Amplification Area

- Abbott *m2000rt* instrument (List No. 9K15-01)
- Abbott RealTime HCV *m2000* ROW System Combined Application CD-ROM (List No. 01L69-040 or higher)
- Abbott *m2000rt* Optical Calibration Kit (List No. 04J71-93)

Other Materials

- Biological safety cabinet approved for working with infectious materials
- Sealable bag for DBS card
- Desiccant pack for DBS card
- Ahlstrom-Munktel TFN paper card and/or Whatman 903
- Sealable plastic bag for disposal of 96-well reaction plate
- RNase-free water (Eppendorf or equivalent)^b
- 1.7 mL molecular biology grade microcentrifuge tubes (Dot Scientific, Inc. or equivalent)^b
- Cotton Tip Applicators (Puritan or equivalent)^b

^bNote: These 3 items are used in the procedure for Monitoring the Laboratory for the Presence of Contamination. Refer to the **QUALITY CONTROL PROCEDURES** section in the Abbott RealTime HCV assay (List No. 4J86) package insert.

Procedural Precautions

- Read the instructions in this IFU carefully before processing specimens.
- Amplification reagents and internal control may be used up to 2 times, as described in this IFU. The Abbott RealTime HCV Calibrators, Negative Control, Low Positive Control, and High Positive Control vials are intended for single-use only and should be discarded after use.
- Use aerosol barrier pipette tips or disposable pipettes only 1 time when pipetting specimens or IC. To prevent contamination to the pipette barrel while pipetting, care should be taken to avoid touching the pipette barrel to the inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.
- Monitoring procedures for the presence of amplification product can be found in the **QUALITY CONTROL PROCEDURES** section.
- To reduce the risk of nucleic acid contamination, clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% sodium hypochlorite or other suitable disinfectant.
- The Abbott RealTime HCV Calibrators and Controls must be prepared in conjunction with the specimens to be tested. The use of the Abbott RealTime HCV Controls and Calibrators is integral to the performance of the Abbott RealTime HCV assay. Refer to the **QUALITY CONTROL PROCEDURES** section.
- The Abbott *m2000sp* Master Mix Addition protocol must be initiated within 1 hour after completion of Sample Preparation.
- If the Abbott *m2000sp* master mix addition protocol is aborted, then discard the amplification reagents that were being used for the second time. New reagent packs being used for the first time may be saved and stored for later use as described in this manual.

- The Abbott *m2000rt* protocol must be started within 50 minutes of the initiation of the Master Mix Addition protocol. If the Abbott *m2000rt* instrument run is not initiated within 50 minutes, or is interrupted or aborted, seal the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.

IMPORTANT: Amplification reagents that will be used a second time must be stored at – 10°C or colder within 50 minutes of the initiation of the master mix addition protocol.

ASSAY PROTOCOL: DBS SAMPLES USING THE ABBOTT *m2000sp* AND THE ABBOTT *m2000rt* INSTRUMENTS

For a detailed description of how to perform an Abbott *m2000sp* instrument and Abbott *m2000rt* instrument protocol, refer to the Abbott *m2000sp* and Abbott *m2000rt* Operations Manuals, Operating Instructions sections. The DBS protocol requires Abbott *m2000sp* Software Version 6.0 or higher. Please follow Abbott *m2000sp* Operations Manual (List 09K20 version 6 or higher).

A total of 96 samples can be processed in each run. An Abbott RealTime HCV Negative Control, Low Positive Control, and High Positive Control are included in each run, therefore allowing a maximum of 93 DBS specimens to be processed per run when calibrators are not included. Process Calibrators and Controls directly as liquid samples; step 1 through step 5 are for DBS only. Do not process plasma or serum specimens on any run where DBS protocol is used. For each DBS sample, a single Abbott Master Mix Tube should be used for the entire sample processing procedure. Ensure DBS samples are labeled throughout processing.

1. Fill the Abbott Master Mix Tube with 1.3 mL of *m*DBS buffer from the Abbott *m*Sample Preparation system DBS Buffer Kit (List No. 09N02-001).

NOTE: Do not use *m*Lysis Buffer or any other reagents for this step

2. Hold perforated DBS paper card above the Abbott Master Mix Tube.
3. Push the DBS circle out of the card using a clean pipette tip, 1 DBS circle per Abbott Master Mix Tube. Each DBS should be approximately one-half-inch (12 millimeters) in diameter.

NOTE: USE A NEW PIPETTE TIP FOR EACH DBS SAMPLE TO PREVENT CROSS-CONTAMINATION.

Alternatively, from non-perforated DBS paper, cut out 1 entire DBS (approximately one-half-inch in diameter) for each specimen and transfer to an Abbott Master Mix Tube.

NOTE: Avoid direct contact of the cutting surface with DBS specimens. Decontaminate the tool used to cut DBS between specimens if necessary, according to good laboratory practices.

4. Ensure that the DBS circle is fully submerged in the *m*DBS Buffer by tapping the tube or using the pipette tip to push the DBS into the buffer.

NOTE: If a pipette tip is used to push the DBS into the buffer, ensure that the pipette tip does not cause DBS buffer volume loss due to liquid containment in the tip and/or absorption of the buffer by the tip filter.

5. Manually shake or swirl the sample tubes. Do not vortex the samples. Incubate for 30 minutes (± 2 minutes) at ambient temperature.
6. Meanwhile, thaw assay controls and internal control (IC) at 15°C to 30°C or at 2°C to 8°C. Thaw calibrators at 15°C to 30°C or at 2°C to 8°C only if performing a calibration run.

NOTE: Once thawed, assay controls, IC, and calibrators can be stored at 2°C to 8°C for up to 24 hours before use.

7. Vortex each assay calibrator and each control 3 times for 2 to 3 seconds before use. Ensure that the contents of each vial are at the bottom by tapping the vials on the bench to bring liquid to the bottom of the vial. Ensure bubbles or foam is not generated; if present, remove with a sterile pipette tip, using a new tip for each vial.
8. Thaw amplification reagents at 15°C to 30°C or at 2°C to 8°C and store at 2°C to 8°C until required for the amplification master mix procedure.
9. Select **new** and/or **partial** amplification reagent packs to be used in the run. Refer to the Abbott *m2000sp* Operations Manual (List No. 09K20 version 6 or higher), Operating Instructions section, for instructions pertaining to amplification reagent pack inventory management. Amplification reagent packs must have the same lot number.

Thaw **new** amplification reagents at 15°C to 30°C or at 2°C to 8°C and store at 2°C to 8°C until required for the amplification master mix procedure.

NOTE: Once thawed, the new amplification reagents can be stored at 2°C to 8°C for up to 24 hours if not used immediately. Partial amplification reagent packs being used a second time should NOT be stored at 2°C to 8°C before use. They should be kept at –10°C or colder until needed for master mix addition. Once removed from the freezer, cumulative room temperature exposure should not exceed 25 minutes, including instances where packs are removed from storage, but not used. If 25 minutes is exceeded, discard the partial amplification reagent packs.

The following table shows the number of sample preparation reagents and internal control vials needed based on the number of reactions.

Sample Preparation Reagents and Internal Control Requirements				
Reagent	1 to 24 Reactions	25 to 48 Reactions	49 to 72 Reactions	73 to 96 Reactions
<i>m</i> Microparticles	1 bottle	2 bottles	2 bottles	2 bottles
<i>m</i> Lysis	1 bottle	2 bottles	3 bottles	4 bottles
<i>m</i> Wash 1	1 bottle	2 bottles	3 bottles	4 bottles
<i>m</i> Wash 2	1 bottle	2 bottles	3 bottles	4 bottles
<i>m</i> Elution Buffer	1 bottle	2 bottles	3 bottles	4 bottles
Internal Control ^a	1 new vial or 1 partial vial	1 new vial or 2 partial vials	2 new vials or 3 partial vials	2 new vials or 4 partial vials

^a A combination of new and partial vials of Internal Control may be used.

- Gently invert the Abbott *m*Sample Preparation bottles to ensure a homogeneous solution. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.
 - Vortex each IC 3 times for 2 to 3 seconds before use.
 - Use a calibrated precision pipette **DEDICATED FOR INTERNAL CONTROL USE ONLY** to add 500 µL of IC to each bottle of *m*Lysis Buffer. Mix by gently inverting the container 5 to 10 times to minimize foaming. Partial vials of IC can be recapped and stored at –10°C or colder for a second use.
 - Place the low and high positive controls, the negative control, and the calibrators, if applicable, into the Abbott *m*2000sp sample racks.
 - After the incubation is complete, manually shake or swirl the DBS sample tubes and then place them into the Abbott *m*2000sp sample racks.
- NOTE:** Ensure that the Abbott *m*2000sp sample racks have been calibrated specifically for the Abbott RealTime HCV procedure for DBS.
- Load the sample racks carefully to avoid splashing. If used, bar codes on tube labels must face right for scanning. Ensure that each tube is placed securely in the sample rack so that the bottom of the tube reaches the inside bottom of the rack.
 - Load filled sample racks onto the Abbott *m*2000sp in consecutive sample rack positions, with the first rack farthest to the right on the worktable, and any additional rack progressively to the left of the first rack.
 - Place the 5 mL Reaction Vessels into the Abbott *m*2000sp 1 mL subsystem carrier.
 - Load the Abbott *m*Sample Preparation System reagents and the Abbott 96 Deep-Well Plate on the Abbott *m*2000sp worktable as described in the Abbott *m*2000sp Operations Manual, Operating Instructions section.
 - From the Protocol screen, select the HCV RNA DBS application file. Initiate the sample extraction protocol as described in the Abbott *m*2000sp Operations Manual, Operating Instruction section.
 - Enter calibrator (needed if a calibration curve has not been stored on the Abbott *m*2000rt) and control lot specific values in the Sample Extraction: Worktable Setup, Calibrator and Control fields. Lot-specific values are specified in each Abbott RealTime HCV Calibrator and Control Kit Card.

- The Abbott *m*2000sp Master Mix Addition protocol (step 28) must be initiated within 1 hour after completion of Sample Preparation.

NOTE: Change gloves before handling the amplification reagents.

- Load the amplification reagents and the master mix tube (if needed) on the Abbott *m*2000sp worktable after sample preparation is completed. The following table shows the number of amplification reagent packs needed based on the number of reactions. If only 1 amplification reagent pack is being used, no master mix tube is required.

Amplification Reagent Pack Requirements ^a			
1 to 24 Reactions	25 to 48 Reactions	49 to 72 Reactions	73 to 96 Reactions
1 if new ; up to 4 with partial packs	2 if new ; up to 4 with partial packs	3 if new ; up to 4 with partial packs	4 new or partial packs

^a Refer to the Abbott *m*2000sp Operations Manual (List No. 09K20 version 6 or higher) for instructions on inventory management to determine the maximum number of reactions that can be tested with the partial packs selected.

- Partial amplification reagent packs can only be used on the same Abbott *m*2000sp instrument used for the reagent pack's initial preparation. Using an amplification reagent pack for a second time on a different instrument will result in an error, which may delay the run.
- Partial and new amplification reagent packs may be used together.

IMPORTANT: Partial amplification reagent packs should be stored at –10°C or colder until immediately before the second use.

Confirm that master mix is thawed before placing partial pack(s) on the Abbott *m*2000sp worktable. Once removed from –10°C or colder storage, partial amplification reagent packs being used a second time must be used within 25 minutes or discarded. This applies to cumulative room temperature exposure, including instances where packs are removed from storage, but not used.

- Ensure that the contents of new amplification reagent packs are at the bottom of the vials prior to opening the amplification reagents by tapping the vials in an upright position on the bench 5 to 10 times.
 - Do not tap partial amplification reagent packs being used a second time. Tapping may result in loss of master mix volume in the cap.
 - Remove caps. If a new amplification reagent pack will be stored for a second use, the vials will need to be recapped for storage. If planning to reuse the original caps to recap the reagent vials, save the original caps. If planning to use fresh caps to recap the reagent vials, original caps may be discarded.
- If the uracil N-glycosylase (UNG) is used for contamination control, follow the instructions below.**
- Use a PIPETTE DEDICATED FOR REAGENT USE ONLY to add specified volume of 1 U/µL UNG (List No. 6L87-02) to the reagent vial in position 3 of new and partial amplification reagent packs. Manually mix by pipetting gently up and down for all partial packs. Do not mix new packs.
 - Use the table below to determine the volume of UNG to add to the reagent vial in position 3 of new and partial amplification reagent packs. The reagent vial in position 3 of a new reagent pack contains the thermostable rTth polymerase enzyme. The reagent vial in position 3 of a partial reagent pack contains the master mix.

Volume of UNG to Add to the Reagent Vial in Position 3 of Each Amplification Reagent Pack	
Tests Remaining in the Pack	Add This Volume of UNG (µL)
1	5
2	6
3	7
4	8
5	9
6	10
7	11
8	12
9	13
10	13
11	14
12	15
13	16
14	17
15	18
16	19
17	20
18	21
19	22
20	23
21	24
22	25
23	26
24 (new pack)	27

NOTE: The volume of UNG added to the reagent vial in position 3 depends upon the number of tests remaining in the reagent pack, and not the number of samples being run. Refer to the Abbott *m2000sp* Operations Manual for instructions pertaining to amplification reagent pack inventory management and how to determine the number of tests remaining in a reagent pack.

- Partial amplification packs are loaded to the left of new amplification reagent packs on the Abbott *m2000sp* worktable.
- Ensure that amplification reagent packs are firmly seated on the instrument.
- Select the appropriate deep-well plate that matches the corresponding sample preparation extraction. Initiate the Abbott *m2000sp* Master Mix Addition protocol. Follow the instructions as described in the Abbott *m2000sp* Operations Manual, Operating Instructions section.

NOTE: The operator should not manually fill any empty/unfilled wells in the Abbott 96-Well Optical Reaction Plate.

- After sample extraction is complete, the Abbott *m2000sp* automatically fills any empty wells in the Abbott 96-Well Optical Reaction Plate with *mElution* buffer when there are greater than 48 samples processed within a run. Plate fill is not performed for runs containing 48 samples or fewer.
- If prompted by the instrument, Reagent Carrier 2 should remain in place, minimally containing the reagent vessel for *mElution* Buffer (Reagent Carrier 2, location 6). If this reagent vessel has been unloaded, place a new reagent vessel with the *mElution* Buffer label into Reagent Carrier 2, location 6. System fluid will be added to the reagent vessel and used to fill empty wells. Once this process is complete, the system will continue with the master mix addition.

NOTE: System instructions for use of the automated plate-filling feature are found in the Abbott *m2000sp* Operations Manual (List No. 09K20 version 6 or higher), section 5, Operating Instructions, Sample Extraction-Closed Mode.

- The Abbott *m2000rt* protocol (step 33) must be started within 50 minutes of the initiation of the Master Mix Addition protocol (step 28).

NOTE: If the run is aborted for any reason subsequent to step 28, a new Abbott 96-Well Optical Reaction Plate must be used if the Abbott *m2000sp* Master Mix Addition Protocol (step 28) will be repeated.

- Switch on and initialize the Abbott *m2000rt* instrument in the Amplification Area.

NOTE: The Abbott *m2000rt* requires 15 minutes to warm-up.

NOTE: Remove gloves before returning to the sample preparation area.

- Seal the Abbott 96-Well Optical Reaction Plate after the Abbott *m2000sp* instrument has completed addition of samples and master mix according to the Abbott *m2000sp* Operations Manual, Operating Instructions section.

- Place the sealed optical reaction plate into the Abbott Splash-Free Support Base for transfer to the Abbott *m2000rt* instrument.

If the uracil N-glycosylase (UNG) is used, follow the instructions below:

- Transfer the sealed optical plate to the Abbott Splash-Free Support Base and incubate at room temperature (15°C to 30°C) for 10 minutes. Following room temperature incubation, transfer the Abbott 96-Well Optical Reaction Plate on the Abbott Splash-Free Support Base to the Abbott *m2000rt* instrument.

- Place the Abbott 96-Well Optical Reaction Plate in the Abbott *m2000rt* instrument. From the Protocol screen, select the HCV RNA DBS application file. Initiate the protocol as described in the Abbott *m2000rt* Operations Manual, Operating Instructions section.

NOTE: Test order transfer through the use of CD-ROM or network connection with export and import features of the *m2000sp* and *m2000rt* software is recommended. If creating the Abbott *m2000rt* test order manually, enter sample IDs in the corresponding PCR tray locations according to the “Wells for Selected Plate” grid, found on the detail screen of the “PCR Plate Results” on the Abbott *m2000sp*. See Section 5 of the Abbott *m2000sp* Operations Manual.

- If a prepared partial amplification reagent pack is to be used a second time, cap the 3 reagent vials with the saved caps or new caps (List No. 3N20-01) and promptly store the reagents at –10°C or colder, protected from light, and in an upright position. Discard any amplification reagent packs that are exhausted or have been used twice.

IMPORTANT: Amplification reagents that will be used a second time must be stored at –10°C or colder within 50 minutes of the initiation of the master mix addition protocol.

POST PROCESSING PROCEDURES

- Remove the Abbott 96 Deep-Well Plate from the worktable and dispose of according to the Abbott *m2000sp* Operations Manual.
- Place the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.
- Clean the Abbott Splash-Free Support Base before next use, according to the Abbott *m2000rt* Operations Manual.

QUALITY CONTROL PROCEDURES

Abbott *m2000rt* Optical Calibration

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Assay Calibration

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert. For DBS testing using Abbott *m2000sp* and Abbott *m2000rt*: Once an Abbott RealTime HCV calibration is accepted and stored, it may be used for 6 months. During this time, all subsequent samples may be tested without further calibration unless:

- An Abbott RealTime HCV Amplification Reagent Kit with a new lot number is used.
- An Abbott *mSample* Preparation System (4 x 24 Preps) with a new lot number is used.
- A new Abbott RealTime HCV RNA DBS application file is installed.

Detection of Inhibition

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Negative and Positive Controls

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

Monitoring the Laboratory for the Presence of Contamination

Refer to the Abbott RealTime HCV assay (List No. 4J86) package insert.

RESULTS FOR DBS SPECIMENS

The reported DBS sample concentration result from the Abbott m2000rt DBS protocol run represents the HCV viral concentration in the plasma of the whole blood specimen from which the DBS specimen is obtained. The Abbott m2000rt instrument automatically reports the results on the Abbott m2000rt workstation. Assay results can be reported in International Units IU/mL, or Log IU/mL.

Interpretation of Results

Result	Interpretation
Not Detected	Target not detected
<2.66 Log IU/mL ^a	Detected
2.66 to 8.00 Log IU/mL	
>8.00 Log IU/mL	>ULQ ^b

^a 462 IU/mL

^b ULQ = upper limit of quantitation

LIMITATIONS OF THE PROCEDURE

- **FOR IN VITRO DIAGNOSTIC USE**
- Optimal performance of this test requires appropriate specimen collection, storage, and transport to the test site (refer to the **SPECIMEN COLLECTION, STORAGE, AND TRANSPORT TO THE TEST SITE** section of this package insert).
- Whole blood specimens collected in EDTA tubes may be used to generate DBS that are tested with the Abbott RealTime HCV assay for DBS. The use of other anticoagulants has not been validated with the Abbott RealTime HCV assay for DBS.
- Use of the Abbott RealTime HCV assay for DBS is limited to personnel who have been trained in the procedures of DBS preparation and pre-processing, a molecular diagnostic assay and the Abbott m2000 system.
- The instruments and assay procedures reduce the risk of contamination by amplification product. However, nucleic acid contamination from the calibrators, positive controls, or specimens must be controlled by good laboratory practices and careful adherence to the procedures specified in this package insert.
- As with any diagnostic test, results from the Abbott RealTime HCV assay for DBS should be interpreted in conjunction with other clinical and laboratory findings. A specimen with a result of "Not Detected" cannot be presumed to be negative for HCV RNA.
- Inappropriate collection, drying, storage, and handling of DBS cards may lead to inaccurate test results.
- Ahlstrom-Munktel TFN and Whatman 903 DBS cards may be used. The use of other cards has not been validated with Abbott RealTime HCV assay.

SPECIFIC PERFORMANCE CHARACTERISTICS FOR DBS SPECIMENS

Limit of Detection (LoD)

The LoD of the Abbott RealTime HCV assay is 462 IU/mL with the DBS sample type. The LoD was determined by testing DBS samples prepared with a dilution series of the WHO International Standard for Hepatitis C Virus RNA diluted in HCV negative whole blood. Testing for each HCV RNA concentration was performed with 4 lots of amplification reagents across 3 days. The results, representative of the analytical sensitivity performance of the Abbott RealTime HCV assay for DBS, are summarized in **Table 1**.

Table 1. Abbott RealTime HCV for DBS Limit of Detection (LoD)


HCV RNA (IU/mL)	Number Tested	Number Detected	Percent Detected
500	84	83	99
400	84	80	95
350	84	75	89
300	84	73	87
150	84	54	64
50	84	26	31

Probit analysis of the data determined that the concentration of HCV RNA detected from dried blood spot with 95% probability was 462 IU/mL (95% CI 381 to 593 IU/mL).

Linear Range

The upper limit of quantitation (ULQ) for the Abbott RealTime HCV assay is 100 million IU/mL and the lower limit of quantitation is equivalent to the LoD (462 IU/mL) for the DBS claim. A 7-member DBS panel was prepared with HCV armored RNA diluted to concentrations ranging from 200 million IU/mL to 316 IU/mL in HCV negative whole blood. The results, representative of the Abbott RealTime HCV assay linearity for DBS, are shown in **Figure 1**.

Figure 1. Abbott RealTime HCV Assay for DBS Linearity

The Abbott RealTime HCV assay for DBS is shown to be linear across the range of HCV RNA concentrations tested (316 to 200 million IU/mL, 2.5 to 8.3 Log IU/mL).

Specificity

The target specificity of the Abbott RealTime HCV assay for DBS is greater than or equal to 99.5%. Specificity was determined by analyzing 119 HCV sero-negative specimens, across 2 lots of amplification reagents. The overall specificity was 100% (119/119).

Precision

Precision was evaluated by testing a 6-member DBS panel targeted to cover the HCV RNA concentration range from 316 IU/mL to 200 million IU/mL. Panel members with concentrations lower than 5 Log IU/mL (100,000 IU/mL) were prepared with HCV positive clinical specimen diluted in HCV negative whole blood, while panel members with concentrations greater than 5 Log IU/mL (100,000 IU/mL) were prepared with HCV armored RNA diluted in HCV negative whole blood. Each panel member was tested in 5 replicates once per day for 5 days, on 3 Abbott m2000sp/rt pairs with 3 amplification reagent lots. Within run, between run, and inter-assay (within-run and between-run) standard deviations (SD) were determined. The Abbott RealTime HCV assay for DBS was designed to achieve an inter-assay standard deviation (SD) of less than or equal to 0.25 Log IU/mL of HCV RNA for samples containing HCV concentrations from 10,000 to 100 million IU/mL. The results, representative of the precision of the Abbott RealTime HCV assay for DBS, are summarized in **Table 2**.

Table 2. Abbott RealTime HCV Assay for DBS Precision

Panel Member	N	Mean Concentration (IU/mL)	Mean Concentration (Log IU/mL)	Within-Run Component SD ^{a,e}	Between-Run Component SD ^a	Inter-Assay SD ^{a,b}	Between-Lot/Instrument Component SD ^a	Total ^{c,e}
1	63 ^d	420	2.62	0.26	0.07	0.27	0.00	0.27
2	75	1,435	3.16	0.23	0.00	0.23	0.09	0.25
3	75	13,211	4.12	0.06	0.04	0.07	0.05	0.09
4	75	395,427	5.60	0.04	0.03	0.05	0.02	0.05
5	75	14,008,770	7.15	0.12	0.04	0.13	0.05	0.14
6	75	316,713,566	8.50	0.04	0.03	0.05	0.05	0.07

^a Standard deviations (SD) are in Log IU/mL.

^b Inter-Assay contains Within-Run and Between-Run components.

^c Total includes Within-Run, Between-Run, and Between-Lot/Instrument components.

^d Mean concentration of panel member 1 is below LLoQ. The precision estimates are based on the replicates that are quantitated.

^e Within-Run SD is also known as Repeatability. Total SD is also known as Reproducibility.

Repeatability and Reproducibility for Panel Member 1 (near LLoQ) were 0.26 Log IU/mL and 0.27 Log IU/mL respectively, and for Panel Member 2 (near LLoQ) were 0.23 Log IU/mL and 0.25 Log IU/mL respectively. Repeatability and reproducibility for Panel Member 5 (near ULQ) were 0.12 Log IU/mL and 0.14 Log IU/mL respectively, and for Panel Member 6 (near ULQ) were 0.04 Log IU/mL and 0.07 Log IU/mL respectively.

Correlation

HCV RNA quantitation was compared between the Abbott RealTime HCV assay using dried blood spots and the CE-marked comparator assay, Abbott RealTime HCV testing using human plasma. Plasma and DBS specimens from a total of 346 subjects (collected and tested in Georgia, Cameroon, and Greece) were included in the analysis. Both DBS prepared from venous blood and DBS prepared from capillary blood (finger prick) were tested by the Abbott RealTime HCV assay for DBS. The quantitative results from specimens that fell within the assays' dynamic ranges were analyzed by the least squares linear regression method and Bland-Altman bias method (DBS capillary versus plasma N=346, DBS venous versus plasma N=346, and DBS capillary versus DBS venous N=345). For the correlation between DBS capillary and plasma, the correlation coefficient was 0.843, the slope was 0.91 (95% CI 0.85 to 0.97), and the intercept was 0.48 Log IU/mL (95% CI 0.12 to 0.84) (Figure 2). For the correlation between DBS venous and plasma, the correlation coefficient was 0.840, the slope was 0.94 (95% CI 0.88 to 1.00), and the intercept was 0.31 Log IU/mL (95% CI -0.07 to 0.69) (Figure 3). For the correlation between DBS capillary and DBS venous, the correlation coefficient was 0.943, the slope was 0.91 (95% CI 0.87 to 0.94), and the intercept was 0.53 Log IU/mL (95% CI 0.33 to 0.72) (Figure 4). Additionally, Bland-Altman plots for these same comparisons are presented in Figure 5, Figure 6, and Figure 7. The mean bias difference between DBS capillary and plasma was -0.05 Log IU/mL. The mean bias difference between DBS venous and plasma was -0.04 Log IU/mL. The mean bias difference between DBS capillary and DBS venous was -0.02 Log IU/mL.

Predicted bias between the DBS Capillary and Plasma was evaluated at LLoQ and ULQ target levels. In addition, predicted bias between DBS Venous and Plasma was evaluated at LLoQ and ULQ target levels. Refer to Table 3 for the predicted bias summary.

Table 3. Predicted Bias Between DBS (Capillary or Venous) and Plasma at LLoQ and ULQ

Comparison	Level	Viral Load Target Level (Log IU/mL)	Predicted Bias (Log IU/mL)
DBS Capillary vs Plasma	LLoQ	2.66	0.24
	ULQ	8.00	-0.26
DBS Venous vs Plasma	LLoQ	2.66	0.15
	ULQ	8.00	-0.17

Figure 2. DBS Capillary and Plasma


Figure 5. DBS Capillary Versus Plasma


Figure 3. DBS Venous and Plasma


Figure 6. DBS Venous Versus Plasma


Figure 4. DBS Capillary and DBS Venous



Figure 7. DBS Capillary Versus DBS Venous

BIBLIOGRAPHY

1. Gower E, Estes C, Blach S, Razavi-Shearer K, Razavi H. Global epidemiology and genotype distribution of the hepatitis C virus infection. *J Hepatol.* 2014;61(1 Suppl):S45-S57. doi:10.1016/j.jhep.2014.07.027
2. Guidelines for the screening, care and treatment of persons with chronic hepatitis C infection. World Health Organization (WHO); updated April 2016. <https://apps.who.int/iris/handle/10665/205035>
3. AASLD/IDSA HCV guidance panel. *Recommendations for testing, managing, and treating hepatitis C.* AASLD/IDSA; updated October 24, 2022. <https://www.hcvguidelines.org/>
4. European Association for the Study of the Liver. EASL Recommendations on Treatment of Hepatitis C 2018. *J Hepatol.* 2018;69(2):461-511. doi:10.1016/j.jhep.2018.03.026
5. Miller MH, Agarwal K, Austin A, et al. Review article: 2014 UK consensus guidelines - hepatitis C management and direct-acting anti-viral therapy. *Aliment Pharmacol Ther.* 2014;39(12):1363-1375. doi:10.1111/apt.12764
6. Davis GL, Esteban-Mur R, Rustgi V, et al. Interferon alfa-2b alone or in combination with ribavirin for the treatment of relapse of chronic hepatitis C. International Hepatitis Interventional Therapy Group. *N Engl J Med.* 1998;339(21):1493-1499. doi:10.1056/NEJM199811193392102
7. Manns MP, McHutchison JG, Gordon SC, et al. Peginterferon alfa-2b plus ribavirin compared with interferon alfa-2b plus ribavirin for initial treatment of chronic hepatitis C: a randomised trial. *Lancet.* 2001;358(9286):958-965. doi:10.1016/s0140-6736(01)06102-5
8. Fried MW, Shiffman ML, Reddy KR, et al. Peginterferon alfa-2a plus ribavirin for chronic hepatitis C virus infection. *N Engl J Med.* 2002;347(13):975-982. doi:10.1056/NEJMoa020047
9. Hadziyannis SJ, Sette H Jr, Morgan TR, et al. Peginterferon-alpha2a and ribavirin combination therapy in chronic hepatitis C: a randomized study of treatment duration and ribavirin dose. *Ann Intern Med.* 2004;140(5):346-355. doi:10.7326/0003-4819-140-5-200403020-00010
10. Orlent H, Deltenre P, Francque S, et al. Update of the Belgian Association for the Study of the Liver guidelines for the treatment of chronic hepatitis C genotype 1 with protease inhibitors. *Acta Gastroenterol Belg.* 2012;75(2):245-259.
11. Halfon P, Bourlière M, Pénaranda G, Khiri H, Ouzan D. Real-time PCR assays for hepatitis C virus (HCV) RNA quantitation are adequate for clinical management of patients with chronic HCV infection. *J Clin Microbiol.* 2006;44(7):2507-2511. doi:10.1128/JCM.00163-06
12. Kowdley KV, Nelson DR, Lalezari JP, et al. On-treatment HCV RNA as a predictor of sustained virological response in HCV genotype 3-infected patients treated with daclatasvir and sofosbuvir. *Liver Int.* 2016;36(11):1611-1618. doi:10.1111/liv.13165
13. Yoshida EM, Sulkowski MS, Gane EJ, et al. Concordance of sustained virological response 4, 12, and 24 weeks post-treatment with sofosbuvir-containing regimens for hepatitis C virus. *Hepatology.* 2015;61(1):41-45. doi:10.1002/hep.27366
14. WHO Guidelines on Hepatitis B and C Testing. World Health Organization (WHO); 2017. <https://www.who.int/publications/i/item/9789241549981>
15. Saldanha J, Heath A, Aberham C, et al. World Health Organization collaborative study to establish a replacement WHO international standard for hepatitis C virus RNA nucleic acid amplification technology assays. *Vox Sang.* 2005;88(3):202-204. doi:10.1111/j.1423-0410.2005.00606.x

TECHNICAL ASSISTANCE

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

SUMMARY OF SAFETY AND PERFORMANCE STATEMENT

A summary of safety and performance (SSP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSP location after the launch of European Database on Medical Devices. Search for device using UDI-DI provided on the outer packaging of the device.

Abbott RealTime HCV is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck- Ring 2, 65205 Wiesbaden, Germany

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

© 2020, 2024 Abbott. All Rights Reserved.

51-608440/R3


August 2024

APPENDIX 1. OVERVIEW OF THE ABBOTT REALTIME HCV AMPLIFICATION REAGENT EXTENDED USE FEATURE

The amplification reagent extended use feature allows for the use of an amplification reagent pack and internal control (IC) a total of 2 times. Amplification reagent packs that have not yet been used to prepare master mix are referred to as **new** amplification reagent packs. Amplification reagent packs that have been used once and contain prepared master mix are referred to as **partial** amplification reagent packs. Refer to the instructions provided in this manual for additional details.

Storage Conditions (Amplification Reagent Pack and Internal Control Vials)		
Pack Type	Storage Temperature	Storage Time
New Packs	-10°C or colder	Until date shown on label
New IC	-10°C or colder	Until date shown on label
Partial Packs	-10°C or colder (protected from light)	Up to 14 days after initial use
Partial IC	-10°C or colder	Up to 14 days after initial use

- If using UNG for contamination control procedure, refer to the Assay Protocol step 25 and 32 for specific instructions.

- Amplification reagent packs eligible for extended use must have a 6-digit serial number above the barcode.
- Partial** amplification reagent packs can only be used a second time on the same instrument as the initial use. Using them on a different instrument will generate a processing error, which may delay the run.
- Partial** and **new** amplification reagent packs may be used together. All amplification reagent packs used on the instrument for a run must have the same lot number.

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1. Gower E, Estes C, Blach S, Razavi-Shearer K, Razavi H. Global epidemiology and genotype distribution of the hepatitis C virus infection. *J Hepatol*. 2014;61(1 Suppl):S45-S57. doi:10.1016/j.jhep.2014.07.027
2. Guidelines for the screening, care and treatment of persons with chronic hepatitis C infection. World Health Organization (WHO); updated April 2016. <https://apps.who.int/iris/handle/10665/205035>
3. AASLD/IDSA HCV guidance panel. *Recommendations for testing, managing, and treating hepatitis C*. AASLD/IDSA; updated October 24, 2022. <https://www.hcvguidelines.org/>
4. European Association for the Study of the Liver. EASL Recommendations on Treatment of Hepatitis C 2018. *J Hepatol*. 2018;69(2):461-511. doi:10.1016/j.jhep.2018.03.026
5. Miller MH, Agarwal K, Austin A, et al. Review article: 2014 UK consensus guidelines - hepatitis C management and direct-acting anti-viral therapy. *Aliment Pharmacol Ther*. 2014;39(12):1363-1375. doi:10.1111/apt.12764
6. Davis GL, Esteban-Mur R, Rustgi V, et al. Interferon alfa-2b alone or in combination with ribavirin for the treatment of relapse of chronic hepatitis C. International Hepatitis Interventional Therapy Group. *N Engl J Med*. 1998;339(21):1493-1499. doi:10.1056/NEJM199811193392102
7. Manns MP, McHutchison JG, Gordon SC, et al. Peginterferon alfa-2b plus ribavirin compared with interferon alfa-2b plus ribavirin for initial treatment of chronic hepatitis C: a randomised trial. *Lancet*. 2001;358(9286):958-965. doi:10.1016/s0140-6736(01)06102-5
8. Fried MW, Shiffman ML, Reddy KR, et al. Peginterferon alfa-2a plus ribavirin for chronic hepatitis C virus infection. *N Engl J Med*. 2002;347(13):975-982. doi:10.1056/NEJMoa020047
9. Hadziyannis SJ, Sette H Jr, Morgan TR, et al. Peginterferon-alpha2a and ribavirin combination therapy in chronic hepatitis C: a randomized study of treatment duration and ribavirin dose. *Ann Intern Med*. 2004;140(5):346-355. doi:10.7326/0003-4819-140-5-200403020-00010
10. Orlent H, Deltenre P, Francque S, et al. Update of the Belgian Association for the Study of the Liver guidelines for the treatment of chronic hepatitis C genotype 1 with protease inhibitors. *Acta Gastroenterol Belg*. 2012;75(2):245-259.
11. Halfon P, Bourlière M, Pénaranda G, Khiri H, Ouzan D. Real-time PCR assays for hepatitis C virus (HCV) RNA quantitation are adequate for clinical management of patients with chronic HCV infection. *J Clin Microbiol*. 2006;44(7):2507-2511. doi:10.1128/JCM.00163-06
12. Kowdley KV, Nelson DR, Lalezari JP, et al. On-treatment HCV RNA as a predictor of sustained virological response in HCV genotype 3-infected patients treated with daclatasvir and sofosbuvir. *Liver Int*. 2016;36(11):1611-1618. doi:10.1111/liv.13165
13. Yoshida EM, Sulkowski MS, Gane EJ, et al. Concordance of sustained virological response 4, 12, and 24 weeks post-treatment with sofosbuvir-containing regimens for hepatitis C virus. *Hepatology*. 2015;61(1):41-45. doi:10.1002/hep.27366
14. WHO Guidelines on Hepatitis B and C Testing. World Health Organization (WHO); 2017. <https://www.who.int/publications/i/item/9789241549981>
15. Saldanha J, Heath A, Aberham C, et al. World Health Organization collaborative study to establish a replacement WHO international standard for hepatitis C virus RNA nucleic acid amplification technology assays. *Vox Sang*. 2005;88(3):202-204. doi:10.1111/j.1423-0410.2005.00606.x

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1. Gower E, Estes C, Blach S, Razavi-Shearer K, Razavi H. Global epidemiology and genotype distribution of the hepatitis C virus infection. *J Hepatol.* 2014;61(1 Suppl):S45-S57. doi:10.1016/j.jhep.2014.07.027
2. Guidelines for the screening, care and treatment of persons with chronic hepatitis C infection. World Health Organization (WHO); updated April 2016. <https://apps.who.int/iris/handle/10665/205035>
3. AASLD/IDSA HCV guidance panel. *Recommendations for testing, managing, and treating hepatitis C.* AASLD/IDSA; updated October 24, 2022. <https://www.hcvguidelines.org/>
4. European Association for the Study of the Liver. EASL Recommendations on Treatment of Hepatitis C 2018. *J Hepatol.* 2018;69(2):461-511. doi:10.1016/j.jhep.2018.03.026
5. Miller MH, Agarwal K, Austin A, et al. Review article: 2014 UK consensus guidelines - hepatitis C management and direct-acting anti-viral therapy. *Aliment Pharmacol Ther.* 2014;39(12):1363-1375. doi:10.1111/apt.12764
6. Davis GL, Esteban-Mur R, Rustgi V, et al. Interferon alfa-2b alone or in combination with ribavirin for the treatment of relapse of chronic hepatitis C. International Hepatitis Interventional Therapy Group. *N Engl J Med.* 1998;339(21):1493-1499. doi:10.1056/NEJM199811193392102
7. Manns MP, McHutchison JG, Gordon SC, et al. Peginterferon alfa-2b plus ribavirin compared with interferon alfa-2b plus ribavirin for initial treatment of chronic hepatitis C: a randomised trial. *Lancet.* 2001;358(9286):958-965. doi:10.1016/s0140-6736(01)06102-5
8. Fried MW, Shiffman ML, Reddy KR, et al. Peginterferon alfa-2a plus ribavirin for chronic hepatitis C virus infection. *N Engl J Med.* 2002;347(13):975-982. doi:10.1056/NEJMoa020047
9. Hadziyannis SJ, Sette H Jr, Morgan TR, et al. Peginterferon-alpha2a and ribavirin combination therapy in chronic hepatitis C: a randomized study of treatment duration and ribavirin dose. *Ann Intern Med.* 2004;140(5):346-355. doi:10.7326/0003-4819-140-5-200403020-00010
10. Orlent H, Deltenre P, Francque S, et al. Update of the Belgian Association for the Study of the Liver guidelines for the treatment of chronic hepatitis C genotype 1 with protease inhibitors. *Acta Gastroenterol Belg.* 2012;75(2):245-259.
11. Halfon P, Bourlière M, Pénaranda G, Khiri H, Ouzan D. Real-time PCR assays for hepatitis C virus (HCV) RNA quantitation are adequate for clinical management of patients with chronic HCV infection. *J Clin Microbiol.* 2006;44(7):2507-2511. doi:10.1128/JCM.00163-06
12. Kowdley KV, Nelson DR, Lalezari JP, et al. On-treatment HCV RNA as a predictor of sustained virological response in HCV genotype 3-infected patients treated with daclatasvir and sofosbuvir. *Liver Int.* 2016;36(11):1611-1618. doi:10.1111/liv.13165
13. Yoshida EM, Sulkowski MS, Gane EJ, et al. Concordance of sustained virological response 4, 12, and 24 weeks post-treatment with sofosbuvir-containing regimens for hepatitis C virus. *Hepatology.* 2015;61(1):41-45. doi:10.1002/hep.27366
14. WHO Guidelines on Hepatitis B and C Testing. World Health Organization (WHO); 2017. <https://www.who.int/publications/item/9789241549981>
15. Saldanha J, Heath A, Aberham C, et al. World Health Organization collaborative study to establish a replacement WHO international standard for hepatitis C virus RNA nucleic acid amplification technology assays. *Vox Sang.* 2005;88(3):202-204. doi:10.1111/j.1423-0410.2005.00606.x

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
1. Gower E, Estes C, Blach S, Razavi-Shearer K, Razavi H. Global epidemiology and genotype distribution of the hepatitis C virus infection. *J Hepatol*. 2014;61(1 Suppl):S45-S57. doi:10.1016/j.jhep.2014.07.027
2. Guidelines for the screening, care and treatment of persons with chronic hepatitis C infection. World Health Organization (WHO); updated April 2016. <https://apps.who.int/iris/handle/10665/205035>
3. AASLD/IDSA HCV guidance panel. *Recommendations for testing, managing, and treating hepatitis C*. AASLD/IDSA; updated October 24, 2022. <https://www.hcvguidelines.org/>
4. European Association for the Study of the Liver. EASL Recommendations on Treatment of Hepatitis C 2018. *J Hepatol*. 2018;69(2):461-511. doi:10.1016/j.jhep.2018.03.026
5. Miller MH, Agarwal K, Austin A, et al. Review article: 2014 UK consensus guidelines - hepatitis C management and direct-acting anti-viral therapy. *Aliment Pharmacol Ther*. 2014;39(12):1363-1375. doi:10.1111/apt.12764
6. Davis GL, Esteban-Mur R, Rustgi V, et al. Interferon alfa-2b alone or in combination with ribavirin for the treatment of relapse of chronic hepatitis C. International Hepatitis Interventional Therapy Group. *N Engl J Med*. 1998;339(21):1493-1499. doi:10.1056/NEJM199811193392102
7. Manns MP, McHutchison JG, Gordon SC, et al. Peginterferon alfa-2b plus ribavirin compared with interferon alfa-2b plus ribavirin for initial treatment of chronic hepatitis C: a randomised trial. *Lancet*. 2001;358(9286):958-965. doi:10.1016/S0140-6736(01)06102-5
8. Fried MW, Shiffman ML, Reddy KR, et al. Peginterferon alfa-2a plus ribavirin for chronic hepatitis C virus infection. *N Engl J Med*. 2002;347(13):975-982. doi:10.1056/NEJMoa020047
9. Hadziyannis SJ, Sette H Jr, Morgan TR, et al. Peginterferon-alpha2a and ribavirin combination therapy in chronic hepatitis C: a randomized study of treatment duration and ribavirin dose. *Ann Intern Med*. 2004;140(5):346-355. doi:10.7326/0003-4819-140-5-200403020-00010
10. Orlent H, Deltenre P, Francque S, et al. Update of the Belgian Association for the Study of the Liver guidelines for the treatment of chronic hepatitis C genotype 1 with protease inhibitors. *Acta Gastroenterol Belg*. 2012;75(2):245-259.
11. Halfon P, Bourlière M, Pénaranda G, Khiri H, Ouzan D. Real-time PCR assays for hepatitis C virus (HCV) RNA quantitation are adequate for clinical management of patients with chronic HCV infection. *J Clin Microbiol*. 2006;44(7):2507-2511. doi:10.1128/JCM.00163-06
12. Kowdley KV, Nelson DR, Lalezari JP, et al. On-treatment HCV RNA as a predictor of sustained virological response in HCV genotype 3-infected patients treated with daclatasvir and sofosbuvir. *Liver Int*. 2016;36(11):1611-1618. doi:10.1111/liv.13165
13. Yoshida EM, Sulkowski MS, Gane EJ, et al. Concordance of sustained virological response 4, 12, and 24 weeks post-treatment with sofosbuvir-containing regimens for hepatitis C virus. *Hepatology*. 2015;61(1):41-45. doi:10.1002/hep.27366
14. WHO Guidelines on Hepatitis B and C Testing. World Health Organization (WHO); 2017. <https://www.who.int/publications/item/9789241549981>
15. Saldanha J, Heath A, Aberham C, et al. World Health Organization collaborative study to establish a replacement WHO international standard for hepatitis C virus RNA nucleic acid amplification technology assays. *Vox Sang*. 2005;88(3):202-204. doi:10.1111/j.1423-0410.2005.00606.x

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