WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

Product: Alinity m HR HPV
WHO reference number: PQDx 09308-027-00

Alinity m HR HPV with product codes 09N15-090, 09N15-080, and 09N15-092, manufactured by Abbott Molecular Inc., CE-marked regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 27 March 2025.

Summary of WHO's Prequalification Assessment for the Alinity m HR HPV

	Date	Outcome
Prequalification listing	27 March 2025	listed
Abridged Dossier assessment	26 March 2025	MR
Product performance evaluation	06 December 2024	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product, for which the WHO has been notified and has undertaken a review. The amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment and change request reference, where applicable.	Date of report amendment
2.0	Incorporate self-collected vaginal specimens to the Alinity m HR HPV assay. The following changes are necessary to implement the specimen type of self-collected vaginal samples: Creation of a new list number (09N15-092) for the Alinity m HR HPV AMP Kit containing the additional self-collected vaginal sample claim in the Instructions for Use. Creation of an application specification file (ASF) for self-collected vaginal sample specimens (List No. 09N15-07A). Creation of a simpli-COLLECT™ HPV Collection Kit (List No. 09R03-030) and Instructions for Use for both laboratory professionals and patients. Creation of additional reagent kits for pre-analytical processing of self-collected vaginal samples; Alinity m Specimen Buffer (List No. 09R03-050 and 09-30-060). The proposed change does not impact the function, usability, or safety of the Alinity m HR-HPV assay. The change does not include any changes to formulations or	30 January 2026

	specifications of the current assay reagents used or changes to the final release testing of the product. The changes do not affect the manufacturing facility location, principles of operation, or device design. The change is not the result of an adverse event for the Alinity m HR HPV assay. Hardware changes are not required to implement the proposed change (PQC-IVD-2025-0047).	
--	--	--

Intended use

For product list (product codes 09N15-090, 09N15-080)

"The Alinity m High Risk

(HR) HPV assay is a qualitative in vitro test for use with the automated Alinity m System for the detection of DNA from 14 high-risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in clinical specimens. The assay specifically identifies HPV genotypes 16, 18, and 45 while reporting the concurrent detection of the other high-risk genotypes (31/ 33/ 52/ 58) and (35/ 39/ 51/ 56/ 59/ 66/ 68) at clinically relevant infection levels.

The Alinity m HR HPV assay is intended for the following uses:

- To screen patients with ASC-US (atypical squamous cells of undetermined significance) cervical cytology test results to determine the need for referral to colposcopy. The results of this test are not intended to prevent women from proceeding to colposcopy.*
- To be used with cervical cytology to adjunctively screen to assess the presence or absence of high-risk HPV genotype.*
- To be used as a first-line primary screening test to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease.*
- To assess the presence or absence of HPV genotypes 16 and 18 to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease with or without cervical cytology.*

The results from the Alinity m HR HPV, together with the physician's assessment of cytology, history, other risk factors, and professional guidelines, may be used to guide patient management.

INTENDED USER

The intended users for Alinity m HR HPV AMP Kit are laboratory professionals."

For product list (product code 09N15-092)

According to the intended use claim from Abbott Molecular Inc., *"The Alinity m High Risk (HR) HPV assay is a qualitative in vitro test for use with the automated Alinity m System for the detection of DNA from 14 high-risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in cervical specimens and self-collected vaginal specimens. The assay specifically identifies HPV genotypes 16, 18, and 45 while reporting the*

concurrent detection of the other high-risk genotypes (31/ 33/ 52/ 58) and (35/ 39/ 51/ 56/ 59/ 66/68) at clinically relevant infection levels.

The Alinity m HR HPV assay is intended for the following uses:

- *To screen patients with ASC-US (atypical squamous cells of undetermined significance) cervical cytology test results to determine the need for referral to colposcopy. The results of this test are not intended to prevent women from proceeding to colposcopy.*
- *To be used with cervical cytology to adjunctively screen to assess the presence or absence of high-risk HPV genotype.*
- *To be used as a first-line primary screening test to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease.*
- *To assess the presence or absence of HPV genotypes 16 and 18 to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease with or without cervical cytology.*

The results from the Alinity m HR HPV, together with the physician's assessment of cytology, history, other risk factors, and professional guidelines, may be used to guide patient management.

INTENDED USER

The intended users for Alinity m HR HPV AMP Kit are laboratory professionals."

Test kit contents:

Component	Number of tests and product codes
Alinity m HR HPV AMP Kit • Alinity m HR HPV AMP TRAY 1 • Alinity m HR HPV ACT TRAY 2	384 tests, product code 09N15-090. 4 trays / 96 tests each 4 trays / 96 tests each
Alinity m HR HPV CTRL Kit • Alinity m HR HPV Negative CTRL • Alinity m HR HPV Positive CTRL	Product code 09N15-080 12 tubes x 0.6 mL 12 tubes x 0.6 mL
Alinity m HR HPV AMP Kit • Alinity m HR HPV AMP TRAY 1 • Alinity m HR HPV ACT TRAY 2	384 tests, product code 09N15-092 4 trays / 96 tests each 4 trays / 96 tests each

Items required but not provided:

- 09N18-001 Alinity m Sample Prep Kit 1
- 09N20-001 Alinity m Lysis Solution

- 09N20-002 Alinity m Ethanol Solution or 09N20-012 Alinity m Bottle for Ethanol Use (filled with customer supplied 190 Proof Ethanol, ACS, Denaturant free)
- 09N20-003 Alinity m Diluent Solution
- 09N20-004 Alinity m Vapor Barrier Solution
- 09N49-010 Alinity m Transport Tube Pierceable Capped
- 09N49-011 Alinity m Transport Tubes
- 09N49-012 Alinity m Pierceable Caps
- 09N65-050 Alinity m Labeled Tube with Pierceable Cap
- Alinity m HR HPV Application Specification File
- Vortex mixer
- Calibrated pipettes capable of delivering 100 to 1000 µL
- Aerosol barrier pipette tips for 100 to 1000 µL pipettes
- Plate adapter for 384 well plates (such as Eppendorf Catalog No. 022638955)
- Centrifuge with swing plate rotor capable of accommodating the plate adapter and capable of ≥ 100 g,

Storage:

The Alinity m HR HPV AMP Kit (both list numbers 09N15-090 and 09N15-092) must be stored between 2 °C and 8 °C.

The Alinity m HR HPV CTRL Kit must be stored between -25 and -15 °C.

Shelf-life upon manufacture¹:

24 months.

Product dossier assessment

Abbott Molecular Inc. submitted an abridged product dossier for the Alinity m HR HPV as per the “Instructions for compilation of a product dossier” (PQDx_018). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO. The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 26 March 2025.

Based on the product dossier screening and assessment findings, the abridged product dossier for Alinity m HR HPV meets WHO prequalification requirements.

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pqweb/vitro-diagnostics/who-public-inspection-reports>

Based on the site inspection and corrective action plan review, the quality management system for the Alinity m HR HPV meets WHO prequalification requirements.

Product performance evaluation

The Alinity m HR HPV assay was evaluated by the Scottish HPV Reference Laboratory, United Kingdom on behalf of WHO in the 3rd quarter of 2024 according to protocol PQDx_255, version 2.0.

Alinity m HR HPV is considered as one of the possible comparator assays for the virologic evaluation of HPV test for prequalification (per protocol PQDx_255, v2.0). As a result, Alinity m HR HPV was waived from virologic evaluation, and the evaluation was limited to the analytical performance evaluation and operational characteristics.

Analytical performance evaluation

Analytical performance characteristics	
Limit of detection (LoD)	The LoD was estimated at 305.2 IU/mL (95% CI: 215.7-431.8) for HPV Type 16; 287.4 IU/mL (95% CI: 180.4-458.0) for HPV Type 18; and 3537 IU/mL (95% CI: 10^{-5} to 10^{12}) for HPV Type 31.
Reproducibility	The hit rates for detection of HPV 16, HPV 18 and HPV 31 at approx. 10^4 IU/mL were 100%, 100%, and 100%, respectively.
Genotype detection	The following genotypes from the genotype detection panel were detected: HPV16, HPV18, HPV31, HPV33, HPV45, HPV52, HPV58, in agreement with the manufacturer's claim. HPV6 and HPV11 were

	not detected, in agreement with manufacturer's claim.
Cross-contamination / carry-over	No cross-contamination was observed when high positive and negative specimens were tested alternatively.

Operational characteristics and ease of use

This assay requires laboratory equipment and cannot be performed in laboratories with limited facilities or in non-laboratory settings. The instrument requires a stable source of electricity and significant physical space. Furthermore, training and implementation of good laboratory practice is essential to obtaining accurate results. Adequate technical support from manufacturer or representative is critical.

The assay was found easy to use by the operators performing the evaluation. It should be noted that the operators were skilled and experienced in molecular laboratory testing, including on the specific platform under evaluation.

Key operational characteristics	
Specimen type(s) and volume	Cervical specimens collected in Alinity m Cervi-Collect Specimen Collection Kit, ThinPrep PreservCyt Solution, or SurePath Preservative Fluid. The minimal volume required in the tube is 0.55 mL.
Time to result for one test	<120 minutes from test initiation to completion
Operator hands-on time for one test/run	Variable as random-access machine and depending on number of specimens tested in a run. For a full rack of specimens (n = 12), hands-on time is approximately 10 minutes to account for vortexing, scanning and addition of assay orders, and loading of specimens.
Level of automation	Fully automated. Once specimens have been loaded, the Alinity m System performs extraction and amplification of nucleic acid, then displays results when complete.
Quality controls	QCs are provided by the manufacturer but should be purchased separately:

	Alinity m HR HPV CTRL Kit contains 12 HPV negative vials and 12 HPV positive vials which should be run at a minimum frequency of every 48 hours.
Operating temperature	15°C - 28°C
Result display and connectivity	Results are displayed on the instrument and can be interfaced. The results can be exported to the laboratory information system and other health information systems, in addition to CSV Excel exports or single result reports.
Power sources	Main power. The use of a UPS is recommended, as stable electricity is required.
Biosafety (<i>outside of infectious specimen handling</i>)	Operators reported biosafety considerations. Some reagents contain hazardous components, including guanidine hydrochloride and guanidine thiocyanate, that may cause skin or eye irritation, are harmful if swallowed, cause skin burns and eye damage, and are harmful to aquatic life. Furthermore, reagents and waste containing, or contaminated with, guanidine hydrochloride and/or guanidine thiocyanate pose an environmental hazard as the reagents react with sodium hypochlorite and can produce toxic fumes.
Waste	The volume of waste generated by the Alinity m is difficult to quantify as the machine is random-access and was shared between different departments at the time of evaluation. Both solid and liquid waste was generated. Waste disposal requires specific measures in addition to usual laboratory biohazard waste disposal procedures. Waste containing, or contaminated with, guanidine hydrochloride and/or guanidine thiocyanate require specialist chemical uplift and should not be discarded down sinks. This this poses an environmental hazard as the reagents react with sodium hypochlorite and can produce toxic fumes.

Calibration	Calibration of the Alinity m System is performed as part of annual preventative maintenance by an Abbott Molecular Inc. engineer.
Maintenance	Weekly and monthly maintenance are required as a minimum and enforced by the Alinity m System. In addition, yearly preventative maintenance is performed by an Abbott Molecular Inc. engineer.
Other specific requirements	The Alinity m System has a considerable footprint (approximately 248 cm (length) x 102 cm (width) x 188 cm (height)) and requires consistent maintenance and service by Abbott-authorized engineers/service providers.

Based on these results, the performance evaluation for Alinity m HR HPV meets the WHO prequalification requirements.

Labelling review

The labelling submitted for the Alinity m HR HPV was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

Controlled Labelling References

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Alinity m HR HPV AMP Kit (09N15-092)	R1	Jan 2024	53-602370
Outer box artwork	simpli-COLLECT HPV Collection Kit, Case Label	R1	Oct 2024	AL1218.1
Outer box artwork	simpli-COLLECT HPV Collection Kit, Box Label	R1	Oct 2024	AL1217.1
Outer box artwork	Alinity m Specimen Buffer Kit II (09R03-050)	R2	Sep 2024	53-602969
Outer box artwork	Alinity m Specimen Buffer Kit III (09R03-060)	R1	Dec 2023	AL1196.1
Pouch / Device label	Alinity m HR HPV AMP Pouch	R1	June 2024	53-602094
Pouch / Device label	Alinity m HR HPV ACT Pouch	R1	Mar 2018	53-602093
Pouch / Device label	simpli-COLLECT HPV Collection Kit, Outer Pouch (Blister Pack)	R1	Oct 2024	AL1216.1
Reagent bottle labels	Alinity m Specimen Buffer Kit II (09R03-050), Bottle Label	R2	Sep 2024	53-602970
Accessory labeling (e.g., pipettes, buffer caps)	Alinity m Sample Tube for Dry Swab	R1	Nov 2023	AL1202.1
Accessory labeling (e.g., pipettes, buffer caps)	Alinity m Specimen Buffer Kit III, Tube Label	R1	Dec 2023	AL1195.1

Instructions for Use (IFU)	Alinity m HR HPV AMP Kit (09N15-092)	R1	April 2024	53-608488
Instructions for Use (IFU)	simpli-COLLECT HPV Collection Kit, Laboratory Personnel (09R03-030)	R1	November 2024	53-608547
Instructions for Use (IFU)	simpli-COLLECT HPV Collection Kit, Patient (09R03-030)	R2	February 2025	53-608546
Instructions for Use (IFU)	Alinity m Specimen Buffer Kit II (09R03-050)	R2	August 2024	53-608399
Instructions for Use (IFU)	Alinity m Specimen Buffer Kit III (09R03-060)	R1	April 2024	53-608469

1. Labels

Alinity m

Abbott

HR HPV

2

UNIT 1

IVD

Σ 96

ACT TRAY

2°C 8°C

REF 9N15M

53-602093/R1

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

LOT

Colors:

■ Pantone 376 C

■ Pantone 185 C

■ Black

Labeling: Renee

Alinity m

HR HPV

1

Abbott

UNIT 1

IVD

Σ 96

53-602094/R1

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

AMP TRAY

2°C 8°C

REF 9N15P

LOT

Colors:

■ Pantone Reflex Blue C

■ Pantone 185 C

■ Black

Labeling: Renee

Alinity^m

HR HPV AMP Kit  **Abbott**

HR HPV AMP Kit

AMP TRAY 4 x  96

ACT TRAY 4 x  96

REF 09N15-092

H302, H316, H317, H370,
H412, P260, P264, P272,
P273, P280, P301+P312,
P302+P352, P308+P311,
P333+P313, P362+P364,
P501


IVD

In Vitro Test

www.molecular.abbott/portal

REF 09N15-092

53-602370/R1

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

EC REP

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

LOT

Alinity^m

Specimen Buffer Kit II


SPECIMEN BUFFER II 4 x 237 mL

EUH032, P501

IVD

In Vitro Test

For In Vitro Diagnostic Use

53-602969/R2

PN 73012

R2

www.molecular.abbott

REF

09R03-050

EC REP

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

LOT

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Alinity m

Specimen Buffer
Kit III

Specimen Buffer Kit III


REF 09R03-060

IVD

In Vitro Test

For In Vitro Diagnostic Use

REF 09R03-060

SPECIMEN BUFFER III 36 x 2.95 mL

EUH032
P501

R1
www.molecular.abbott/portal

AL1196.1

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

EC REP

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

PRODUCED BY MML Diagnostics Packaging Troutdale, OR 97060 USA

LOT

simpli-COLLECT™ HPV Collection Kit

SAMPLE TUBE 1

SWAB 1

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

PRODUCED BY

MML Diagnostics Packaging, Inc., Troutdale, OR 97060 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

R1

AL1216.1

OTHER

PAP

Abbott

REF

09R03-030

IVD

In Vitro Test

STORE AT

simpli-COLLECT™

HPV Collection Kit

UNIT 40

SAMPLE TUBE 1


SWAB 1

REF 09R03-030

IVD

In Vitro Test

STORE AT

www.molecular.abbott

AL1217.1

LOT

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

EC REP

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

PRODUCED BY MML Diagnostics Packaging, Inc., Troutdale, OR 97060 USA

2. Instructions for use²

²WHO assessed the English versions of the IFU. The manufacturer is responsible for ensuring the correct translation into other languages.

Created April 2024

REF 09N15-092
53-608488/R1

CUSTOMER SERVICE: 1-800-553-7042

CUSTOMER SERVICE INTERNATIONAL: CALL YOUR ABBOTT REPRESENTATIVE

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NOTICE TO USER

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the customer service or technical assistance section of these instructions.

NAME

Alinity m HR HPV AMP Kit

INTENDED USE

The Alinity m High Risk (HR) HPV assay is a qualitative in vitro test for use with the automated Alinity m System for the detection of DNA from 14 high-risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in cervical specimens and self-collected vaginal specimens. The assay specifically identifies HPV genotypes 16, 18, and 45 while reporting the concurrent detection of the other high-risk genotypes (31/ 33/ 52/ 58) and (35/ 39/ 51/ 56/ 59/ 66/ 68) at clinically relevant infection levels.

The Alinity m HR HPV assay is intended for the following uses:

- To screen patients with ASC-US (atypical squamous cells of undetermined significance) cervical cytology test results to determine the need for referral to colposcopy. The results of this test are not intended to prevent women from proceeding to colposcopy.
- To be used with cervical cytology to adjunctively screen to assess the presence or absence of high-risk HPV genotype.
- To be used as a first-line primary screening test to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease.
- To assess the presence or absence of HPV genotypes 16 and 18 to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease with or without cervical cytology.

The results from the Alinity m HR HPV, together with the physician's assessment of cytology, history, other risk factors, and professional guidelines, may be used to guide patient management.

INTENDED USER

The intended users for Alinity m HR HPV AMP Kit are laboratory professionals.

SUMMARY AND EXPLANATION OF THE TEST

HPV is a small, non-enveloped, double-stranded DNA virus (approximately 8,000 base pairs) that replicates in the nucleus of squamous epithelial cells and induces hyperproliferative lesions.¹ HPV infections are among the most common sexually transmitted infections.² Most HPV infections have a benign clinical consequence and are cleared spontaneously.³ However, persistent HPV infection may result in progression to cervical cancer.⁴⁻⁷ More than two hundred different HPV genotypes have been identified,⁸ among which over forty infect mucosal and genital epithelia.⁹ Genital HPV genotypes are generally classified into high risk (HR) and low risk (LR) groups based on their carcinogenic potential. HR HPV genotypes (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) are associated with invasive cervical cancer (squamous cell carcinoma or adenocarcinoma) or its immediate precursor (high-grade squamous intraepithelial lesion, cervical intraepithelial neoplasia, carcinoma in situ, or adenocarcinoma in situ), whereas LR HPV genotypes induce benign lesions and are not associated with cervical cancer.¹⁰⁻¹³ Three of the 14 HR HPV genotypes, 16, 18 and 45, are associated with approximately 70-80% of invasive cervical cancer cases worldwide.¹⁴⁻¹⁶

For adenocarcinoma, which is more difficult to detect, HPV 16, 18, and 45 have been observed in 78-94% of cases.^{14,15,17} Infection by HPV 16 or HPV 18 is associated with higher risk of disease progression compared to other HR HPV genotypes.¹⁸

Compared with cervical screening methods identifying cytological abnormalities, molecular tests that specifically detect the presence of HR HPV DNA can potentially increase sensitivity and cost-effectiveness of cervical cancer screening programs.¹⁹⁻²⁵ Furthermore, HPV DNA tests can be effectively used in triaging patients with equivocal cytology, in post therapeutic follow-up and in monitoring vaccine efficacy.²⁶⁻²⁸

The Alinity m HR HPV assay is a qualitative in vitro test for use with the automated Alinity m System for the detection of DNA from 14 high risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in cervical specimens and self-collected vaginal specimens. The assay specifically identifies HPV genotypes 16, 18, and 45 while reporting the concurrent detection of the other high risk genotypes (31/ 33/ 52/ 58 and 35/ 39/ 51/ 56/ 59/ 66/ 68) at clinically relevant infection levels.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The Alinity m HR HPV assay requires 2 separate assay specific kits:

- Alinity m HR HPV AMP Kit (09N15-092) consisting of 2 types of multi-well assay trays. The amplification trays (AMP Trays) contain lyophilized, unit-dose PCR amplification/detection reagents. The activation trays (ACT Trays) contain liquid unit-dose activation reagent. The intended storage condition for the Alinity m HR HPV AMP Kit is 2°C to 8°C.
- Alinity m HR HPV CTRL Kit (09N15-080) consisting of negative controls and positive controls, each supplied as liquid in single-use tubes. The intended storage condition for the Alinity m HR HPV CTRL Kit is -25°C to -15°C.

The Alinity m HR HPV assay utilizes real-time polymerase chain reaction (PCR) to amplify and detect DNA sequences from 14 HR HPV genotypes (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) and human genomic DNA sequences that have been extracted from cervical and self-collected vaginal specimens.

Cervical specimens collected with the Alinity m Cervi-Collect Specimen Collection Kit or specimens collected in ThinPrep PreservCyt Solution or SurePath Preservative Fluid can be used with the Alinity m HR HPV assay. Cervical specimens collected in ThinPrep PreservCyt Solution or SurePath Preservative Fluid are transferred to an Alinity m Transport Tube for processing on the Alinity m System.

Self-collected vaginal specimens collected in a clinical or non-clinical setting using the simpli-COLLECT™ HPV Collection Kit or the Evalyn® Brush (Rovers Medical Devices) device can be used with the Alinity m HR HPV assay.

Self-collected vaginal specimens collected in a clinical or non-clinical setting using the simpli-COLLECT HPV Collection Kit are prepared using the Alinity m Specimen Buffer Kit II for processing on the Alinity m System. Self-collected vaginal specimens collected in a clinical or non-clinical setting using the Evalyn Brush are transferred to an Alinity m Specimen Buffer III tube from the Alinity m Specimen Buffer Kit III for processing on the Alinity m System.

The steps of the Alinity m HR HPV assay consist of sample preparation, PCR assembly, amplification/detection, and result calculation and reporting. All steps of the Alinity m HR HPV assay procedure are executed automatically by the Alinity m System.

The Alinity m System is designed to be a random access analyzer that can perform the Alinity m HR HPV assay in parallel with other Alinity m assays on the same instrument.

Nucleic acids from specimens are extracted using the Alinity m Sample Prep Kit 1, Alinity m Lysis Solution, Alinity m Ethanol Solution or Alinity m Bottle for Ethanol Use (filled with customer supplied 190 Proof Ethanol, ACS, Denaturant free) and Alinity m Diluent Solution. The Alinity m System employs magnetic microparticle technology to facilitate nucleic acid capture, wash and elution.

The resulting purified nucleic acids are then combined with liquid unit-dose Alinity m HR HPV activation reagent and lyophilized unit-dose Alinity m HR HPV amplification/detection reagents and transferred into a reaction vessel. Alinity m Vapor Barrier Solution is then added to the reaction vessel which is then transferred to an amplification/detection unit for PCR amplification, and real-time fluorescence detection of HR HPV.

The Alinity m HR HPV amplification/detection reagents include primers and probes that amplify and detect an endogenous human beta globin (BG) sequence as sample validity control for cell adequacy, sample extraction and amplification efficiency. The Alinity m HR HPV amplification/detection reagent also contains Uracil-DNA Glycosylase (UDG) as a contamination control for amplicons containing uracil, which may be present in molecular laboratories.

Assay controls are tested at or above an established minimum frequency to help ensure that instrument and reagent performance remains satisfactory. During each control event, a negative control and a positive control are processed through sample preparation and PCR procedures that are identical to those used for specimens.

The possibility of nucleic acid contamination on the Alinity m System is minimized because:

- Aerosol barrier pipette tips are used for all pipetting. The pipette tips are discarded after use.
- PCR amplification and detection is carried out automatically in a sealed reaction vessel.
- Disposal of the reaction vessel is performed automatically by the Alinity m System.

For additional information on system and assay technology, refer to the Alinity m System Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity m HR HPV AMP Kit List No. 09N15-092

The Alinity m HR HPV AMP Kit is comprised of 2 types of multi-well trays: Alinity m HR HPV AMP TRAY 1 and Alinity m HR HPV ACT TRAY 2.

Each Alinity m HR HPV AMP TRAY 1 (individually packed in a foil pouch with a desiccant bag) contains 96 unit-dose lyophilized amplification reagent wells. One reagent well is used per test.

- Amplification reagent wells consist of synthetic oligonucleotides, DNA Polymerase, excipient, dNTPs, Uracil-DNA Glycosylase in a buffered solution.

Each Alinity m HR HPV ACT TRAY 2 (individually packed in a foil pouch without a desiccant bag) contains 96 unit-dose liquid activation reagent wells. One reagent well is used per test.

- Activation reagent wells consist of magnesium chloride, potassium chloride, and tetramethyl ammonium chloride. Preservative: 0.15% ProClin® 950.

	Quantity
	384 tests
Alinity m HR HPV AMP TRAY 1	4 trays / 96 tests each
Alinity m HR HPV ACT TRAY 2	4 trays / 96 tests each

WARNINGS AND PRECAUTIONS

IVD

- For In Vitro Diagnostic Use

Safety Precautions

Human specimens should be handled as if infectious using safe laboratory procedures, such as those outlined in Biosafety in Microbiological and Biomedical Laboratories,²⁹ OSHA Standard on Bloodborne Pathogens,³⁰ CLSI Document M29-A4,³¹ and other appropriate biosafety practices.³² Therefore all human sourced materials should be considered infectious.

These precautions include, but are not limited to, the following:

- Wear gloves when handling specimens or reagents.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% sodium hypochlorite or other suitable disinfectant.²⁹

Decontaminate and dispose of all potentially infectious materials in accordance with local, state and federal regulations.³²

The following warnings and precautions apply to:
Alinity m HR HPV ACT Tray 2.

DANGER	Contains Tetramethylammonium chloride, and methylisothiazolone
H302	Harmful if swallowed.
H316	Causes mild skin irritation. ^a
H317	May cause an allergic skin reaction.
H370	Causes damage to organs.
H412	Harmful to aquatic life with long lasting effects.
Prevention	
P260	Do not breathe mist / vapours / spray.
P264	Wash hands thoroughly after handling.
P272	Contaminated work clothing should not be allowed out of the workplace.
P273	Avoid release to the environment.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P301+P312	IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P302+P352	IF ON SKIN: Wash with plenty of water.
P308+P311	IF exposed or concerned: Call a POISON CENTER/doctor.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

^a Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29CFR1910.1200 (HCS) 2012 have been implemented.

Important information regarding the safe handling, transport and disposal of this product is contained in the Safety Data Sheet.

Safety Data Sheets are available from your Abbott Representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity m System Operations Manual, Section 7 and Section 8.

Reagent Shipment

	Shipment Condition
Alinity m HR HPV AMP Kit	On dry ice

Reagent Storage

In order to minimize damage to foil pouches, it is recommended that the Alinity m HR HPV AMP TRAY 1 (AMP TRAY 1) and Alinity m HR HPV ACT TRAY 2 (ACT TRAY 2) are stored in the original kit packaging. Open the foil pouch for the reagent trays just prior to loading onto the instrument. Onboard storage time begins when reagents are loaded on the Alinity m System.

	Storage Temperature	Maximum Storage Time
Unopened	2°C to 8°C	Until expiration date
Onboard	System Temperature	30 days (not to exceed expiration date)

Reagent Handling

- Do not use reagents that have been damaged.
- Minimize contact with the surface of reagent trays during handling.
- Only load AMP TRAY 1 and ACT TRAY 2 from the same AMP Kit lot on the same Alinity m Assay Tray Carrier. Do not load AMP TRAY 1 and ACT TRAY 2 from different AMP Kit lots on the same Alinity m Assay Tray Carrier.
- The Alinity m System will track the onboard storage time of AMP TRAY 1 and ACT TRAY 2 while on the instrument. The Alinity m System will not allow the use of AMP TRAY 1 and ACT TRAY 2 if the maximum onboard storage time has been exceeded.
- For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity m System Operations Manual, Section 8.

Indications of Reagent Deterioration

- Deterioration of the reagents may be indicated when a control error occurs or control values are repeatedly out of the specified ranges.
- Reagents are shipped on dry ice and are stored at 2°C to 8°C upon arrival. If reagents arrive in a condition contrary to this recommendation or are damaged, immediately contact your Abbott Representative.
- For troubleshooting information, refer to the Alinity m System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity m HR HPV assay application specification file must be installed on the Alinity m System prior to performing the assay.

For detailed information on viewing and editing the customizable assay parameters, refer to the Alinity m System Operations Manual, Section 2. For information on printing assay parameters, refer to the Alinity m System Operations Manual, Section 5.

For a detailed description of system operating instructions, refer to the Alinity m System Operations Manual, Section 5.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Collection

Users must follow the instructions in the Alinity m Cervi-Collect Specimen Collection Kit (List No. 09N29-001) for collecting cervical specimens with the Alinity m Cervi-Collect Specimen Collection Kit.

Users must follow the respective manufacturer's instructions for collecting cervical specimens in ThinPrep PreservCyt Solution (Hologic Inc.) or SurePath Preservative Fluid (BD).

Patients must follow the instructions in the simpli-COLLECT HPV Collection Kit (List No. 09R03) or Evalyn Brush for self-collecting vaginal specimens in a clinical or non-clinical setting.

For simpli-COLLECT HPV specimens, only the orange-shaft swab provided in the kit should be used. An Alinity m Sample Tube from the simpli-COLLECT HPV Collection Kit containing no swab, a different type of swab, multiple swabs, or other foreign material cannot be used with the Alinity m HR HPV assay, and a new specimen should be collected.

Specimen Types

Specimens collected in ThinPrep PreservCyt Solution, SurePath Preservative Fluid, or using Alinity m Cervi-Collect Specimen Collection Kit (List No. 09N29-001), simpli-COLLECT HPV Collection Kit (List No. 09R03), or Evalyn Brush can be tested with the Alinity m HR HPV assay. For ThinPrep PreservCyt and SurePath specimens, the aliquot that is removed either prior to or after cytological processing from the collection vial can be used.

Collection Device	Specimen Type ^a
Alinity m Cervi-Collect Specimen Collection Kit	Cervical Specimens
ThinPrep PreservCyt Solution	Cervical Specimens
SurePath Preservative Fluid	Cervical Specimens
simpli-COLLECT HPV Collection Kit	Vaginal Specimens
Evalyn Brush	Vaginal Specimens

^a The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to use the correct specimen types in the assay.

Specimen Storage

Collection Device	Temperature	Maximum Storage Time from Date of Collection
Alinity m Cervi-Collect Specimen Collection Kit	2°C to 30°C	6 months
	-25°C to -15°C	6 months ^a
ThinPrep PreservCyt Solution	2°C to 30°C	6 months
	-25°C to -15°C	6 months
SurePath Preservative Fluid	15°C to 30°C	2 weeks
	2°C to 8°C	6 months
	-25°C to -15°C	6 months
simpli-COLLECT HPV Collection Kit	2°C to 30°C	30 days (dry) ^b
	2°C to 30°C	30 days (resuspended) ^{b,c}
	-25°C to -15°C	6 months (resuspended) ^{a,c}
Evalyn Brush	2°C to 30°C	30 days (dry) ^b
	2°C to 30°C	30 days (resuspended) ^{b,d}
	-25°C to -15°C	6 months (resuspended) ^{a,d}

^a Avoid more than 4 freeze-thaw cycles.

^b Total storage time for dry and resuspended specimens at 2°C to 30°C must not exceed a combined total of 30 days from date of collection.

^c Resuspend the specimen in Alinity m Specimen Buffer II.

^d Resuspend the specimen in an Alinity m Specimen Buffer III tube pre-filled with buffer.

Specimen Shipping

Ship cervical specimens according to the recommended storage temperature listed in the **Specimen Storage** section. Ship simpli-COLLECT HPV specimens and Evalyn Brush specimens at ambient temperature. Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical, diagnostic, or biological specimens.

Preparation for Analysis

Cervical Specimens

Alinity m Cervi-Collect Specimen Collection Kit specimens arriving in the laboratory containing a cervical brush cannot be tested with the Alinity m HR HPV assay. The collection of a new specimen is required. For ThinPrep PreservCyt specimens, vortex each specimen for 15 to 20 seconds and immediately transfer a volume between 0.55 mL and 2.00 mL to an Alinity m Transport Tube (List No. 09N49-010 or 09N49-011) or Alinity m Labeled Tube with Pierceable Cap (List No. 09N65-050) prior to loading on the Alinity m System. The Alinity m Transport Tube may be capped with an Alinity m Pierceable Cap (List No. 09N49-012).

For SurePath specimens from the original SurePath vial, vortex each specimen for 15 to 20 seconds and immediately transfer 0.5 mL of each specimen to an Alinity m Transport Tube (List No. 09N49-010 or 09N49-011) or Alinity m Labeled Tube with Pierceable Cap (List No. 09N65-050) prior to loading on the Alinity m System. The Alinity m Transport Tube may be capped with an Alinity m Pierceable Cap (List No. 09N49-012).

NOTE: Do not use specimens which appear bloody or have a dark brown color.

NOTE: Do not use Alinity m Aliquot Tube (List No. 09N49-013).

Self-collected Vaginal Specimens

Self-collected vaginal specimens using the simpli-COLLECT HPV Collection Kit:

1. Uncap the Alinity m Sample Tube for Dry Swab and place in a sample rack or a tube holder.
2. Transfer 2.5 mL of Alinity m Specimen Buffer II from the Alinity m Specimen Buffer Kit II (List No. 09R03-050) into the Alinity m Sample Tube for Dry Swab.
3. Replace the pierceable cap on the tube.
NOTE: The swab should remain in the tube while being processed on the Alinity m System.

Self-collected vaginal specimens using the Evalyn Brush:

1. Uncap an Alinity m Specimen Buffer Kit III tube. Place uncapped tube in a sample rack or a tube holder to avoid spilling buffer out of the tube; place the pierceable cap on a clean table liner to avoid contamination.
2. Hold the Evalyn Brush device and remove the pink cap by squeezing the cap sides.
3. Push the pink plunger to expose the brush head and the base of the brush head. Do not touch the Evalyn Brush white bristles.
4. If using Evalyn Tweezers to detach the brush head, position the tweezers around the base of the brush head. If using a ReMover to transfer the brush head, guide the Evalyn Brush into the ReMover until the casing hits the ReMover.
NOTE: A new pair of Evalyn Tweezers or a new ReMover must be used for each brush.
5. Squeeze the Evalyn Tweezers or the ReMover firmly together and pull the brush head off. Keep the Evalyn Tweezers or ReMover compressed.
6. Hold the Evalyn Tweezers or ReMover over the opening of the Alinity m Specimen Buffer Kit III tube and release the brush head into the tube with the bristles oriented up.
7. Replace the pierceable cap on the tube.

PROCEDURE

Materials Provided

09N15-092 Alinity m HR HPV AMP Kit

Materials Required but not Provided

- 09N15-080 Alinity m HR HPV CTRL Kit
- 09N18-001 Alinity m Sample Prep Kit 1
- 09N20-001 Alinity m Lysis Solution
- 09N20-002 Alinity m Ethanol Solution or 09N20-012 Alinity m Bottle for Ethanol Use (filled with customer supplied 190 Proof Ethanol, ACS, Denaturant free)
- 09N20-003 Alinity m Diluent Solution
- 09N20-004 Alinity m Vapor Barrier Solution
- 09N49-010 Alinity m Transport Tube Pierceable Capped
- 09N49-011 Alinity m Transport Tubes
- 09N49-012 Alinity m Pierceable Caps
- 09N65-050 Alinity m Labeled Tube with Pierceable Cap
- 09R03-050 Alinity m Specimen Buffer Kit II
- 09R03-060 Alinity m Specimen Buffer Kit III
- Alinity m HR HPV Application Specification File
- Vortex mixer
- Calibrated pipettes capable of delivering 100 to 1000 µL
- Aerosol barrier pipette tips for 100 to 1000 µL pipettes
- Plate adapter for 384 well plates (such as Eppendorf Catalog No. 022638955)
- Centrifuge with swing plate rotor capable of accommodating the plate adapter and capable of ≥ 100g

For information on materials required for operation of the instrument, refer to the Alinity m System Operations Manual, Section 1.

For general operating procedures, refer to the Alinity m System Operations Manual, Section 5.

For optimal performance, it is important to perform routine maintenance as described in the Alinity m System Operations Manual, Section 9.

Procedural Precautions

- Read the instructions in this package insert carefully before processing samples.
- Use aerosol barrier pipette tips or disposable pipettes only one time when pipetting specimens. To prevent contamination to the pipette barrel while pipetting, care should be taken to avoid touching the pipette barrel to the inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.
- Do not use specimens if the specimen tube is damaged or if buffer has leaked from the specimen tube. Discard unused, damaged or leaking specimen tubes in accordance with local, state, and federal regulations.
- Work area and instrument platforms must be considered potential sources of contamination.
- Ensure the Alinity m HR HPV AMP TRAY 1 is tapped prior to loading on the Alinity m System per instructions in the **Assay Procedure** section.
- Ensure the Alinity m HR HPV ACT TRAY 2 is centrifuged prior to loading on the Alinity m System per instructions in the **Assay Procedure** section.
- Monitoring procedures for the presence of amplification product can be found in the Alinity m System Operations Manual, Section 9.
- To reduce the risk of nucleic acid contamination, clean and disinfect spills of specimens by including the use of a tuberculocidal disinfectant such as 1.0% (v/v) sodium hypochlorite or other suitable disinfectant.
- To prevent contamination, change to new gloves before handling the Alinity m Sample Prep Kit 1, assay trays, system solutions, Integrated Reaction Unit (IRU) sleeves, and pipette tips. Also change to new gloves whenever they are contaminated by a specimen, a control, or a reagent. Always use powder-free gloves.
- The use of the Alinity m HR HPV CTRL Kit is integral to the performance of the Alinity m HR HPV assay. Refer to the **QUALITY CONTROL PROCEDURES** section of this package insert for details. Refer to the Alinity m HR HPV CTRL Kit package insert for preparation and usage.
- The Alinity m HR HPV control reagents are contained in single-use tubes with pierceable caps. Avoid contamination or damage to the caps after removal from their original packaging. Discard tubes after use.

Assay Procedure

Prior to loading on the Alinity m System, hold the AMP TRAY 1 by the edges with the label facing up and tap three times on the bench.

Prior to loading on the Alinity m System, ACT TRAY 2 must be centrifuged as follows:

1. Load ACT TRAY 2 onto the plate adapter (Eppendorf Catalog No. 022638955).
2. Load the plate adapter (with ACT TRAY 2) on a swing plate centrifuge capable of accommodating the plate adapter. Spin at 100 to 800g for 1 to 5 minutes to remove potential bubbles.
3. Immediately following centrifugation, carefully transfer the ACT TRAY 2 to the Alinity m Assay Tray Carriers. Take care to minimize disturbance to the ACT TRAY 2. Load the tray carriers per the Alinity m System Operations Manual, Section 5.
4. If disturbance occurs during transfer that could potentially introduce bubbles (eg, dropping, bumping, inversion of the ACT TRAY 2), re-centrifuge the ACT TRAY 2.
5. Proceed with the **Reagent and sample inventory management** procedure per the Alinity m System Operations Manual, Section 5.

For a detailed description of how to run an assay, refer to the Alinity m System Operations Manual, Section 5.

Prior to testing specimens, check control status. If control testing is required, refer to **QUALITY CONTROL PROCEDURES** section. Controls may be tested separately or with specimens.

For preparation of samples, refer to the instructions under **SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS: Preparation for Analysis** section. Vortex specimens for 15 to 20 seconds prior to loading on the Alinity m System. If the specimen is stored frozen, it must be completely thawed prior to sample preparation.

From the Create Order screen, select the assay (HR HPV 07A) being tested.

The Alinity m System will track the onboard storage time of amplification reagents, controls, and specimens while on the instrument. The Alinity m System will not allow the use of amplification reagents, controls, or process specimens that have exceeded the allowable onboard storage time. Specimen tubes need to meet the requirements for minimum sample volume and use of caps when loaded on the Alinity m System. Specimen tubes may be placed on the Alinity m Universal Sample Rack (sample rack) onboard the system for up to 4 hours prior to processing. Specimens collected with the simpli-COLLECT HPV Collection Kit and Evalyn Brush must be loaded onto the sample rack capped and with the retention bar on. Clean the retention bar after each use.

Tube Type ^a	List No.	Minimum Specimen Volume	Cap Requirement on Instrument
Alinity m Cervi-Collect Specimen Collection Transport Tube	09N29-001	0.55 mL	Capped or Uncapped
Alinity m Transport Tube Pierceable Capped	09N49-010	ThinPrep: 0.55 mL SurePath: 0.50 mL	Capped or Uncapped
Alinity m Transport Tube Alinity m Pierceable Cap	09N49-011 09N49-012	ThinPrep: 0.55 mL SurePath: 0.50 mL	Capped or Uncapped
Alinity m Labeled Tube with Pierceable Cap	09N65-050	ThinPrep: 0.55 mL SurePath: 0.50 mL	Capped or Uncapped
Alinity m Sample Tube for Dry Swab (for simpli-COLLECT HPV Collection Kit)	AB887	Fill volume allows initial test and one retest	Capped
Alinity m Specimen Buffer III (for Evalyn Brush)	09R03-060	Fill volume allows initial test and one retest	Capped

^a Do not use Alinity m Aliquot Tube (List No. 09N49-013)

Post Processing Procedure

Upon completion of sample processing, specimens collected with the Alinity m Cervi-Collect Specimen Collection Kit and ThinPrep PreservCyt and SurePath specimens that have been transferred in Alinity m Transport Tubes can be recapped using new, unused Alinity m Pierceable Caps (List No. 09N49-012). Store specimens according to the **Specimen Storage** section.

Upon completion of sample processing, specimens collected with the simpli-COLLECT HPV Collection Kit and Evalyn Brush can be recapped using new, unused Alinity m Pierceable Caps (List No. 09N49-012). Store specimens according to the **Specimen Storage** section.

QUALITY CONTROL PROCEDURES

Negative and Positive Controls

One Alinity m HR HPV Negative CTRL and one Alinity m HR HPV Positive CTRL are recommended to be tested, at or above the minimum frequency of once every 48 hours, to monitor the performance of the assay and Alinity m System. Valid results for all control levels must be obtained before specimen results are reported. Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

If quality control results do not meet the acceptance criteria, refer to the Alinity m System Operations Manual, Section 10, for troubleshooting information. A flag is displayed for specimens when a control result is invalid. All of the specimens processed following an invalid assay control must be retested.

If control results are invalid, refer to the Alinity m System Operations Manual, Section 5 for a description of quality control flags, and Section 10 for troubleshooting information.

The presence of HPV must not be detected in the negative control. HPV detected in the negative control is indicative of contamination by other samples or by amplified product. To avoid contamination, clean the Alinity m System and repeat sample processing for controls and specimens following the Procedural Precautions in this package insert. Monitoring procedures for the presence of amplification product can be found in the Alinity m System Operations Manual, Section 9.

If negative controls are persistently reactive, contact your Abbott Representative.

Detection of Inhibition and/or Cell Inadequacy

The Alinity m HR HPV assay detects the endogenous human beta globin (BG) sequence as Cellular Control (CC) signal to evaluate cell adequacy, sample extraction and amplification efficiency. The CC cycle number (CN) is assessed against a defined range to determine specimen validity.

A Flag or Message Code is displayed when CC CN value of a specimen is out of the established range.

- If the HPV signal(s) in a specimen is detected and CC CN is out of range, the specimen will yield an HR HPV Detected interpretation. A Flag will be displayed.
- If the HPV signal(s) in a specimen is not detected and CC CN is out of range, no HPV result will be reported and a Message Code will be displayed.

Refer to the Alinity m System Operations Manual, Section 5 for an explanation of the corrective actions for Flags.

Refer to the Alinity m System Operations Manual, Section 10 for an explanation of the corrective actions for Message Codes.

RESULTS

Results and interpretations for the Alinity m HR HPV assay are determined based on the comparison of specimen's cycle number (CN) values in each HPV signal against signal-specific, established cutoff values. Five HPV signals measured in separate fluorescence channels corresponding to HPV 16, HPV 18, HPV 45, HPV 31/33/52/58 (pooled signal: a single signal corresponding to any individual genotype or combination of genotypes within the group) and HPV 35/39/51/56/59/66/68 (pooled signal) are evaluated for each sample. Each signal is either determined as "Detected" if the CN is less than or equal to a fixed assay cutoff cycle for that signal or is determined as "Not Detected" if either the CN is not generated or the CN is greater than the assay cutoff cycle. All the detected signals (HPV 16, HPV 18, HPV 45, Other HR HPV A or Other HR HPV B) are reported in the sample result. The detection of Other HR HPV A is reported when HPV 31/33/52/58 signal is detected. The detection of Other HR HPV B is reported when HPV 35/39/51/56/59/66/68 signal is detected. Samples with any of the HR HPV signals detected have an interpretation of "HR HPV Detected." Samples with all HR HPV signals not detected will have an interpretation of "Not Detected."

Interpretation of Results

Results	Interpretation	Flags
Not Detected	Not Detected	
HPV 16	HR HPV Detected	
HPV 18	HR HPV Detected	
HPV 45	HR HPV Detected	
Other HR HPV A	HR HPV Detected	
Other HR HPV B	HR HPV Detected	
HPV16; HPV18	HR HPV Detected	
HPV16; HPV45	HR HPV Detected	
HPV16; Other HR HPV A	HR HPV Detected	
HPV16; Other HR HPV B	HR HPV Detected	
HPV18; HPV 45	HR HPV Detected	
HPV18; Other HR HPV A	HR HPV Detected	
HPV18; Other HR HPV B	HR HPV Detected	
HPV45; Other HR HPV A	HR HPV Detected	
HPV45; Other HR HPV B	HR HPV Detected	
Other HR HPV A; Other HR HPV B	HR HPV Detected	
HPV16; HPV18; HPV45	HR HPV Detected	
HPV16; HPV18; Other HR HPV A	HR HPV Detected	
HPV16; HPV18; Other HR HPV B	HR HPV Detected	
HPV16; HPV45; Other HR HPV A	HR HPV Detected	
HPV16; HPV45; Other HR HPV B	HR HPV Detected	
HPV18; Other HR HPV A; Other HR HPV B	HR HPV Detected	
HPV45; Other HR HPV A; Other HR HPV B	HR HPV Detected	
HPV16; HPV18; HPV45; Other HR HPV A	HR HPV Detected	
HPV16; HPV18; HPV45; Other HR HPV B	HR HPV Detected	
HPV16; HPV18; Other HR HPV A; Other HR HPV B	HR HPV Detected	
HPV16; HPV45; Other HR HPV A; Other HR HPV B	HR HPV Detected	
HPV18; HPV45; Other HR HPV A; Other HR HPV B	HR HPV Detected	
HPV16; HPV18; HPV45; Other HR HPV A; Other HR HPV B	HR HPV Detected	

Flags, Results Codes, and Message Codes

Some results may contain information in the Flags and Codes fields. For a description of the flags and result codes that may appear in these fields, refer to the Alinity m System Operations Manual, Section 5.

For a description of message codes, refer to the Alinity m System Operations Manual, Section 10.

LIMITATIONS OF THE PROCEDURE

- Optimal performance of this test requires appropriate specimen collection and handling (refer to the **SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS** section of this package insert).
- Clinician-collected cervical specimens collected with the Alinity m Cervi-Collect Specimen Collection Kit, collected in ThinPrep PreservCyt Solution, or SurePath Preservative Fluid may be used with the Alinity m HR HPV assay. Self-collected vaginal specimens collected with the simpli-COLLECT HPV Collection Kit and Evalyn Brush in a clinical or non-clinical setting may be used with the Alinity m HR HPV assay. Performance with other specimen types or collection devices has not been evaluated.
- The instruments and assay procedures reduce the risk of contamination by amplification product. However, nucleic acid contamination from the positive controls or specimens must be controlled by good laboratory practice and careful adherence to the procedures specified in this package insert.
- If the HPV results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- As with any diagnostic test, results from the Alinity m HR HPV assay should be interpreted in conjunction with other clinical and laboratory findings.

SPECIFIC PERFORMANCE CHARACTERISTICS

Clinical Sensitivity and Specificity in Screening Population

The clinical sensitivity and specificity for detection of $\geq$ CIN2 among a screening population (age $\geq$ 30 years) were determined for Alinity m HR HPV in comparison with hc2 High-Risk HPV DNA Test (hc2) and Abbott RealTime High Risk (HR) HPV (Table 1). All specimens were collected in ThinPrep PreservCyt Solution.

Assay	Sensitivity (%)		Specificity (%)	
	Estimate (95% CI)	n/N ^a	Estimate (95% CI)	n/N ^b
Alinity m HR HPV	100.0 (94.7, 100.0)	68/68	93.2 (92.2, 94.0)	2768/2971
hc2	95.6 (87.8, 98.5)	65/68	92.9 (91.9, 93.7)	2759/2971
Abbott RealTime HR HPV	100.0 (94.7, 100.0)	68/68	94.2 (93.3, 95.0)	2794/2967

^a n represents the number of $\geq$ CIN2 specimens detected. N represents the number of $\geq$ CIN2 specimens tested.

^b n represents the number of disease negative specimens that were not detected. N represents the number of disease negative specimens that were tested. Negative disease status is defined as either $<$ CIN2 histology or negative cytology when histology result is unknown.

Among the 68 $\geq$ CIN2 specimens, 34 were $\geq$ CIN3. The clinical sensitivity for detection of $\geq$ CIN3 was 100.0% (34/34; 95% CI 89.8% to 100.0%) for Alinity m HR HPV, 97.1% (33/34; 95% CI 85.1% to 99.5%) for hc2, and 100.0% (34/34; 95% CI 89.8% to 100.0%) for Abbott RealTime HR HPV.

Clinical Sensitivity in ASC-US Population

The clinical sensitivity for detection of $\geq$ CIN2 among patients with ASC-US cytology was determined for Alinity m HR HPV in comparison with hc2 and Abbott RealTime HR HPV (Table 2). All specimens were collected in ThinPrep PreservCyt Solution.

Assay	Estimate (95% CI)	n/N ^a
Alinity m HR HPV	96.8 (83.8, 99.4)	30/31
hc2	93.5 (79.3, 98.2)	29/31
Abbott RealTime HR HPV	96.8 (83.8, 99.4)	30/31

^a n represents the number of $\geq$ CIN2 specimens detected. N represents the number of $\geq$ CIN2 specimens tested.

Clinical Specificity in Screening Population with Normal Cytology

The clinical specificity in a screening population (age $\geq$ 30 years) with normal cytology was determined for Alinity m HR HPV in comparison with hc2 and Abbott RealTime HR HPV (Table 3). All specimens were collected in ThinPrep PreservCyt Solution.

Assay	Estimate (95% CI)	n/N ^a
Alinity m HR HPV	92.8 (91.8, 93.7)	2762/2976
hc2	92.5 (91.5, 93.4)	2753/2976
Abbott RealTime HR HPV	93.8 (92.9, 94.6)	2788/2972

^a n represents the number of specimens that were not detected. N represents the number of specimens tested.

Accuracy in Identification of HPV 16 and/or HPV 18 in Women with Cervical Disease

The performance of Alinity m HR HPV in identification of HPV 16 and/or HPV 18 in $\geq$ CIN2 was evaluated based on the results from a screening population (age $\geq$ 30 years) and an ASC-US population (Table 4). Out of 68 $\geq$ CIN2 specimens from the screening population, 68 had an interpretation of "HR HPV Detected" from both Alinity m HR HPV and Abbott RealTime HR HPV. Thirty-five specimens were detected as HPV 16 and/or HPV 18 by both assays. Thirty-three specimens were detected as non-HPV 16/18 by both assays. The overall agreement for detection of HPV 16 and/or HPV 18 between Alinity m HR HPV and Abbott RealTime HR HPV was 100.0% (68/68).

Table 4. Genotyping Accuracy for HPV 16 and/or HPV 18 in Screening Population

Alinity m HR HPV	Abbott RealTime HR HPV	
	HPV 16 and/or HPV 18 Detected ^a	Non-HPV 16/18 High Risk Genotype(s) Detected ^b
HPV 16 and/or HPV 18 Detected ^a	35	0
Non-HPV 16/18 High Risk Genotype(s) Detected ^b	0	33

^a These specimens were detected for HPV 16 and/or HPV 18 signal(s) with or without non-HPV 16/18 HR HPV signal detected.

^b These specimens were not detected for HPV 16 or HPV 18 signal and detected for non-HPV 16/18 HR HPV signal.

Out of 31 $\geq$ CIN2 specimens from the ASC-US population, 30 had an interpretation of "HR HPV Detected" from both Alinity m HR HPV and Abbott RealTime HR HPV (Table 5). Twenty-two specimens were detected as HPV 16 and/or HPV 18 by both assays. Eight specimens were detected as non-HPV 16/18 by both assays. The overall agreement for detection of HPV 16 and/or HPV 18 between Alinity m HR HPV and Abbott RealTime HR HPV was 100.0% (30/30).

Table 5. Genotyping Accuracy for HPV 16 and/or HPV 18 in ASC-US Population

Alinity m HR HPV	Abbott RealTime HR HPV	
	HPV 16 and/or HPV 18 Detected ^a	Non-HPV 16/18 High Risk Genotype(s) Detected ^b
HPV 16 and/or HPV 18 Detected ^a	22	0
Non-HPV 16/18 High Risk Genotype(s) Detected ^b	0	8

^a These specimens were detected for HPV 16 and/or HPV 18 signal(s) with or without non-HPV 16/18 HR HPV signal detected.

^b These specimens were not detected for HPV 16 or HPV 18 signal and detected for non-HPV 16/18 HR HPV signal.

Estimate of Relative Disease Risk Associated with Different Genotype Results

The relative risks of having cervical disease ($\geq$ CIN2) were estimated, by calculating the ratio of absolute risks, for HPV 16 and/or HPV 18 Detected vs. Non-HPV 16/18 HR HPV Detected results in a screening population (age $\geq$ 30 years) and an ASC-US Population (Table 6). The relative risk was 2.3 in screening population, and 2.1 in ASC-US population.

Table 6. Relative Risk of Cervical Disease Associated with Different Genotype Results (HPV 16 and/or HPV 18 Detected vs Non-HPV 16/18 HR HPV Detected)

Population	Relative Risk	95% CI
Screening	2.3	(1.6, 3.5)
ASC-US	2.1	(1.1, 4.1)

Genotype Inclusivity and Partial Genotyping

The ability of Alinity m HR HPV to detect 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) and to specifically identify HPV genotypes 16, 18 and 45 while reporting the concurrent detection of the other high-risk genotypes (Other HR HPV A: 31/33/52/58; Other HR HPV B: 35/39/51/56/59/66/68) was evaluated. Forty samples containing HPV DNA targets from 14 genotypes, individually and in combinations, were tested as listed in Table 7. Results from all samples including 14 with single genotype, 10 with 2 genotypes, 10 with 3 genotypes, 5 with 4 genotypes and 1 with 5 genotypes were reported accurately. The presence or absence of HPV DNA for each designated signal, HPV 16, HPV 18, HPV 45, Other HR HPV A and Other HR HPV B, was accurately determined in each case.

Table 7. Genotype Detection and Partial Genotyping Capability		
Sample No.	Genotype(s)	Reported Result
1	HPV 16	HPV 16
2	HPV 18	HPV 18
3	HPV 31	Other HR HPV A
4	HPV 33	Other HR HPV A
5	HPV 35	Other HR HPV B
6	HPV 39	Other HR HPV B
7	HPV 45	HPV 45
8	HPV 51	Other HR HPV B
9	HPV 52	Other HR HPV A
10	HPV 56	Other HR HPV B
11	HPV 58	Other HR HPV A
12	HPV 59	Other HR HPV B
13	HPV 66	Other HR HPV B
14	HPV 68	Other HR HPV B
15	HPV 16 and HPV 18	HPV 16; HPV 18
16	HPV 16 and HPV 45	HPV 16; HPV 45
17	HPV 16 and HPV 58	HPV 16; Other HR HPV A
18	HPV 16 and HPV 39	HPV 16; Other HR HPV B
19	HPV 18 and HPV 45	HPV 18; HPV 45
20	HPV 18 and HPV 58	HPV 18; Other HR HPV A
21	HPV 18 and HPV 39	HPV 18; Other HR HPV B
22	HPV 45 and HPV 58	HPV 45; Other HR HPV A
23	HPV 45 and HPV 39	HPV 45; Other HR HPV B
24	HPV 58 and HPV 39	Other HR HPV A; Other HR HPV B
25	HPV 16, HPV 18, and HPV 45	HPV 16; HPV 18; HPV 45
26	HPV 16, HPV 18, and HPV 58	HPV 16; HPV 18; Other HR HPV A
27	HPV 16, HPV 18, and HPV 39	HPV 16; HPV 18; Other HR HPV B
28	HPV 16, HPV 45, and HPV 58	HPV 16; HPV 45; Other HR HPV A
29	HPV 16, HPV 45, and HPV 39	HPV 16; HPV 45; Other HR HPV B
30	HPV 16, HPV 58, and HPV 39	HPV 16; Other HR HPV A; Other HR HPV B
31	HPV 18, HPV 45, and HPV 58	HPV 18; HPV 45; Other HR HPV A
32	HPV 18, HPV 45, and HPV 39	HPV 18; HPV 45; Other HR HPV B
33	HPV 18, HPV 58, and HPV 39	HPV 18; Other HR HPV A; Other HR HPV B
34	HPV 45, HPV 58, and HPV 39	HPV 45; Other HR HPV A; Other HR HPV B
35	HPV 16, HPV 18, HPV 45, and HPV 58	HPV 16; HPV 18; HPV 45; Other HR HPV A
36	HPV 16, HPV 18, HPV 45, and HPV 39	HPV 16; HPV 18; HPV 45; Other HR HPV B
37	HPV 16, HPV 18, HPV 58, and HPV 39	HPV 16; HPV 18; Other HR HPV A; Other HR HPV B
38	HPV 16, HPV 45, HPV 58, and HPV 39	HPV 16; HPV 45; Other HR HPV A; Other HR HPV B
39	HPV 18, HPV 45, HPV 58, and HPV 39	HPV 18; HPV 45; Other HR HPV A; Other HR HPV B
40	HPV 16, HPV 18, HPV 45, HPV 58, and HPV 39	HPV 16; HPV 18; HPV 45; Other HR HPV A; Other HR HPV B

Limit of Detection for High Risk HPV Genotypes

Limit of detection (LoD) of Alinity m HR HPV was determined by testing plasmids for 14 HR HPV genotype sequences (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) in ThinPrep PreservCyt Solution containing an HPV negative human cell line (C33A) background. In addition, LoD was determined by testing HPV positive cell lines, SiHa (HPV 16) and HeLa (HPV18), in ThinPrep PreservCyt Solution containing C33A cell line background. Four hundred microliters of sample is used per assay. Each plasmid or cell line was diluted to 4 concentrations, including at, below and above the expected LoD levels, and was tested with 2 lots of amplification reagents. For each plasmid or cell line, a total of 72 replicates were tested across all concentrations (ie, 4 concentrations, 6 replicates per day over 3 days per concentration) with each amplification reagent lot. The LoD was defined as a concentration having a ≥95% detection rate with all higher concentrations having a ≥95% detection rate. The LoD was 40 cells per assay for SiHa (HPV 16) and HeLa (HPV 18); 240 copies per assay for HPV 16 and 18; 500 copies per assay for HPV 45; 2000 copies per assay for HPV 33, 35, 51, 52, and 59; and 5000 copies per assay for HPV 31, 39, 56, 58, 66, and 68.

To evaluate the LoD for vaginal collection devices, plasmids for the 14 HR HPV genotype sequences and HPV positive cell lines SiHa (HPV 16) and HeLa (HPV 18) were prepared at LoD in Specimen Buffer containing an HPV negative human cell line background. Testing was performed using 1 lot of amplification reagent.

100% detection was observed at 40 cells per assay for SiHa (HPV16) and HeLa (HPV18) cells, 240 copies per assay for HPV 16 and 18, 500 copies per assay for HPV 45, 2000 copies per assay for HPV 33, 35, 51, and 59, 4000 copies per assay for HPV 52, and 5000 copies per assay for HPV 31, 39, 56, 58, 66, and 68.

Reproducibility

A total of 581 specimens were tested to evaluate the reproducibility of Alinity m HR HPV within the same laboratory. Each of the specimens was tested twice by two operators in the same laboratory. The overall percent agreement between the two tests (**Table 8**) was 97.9% (569/581; 95% CI: 96.4% to 98.8%), with average positive agreement of 96.9% (370/382; 95% CI: 94.9% to 98.5%) and average negative agreement of 98.5% (768/780; 95% CI: 97.5% to 99.2%).

Table 8. Intra-laboratory Reproducibility

First Test	Second Test	
	HR HPV Detected	Not Detected
	185	5
Not Detected	7	384

A total of 560 specimens were tested to evaluate the reproducibility of Alinity m HR HPV between two laboratories at Abbott using different Alinity m Systems. The overall percent agreement between the two tests (**Table 9**) was 97.5% (546/560; 95% CI: 95.8% to 98.5%), with average positive agreement of 96.1% (342/356; 95% CI: 93.7% to 97.9%) and average negative agreement of 98.2% (750/764; 95% CI: 97.1% to 99.0%).

Table 9. Inter-laboratory Reproducibility

First Laboratory	Second Laboratory	
	HR HPV Detected	Not Detected
	171	9
Not Detected	5	375

Analytical Specificity – Potential Cross-Reactants

The analytical specificity of Alinity m HR HPV was evaluated with a panel of microorganisms (**Table 10**) in HPV negative samples and HR HPV positive samples (containing HPV at 3 times the limit of detection). The panel included low-risk HPV, bacteria, viruses, protozoan, yeast, and human cellular DNA. No cross-reactivity or interference in Alinity m HR HPV performance was observed in the presence of the tested microorganisms.

Table 10. Microorganisms Tested

Low Risk HPV	Bacteria	Virus
HPV 6	<i>Bacteroides fragilis</i>	Adenovirus
HPV 11	<i>Bacteroides ureolyticus</i>	Cytomegalovirus (CMV)
HPV 13	<i>Bifidobacterium adolescentis</i>	Epstein Barr Virus (EBV)
HPV 26	<i>Chlamydia trachomatis</i>	Hepatitis B virus (HBV)
HPV 30	<i>Clostridium perfringens</i>	Herpes simplex virus I
HPV 32	<i>Corynebacterium genitalium</i>	Herpes simplex virus II
HPV 34	<i>Enterococcus faecalis</i>	Human Herpes virus 3
HPV 40	<i>Enterobacter cloacae</i>	Hepatitis C virus (HCV)
HPV 42	<i>Escherichia coli</i>	Human immunodeficiency virus (HIV-1)
HPV 43	<i>Fusobacterium necrophorum</i>	
HPV 44	<i>Gardnerella vaginalis</i>	
HPV 53	<i>Haemophilus ducreyi</i>	Protozoan
HPV 54	<i>Klebsiella pneumoniae ss ozaenae</i>	Trichomonas vaginalis
HPV 55B	<i>Lactobacillus acidophilus</i>	
HPV 57	<i>Mycoplasma genitalium</i>	Yeast
HPV 61	<i>Mycoplasma hominis</i>	
HPV 67	<i>Neisseria gonorrhoeae</i>	Candida albicans
HPV 69	<i>Neisseria meningitidis Serogroup A</i>	
HPV 70	<i>Peptostreptococcus anaerobius</i>	Other
HPV 73	<i>Proteus mirabilis</i>	Human Cellular DNA
HPV 82	<i>Pseudomonas aeruginosa</i>	
HPV 85	<i>Staphylococcus aureus</i>	
	<i>Staphylococcus epidermidis</i>	
	<i>Streptococcus agalactiae</i>	
	<i>Streptococcus pneumoniae</i>	
	<i>Streptococcus pyogenes</i>	
	<i>Treponema pallidum</i>	
	<i>Ureaplasma urealyticum</i>	

Analytical Specificity – Potentially Interfering Substances

The effect of potentially interfering endogenous and exogenous substances that may be present in cervical specimens on Alinity m HR HPV performance was assessed by testing HPV negative samples and HPV positive samples (SiHa and HeLa cell lines at 3 times the limit of detection) in ThinPrep PreservCyt Solution containing C33A cell line background. Blood was evaluated at a concentration of 10% (v/v), mucus at 5% (v/v), and peripheral blood mononuclear cells (PBMC) at approximately 1×10^6 cells/mL. The exogenous interference substances were evaluated at 0.5% (w/v) by spiking into the sample. No interference in Alinity m HR HPV performance was observed in the presence of the tested substances (Table 11).

Table 11. Substances Tested

Mucus	Up & Up Lubricating Liquid
PBMC	KY Jelly
Whole Blood	VCF Contraceptive Gel
Norforms Deodorant Suppositories	Conceptrol Vaginal Contraceptive Gel
Clotrimazole Vaginal Cream	Hydrocortisone
Terazol-3 Vaginal Cream	Zovirax
Monistat 3	Up & Up Povidone - Iodine
Metrogel-Vaginal	Sigma Acetic Acid
KY Warming Liquid	

Carryover

The carryover rate for Alinity m HR HPV assay was evaluated in 3 studies. In study 1, the sample input rack and sample processing unit carryover rate was determined by analyzing 265 replicates of HPV negative samples processed from alternating positions in the sample rack with high concentration HPV positive samples (representative of cervical specimens) at 10,000,000 copies/mL, across multiple runs. HPV was not detected in any HPV negative sample, resulting in an overall carryover rate of 0.0% (95% CI: 0.0% to 1.4%).

In study 2, the sample input rack and sample processing unit carryover rate for vaginal collection devices was determined by analyzing 240 replicates of HPV negative samples processed from alternating positions in the sample rack with high concentration HPV positive sample prepared in Specimen Buffer at 10,000,000 copies/mL, across multiple runs. HPV was detected in 2 HPV negative samples, resulting in an overall carryover rate of 0.8% (95% CI: 0.2% to 3.0%).

In study 3, the AMP tray carryover rate was determined by analyzing 285 replicates of HPV negative samples processed from alternating positions on the AMP tray with high concentration HPV positive samples at >10,000,000 copies/mL, across multiple runs. HPV was not detected in any HPV negative sample, resulting in an overall carryover rate of 0.0% (95% CI: 0.0% to 1.3%).

Agreement Between Specimen Types

Specimens collected in ThinPrep PreservCyt Solution and with the Alinity m Cervi-Collect Specimen Collection Kit from the same 277 patients were tested with Alinity m HR HPV. The overall percent agreement between the Alinity m HR HPV interpretations of the two specimen types (Table 12) was 94.6% (262/277; 95% CI: 91.3% to 96.7%), with average positive agreement of 93.4% (214/229; 95% CI: 89.7% to 96.5%) and average negative agreement of 95.4% (310/325; 95% CI: 92.7% to 97.6%).

Table 12. Agreement between Cervi-Collect and ThinPrep PreservCyt Specimens

	Cervi-Collect	
	HR HPV Detected	Not Detected
ThinPrep PreservCyt	HR HPV Detected	107
	Not Detected	7
		8
		155

Specimens collected in SurePath Preservative Fluid and ThinPrep PreservCyt Solution from the same 276 patients were tested with Alinity m HR HPV. The overall percent agreement between the Alinity m HR HPV interpretations of the two specimen types (Table 13) was 93.8% (259/276; 95% CI: 90.4% to 96.1%), with average positive agreement of 91.5% (184/201; 95% CI: 87.2% to 95.1%) and average negative agreement of 95.2% (334/351; 95% CI: 92.6% to 97.3%).

Table 13. Agreement between SurePath and ThinPrep PreservCyt Specimens

	SurePath	
	HR HPV Detected	Not Detected
ThinPrep PreservCyt	HR HPV Detected	92
	Not Detected	7
		10
		167

Agreement Between ThinPrep PreservCyt Pre-Cytology and Post-Cytology Samples

The aliquots of 119 ThinPrep PreservCyt specimens that were removed both prior to and after cytological processing (ie, pre-cytology and post-cytology samples) were tested with Alinity m HR HPV. The overall percent agreement between the Alinity m HR HPV interpretations of the pre-cytology and post-cytology samples (Table 14) was 97.5% (116/119; 95% CI: 92.8% to 99.1%), with average positive agreement of 97.1% (102/105; 95% CI: 93.2% to 100.0%) and average negative agreement of 97.7% (130/133; 95% CI: 94.7% to 100.0%).

Table 14. Agreement between Pre-Cytology and Post-Cytology ThinPrep PreservCyt Samples

	Post-Cytology	
	HR HPV Detected	Not Detected
Pre-Cytology	HR HPV Detected	51
	Not Detected	0
		3
		65

Agreement Between SurePath Pre-Cytology and Post-Cytology Samples

The aliquots of 112 SurePath specimens that were removed both prior to and after cytological processing (ie, pre-cytology and post-cytology samples) were tested with Alinity m HR HPV.³³ The overall percent agreement between the Alinity m HR HPV interpretations of the pre-cytology and post-cytology samples (Table 15) was 92.9% (104/112; 95% CI: 86.5% to 96.3%), with average positive agreement of 90.5% (76/84; 95% CI: 83.4% to 96.3%) and average negative agreement of 94.3% (132/140; 95% CI: 89.6% to 97.8%).

Table 15. Agreement between Pre-Cytology and Post-Cytology SurePath Samples

	Post-Cytology	
	HR HPV Detected	Not Detected
Pre-Cytology	HR HPV Detected	38
	Not Detected	1
		7
		66

Agreement Between Self-Collected simpli-COLLECT Vaginal Specimens and Clinician-Collected Cervical Specimens

A comparison of Alinity m HR HPV assay results from self-collected vaginal specimens and clinician-collected cervical specimens was performed using paired specimens from study participants who were 25 years or older presenting for routine cervical cancer screening. Each participant collected a vaginal specimen using the simpli-COLLECT HPV Collection Kit (ie, simpli-COLLECT) in a simulated home setting. Following the self-collection of the vaginal specimen, a cervical specimen was collected from each participant by the clinician and placed in ThinPrep PreservCyt Solution (ie, PreservCyt) per standard of care.

968 participants with valid results from paired clinician-collected cervical specimen and self-collected vaginal specimen using the simpli-COLLECT Kit were evaluated. For 85 subjects with positive cervical results, 81 had a positive result for the self-collected vaginal specimen with the simpli-COLLECT Kit (ie, PPA 95.3%). For 883 subjects with negative cervical results, 794 had a negative result for the self-collected vaginal specimen with simpli-COLLECT Kit (ie, NPA 89.9%). Refer to Table 16 for the positive agreement, negative agreement, and overall percent agreement.

Table 16. Agreement between Self-Collected Vaginal Specimens Using simpli-COLLECT and Clinician-Collected Cervical PreservCyt Specimens		
	Estimate (95% CI)	n/N ^a
Positive Percent Agreement	95.3 (88.5, 98.2)	81/85
Negative Percent Agreement	89.9 (87.8, 91.7)	794/883
Overall Percent Agreement	90.4 (88.4, 92.1)	875/968

^a n represents the number of subjects with self-collected vaginal specimen results in agreement with clinician-collected cervical specimen results. N represents the number of subjects.

Agreement Between Self-Collected Evalyn Brush Vaginal Specimens and Clinician-Collected Cervical Specimens

A comparison of Alinity m HR HPV assay results from self-collected vaginal specimens and clinician-collected cervical specimens was performed using paired specimens from study participants who were 25 years or older presenting for routine cervical cancer screening. Each participant collected a vaginal specimen using an Evalyn Brush (Rovers Medical Devices) in a simulated home setting. Following the self-collection of the vaginal specimen, a cervical specimen was collected from each participant by the clinician and placed into PreservCyt per standard of care.

960 participants with valid results from paired clinician-collected cervical specimen and self-collected vaginal specimen using Evalyn Brush were evaluated. For 84 subjects with positive cervical results, 80 had a positive result for the self-collected vaginal specimen with Evalyn Brush (ie, PPA 95.2%). For 876 subjects with negative cervical results, 808 had a negative result for the self-collected vaginal specimen with Evalyn Brush (ie, NPA 92.2%). Refer to **Table 17** for the positive agreement, negative agreement, and overall percent agreement.

Table 17. Agreement between Self-Collected Vaginal Specimens Using Evalyn Brush and Clinician-Collected Cervical PreservCyt Specimens		
	Estimate (95% CI)	n/N ^a
Positive Percent Agreement	95.2 (88.4, 98.1)	80/84
Negative Percent Agreement	92.2 (90.3, 93.8)	808/876
Overall Percent Agreement	92.5 (90.7, 94.0)	888/960

^a n represents the number of subjects with self-collected vaginal specimen results in agreement with clinician-collected cervical specimen results. N represents the number of subjects.

BIBLIOGRAPHY

- Howley PM. Papillomaviridae: the viruses and their replication. In: Fields BN, Knipe DM, Howley PM, eds. *Fields Virology*. 3rd ed. Lippincott-Raven Publishers; 1996:947-78.
- CDC. Genital HPV Infection - CDC Fact Sheet. 2022; <https://www.cdc.gov/std/hpv/stdfact-hpv.htm>.
- zur Hausen H. Papillomaviruses and cancer: from basic studies to clinical application. *Nat Rev Cancer*. 2002;2(5):342-350. doi:10.1038/nrc798
- Walboomers JM, Jacobs MV, Manos MM, et al. Human papillomavirus is a necessary cause of invasive cervical cancer worldwide. *J Pathol*. 1999;189(1):12-19. doi:10.1002/(SICI)1096-9896(199909)189:1<12::AID-PATH431>3.0.CO;2-F
- Snijders PJ, Steenbergen RD, Heideman DA, Meijer CJ. HPV-mediated cervical carcinogenesis: concepts and clinical implications. *J Pathol*. 2006;208(2):152-164. doi:10.1002/path.1866
- Kjaer SK, van den Brule AJ, Paull G, et al. Type specific persistence of high risk human papillomavirus (HPV) as indicator of high grade cervical squamous intraepithelial lesions in young women: population based prospective follow up study. *BMJ*. 2002;325(7364):572. doi:10.1136/bmj.325.7364.572
- Cuschieri KS, Cubie HA, Whitley MW, et al. Persistent high risk HPV infection associated with development of cervical neoplasia in a prospective population study. *J Clin Pathol*. 2005;58(9):946-950. doi:10.1136/jcp.2004.022863
- Bravo IG, Féliz-Sánchez M. Papillomaviruses: Viral evolution, cancer and evolutionary medicine. *Evol Med Public Health*. 2015;2015(1):32-51. doi:10.1093/emph/eov003
- de Villiers EM, Fauquet C, Broker TR, Bernard HU, zur Hausen H. Classification of papillomaviruses. *Virology*. 2004;324(1):17-27. doi:10.1016/j.virol.2004.03.033
- IARC Monographs on the evaluation of carcinogenic risks to humans. Human Papillomaviruses. Lyon: International Agency for Research on Cancer 2007; Volume 90.
- Muñoz N, Bosch FX, de Sanjosé S, et al. Epidemiologic classification of human papillomavirus types associated with cervical cancer. *N Engl J Med*. 2003;348(6):518-527. doi:10.1056/NEJMoa021641
- Clifford GM, Smith JS, Plummer M, Muñoz N, Franceschi S. Human papillomavirus types in invasive cervical cancer worldwide: a meta-analysis. *Br J Cancer*. 2003;88(1):63-73. doi:10.1038/sj.bjc.6600688
- Muñoz N, Castellsagué X, Berrington de González A, Gissmann L. Chapter 1: HPV in the etiology of human cancer. *Vaccine*. 2006;24 Suppl 3:S3/1-S3/10. doi:10.1016/j.vaccine.2006.05.115
- de Sanjose S, Quint WG, Alemany L, et al. Human papillomavirus genotype attribution in invasive cervical cancer: a retrospective cross-sectional worldwide study. *Lancet Oncol*. 2010;11(11):1048-1056. doi:10.1016/S1470-2045(10)70230-8
- Li N, Franceschi S, Howell-Jones R, Snijders PJ, Clifford GM. Human papillomavirus type distribution in 30,848 invasive cervical cancers worldwide: Variation by geographical region, histological type and year of publication. *Int J Cancer*. 2011;128(4):927-935. doi:10.1002/ijc.25396
- Smith JS, Lindsay L, Hoots B, et al. Human papillomavirus type distribution in invasive cervical cancer and high-grade cervical lesions: a meta-analysis update. *Int J Cancer*. 2007;121(3):621-632. doi:10.1002/ijc.22527
- Duska LR. Can we improve the detection of glandular cervical lesions: the role and limitations of the Pap smear diagnosis atypical glandular cells (AGC). *Gynecol Oncol*. 2009;114(3):381-382. doi:10.1016/j.ygyno.2009.07.008
- Khan MJ, Castle PE, Lorincz AT, et al. The elevated 10-year risk of cervical precancer and cancer in women with human papillomavirus (HPV) type 16 or 18 and the possible utility of type-specific HPV testing in clinical practice. *J Natl Cancer Inst*. 2005;97(14):1072-1079. doi:10.1093/jnci/dji187
- Davies P, Arbyn M, Dillner J, et al. A report on the current status of European research on the use of human papillomavirus testing for primary cervical cancer screening. *Int J Cancer*. 2006;118(4):791-796. doi:10.1002/ijc.21611
- Cuzick J, Clavel C, Petry KU, et al. Overview of the European and North American studies on HPV testing in primary cervical cancer screening. *Int J Cancer*. 2006;119(5):1095-1101. doi:10.1002/ijc.21955
- Mayrand MH, Duarte-Franco E, Rodrigues I, et al. Human papillomavirus DNA versus Papanicolaou screening tests for cervical cancer. *N Engl J Med*. 2007;357(16):1579-1588. doi:10.1056/NEJMoa071430

22. Goldie SJ, Gaffikin L, Goldhaber-Fiebert JD, et al. Cost-effectiveness of cervical-cancer screening in five developing countries. *N Engl J Med*. 2005;353(20):2158-2168. doi:10.1056/NEJMsa044278
23. Kim JJ, Wright TC, Goldie SJ. Cost-effectiveness of human papillomavirus DNA testing in the United Kingdom, The Netherlands, France, and Italy. *J Natl Cancer Inst*. 2005;97(12):888-895. doi:10.1093/jnci/dji162
24. Goldie SJ, Kim JJ, Wright TC. Cost-effectiveness of human papillomavirus DNA testing for cervical cancer screening in women aged 30 years or more. *Obstet Gynecol*. 2004;103(4):619-631. doi:10.1097/01.AOG.0000120143.50098.c7
25. Meijer CJ, Berkhof J, Castle PE, et al. Guidelines for human papillomavirus DNA testing for primary cervical cancer screening in women 30 years and older. *Int J Cancer*. 2009;124(3):516-520. doi:10.1002/ijc.24010
26. Cuschieri KS, Cubie HA. The role of human papillomavirus testing in cervical screening. *J Clin Virol*. 2005;32 Suppl 1:S34-S42. doi:10.1016/j.jcv.2004.11.020
27. Franco EL, Cuzick J. Cervical cancer screening following prophylactic human papillomavirus vaccination. *Vaccine*. 2008;26 Suppl 1:A16-A23. doi:10.1016/j.vaccine.2007.11.069
28. Stanley M, Villa LL. Monitoring HPV vaccination. *Vaccine*. 2008;26 Suppl 1:A24-A27. doi:10.1016/j.vaccine.2007.11.059
29. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 6th ed. Washington, DC: US Government Printing Office; December 2020. [Also available online. Type> www.cdc.gov, search>BMBL>look up sections III and IV.]
30. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, *Bloodborne pathogens*.
31. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
32. World Health Organization. *Laboratory Biosafety Manual*. 4th ed. Geneva: World Health Organization; 2020.
33. Jang D, Ratnam S, Smieja M, et al. Comparison of Alinity m HPV and cobas HPV assays on cervical specimens in diverse storage media. *Tumour Virus Res*. 2021;12:200224. doi:10.1016/j.tvr.2021.200224

KEY TO SYMBOLS

	Reference Number
	In Vitro Diagnostic Medical Device
	Lot Number
	In Vitro Test
	AMP TRAY
	ACT TRAY
	Unit
	Health Hazard
	Exclamation Mark
	Consult Instructions for Use
	Temperature Limit
	Use By Date
	Contains sufficient for <n> tests
	Authorized Representative in the European Community
	Manufacturer

TECHNICAL ASSISTANCE

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

Abbott Molecular Inc. is the legal manufacturer of the Alinity m HR HPV AMP Kit.

SUMMARY OF SAFETY AND PERFORMANCE STATEMENT

A summary of safety and performance (SSP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSP location after the launch of European Database on Medical Devices. Search for device using UDI-DI provided on the outer packaging of the device.

The Alinity m HR HPV AMP Kit is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

©2024 Abbott. All Rights Reserved.

Alinity is a trademark of Abbott. simpli-COLLECT is a trademark of Abbott. All other trademarks are property of their respective owners.
53-608488/R1

April 2024

Specimen Buffer Kit II

REF 09R03-050
53-608399/R2

Created August 2024

REF 09R03-050

53-608399/R2

PN MD2361

9PIMD2361rev00

CUSTOMER SERVICE: 1-800-553-7042

CUSTOMER SERVICE INTERNATIONAL:

CALL YOUR ABBOTT REPRESENTATIVE

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NOTICE TO USER

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the customer service or technical assistance section of these instructions.

NAME

Alinity m Specimen Buffer Kit II

INTENDED USE

The Alinity m Specimen Buffer Kit II contains bottles of bulk buffer solution intended for the resuspension of specimens for testing on the automated Alinity m System; refer to the appropriate Alinity m assay package inserts for additional information.

INTENDED USER

The intended users for the Alinity m Specimen Buffer Kit II are laboratory and healthcare professionals.

REAGENTS

Kit Contents

Alinity m Specimen Buffer Kit II List No. 09R03-050

The Alinity m Specimen Buffer Kit II contains the following:

- Alinity m Specimen Buffer II bottles
- Each bottle contains guanidine thiocyanate in Tris buffer.

Reagent	Quantity
Alinity m Specimen Buffer II	4 bottles × 237 mL

WARNINGS AND PRECAUTIONS

IVD

For In Vitro Diagnostic Use

Safety Precautions

Wear disposable gloves while handling specimens and wash hands thoroughly afterward. Use of protective eyewear is recommended.

CAUTION: This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled with appropriate biosafety practices.^{1,2}

The following warnings and precautions apply to:

Alinity m Specimen Buffer II

Hazard-determining components of labeling:

Guanidine thiocyanate

EUH032 Contact with acids liberates very toxic gas.

Disposal

P501 Dispose of contents / container in accordance with local regulations.

Important information regarding the safe handling, transport and disposal of this product is contained in the Safety Data Sheet.

Safety Data Sheets are available from your Abbott Representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity m System Operations Manual; Section 8.

Reagent Shipping

	Shipment Condition
Alinity m Specimen Buffer Kit II	15°C to 30°C

Reagent Storage

	Storage Temperature	Maximum Storage Time
Unopened	15°C to 30°C	Until expiration date

Reagent Handling

- Do not use reagents that have been damaged.
- Do not use past expiration date.
- Alinity m Specimen Buffer Kit II is shipped at 15°C to 30°C and is stored at 15°C to 30°C upon arrival. If reagents arrive in a condition contrary to this recommendation or are damaged, immediately contact your Abbott Representative.

PROCEDURE

Materials Provided

09R03-050 Alinity m Specimen Buffer Kit II

Instructions for Use

Refer to the appropriate Alinity m assay Package Insert for the Alinity m Specimen Buffer Kit II instructions for use.

QUALITY CONTROL PROCEDURES

Refer to the **QUALITY CONTROL PROCEDURES** section of the appropriate Alinity m assay Package Insert.

BIBLIOGRAPHY

- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. CLSI; 2014.
- US Environmental Protection Agency. *EPA Guide for Infectious Waste Management*. US Environmental Protection Agency; 1986:1-1-5-5, R1-R3, A1-A24. Publication No. EPA/530-SW-86-014.

KEY TO SYMBOLS

REF	Reference Number
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PN	Part Number
In Vitro Test	In Vitro Test
SPECIMEN BUFFER II	Specimen Buffer II
	Consult Instructions for Use or Consult Electronic Instructions for Use
	Temperature Limit
For In Vitro Diagnostic Use	For In Vitro Diagnostic Use
	Use By Date
EC REP	Authorized Representative in the European Community/ European Union
	Manufacturer

TECHNICAL ASSISTANCE

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

The Alinity m Specimen Buffer Kit II is manufactured for Abbott Molecular Inc. by Promega Corporation, Madison, WI 53711 USA.

Abbott Molecular Inc. is the legal manufacturer of the Alinity m Specimen Buffer Kit II.

The Alinity m Specimen Buffer Kit II is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

©2024 Abbott. All Rights Reserved.

Alinity is a trademark of Abbott. All other trademarks are property of their respective owners.

53-608399/R2

August 2024

Specimen Buffer Kit III

REF 09R03-060
53-608469/R1

Created April 2024

REF 09R03-060
53-608469/R1**CUSTOMER SERVICE: 1-800-553-7042**
CUSTOMER SERVICE INTERNATIONAL:
CALL YOUR ABBOTT REPRESENTATIVE

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NOTICE TO USER

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the customer service or technical assistance section of these instructions.

NAME

Alinity m Specimen Buffer Kit III

INTENDED USE

The Alinity m Specimen Buffer Kit III contains transport tubes filled with buffer solution intended for the resuspension of specimens for testing on the automated Alinity m System; refer to the appropriate Alinity m assay package inserts for additional information.

INTENDED USER

The intended users for the Alinity m Specimen Buffer Kit III are laboratory and healthcare professionals.

REAGENTS**Kit Contents**

Alinity m Specimen Buffer Kit III List No. 09R03-060

The Alinity m Specimen Buffer Kit III contains the following:

- Transport tubes with pierceable cap containing 2.95 mL Alinity m Specimen Buffer III (guanidine thiocyanate in Tris buffer)

Reagent	Quantity
Alinity m Specimen Buffer III	36 tubes

WARNINGS AND PRECAUTIONS**IVD****For In Vitro Diagnostic Use****Safety Precautions**

Wear disposable gloves while handling specimens and wash hands thoroughly afterward. Use of protective eyewear is recommended.

CAUTION: This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled with appropriate biosafety practices.^{1,2}

The following warnings and precautions apply to:

Alinity m Specimen Buffer III

Hazard-determining components of labeling:

Guanidine thiocyanate

EUH032 Contact with acids liberates very toxic gas.

Disposal

P501 Dispose of contents / container in accordance with local regulations.

Important information regarding the safe handling, transport and disposal of this product is contained in the Safety Data Sheet.

Safety Data Sheets are available from your Abbott Representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity m System Operations Manual; Section 8.

Reagent Shipping

	Shipment Condition
Alinity m Specimen Buffer Kit III	15°C to 30°C

Reagent Storage

	Storage Temperature	Maximum Storage Time
Unopened	15°C to 30°C	Until expiration date

Reagent Handling

- Do not use reagents that have been damaged.
- Do not use past expiration date.
- Alinity m Specimen Buffer Kit III is shipped at 15°C to 30°C and is stored at 15°C to 30°C upon arrival. If reagents arrive in a condition contrary to this recommendation or are damaged, immediately contact your Abbott Representative.
- Alinity m Specimen Buffer Kit III contains single-use tubes with pierceable caps. Avoid contamination or damage to the caps after removal from their original packaging.
- Specimen labels, including barcodes, must be placed below the blue ring present on the tube. Blue ring must remain entirely visible after label placement.

PROCEDURE**Materials Provided**

09R03-060 Alinity m Specimen Buffer Kit III

Instructions for Use

Refer to the appropriate Alinity m assay Package Insert for the Alinity m Specimen Buffer Kit III instructions for use.

QUALITY CONTROL PROCEDURES

Refer to the **QUALITY CONTROL PROCEDURES** section of the appropriate Alinity m assay Package Insert.

BIBLIOGRAPHY

- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. CLSI; 2014.
- US Environmental Protection Agency. *EPA Guide for Infectious Waste Management*. US Environmental Protection Agency; 1986:1-1-5-5, R1-R3, A1-A24. Publication No. EPA/530-SW-86-014.

KEY TO SYMBOLS

REF	Reference Number
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
In Vitro Test	In Vitro Test
For In Vitro Diagnostic Use	For In Vitro Diagnostic Use
SPECIMEN BUFFER III	Specimen Buffer III
PRODUCED BY	Produced by
QTY	Quantity
	Consult Instructions for Use
	Temperature Limit
	Use By Date
EC REP	Authorized Representative in the European Community
	Manufacturer

TECHNICAL ASSISTANCE

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

The Alinity m Specimen Buffer Kit III is manufactured for Abbott Molecular Inc. by MML Diagnostics Packaging, Inc., Troutdale, OR 97060 USA.

Abbott Molecular Inc. is the legal manufacturer of the Alinity m Specimen Buffer Kit III.

The Alinity m Specimen Buffer Kit III is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

 Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

EC REP Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

©2024 Abbott. All Rights Reserved.

Alinity is a trademark of Abbott. All other trademarks are property of their respective owners.

53-608469/R1

April 2024

simpli-COLLECT™ HPV Collection Kit

en
REF 09R03-030
53-608547/R1

Created November 2024

REF 09R03-030
53-608547/R1

CUSTOMER SERVICE: 1-800-553-7042 CUSTOMER SERVICE INTERNATIONAL: CALL YOUR ABBOTT REPRESENTATIVE

Instructions must be carefully followed. Reliability of results cannot be guaranteed if there are any deviations from these instructions.

NOTICE TO USER

If a serious incident occurs in relation to this device, the incident should be reported to the manufacturer and to the appropriate competent authority of the member state in which the user and/or the patient is established. To report to the manufacturer, see the contact information provided in the Customer Service section or Technical Assistance section of these instructions.

NAME

simpli-COLLECT™ HPV Collection Kit

INTENDED USE

The simpli-Collect HPV Collection Kit is intended for the self-collection and transport of vaginal swab specimens for testing with the Alinity m HR HPV assay. Swab specimens may be collected and prepared for transport in a healthcare or non-healthcare environment by a lay user. The collected specimens are intended to be tested on the automated Alinity m System; refer to the Alinity m HR HPV assay package insert for additional information.

INTENDED USER

The intended users for the simpli-COLLECT HPV Collection Kit are lay users, and laboratory and healthcare professionals.

SUMMARY AND EXPLANATION OF THE TEST

The simpli-COLLECT HPV Collection Kit allows patients to self-collect vaginal swab specimens without supervision. The kit and patient instructions for use are given or sent to the patient. The patient collects and transfers the specimen to the Alinity m Sample Tube before returning it to the laboratory for analysis.

CONTENTS

Case Contents

Component	Quantity
simpli-COLLECT HPV Collection Kit	40 individually-wrapped simpli-COLLECT HPV Collection Kits and 40 patient instructions for use per box, 10 boxes per case

simpli-COLLECT HPV Collection Kit Contents

Each non-reusable simpli-COLLECT HPV Collection Kit includes the following:

- 1 Alinity m Sample Tube
- 1 Specimen Collection Swab

WARNINGS AND PRECAUTIONS

IVD

For In Vitro Diagnostic Use

Safety Precautions

- Follow all instructions in this package insert.
- Do not send the simpli-COLLECT HPV Collection Kit to the patient if the kit is damaged. Discard unused or damaged kits in accordance with local, regional, and national regulations.^{1,2}
- Do not send the simpli-COLLECT HPV Collection Kit to the patient if the kit is beyond its expiration date.
- Wear disposable gloves while handling specimens and wash hands thoroughly afterward. Use of protective eyewear is recommended.
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local, regional, and national regulations.^{1,2}

CAUTION: This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled with appropriate biosafety practices.^{1,2}

Storage

	Storage Temperature	Maximum Storage Time
Unopened	15°C to 30°C	Until expiration date

PROCEDURE

Materials Provided (per case)

- 10 boxes, each box contains 40 simpli-COLLECT HPV Collection Kits and 40 patient instructions for use

Specimen Collection, Storage, and Transport

Provide the simpli-COLLECT HPV Collection Kit along with the Patient Instructions for Use to the patient. Specimens must include identification information, including the date of collection. Specimens should be transported at ambient temperature. Refer to the Alinity m HR HPV (List No. 09N15-092) package insert for information on specimen storage. For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable local, national, and international regulations covering the transport of clinical, diagnostic, or biological specimens.

LIMITATIONS OF THE PROCEDURE

- Optimal performance of this kit requires appropriate specimen collection, handling, preparation, and storage. Refer to the Alinity m HR HPV package insert for additional information.
- This kit should only be used to collect vaginal swab samples. Refer to the patient instructions for use and the Alinity m HR HPV (List No. 09N15-092) package insert for additional information.
- Pregnant women should not use the simpli-COLLECT HPV Collection Kit.

BIBLIOGRAPHY

1. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
2. US Environmental Protection Agency. *EPA Guide for Infectious Waste Management*. US Environmental Protection Agency, 1986:1-1–5-5, R1-R3, A1-A24. Publication No. EPA/530-SW-86-014.

KEY TO SYMBOLS

REF	Reference Number
LOT	Lot Number
IVD	In Vitro Diagnostic Medical Device
In Vitro Test	In Vitro Test
STORE AT	Store At
UNIT	Unit
SWAB	Swab
SAMPLE TUBE	Sample Tube
GTIN	Global Trade Item Number
PRODUCED BY	Produced By
	Do Not Use if Package is Damaged
	Do Not Reuse
	Temperature Limit
	Expiration Date
	Plain Paper
	Ethylene-vinyl acetate (EVA)
	Consult Instructions for Use or Electronic Instructions for Use
EC REP	Authorized Representative in the European Community/European Union
	Manufacturer

TECHNICAL ASSISTANCE

For technical assistance, call Abbott Technical Services at 1-800-553-7042 (within the US) or +49-6122-580 (outside the US), or visit the Abbott website at www.molecular.abbott.

Abbott Molecular Inc. is the legal manufacturer of the simpli-COLLECT™ HPV Collection Kit.

simpli-COLLECT™ HPV Collection Kit is manufactured for Abbott Molecular Inc. by MML Diagnostics Packaging, Inc., Troutdale, OR 97060 USA.

The simpli-COLLECT™ HPV Collection Kit is imported into the European Union by Abbott Diagnostics GmbH, located at Max-Planck-Ring 2, 65205 Wiesbaden, Germany.

Abbott Molecular Inc.
1300 East Touhy Avenue
Des Plaines, IL 60018 USA

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

©2024 Abbott. All Rights Reserved.

Alinity is a trademark of Abbott. simpli-COLLECT is a trademark of Abbott. All other trademarks are property of their respective owners.

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You will receive a kit containing the following:

- 1 Tube (Sample Tube)
- 1 Swab

If you do not follow the instructions, your result may be incorrect.

- Store kit before and after use at room temperature (between 15°C to 30°C).
 - Return your sample within 24 hours.
 - Do not use if you are pregnant.
 - Do not use genital area immediately before collecting sample.
 - Do not use the kit beyond its **EXPIRATION DATE** printed on kit packaging.
- ¹ Hours
- ¹ Stunden

DO NOT touch the white tip of the swab or lay the swab down. If the white tip is touched or the swab is laid down or dropped, your results may not be accurate. You need to request a new **simpli-COLLECT HPV** Collection Kit.

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Contact your healthcare provider with questions regarding test results/ diagnosis or assistance with obtaining a new kit.

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For additional product support please visit: www.molecular.abbottni/en/simpli-collect-hpv

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EN simpli-COLLECT is a trademark of Abbott. All other trademarks are property of their respective owners. © [2025 Abbott](#). All rights reserved. [February 2025](#)