WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

Product: Bioline HCV

WHO reference number: PQDx 0257-012-00

Bioline HCV with product codes 02FK10, 02FK16 and 02FK17, manufactured by Abbott Diagnostics Korea Inc., Rest-of-World (RoW) regulatory version was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 29 November 2016.

Summary of the WHO prequalification assessment for the Bioline HCV

	Date	Outcome
Prequalification listing	29 November 2016	listed
Dossier review	18 October 2016	MR
Product performance evaluation	8 August 2016	MR

MR: Meets requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendments and change request reference where applicable.	Date of report amendment
1.0 to 4.0	Inclusion of Instructions for Use.	8 March 2016
5.0	Correction of a typographical error	20 December 2016
6.0	Addition of specimen collected from a finger prick and subsequently two additional product codes (02FK16 and 02FK17)	26 March 2018
7.0	Addition of supplier for safety lancet	10 September 2018
8.0	The product name was changed from SD BIOLINE HCV to Bioline HCV.	3 March 2020

	The manufacturer's name changed from Standard Diagnostics Inc to Abbott Diagnostics Korea Inc	
9.0	Accessory labelling changes according to accessory manufacturer information change (Manufacturer name change and EU representative address change) (PQC-IVD-2024-0033).	24 March 2026

Intended use:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“the BiolineHCV is an in vitro immunochromatographic, rapid assay designed for the qualitative detection of antibodies specific to HCV, in human serum, plasma (heparin, EDTA and sodium citrate) or whole blood. Bioline HCV is intended only for professional use as the initial test, as an aid to diagnosis. Reactive specimens should be reflexed for additional testing, either by nucleic acid testing (NAT) technologies for the detection of HCV RNA or HCV core antigen testing, to identify current HCV infection. This product is intended for use in a population with high HCV prevalence or who have a history of HCV risk exposure/behaviour including pregnant women. This test may not be suitable for diagnosis of early infection or blood donation screening. The performance of the assay has not been established for populations of infants or children.”*

Product test kit contents:

Component	Product code(s)		
	02FK10	02FK16	02FK17
Specimen procedure(s)	Serum, plasma and whole blood (Venous whole blood and finger-prick blood)		
Test devices with desiccant in an individual foil pouch	30 T/kit	25 T/kit	25 T/kit
Assay diluent	1 x 5 mL/vial	1 x 5 mL/vial	1 x 5 mL/vial
Capillary pipette(s)	n/a	25 (10 µL each)	25 (10 µL each)
Sterile lancet(s)	n/a	25 (Sterile lancet)	25 (Safety lancet)
Alcohol swap(s)	n/a	25	25
Instructions for use	1	1	1

Items required but not provided:

- Micropipette
- Protective gloves
- Timer
- Biohazard container

Storage:

The test kit must be stored between 1 °C and 30 °C.

Shelf-life upon manufacture¹:

24 months.

Product dossier assessment

Abbott Diagnostics Korea Inc. (Formerly Standard Diagnostics, Inc.) submitted a product dossier for Bioline HCV as per the “*Instructions for compilation of a product dossier*” (PQDx_018 v1). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO. The manufacturer's responses to the nonconformities found during dossier assessment findings were accepted on 18 October 2016.

Based on the product dossier assessment findings, the product dossier for Bioline HCV meets WHO prequalification requirements.

Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

(iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pqweb/vitro-diagnostics/who-public-inspection-reports>

Product performance evaluation

Bioline HCV was evaluated by WHO in the second quarter of 2016 using plasma specimens. From this evaluation, we drew the following conclusions:

- Bioline HCV assay is an immunochromatographic assay for the detection of antibodies to HCV in human serum, plasma and venous whole blood.
- A volume of 10 µL of a specimen is needed to perform the assay.
- This type of assay does require laboratory equipment (i.e. precision pipette and tips) and can be performed in laboratories with limited facilities.

In this limited performance evaluation on a panel of 483 specimens, we found the performance summarized below:

Performance characteristics in comparison with an agreed reference standard		
	Initial (95% CI)	Final (95% CI)
Sensitivity %	98.8% (95.6 – 99.7%)	100% (97.76 – 100%)
Specificity %	100% (98.85 – 100%)	100% (98.85 – 100%)
Invalid rate %	0%	

Additional performance characteristics	
Sensitivity during seroconversion on 4 seroconversion panels in comparison with a benchmark assay; DiaSorin Anti-HCV Murex EIA (Version 4.0)	Seroconversion sensitivity index of +2.0, therefore detection is 2 specimens (i.e. bleedings) later than the benchmark assay
Analytical sensitivity on a mixed titer panel in comparison with an agreed reference standard	15 of 15 specimens were correctly classified
Lot to lot variation on a dilution panel in comparison with an agreed reference standard	Acceptable
Key operational characteristics	
Validated specimen types	Serum, plasma (heparin, EDTA and sodium citrate), venous whole blood
Number of steps	3 without precision required
Time to result	5 minutes
Endpoint stability	20 minutes
Internal QC	Yes, the control line on the test device is an internal procedure control. The absence of the control line indicates that insufficient or improper assay diluent was added to the device.
In-use stability of reagents	Until expiry date

Labelling review

The labelling submitted for Bioline HCV was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

02FK10

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-11-26	02FK10-40-A0
Pouch / Device label	Device pouch label	A0	2019-11-01	FP26-A0
Reagent bottle labels	Assay diluent labels	A0	2019-10-25	FB07-A0
Instructions for Use (IFU)	IFU	A1	2022-12-23	02FK10/02FK16/02FK17-04-A1

02FK16

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-11-26	02FK16-40-A0
Pouch / Device label	Device pouch label	A0	2019-11-01	FP26-A0
Reagent bottle labels	Assay diluent labels	A0	2019-10-25	FB07-A0
Capillary Tube	40005439	0	2024-05-21	D10151501
Sterile non-safety lancet	P65-008	0	2023-11-10	D10151501
Sterile non-safety lancet	40006487	0	2024-10-01	D10151501
IFU	IFU	A1	2022-12-23	02FK10/02FK16/02FK17-04-A1

02FK17

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-11-26	02FK17-40-A0
Pouch / Device label	Device pouch label	A0	2019-11-01	FP26-A0
Reagent bottle labels	Assay diluent labels	A0	2019-10-25	FB07-A0
Capillary tube	40005439	0	2024-05-21	D10151501
Sterile safety lancet	P65-013/000/001	0	2023-11-10	D10151501
Sterile safety lancet	P65-013/003	0	2023-11-10	D10151501
Sterile safety lancet	P65-013/006	0	2023-11-10	D10151501
IFU	IFU	A1	2022-12-23	02FK10/02FK16/02FK17-04-A1

Labels

Bioline[®] HCV

Abbott

Abbott **Bioline[®] HCV**

EN Rapid, immunochromatographic test for the detection of antibody against Hepatitis C virus in human serum, plasma or whole blood	FR Test immunochromatographique rapide pour la détection des anticorps dirigés contre le virus de l'hépatite c dans le sérum, le plasma ou le sang total humains.	ES Rápida, prueba inmunocromatográfica para la detección de anticuerpos contra el virus de hepatitis C en el suero, plasma o sangre total.	PT TESTE imunocromatográfico rápido para a detecção do anticorpo anti-vírus de hepatite c em soro, plasma ou sangue total humanos.
CONTENTS: • 30 Test devices with desiccant in individual foil pouches • Assay diluent (1 x 5ml/vial) • 1 Instructions for use	CONTENU : • 30 cassettes-tests emballées individuellement avec un agent déshydratant • Diluant de test (1 x 5ml/flacon) • 1 Mode d'emploi	CONTENIDO : • 30 Dispositivos individuales de la prueba con un agente desecante • Diluyente del ensayo (1 x 5ml/vial) • 1 Instrucciones de uso	CONTEÚDO : • 30 dispositivos de teste individualmente embalados em bolsas de alumínio com um dessecante • Diluente de ensaio (1 x 5ml/frasco) • 1 Instruções de utilização

02FK10-40-A0

PRODUCT - Bioline™ HCV 30T

SIZE - 152 X 115 X 70 MM

 PMS 303 C
 PMS 143 C

PN. - 02FK10-40-A0
 REV. DATE - 2019.09.26
 CREATOR - THINKBOX

Abbott

Abbott

Abbott

Bioline[®] HCV

PRINT AREA

Bioline[®] HCV


REF 02FK10

Bioline[®] HCV

Abbott Diagnostics Korea Inc.
65, Boraehallim, Cheongju, Yongsin-si,
Gyeonggi-do, 17099, Republic of Korea
abbott.com/koct

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

HCV
Bioline™

02FK16-40-A0

EN
Rapid, immunochromatographic test for the detection of antibody against Hepatitis C virus in human serum, plasma or whole blood.

CONTENTS:
• 25 Test devices with desiccant in individual foil pouches
• Assay diluent (1 x 5 mL/vial)
• 25 Capillary pipettes (10 µL)
• 25 Sterile lancets, 25 Alcohol swabs
• 1 Instructions for use

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

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Bioline[®]
HCV

Abbott

Bioline[®]
HCV

Abbott

Abbott

Abbott

EN
Rapid, immunochromatographic test for the detection of antibody against Hepatitis C virus in human serum, plasma or whole blood

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CONTENTS :

- 25 Test devices with desiccant in individual foil pouches
- Assay diluent (1 x 5 ml/vial)
- 25 Capillary pipettes (10 µl)
- 25 Safety lancets, 25 Alcohol swabs
- 1 Instructions for use

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02FK17-40-A0

PRODUCT - Bioline™ HCV
25T, Safety Lancet

PN. - 02FK17-40-A0
REV. DATE - 2019_09_26
CREATOR - THINKBOX


SIZE - 152 X 115 X 75 MM

PMS 303 C
PMS 143 C

PRINT AREA

Bioline[®]
HCV

1°C 34°F | 30°C 86°F | 25 | IVD | REF 02FK17

Bioline[®]
HCV

Abbott Diagnostics Korea Inc.
65, Borahagalli-ro, Gilheung-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea
abbott.com/pcot

© 2019 Abbott. All rights reserved. All trademarks
referenced are trademarks of either the Abbott group of
companies or their respective owners.

PRODUCT - Bioline™ AL pouch 1-30 °C
Non-CE

SIZE - 108 X 32 MM

 PMS 303 C
 PMS 2925 C

CAT NO. - FP26-A0
REV. DATE - 2019.10.28
CREATOR - SCOTT KIM

 **Abbott**

Bioline™
RAPID DIAGNOSTIC TEST

 **Abbott Diagnostics Korea Inc.**
65, Borahagal-ro, Giheung-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea

1°C - 30°C
34°F - 86°F

FP26-A0

ITEM ITEM: Bioline™

LOT LOT NO.:

MFG.DATE (YYYY.MM.DD)
EXP.DATE (YYYY.MM.DD)

11mm

PRODUCT - Biline™

Assay diluent label

SIZE - 50 X 16 MM

PMS 303 C

PMS 2925 C

CAT NO. - FB07-A0

REV. DATE - 2019.10.21

CREATOR - SCOTT KIM

	ITEM HCV
Biline™ ASSAY DILUENT	VOL 5ml/vial
IVD 1°C → 30°C	LOT XXXXXXXXX
<small>Abbott Diagnostics Korea Inc. FB07-A0</small>	Y YYY.MM.DD

XXXXXXXX_Capillary Tube (XXul)XXPCS_60×60

XXµl Capillary Tube

Code No. : XXXXXXXX
LOT No. : XXXXXXXX
Quantity : XXPCS
MFG Date : XXXX.XX.XX

(01)0880989962XXXX

HLB Life Science Co.,Ltd
101-908, Digital Empirell, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16681, Republic Korea

EC REP Cmc Medical Devices & Drugs S.L.
C/ Horacio Lengo n18 C.P 29006 Málaga-Spain

XXµl Capillary Tube

Code No. : XXXXXXXX
LOT No. : XXXXXXXX
Quantity : XXPCS
MFG Date : XXXX.XX.XX

(01)0880989962XXXX

HLB Life Science Co.,Ltd
101-908, Digital Empirell, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16681, Republic Korea

EC REP Cmc Medical Devices & Drugs S.L.
C/ Horacio Lengo n18 C.P 29006 Málaga-Spain

PANTONE Reflex Blue C
Shandong
Size:50X40


Sterile Lancet for single use

LOT LOT NO./Type : 2305267303 / 28G
EXP. DATE : 2028-06-04
QTY : 25pcs

STERILE R **CE** 0123 
010004451700443617288004193285287101

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 60, 20537 Hamburg, Germany

Shandong Lianfa Medical Plastic Products Co., Ltd.
No. 3 Shuangshihai Sanjian Road, 250200 Zhangqiu City, Jinan, Shandong
PEOPLE'S REPUBLIC OF CHINA

(01)16957138602815
 (11)YYMMDD
 (17)YYMMDD
 (10)XXXXXXXX

Shanghai International Holding
 Corp. GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel:0049-40-2513175
 E-mail: shholding@hotmail.com

**Twist-off lancet (use with lancing device)
 25 pcs (28G I)**

[Instructions]

1. Unscrew the lancing device cap from the body of the lancing device.
2. Insert a sterile lancet into the lancing device firmly.
3. Twist off the protective cap of the lancet.
4. Carefully screw the lancing device cap back to the lancing device.
5. Adjust the puncture depth.
6. Put the lancing device against the side of the finger and press the release button to prick your fingertip.
7. Apply a drop of blood to test strip.
8. Recap lancet and discard it in appropriate container.
9. Place the lancing device cap back to the lancing device.

[Structure] The Twist-off Lancets is consist of needle, main body and protective cap.

[Performance] The device is sterilized by irradiation.

[Intended use] The twist-off lancet (use with lancing device) is a single use, sterile medical device, which is designed for use of micro blood sampling puncture to obtain capillary blood samples from the fingertip.

[Contraindications] None.

[Shelf life] 5 years.

[Storage] The products should be stored in a well-ventilated room and kept away from light, moisture and corrosive gases.

Tianjin Huahong Technology Co., Ltd.


A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone
 (Air Port Industrial Park), 300308 Tianjin, China
 Tel: 0086 22 58775333 400 8587753 Fax: 0086 22 58775353
 E-mail: huahongsales@163.com

[Cautions]

1. Read the instructions carefully before use.
2. Use the lancets within the expiry date.
3. Do not use if protective cap is removed or broken previously.
4. Do not store lancet in lancing device, discard lancet after use, do not re-use.
5. For your safety, please do not use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals.
6. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

 The lancet and lancing device must be compatible with each other. When you use Huahong's sterile lancets with other manufacturers' lancing devices, please be sure that these products are compatible. We recommend that you shall consult a health care professional before use.


Instructions for use :

Twist protective cap, pull off and discard it.

Apply blood collection to the disinfected sampling site and press it to the bottom quickly.

Squeeze the blood collection site to obtain the corresponding blood collection volume.

The sterile lancets after use should be placed in a special recovery container.

(01)16957138604734
 (11)YYMMDD
 (17)YYMMDD
 (10)XXXXXX

Sterile lancets (Safety lancets)

Type : 28G- XXIII

QTY : 25pcs

STERILE R Irradiation Sterilization CE 0197

Tianjin Huahong Technology Co., Ltd.
 A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone
 (Air Port Industrial Park), 300308 Tianjin, China
 Tel: 0086 22 58775333 400 8587753 Fax: 0086 22 58775353
 E-mail: huahongsales@163.com

EC/REP Shanghai International Holding Corp. GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel:0049-40-2513175
 E-mail: shholding@hotmail.com

[Structure] The safety Lancets consist of needle core, button, housing, bottom, spring and protective cap.

[Performance] The device is sterilized by irradiation and the needle tip is sterile.

[Intended use] The safety lancets are single use medical devices, which are designed to collect capillary blood sample.

[Contraindications] None.

[Cautions]

1. Read the instructions carefully before use.
2. Use the lancets within the expiry date.
3. Do not use if protective cap is removed or broken previously.
4. Discard the lancet after use, do not re-use.
5. For your safety, never use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals.
6. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

[Storage]

The products should be stored in a well-ventilated room and kept away from light , moisture and corrosive gases.

We recommend that you shall consult a health care professional before use.

LOT

PANTONE Reflex Blue C
Shandong
Size:50X40

Safety Lancet

LOT LOT NO./Type : 23032692 / 28G
EXP. DATE : 2028-05-09
QTY : 25pcs

STERILE R **CE** 0123 


(010004851700104720000001023032692)

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 60, 20537 Hamburg, Germany

Shandong Lianfa Medical Plastic Products Co., Ltd.
No. 3 Shuangshihai Sanjian Road 250200 Zhangou City, Jinan, Shandong
PEOPLE'S REPUBLIC OF CHINA


FRONT SIDE OF THE BAG

BACK SIDE OF THE BAG

WHITE BACKGROUND

WHITE BACKGROUND

- transparent window where the lancets are visible

UDI code printed directly on the bag during packing process

Instructions for use²

² English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

- 2) **Precise assay**
 Bioline™ HCV may exhibit precise effect. (False non-reactive result) in specimens which have higher than S/COD of approx. 11.0 in Abbott Architect and S/COD approx. 3.0 in Orbia HCV version 3.0 ELISA and HCV assays.
- 3) **Specimen matrix**
 Validation on whole blood was done by testing 500 negative and 100 positive anti-HCV specimens. The performance of Bioline™ HCV on whole blood was comparable to the performance on plasma specimens. (See table 1 and table 3 presented in study 2 above.) Validation on serum and different plasma specimen types (EDTA / Heparin / Sodium citrate) was performed by Geneva Red Cross. The results obtained on negative and positive specimens are identical on serum, EDTA plasma, heparin plasma and citrate plasma.
- | Specimen type | No. of positive specimens | Bioline™ HCV reactive/No. of positive specimens | No. of Bioline™ HCV non-reactive/No. of negative specimens |
|-------------------|---------------------------|---|--|
| Serum | 25/25 | 25/25 | 25/25 |
| EDTA plasma | 25/25 | 25/25 | 25/25 |
| Heparin plasma | 25/25 | 25/25 | 25/25 |
| No-Citrate plasma | 25/25 | 25/25 | 25/25 |
- 4) **Complement factors interference in fresh serum specimens**
 In total 25 negative specimens, spiked with an anti-HCV positive specimen were tested within 24 hours after collection and retested after being stored at 4 °C for 1, 2, 3 and 4 days. No differences were observed on the results obtained on the fresh specimens and the same specimens stored for 1 to 4 days at 4 °C.

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Product Disclaimer:
 While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of this diagnosis should consult a doctor for further confirmation of the test result.

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Bibliography of suggested reading / Bibliographie			
1. Park C, Condon W, Chu L, Spector AB, Fromont-Ravot C, Charlier M, Rice E. Expression and Identification of Hepatitis C Virus Polyprotein Cleavage Products. <i>Journal of Virology</i> , March 1998, p.1385-1395.	4. A. Yoshikawa, K. Takahashi, S. Kishimoto : Serodiagnosis of hepatitis C virus infection by ELISA for antibodies against the putative core protein (p20C) expressed in Escherichia coli. <i>Journal of Immunological Methods</i> , 148 (1992) 143-150.	7. Young Gyu Cho, Min Kyung Yi, Young Eui Jung, Chang Wuk Kim and Young Chul Sung : Cloning and Characterization of the Highly Immunogenic Region of HCV Genome from Korean Patients. <i>Mol. Cells</i> , Vol. 3, 4: 477-480.	10. S. Ohnawa, E. Caccamo, S. Griva, F. Gambetta, R. Calogeri, C. Ross and F. Bonelli : Expression in E. coli and purification of a chimeric p22-N53 recombinant antigen of Hepatitis C Virus(HCV). <i>Federation of European Biochemical Societies</i> , Volume 324, number 3, 253-257.

Bibliography of suggested reading / Bibliographie			
2. Young Gyu Cho, Min Kyung Yi, Young Eui Jung, Chang Wuk Kim and Young Chul Sung : Cloning and Characterization of the Highly Immunogenic Region of HCV Genome from Korean Patients. <i>Mol. Cells</i> , Vol. 3, 4: 477-480.	5. Mitchell BL, et al. Impact freeze-thaw cycles and storage time on plasma samples used in mass spectrometry based biomarker discovery projects. <i>Cancer Informatics</i> , 2007;6:168-184.	8. Cuhadar S, et al. The effect of storage time and freeze-thaw cycles on the stability of serum samples. <i>Biochimica Medica</i> , 2013;23(1):70-77.	9. A. Yoshikawa, K. Takahashi, S. Kishimoto : Serodiagnosis of hepatitis C virus infection by ELISA for antibodies against the putative core protein (p20C) expressed in Escherichia coli. <i>Journal of Immunological Methods</i> , 148 (1992) 143-150.

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