WHO Prequalification of In Vitro Diagnostics Programme PUBLIC ASSESSMENT REPORT

Product: Bioline HIV/Syphilis Duo
Number: PQDx 0179-012-00

Bioline HIV/Syphilis Duo with product codes 06FK30, 06FK35 and 06FK37, manufactured by Abbott Diagnostics Korea Inc., Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 28 October 2015.

Summary of prequalification status for the Bioline HIV/Syphilis Duo

	Date	Outcome
Status on PQ list	28 October 2015	listed
Dossier assessment	14 July 2015	MR
Product performance evaluation	Quarter (Q) 3 2014-Q1 2015	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendments and change request reference where applicable.	Date of report amendment
4.0	Updated labelling (IFU and labels). The product was rebranded from SD BIOLINE HIV/Syphilis Duo to Alere HIV/Syphilis Duo. Corrections/clarifications on IFU. Alere logo on IFU and boxes.	15 June 2017
5.0	Rebranding of Alere HIV/Syphilis Duo to SD BIOLINE HIV/Syphilis Duo. Standard Diagnostics has been merged to Abbott during the change process of the brand name of this product. The rebranding from SD BIOLINE HIV/Syphilis Duo to Alere HIV/Syphilis Duo (as reported in 2 December 2016) did not proceed as planned.	07 December 2018.
6.0	Changes to the manufacturer name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc and	20 August 2020

	product name from SD BIOLINE HIV/Syphilis Duo to Bioline HIV/Syphilis Duo.	
7.0	Addition of safety lancet resulting in a new product configuration with product code 06FK37.	22 March 2022
8.0	Accessory labelling changes according to accessory manufacturer information change (Manufacturer name change and EU representative address change) (PQC-IVD-2024-0033).	19 March 2026
9.0	Addition of a capillary tube label (PQC-IVD-2025-0118).	15 May 2026

Intended use:

According to Abbott Diagnostics Korea Inc, “*Bioline HIV/Syphilis Duo test is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 including subtype-O, HIV-2 and Syphilis (Treponema pallidum) in human serum, plasma or whole blood. The Bioline HIV/Syphilis Duo is intended only for professional use, and for in vitro diagnostic use. This test may not be suitable for diagnosis of early infection or blood donation screening. Reactive specimens should be confirmed by a supplemental assay such as ELISA or Western Blot for antibodies to HIV-1/HIV-2 and the treponemal test (e.g. TPPA, TPHA) and non-treponemal test (e.g. RPR and VDRL) for antibodies to T. pallidum*”.

Test kit contents

	25 tests/kit Product code 06FK30	25 tests/kit Product code 06FK35	25 tests/kit Product code 06FK37
Test devices Individually foil pouched with a desiccant	25	25	25
Assay diluent Plastic bottle	1 vial of 4mL	1 vial of 4mL	1 vial of 4mL
Specimen transfer devices (capillary pipettes) Disposable, 20µl	N/A	25 units	25 units
Lancets	N/A	25 units of sterile lancets	25 units of sterile safety lancets
Alcohol swabs Disposable	N/A	25 units	25 units
Instructions for Use	1 unit	1 unit	1 unit

Items required but not provided

- Micropipette,
- Protective gloves
- Timer
- Biohazard container

Storage

The test kit must be stored between 1 - 30 °C.

Shelf-life¹:

24 months.

Product dossier assessment

Abbott Diagnostics Korea Inc submitted a product dossier for BiolineHIV/Syphilis Duo as per the "Instructions for compilation of a product dossier" (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2).

Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for Bioline HIV/Syphilis Duo for prequalification.

Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pqweb/vitro-diagnostics/who-public-inspection-reports>

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

Product performance evaluation

Bioline HIV/Syphilis Duo was evaluated at the Institute of Tropical Medicine, Antwerp, Belgium, on behalf of WHO in the last quarter of 2014/first quarter of 2015 using serum/plasma specimens. From this evaluation, we drew the following conclusions:

Bioline HIV/Syphilis Duo is an immunochromatographic assay for the detection of HIV-1/2 and syphilis antibodies in human serum/plasma, fingerstick/venous whole blood specimens. A volume of 10 µL of serum/plasma is needed to perform the assay. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities and non-laboratory settings. Reading of the results is done visually i.e. subjectively read.

In this limited evaluation a panel of 400 clinically-derived specimens was used. For antibodies to HIV, we observed an initial sensitivity (N=200) of 100% (95% CI: 98.2% - 100%) and an initial specificity (N=200) of 99.5% (95% CI: 97.2% - 100%) compared to the reference assays (Vironostika HIV Ag/Ab [bioMérieux] and Enzygnost anti-HIV 1/2 [Siemens Healthcare Diagnostics] in parallel; followed by INNO-LIA HIV I/II Score [Fujirebio]). The final sensitivity was 100% (95% CI: 98.2% - 100%) and the final specificity was 99.5% (95% CI: 97.2% - 100%) compared to the reference assays. Lot to lot variation was acceptable.

For the case of antibodies to *Treponema pallidum*, we observed an initial sensitivity (N=200) of 86.5% (95% CI: 81.0% - 90.9%) and an initial specificity (N=200) of 99.5% (95% CI: 97.2% - 100%) compared to the reference assays (Vitros Syphilis [Ortho clinical Diagnostics] followed by Serodia TPPA [Fujirebio] to confirm initial reactive results). The final sensitivity was 87.0% (95% CI: 81.5% - 91.3%) and the final specificity was 99.5% (95% CI: 97.2% - 100%) compared to the reference assays. Lot to lot variation was acceptable.

For HIV, eight seroconversion panels were evaluated. Bioline HIV/Syphilis Duo detected on average 0.125 specimens earlier than the benchmark assay (Enzygnost Anti-HIV 1/2 Plus).

For syphilis one seroconversion panel was evaluated. Bioline HIV/Syphilis Duo detected 2 specimens earlier than the benchmark assay (Vitros Syphilis).

For the HIV mixed titer panel, Bioline HIV/Syphilis Duo correctly identified all specimens compared to the reference result.

For the syphilis mixed titer panel, Bioline HIV/Syphilis Duo correctly identified all specimens compared to the reference result.

For the 1st International Reference Panel for anti-HIV [NIBSC code 02/210] Bioline HIV/Syphilis Duo gave an indeterminate result for the HIV-1 group O specimen. Bioline HIV/Syphilis Duo correctly identified all HIV-1 group M subtypes and HIV-2 specimens.

For the 1st International Reference Panel for anti-treponemal [NIBSC code 05/132], Bioline HIV/Syphilis Duo correctly identified all specimens.

In this study, 2.5% of the results for syphilis were recorded as indeterminate, 0% for HIV. Results were interpreted independently by three technicians; the inter-reader variability was 0% for HIV and 4% for syphilis. The invalid rate was 0%.

Labelling review

The labelling submitted for Bioline HIV/Syphilis Duo was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

Controlled Labelling References

06FK30

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-11-20	06FK30-40-A0
Pouch / Device label	Device pouch label	A0	2019-11-01	FP26-A0
Reagent bottle labels	Assay diluent labels	A0	2019-11-25	FB09-A0
Instructions for Use (IFU)	IFU	A0	2019-11-20	06FK30/06FK35-04-A0

06FK35

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-11-20	06FK35-40-A0
Pouch / Device label	Device pouch label	A0	2019-11-01	FP26-A0
Reagent bottle labels	Assay diluent labels	A0	2019-11-25	FB09-A0
Capillary Tube	40005438	2	2024-05-21	D10151501
Capillary Tube	40001927	0	2026-03-23	D10167079
Sterile non-safety lancet	P65-008	0	2023-11-10	D10151501
Sterile non-safety lancet	40006487	0	2024-10-01	D10151501
IFU	IFU	A0	2019-11-20	06FK30/06FK35-04-A0

06FK37

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2021-04-08	06FK37-40-A0
Pouch / Device label	Device pouch label	A0	2019-11-01	FP26-A0
Reagent bottle labels	Assay diluent labels	A0	2019-11-25	FB09-A0
Capillary Tube	40005438	2	2024-05-21	D10151501
Capillary Tube	40001927	0	2026-03-23	D10167079
Sterile safety lancet	P65-013/000/001	0	2023-11-10	D10151501
Sterile safety lancet	P65-013/003	0	2023-11-10	D10151501
Sterile safety lancet	P65-013/006	0	2023-11-10	D10151501
IFU	IFU	A0	2021-04-07	06FK37-04-A0

Labels

PRODUCT - Bioline™ HIV/Syphilis Duo 25T
 SIZE - 152 X 115 X 58 MM

 PMS 303 C
 PMS 3025 C
 PMS 185 C

PN. - 06FK30-40-A0
 REV. DATE - 2019.10.22
 CREATOR - THINKBOX

HIV/Syphilis Duo
 Bioline™

HIV/Syphilis Duo
 Bioline™

Abbott **HIV/Syphilis Duo**
 Bioline™

HIV/Syphilis Duo
 Bioline™

EN
 Rapid, immunochromatographic assay for the qualitative detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or Treponema pallidum (TP) simultaneously in human serum, plasma or whole blood

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CONTENTS:
 • 25 Test devices with desiccant in individual foil pouches
 • Assay diluent (1 x 4 ml/vial) 
 • 1 Instructions for use

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06FK30-40-A0 © 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Abbott Diagnostics Korea Inc.
 65, Borshagal-ro, Gyeong-gu, Yongin-si,
 Gyeonggi-do, 17099, Republic of Korea
 abbott.com/ko-kr

PRINT AREA

Bioline™ HIV/Syphilis Duo

1°C 30°C / 34°F 86°F  25       **REF 06FK30**

Bioline™ HIV/Syphilis Duo Assay diluent
Proclin™ 300 (5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)
Skin sensitization Category 1

Warning
 H317: May cause an allergic skin reaction
 P272: Avoid breathing dust/fume/gas/mist/vapours/spray
 P273: Contaminated work clothing should not be allowed out of the workplace
 P280: Wear protective gloves/protective clothing/eye protection/face protection
 P302+P352: IF ON SKIN: Wash with plenty of soap and water
 P321: Specific treatment
 P333+P313: If skin irritation or rash occurs: Get medical advice/attention
 P363: Wash contaminated clothing before reuse
 P501: Dispose of contents/container in accordance with local/regional/national/international regulation

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06FK35-40-A0

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

PRINT AREA

Abbott Diagnostics Korea Inc.
 65, Borahae-gil, Gyeongju, Gyeongsang-do, 71099, Republic of Korea
 abbott.com/ko/kr

Bioline™
HIV/Syphilis Duo

1°C - 30°C / 34°F - 86°F

Bioline™ HIV/Syphilis Duo Assay diluent
Proclin™ 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)
Skin sensitization Category 1

Warning

H317: May cause an allergic skin reaction
 P272: Avoid breathing dust/fume/gas/mist/vapours/spray
 P273: Contaminated work clothing should not be allowed out of the workplace
 P280: Wear protective gloves/protective clothing/eye protection/face protection
 P302+P352: IF ON SKIN: Wash with plenty of soap and water
 P321: Specific treatment
 P333+P313: If skin irritation or rash occurs: Get medical advice/attention
 P363: Wash contaminated clothing before reuse
 P501: Dispose of contents/container in accordance with local/regional/national/international regulation

PRODUCT - Bioline™ HIV/Syphilis Duo
25T, Safety Lancet

PN. - 06FK37-40-A0
REV. DATE - 2021. 03. 18
CREATOR - THINKBOX

SIZE - 152 X 115 X 75 MM

-  PMS 303 C
-  PMS 3025 C
-  PMS 185 C

Bioline™
HIV/Syphilis Duo

Bioline™
HIV/Syphilis Duo

Abbott Bioline™
HIV/Syphilis Duo

Abbott

Abbott Bioline™
HIV/Syphilis Duo

EN
Rapid, immunochromatographic assay for the qualitative detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or Treponema pallidum (TP) simultaneously in human serum, plasma or whole blood

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06FK37-40-A0 © 2021 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

PRINT AREA

Abbott Diagnostics Korea Inc.
65, Borahgal-ro, Gyeong-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea
globalpointofcare.abbott

Bioline™
HIV/Syphilis Duo

1°C 34°F | 30°C 86°F | 25 | 1 | 2 | IVD | REF 06FK37

Bioline™ HIV/Syphilis Duo Assay diluent
Prodn™ 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)
Skin sensitization Category 1

- Warning**
- H317: May cause an allergic skin reaction
 - P261: Avoid breathing dust/fume/gas/mist/vapours/spray
 - P273: Contaminated work clothing should not be allowed out of the workplace
 - P280: Wear protective gloves/protective clothing/eye protection/face protection
 - P302+P352: IF ON SKIN: Wash with plenty of soap and water
 - P321: Specific treatment
 - P333+P313: If skin irritation or rash occurs: Get medical advice/attention
 - P363: Wash contaminated clothing before reuse
 - P501: Dispose of contents/container in accordance with local/regional/national/international regulation

PRODUCT - Bioline™ AL pouch 1-30 °C
Non-CE

SIZE - 108 X 32 MM

 PMS 303 C
 PMS 2925 C

CAT NO. - FP26-A0
REV. DATE - 2019.10.28
CREATOR - SCOTT KIM

 **Abbott**

Bioline™
RAPID DIAGNOSTIC TEST

 **Abbott Diagnostics Korea Inc.**
65, Borahagal-ro, Giheung-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea

1°C - 30°C
34°F - 86°F

FP26-A0

ITEM ITEM: Bioline™

LOT LOT NO.:

MFG.DATE (YYYY.MM.DD)
EXP.DATE (YYYY.MM.DD)


11mm

PRODUCT - Bioline™

Assay diluent label

(CHEMICAL HAZARD SYMBOLS)

SIZE - 50 X 16 MM

PMS 303 C

PMS 2925 C

CAT NO. - FB09-A0

REV. DATE - 2019. 11. 12

CREATOR - SCOTT KIM

Abbott **ITEM** HIV/Syphilis Duo

Bioline™
ASSAY DILUENT **VOL** 4ml/vial

IVD 1°C - 30°C **LOT** XXXXXXXXXX

Abbott Diagnostics Korea Inc. FB09-A0 **YYYY.MM.DD**

40005438_Capillary Tube (20ul)25PCS_60×60

20µl Capillary Tube

Code No. : 40005438
LOT No. :
Quantity : 25PCS
MFG Date :

 HLB Life Science Co.,Ltd
101-908, Digital EmpireII, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16681, Republic Korea
 Cmc Medical Devices & Drugs S.L.
C/ Horacio Lengo n18 C.P 29006 Málaga-Spain

(01) 08809899620106

20µl Capillary Tube

Code No. : 40005438
LOT No. :
Quantity : 25PCS
MFG Date :

 HLB Life Science Co.,Ltd
101-908, Digital EmpireII, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16681, Republic Korea
 Cmc Medical Devices & Drugs S.L.
C/ Horacio Lengo n18 C.P 29006 Málaga-Spain

(01) 08809899620106

Capillary tube

Model: Single marked line

Specification: 20µl

FOR PROFESSIONAL USE ONLY

Jiangsu Changfeng Medical Industry Co.,ltd

Address: 189 Tongda Road, Touqiao Town, Yangzhou, China

Tel:+86-13816762056

E-mail: Jacky. cao@changfeng-medical.com

EC

REP

Llins Service & Consulting GmbH

Address: Obere Seegasse 34/2, 69124, Heidelberg, Germany

Email: info@llins-service.com

UDI

(01) XXXXXXXXXXXXXXXX

(11) YYMMDD

(17) XXXXXX

(10) XXXXXXXX

YYYY-MM-DD

YYYY-MM-DD

LOT

XXXXXXXXXX

Made in China

A/0

IVD

PANTONE Reflex Blue C
Shandong
Size:50X40

Sterile Lancet for single use


LOT LOT NO./Type : 2305267303 / 28G
EXP. DATE : 2028-06-04
QTY : 25pcs

STERILE R **CE** 0123 

0100040517004436172880040103205207101

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 60, 20537 Hamburg, Germany

Shandong Lianfa Medical Plastic Products Co., Ltd.
No. 3 Shuangshihai Sanjian Road, 250200 Zhangqiu City, Jinan, Shandong
PEOPLE'S REPUBLIC OF CHINA

(01)16957138602815
 (11)YYMMDD
 (17)YYMMDD
 (10)XXXXXXXX

LOT

EC REP

Shanghai International Holding
 Corp. GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel: 0049-40-2513175
 E-mail: shholding@hotmail.com

**Twist-off lancet (use with lancing device)
 25 pcs (28G I)**

[Instructions]

1. Unscrew the lancing device cap from the body of the lancing device.
2. Insert a sterile lancet into the lancing device firmly.
3. Twist off the protective cap of the lancet.
4. Carefully screw the lancing device cap back to the lancing device.
5. Adjust the puncture depth.
6. Put the lancing device against the side of the finger and press the release button to prick your fingertip.
7. Apply a drop of blood to test strip.
8. Recap lancet and discard it in appropriate container.
9. Place the lancing device cap back to the lancing device.

[Structure] The Twist-off Lancets is consist of needle, main body and protective cap.

[Performance] The device is sterilized by irradiation.

[Intended use] The twist-off lancet (use with lancing device) is a single use, sterile medical device, which is designed for use of micro blood sampling puncture to obtain capillary blood samples from the fingertip.

[Contraindications] None.

[Shelf life] 5 years.

[Storage] The products should be stored in a well-ventilated room and kept away from light, moisture and corrosive gases.

Tianjin Huahong Technology Co., Ltd.

A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone
 (Air Port Industrial Park), 300308 Tianjin, China
 Tel: 0086 22 58775333 400 8587753 Fax: 0086 22 58775353
 E-mail: huahongsales@163.com

[Cautions]


1. Read the instructions carefully before use.
2. Use the lancets within the expiry date.
3. Do not use if protective cap is removed or broken previously.
4. Do not store lancet in lancing device, discard lancet after use, do not re-use.
5. For your safety, please do not use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals.
6. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

 The lancet and lancing device must be compatible with each other. When you use Huahong's sterile lancets with other manufacturers' lancing devices, please be sure that these products are compatible. We recommend that you shall consult a health care professional before use.

STERILE R


CE 0197

Instructions for use :

Twist protective cap, pull off and discard it.

Apply blood collection to the disinfected sampling site and press it to the bottom quickly.

Squeeze the blood collection site to obtain the corresponding blood collection volume.

The sterile lancets after use should be placed in a special recovery container.

(01)16957138604734
 (11)YYMMDD
 (17)YYMMDD
 (10)XXXXXX

Sterile lancets (Safety lancets)

Type : 28G- XXIII Σ_{25T} QTY : 25pcs

STERILE R Irradiation Sterilization **CE 0197**

Tianjin Huahong Technology Co., Ltd.
 A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone
 (Air Port Industrial Park), 300308 Tianjin, China
 Tel: 0086 22 58775333 400 8587753 Fax: 0086 22 58775353
 E-mail: huahongsales@163.com

EC/REP Shanghai International Holding Corp. GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel:0049-40-2513175
 E-mail: shholding@hotmail.com

[Structure] The safety Lancets consist of needle core, button, housing, bottom, spring and protective cap.

[Performance] The device is sterilized by irradiation and the needle tip is sterile.

[Intended use] The safety lancets are single use medical devices, which are designed to collect capillary blood sample.

[Contraindications] None.

[Cautions]

1. Read the instructions carefully before use.
2. Use the lancets within the expiry date.
3. Do not use if protective cap is removed or broken previously.
4. Discard the lancet after use, do not re-use.
5. For your safety, never use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals.
6. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

[Storage]

The products should be stored in a well-ventilated room and kept away from light , moisture and corrosive gases.

We recommend that you shall consult a health care professional before use.

LOT

PANTONE Reflex Blue C
Shandong
Size:50X40

Safety Lancet

LOT LOT NO./Type : 23032692 / 28G
 EXP. DATE : 2028-05-09
 QTY : 25pcs

STERILE R  **0123** 


(010004851700104720000001023032692)

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 60, 20537 Hamburg, Germany

 **Shandong Lianfa Medical Plastic Products Co., Ltd.**
No. 3 Shuangshihai Sanjian Road 250200 Zhangou City, Jinan, Shandong
PEOPLE'S REPUBLIC OF CHINA


FRONT SIDE OF THE BAG

BACK SIDE OF THE BAG

WHITE BACKGROUND

WHITE BACKGROUND

- transparent window where the lancets are visible

UDI code printed directly on the bag during packing process

Instructions for use²

² English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

Abbott

REF 06FK30, 06FK35

Bioline™

HIV/Syphilis Duo

Antibodies to HIV-1/2 and *Treponema pallidum* Test
Test des anticorps du VIH-1/2 et du *Treponema pallidum*
Prueba para anticuerpos de HIV-1/2 y *Treponema pallidum*
Anticorpos para teste de HIV-1/2 e *Treponema pallidum*

INTENDED USE

The Bioline™ HIV/Syphilis Duo test is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 including subtype-O, HIV-2 and Syphilis (*Treponema pallidum*) in human serum, plasma or whole blood. The Bioline™ HIV/Syphilis Duo is intended only for professional use, and for *in vitro* diagnostic use. This test may not be suitable for diagnosis of early infection or blood donation screening. Reactive specimens should be confirmed by a supplemental assay such as ELISA or Western Blot for antibodies to HIV-1/HIV-2 and the treponemal test (e.g. TPPA, TPHA) and non-treponemal test (e.g. RPR and VDRL) for antibodies to *T. pallidum*.

INTRODUCTION

HIV (Human Immunodeficiency Virus) is recognized as the etiologic agent of Acquired Immune Deficiency Syndrome (AIDS). The virus is transmitted by sexual contact, exposure to infected blood, certain body fluids or tissues, and from mother to fetus or child during the perinatal period. HIV-1 has been isolated from patients with AIDS and AIDS related complex, and from healthy persons with high potential risk of developing AIDS. Patients with HIV-2 are found primarily in parts of West Africa. Its course is marked by increasing levels of viral replication and the emergence of more virulent viral strains. HIV-1 and HIV-2 are similar in their morphology, cell tropism, host interaction and generic structure. Serological studies have determined that HIV-1 and HIV-2 have multiple common epitopes in core antigens but much less so in the envelope antigens. This clinical diagnosis of HIV may include the detection of antibodies to HIV-1/2 in human plasma or serum by immunoassay. The presence of HIV can be identified by detection of antibodies to HIV 1/2 in human serum, plasma and whole blood by immunoassay. This advanced assay utilizes recombinant antigens targeted against immunogenic proteins. The major immunoreactive antigens of these proteins are HIV-1 gp41, p24 and HIV-2 gp36.

The spirochetal bacillus *Treponema pallidum* (TP) is the causative agent of the venereal disease syphilis. The presence of TP can be identified by detection of antibodies to TP in human blood by immunoassay. This advanced assay utilizes recombinant antigens targeted against immunogenic TP membrane proteins.⁷

TEST PRINCIPLE

The Bioline™ HIV/Syphilis Duo test is a solid phase immunochromatographic assay for the qualitative, simultaneous detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or *Treponema pallidum* in human serum, plasma or whole blood. The Bioline™ HIV/Syphilis Duo contains a nitrocellulose membrane strip with two pre-coated test line regions. Test line region HIV is pre-coated with recombinant HIV-1 capture antigen (gp41), recombinant HIV-2 capture antigen (gp36) and recombinant HIV-sub O antigen; test line region SYP is pre-coated with recombinant *Treponema pallidum* antigens (17KDa). When specimen and assay diluent is added to

the specimen pad, it moves across the antigen-colloid gold conjugate pad. This mixture continues to migrate through the solid phase to the immobilized antigen on the membrane and forms a visible line. This antigen-antibody-antigen gold particle complex forms with high degree of sensitivity and specificity. If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible HIV and C lines in the window. If antibodies to syphilis are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible SYP and C lines in the window. If antibodies to HIV-1 and/or HIV-2 and TP are undetectable, the antigen-colloid gold flows past the window and one control line appears in the window.

MATERIALS PROVIDED AND ACTIVE INGREDIENTS OF MAIN COMPONENTS

1. The Bioline™ HIV/Syphilis Duo test kit contains the following items to perform the assay:

Catalog No.	06FK30	06FK35
Contents	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 4 ml/vial)  3. 1 Instructions for use	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 4 ml/vial)  3. 25 Capillary pipettes (20 µl), 25 Sterile lancets, 25 Alcohol swabs 4. 1 Instructions for use

2. Active ingredients of main components:

1 test strip includes:	Assay diluent includes:
<ul style="list-style-type: none"> - Gold conjugates: Recombinant HIV-1 antigen – gold colloid (0.029±0.006 µg), Recombinant HIV-2 antigen – gold colloid (0.043±0.009 µg), Recombinant Treponema pallidum antigen – gold colloid (0.017±0.003 µg), Chicken IgY– gold colloid (0.038±0.008 µg) - Test line HIV: Recombinant HIV-1 antigen (gp41) (0.384±0.077 µg), Recombinant HIV-2 antigen (gp36) (0.640±0.128 µg), Recombinant HIV - sub O antigen (0.224±0.045 µg) - Test line SYP: Recombinant Treponema pallidum antigen (17KDa) (0.480±0.096 µg) - Control line: Mouse monoclonal anti-chicken IgY (0.640±0.128 µg) 	<ul style="list-style-type: none"> - 20 mM Tris-HCl buffer (100 mM) - Tricine (100 mM) - Casein (0.2 w/v%) - Proclin™ 300 (0.05 w/v%) - Sucrose (5 w/v%) - Sodium azide (0.01 w/v%)

- Classification according to Regulation (EC) No 1272/2008:

Product identifier	Trade name	Bioline™ HIV/Syphilis Duo Assay diluent		
	Hazardous substance	Proclin™ 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)		
Classification	Skin sensitization Category 1		Signal Word	Warning
Hazard Pictogram			Hazard Statement	H317: May cause an allergic skin reaction
Precautionary statements				
Prevention	P261: Avoid breathing dust/fume/gas/mist/vapours/spray P272: Contaminated work clothing should not be allowed out of the workplace P280: Wear protective gloves/protective clothing/eye protection/face protection			
Response	P302+P352: IF ON SKIN: Wash with plenty of soap and water P321: Specific treatment P333+P313: If skin irritation or rash occurs: Get medical advice/attention P363: Wash contaminated clothing before reuse			
Disposal	P501: Dispose of contents/container in accordance with local/regional/national/international regulation			

MATERIALS REQUIRED BUT NOT PROVIDED

- Micropipette, Protective gloves, Timer, Biohazard container

KIT STORAGE AND STABILITY

- The test kit should be stored at a temperature between 1 °C and 30 °C. Do not freeze the kit or its components. Assay diluent cap should be kept firmly sealed between each use.
- Assay diluent may be opened and resealed for each assay. It is stable until expiration date if kept at 1 - 30 °C.
- The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.

4. Do not use the test kit beyond its expiration date. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

WARNINGS

1. The test devices are for *in vitro* diagnostic use only. Do not reuse the test device.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be proficient.
3. Do not use the pipette by mouth, smoke, drink, eat, apply cosmetics, or handle contact lenses in areas where specimens or kit components are being handled.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
7. Do not mix or interchange different specimens.
8. Do not eat the desiccant from the foil pouch.
9. Avoid splashing or aerosol formation of specimen and assay diluent.
10. Do not mix or interchange components among different lots or those for other products.
11. Do not drink assay diluent.
12. Care should be taken to avoid contamination of the bottle nozzle when dropping assay diluent into the specimen well.
13. The assay diluent contains a proprietary antimicrobial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.
14. The assay diluent contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of these reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
15. Safety data sheet available for professional user on request.

SPECIMEN COLLECTION AND HANDLING

1. Whole blood

[Collection by venipuncture]

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
- If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the blood specimen must be tested within 3 days of refrigeration.
- Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.
- Bring blood specimens to room temperature (15 - 30 °C) prior to use.

[Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip then prick the lateral side of the finger with the sterile lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new 20 µl capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

2. Plasma or serum

- If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8 °C. For storage period longer than 2 weeks, freezing (below -20 °C) is required. Bring plasma or serum specimens to room temperature (15 - 30 °C) prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

[Plasma]

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) then centrifuge the tube to generate a plasma specimen.

[Serum]

- Using venipuncture, draw whole blood into the collection tube (NOT containing anticoagulants) then leave for 30 minutes to allow blood coagulation to occur. Centrifuge the tube to generate a serum specimen.

*Anticoagulants including heparin, EDTA and sodium citrate do not affect the test result. Use of other anticoagulants has not been validated. Their use may affect the test result.


TEST PROCEDURE (REFER TO FIGURE)

1. Allow all kit components and specimens to reach a temperature between 15 °C and 30 °C prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. [Using a micropipette]
Dispense 10 µl of plasma or serum specimen or dispense 20 µl of whole blood specimen into the specimen well marked “S”.


Or,

[Using a capillary pipette]

Dispense 20 µl of drawn whole blood specimen into the specimen well marked “S”.

 **Caution:** Lightly touch the capillary pipette to the specimen pad while dispensing.

4. Dispense 3 drops (approximately 100 µl) of assay diluent into the specimen well marked “S”.

 **Caution:** Do not let bottle nozzle touch device in order to avoid cross-contamination. Hold bottle vertically while dispensing. If you do not hold the bottle vertically, it can lead to inaccurate results. Exactly, 3 drops should be added. Adding more than 3 drops may result in reddish color background or an invalid result.

5. As the test begins to work, you will see purple color move across the result window in the center of the test device.
6. Interpret test results 15 - 20 minutes after adding assay diluent. Do not read after 20 minutes.

 **Caution:** If the test result is not legible after 15 minutes due to high background color, read again later but within 20 minutes of adding the diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.

TEST INTERPRETATION (REFER TO FIGURE)

The characters “HIV”, “SYP” and “C” are printed on the test device. They correspond to the positions of the HIV, Syphilis and control lines in the test window.

- **Nonreactive result:**

The presence of only one purple colored line (C) within the result window indicates a nonreactive result.

- **Reactive result:**

 **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

- 1) HIV-1/2 reactive: The presence of both test line HIV and the control line (C) indicates a reactive result for anti-HIV-1/2.
- 2) Syphilis reactive: The presence of both test line SYP and the control line (C) indicates a reactive result for anti-*Treponema pallidum*.

3) HIV-1/2 and Syphilis reactive: The presence of test line HIV, test line SYP and the control line (C) indicates reactive result for both anti-HIV-1/2 and anti-*Treponema pallidum*.

- **Invalid result:**

Absence of the control line (C) and/or presence of a pink/purple smear in the result window indicate an invalid result. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

TEST LIMITATIONS

1. A reactive result with the Bioline™ HIV/Syphilis Duo test indicates presence of antibodies to *Treponema pallidum*. A reactive result alone is insufficient for diagnosis of Syphilis infection. Results must be interpreted along with other clinical information available to the physician.
2. A reactive result may indicate infection with HIV-1/2. Immunochromatographic testing alone cannot be used to diagnose AIDS. An AIDS diagnosis can only be made in a clinical setting if an individual meets the case definition for AIDS established by the Centers for Disease Control (CDC). Reactive specimens should be confirmed by a supplemental assay such as ELISA or Western Blot test.
3. A nonreactive result does not eliminate the possibility of infection with HIV-1/2 and/or *Treponema pallidum*. The specimen may contain low levels of antibodies that cannot be detected by Bioline™ HIV/Syphilis Duo. If a test result is nonreactive and clinical symptoms persist, additional testing using other clinical methods is recommended.
4. Some known HIV-infected persons taking antiretroviral medication have been shown to produce false negative results when tested by rapid diagnostic tests.^{9,10}
5. Where clinical presentation or other data would suggest an inconsistent test result then the individual should be tested by nucleic acid testing (NAT) technologies immediately and/or retested for antibodies to HIV after more than 21 days since the original testing.

INTERNAL QUALITY CONTROL

The Bioline™ HIV/Syphilis Duo test device has test line HIV, test line SYP and the control line C on the surface of the device. The lines are not visible before applying a specimen. The control line is used for procedural control and shows that the diluent has been applied successfully and that the active ingredients of main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a reactive specimen control.

PERFORMANCE CHARACTERISTICS

Performance Summary

Total	HIV	Syphilis
Diagnostic Sensitivity	99.80 %	90.00 %*
Clinical Sensitivity	99.91 %	99.67 %
Diagnostic Specificity	100.00 %	99.90 %
Clinical Specificity	99.67 %	99.72 %

*Diagnostic Sensitivity on HIV/Syphilis co-infected patient specimens was lower than on clinical specimens not co-infected.

1. Diagnostic sensitivity
Clinical performance studies were performed by two external evaluation centres, the Institute of Tropical Medicine in Belgium and the German Red Cross in Germany.
- 1) Anti-HIV-1 and anti-HIV-2 reactive specimens
400 anti-HIV-1 reactive specimens, among which 40 non-B subtype and 100 anti-HIV-2 reactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 500 reactive specimens, is 99.8%.

	Bioline™ HIV/Syphilis Duo		Sensitivity (95 % CI)
	HIV nonreactive	HIV reactive	
Anti-HIV-1 reactive	1	259	
Anti-HIV-1 reactive non-B subtype	0	40	
Anti-HIV-1 and anti-TP reactive	0	100	
Anti-HIV-2 reactive	0	100	
Total	1	499	

- 2) Anti-*Treponema pallidum* reactive specimens
150 anti-*Treponema pallidum* reactive specimens, among which 100 anti-HIV-1 reactive and 50 anti-HIV-1 nonreactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-*Treponema pallidum* antibody detection, calculated on 150 confirmed reactive specimens is 90.0%.

	Bioline™ HIV/Syphilis Duo		
	Syphilis nonreactive	Syphilis reactive	Sensitivity (95 % CI)
Anti-HIV-1 reactive	13	87	
Anti-HIV-1 nonreactive	2	48	
Total	15	135	90.0 % (84.0 - 94.3 %)

2. Diagnostic specificity

In total 1,000 EDTA plasma and 500 whole blood specimens from blood donors were tested. The specimens originated from 2 collection sites in Germany; Frankfurt and Kassel.

		HIV		Syphilis	
		Nonreactive	Reactive	Nonreactive	Reactive
Blood donors	EDTA Plasma	1,000	0	1,000	0
	Whole blood	500	0	499	1
Specificity (95 % CI)		100 % (99.7 - 100 %)		99.9 % (99.6 - 100 %)	

3. Laboratories in Ghana, Mexico, Laos, Togo, Kenya, and Myanmar participated in the evaluation during 2012–2013. Each site characterized sera using *T. pallidum* particle agglutination assay or *T. pallidum* hemagglutination assay and HIV enzyme immunoassay, Western Blot, and/or HIV antibody rapid tests.

1) Laboratory Performance for Detection of Anti-HIV

Country	Year	n	True Reactive	False Reactive	False Nonreactive	True Nonreactive	Sensitivity Estimate (95 % CI)	Specificity Estimate (95 % CI)
Ghana	2012	400	250	0	0	150	100 % (98.54 - 100 %)	100 % (97.57 - 100 %)
Togo	2013	310	203	0	0	107	100 % (98.20 - 100 %)	100 % (96.61 - 100 %)
Myanmar	2013	245	114	1	0	130	100 % (96.82 - 100 %)	99.24 % (95.82 - 99.98 %)
Kenya	2013	698	345	0	1	352	99.71 % (98.40 - 99.99 %)	100 % (98.96 - 100 %)
Mexico	2013	527	158	0	0	369	100 % (97.69 - 100 %)	100 % (99.01 - 100 %)
Laos	2013	156	53	3	0	100	100 % (93.28 - 100 %)	97.09 % (91.72 - 99.40 %)
Total	-	2,336	1,123	4	1	1,208	99.91 % (99.51 - 100 %)	99.67 % (99.16 - 99.91 %)

2) Laboratory performance for Detection of Anti-*Treponema pallidum*

Country	Year	n	True Reactive	False Reactive	False Nonreactive	True Nonreactive	Sensitivity Estimate (95 % CI)	Specificity Estimate (95 % CI)
Ghana	2012	400	250	1	0	149	100 % (98.54 - 100 %)	99.33 % (96.34 - 99.98 %)
Togo	2013	241	88	1	0	152	100 % (95.89 - 100 %)	99.35 % (96.41 - 99.98 %)
Myanmar	2013	200	74	1	1	124	98.67 % (92.79 - 99.97 %)	99.20 % (95.62 - 99.98 %)
Kenya	2013	698	85	0	0	613	100 % (95.75 - 100 %)	100 % (99.40 - 100 %)
Mexico	2013	414	106	1	1	306	99.07 % (94.90 - 99.98 %)	99.67 % (98.20 - 99.99 %)
Laos	2013	106	6	0	0	100	100 % (54.07 - 100 %)	100 % (96.38 - 100 %)
Total	-	2,059	609	4	2	1,444	99.67 % (98.82 - 99.96 %)	99.72 % (99.29 - 99.92 %)

4. Sensitivity on seroconversion panels

30 commercially available HIV seroconversion panels were tested showing how early detection can be obtained with the Bioline™ HIV/Syphilis Duo test. 20 HIV Ab panels (pertaining 62 early seroconversion samples) were from the Institute of Tropical Medicine and 10 panels (pertaining 12 early seroconversion samples) were from the Paul-Ehrlich-Institute. The performance of the Bioline™ HIV/Syphilis Duo test on seroconversion panels is comparable to other CE marked Anti-HIV-1/2 screening and rapid tests.

5. Analytical specificity

• Cross reactivity

The following 16 potential cross-reacting pathogens had no effect on test results of Bioline™ HIV/Syphilis Duo.

Anti-HBs	EBV	HIV (to SYP line)	Syphilis (to HIV line)
<i>Borrelia burgdorferi</i>	HAV	HTLV	Toxoplasma
Chlamydia	HBsAg	Influenza	<i>Trypanosoma cruzi</i> I / II
CMV	HCV	<i>Plasmodium vivax</i>	<i>Trypanosoma gambiense</i>

• Interfering substances

The following 9 substances had no effect on test results of Bioline™ HIV/Syphilis Duo.

: Pregnant women, High cholesterol (245 mg/dL), High bilirubin (1.4 mg/dL), Rheumatoid factor (28 IU/ml), Lipemic, Hemolyzed, Autoimmune, Alcoholic cirrhosis, Multiparous pregnancy


• Pharmaceutical substances

There was no significant interference with the following 25 drugs on Bioline™ HIV/Syphilis Duo. All drugs were tested at concentrations of 250 µg/ml.

Abacavir	Cholecalciferol	Folic acid	L-ascorbic acid	Pantoprazole
Acetaminophen	Cyclobenzaprine	Hydrochlorothiazide	Magnesium sulfate	Pyrazinamide
Acetylsalicylic acid	Darunavir	Ibuprofen	Metformin	Rifampicin
Amoxicillin	Diclofenac	Iron Chloride	Naproxen	Ritonavir
Aspirin	Ergocalciferol	Isoniazid	Nevirapine	Salicylic acid

6. Reproducibility of the Bioline™ HIV/Syphilis Duo has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

PREPARATION

Open the foil pouch and look for the following:

1. Test device
2. Desiccant

Then, label the device with the patient identifier.


● SP10 / WB20 μ l : Serum 10 μ l or Plasma 10 μ l or Whole blood 20 μ l

□▶●●● : Assay diluent 3 drops


TEST PROCEDURE

I. Blood (by venipuncture), Plasma or Serum specimen


- 1 Using a micropipette draw plasma or serum: **10 μ l** specimen
or
whole blood: **20 μ l** specimen.

- 2** Dispense
10 μl of plasma or serum specimen
or
20 μl of whole blood specimen
into the specimen well marked “S”.

- 3** Dispense **3 drops** (approximately 100 μl) of assay diluent into the specimen well marked “S”. Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.


- 4** Interpret test results 15 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.


II. Blood specimen (with a lancet)


1 Clean the area to be lanced with an alcohol swab.


2 Squeeze the fingertip then prick the lateral side of the finger with the sterile lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.


3 Immerse the open end of a new **20 µl** capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

- 4** Dispense **20 μ l** of drawn whole blood specimen in the specimen well marked “S”. Lightly touch the capillary pipette to the specimen pad while dispensing.

- 5** Dispense **3 drops** (approximately 100 μ l) of assay diluent into the specimen well marked “S”. Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.


- 6** Interpret test results 15 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.

INTERPRETATION

Nonreactive

The presence of only one purple colored line (C) within the result window indicates a nonreactive result.



Reactive

⚠ Caution: The presence of any test line, no matter how faint, the result is considered reactive.



HIV-1/2 Reactive

Presence of both “C” and “HIV” lines





Syphilis Reactive

Presence of both “C” and “SYP” lines



HIV-1/2 and Syphilis Reactive

Presence of “C”, “HIV” and “SYP” lines

Invalid

Absence of the control line (C) and/or presence of a pink/purple smear in the result window indicate an invalid result. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

Product. Disclaimer:

While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:

The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive diagnosis using this product.

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1. Miller JN : Value and limitation of nontreponemal and treponemal tests in the laboratory diagnosis of syphilis. Clin. Obstet Gynecol 18: 191-203, 1975.
2. Lefevre J-C, Bertrand M-A, Baogriaud R : Evaluation of Captia enzyme immunoassays for detection of immunoglobulins G and M to *Treponema pallidum* in syphilis. J Clin Microbiol 28 : 1704-1707, 1990.
3. V.S. Ivanov, Z.K. Suvorova, L.D. Tchikin, A.T. Kozitch and V.T. Ivanov: Effective method for synthetic peptide immobilization that increases the sensitivity and specificity of ELISA procedures Journal of Immunological Methods, 153 (1992) 229-233.
4. Mi Jin Sohn, Young Hae Chong, Ji Eun Chang, Young Ik Lee: Overexpression and simple purification of human immunodeficiency virus-1 gag epitope derived from a recombinant antigen in *E. coli* and its use in ELISA. Journal of Biotechnology 34 (1994) 149-155.
5. Laboratory diagnostics for HIV Infection. Butto S. et al. Centro Nazionale AIDS, Istituto Superiore di Sanita, Rome, Italy. Ann Ist Super Sanita. 2010 Vol. 46 no.1: 24-33.
6. HIV Diagnostic Testing: Evolving Technology and Testing Strategies. Daskalakis D. et al. 2011. New York University Medical Center, USA. Top Antivir Med. 2011 Vol. 19 (1): 18-22.
7. Laboratorial Diagnosis of Syphilis, Syphilis – Recognition, Description and Diagnosis. 2011 Neuzha Satomi Sato et al. (Ed.), ISBN: 978-953-307-554-9, InTech.
8. The laboratory diagnosis of Syphilis. 2005 Ratnam et al. Can J Infect Dis Med Microbiol Vol 16 No 1 January/February 2005.
9. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
10. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.
11. Websites For Peer reviews
- [http://www.ijgo.org/article/S0020-7292\(15\)00201-5/references](http://www.ijgo.org/article/S0020-7292(15)00201-5/references)
- http://sti.bmj.com/content/89/Suppl_1/A56.1
- <http://www.ncbi.nlm.nih.gov/pubmed/25908677>
- <http://www.ncbi.nlm.nih.gov/pubmed/25734088>
- <http://www.ncbi.nlm.nih.gov/pubmed/25339379>
- <http://www.ncbi.nlm.nih.gov/pubmed/25391806>

Abbott Diagnostics Korea Inc.

65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, 17099,
Republic of Korea

Tel : 82-31-899-2800 | Fax : 82-31-899-2840

abbott.com/poct

Date issued : 2019. 10
06FK30/06FK35-04-A0

© 2019 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Abbott

REF 06FK37

Bioline™

HIV/Syphilis Duo

Antibodies to HIV-1/2 and *Treponema pallidum* Test

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TEST PRINCIPLE

The Bioline™ HIV/Syphilis Duo test is a solid phase immunochromatographic assay for the qualitative, simultaneous detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or *Treponema pallidum* in human serum, plasma or whole blood. The Bioline™ HIV/Syphilis Duo contains a nitrocellulose membrane strip with two pre-coated test line regions. Test line region HIV is pre-coated with recombinant HIV-1 capture antigen (gp41), recombinant HIV-2 capture antigen (gp36) and recombinant HIV-sub O antigen; test line region SYP is pre-coated with recombinant *Treponema pallidum* antigens (17KDa). When specimen and assay diluent is added to the specimen pad, it moves across the antigen-colloid gold conjugate pad. This mixture continues to migrate through the solid phase to the immobilized antigen on the membrane and forms a visible line. This antigen-antibody-antigen gold particle complex forms with high degree

of sensitivity and specificity. If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible HIV and C lines in the window. If antibodies to syphilis are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible SYP and C lines in the window. If antibodies to HIV-1 and/or HIV-2 and TP are undetectable, the antigen-colloid gold flows past the window and one control line appears in the window.

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- The Bioline™ HIV/Syphilis Duo test kit contains the following items to perform the assay:
 - 25 Test devices with desiccant in individual foil pouches
 - Assay diluent (1 x 4 ml/vial) ⚠
 - 25 Capillary pipettes (20 µl), 25 Safety lancets, 25 Alcohol swabs
 - 1 Instructions for use
- Active ingredients of main components:

1 test strip includes:	Assay diluent includes:
<ul style="list-style-type: none"> - Gold conjugates: Recombinant HIV-1 antigen – gold colloid (0.029±0.006 µg), Recombinant HIV-2 antigen – gold colloid (0.043±0.009 µg), Recombinant Treponema pallidum antigen – gold colloid (0.017±0.003 µg), Chicken IgY– gold colloid (0.038±0.008 µg) - Test line HIV: Recombinant HIV-1 antigen (gp41) (0.384±0.077 µg), Recombinant HIV-2 antigen (gp36) (0.640±0.128 µg), Recombinant HIV - sub O antigen (0.224±0.045 µg) - Test line SYP: Recombinant Treponema pallidum antigen (17kDa) (0.480±0.096 µg) - Control line: Mouse monoclonal anti-chicken IgY (0.640±0.128 µg) 	<ul style="list-style-type: none"> - 20 mM Tris-HCl buffer (100 mM) - Tricine (100 mM) - Casein (0.2 w/v%) - Proclin™ 300 (0.05 w/v%) - Sucrose (5 w/v%) - Sodium azide (0.01 w/v%)

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[Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new 20 µl capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

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TEST PROCEDURE (REFER TO FIGURE)

1. Allow all kit components and specimens to reach a temperature between 15 °C and 30 °C prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. [Using a micropipette]
Dispense 10 µl of plasma or serum specimen or dispense 20 µl of whole blood specimen into the specimen well marked “S”.
Or,
[Using a capillary pipette]
Dispense 20 µl of drawn whole blood specimen into the specimen well marked “S”.
⚠ **Caution:** Lightly touch the capillary pipette to the specimen pad while dispensing.
4. Dispense 3 drops (approximately 100 µl) of assay diluent into the specimen well marked “S”.
⚠ **Caution:** Do not let bottle nozzle touch device in order to avoid cross-contamination. Hold bottle vertically while dispensing. If you do not hold the bottle vertically, it can lead to inaccurate results. Exactly, 3 drops should be added. Adding more than 3 drops may result in reddish color background or an invalid result.
5. As the test begins to work, you will see purple color move across the result window in the center of the test device.
6. Interpret test results 15 - 20 minutes after adding assay diluent. Do not read after 20 minutes.
⚠ **Caution:** If the test result is not legible after 15 minutes due to high background color, read again later but within 20 minutes of adding the diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.

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- **Reactive result:**

⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

- 1) HIV-1/2 reactive: The presence of both test line HIV and the control line (C) indicates a reactive result for anti-HIV-1/2.
- 2) Syphilis reactive: The presence of both test line SYP and the control line (C) indicates a reactive result for anti-*Treponema pallidum*.
- 3) HIV-1/2 and Syphilis reactive: The presence of test line HIV, test line SYP and the control line (C) indicates reactive result for both anti-HIV-1/2 and anti-*Treponema pallidum*.

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1. Diagnostic sensitivity

Clinical performance studies were performed by two external evaluation centres, the Institute of Tropical Medicine in Belgium and the German Red Cross in Germany.

1) Anti-HIV-1 and anti-HIV-2 reactive specimens

400 anti-HIV-1 reactive specimens, among which 40 non-B subtype and 100 anti-HIV-2 reactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 500 reactive specimens, is 99.8%.

	Bioline™ HIV/Syphilis Duo		Sensitivity (95 % CI)
	HIV nonreactive	HIV reactive	
Anti-HIV-1 reactive	1	259	
Anti-HIV-1 reactive non-B subtype	0	40	
Anti-HIV-1 and anti-TP reactive	0	100	
Anti-HIV-2 reactive	0	100	
Total	1	499	99.8 % (98.9 - 100 %)

2) Anti-*Treponema pallidum* reactive specimens

150 anti-*Treponema pallidum* reactive specimens, among which 100 anti-HIV-1 reactive and 50 anti-HIV-1 nonreactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-*Treponema pallidum* antibody detection, calculated on 150 confirmed reactive specimens is 90.0%.

	Bioline™ HIV/Syphilis Duo		
	Syphilis nonreactive	Syphilis reactive	Sensitivity (95 % CI)
Anti-HIV-1 reactive	13	87	
Anti-HIV-1 nonreactive	2	48	
Total	15	135	90.0 % (84.0 - 94.3 %)

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
- Reproducibility of the Bioline™ HIV/Syphilis Duo has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

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


- 2** Dispense **10 μl** of plasma or serum specimen or **20 μl** of whole blood specimen into the specimen well marked “S”.

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


2 Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.


OR


OR


- 3** Immerse the open end of a new **20 µl** capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

- 4** Dispense **20 µl** of drawn whole blood specimen in the specimen well marked "S". Lightly touch the capillary pipette to the specimen pad while dispensing.

- 5** Dispense **3 drops** (approximately 100 µl) of assay diluent into the specimen well marked "S". Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.

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