WHO Prequalification of Diagnostics Programme PUBLIC ASSESSMENT REPORT

Product: Bioline Malaria Ag P.f
PQDx Reference Number: PQDx 0031-012-01

The Bioline Malaria Ag P.f with product code 05FK50, 05FK53, 05FK51 and 05FK52, manufactured by Abbott Diagnostics Korea Inc, CE-mark regulatory version was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 06 December 2010.

Summary of WHO Prequalification Assessment for the Bioline Malaria Ag P.f

	Date	Outcome
Prequalification listing	06 December 2010	listed
Dossier assessment	02 December 2010	MR
Product performance evaluation	Round 1 (2009) and Round 5 (2014)	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product, for which the WHO has been notified and has undertaken a review. The amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment and change request reference, where applicable.	Date of report amendment
1.0-3.0	Editorial amendments and addition of a product code 05FK53	2 February 2016
4.0	Changes to the labelling of the cassette	2 May 2016
5.0	Addition of product codes 05FK51 and 05FK52	22 December 2016
6.0	Changes to the manufacturer name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc and product name from SD BIOLINE Malaria Ag P.f and SD BIOLINE Malaria Ag P.f POCT to Bioline Malaria P.f	20 August 2020
7.0	Addition of the affiliate company of Abbott Group (Rapid test manufacturing site in China) for device assembly, pouch sealing and packing process.	15 January 2021
8.0	Accessory labelling changes according to accessory manufacturer information change (Manufacturer name	23 March 2026

	<p>change and EU representative address change) (PQC-IVD-2024-0034). Editorial and Instructional updates to the IFU for Bioline Malaria Ag P.f for alignment with the globally distributed product. Changes are limited to administrative changes, changes to correct errors, and clarification of existing information without changing the procedure and do not impact quality, safety or performance (PQC-IVD-2025-0043).</p>	
--	---	--

Intended use:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“The Bioline Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of Malaria Plasmodium falciparum in human whole blood. Bioline Malaria Ag P.f test is intended for professional use, only for an initial screening test and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear”*

Test kit contents:

	25T/kit (product code 05FK50)	1T/kit x 25ea (product code 05FK53)	25T/kit (product code 05FK51)	1T/kit x 25ea (product code 05FK52)
Test cassettes individually packed in foil pouch with a desiccant	25 test devices	25x 1 test device	25 test devices	25x 1 test device
Assay diluent dispensed in plastic bottle	1 x 5 mL/bottle	25x 180 µL/vial	1 x 5 mL/bottle	25 x 180 µL/vial
Specimen transfer devices disposable (5 µL)	25 units of 5 µL	25 pouches each containing	25 units of 5 µL	25 pouches each containing

Lancets (optional) disposable, sterile	25 units	one specimen transfer device (5µL), one lancet, one alcohol swab	25 units (safety lancets)	one specimen transfer device (5 µL), one safety lancet , one alcohol swab
Alcohol swabs (optional) disposable	25 units		25 units	
Instructions for use	1 copy	1 copy	1 copy	1 copy
Summarized instructions for use	N/A	25 copies	N/A	25 copies

Storage:

The test kit must be stored between 1 and 40 °C.

Shelf-life upon manufacture¹:

24 months.

Dossier assessment

Abbott Diagnostics Korea Inc submitted a product dossier for Bioline Malaria Ag P.f as per the *“Instructions for compilation of a product dossier”* (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2).

Based on the product dossier screening and assessment findings, the Bioline Malaria Ag P.f product dossier meets WHO prequalification requirements.

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pqweb/vitro-diagnostics/who-public-inspection-reports>

Product performance evaluation

Bioline Malaria Ag P.f was submitted to Round 1 and Round 5 of WHO Product Testing Programme in 2009 and 2014². The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild type parasite panel, and a mixed panel of *P. falciparum* and non- *P. falciparum* to determine specificity. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

In Round 1, the following results were observed: *P. falciparum* panel detection score (97.5% at 200 parasites/μL), false-positive rates (2.4% for clean negatives, 0.0% for *P. vivax* at 200 parasites/μL, 0.0 % for *P. vivax* at 2000 or 5000 parasites/μL) and invalid rate (0.0%).

In Round 5, the following results were observed: *P. falciparum* panel detection score (95.0% at 200 parasites/μL), false-positive rates (0.0% for clean negatives, 0.0% for *P. vivax* at 200 parasites/μL, 2.9% for *P. vivax* at 2000 or 5000 parasites/μL) and invalid rate (0%), for Bioline Malaria Ag P.f meets the current laboratory evaluation requirements for prequalification.

Summary performance characteristics	Panel detection score	False positive rate (%)		Invalid rate (%)
	200 parasites/μL	200 parasites/μL	Clean negatives	
	Pf	Pv		
Round 1	97.47	0.0	2.38	0.0
Round 5	95.0	0.0	0.0	0.0

² The same product was also submitted to Round 1 Product Testing.

Based on these results, the performance evaluation for Bioline Malaria Ag P.f meets the WHO prequalification requirements.

Labelling review

The labelling submitted for Bioline Malaria Ag P.f was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

Controlled Labelling References

05FK50

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-09-06	05FK50-40-A0
Pouch / Device label	Device pouch label	A0	2025-06-12	FP14-A0
Assay diluent labels	Assay diluent labels	A0	2019-10-25	FB07-A0
Inverted cup	40008361	0	2025-09-25	D10151501
Sterile lancet	P65-008	0	2023-11-10	D10151501
Sterile lancet	40006487	0	2024-10-01	D10151501
Instructions for Use (IFU)	IFU	A2	2020-09-23	05FK50-04-A2

05FK51

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Package	A0	2019-09-06	05FK51-40-A0
Pouch / Device label	Device pouch label	A0	2025-06-12	FP14-A0
Assay diluent labels	Assay diluent labels	A0	2019-10-25	FB07-A0
Inverted cup	40008361	0	2025-09-25	D10151501
Sterile safety lancet	P65-013/000/001	0	2023-11-10	D10151501
IFU	IFU	A2	2020-09-23	05FK51-04-A2

05FK52

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Outer box artwork	A0	2019-09-06	05FK52-40-A0
Pouch / Device label	Device pouch label	A0	2025-06-12	FP14-A0
POCT Pouch	Pouch label	A0	2019-09-09	05FK52-P-A0
Inverted cup*	40005865	N/A	N/A	N/A
Sterile safety lancet*	40005867	N/A	N/A	N/A
Ampule*	P60-009/009	N/A	N/A	N/A
IFU*	Summarized IFU	A2	2021-04-09	05FK52-02-C-A2
IFU	IFU	A2	2020-09-23	05FK52-04-A2
Reference Buffer guide	Malaria Buffer guide	A0	2022-02-07	05FK00-IC-A0


* These items are packaged together in a POCT pouch

05FK53

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	Outer box artwork	A0	9/6/2019	05FK53-40-A0
Pouch / Device label	Device pouch label	A0	2025-06-12	FP14-A0
POCT Pouch	Pouch label	A0	9/9/2019	05FK53-P-A0
Inverted cup*	40005865	N/A	N/A	N/A
Sterile lancet*	40005864	N/A	N/A	N/A
Sterile lancet*	40009512	N/A	N/A	N/A
Ampule*	P60-009/009	N/A	N/A	N/A
IFU	Summarized IFU	A2	4/9/2021	05FK53-02-C-A2
IFU	IFU	A2	9/23/2020	05FK53-04-A2
Reference Buffer guide	Malaria Buffer guide	A0	2/7/2022	05FK00-IC-A0

* These items are packaged together in a POCT pouch

Labels

Malaria Ag P.f
Bioline™

Malaria Ag P.f

**Bioline™
Malaria Ag P.f**


PRINT AREA

**Bioline™
Malaria Ag P.f**

1°C - 40°C / 34°F - 104°F | 25 | 1 | 2 | | 05FK50

**Bioline™
Malaria Ag P.f**

EC REP

MT Promedi Consulting GmbH
Altenhofstrasse 80,
66386 St. Ingbert, Germany


© 2019 Abbott. All rights reserved. All trademarks
referenced are trademarks of either the Abbott group of
companies or their respective owners.

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CONTENTS:
• 25 Test devices with desiccant in individual foil pouch
• Assay diluent (1 x 5 ml/vial)
• 25 Disposable specimen applicators (5 µl), 25 Safety lancets, 25 Alcohol swabs
• 1 Instructions for use

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

EC REP

MT Promed Consulting GmbH
Altenhofstrasse 80,
66396 St. Ingbert, Germany

© 2019 Abbott. All rights reserved. All trademarks referred are trademarks of either the Abbott group of companies or their respective owners.

Abbott Diagnostics Korea Inc.
65, Borahagall-ro, Gyeong-gu, Yongin-si,
Gyeonggi-do, 31099, Republic of Korea
abbott.com/peet

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CONTENTS :
• 25 Test devices with desiccant in individual foil pouch
• 25 Assay diluents
• 25 Disposable specimen applicators (5 µl),
• 25 Safety lancets, 25 Alcohol swabs
• 25 Summarized instructions for use
• 1 Instructions for use

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


© 2019 Abbott. All rights reserved. All trademarks
 referenced are trademarks of either the Abbott group of
 companies or their respective owners.

PRODUCT - Bioline™ Malaria Ag P.f
SIZE - 220 X 130 X 70MM
PN - 05FK53-40-A0
REV. DATE - 2019.08.29
CREATOR - THINKBOX
PMS 303 C
PMS 2925 C

PRODUCT - Bioline™ Malaria Ag P.f
(Safety lancet, POCT spec) Outer pouch

SIZE - 145 x 50mm

 PMS 303 C
 PMS 2925 C

CAT NO. - 05FK52-P-A0
REV. DATE - 2019. 08. 29
CREATOR - THINKBOX

Bioline™
Malaria Ag P.f

Malaria Antigen P.f (HRP2) Rapid Test

CONTENTS : Test device, Assay diluent, Disposable specimen applicator (5 µl), Safety lancet, Alcohol swab, Summarized instructions for use

 **Abbott Diagnostics Korea Inc.**
65, Borahagel-ro, Giheung-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea

1°C / 34°F  40°C / 104°F

REF 05FK52

05FK52-P-A0

PRODUCT - Bioline™ Malaria Ag P.f
(POCT spec) Outer pouch

SIZE - 145 x 50mm

MATERIAL -

 PMS 303 C
 PMS 2925 C

CAT NO. - 05FK53-P-A0
REV. DATE - 2019. 08. 29
CREATOR - THINKBOX

Bioline™
Malaria Ag P.f

Malaria Antigen P.f (HRP2) Rapid Test

CONTENTS : Test device, Assay diluent, Disposable specimen applicator (5 µl), Sterile lancet, Alcohol swab, Summarized instructions for use

 **Abbott Diagnostics Korea Inc.**
65, Borahagel-ro, Giheung-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea

1°C / 34°F
40°C / 104°F

REF 05FK53

05FK53-P-A0

PRODUCT - Bioline™ AL pouch 1-40 °C

SIZE - 108 X 32 MM

PMS 303 C

PMS 2925 C

CAT NO. - FP14-A0

REV. DATE - 2019.10.28

CREATOR - SCOTT KIM

Bioline™
RAPID DIAGNOSTIC TEST

 Abbott Diagnostics Korea Inc.
65, Borahagal-ro, Giheung-gu, Yongin-si,
Gyeonggi-do, 17099, Republic of Korea

1°C - 40°C
34°F - 104°F

FP14-A0

ITEM ITEM: Bioline™

LOT LOT NO.:

MFG.DATE (YYYY.MM.DD)
EXP.DATE (YYYY.MM.DD)

11mm

PANTONE Reflex Blue C
Shandong
Size:50X40

Sterile Lancet for single use

LOT LOT NO./Type : 2305267303 / 28G
EXP. DATE : 2028-06-04
QTY : 25pcs

STERILE R **CE** 0123 

0100040517004436172880040103205207101

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 60, 20537 Hamburg, Germany

Shandong Lianfa Medical Plastic Products Co., Ltd.
No. 3 Shuangshihai Sanjian Road 250200 Zhangqiu City, Jinan, Shandong
PEOPLE'S REPUBLIC OF CHINA

PANTONE Reflex Blue C
Shandong
Size:50X40

Sterile Lancet for single use

LOT LOT NO./Type : 2305267303 / 28G
EXP. DATE : 2028-06-04
QTY : 25pcs

STERILE R **CE** 0123 

0100040517004436172880040103205207101

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 60, 20537 Hamburg, Germany

Shandong Lianfa Medical Plastic Products Co., Ltd.
No. 3 Shuangshihai Sanjian Road, 250200 Zhangqiu City, Jinan, Shandong
PEOPLE'S REPUBLIC OF CHINA

XXXXXXXXX_Inverted Cup(80mm)XXPCS_55×60

Inverted Cup 80mm (5µl)

Code No. : XXXXXXXX
LOT No. : XXXXXXXXXX
Quantity : XXPCS
MFG Date : XXXX.XX.XX

HLB Life Science Co.,Ltd
101-908, Digital EmpireII, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16681, Republic Korea

(0 1) 0 8 8 0 9 8 9 9 6 2 X X X X

EC REP Cmc Medical Devices & Drugs S.L.
C/ Horacio Lengo n18 C.P 29006 Málaga-Spain

Inverted Cup 80mm (5µl)

Code No. : XXXXXXXX
LOT No. : XXXXXXXXXX
Quantity : XXPCS
MFG Date : XXXX.XX.XX


HLB Life Science Co.,Ltd
101-908, Digital EmpireII, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16681, Republic Korea

(0 1) 0 8 8 0 9 8 9 9 6 2 X X X X


EC REP Cmc Medical Devices & Drugs S.L.
C/ Horacio Lengo n18 C.P 29006 Málaga-Spain


	10mm	33.75mm	67.5mm	33.75mm	10mm
10mm	<div style="display: flex; justify-content: space-between;"> <div style="width: 20%;">   <p>(01)16957138602815 (11)YYMMDD (17)YYMMDD (10)XXXXXXXX</p> <p>LOT</p>   <p>EC REP Shanghai International Holding Corp. GmbH (Europe) Eiffestrasse 80, 20537 Hamburg, Germany Tel: 0049-40-2513175 E-mail: shholding@hotmail.com</p> </div> <div style="width: 50%; text-align: center;"> <p>Twist-off lancet (use with lancet device) 25 pcs (28G I)</p> <p>[Instructions]</p> <ol style="list-style-type: none"> Unscrew the lancet device cap from the body of the lancet device. Insert a sterile lancet into the lancet device firmly. Twist off the protective cap of the lancet. Carefully screw the lancet device cap back to the lancet device. Adjust the puncture depth. Put the lancet device against the side of the finger and press the release button to prick your fingertip. Apply a drop of blood to test strip. Recap lancet and discard it in appropriate container. Place the lancet device cap back to the lancet device. <p>[Structure] The Twist-off Lancets is consist of needle, main body and protective cap.</p> <p>[Performance] The device is sterilized by irradiation.</p> <p>[Intended use] The twist-off lancet (use with lancet device) is a single use, sterile medical device, which is designed for use of micro blood sampling puncture to obtain capillary blood samples from the fingertip.</p> <p>[Contraindications] None.</p> <p>[Shelf life] 5 years.</p> <p>[Storage] The products should be stored in a well-ventilated room and kept away from light, moisture and corrosive gases.</p> <p>Tianjin Huahong Technology Co., Ltd. A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone (Air Port Industrial Park), 300308 Tianjin, China Tel: 0086 22 58775333 400 8587753 Fax: 0086 22 58775353 E-mail: huahongsales@163.com</p> </div> <div style="width: 25%;"> <p>[Cautions]</p> <ol style="list-style-type: none"> Read the instructions carefully before use. Use the lancets within the expiry date. Do not use if protective cap is removed or broken previously. Do not store lancet in lancet device, discard lancet after use, do not re-use. For your safety, please do not use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals. In the event of any serious adverse event, please report to the manufacturer and the competent authority. <p> The lancet and lancet device must be compatible with each other. When you use Huahong's sterile lancets with other manufacturers' lancet devices, please be sure that these products are compatible. We recommend that you shall consult a health care professional before use.</p> <p>STERILE R</p>    <p> CE 0197</p> </div> </div>				
75mm					
10mm					

Instructions for use :

Twist protective cap, pull off and discard it.

Apply blood collection to the disinfected sampling site and press it to the bottom quickly.

Squeeze the blood collection site to obtain the corresponding blood collection volume.

The sterile lancets after use should be placed in a special recovery container.

(01)16957138604734
 (11)YYMMDD
 (17)YYMMDD
 (10)XXXXXX

Sterile lancets (Safety lancets)

Type : 28G- XXIII Σ_{25T} QTY : 25pcs

STERILE R Irradiation Sterilization **CE 0197**

Tianjin Huahong Technology Co., Ltd.
 A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone
 (Air Port Industrial Park), 300308 Tianjin, China
 Tel: 0086 22 58775333 400 8587753 Fax: 0086 22 58775353
 E-mail: huahongsales@163.com

EC/REP Shanghai International Holding Corp. GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel:0049-40-2513175
 E-mail: shholding@hotmail.com

[Structure] The safety Lancets consist of needle core, button, housing, bottom, spring and protective cap.

[Performance] The device is sterilized by irradiation and the needle tip is sterile.

[Intended use] The safety lancets are single use medical devices, which are designed to collect capillary blood sample.

[Contraindications] None.

[Cautions]

1. Read the instructions carefully before use.
2. Use the lancets within the expiry date.
3. Do not use if protective cap is removed or broken previously.
4. Discard the lancet after use, do not re-use.
5. For your safety, never use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals.
6. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

[Storage]

The products should be stored in a well-ventilated room and kept away from light , moisture and corrosive gases.

We recommend that you shall consult a health care professional before use.

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
Instructions for use³

³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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ENGLISH

About the test

[Introduction]
Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is caused by a parasite that is transmitted from one human to another through the bite of infected *Anopheles* mosquitoes. There are four types of malaria that can infect humans: *Plasmodium falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and are released in another form: merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 200 million clinical cases and 0.5 million malaria-caused deaths per year. At present, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

[Test principle]
Bioline™ Malaria Ag P.f test device contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to HRP2 of *P. falciparum* on test line "P.f". The mouse monoclonal antibodies specific to HRP2 of *P. falciparum* colloidal gold conjugate reacts with the Malaria *Plasmodium falciparum* antigen in the specimen. They move along the membrane chromatographically to the test region "P.f" and form a visible line as the antibody-antigen-antibody gold particle complex with high degree of sensitivity and specificity.

[Intended Use]
Bioline™ Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of *Malaria Plasmodium falciparum* in human whole blood. Bioline™ Malaria Ag P.f test is intended only for professional use and for in vitro diagnostic use. Reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear.

Materials provided and active ingredients of main component

- Bioline™ Malaria Ag P.f test kit contains the following items to perform the assay:
 - 25 Test devices with desiccant in individual foil pouches
 - Assay diluent (1 x 5 ml/vial)
 - 25 Disposable specimen applicators (5 µl), 25 Sterile lancets, 25 Alcohol swabs
 - 1 Instructions for use
 - Active ingredients of main component
- Test strip ingredients : Gold conjugate : Mouse monoclonal antibodies specific to P.f HRP2 - gold colloid (0.10±0.02 µg), Test line P.f: Mouse monoclonal antibodies specific to P.f HRP2 (0.5±0.1 µg), Control line : Goat anti- mouse IgG (0.10±0.2 µg)
- Assay diluent includes: Bovine serum albumin (q.s.), 200 mM Phosphate Saline Buffer (q.s.), TRITON™ X-100 (q.s.), Sodium azide (q.s.)

Materials required but not provided

- Micropipette, Protective gloves, Timer, Biohazard container

Kit storage and stability

- The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the kit or its components. Assay diluent cap should be kept firmly sealed between each use.
 - 25 Test devices with desiccant in individual foil pouches
 - Assay diluent (1 x 5 ml/vial)
 - 25 Disposable specimen applicators (5 µl), 25 Sterile lancets, 25 Alcohol swabs
 - 1 Instructions for use
 - Active ingredients of main component
- Test strip ingredients : Gold conjugate : Mouse monoclonal antibodies specific to P.f HRP2 - gold colloid (0.10±0.02 µg), Test line P.f: Mouse monoclonal antibodies specific to P.f HRP2 (0.5±0.1 µg), Control line : Goat anti- mouse IgG (0.10±0.2 µg)
- Assay diluent includes: Bovine serum albumin (q.s.), 200 mM Phosphate Saline Buffer (q.s.), TRITON™ X-100 (q.s.), Sodium azide (q.s.)

Warnings

- For in vitro diagnostic use only. Do not reuse the test device and kit components.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
- Avoid splashing or aerosol formation of specimen and assay diluent. Clean up spills thoroughly using an appropriate disinfectant.
- Decomtaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. lancet, specimen applicator, test device) in a biohazard container as if they were infectious waste.
- Do not mix or interchange different lot's specimens.
- Do not mix reagent of different lots or those for other products.
- Do not pipette by mouth.
- Do not eat the desiccant in the foil pouch.
- Do not drink assay diluent.
- The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety procedures are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

Specimen collection and storage
[Collection by venipuncture]

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
- If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the blood specimen must be tested within 3 days.
- Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.

No. of parasites / µl of blood	Microscopy (No. of positive)	Bioline™ Malaria Ag P.f (No. of positive)	Sensitivity (%)
1-50	16	15	93.8 %
51-100	35	35	100 %
101-500	96	96	100 %
501-1,000	81	81	100 %
1,001-5,000	79	79	100 %
>5,000	61	61	100 %
Total	368	367	99.7 %

Performance characteristics

- Internal evaluation of Bioline™ Malaria Ag P.f test
 - A. Bioline™ Malaria Ag P.f test for the detection of *P. falciparum* malaria according to the level of parasite (Table 1)
 - B. Sensitivity and specificity of Bioline™ Malaria Ag P.f test. (Table 2)
- Reproducibility of Bioline™ Malaria Ag P.f has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

Product Disclaimer :
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning :
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.

Reference method (Microscopic examination)	Positive	Negative	Total
P.f Positive	367	1	368
Negative	1	199	200
Sensitivity (95 % CI)	99.7 % (98.5 - 100 %)		
Specificity (95 % CI)	99.5 % (97.2 - 99.9 %)		

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
PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

- EN** Open the package and look for the following:
 1. Test device with desiccant in individual foil pouch
 2. Assay diluent
 3. Disposable specimen applicator (5 µl)
 (Capillary pipette or Disposable inverted cup)
 4. Sterile lancet
 5. Alcohol swab
 6. Summarized instructions for use
 7. Instructions for use

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








ENGLISH

Intended use

The Bioline™ Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of Malaria *Plasmodium falciparum* in human whole blood.

REF 05FK52

Materials provided

- Test device
- Assay diluent
- Inverted cup (5 µl)
- Safety lancet
- Alcohol swab

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the kit or its components.
Note : When stored at refrigerator, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to the test. Do not open the pouch whilst components come to room temperature.
2. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
3. Do not use test kit beyond its expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.
6. When transporting or storing the kit, avoid exposure to high temperature (above 45 °C) for a period longer than 1 week.

Specimen collection using a lancet

1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the end of the fingertip and pierce with a safety lancet provided.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. With a disposable inverted cup (5 µl) provided, dip the circular end of an inverted cup into the blood specimen and carefully place the circular end of the inverted cup into the round specimen well.

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© 2021 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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



- 4 **EN** Read the test result after 15 minutes.
CAUTION: Do not read test results after 30 minutes.; reading after 30 minutes can yield erroneous results.

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
Positive / Positif

- EN** **CAUTION:** The presence of any test line, no matter how faint, the result is considered positive.

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
Invalid / Non valide

- EN** If the control line "C" is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

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ENGLISH

Intended Use

The Bioline™ Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of Malaria *Plasmodium falciparum* in human whole blood.

REF 05FK53

Materials Provided

- Test device
- Assay diluent
- Inverted cup (5 µl)
- Sterile lancet
- Alcohol swab

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the kit or its components.
Note : When stored at refrigerator, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to the test. Do not open the pouch whilst components come to room temperature.
2. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
3. Do not use test kit beyond its expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.
6. When transporting or storing the kit, avoid exposure to high temperature (above 45 °C) for a period longer than 1 week.

Specimen collection using a lancet

1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the end of the fingertip and pierce with a sterile lancet provided.
3. Wipe away the first drop of blood with alcohol swab.
4. With a disposable inverted cup (5 µl) provided, dip the circular end of an inverted cup into the blood specimen.

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