WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Determine Antenatal Care Panel WHO reference number: PQDx 13226-013-00

The Determine Antenatal Care Panel, with product code CPK-038, manufactured by Abbott Diagnostics Medical Co., Ltd, Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 10 July 2025.

Summary of the WHO prequalification assessment for the Determine Antenatal Care Panel

	Date	Outcome
Prequalification listing	10 July 2025	listed
Dossier assessment	2 July 2025	MR
Site inspection of the quality management system	31 March 2023	MR
Product performance evaluation	Waived ¹	MR

MR: Meets Requirements

Intended use

According to the intended use claim from Abbott Diagnostics Medical Co., Ltd., “*Determine Antenatal Care Panel is composed of three in vitro tests to be used for pregnant women (12 years of age or older) for the qualitative detection of HIV, Hepatitis B Surface Antigen, and Treponema pallidum Syphilis antibodies in human capillary whole blood. The Determine Antenatal Care Panel is intended for use by health-care professionals in a laboratory and non-laboratory (point of care) settings. These tests are not intended for screening blood donors. These tests are not automated.*

Determine Syphilis TP (Test Strip 1)

Determine Syphilis TP is an in vitro visually read, qualitative immunoassay to be used for pregnant women (12 years of age or older) for the detection of antibodies (IgG and IgM) to Treponema pallidum, which is the bacteria that causes syphilis infection, in human capillary whole blood. The test is intended as an aid in the diagnosis of syphilis in pregnant women.

¹ The performance evaluation of the Determine Antenatal Care Panel was waived as the individual test strips included in the panel—Determine Syphilis TP, Determine HIV Early Detect, and Determine HBsAg 2—have each undergone separate performance evaluations and are WHO prequalified.

Determine HIV Early Detect (Test Strip 2) Determine HIV Early Detect is an in vitro, visually read, qualitative immunoassay to be used for pregnant women (12 years of age or older) for the detection of antibodies (Ab) to HIV-1 and HIV-2 and the detection of non-immunocomplexed (free) HIV-1 p24 antigen (Ag) in human capillary whole blood. Therefore, any reactivity on 1) the Ab bar alone 2) the Ag bar alone or 3) both Ab and Ag bars simultaneously, is considered a reactive result suggestive of infection with HIV. This test is intended as an aid in the diagnosis of HIV-1/HIV-2 infection in pregnant women.

Determine HBsAg 2 (Test Strip 3) Determine HBsAg 2 is an in vitro visually read, qualitative immunoassay to be used for pregnant women (12 years of age or older) for the detection of Hepatitis B Surface Antigen (HBsAg) in human capillary whole blood. The test is intended as an aid in the diagnosis of HBV infection through detection of HBsAg in pregnant women."

Assay description

According to the claim of the assay description from Abbott Diagnostics Medical Co., Ltd.,

"Determine Syphilis TP (Test Strip 1)

This test is an immunochromatographic test for the qualitative detection of antibodies to Treponema pallidum antigens.

A sample is added to the sample pad. The sample migrates through the conjugate pad and mixes with the Treponema pallidum antigen-selenium colloid conjugate. This mixture continues to migrate through the solid phase to the immobilized Treponema pallidum antigens at the patient window site.

If antibodies to Treponema pallidum are present in the sample, the antibodies bind to the Treponema pallidum antigen-selenium colloid and the Treponema pallidum antigen at the patient window, forming a red bar at the patient window site.

If antibodies to Treponema pallidum are not present, the Treponema pallidum antigen-selenium colloid flows past the patient window site, and no red bar is formed at the patient window site.

To ensure assay validity, a procedural control bar is incorporated in the assay device at the control window.

Determine HIV Early Detect (Test Strip 2)

This is an immunochromatographic test for the qualitative detection of free HIV-1 p24 antigen and antibodies to HIV-1 and HIV-2. The sample mixes with the biotinylated anti-p24 antibodies and selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens, synthetic HIV-2 peptide and anti-p24 mouse monoclonal antibody. This mixture migrates through the solid phase to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides at the Antibody (Ab) window, immobilized avidin at the Antigen (Ag) window.

If antibodies to HIV-1, HIV-2 or HIV-1 group O are present in the sample, the antibodies bind to the selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides, forming a red bar at the Ab window site. If

antibodies to HIV-1 and HIV-2 are absent, the selenium colloid-conjugates flow past the Ab window and no red bar is formed at the Ab window site.

If free HIV-1 p24 antigen is present in the sample, the antigen binds to the biotinylated anti-p24 antibodies and the selenium colloid-conjugate coated with anti-p24 mouse monoclonal antibody. This complex binds to the immobilized avidin forming a red bar at the Ag window site. If HIV-1 p24 antigen is not present, both the biotinylated anti-p24 antibodies and selenium colloid-conjugate flow past the Ag window and no red bar is formed at the Ag window site. To ensure assay validity, a procedural control bar is incorporated in the assay device at the control window.

Determine HBsAg 2 (Test Strip 3)

This test is an immunochromatographic test for the qualitative detection of Hepatitis B Surface Antigen (HBsAg). The sample mixes with the biotinylated anti-HBsAg mouse monoclonal antibodies and black particles coated with anti-HBsAg mouse monoclonal antibodies. This mixture migrates along the solid phase to the immobilized avidin at the patient bar. If HBsAg is present in the sample, the antigen binds to the biotinylated anti-HBsAg mouse monoclonal antibodies and the black particles coated with anti-HBsAg mouse monoclonal antibodies. This complex binds to the immobilized avidin forming a black bar on the test strip. If HBsAg is not present, black particles flow past the patient bar and no black bar is formed on the test strip. To ensure the assay validity, a procedural control bar is incorporated in the assay device. The red control bar will appear.”

Test kit contents:

Each Kit contains:

- 10 Determine Antenatal Care Panel foil pouches. Each foil pouch contains:
 - 1 Test Card with 3 test strips attached (1 Determine Syphilis TP Test Strip, 1 Determine HIV Early Detect Test Strip and, 1 Determine HBsAg 2 Test Strip)
 - 1 Sample Collection Tube
 - 1 Dropper
 - 1 Safety Lancet (sterile)
 - 1 Alcohol Pad (sterile)
 - 1 Gauze (sterile)
 - 1 Plaster (sterile)
 - 1 Desiccant
- 1 Instructions for Use (IFU)
- 1 Quick Reference Instructions (QRI)
- 1 Chase Buffer (M2683)

NOTE: The Determine Antenatal Care Panel includes one bottle (2.5 mL) of Chase Buffer (M2683), sufficient for the 10 panels provided. If additional Chase Buffer is needed, it can be ordered using catalogue number 7D2243 (white capped bottle). The 7D2243 Chase Buffer is compatible with Determine Antenatal Care Panel, Determine HIV Early Detect, Determine HBsAg 2, and Determine Syphilis TP products.

Items required but not provided:

- Timer
- Sharps container
- Biohazard waste container
- Protective gloves
- Pen or pencil

Storage:

Determine Antenatal Care Panel foil pouches and Chase Buffer must be stored at 2 °C to 30 °C.

Shelf-life upon manufacture:

14 months.

Warnings/limitations:

Refer to the current version of the manufacturer's instructions for use attached to this public report.

Prioritisation for prequalification

Based on the established criteria, the Determine Antenatal Care Panel was given priority for the WHO prequalification assessment.

Product dossier assessment

Abbott Diagnostics Medical Co., Ltd. submitted a product dossier for the Determine Antenatal Care Panel as per the "Instructions for compilation of a product dossier" (PQDx_018). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO. The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 2 July 2025.

Based on the product dossier screening and assessment findings, the Determine Antenatal Care Panel meets WHO prequalification requirements.

Manufacturing site inspection

At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality. Routine inspections of the Manufacturing site will be conducted with copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the

information on the most current assessment performed at a manufacturing site for in vitro diagnostic products and summarises the assessment findings.

<https://extranet.who.int/pgweb/vitro-diagnostics/who-public-inspection-reports>

All published WHOPIRs are with the agreement of the manufacturer.

Based on the site inspection and corrective action plan review, the quality management system for the Determine Antenatal Care Panel meets WHO prequalification requirements.

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




4. Instructions for use²

²WHO assessed the English versions of the IFU. The manufacturer is responsible for ensuring the correct translation into other languages.

Abbott

Determine™
Antenatal Care Panel

Determine™ Syphilis TP

Determine™ HIV Early Detect

Determine™ HBsAg 2

INSTRUCTIONS FOR USE

MODE D'EMPLOI

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REF **CPK-038**

Scan for video instructions.

NE

DETERMINE™ ANTENATAL CARE PANEL

Read the Instructions for Use in its entirety prior to testing.

Table of Contents

. PRODUCT NAMES AND INTENDED USE.....	4
. SUMMARY AND EXPLANATION OF THE TESTS.....	4
. BIOLOGICAL PRINCIPLES OF THE TESTS	5
. MATERIALS.....	6
. KIT STORAGE AND STABILITY	7
. WARNINGS.....	7
. TEST PREPARATION	8
. SAMPLE COLLECTION HANDLING.....	8
. TEST PROCEDURE	9
. TEST INTERPRETATION	16
. TEST LIMITATIONS	21
. INTERNAL QUALITY CONTROL.....	22
. PERFORMANCE CHARACTERISTICS.....	23
. REFERENCE DOCUMENTS.....	34
. ADVICE LINE.....	35
. GLOSSARY OF SYMBOLS.....	35

PRODUCT NAMES AND INTENDED USE

Determine™ Antenatal Care Panel is composed of three *in vitro* tests to be used for pregnant women (12 years of age or older) for the qualitative detection of HIV, Hepatitis B Surface Antigen, and *Treponema pallidum* Syphilis antibodies in human capillary whole blood. The Determine™ Antenatal Care Panel is intended for use by health-care professionals in a laboratory and non-laboratory (point of care) settings. These tests are not intended for screening blood donors. These tests are not automated.

Determine™ Syphilis TP (Test Strip 1)

Determine™ Syphilis TP is an *in vitro* visually read, qualitative immunoassay to be used for pregnant women (12 years of age or older) for the detection of antibodies (IgG and IgM) to *Treponema pallidum*, which is the bacteria that causes syphilis infection, in human capillary whole blood. The test is intended as an aid in the diagnosis of syphilis in pregnant women.

Determine™ HIV Early Detect (Test Strip 2)

Determine™ HIV Early Detect is an *in vitro*, visually read, qualitative immunoassay to be used for pregnant women (12 years of age or older) for the detection of antibodies (Ab) to HIV-1 and HIV-2 and the detection of non-immunocomplexed (free) HIV-1 p24 antigen (Ag) in human capillary whole blood. Therefore, any reactivity on 1) the Ab bar alone 2) the Ag bar alone or 3) both Ab and Ag bars simultaneously, is considered a reactive result suggestive of infection with HIV. This test is intended as an aid in the diagnosis of HIV-1/HIV-2 infection in pregnant women.

Determine™ HBsAg 2 (Test Strip 3)

Determine™ HBsAg 2 is an *in vitro* visually read, qualitative immunoassay to be used for pregnant women (12 years of age or older) for the detection of Hepatitis B Surface Antigen (HBsAg) in human capillary whole blood. The test is intended as an aid in the diagnosis of HBV infection through detection of HBsAg in pregnant women.

SUMMARY AND EXPLANATION OF THE TESTS

Determine™ Syphilis TP (Test Strip 1)

Syphilis is caused by infection with the bacterium *Treponema pallidum*¹ which can be transmitted congenitally or by sexual contact. The disease can evolve into a latent phase in which syphilis is clinically inapparent. Serologic tests (nontreponemal specific and treponemal specific) are currently the primary method for syphilis diagnosis and management. Nontreponemal tests (VDRL, RPR, etc.) are generally used for screening, and treponemal tests (TPHA/TPPA, FTA-ABS, etc.) are used as confirmatory tests.

Determine™ HIV Early Detect (Test Strip 2)

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes a depletion of CD4 helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The multiplication of the HIV in the infected cells leads to cell rupture and the release of HIV virus particles, which are first detected in the form of HIV RNA and next in the form of HIV antigen^{2,3}. This is followed by production of specific antibodies to either HIV-1 or HIV-2^{4,5,6}. HIV antigen may become undetectable at this time because of the formation of antibody-antigen complexes.

Determine™ HBsAg 2 (Test Strip 3)

Hepatitis B virus (HBV), is a DNA virus transmitted percutaneously, sexually and perinatally. Worldwide it causes deaths from cirrhosis, liver failure and hepatocellular carcinoma⁷. HBsAg seropositivity is the first serological marker to appear in acute HBV infection⁸.

BIOLOGICAL PRINCIPLES OF THE TESTS

Determine™ Syphilis TP (Test Strip 1)

This test is an immunochromatographic test for the qualitative detection of antibodies to *Treponema pallidum* antigens.

A sample is added to the sample pad. The sample migrates through the conjugate pad and mixes with the *Treponema pallidum* antigen-selenium colloid conjugate. This mixture continues to migrate through the solid phase to the immobilized *Treponema pallidum* antigens at the patient window site. If antibodies to *Treponema pallidum* are present in the sample, the antibodies bind to the *Treponema pallidum* antigen-selenium colloid and the *Treponema pallidum* antigen at the patient window, forming a red bar at the patient window site.

If antibodies to *Treponema pallidum* are not present, the *Treponema pallidum* antigen-selenium colloid flows past the patient window site, and no red bar is formed at the patient window site. To ensure assay validity, a procedural control bar is incorporated in the assay device at the control window.

Determine™ HIV Early Detect (Test Strip 2)

This is an immunochromatographic test for the qualitative detection of free HIV-1 p24 antigen and antibodies to HIV-1 and HIV-2. The sample mixes with the biotinylated anti-p24 antibodies and selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens, synthetic HIV-2 peptide and anti-p24 mouse monoclonal antibody. This mixture migrates through the solid phase to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides at the Antibody (Ab) window, immobilized avidin at the Antigen (Ag) window.

If antibodies to HIV-1, HIV-2 or HIV-1 group O are present in the sample, the antibodies bind to the selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides, forming a red bar at the Ab window site. If antibodies to HIV-1 and HIV-2 are absent, the selenium colloid-conjugates flow past the Ab window and no red bar is formed at the Ab window site.

If free HIV-1 p24 antigen is present in the sample, the antigen binds to the biotinylated anti-p24 antibodies and the selenium colloid-conjugate coated with anti-p24 mouse monoclonal antibody. This complex binds to the immobilized avidin forming a red bar at the Ag window site. If HIV-1 p24 antigen is not present, both the biotinylated anti-p24 antibodies and selenium colloid-conjugate flow past the Ag window and no red bar is formed at the Ag window site. To ensure assay validity, a procedural control bar is incorporated in the assay device at the control window.

Determine™ HBsAg 2 (Test Strip 3)

This test is an immunochromatographic test for the qualitative detection of Hepatitis B Surface Antigen (HBsAg). The sample mixes with the biotinylated anti-HBsAg mouse monoclonal antibodies and black particles coated with anti-HBsAg mouse monoclonal antibodies. This mixture migrates along the solid phase to the immobilized avidin at the patient bar. If HBsAg is present in the sample, the antigen binds to the biotinylated anti-HBsAg mouse monoclonal antibodies and the black particles coated with anti-HBsAg mouse monoclonal antibodies. This complex binds to the immobilized avidin forming a black bar on the test strip. If HBsAg is not present, black particles flow past the patient bar and no black bar is formed on the test strip. To ensure the assay validity, a procedural control bar is incorporated in the assay device. The red control bar will appear.

MATERIALS

MATERIALS PROVIDED

Each Kit contains:

- 10 Determine™ Antenatal Care Panel foil pouches. Each foil pouch contains:
 - 1 Test Card with 3 test strips attached (1 Determine™ Syphilis TP Test Strip, 1 Determine™ HIV Early Detect Test Strip and 1 Determine™ HBsAg 2 Test Strip)
 - 1 Sample Collection Tube
 - 1 Dropper
 - 1 Safety Lancet (sterile)
 - 1 Alcohol Pad (sterile)
 - 1 Gauze (sterile)
 - 1 Plaster (sterile)
 - 1 Desiccant
- 1 Instructions for Use (IFU)
- 1 Quick Reference Instructions (QRI)
- 1 Chase Buffer (M2683)

NOTE 1:

The Determine™ Antenatal Care Panel includes one bottle (2.5 mL) of Chase Buffer (M2683), sufficient for the 10 panels provided. If additional Chase Buffer is needed, it can be ordered using catalogue number 7D2243 (white capped bottle). The 7D2243 Chase Buffer is compatible with Determine™ Antenatal Care Panel, Determine™ HIV Early Detect, Determine™ HBsAg 2, and Determine™ Syphilis TP products.

NOTE 2:

Desiccant absorbs moisture during storage. Discard the Desiccant after opening foil pouch. All components in the foil pouch are for single use. Do not reuse. Check sterile components. Do not use if damaged.

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Sharps container
- Biohazard waste container
- Protective gloves
- Pen or pencil

Active Ingredients of Main Components

Determine™ HIV Early Detect: 1 test strip coated with HIV-1/2 recombinant antigen, synthetic peptides, anti-p24 antibodies and avidin

Determine™ Syphilis TP: 1 test strip coated with *Treponema pallidum* antigen

Determine™ HBsAg 2: 1 test strip coated with anti-HBsAg mouse monoclonal antibodies

Chase Buffer is a phosphate buffer preserved with antimicrobial agents

KIT STORAGE AND STABILITY

1. Determine™ Antenatal Care Panel foil pouches and Chase Buffer should be stored at a temperature between 2 °C to 30 °C until expiration date. Do not expose any of the kit contents to freezing or high temperatures. The operating temperature range of the kit is 18 °C to 30 °C.
NOTE: When stored in a refrigerator, allow all kit components to reach room temperature (18 °C to 30 °C) for a minimum of 30 minutes prior to performing the test. **DO NOT OPEN** the foil pouch while components come to room temperature.
2. Chase Buffer cap should be kept firmly sealed between each use to avoid evaporation or spillage.
3. Perform the tests immediately after removing the foil cover, to avoid exposure to humidity/ moisture.
4. The shelf life of the kit is indicated on the outer packaging. Do not use the test kit beyond the expiration date.
5. Do not use the test kit if the foil pouch is damaged or the seal is broken.
6. Determine™ Antenatal Care Panel foil pouch is not re-sealable. Testing materials cannot be stored in an opened Determine™ Antenatal Care Panel foil pouch. Once a foil pouch is opened, the test materials should either be used immediately or discarded.

WARNINGS

1. For *in vitro* diagnostic use only.
2. All components in the foil pouch are for single use. Do not reuse.
3. The instructions must be followed exactly to achieve accurate results.
4. Running the test in high temperature/low humidity may affect the appearance of the patient bar. If the test strip is partially dried and it is difficult to read the result at 20 - 25 minutes, the test should be repeated using a new sample and a new test and result read at 20 minutes. When the test strip is partially dried, it appears as mixed white spot and grayish area.
5. This product contains *Treponema pallidum* sourced from potentially infectious components. No known test method can offer complete assurance that products derived from microorganisms will not transmit infection. Therefore, it is recommended that the test device be considered potentially infectious and handled with appropriate biosafety practices.^{9,10}
6. Do not eat, smoke, apply cosmetics or handle contact lenses where these materials are handled.
7. Wear protective gloves while handling the samples and wash hands thoroughly afterwards.
8. Avoid splashing or aerosol formation of the sample and Chase Buffer.
9. Use new gloves for each patient.
10. When collecting the blood sample, do not squeeze or apply strong repetitive pressure (milking) to the site directly; this may result in hemolysis or tissue-fluid contamination of the specimen.
11. Ensure sample is applied to all 3 tests within 10 minutes of sample collection, to avoid coagulation. A delay of sample application between tests may affect the accuracy of the results.
12. Dispose of all samples, tested devices and potentially contaminated materials (i.e., Sample Collection Tube, Gauze) in a biohazard container.
13. Clean and disinfect all spills of samples or reagents using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.¹¹
14. Dispose of Safety Lancet in a sharps container immediately after use.^{9,10,11,12,13}
15. Do not mix or interchange different samples.
16. Use only the disposable Sample Collection Tube and Dropper provided.

TEST PREPARATION

- Please read the Instructions for Use in its entirety, prior to use.
- All kit components must be brought to between 18°C to 30°C prior to testing.

PREPARE WORK AREA

1. Ensure work area is clean.
2. Remove Test Card from the foil pouch.
3. Remove the following items from the foil pouch and place in the work area:
 - Sample Collection Tube
 - Dropper
 - Safety Lancet (sterile)
 - Alcohol Pad (sterile)
 - Gauze (sterile)
 - Plaster (sterile)
 - Desiccant
4. Place the Test Card on a clean, flat surface.
5. Label or write the patient identification on the top white area of the Test Card.
6. Place the Chase Buffer in the work area.
7. Place a timing device beside the work area.

NOTE:

Desiccant absorbs moisture during storage. Discard the Desiccant after opening foil pouch. All components in the foil pouch are for single use. Do not reuse. Check sterile components. Do not use if damaged.

SAMPLE COLLECTION HANDLING


Precautions^{11,12}

1. All blood samples are potentially infectious. Respect universal precautions.^{11,13}
2. Wear disposable gloves and protective clothing during the entire procedure.
3. Discard used Safety Lancets into a biohazard sharps container.
4. Open the Sample Collection Tube for ease of sample collection, prior to using the Safety Lancet.
5. Discard lab consumables and used test strips into biohazard waste container in accordance with local regulations.
6. Collect biohazardous waste in appropriate container. Do not mix sharps with other waste.


TEST PROCEDURE

Perform Fingerstick^{13,14}

- 1 To ensure that 250 μ L of blood can be collected, ask the patient to warm their hands by massaging to increase blood flow.


- 2 Choose the fingertip of the middle or ring finger (whichever is the least callused).

- 3 Disinfect fingertip with Alcohol Pad; allow to air dry.

- 4 Position the hand palm-side up and place clean finger on the table.


TEST PROCEDURE

Perform Fingersick^{13,14}

- 5** **TWIST** off and discard Safety Lancet cap.


NOTE: The cap is the sterile barrier, do not use if cap is damaged or previously removed. Use a new Safety Lancet for each person.


- 6** Place the Safety Lancet off-center on the fingertip.

Firmly press the Safety Lancet against the finger to activate the Safety Lancet to prick finger.


Dispose of the Safety Lancet in an appropriate biohazard sharps container after use.

- 7** Wipe away the first drop of blood with the Gauze.


- 8** Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times while avoiding the site of the puncture.

TEST PROCEDURE


Perform Fingerstick^{13,14}

- 9** Collect the blood into the Sample Collection Tube.
Fill the Sample Collection Tube with blood up to the fill line (250 μ L).


CAUTION: Proceed promptly through steps 10 to 12.

- 10** Ask the patient to cover the puncture with Gauze.


- 11** Close the cap on the Sample Collection Tube.

NOTE: Replace the cap on the Sample Collection Tube and press down firmly until you hear a click. Ensure Sample Collection Tube is fully closed.

- 12** Immediately invert the Sample Collection Tube 10 times to thoroughly mix blood.

CAUTION: Sample must be applied to all 3 test strips within 10 minutes of sample collection as outlined in section 'Running the Test', to avoid coagulation. A delay of sample application between tests may affect the accuracy of the results.

TEST PROCEDURE


How to use Dropper

⚠ CAUTION: Do not perform these steps immediately. This section is an instruction on 'How to use Dropper' for carrying out steps 3, 4 and 5 in the section 'Running the Test'.

- A** Confirm the location of the line on the Dropper.


- B** Gently squeeze the Dropper and insert into the Sample Collection Tube.

- C** Slowly release pressure to draw the blood up to the fill line of the Dropper.


NOTE: If blood is collected above or below the line, it may cause incorrect or invalid results.

- D** Gently squeeze the Dropper to add the **full blood sample** (50 μ L) to the test strip.

Touch the Dropper to the blood on the test strip to **extract all the blood** from the Dropper.


Keep squeezing the Dropper as you remove the Dropper away from the test strip.


TEST PROCEDURE

Running the Test


- 1 Carefully remove the protective foil covers from all 3 test strips (Syphilis, HIV and HBsAg2).
NOTE: Peel the foil from the top down.
 Avoid pressing down on the individual test strips.
⚠ CAUTION: After removing the protective foils, start the tests immediately, no later than 30 minutes.

- 2 Remove the cap from the Sample Collection Tube.

- 3 Using the Dropper, add the **full blood sample** (50 µL) to Syphilis test strip.
NOTE: Reference section 'How to use Dropper'.


- 4 Using the same Dropper, add the **full blood sample** (50 µL) to HIV test strip.
NOTE: Reference section 'How to use Dropper'.


TEST PROCEDURE

Running the Test

- 5** Using the same Dropper, add the **full blood sample** (50 μL) to HBsAg 2 test strip.
NOTE: Reference section 'How to use Dropper'.


- 6** Start timer for 20 minutes.

- 7** Wait until the blood is absorbed and there is no droplet on the sample pad (approximately 1 minute).


- 8** Start with the Syphilis test, add **1 drop of Chase Buffer** (25 μL) on the sample pad of **each Test Strip** (Syphilis, HIV & HBsAg 2). Hold the bottle vertically to dispense the correct drop volume.

⚠ CAUTION: Keep tip of Chase Buffer a **minimum of 1 cm** away from the sample pad to avoid contamination.


If blood sample splashes the Chase Buffer nozzle, discard the bottle of Chase Buffer.


TEST PROCEDURE

Running the Test


- 9 Close the cap after using the Chase Buffer every time.

- 10 Discard the Sample Collection Tube and Dropper in the biohazard container.

- 11 After the finger has stopped bleeding, remove Gauze and apply Plaster.

TEST INTERPRETATION


WAIT 20 minutes before reading results.
Read results between 20 and 25 minutes.

Do NOT read results after 25 minutes
Results after 25 minutes may not be accurate.

How to read the test results:


If a line is visible, **no matter how faint**, interpret the line as visible.

LOOK very carefully for ANY hint of a line.

NOTE: If it is difficult to read the test result due to reddish background, it is recommended to repeat the test using a new sample and a new test.

INTERPRETATION OF SYPHILIS RESULTS (Test Strip 1)

Control Bar

The Control Bar is to confirm that the test is working.

Patient Bar

The Patient Bar indicates if the Syphilis antibody has been found.


NEGATIVE (Control Bar ONLY)

NEG
⊖

If the Control Bar is present, **no matter how faint**, and no Patient Bar is visible, the test is showing a **NEGATIVE** result.

INVALID (No Control Bar)

Invalid
⊗

If no Control Bar is present, the test did not work and is **INVALID**.

Repeat the test using a new sample and a new test.

If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call the Advice Line.

POSITIVE (Control Bar AND Patient Bar)


POS
⊕


A Control Bar and a Patient Bar, **no matter how faint**, is a **POSITIVE** result.

TEST INTERPRETATION

INTERPRETATION OF HIV Early Detect RESULTS (Test Strip 2)


NEGATIVE (Control Bar ONLY)

If the Control Bar is present, **no matter how faint**, and no Patient Bars are visible, the test is showing a **NEGATIVE** result.

INVALID (No Control Bar)

If no Control Bar is present, the test did not work and is **INVALID**.

Repeat the test using a new sample and a new test.

If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call the Advice Line.

POSITIVE (Control Bar AND One or Two Patient Bar/s)

POS
+

ANTIGEN (p24) POSITIVE (Control and AG Bars)

A Control Bar and an AG Bar, **no matter how faint**, is an **ANTIGEN POSITIVE** result.

The presence of **ONLY** an AG Bar (and Control Bar) suggests that the infection is at an early stage.

Follow-up testing may be suggested in order to track the expected future detection of antibodies.

POS
+

ANTIBODY POSITIVE (Control and AB Bars)

A Control Bar and an AB Bar, **no matter how faint**, is an **ANTIBODY POSITIVE** result.

POS
+

ANTIGEN (p24) AND ANTIBODY POSITIVE (Control Bar, AG Bar and AB Bar)

A Control Bar, an AG Bar and an AB Bar, **no matter how faint**, is a **POSITIVE** result.

TEST INTERPRETATION

INTERPRETATION OF HBsAg 2 RESULTS (Test Strip 3)


Control Bar

The Control Bar is to confirm that the test is working.

Patient Bar


The Patient Bar indicates if the Hepatitis B Surface Antigen has been found.

NEGATIVE (Control Bar ONLY)

If the red Control Bar is present, **no matter how faint**, and no Patient Bars are visible, the test is showing a **NEGATIVE** result.

INVALID (No Control Bar)


If no Control Bar is present, the test did not work and is **INVALID**.

Repeat the test using a new sample and a new test.

If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call the Advice Line.

POSITIVE (Control Bar AND Patient Bar)

A **red** Control Bar and a **black** Patient Bar, **no matter how faint**, is a **POSITIVE** result.

TEST LIMITATIONS

The test procedure, precautions and interpretation of results must be followed precisely.

1. Capillary whole blood samples collected using a sample collection method other than what is provided in this kit have not been validated for use with Determine™ Antenatal Care Panel and may give incorrect results.
2. Extreme environmental conditions, such as high temperatures or freezing may adversely affect test performance. Refer to KIT STORAGE AND STABILITY section for details.

a. Determine™ Syphilis TP (Test Strip 1)

1. The test is designed to detect antibodies to *Treponema pallidum* in human capillary whole blood. Other body fluids or pooled samples may not give accurate results and should not be used.
2. The intensity of the patient bar does not necessarily correlate to the titer of antibody in the sample.
3. Positive results for Syphilis TP may require confirmatory or supplemental testing before a clinical diagnosis or treatment decision is made.
4. The test detects antibodies to *Treponema pallidum*, which may persist for life and does not differentiate between current, past, or treated infection; therefore, follow-up testing with a non-treponemal assay (e.g., RPR or VDRL) is recommended.
5. No test provides absolute assurance that a sample does not contain low levels of antibodies to *Treponema pallidum*, such as those present at a very early stage of infection, or antibodies with low reactivity to the *Treponema pallidum* antigens. Therefore, a negative result at any time does not preclude the possibility of exposure to infection with syphilis.
6. Biotin treatment higher than 200 mg per day may lead to decreased control bar intensity. The concentration of biotin up to 750 ng/mL in serum did not impact the control bar. There is no impact to the patient bar by biotin.

b. Determine™ HIV Early Detect (Test Strip 2)

1. This test is designed to detect antibodies (Ab) to HIV-1 and/or HIV-2 and non-immunocomplexed (free) HIV-1 p24 antigen (Ag) in human capillary whole blood samples. Other body fluids or pooled samples may not give accurate results and should not be used.
2. The intensity of the Ab and Ag bars does not necessarily correlate to the titer of antibody and antigen in the sample, respectively.
3. A reactive result for antibodies to HIV 1/2 combined with a non-reactive result for HIV-1 p24 antigen does not preclude the possibility of acute HIV infection.
4. Reactivity for both antibody and antigen bars or for antigen reactivity only may be confirmed by 4th generation enzyme immunoassay (EIA) or nucleic acid test (NAT) where the IVD is indicated as an aid to diagnosis. Antibody reactivity can also be confirmed by using another HIV Ab Rapid Diagnostic Test. If non-reactive in subsequent testing, retest at least 14 days after the original test is recommended. Follow local guideline for confirmatory methods, if applicable.
5. No test provides absolute assurance that a sample does not contain low levels of HIV-1 p24 antigen and/or antibodies to HIV-1 and HIV-2 such as those present at a very early stage or late stage of HIV infection.

6. A non-reactive result for both antibodies to HIV-1/2 and HIV-1 p24 antigen does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2 viruses.
7. The absence of Ag bar may occur when all p24 antigen is bound by antibodies. When high levels of antibodies against the p24 antigen are present in the blood after seroconversion, the antibodies tend to bind to the antigens, forming immunocomplexes. Determine™ HIV Early Detect detects only non-immunocomplexed (free) antigens; it does not detect immunocomplexed (bound) antigens.
8. Some known HIV-infected persons taking antiretroviral medication have been shown to produce false negative results when tested by rapid diagnostic tests.^{15,16,17}
9. Where clinical presentation is inconsistent with the non-reactive Determine™ HIV Early Detect test result, the individual should be re-tested for antibodies to HIV with Determine™ HIV Early Detect at least 14 days after the original test, or with a NAT assay that is indicated for diagnosis.
10. Biotin treatment higher than 20 mg per day may lead to decreased Ag bar intensity. Biotin concentrations up to 200 ng/mL in serum or whole blood did not impact the Ag sensitivity. There is no impact to Ab bar by biotin.

c. Determine™ HBsAg 2 (Test Strip 3)

1. The test is designed to detect Hepatitis B Surface Antigen (HBsAg) in human capillary whole blood. Other body fluids or pooled samples may not give accurate results and should not be used.
2. The intensity of the patient bar does not necessarily correlate with the titer of antigen in the sample.
3. Positive results for HBsAg may require confirmatory or supplemental testing before a clinical diagnosis or treatment decision is made.
4. For HBsAg 2, while a positive result reflects current hepatitis B infection, further testing (e.g., HBV DNA, HBeAg, liver function tests) is necessary to determine the stage and infectivity of the disease, particularly for pregnant women, where timely intervention is critical.
5. No test provides absolute assurance that a sample does not contain low levels of HBsAg such as those present at a very early stage of infection.
6. A non-reactive result at any time does not preclude the possibility of exposure to or infection with HBV.
7. The use of anticoagulants other than EDTA for human capillary whole blood samples has not been validated for use with Determine™ HBsAg 2 and may give incorrect results.
8. Biotin treatment higher than 5 mg per day may lead to decreased patient bar intensity. Biotin concentrations up to 50 ng/mL in serum did not impact the sensitivity.
9. If Hepatitis B antiviral prophylaxis is initiated in pregnancy following a positive result on Determine™ HBsAg 2 test, confirmatory testing is recommended after delivery to assess long-term treatment eligibility.¹⁸

INTERNAL QUALITY CONTROL

Determine™ Syphilis TP (Test Strip 1)

To ensure assay validity, a procedural control is incorporated in the device and is labeled “Control”. If the control bar does not appear by assay completion, the test result is invalid. Repeat the test using a new sample and a new test. A visible control bar confirms a lateral flow through the membrane.

Determine™ HIV Early Detect (Test Strip 2)

To ensure assay validity, a procedural control is incorporated in the device and is labeled “C”. If the control bar does not appear by assay completion, the test result is invalid. Repeat the test using a new sample and a new test. A visible control bar confirms a lateral flow through the membrane.

Determine™ HBsAg 2 (Test Strip 3)

To ensure assay validity, a procedural control system is incorporated in the device. This red control bar will appear in the window. If the control bar does not appear by assay completion, the test result is invalid. Repeat the test using a new sample and a new test. A visible control bar confirms a lateral flow through the membrane.

PERFORMANCE CHARACTERISTICS

Performance data is an aggregate of whole blood (venipuncture and fingerstick), serum and plasma specimens data.

Clinical Performance

Determine™ Syphilis TP (Test Strip 1)

SENSITIVITY AND SPECIFICITY

Serum specimens from 266 cases of syphilis and 1133 cases of non-syphilis from three clinical sites in Japan were tested by Determine Syphilis TP (Table 1).

Table 1: Sensitivity and Specificity of Determine™ Syphilis TP for Detection of Syphilis

	Population	Number of Specimens Tested	Determine™ Syphilis TP
Sensitivity (95% CI)	Syphilis	266	99.62% (265/266) (97.92% - 99.99%)
Specificity (95% CI)	Non-syphilis	1133	99.74% (1130/1133) (99.23% - 99.95%)

Whole blood specimens with paired serum and plasma from 108 cases of syphilis and 101 cases of non-syphilis from three clinical sites in Uganda were tested with Determine Syphilis TP (Table 2).

Table 2: Sensitivity and Specificity of Determine™ Syphilis TP in Whole Blood (fingerstick and venipuncture), Serum and Plasma matched specimens

Population		Number of Specimens Tested	Determine™ Syphilis TP			
			Whole Blood (fingerstick)	Whole Blood (venipuncture)	Plasma	Serum
Sensitivity (95% CI)	Syphilis	108	93.5% (101/108) (87.1-97.4%)	98.1% (106/108) (93.5-99.8%)	98.1% (106/108) (93.5-99.8%)	98.1% (106/108) (93.5-99.8%)
Specificity (95% CI)	Non-syphilis	101	99.0% (100/101) (94.6-100%)	99.0% (100/101) (94.6-100%)	99.0% (100/101) (94.6-100%)	99.0% (100/101) (94.6-100%)

Determine™ HIV Early Detect (Test Strip 2)

SENSITIVITY

Sensitivity was evaluated by testing confirmed HIV antibody positive specimens, commercial seroconversion panels and HIV-1 subtype virus panels.

HIV Antibody Positive Specimens

A total of 663 confirmed HIV positive serum or plasma specimens were tested (Table 3). The antibody diagnostic sensitivity of Determine™ HIV Early Detect on this population of specimens was calculated to be 100%.

Table 3: HIV Antibody Positive Specimens

Types	Number of Specimens Tested	Reactive by Determine™ HIV Early Detect	Sensitivity (95% CI)
HIV-1	421	421	100% (99.1-100%)
HIV-1 non B subtypes	42	42	100% (91.6-100%)
HIV-2	100	100	100% (96.3-100%)
HIV*	100	100	100% (96.3-100%)
Total	663	663	100% (99.4-100%)

* Specimens collected in Europe

SPECIFICITY

A total of 2509 confirmed HIV negative serum or plasma specimens were tested with the Determine™ HIV Early Detect and the specificity was determined for the antibody (Ab) bar and for the antigen (Ag) bar (Table 4). The specificity was 99.92% for the antibody (Ab) bar and 99.72% for the antigen (Ag) bar.

Table 4: Specificity of Determine HIV Early Detect

Population	Number of Tested Specimens	Determine™ HIV Early Detect					
		Ab bar		Ag bar		Combined Ab and Ag bar	
		Reactive	Non-reactive	Reactive	Non-reactive	Reactive	Non-reactive
Seronegative specimens*	1749	0	1749	2	1747	2	1747
Hospitalized specimens	218	1	217	2	216	3	215
Pregnant women	206	0	206	0	206	0	206
Potentially cross-reacting specimens**	336	1	335	3	333	4	332
Total	2509	2	2507	7	2502	9	2500
Specificity (95% CI)		99.92% (99.71-99.98%)		99.72% (99.43-99.86%)		99.64% (99.32-99.81%)	

* Including specimens collected in Europe (300), USA (1299) and Africa (150)

** Disease states other than HIV, and potentially interfering substances: Rheumatoid factor (1), antinuclear antibody, systemic lupus erythematosus, high cholesterol, high total protein, high IgM (1), human anti mouse IgG, other infections (HBV, HCV, acute HAV (1), HTLV I or II, CMV, Toxo

IgG, Syphilis, HSV 1/2, varicella zoster virus, influenza A and B, malaria, visceral leishmaniasis, and acute or chronic EBV), Flu vaccinated patients, measles (1) and dialysis patient specimen. Numbers in () indicates number of specimens that showed false positive results.

SAMPLE TYPE

1. Antibody sensitivity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

Multiple (matched) specimens: Seropositive specimens from a total of 50 individuals from Africa were tested with the Determine™ HIV Early Detect (Table 5). Multiple (matched) specimens, 50 serum specimens, 50 plasma specimens, 50 whole blood (venipuncture) specimens and 50 whole blood (fingerstick) specimens were obtained from the donors in Africa. The results obtained from all specimen matrices showed 100% correlation, demonstrating that Determine™ HIV Early Detect gives identical results for these types of specimen matrices.

Table 5: Antibody sensitivity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

No. of matched specimens tested	Type of Specimens and No. of Reactive by Determine™ HIV Early Detect (Ab detection)				Correlation between matrices (95% CI)
	Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick	
50	50	50	50	50	100% (92.9-100%)

2. Specificity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

Multiple (matched) specimens: Seronegative specimens from a total of 50 individuals from Africa were tested with the Determine™ HIV Early Detect (Table 6). Multiple (matched) specimens, 50 serum specimens, 50 plasma specimens, 50 whole blood (venipuncture) specimens and 50 whole blood (fingerstick) specimens were obtained from the donors in Africa. The results obtained from all specimen matrices showed 100% correlation, demonstrating that Determine™ HIV Early Detect gives identical results for these types of specimen matrices.

Table 6: Specificity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

No. of matched specimens tested	Type of Specimens and No. of Non-reactive by Determine™ HIV Early Detect				Correlation between matrices (95% CI)
	Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick	
50	50	50	50	50	100% (92.9-100%)

Determine™ HBsAg 2 (Test Strip 3)

SENSITIVITY

Sensitivity was evaluated by testing confirmed HBsAg positive specimens, commercial seroconversion panels, HBV genotype and serotype panels.

HBsAg Positive Specimens

A total of 437 confirmed HBsAg positive specimens were tested with Determine™ HBsAg 2 and a commercially available rapid HBsAg test (Table 7). The sensitivity (95% CI) of Determine™ HBsAg 2 on this population of specimens was calculated to be 98.4% (96.7-99.4%). The non-reactive specimens with Determine™ HBsAg 2 were 0.06, 0.06, 0.09, 0.10, 0.13, 0.33 and 1.29 IU/mL by a quantitative HBsAg test kit.

Table 7: HBsAg Positive Specimens

Types	Number of Specimens Tested	Reactive by Determine™ HBsAg 2	Reactive by alternative rapid HBsAg Test
HBsAg positive specimens	369	362	294
HBV genotype determined specimens*	51	51	43
HBsAg serotype determined specimens**	17	17	17
Plasma sensitivity (95% CI)	369	98.6% (364/369) (96.9-99.6%)	81.3% (300/369) (76.9-85.2%)
Serum sensitivity (95% CI)	68	97.1% (66/68) (89.8-99.6%)	79.4% (54/68) (67.9-88.3%)
Total Sensitivity (95% CI)	437	98.4% (430/437) (96.7-99.4%)	81.0% (354/437) (77.0-84.6%)

* Genotypes: A, B, C, D, E, F, G, H, A/E, A/F, D/E, D/F and D/G

** Serotypes: adw2, adr, ayw2, ayw3 and ayw4 including 15 specimens of # 6100/09 PEI Hepatitis B virus genotype panel

SPECIFICITY

A total of 1,736 confirmed HBsAg negative serum or plasma specimens were tested with the Determine™ HBsAg 2 and the specificity (95% CI) was determined (Table 8). The specificity was 99.5% (99.1-99.8%).

Table 8: HBsAg Negative Specimens

Population*	Determine™ HBsAg 2		Rapid HBsAg Test	
	total	Non-reactive	total	Non-reactive
Seronegative specimens	1027	1027	1027	1027
Clinical specimens	213	209	213	210
Pregnant women	211	210	206	205
Potentially cross-reacting specimens**	285	282	204	202
Plasma specificity (95% CI)	542	100% (542/542) (99.3-100%)	522	100% (522/522) (99.3-100%)
Serum specificity (95% CI)	1194	99.3% (1186/1194) (98.7-99.7%)	1128	99.5% (1122/1128) (98.8-99.8%)

Population*	Determine™ HBsAg 2		Rapid HBsAg Test	
	total	Non-reactive	total	Non-reactive
Total Specificity (95% CI)	1736	99.5% (1728/1736) (99.1-99.8%)	1650	99.6% (1644/1650) (99.2-99.8%)

* Including specimens collected in USA (1240) and Africa (20)

**Disease states other than HBV and potentially interfering substances: Rheumatoid factor (20/20), antinuclear antibody (15/15), systemic lupus erythematosus (15/15), high cholesterol (9/9), high bilirubin (10/10), high total protein (3/3), high IgM (10/10), high IgG (4/4), human anti-mouse IgG (20/20), other infections - Hepatitis A virus (6/6), Hepatitis C virus (15/15), Hepatitis E virus (5/5), HTLV I (9/9), Cytomegalovirus (10/10), Toxoplasma IgG (10/10), Syphilis (8/8), Herpes simplex virus ½ (10/10), Epstein-Barr virus (9/10), Flu vaccinated patients (10/10), non-virus liver disease (19/20), Measles (5/6), Influenza A (5/5), Influenza B (5/5), Varicella zoster virus (5/5), HIV-1 (5/5), HIV-2 (5/5), Malaria (5/5), Visceral Leishmaniasis (5/5), Gonorrhoea (5/5), Trichomonas (5/5), Schistosomiasis (5/5), and dialysis patient specimens (10/10). The following diseases have not been tested: multiple blood transfusions, sickle-cell disease, yellow fever, human African trypanosomiasis and yellow fever vaccination.

SPECIMEN TYPE

All specimen matrices (serum, plasma, whole blood venipuncture and whole blood fingerstick) were tested.

Table 9: HBsAg sensitivity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

Population	No. of Matched Specimens Tested		Type of Specimens and No. of reactive or non-reactive by Determine™ HBsAg 2			
			Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick**
HBsAg positive	145	Reactive	142	143	141	141
		Non-Reactive	3	2	4	4
		Sensitivity (95% CI)	97.9% (94.1-99.6%)	98.6% (95.1-99.8%)	97.2%*** (93.1-99.2%)	97.2%*** (93.1-99.2%)
HBsAg negative	203	Reactive	1	0	0 <3**>	0 <3**>
		Non-Reactive	202	202*	203 <200**>	203 <200**>
		Specificity (95% CI)	99.5% (97.3-100%)	100% (98.2-100%)	100% <98.5**> (98.2-100%)	100% <98.5**> (98.2-100%)

The results obtained from all specimen matrices showed correlation, except for the following. One positive specimen (0.27 IU/mL) was reactive with plasma but non-reactive with serum, venipuncture whole blood and fingerstick. One positive specimen (0.14 IU/mL) was reactive with serum and plasma but non-reactive with venipuncture whole blood and fingerstick. Two seropositive specimens (0.08 IU/mL and 0.14 IU/mL) were non-reactive with all specimen matrices at 15 minutes. At 30 minutes, results for serum and plasma were reactive.

One negative specimen was non-reactive with plasma, venipuncture whole blood and fingerstick but reactive with serum. Three negative venipuncture whole blood specimens were reactive. Three negative fingerstick specimens were reactive.

EDTA Capillary Tubes were used for 225 specimens. MICROSAFE® Tubes were used for 123 specimens. Prospective multiple (matched) specimens of whole blood, serum and plasma from 348 individuals (145 HBsAg positive) in European countries were tested with Determine™ HBsAg 2 (Table 9).

* There was no plasma result for 1 study subject.

** Whole blood initially reactive results were adjudicated utilizing the serum or plasma results; all were non-reactive on serum and plasma testing. The final assessment is that the samples were nonreactive. The true negative but initially reactive whole blood results are listed in parentheses <> in Table 9.

*** When using the cut-off 0.13 IU/mL, the sensitivity at 15 and at 30 minutes for fingerstick and venous whole blood samples was 97.9%.

Analytical Performance

Determine™ Syphilis TP (Test Strip 1)

1. Analytical sensitivity

WHO International Standard (NIBSC code: 05/122) for Syphilis was tested with Determine Syphilis TP. Determine Syphilis TP detected 18.75 mIU/mL.

2. Seroconversion Panel Testing

Commercially available seroconversion panels were evaluated.

Determine Syphilis TP detected earlier than the CE marked CLIA test and at the same time as the RPR test.

Table 10: Seroconversion Panel

Seroconversion Panel	Days from 1 st bleeding date	Non-treponema RPR test	CE marked CLIA test	Determine Syphilis TP
PSS901-01	0	-	-	-
PSS901-02	5	-	-	-
PSS901-03	10	-	-	-
PSS901-04	13	-	-	-
PSS901-05	31	+	-	+
PSS901-06	45	+	+	+
PSS901-07	48	+	+	+
PSS901-08	52	+	+	+
PSS901-09	59	+	+	+

3. Hook Effect

Twenty high titer positive specimens (TPPA titer ≥ 10240) were tested with the Determine Syphilis TP and all were positive. However, 12/20 showed a lower titer when the specimen was tested neat than when it was tested diluted, suggesting that the Determine Syphilis TP may be subject to high dose hook effect. However, in no case was the lower band score negative or weakly positive indicating that a false negative result in a high titer specimen is unlikely.

4. Cross reactivity and endogenous Interference

Determine Syphilis TP was evaluated for potential cross-reactivity and endogenous interference with 381 samples of patients which include infections other than syphilis, immunization and endogenous interfering substances.

Table 11: Cross Reactivity

Category	Sample group	Number of specimens tested	Determine™ Syphilis TP	
			Positive	Negative
Viral infections	HIV	10	0	10
	Hepatitis A virus	5	0	5
	Hepatitis B virus (HBsAg)	16	0	16
	Hepatitis C virus	16	0	16
	Cytomegalovirus (CMV) IgG	17	0	17
	Cytomegalovirus (CMV) IgM	14	0	14
	Epstein-Barr virus (EBV) IgG	14	0	14
	Epstein-Barr virus (EBV) IgM	13	1**	12
	Herpes simplex virus (HSV)1/2	12	0	12
Immunization	Flu Vaccine Specimens	10	0	10
Interfering substances	Antiretroviral treated specimens	3	0	3
	Human anti mouse antibody (HAMA)	10	0	10
	Anti-nuclear antibody (ANA)	10	0	10
	Rheumatoid Factor (RF)	10	0	10
	Systemic lupus erythematosus (SLE)	5	0	5
	Triglyceride	5	0	5
	High IgM	5	0	5
	Cholesterol	5	0	5
	High total protein	4*	0	4*
	Multiparous pregnancy (MP)	197	2***	195
Total of cross reactivity and interfering substances		381	3	378

*: 1 out of 5 high total protein samples was assay invalid due to no red color of control bar.

**: 1 out of 13 EBV IgM samples was false positive.

***: 2 out of 197 Multiparous pregnancy (MP) samples showed false positive.

Total specificity in the potential cross-reactivity and endogenous interference samples is 99.21% (95% CI 97.72-99.84%) representing no significant difference from clinical study.

5. Interference of exogenous substances

The interference of Bilirubin, Hemoglobin, Triglycerides, and EDTA was tested by spiking the interference materials into serum samples. 20 mg/dL Bilirubin, 500 mg/dL Hemoglobin, 1 g/dL Triglyceride, and 200 mg/dL EDTA showed no interference.

The interference of relevant medicines were evaluated.

There is no interference for spiked positive and negative for the medicines in Table 12.

Table 12: Interference

Type	Medicine	Concentration
Anti-parasitic medications	Eflornithine	16.8 mg/ml
	Metronidazole	0.123 mg/ml
Anti-bacterial medications	Amoxicillin	54.0 µg/ml
	Ciprofloxacin	12.0 µg/ml
	Azithromycin	11.1 µg/ml
	Levofloxacin	36.0 µg/ml

Type	Medicine	Concentration
Anti-malarial medications	Mefloquine	9.84 µg/ml
	Quinine	53.7 µg/ml
	Pyrimethamine	2.58 µg/ml
	Primaquine	7.23 µg/ml
Anti-retroviral medication	Ritonavir	44.4 µg/ml
Anti-tuberculosis medication	Isoniazid (INH)	60.0 µg/ml
Antiviral	Ribavirin	11.0 µg/ml
	Entecavir	24.6 ng/ml
	Valganciclovir hydrochloride	21.4 µg/ml
Common over-the-counter analgesic medications	Aspirin	0.782 mg/ml
	Paracetamol	0.156 mg/ml

6. Reproducibility

Reproducibility of Determine Syphilis TP has been determined using in-house reference panels. There were no differences observed:

- within-run (replicates tested by one operator)
- day to day (five days)
- between runs (three different test operators)
- between sites (three sites)
- lot-to-lot (3 different lots).

Determine™ HIV Early Detect (Test Strip 2)

1. Analytical Sensitivity of HIV-1 p24 antigen

The analytical sensitivity of the Determine™ HIV Early Detect was evaluated by testing WHO International Standard HIV-1 p24 Antigen (NIBSC code 90/636). The results demonstrated that the test could detect a concentration of 2 IU/mL HIV-1 p24 antigen.

2. HIV Seroconversion Panels

The sensitivity of the Determine™ HIV Early Detect was evaluated using 31 sets of seroconversion panels; each including early seroconversion panel members (98 early samples in total). The results were compared with the results of the commercially available CE marked HIV 4th generation rapid immunochromatographic test kit, CE marked HIV 3rd generation rapid immunochromatographic test kit, HIV antibody EIA kit, and HIV antigen EIA kit. Additionally, the results of 20 seroconversion panels were compared with the data of the CE marked HIV 4th generation kit (Chemiluminescent microparticle immunoassay (CMIA)). On 28 of 31 panels, Determine™ HIV Early Detect detected HIV infection earlier than a 3rd generation rapid immunochromatographic test kit. On 11 of 31 panels, Determine™ HIV Early Detect detected HIV infection earlier than a rapid 4th generation immunochromatographic test kit. Detection of 6 of 20 panels by Determine™ HIV Early Detect was delayed by 1 bleed date when compared to a 4th generation kit (CMIA). Detection of 8 of 31 panels by Determine™ HIV Early Detect was delayed by 1 bleed date when compared to the HIV Antigen Kit (EIA).

Table 13: HIV Seroconversion Panels

No.	Panel	Days to First Reactive Results				
		Determine™ HIV Early Detect	4 th gen rapid test	3 rd gen rapid test	4 th gen kit (CMIA)*	HIV antigen EIA*
		Ag/Ab	Ag/Ab	Ab	Ag/Ab	Ag
1	PRB943	12	12	19	7	7
2	PRB956	47	50	>50	47	47
3	PRB958	9	15	15	7	7
4	PRB961	27	29	>29	27	27
5	PRB962	14	14	>17	14	14
6	PRB963	17	17	>21	17	17
7	PRB966	48	51	51	44	44
8	PRB967	17	17	17	17	17
9	PRB968	26	26	28	26	26
10	PRB969	70	70	70	63	63
11	PRB970	7	7	14	0	0
12	PRB972	18	18	18	18	18
13	PRB973	11	11	>11	7	7
14	PRB974	9	16	>16	9	9
15	PRB975	14	>14	>14	14	14
16	PRB976	7	7	>9	7	7
17	PRB977	13	13	15	13	13
18	PRB924	26	26	33	NT	26
19	PRB926	7	7	27	NT	7
20	PRB934	0	7	7	NT	0
21	PRB935	28	28	43	NT	28
22	PRB940	7	7	15	NT	7
23	PRB945	13	13	15	NT	13
24	PRB947	9	9	11	NT	9
25	PRB948	23	23	>23	NT	23
26	PRB949 **	18	18	>18	NT	18
27	PRB951	11	11	19	NT	8
28	PRB955	7	12	12	NT	3
29	9018	28	35	35	28	28
30	9022	25	32	32	25	25

No.	Panel	Days to First Reactive Results				
		Determine™ HIV Early Detect	4 th gen rapid test	3 rd gen rapid test	4 th gen kit (CMIA)*	HIV antigen EIA*
		Ag/Ab	Ag/Ab	Ab	Ag/Ab	Ag
31	9079	40	42	55	40	40
Total		31	31	31	20	31
Detected earlier than Determine™ HIV Early Detect			0	0	6	8
Detected same day as Determine™ HIV Early Detect			20	3	14	23
Detected later than Determine™ HIV Early Detect			11	28	0	0

* Vendor certificate data, **PRB949 No.5 was not tested. NT: not tested

3. HIV-1 Antibody subtype

171 HIV-1 Antibody subtype specimens and diluted samples were tested: (A(15), CRF03_AB(1), CRF01_AE(15), CRF02_AG(16), CRF09_A/U(1), B(16), C(63), D(6), CRF05_D/F(1), F(7), G(9), H(4), H/A1(1), H/U(1), J(5), K(4), K/CRF09(1), CRF11_cpx(1), O(4)). The Determine™ HIV Early Detect detected all the specimens and diluted samples.

4. HIV-1 p24 Subtype Antigens

50 HIV virus panels were prepared by diluting 50 different cultured HIV isolates including HIV-1 group M subtypes: A(4), B(8), C(7), D(5), F(5), G(2), H(1), and CRF01-AE(10), CRF02-AG(2), HIV-1 group O(4) and HIV-2(2) with HIV negative human serum. The Determine™ HIV Early Detect detected all HIV-1 panels and did not detect HIV-2 panels.

5. Interference

- 20 mg/dL bilirubin, 500 mg/dL hemoglobin, 3300 mg/dL triglyceride, and 12 g/dL protein did not interfere with test results.
- At 15 g/dL protein, the test result presented false positive on Ab bar with Ab negative specimens. No interference was observed on Ag bar.
- When specimens diagnosed as positive to rheumatoid factor were spiked with p24 antigen, the test result of 2 out of 5 specimens presented false negative on Ag bar.

Determine™ HBsAg 2 (Test Strip 3)

1. Analytical sensitivity of HBsAg

The analytical sensitivity of Determine™ HBsAg 2 was evaluated by testing the WHO International Standard for HBsAg (NIBSC code 12/226). The results demonstrated that the test could detect a concentration of 0.1 IU/mL HBsAg.

2. HBV Seroconversion Panels

The sensitivity of Determine™ HBsAg 2 was evaluated using 32 sets of seroconversion panels; each including early seroconversion panel members. The results were compared with the results of a commercially available rapid immunochromatographic HBsAg test kit and a quantitative HBsAg test kit (Chemiluminescent microparticle immunoassay (CMIA)). On all sets, Determine™ HBsAg 2 detected HBsAg earlier than a rapid immunochromatographic test kit. On 4 of 32 sets, Determine™ HBsAg 2 detected HBsAg earlier than the CMIA kit. Detection of 4 of 32 sets by Determine™ HBsAg 2 was delayed by 1 bleed date when compared to the CMIA kit.

The first reactive date followed by continuously reactive results was compared with the first reactive date with the Determine™ HBsAg 2.

3. HBsAg Mutant Detection

HBsAg mutant susceptibility was evaluated with Determine™ HBsAg 2. A panel consisting of 14 different recombinant HBsAg mutant panels was tested. The mutant panel consisted of the following mutations, P120Q, T123A, T126N, T126S, Q129R, Q129H, Q129L, M133H, M133L, K141E, P142S, T143K, D144A and G145R. All 14 specimens were detected with Determine™ HBsAg 2 at 0.25 IU/mL, and 11 out of 14 specimens at 0.125 IU/mL.

4. Interference

- 20 mg/dL bilirubin, 500 mg/dL hemoglobin, 3300 mg/dL triglyceride, and 12 g/dL protein did not interfere with test results.
- At 15 g/dL protein, no interference was observed with HBsAg negative sample, and HBsAg positive sample at 0.2 IU/mL or more. With HBsAg positive sample at 0.1 IU/mL, the test result presented false negative at 15 minutes reading time, but showed positive results at 30 minutes reading time.

NOTE ON LOT NUMBER

The manufacturing process produces different lot numbers for the kit and test cards; these lot numbers are traceable.

REFERENCE DOCUMENTS

1. Meyer JC. Laboratory Diagnosis of Syphilis. *Curr Probl Dermatol*. 1996;24:1-11.
2. Pilcher C, Eron JJ, Galvin S, et al. (2004) Acute HIV revisited: new opportunities for treatment and prevention. *The Journal of Clinical Investigations* 113(7): 937-945.
3. Respass RA, Rayfield MA and Dondero TJ (2001) Laboratory testing and rapid HIV assays: applications for HIV surveillance in hard-to-reach populations. *AIDS* 15 Supplement 3: S49-S59.
4. Louie B, Wong E, Klausner JD, et al. (2008) Assessment of rapid tests for detection of Human Immunodeficiency Virus-specific antibodies in recently infected individuals. *Journal of Clinical Microbiology* 46(4): 1494-1497.
5. Lyons MS, Lindsell CJ, Hawkins DA, et al. (2008) Contributions to early HIV diagnosis among patients linked to care vary by testing venue. *BMC Public Health* 8:220.
6. Gurtler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP- 5180) from Cameroon. *Journal of Virology*. 1994; 68(3):1581-1585.
7. Yun-Fan Liaw, Chia-Ming Chu (2009) Hepatitis B virus infection. *Safety lancet* 373:582-592.
8. Jules L. Dienstag (2008) Hepatitis B Virus Infection. *The New England Journal of Medicine* 359:1486-1500.
9. World Health Organization. *Laboratory Biosafety Manual* Fourth edition. Geneva: World Health Organization, 2020.
10. CLINICAL AND LABORATORY STANDARDS INSTITUTE. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline-Fourth Edition* M29-A3 Vol. 25 No. 10.
11. Clinical and Laboratory Standards Institute. *Clinical Laboratory Waste Management; Approved Guideline – Third Edition*. GP05-A3 Vol.31 No.3 January 2011.
12. EPA Guide for Infectious Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
13. Clinical and Laboratory Standards Institute. *Collection of Capillary Blood Specimens; Approved Standard - 7th GP42-Ed7* September 2020.
14. CLINICAL AND LABORATORY STANDARDS INSTITUTE. *Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard - Fifth Edition* H4-A5 Vol. 24 No. 21.
15. Delaney KP, Branson BM, Uniyal A, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
16. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
17. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.
18. World Health Organization. *Guidelines for the prevention, diagnosis, care and treatment for people with chronic hepatitis B infection*, March 2024.

ADVICE LINE

For further information, please contact your distributor, or call one of the following Abbott Technical Support Care Centers:

Region	Phone	E-Mail Address
Europe	+ (44) 161 483 9032	EME.TechSupport@abbott.com
Middle East	+ (965) 2202 2828	EME.TechSupport@abbott.com
Asia Pacific	+ (61) 7 3363 7711	AP.TechSupport@abbott.com
Africa	+ (27) 10 500 9700	arcis.techsupport@abbott.com
Russia & Commonwealth of Independent States (CIS)	+ (7) 499 403 9512	arcis.techsupport@abbott.com
Latin America	+ (57) 601 482 4033	LA.TechSupport@abbott.com

GLOSSARY OF SYMBOLS

	Catalogue number		<i>In Vitro</i> Diagnostic Medical Device
	Consult instructions for use		Temperature limit
	Keep dry		Contains sufficient for n tests
	Use by date		Do not use if package is damaged and consult instructions for use
	Batch code		Manufacturer
	Keep away from sunlight		Do not re-use
	Unique Device Identification		Caution
	Patient Information Website		Fragile, handle with care
	Biological risks		Patient Identification
	Medical device		Authorized representative in the European Community/ European Union
	Sterilized using irradiation		Sterilized using ethylene oxide

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Determine™ HIV Early Detect (bandelette de test 2)

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Assembled by : ABON Biopharm (Hangzhou) Co., Ltd.