WHO Prequalification of Diagnostics Programme PUBLIC REPORT

Product: Determine HIV-1/2¹
Number: PQDx 0033-013-00

Determine HIV-1/2 with product codes **7D2342**, **7D2343**, **7D2343SET**, and **7D2343SETS** manufactured by **Abbott Diagnostics Medical Co., Ltd²**, **Rest-of-World regulatory version** (non-CE-marked regulatory version) was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 25 November 2011.

Summary of prequalification status for Determine HIV-1/2

	Date	Outcome
Prequalification listing	25 November 2011	listed
Dossier review	22 November 2011	MR
Site inspection (s) of the quality management system	24 May 2019	MR
Product performance evaluation	11 November 2011	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Change of shelf life, extended from 14 months to 18 months	5 July 2016
3.0	Change of physical address of the Chiba Logistics Center to 483-2 Matsuhidai, Matsudo-shi, Chiba, 270-2214 Japan	18 October 2016
4.0	Change of HIV-1 and HIV-2 antibodies, and their reagent formulation, in the Control bar solution to address customer complaints regarding the weakness of control line intensity.	15 February 2017

¹ Product name was changed from Alere Determine HIV-1/2 to Determine HIV 1/2.

² Manufacturer name changed from Alere Medical Co.Ltd to Abbott Diagnostics Medical Co.Ltd.

5.0	Introduction of safety blood lancets as a component of whole blood assay SET. Addition of new configurations of its IVDs, where lancets currently provided with the IVDs are replaced with contact-activated safety lancets. Add a product code corresponding to the new safety lancet (7D2343SETS (Alere Determine HIV-1/2)).	1 March 2017
6.0	Shelf life extension of the Chase Buffer (7D2243) from 18 to 24 months. The Chase Buffer is included in the following SET products: PQDx 0033-013-00: Alere Determine HIV-1/2 SET (7D2343SET) (100 tests/kit set for whole blood assay). Alere Determine HIV-1/2 SETS (7D2343SETS) (100 tests/kit set for whole blood assay with safety lancet).	30 July 2018
7.0	Change the manufacturer's name from Alere Medical Co., Ltd. to Abbott Diagnostics Medical Co., Ltd. Along with this, product names were changed. Product name changed from Alere Determine HIV-1/2 to Determine HIV 1/2. Products themselves are completely the same.	11 December 2019
8.0	IFU updates and stating the status of post-prequalification commitments.	13 January 2022
9.0	Downsizing of Determine HIV-1/2 SET for the benefits of shipments and storage. This SET Downsizing includes the changes of Capillary Tubes from glass to plastic (plastic is safer than glass).	14 April 2022
10.0	Updates to the labelling of the Chase Buffer (7D2243) included in each SET product.	13 July 2023
11.0	Removed the CE mark from the included lancets from the external supplier; Health Canada registration of the lancets is retained.	18 September 2025

Intended use:³

According to the claim of Abbott Diagnostics Medical Co., Ltd, "*the DetermineHIV-1/2 is an In Vitro, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals. The test is for professional use only*".

³ This product is one that uses Protein A to detect human IgG antibodies. Protein A is also able to detect other classes of human antibody (IgA, IgD, IgE and IgM) but not as reliably as it does IgG. This product has been prequalified with respect to its ability to detect human IgG antibodies. Any claim to detect other types of antibodies on this kind of product has not been validated based on WHO prequalification requirements.

Assay principle:

According to the claim of Abbott Diagnostics Medical Co., Ltd, "*Determine HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2. Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site.*

If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site.

If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site.

To ensure assay validity, a procedural control bar is incorporated in the assay device."

Test kit contents:

Component	20 tests (7D2342)	100 tests (7D2343)	100 tests (7D2343SET)	100 tests (7D2343SETS)
Test cards	2 cards of 10 tests/cards	10 cards of 10 tests/card	20 cards of 5 tests/card	10 cards of 10 tests/card
Instructions for use (IFU)	1	1	1	1
Chase buffer (7D2243)	Not provided	Not provided	1 bottle x 2.5ml and IFU	1 bottle x 2.5ml and IFU
Capillary tubes, plastic	Not provided	Not provided	100	Not provided
EDTA capillary tubes, glass (7D2222)	Not provided	Not provided	Not provided	100
Blood lancets, sterile	Not provided	Not provided	100	Not provided
Blood lancets, sterile (safety)	Not provided	Not provided	Not provided	100

Items required but not provided within the test kit:

- Disposable gloves, timing device
- Micropipette capable of delivering 50 µL (other than fingerstick)
- Alcohol swab, gauze pad, lancet (for fingerstick).

For testing Whole Blood samples

- Chase Buffer (7D2243) 1 Bottle (2.5 mL) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

For testing Whole Blood samples (fingerstick assay)

- EDTA Capillary Tubes (7D2222). Capillary tubes are required for 7D2342 and 7D2343.
- For 7D2343SET, Plastic capillary tubes are provided within the SET.

Storage:

The test kit should be stored at 2-30 °C.

Shelf life upon manufacture:

18 months Determine HIV-1/2 test card

Prioritization for prequalification

Based on the established eligibility criteria, Determine HIV-1/2 was given priority for the WHO prequalification assessment.

Product dossier assessment

Abbott Diagnostics Medical Co., Ltd (formerly called Alere Medical Co. Ltd) submitted a product dossier for Determine HIV-1/2 as per the "Instructions for compilation of a product dossier" (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO.

Based on the product dossier screening and assessment findings, the product dossier for Determine HIV-1/2 meets WHO prequalification requirements.

Manufacturing site inspection

A comprehensive inspection was performed at the site of manufacture of Determine HIV-1/2 manufactured by Abbott Diagnostics Medical Co., Ltd. (formerly called Alere Medical Co. Ltd) at 357 Matsuhidai Matsudo-shi, Chiba-ken 270-2214, Japan, between 26-27 November 2018⁴. The inspection procedure is described in "Information for manufacturers on WHO

⁴ Initial inspection took place on 28 September to 1 October 2010

prequalification inspection procedures for the sites of manufacture of diagnostics" (PQDx_014 v4).

The inspection found that Abbott Diagnostics Medical Co. Ltd had an established quality management system and manufacturing practices in place that should ensure the manufacture of a product of consistent quality.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 24 May 2019.

Based on the site inspection and corrective action plan review, the quality management system for Determine HIV-1/2 meets WHO prequalification requirements.

Product performance evaluation

Determine HIV-1/2 was evaluated by WHO in the third quarter of 2011 at the Institute of Tropical Medicine, Belgium – a WHO Collaborating Centre for HIV/AIDS Diagnostics and Laboratory Support. The laboratory evaluation was conducted according to the "WHO Protocol for the laboratory evaluation of HIV serology assays" (PQDx_030 V1.0) and drew the following conclusions:

Determine HIV-1/2 is an immunochromatographic rapid diagnostic test for the detection of antibodies to HIV-1/2 in human serum, plasma, and whole blood. A volume of 50µl of serum, plasma or venous/capillary whole blood is required to perform the test procedure. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities. Reading of the results is performed visually, i.e. subjective reading.

In this limited performance evaluation using a panel of 1079 biological specimens, we observed an initial sensitivity (95% CI) of 100% (99.1% - 100%) and an initial specificity (95% CI) of 97.87% (96.4% - 98.8%) compared to the reference assays. The final sensitivity (95% CI) was 100% (99.1% - 100%), and the final specificity (95% CI) was 98.93% (97.8% - 99.6%) compared to the reference assays. In this study, 0.3% of the overall results were recorded as indeterminate. Results were interpreted independently by three technicians; the inter-reader variability was 1.4%. The invalid rate was 0.3% for initial testing and 0.1% for repeat testing.

Labelling

1. 7D2342, 7D2343

1.1 Labels

1.1.1 7D2342 Outer Label

DETERMINE™
HIV-1/2

EN

For *In Vitro* Diagnostic Use
Determine™ HIV-1/2 is a visual read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2.
For professional use only.
7D2243 Chase buffer is required for whole blood testing.
Kit contains:
2 HIV-1/2 recombinant antigen and synthetic peptide coated test cards.

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© 2019 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. July, 2019

241563/R8

Determine HIV-1/2
ROW 20 Test Pouch Label

PN: 241563/R8

Pouch Size:
544mm(w) x 160mm(h)

Label Artwork Size:
204mm(w) x 132mm(h)

PMS 303

PMS 2269

White

1.1.2 7D2343 Outer Label

DETERMINE™
HIV-1/2

EN

For *In Vitro* Diagnostic Use
Determine™ HIV-1/2 is a visual read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2.
For professional use only.
7D2243 Chase buffer is required for whole blood testing.

Kit contains:
10 HIV-1/2 recombinant antigen and synthetic peptide coated test cards.

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White

1.2 Instructions for Use⁵

⁵ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

REF

Catalogue Number

IVD

In Vitro Diagnostic Medical Device

2°C ~ 30°C

Store at 2-30°C

Key to symbols used

Σ 20

Contains Sufficient for 20 tests

Σ 100

Contains Sufficient for 100 tests

☀

Keep away from sunlight

⌛

Do not reuse

⊗

Do not use if package is damaged

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE
The Determine™ HIV-1/2 is an *In Vitro*, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals. The test is for professional use only.

SUMMARY AND EXPLANATION OF THE TEST
AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of the AIDS virus elicits the production of specific antibodies to either HIV-1 or HIV-2.^{1,2,3}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE
Determine™ HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2. Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site. If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site. If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site. To insure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS
Determine™ HIV-1/2 Serum/Plasma Assay, 20 Tests (7D2342) or 100 Tests (7D2343)
• Determine™ HIV-1/2 Test Card, 2 cards or 10 cards (10 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.

ACCESSORIES (required but not provided)
For testing Whole Blood samples
• 1 Bottle (2.5 mL) Chase Buffer (7D2243) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

Whole Blood (fingerstick assay)
• EDTA Capillary Tubes (7D2222)

Materials Required But Not Provided
Disposable gloves, timing device
Micropipette capable of delivering 50 µL (other than fingerstick)
Alcohol swab, gauze pad, Lancet (for fingerstick)

WARNINGS AND PRECAUTIONS
For In Vitro Diagnostic Use.

CAUTION:
Appropriate biosafety practices^{4,5} should be used when handling specimens and reagents. These precautions include, but are not limited to the following:
• **Wear gloves.**
• **Do not pipette by mouth.**
• **Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.**
• **Clean and disinfect all spills of specimens or reagents using a suitable disinfectant, such as 0.5% sodium hypochlorite.^{6,7}**
• **Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.^{8,9}**

STORAGE
The Determine™ HIV-1/2 Test Cards and Chase Buffer must be stored at 2-30°C until expiration date. Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date. Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.

SPECIMEN COLLECTION
Serum, Plasma, and Whole Blood Collection by Venipuncture
Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.
NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingerstick¹⁰
Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

- Choose the fingertip of the middle, ring, or index finger (whichever is the least callused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
- Clean fingertip with alcohol; allow to air dry. Position the hand palm-side up.
- Use a new lancet for each person. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
- Wipe away the first drop of blood with a sterile gauze pad.
- Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood*. Avoid air bubbles.

*If EDTA Capillary Tubes (7D2222) will be used, fill the tube with blood between the 2 marked lines (50 µL).

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole serum collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE
The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.
This test should be performed at 15 to 40°C.

NOTES:

- Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.**
- Assay should be initiated within 2 hours after removing the protective foil cover from each test.**

- Remove the protective foil cover from each test.
- For serum or plasma samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
- For whole blood (venipuncture) samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait one minute, then apply one drop of Chase Buffer to the sample pad.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
- For whole blood (fingerstick) samples:
 - Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
 - Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.

QUALITY CONTROL
To insure assay validity, a procedural control is incorporated in the device and is labeled “Control”. If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS
POSITIVE (Two Bars)
Red bars appear in both the control window (labeled “Control”) and the patient window (labeled “Patient”) of the strip. Any visible red bar in the patient window should be interpreted as positive.

NEGATIVE (One Bar)
One red bar appears in the control window of the strip (labeled “Control”), and no red bar appears in the patient window of the strip (labeled “Patient”).

INVALID (No Bar)
If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

NOTES:

- The test result is positive even if the patient bar appears lighter or darker than the control bar.
- The control bar may exhibit a weak intensity for some patient samples, particularly those with high titer HIV.
- Upon appearance of a red bar in the control window, no matter how faint, the test result is considered valid.
- If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

LIMITATIONS OF THE PROCEDURE

- The Determine™ HIV - 1/2 test is designed to detect antibodies to HIV - 1 and HIV - 2 in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- A negative result with Determine™ HIV-1/2 does not exclude the possibility of infection with HIV. A false negative result can occur in the following circumstances:
 - low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
 - infection with a variant of the virus that is less detectable by the Determine™ HIV-1/2 assay configuration
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration
 - specimen handling conditions which result in loss of HIV antibody multivalency
 - HIV-infected persons taking antiretroviral medication^{11,12,13}

For these reasons care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA may give incorrect results.
- Neonates of HIV-infected mothers may carry maternal antibodies to HIV for up to around eighteen months, which may not necessarily indicate the true infection status of the new born.

PERFORMANCE CHARACTERISTICS
SPECIFICITY
A total of 1594 serum and plasma specimens from Asia, West Africa, and North America were tested by Determine™ HIV-1/2 and a commercially available test (Table I).

Table I Specificity of Determine™ HIV-1/2			
Population	Number of Specimens Tested	Negative by Determine™ HIV-1/2	Negative by a Commercially Available Test***
Seronegative Serum	908	907/908 (99.89%)	908/908 (100.00%)
Plasma	403	403/403 (100.00%)	403/403 (100.00%)
Pregnant Females	58*	57/57 (100.00%)	57/57 (100.00%)
West Africans	49	48/49 (97.96%)	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	176*	173/175 (98.86%)	174/175 (99.45%)
Total	1594**	1588/1592 (99.75%)	1590/1592 (99.87%)


* One specimen from a pregnant female and an HCV positive patient were positive by both Determine™ HIV-1/2 and the commercially available test. Both specimens confirmed positive by HIV-1 Western Blot.
** 456 specimens were from North America, 1089 specimens were from Asia, and 49 specimens were from Africa.
***The reference method of a commercially available test is particle agglutination.

A total of 3663 seronegative serum and plasma specimens from North America, Asia, and Africa were tested by Determine™ HIV-1/2 and commercially available tests (Table II). The specimens from North America, Asia, and 49 of 2118 specimens from Africa (referred to as “West Africans” in Table I) were included in Table I. Discordant specimens were confirmed negative by either Western blot or HIV-1 PCR assays.


Table II A Comparison of Determine™ HIV-1/2 Specificity by Geographic Area			
Area	Number of Specimens Tested	Negative by Determine™ HIV-1/2	Negative by Commercially Available Tests*
North America	456	451/454 (99.34%)	453/454 (99.78%)
Asia	1089	1089/1089 (100.00%)	1089/1089 (100.00%)
Africa	2118	2079/2118 (98.16%)	2100/2118 (99.15%)

* The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.


Pipette

EDTA Capillary Tubes

Pipette

A total of 368 seronegative whole blood specimens from Thailand were tested with paired serum and plasma by Determine™ HIV-1/2. Thirty-nine of the whole blood specimens were collected by both venipuncture and fingerstick (Table III).

Table III A Comparison of Determine™ HIV-1/2 Specificity in Seronegative Whole Blood and Paired Serum and Plasma Specimens		
Specimen Type	Number of Specimens Tested	Negative by Determine™ HIV-1/2
Serum	368	368/368 (100.00%)
Plasma	368	368/368 (100.00%)
Whole Blood (venipuncture)	368	368/368 (100.00%)
Whole Blood (fingerstick)	39	39/39 (100.00%)

SENSITIVITY
A total of 869 HIV-1 and HIV-2 antibody positive serum and plasma specimens from Asia, Africa, North and South America were tested by Determine™ HIV-1/2 and a commercially available test (Table IV).

Table IV Sensitivity of Determine™ HIV-1/2			
Population	Number of Specimens Tested	Positive by Determine™ HIV-1/2	Positive by a Commercially Available Test**
HIV-1 Positive	521*	521/521 (100.00%)	521/521 (100.00%)
HIV-2 Positive	114*	114/114 (100.00%)	114/114 (100.00%)
HIV-1 Subtypes A-G	222	222/222 (100.00%)	Not Tested Not Tested
HIV-1 Group O	12	12/12 (100.00%)	Not Tested Not Tested
Total	869	869/869 (100.00%)	635/635 (100.00%)

* 228 specimens were from North America, 296 specimens were from Asia, and 111 specimens were from Africa.
**The reference method of a commercially available test is particle agglutination.

A total of 1653 seropositive serum and plasma specimens from North America, Asia, and Africa were tested by Determine™ HIV-1/2 and commercially available tests (Table V). The specimens from North America, Asia, and 111 of 1129 specimens from Africa (referred to as “HIV-2 Positive” in Table IV) were included in Table IV. Discordant specimens were confirmed HIV-1 positive by either Western blot or HIV-1 PCR assays.

Table V A Comparison of Determine™ HIV-1/2 Sensitivity by Geographic Area			
Area	Number of Specimens Tested	Positive by Determine™ HIV-1/2	Positive by Commercially Available Tests**
North America	228	228/228 (100.00%)	228/228 (100.00%)
Asia	296	296/296 (100.00%)	296/296 (100.00%)
Africa	1129	1128*/1129 (99.91%)	1129/1129 (100.00%)

* One negative specimen by Determine™ HIV-1/2 confirmed positive by HIV-1 PCR.
**The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.

A total of 102 seropositive whole blood specimens from Thailand were tested with paired serum and plasma by Determine™ HIV-1/2. Thirty-two of the whole blood specimens were collected by both venipuncture and fingerstick (Table VI).


Table VI A Comparison of Determine™ HIV-1/2 Sensitivity in Seropositive Whole Blood and Paired Serum and Plasma Specimens		
Specimen Type	Number of Specimens Tested	Positive by Determine™ HIV-1/2
Serum	102	102/102 (100.00%)
Plasma	102	102/102 (100.00%)
Whole Blood (venipuncture)	102	102/102 (100.00%)
Whole Blood (fingerstick)	32	32/32 (100.00%)

The manufacturing process produces different lot numbers for the kit and test cards; these lot numbers are traceable.

BIBLIOGRAPHY


- Piot P, Plummer FA, Mhalu FS, Lambarory JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988; 239: 573-579.
- Weniger BG, Takebe Y, Ou CY, *et al*. The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994; 8(S2): S13-S28.
- Gürtler LG, Hauser PH, Eberle J, *et al*. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3): 1581-1585.
- World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization, 1993.
- National Committee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991: 1-43.
- CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, MMWR 1987; 36(2S): 3S-18S.
- Sehustler LM, Hollinger FB, Dreesman GR, *et al*. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Applied and Environmental Microbiology*. 1981; 42(5): 762-7.
- Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
- EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
- Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
- Delaney KP, Branson BM, *et al*. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
- O’Connell RJ, Merritt TM, Malia JA, *et al*. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
- O’Connell RJ, Agan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

Positive


+

Negative

–

Invalid

Advice Line		
For further information, please contact your distributor, or call to one of the following Abbott Technical Support Care Centers:		
Region	Phone	E-Mail Address
Europe	+ (44) 161 483 9032	EME.TechSupport@abbott.com
Middle East	+ (965) 2202 2828	EME.TechSupport@abbott.com
Asia Pacific	+ (61) 7 3363 7100	AP.TechSupport@abbott.com
Africa	+ (27) 10 500 9700	arcis.techsupport@abbott.com
Russia & CIS	+ (7) 499 403 9512	arcis.techsupport@abbott.com
Latin America	+ (57) 1 482 4033	LA.TechSupport@Abbott.com

357 Matsuhidai, Matsudo-shi,
Chiba, 270-2214 Japan
Tel +81 47 311 5750

© 2021 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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- Piot P, Plummer FA, Mhalu FS, Lambarory JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988; 239: 573-579.
- Weniger BG, Takebe Y, Ou CY, *et al*. The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994; 8(52): S13-S28.
- Gürtler LG, Hauser PH, Eberle J, *et al*. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3): 1581-1585.
- World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization, 1993.
- National Commitee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991: 1-43.
- CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, MMWR 1987; 36(25): 3S-18S.
- Sehustler LM, Hollinger FB, Dreesman GR, *et al*. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Applied and Environmental Microbiology*. 1981; 42(5): 762-7.
- Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
- EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
- Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
- Delaney KP, Branson BM, *et al*. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
- O'Connell RJ, Merritt TM, Malia JA, *et al*. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
- O'Connell RJ, Ajan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

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BIBLIOGRAPHIE

- Piot P, Plummer FA, Mhalu FS, Lambarary JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988; 239: 573-579.
- Weniger BG, Takebe Y, Ou CY, *et al*. The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994; 8(S2): S13-S28.
- Gürtler LG, Hauser PH, Eberle J, *et al*. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3): 1581-1585.
- World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization, 1993.
- National Committee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991: 1-43.
- CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, MMWR 1987; 36(2S): 3S-18S.
- Sehustler LM, Hollinger FB, Dreesman GR, *et al*. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Applied and Environmental Microbiology*. 1981; 42(5): 762-7.
- Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
- EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
- Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
- Delaney KP, Branson BM, *et al*. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
- O'Connell RJ, Merritt TM, Malia JA, *et al*. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
- O'Connell RJ, Agan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

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BIBLIOGRAFIA

- Piot P, Plummer FA, Mhalu FS, Lamborary JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988; 239: 573-579.
- Weniger BG, Takebe Y, Ou CY, *et al.* The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994; 8(S2): S13-S28.
- Gürtler LG, Hauser PH, Eberle J, *et al.* A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3): 1581-1585.
- World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization, 1993.
- National Committee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991: 1-43.
- CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, MMWR 1987; 36(2S): 3S-18S.
- Sehstler LM, Hollinger FB, Dreesman GR, *et al.* Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Applied and Environmental Microbiology*. 1981; 42(5): 762-7.
- Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
- EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
- Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
- Delaney KP, Branson BM, *et al.* Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
- O’Connell RJ, Merritt TM, Malia JA, *et al.* Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
- O’Connell RJ, Agan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

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
EN
For In Vitro Diagnostic Use.
Determine® HIV-1/2 is a visual read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2.
For professional use only.

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



© 2021 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
241975/R1 May, 2021

2.1.2 100 Test Pouch Label (7D2343SET)

2.1.3 Chase Buffer Bottle Label (7D2343SET)

Chase Buffer

Black

Box

Size: (W)22mm x (L)44mm x (H)56mm

PN: 1135945401

1156253101

REF

7D2243

1156253101

April 2022

DETERMINE™

CHASE BUFFER

Abbott

EN

Intended Purpose

Chase Buffer is intended for use with Determine™ in vitro diagnostic rapid test devices. Chase Buffer enables the migration of the specimen and conjugate through the test device when testing whole blood specimens. It is specially formulated and optimised to ensure the Determine™ test performs correctly and controls assay flow. It is for professional use and near-patient testing, and not for self-testing.

For full procedure refer to the diagnostic assay IFUs:

- Determine™ HIV-1/2 (7D23)
- Determine™ HIV Early Detect (7D28)
- Determine™ HBsAg 2 (7D29)
- Determine™ Syphilis TP (7D24)
- Determine™ HBsAg (7D25)

When adding Chase Buffer to the sample pad, hold the bottle vertically. One drop to be used per test (one bottle of Chase Buffer can be used for 100 tests). Avoid drop with bubble, do not add buffer to specimen.

USER TRAINING, QUALIFICATION AND EXPERIENCE REQUIREMENTS

This Chase Buffer is designed for use by professionals in both laboratory and near-patient settings including laboratory technologists. Suitable personnel shall have a consummate level of training and experience in specimen preparation, the running and interpretation of specimens, management of all supplies and waste disposal as per local regulations and generation of reports.

Caution: Do not eat, drink, smoke or apply make-up while using this product. For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Contents	1 Bottle (2.5 mL) Chase Buffer contains: CHASE BUFFER phosphate buffer solution (20mM pH7.4) containing preservatives Nipasept™ and A56620 (saraflaxacin).
Storage Instruction	Recap the Chase Buffer to avoid evaporation or spillage and store at 2-30 °C until expiration date.

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Abbott Rapid Dx International Limited
Parkmore East Business Park, Ballybri,
Galway, Ireland, H91 VK7E

UK Responsible person

Unipath Limited (trading as Alere International)
Clearblue Innovation Centre,
Priory Business Park, Stannard Way, Bedford,
Bedfordshire, England, MK443UP, United Kingdom

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


19mm

Model: I

- 5. Discard used lancets in a suitable container and/or in accordance with local requirements.
- 4. After procedure, help minimize accidental needle-stick by pressing needle-tip into flat side of the saved cap.
- 3. Perform lancing procedure using medically accepted aseptic techniques.
- 2. Pull off the cap to expose needle tip. Save cap for later use (step 04).
- 1. Twist lancet-cap to separate. **NOTE:** Do not use if cap appears damaged or loose before twisting.

NOTE: If using a lancing device, follow instructions for use of that device.

Do not use lancet if protective cap has been previously removed.

Instructions for Use

Litetouch™
Sterile Lancet
medicore™
Austell, GA, USA
1-855-867-2190

100t qty : 100pcs
Item No. : M2721

LOT : 220802T
LOT : 2027-08-01
LOT : 2022-08-02

(01) 06957138607882
(10) 220802T
(17) 270801
(11) 220802
(240) IF28100100

Tianjin Huahong Technology Co., Ltd.
Plant B, No.276, Huagang Road, Tianjin
Pilot Free Trade Zone, China
Tianjin 300308, China
Tel: +86 022-58775333 +86 4008-587-753
Fax: +86 022-58775353
E-mail: huahong@163.com

STERILE R
Irradiation Sterilization
Single Use Only
Rev# 2024-04-22
Made in China

28 gauge

2.2 Instruction for Use⁶

⁶ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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






NAME AND INTENDED USE

The Determine™ HIV-1/2 is an *In Vitro*, visually read, qualitative immuno-assay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals. The test is for professional use only.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of the AIDS virus elicits the production of specific antibodies to either HIV-1 or HIV-2.^{1,2,3}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Determine HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2. Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site. If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site. If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site. To insure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS

- Determine HIV-1/2 SET (7D2343SET) (100 Test for testing whole blood samples)
- Determine HIV-1/2 Test Card 20 cards (5 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.
 - 1 Bottle (2.5 mL) Chase Buffer (7D2243) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.
 - Capillary Tubes
 - Blood Lancet

Materials Required But Not Provided

Disposable gloves, timing device
Micropipette capable of delivering 50 µL (other than fingerstick)
Alcohol swab, gauze pad

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

CAUTION:

The package insert is placed in the kit box. DO NOT put the package insert inside the foil pouch.

Appropriate biosafety practices^{4,5} should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using a suitable disinfectant, such as 0.5% sodium hypochlorite.^{6,7}
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.^{8,9}

STORAGE

The Determine HIV-1/2 Test Cards and Chase Buffer must be stored at 2-30°C until expiration date. Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date. Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture

Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingerstick¹⁰

- Before collecting a finger-stick specimen, place a capillary tube and a lancet on a clean dry surface. After taking a capillary tube, seal to close the plastic bag and store the remaining capillary tubes in the kit box to avoid the sunlight. Do not use a capillary tube that is dirty or bent. Do not touch the open end of the capillary tubes before use.
- Ask the patient to increase blood circulation by rubbing the hands together.
- Choose the fingertip of the middle, ring, or index finger (whichever is the least callused). Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
- Clean the fingertip with alcohol and allow to air-dry. Position the hand, palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
- Wipe away the first drop of blood with a sterile gauze pad.
- Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times.
- Confirm the location of the line on the capillary tube. Use your index finger and thumb to gently squeeze the body of capillary tube, let the open end of capillary tube touch the blood bead, gently release the finger - blood will be sucked into the capillary tube.
- Carefully collect the blood up to the line and do not collect the blood above the line. Reliability of assay results cannot be guaranteed if blood is collected below or above the line.
- If the sample is not enough to fill up to the line of the capillary tube apply, gentle pressure to the finger again to make a drop of blood and take the blood until the line is reached.
- Confirm that no air bubble is present in the sample. If an air bubble is present, or if blood is collected below or above the line after completion of collection, discard the capillary tube with the sample and do not use it.
- Immediately proceed to Testing Procedure 4 for whole blood (finger-stick) samples.

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 5-test card can be removed by bending and tearing at the perforation. This test should be performed at 15 to 40°C.

NOTES:

- The lot number and expiration date are printed on each test card and on the foil pouch.
 - Assay should be initiated within 2 hours after removing the protective foil cover from each test.
 - Recap and store the Chase Buffer at 2-30°C to avoid evaporation or spillage.
- Remove the protective foil cover from each test.
 - For serum or plasma samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
 - For whole blood (venipuncture) samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait one minute, then apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
 - For whole blood (fingerstick) samples:
 - Squeeze the body of capillary tube expelling the blood onto the sample pad of the test device. Apply all blood in the capillary tube to the sample pad. Do not splash the blood.
 - Dispose of the used capillary tubes as biohazardous material according to local regulations.
 - Wait exactly one minute to allow the sample to be absorbed, then apply one drop of Chase Buffer to the sample pad, holding the bottle vertically
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.


QUALITY CONTROL

To insure assay validity, a procedural control is incorporated in the device and is labeled “Control”. If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS

POSITIVE (Two Bars)

Red bars appear in both the control window (labeled “Control”) and the patient window (labeled “Patient”) of the strip. Any visible red bar in the patient window should be interpreted as positive.


NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled “Control”), and no red bar appears in the patient window of the strip (labeled “Patient”).

INVALID (No Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

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LIMITATIONS OF THE PROCEDURE

- The Determine HIV - 1/2 test is designed to detect antibodies to HIV - 1 and HIV - 2 in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- A negative result with Determine HIV-1/2 does not exclude the possibility of infection with HIV. A false negative result can occur in the following circumstances:
 - low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
 - infection with a variant of the virus that is less detectable by the Determine HIV-1/2 assay configuration
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration
 - specimen handling conditions which result in loss of HIV antibody multivalency
 - HIV-infected persons taking antiretroviral medication^{11,12,13}

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PERFORMANCE CHARACTERISTICS
SPECIFICITY

A total of 1594 serum and plasma specimens from Asia, West Africa, and North America were tested by Determine HIV-1/2 and a commercially available test (Table I).

Table I Specificity of Determine HIV-1/2			
Population	Number of Specimens Tested	Negative by Determine HIV-1/2	Negative by a Commercially Available Test***
Seronegative Serum	908	907/908 (99.89%)	908/908 (100%)
Plasma	403	403/403 (100%)	403/403 (100%)
Pregnant Females	58*	57/57 (100%)	57/57 (100%)
West Africans	49	48/49 (97.96%)	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	176*	173/175 (98.86%)	174/175 (99.45%)
Total	1594**	1588/1592 (99.75%)	1590/1592 (99.87%)

*One specimen from a pregnant female and an HCV positive patient were positive by both Determine HIV-1/2 and the commercially available test. Both specimens confirmed positive by HIV-1 Western Blot.

**456 specimens were from North America, 1089 specimens were from Asia, and 49 specimens were from Africa.
***The reference method of a commercially available test is particle agglutination.

A total of 3663 seronegative serum and plasma specimens from North America, Asia, and Africa were tested by Determine HIV-1/2 and commercially available tests (Table II). The specimens from North America, Asia, and 49 of 2118 specimens from Africa (referred to as “West Africans” in Table I) were included in Table I. Discordant specimens were confirmed negative by either Western blot or HIV-1 PCR assays.

Table II A Comparison of Determine HIV-1/2 Specificity by Geographic Area			
Area	Number of Specimens Tested	Negative by Determine HIV-1/2	Negative by a Commercially Available Test*
North America	456	451/454 (99.34%)	453/454 (99.78%)
Asia	1089	1089/1089 (100%)	1089/1089 (100%)
Africa	2118	2079/2118 (98.16%)	2100/2118 (99.15%)

*The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.

A total of 368 seronegative whole blood specimens from Thailand were tested with paired serum and plasma by Determine HIV-1/2. Thirty-nine of the whole blood specimens were collected by both venipuncture and fingerstick (Table III).

Table III A Comparison of Determine HIV-1/2 Specificity in Seronegative Whole Blood and Paired Serum and Plasma Specimens		
Specimen Type	Number of Specimens Tested	Negative by Determine HIV-1/2
Serum	368	368/368 (100%)
Plasma	368	368/368 (100%)
Whole Blood (venipuncture)	368	368/368 (100%)
Whole Blood (fingerstick)	39	39/39 (100%)

SENSITIVITY

A total of 869 HIV-1 and HIV-2 antibody positive serum and plasma specimens from Asia, Africa, North and South America were tested by Determine HIV-1/2 and a commercially available test (Table IV).

Table IV Sensitivity of Determine HIV-1/2			
Population	Number of Specimens Tested	Positive by Determine HIV-1/2	Positive by a Commercially Available Test**
HIV-1 Positive	521*	521/521 (100%)	521/521 (100%)
HIV-2 Positive	114*	114/114 (100%)	114/114 (100%)
HIV-1 Subtypes A-G	222	222/222 (100%)	Not Tested
HIV-1 Group O	12	12/12 (100%)	Not Tested
Total	869	869/869 (100%)	635/635 (100%)

*228 specimens were from North America, 296 specimens were from Asia, and 111 specimens were from Africa.

**The reference method of a commercially available test is particle agglutination.

A total of 1653 seropositive serum and plasma specimens from North America, Asia, and Africa were tested by Determine HIV-1/2 and commercially available tests (Table V). The specimens from North America, Asia, and 111 of 1129 specimens from Africa (referred to as “HIV-2 Positive” in Table IV) were included in Table IV. Discordant specimens were confirmed HIV-1 positive by either Western blot or HIV-1 PCR assays.

Table V A Comparison of Determine HIV-1/2 Sensitivity by Geographic Area			
Area	Number of Specimens Tested	Positive by Determine HIV-1/2	Positive by a Commercially Available Test**
North America	228	228/228 (100%)	228/228 (100%)
Asia	296	296/296 (100%)	296/296 (100%)
Africa	1129	1128*/1129 (99.91%)	1129/1129 (100%)

*One negative specimen by Determine HIV-1/2 confirmed positive by HIV-1 PCR.
**The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.

A total of 102 seropositive whole blood specimens from Thailand were tested with paired serum and plasma by Determine HIV-1/2. Thirty-two of the whole blood specimens were collected by both venipuncture and fingerstick (Table VI).

Table VI A Comparison of Determine HIV-1/2 Sensitivity in Seropositive Whole Blood and Paired Serum and Plasma Specimens		
Specimen Type	Number of Specimens Tested	Positive by Determine HIV-1/2
Serum	102	102/102 (100%)
Plasma	102	102/102 (100%)
Whole Blood (venipuncture)	102	102/102 (100%)
Whole Blood (fingerstick)	32	32/32 (100%)

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BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/BIBLIOGRAFIA

1. Piot P, Plummer FA, Mhalu FS, Lamborary JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988; 239: 573-579.
2. Weniger BG, Takebe Y, Ou CY, *et al*. The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994; 8(S2): S13-S28.
3. Gürtler LG, Hauser PH, Eberle J, *et al*. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3): 1581-1585.
4. World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization, 1993.
5. National Commitee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991: 1-43.
6. CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, *MMWR* 1987; 36(2S): 3S-18S.
7. Sehustler LM, Hollinger FB, Dreesman GR, *et al*. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Applied and Environmental Microbiology*. 1981; 42(5): 762-7.
8. Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
9. EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
- 10.Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
11. Delaney KP, Branson BM, *et al*. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
12. O’Connell RJ, Merritt TM, Malia JA, *et al*. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
13. O’Connell RJ, Agan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

Advice Line

For further information, please contact your distributor, or call to one of the following Abbott Technical Support Care Centers:

Region	Phone	E-Mail Address
Europe	+ (44) 161 483 9032	EME.TechSupport@abbott.com
Middle East	+ (965) 2202 2828	EME.TechSupport@abbott.com
Asia Pacific	+ (61) 7 3363 7100	AP.TechSupport@abbott.com
Africa	+ (27) 10 500 9700	arcis.techsupport@abbott.com
Russia & CIS	+ (7) 499 403 9512	arcis.techsupport@abbott.com
Latin America	+ (57) 1 482 4033	LA.TechSupport@Abbott.com

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