WHO Prequalification of Diagnostics Programme PUBLIC REPORT

Product: HIV 1/2 STAT-PAK
Number: PQDx 0007-006-00

The HIV 1/2 STAT-PAK with product code HIV101, manufactured by Chembio Diagnostic Systems, Inc., 3661 Horseblock Road, Medford, NY 11763, United States (US) of America, regulatory version “solely for export outside of the United States of America, except Nigeria”, was accepted for the WHO list of prequalified diagnostics and was listed on 16 January 2012.

The manufacturer also holds a US Food and Drug Administration (FDA) approved regulatory version of the HIV 1/2 STAT-PAK with product code HIV102.

Summary of prequalification status for the HIV 1/2 STAT-PAK

	Date	Outcome
Prequalification listing	16 January 2012	listed
Dossier assessment	9 January 2012	MR
Site inspection(s) of the quality management system	Chembio Diagnostic Systems, Inc.: 8-10 April 2024 Chembio Diagnostics Malaysia Sdn. Bhd.: 8-10 March 2023	MR MR
Product performance evaluation	5 January 2012	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen due to changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendments	Date of report amendment
2.0	Addition of the manufacturing site, Chembio Diagnostics Malaysia (CDM), based in Kuala Lumpur, Malaysia, for the final manufacturing process of the HIV 1/2 STAT-PAK Assay (Product Code HIV101/60-9500-00) as manufactured at the CDS facility based in Medford, NY. The scope of the proposed change included the final assembly process for HIV STAT-PAK, which consists of processes for buffer filling and labelling, card cutting, cassette assembly, packaging, final kitting, quality control (QC), and product release.	7 December 2022
3.0	The extension of Shelf-Life from 24 months to 33 months for the HIV 1/2 STAT-PAK (Product Code HIV101 / 60-9500-0).	25 June 2024

Intended use

According to the claim of intended use from Chembio Diagnostic Systems, *“The HIV 1/2 STAT-PAK is a single-use immunochromatographic rapid screening test for the detection of antibodies to Human Immunodeficiency Virus Types 1 and 2 (HIV 1/2) in fingerstick whole blood, venous whole blood, serum or plasma specimens. The HIV 1/2 STAT-PAK is intended for use as a point-of-care test to aid in the diagnosis of infection with HIV 1 and HIV 2. The HIV 1/2 STAT-PAK is not approved for use to screen blood, plasma, cell or tissue donors.”*

Assay description

According to the claim of intended use from Chembio Diagnostic Systems, *“the HIV 1/2 STAT-PAK employs a unique combination of a specific antibody binding protein, which is conjugated to colloidal gold dye particles, and HIV 1/2 antigens, which are bound to the membrane solid phase. The sample is applied to the sample (S) well, followed by the addition of running buffer. The buffer facilitates the lateral flow of the released products and promotes the binding of antibodies to the antigens. If present, the antibodies bind to the gold conjugated antibody binding protein.*

In a reactive sample, the dye conjugated-immune complex migrates on the nitrocellulose membrane and is captured by the antigens immobilized in the TEST (T) area producing a pink/purple line. In the absence of HIV antibodies, there is no pink/purple line in the TEST (T) area. The sample continues to migrate along the membrane and produces a pink/purple line in the CONTROL (C) area containing immunoglobulin G antigens. This procedural control serves to demonstrate that the specimen and reagents have been properly applied and have migrated through the device.

A reactive result using the HIV 1/2 STAT-PAK suggests the presence of antibodies to HIV 1 or HIV 2 in the specimen. The HIV 1/2 STAT-PAK is intended as an aid in the diagnosis of infection with HIV-1/2. HIV and AIDS-related conditions are clinical syndromes, and their diagnosis can only be established clinically.

A non-reactive result at any time does not preclude the possibility of exposure to HIV or infection with HIV. An antibody response to recent exposure may take several months to reach detectable levels.”

The test kit contains

- 20 individually pouched test devices
- 1 HIV running buffer (3.5 ml)
- 20 disposable 5µl sample loops
- 1 Instructions for use.

Storage

The test kit must be stored at 8 - 30 °C.

Shelf-life upon manufacture

33 months.

Warnings/limitations

Refer to the current version of the manufacturer’s instructions for use attached to this report.

Prioritization for prequalification

Based on the established eligibility criteria, the HIV 1/2 STAT-PAK was given priority for the WHO prequalification assessment.

Dossier assessment

Chembio Diagnostic Systems, Inc. submitted a product dossier for the HIV 1/2 STAT-PAK as per the Instructions for compilation of a product dossier (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by the WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2).

Based on the product dossier screening and assessment findings, the product dossier for the HIV 1/2 STAT-PAK meets the WHO prequalification requirements.

Manufacturing site inspection

An inspection of Chembio Diagnostics Malaysia Sdn. Bhd. located at N-TATT Building 2nd floor, No. 2-1 Jalan TP5, Taman Perindustrian, UEP, 47600, Subang Jaya, Selangor, Malaysia, was conducted from 8-10 March 2023. An inspection of Chembio Diagnostics, Inc. located at 3661 Horseblock Road, Medford, NY 11763, USA and 1560 5th Avenue, Bay Shore, NY 11706, USA was conducted from 8-10 April 2024. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality. Routine inspections of the Manufacturing sites will be conducted with copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current inspection performed at a manufacturing site for *in vitro* diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection can be found at:

https://www.who.int/diagnostics_laboratory/evaluations/PQDxSiteInspection/en/

All published WHOPIRs are with the agreement of the manufacturer.

Product performance evaluation

HIV 1/2 STAT-PAK was evaluated by WHO in the fourth quarter of 2011 at the Institute of Tropical Medicine, Antwerp, Belgium – a WHO Collaborating Centre for HIV/AIDS Diagnostics and Laboratory Support. The laboratory evaluation was conducted according to the “WHO Protocol for the laboratory evaluation of HIV serology assays” (PQDx_030 V1.0) and drew the following conclusions:

HIV 1/2 STAT-PAK is an immunochromatographic rapid diagnostic test for the detection of antibodies to HIV-1/2 in human serum, plasma, and whole blood. A volume of 5µl of serum, plasma or venous/capillary whole blood is required to perform the test procedure. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities. Reading of the results is performed visually, i.e. subjective reading.

In this limited performance evaluation using a panel of 1079 biological specimens, we observed an initial (sensitivity (95% CI) of 99.29% (97.9% - 100%) and an initial specificity (95% CI) of 100% (99.4% - 100%) compared to the reference assays. The final sensitivity (95% CI) was 99.52% (98.3% - 100%), and the final specificity (95% CI) was 100% (99.4% - 100%) compared to the reference assays. In this study, 0% of the overall results were

recorded as indeterminate. Results were interpreted independently by three technicians; the inter-reader variability was 0.18%. The invalid rate was 0%.

Labelling

- 1. Labels**
- 2. Instructions for use**

1. Labels

Kit Box

HIV 1/2 STAT-PAK®

A Qualitative Screening Test Kit for the Detection of Antibodies to HIV 1/2 in Serum, Plasma, or Whole Blood

FOR IN VITRO DIAGNOSTIC USE

FOR PROFESSIONAL USE ONLY

SEE PRODUCT INSERT FOR COMPLETE DIRECTIONS FOR USE

KIT CONTENTS:

- 20 HIV 1/2 Test Devices
- 20 Sample Loops
- 1 Vial Buffer
- 1 Product Insert

CHEMBIO DIAGNOSTIC SYSTEMS, INC.

Manufactured by:

CHEMBIO DIAGNOSTIC SYSTEMS, INC.

3661 Horseblock Road
Medford, New York 11763 USA
info@chembio.com • www.chembio.com

10-6025-0 Rev. 1

Product: HIV101- 20 TESTS

DO NOT FREEZE

STORE AT 8-30°C

The Mid-York Press, Inc.

Since 1966

PROOF

Job# **41141**.....P.O.#.....

Company...**Chembio Diagnostics**.....

Date Submitted.....Date Returned.....

Approved By.....

Further Changes Necessary.....

Remarks:.....

IMPORTANT NOTICE

Final approval is your (you the customer's) responsibility. Please check for all possible errors before signing off on this proof. The Mid-York Press, Inc. will not be responsible for undetected production errors if the work is printed per the customer's approval or if requests for changes are communicated orally.

The Mid-York Press, Inc.

Note: PMS spot colors on color proof may or may not be accurate to Pantone Matching System.

5/18/05

MY-3127 (6-1/2 x 5 x 2-3/4)


Prints:

293 Blue 

and

Spot Aqueous 

(void for glue only)

Running Buffer

Note: There are no labels for the disposable sample loops (unlabelled). The device pouches are not formally labelled but rather include a hot seal printed lot number and expiry dating information as follows:

Lot: 44MMDDYY

Exp: DD MMM YY

2. Instructions for use¹

¹ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

HIV 1/2 STAT-PAK®

FOR *IN VITRO* DIAGNOSTIC USE

FOR PROFESSIONAL USE ONLY

A Qualitative Screening Test Kit for the Detection of Antibodies to HIV-1/2 in
Human Fingertick and Venous Whole Blood, Serum and Plasma

STORAGE: Store at 8 to 30°C (46 to 86°F)

INTENDED USE

The Chembio HIV 1/2 STAT-PAK is a single-use immunochromatographic, rapid screening test for the detection of antibodies to Human Immunodeficiency Virus Types 1 and 2 (HIV 1/2) in fingerstick whole blood, venous whole blood, serum or plasma specimens. The Chembio HIV 1/2 STAT-PAK is intended for use as a point-of-care test to aid in the diagnosis of infection with HIV 1 and HIV 2. This test is suitable for use in multi-test algorithms designed for the statistical validation of rapid HIV test results. When multiple rapid HIV tests are available, this test should be used in appropriate multi-test algorithms.

RESTRICTIONS — The HIV 1/2 STAT-PAK is not approved for use to screen blood, plasma, cell or tissue donors.

SUMMARY AND EXPLANATION

Discovered in 1983, the Human Immunodeficiency Virus (HIV) is a retrovirus identified as the etiologic agent for Acquired Immunodeficiency Syndrome (AIDS) and AIDS related complex [1]. AIDS is characterized by changes in the population of T-cell lymphocytes that play a key role in the immune defense system. In the infected individual the virus causes a depletion of a subpopulation of T-cells, called T-helper cells, which leaves these patients susceptible to opportunistic infections and certain malignancies. The major routes of transmission are sexual contact, contamination by blood or blood products and mother-to-newborn transmission [2-4].

Although there has been a decrease in the rate of infection in certain countries, the number of persons infected with HIV globally has continued to increase. By the end of 2005 there were approximately 40.3 million people living with HIV/AIDS, an increase from nearly 37.5 million in 2003. An estimated 5 million people were newly infected with HIV/AIDS in 2005. In the same year more than 3 million died of AIDS-related illness; more than 500,000 of these were children [5].

HIV infection, AIDS and AIDS related complex have become a leading cause of illness and death, globally, for the past two decades. As of 2009, a total of 1.8 million persons had died of AIDS or AIDS related complex. Approximately 33.3 million persons across the world are infected with HIV and approximately 20% of these persons may not be aware of their infected status. [6]

The HIV virus consists of a genomic RNA molecule protected by a capsid and an envelope. The HIV envelope is the major target for humoral antibody response. The presence of the virus in patients causes the immune system to elicit the production of antibodies to HIV. The detection of these antibodies can be used as a diagnostic tool.

ELISAs, Western Blots, PCR-based and various other test systems are currently available for HIV 1/2 detection [7-11]. The Chembio HIV 1/2 STAT-PAK is a rapid immunochromatographic test, which is simple and easy to use. The Chembio HIV 1/2 STAT-PAK utilizes immobilized antigens for the detection of antibodies to HIV 1/2 in serum, plasma or whole blood.

PRINCIPLE OF THE TEST

The Chembio HIV 1/2 STAT-PAK employs a unique combination of a specific antibody binding protein, which is conjugated to colloidal gold dye particles, and HIV 1/2 antigens, which are bound to the membrane solid phase. The sample is applied to the sample (S) well followed by the addition of running buffer. The buffer facilitates the lateral flow of the released products and promotes the binding of antibodies to the antigens. If present, the antibodies bind to the gold conjugated antibody binding protein. In a reactive sample, the dye conjugated-immune complex migrates on the nitrocellulose membrane and is captured by the antigens immobilized in the TEST (T) area producing a pink/purple line. In the absence of HIV antibodies, there is no pink/purple line in the TEST (T) area. The sample continues to migrate along the membrane and produces a pink/purple line in the CONTROL (C) area containing immunoglobulin G antigens. This procedural control serves to demonstrate that specimen and reagents have been properly applied and have migrated through the device.

MATERIALS PROVIDED

Each kit contains the items to perform 20 tests:

- 20 STAT-PAK Individually Pouched Test Devices
- 1 HIV Running Buffer (3.5mL)
- 20 Disposable 5µL Sample Loops
- 1 Product Insert

MATERIALS AVAILABLE AS AN ACCESSORY TO THE KIT

Chembio Rapid HIV 1/2 Antibody Test Kit Controls (Catalogue # HIV104)

Each package contains:

- HIV 1 Reactive Control
- HIV 2 Reactive Control
- Nonreactive Control
- Product Insert for the HIV104 Kit Controls

MATERIALS REQUIRED BUT NOT PROVIDED

- Clock, watch or other timing device
- Automatic pipettor capable of delivering 5 µL of sample may be used in lieu of the disposable 5 µL sample loop supplied with the kit, for other than fingerstick specimens
- Disposable gloves
- Sterile gauze (for fingerstick whole blood specimens)
- Sterile alcohol swab (for fingerstick samples only)
- Biohazard disposal container
- Sterile lancet (for fingerstick samples only)
- Collection devices for samples other than fingerstick

WARNINGS

For *in-vitro* diagnostic use

1. Read the Product Insert completely before using this assay. Following the instructions carefully as not doing so may result in inaccurate test results.
2. Use of this test kit with specimen types other than those specifically approved for use with this device may result in inaccurate test results.
3. This test should be performed at 18-30°C (64 to 86°F). If stored refrigerated, ensure that the pouch is brought to operating temperature before performing testing.
4. If the test kit is stored at temperatures outside the storage temperature 8-30°C (46 to 86°F), or used outside the operating temperature 18-30°C (64 to 86°F), use the Kit Controls (Catalogue HIV104) to ensure proper performance of the test.
5. Individual(s) infected with HIV-1 and/or HIV-2 who is receiving highly active antiretroviral therapy (HAART) may produce negative results.

PRECAUTIONS

Safety Precautions

1. Do not eat, drink or smoke in the area where specimens and kit reagents are handled. Avoid any contact between hands, eyes or mouth during specimen collection and testing.
2. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when handling patient specimens.
3. Dispose of all specimens and materials used in the test procedure in a biohazard waste container. Lancets should be placed in a puncture-resistant container prior to disposal. The recommended method of disposal of biohazard waste is autoclaving for a minimum of 1 hour at 121°C. Disposable materials may be incinerated. Liquid wastes may be mixed with appropriate chemical disinfectants. A freshly prepared solution of 10% bleach (0.5% solution of sodium hypochlorite) is recommended. Allow 60 minutes for effective decontamination.

NOTE: Do not autoclave solutions that contain bleach.

4. For additional information refer to: Centers for Disease Control (CDC): Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HIV and Recommendations for Postexposure Prophylaxis [12].

Handling Precautions

1. Do not use any device if the pouch has been perforated.
2. If desiccant packet is missing, DO NOT USE, discard test device, and a new test device should be used.
3. Each device is for single use only.
4. Always check expiration date prior to testing. Do not use the test beyond the expiration date printed on the pouch.
5. Do not mix reagents from different lot numbers of kits.
6. Adequate lighting is required to read the test results.

STORAGE AND STABILITY

The HIV 1/2 STAT-PAK test devices should be stored in unopened pouches at 8 to 30°C (46 to 86°F). Do not freeze. Do not use beyond the indicated expiration date. Test devices are stable until the expiration date marked on the pouch, when stored as indicated. Running Buffer should be stored at 8 to 30°C (46 to 86°F) in the original vial.

SPECIMEN COLLECTION

The Chembio HIV 1/2 STAT-PAK test can be performed on fingerstick whole blood, venous whole blood, serum or plasma specimens. All specimens should be collected following local clinical or laboratory procedures.

Fingerstick Whole Blood

Prepare to perform the fingerstick blood collection procedure. Clean the finger of the person being tested with an antiseptic wipe. Allow the finger to dry thoroughly or wipe dry with a sterile gauze pad.

Using sterile lancet, following laboratory procedure; prick the finger just off the center of the finger and wipe away the first drop with sterile gauze. Avoid squeezing the fingertip to accelerate bleeding as this may dilute the blood with excess tissue fluid.

Collect the sample from the second drop by touching the disposable sample loop provided to the drop of blood until the sample loop is full. Test immediately, following test procedure instructions.

Venous Whole Blood

Draw blood following laboratory procedure for obtaining venous blood. Depending on use, collect sample in a tube containing citrate, heparin or EDTA. Be sure the tube of blood is well mixed before sampling. Dip the sample loop into the blood and allow it to fill. Follow test procedure instructions.

If tested the same day, venous whole blood may be kept at room temperature. Venous whole blood may be stored for up to 3 days between 2 and 8°C (36 to 46°F) before testing.

DO NOT FREEZE WHOLE BLOOD! Allow for refrigerated sample to reach room temperature and mix gently before testing.

Serum or Plasma

Draw blood following laboratory procedures for obtaining serum or plasma specimens. Collect specimen in a clean tube not containing any anticoagulant (serum) or in a tube containing citrate, heparin or EDTA (plasma), following standard laboratory procedures. All specimens should be centrifuged following local clinical or laboratory procedures.

Patient samples perform best when tested immediately after collection. If not tested immediately, specimens should be refrigerated at 2 to 8°C (36 to 46°F) and can be used up to 3 days after collection. If testing within 3 days is not possible, serum or plasma specimens should be frozen at -20°C (-4°F) or colder.

NOTE: If specimens are to be shipped, they should be packed in compliance with regulations covering the transportation of etiologic agents. Venous whole blood, serum and plasma specimens should be shipped refrigerated with cold packs or wet ice.

TEST PROCEDURE


If the specimen to be tested is refrigerated, remove it from the refrigerator and allow it to come to a temperature of 18 to 30°C (64 to 86°F) prior to testing.

1. Remove the Chembio HIV 1/2 STAT-PAK test device from its pouch and place it on a flat surface (it is not necessary to remove the desiccant from the pouch).

NOTE: If desiccant packet is missing, DO NOT USE. Discard the test device and a new test device should be used.

2. Label the test device with patient name or identification number. (See Figure 1 below)

Figure 1

3. Touch the 5 μ L sample loop provided to the specimen, allowing the opening of the loop to fill with the liquid.
4. Holding the sample loop vertically, touch it to the sample pad in the center of the SAMPLE (S) well of the device to dispense $\sim$ 5 μ L of sample (serum, plasma or whole blood) onto the sample pad. (See Figure 2 below)

Figure 2

5. Invert the Running Buffer bottle and hold it vertically (not at an angle) over the sample well. Add 3 drops ($\sim$ 105 μ L) of buffer slowly, dropwise, into the SAMPLE (S) well. (See Figure 3 below)

Figure 3

6. Read the test result 15 minutes after the addition of the Running Buffer. In some cases a test line may appear in less than 15 minutes however, 15 minutes are needed to

report a nonreactive result. Read results in a well-lit area. **Do not read results after 20 minutes.**

NOTE: Discard the used sample loop, test device and any other test materials into a biohazard waste container. Do not re-use loops.

QUALITY CONTROL

Built-in Control Feature

When the test is complete, you will see a pink/purple line in the CONTROL (C) area of the test device on nonreactive as well as reactive samples. This control line serves as an internal control and gives confirmation of sample addition and proper test performance. A pink/purple line will appear in the CONTROL (C) area if the test has been performed correctly and the device is working properly (Please see section: INTERPRETATION OF RESULTS).

External Quality Control

Good Laboratory Practices necessitates testing external control material along with the test samples to ensure proper performance of the test kit. Chembio HIV Reactive and Nonreactive Controls (Catalogue # HIV104) are available separately for use with the Chembio HIV 1/2 STAT-PAK test. The HIV Controls are used to verify the operator's ability to properly perform the test and to interpret the results. The Reactive Control will produce a reactive test result and has been manufactured to produce a faint line in the TEST (T) area. The Nonreactive Control will produce a nonreactive test result. Run the controls as per the TEST PROCEDURE and INTERPRETATION OF RESULTS sections of this insert. It is the responsibility of each facility using the HIV 1/2 STAT-PAK to establish an adequate quality assurance program to ensure the performance of the device under specific locations and conditions.

RUN THE KIT CONTROLS UNDER THE FOLLOWING CIRCUMSTANCES:

- Each new operator prior to performing tests on patient specimens,
- When opening a new test Kit lot,
- Whenever a new shipment of test Kits is received,
- If the temperature of the test storage area falls outside of 8 to 30°C (46 to 86°F),
- If the temperature of the testing area falls outside of 18 to 30°C (64 to 86°F),
- At periodic intervals as indicated by the user facility.

If the HIV Control reagents do not produce the expected results, contact Chembio Diagnostic Systems' Customer Service (+1-631-924-1135).

INTERPRETATION OF RESULTS

Nonreactive Results

Control Line

No Test Line

Nonreactive Result:

One pink/purple line in the CONTROL (C) area, with no line in the TEST (T) area indicates a NONREACTIVE test result. A NONREACTIVE test result at 15 minutes indicates that HIV-1 and HIV-2 antibodies were not detected in the specimen. Test result is interpreted as NEGATIVE for HIV-1 and HIV-2 antibodies. However, this result does not exclude the possibility of HIV infection.

Reactive Results

Control Line

Reactive Test Line


Reactive Result:

Two pink/purple lines, one in the TEST (T) area and one in the CONTROL (C) area indicate a REACTIVE test result. The line in TEST (T) area may look different from the line in the CONTROL (C) area. Intensities of the Test and Control Lines may vary. Test result with visible lines in both TEST (T) and CONTROL (C) areas, regardless of intensity, is considered REACTIVE. A REACTIVE result means that HIV-1 and/or HIV-2 antibodies have been detected in the specimen. The test result is interpreted as Preliminary POSITIVE for HIV-1 and/or HIV-2 antibodies.

NOTE: Intensities of the TEST (T) and CONTROL (C) lines may vary. If any visible line appears in the TEST (T) area and in the CONTROL (C) area, the result is reactive.

Invalid Results

Invalid Result:

A pink/purple line should always appear in the CONTROL (C) area, whether or not a line appears in the TEST (T) area. If there is no distinct pink/purple line visible in the CONTROL (C) area, then the test is INVALID. Any line that appears outside of the CONTROL (C) area or TEST (T) area is an INVALID test. An INVALID test indicates a problem with running the test, either related to the specimen or to the device and cannot be interpreted. It is recommended that the INVALID test be repeated with a new device. Contact Chembio Diagnostic Systems at +1-631-924-1135 ext 114 or info@chembio.com if you are unable to obtain a valid result upon repeat testing.

LIMITATIONS OF THE PROCEDURE

1. The Chembio HIV 1/2 STAT-PAK must be used with capillary (fingerstick) or venous whole blood, serum or plasma only. Using other types of specimens or testing of venipuncture whole blood specimens collected using a tube containing an anticoagulant other than citrate, heparin or EDTA may not yield accurate results. For serum samples, collect blood without anticoagulant.
2. The Chembio HIV 1/2 STAT-PAK must be used in accordance with the instructions in this product insert to obtain accurate results.
3. Be careful not to add more than 3 drops of Running Buffer with blood sample as it may lead to the appearance of a red line at the base of the window.
4. Reading test results earlier than 15 minutes or later than 20 minutes may yield erroneous results.
5. Test devices should be used immediately upon opening the sealed foil pouch.
6. Do not use kit contents beyond labeled expiration date.
7. Ensure finger is completely dry before performing fingerstick.
8. Read results in a well-lit area.
9. A reactive result using the Chembio HIV 1/2 STAT-PAK suggests the presence of antibodies to HIV-1 or HIV-2 in the specimen. The Chembio HIV 1/2 STAT-PAK is intended as an aid in the diagnosis of infection with HIV-1/2. HIV and AIDS-related conditions are clinical syndromes and their diagnosis can only be established clinically.
10. For a reactive result, the intensity of the test line does not necessarily correlate with the titer of antibody in the specimen.
11. A nonreactive result does not preclude the possibility of exposure to HIV or infection with HIV. An antibody response to recent exposure may take several months to reach detectable levels.
12. An individual infected with HIV-1 and/or HIV-2 who is receiving highly active antiretroviral therapy (HAART) may produce false negative results.
13. This assay has not been evaluated for newborn screening, cord blood specimens, or for individuals less than 13 years of age or older than 64 years of age.

EXPECTED RESULTS

This is a qualitative test for the detection of antibodies to HIV 1/2 in whole blood, serum or plasma. As described in the PERFORMANCE CHARACTERISTICS section below, the sensitivity of the Chembio HIV 1/2 STAT-PAK was found to be substantially equivalent to EIA and Western Blot tests when tested on selected performance panels.

PERFORMANCE CHARACTERISTICS

In-house studies demonstrate that the sensitivity of the Chembio HIV 1/2 STAT-PAK is substantially equivalent to the EIA and Western Blot test when tested on BBI performance panels PRZ204 Anti-HIV 1 / 2 Combo Performance Panel, PRF202 Anti-HIV 2 Performance Panel and PRB203 Anti-HIV 1 Mixed Titer Performance Panel as well as BBI Seroconversion panels PRB904 and PRB909.

In an external evaluation of the performance of the HIV 1/2 STAT-PAK using 336 confirmed nonreactive and reactive sera, plasma and whole blood samples, sensitivity was 100% (129/129) and specificity, 100% (207/207). This study also included 34 fingerstick blood samples; 29 reactive and 5 nonreactive. In this same study, excellent analytical sensitivity relative to EIA was demonstrated using BBI seroconversion panels, PRB940 and PRB931.

EFFECT OF POTENTIALLY INTERFERING SUBSTANCES AND UNRELATED MEDICAL CONDITIONS

To evaluate the influence of unrelated medical conditions or interfering substance on the specificity and sensitivity of the Chembio HIV 1/2 STAT-PAK; 208 specimens representing unrelated medical conditions and 110 specimens representing potential interfering substances were tested (Table 1). The specimens were spiked with either saline (Nonreactive) or an HIV-1 reactive serum specimen to a low level of reactivity. All HIV-1 spiked specimens gave reactive results while all unspiked samples, with the exception of one elevated albumin specimen and 14 syphilis specimens, gave nonreactive results. The one elevated albumin specimen and all of the 14 unspiked syphilis specimens with reactive results were subsequently confirmed as infected with HIV-1 using a licensed Western Blot assay. An additional ten known HIV-1 nonreactive, syphilis reactive specimens were tested and yielded expected results.

Table 1

Chembio HIV 1/2 STAT-PAK reactivity against specimens from unrelated medical conditions or containing potential interfering substances

Chembio HIV 1/2 STAT-PAK		
Description	Saline (Nonreactive)	HIV-1/2 (Weak Reactive)
Cirrhosis	20/20	20/20
CMV IgM	20/20	20/20
Recent flu vaccination ¹	11/11	11/11
HBV	21/21	21/21
HCV	19/19	19/19
HTLV-I	11/11	11/11
HTLV-II	10/10	10/10
Multiparous	9/9	9/9
Myeloma	10/10	10/10
Rheumatoid Factor	10/10	10/10
Syphilis ²	15/29	29/29
Tuberculosis	38/38	38/38
Elevated Albumin ³	9/10	10/10
Elevated Bilirubin	10/10	10/10
Citrate	10/10	10/10
DNA	10/10	10/10
EDTA	10/10	10/10
Hemolyzed	10/10	10/10
Heparin	10/10	10/10
Icteric	10/10	10/10
Lipemic	10/10	10/10
Elevated Protein	10/10	10/10
Elevated Triglycerides	10/10	10/10

¹Collected within 6 months of vaccination

²Fourteen samples were confirmed reactive, using a licensed WB assay

³One sample was confirmed as containing HIV antibodies by using a licensed WB assay

PRECISION

Intraassay

Within run precision was determined by using 10 replicates of two specimens containing different levels of HIV 1/2 antibodies. The nonreactive and reactive results were correctly identified 100% of the time.

Interassay

Between run precision was determined by using the same two specimens in 10 different replicates from three different lots of test devices. Again nonreactive and reactive results were observed 100% of the time.

Cross Reactivity and Interference

No cross reactivity was observed from Hepatitis B, Rheumatoid Factor (~80 IU/ml) and hCG (500 mIU/ml) with Chembio HIV 1/2 STAT-PAK. In addition, no interference from bilirubin, hemoglobin and triglycerides was observed.

REFERENCES

1. Essex, M. (1999) Human immunodeficiency viruses in the developing world. *Adv Virus Res* 53 : 71-88.
2. Kanki, P.J., Hopper, J.R. and Essex, M. (1987) The origins of HIV-1 and HTLV-4/HIV-2. *Ann N Y Acad Sci* 511 : 370-375.
3. Nicoll, A., Gill, O.N. (1999) The global impact of HIV infection and disease. *Commun Dis Publ Health* 2 : 85-95.
4. Valdiserri R.O., Holtgrave, D.R., West, G.R. (1999) Promoting early diagnosis and entry into care. *AIDS* 13 : 2317-2330.
5. UNAIDS / WHO Press Release (2005) HIV Infection Rates Decreasing in Several Countries but Globally Number of People Living with HIV Continue to Rise. Available at www.unaids.org/epi/2005/doc/docs/PR_EpiUpdateNov05en.pdf: 1-3.
6. World Health Organization (WHO): Updated HIV surveillance, estimates, monitoring and evaluation. <http://www.who.int/hiv/topics/me/en/index.html>
7. Essex, M., Kanki, P. J., Marlink, R., et al. (1990) Antigenic characterization of the human immunodeficiency viruses. *J Am Acad Dermatol* 22 : 1206-1210.
8. Essex, M., McLane, M.F., Lee, T.H. et al. (1983) Antibodies to cell membrane antigens associated with human T-cell leukemia virus in patients with AIDS. *Science* 220 : 859-862.
9. Gallo, R.C., Saluahuudin, S.Z., Popovic, M., et al. (1984) Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 224 : 500-503.
10. Kenealy, W., Reed, D., Cybulsky, R., et.al. (1987) Analysis of human serum antibodies to human immunodeficiency virus (HIV) using recombinant ENV and GAG antigens. *AIDS Res Human Retrovir* 3 : 95-105.
11. Kovacs, A., Xu, J., Rasheed, S., et al. (1995) Comparison of a rapid non-isotopic polymerase chain reaction assay with four commonly used methods for the early diagnosis of human Immunodeficiency Virus type 1 infection in neonates and children. *Pediatr Infect Dis J* 14 : 948-954.
12. Centers for Disease Control (CDC): Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HIV and Recommendations for Postexposure Prophylaxis. *MMWR* 2005; 54(RR-9): 1-17.

Available Accessories:

Catalogue # HIV104: HIV Reactive and Nonreactive Controls

***This product insert is similar to the FDA approved Product Insert.
FDA version of Product Insert is only available in the United States.***

***Certain language requirements specific to clinical laboratories in the
United States have been removed.***

Manufactured by:
CHEMBIO DIAGNOSTIC SYSTEMS, INC.
3661 HORSEBLOCK ROAD
MEDFORD, NY 11763, USA

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1. Essex, M. (1999) Human immunodeficiency viruses in the developing world. *Adv Virus Res* 53 : 71-88.
2. Kanki, P.J., Hopper, J.R. y Essex, M. (1987) The origins of HIV-1 and HTLV-4/HIV/2. *Ann N Y Acad Sci* 511 : 370-375.
3. Nicoll, A., Gill, O.N. (1999) The global impact of HIV infection and disease. *Commun Dis Publ Health* 2 : 85-95.
4. Valdiserri R.O., Holtgrave, D.R., West, G.R. (1999) Promoting early diagnosis and entry into care. *AIDS* 13 : 2317-2330.
5. Artículo de prensa de ONUSIDA/ OMS (2005) HIV Infection Rates Decreasing in Several Countries but Globally Number of People Living with HIV Continue to Rise. Disponible en www.unaids.org/epi/2005/doc/docs/PR_EpiUpdateNov05en.pdf: 1-3.
6. ORGANIZACIÓN MUNDIAL DE LA SALUD (OMS): Updated HIV surveillance, estimates, monitoring and evaluation. <http://www.who.int/hiv/topics/me/en/index.html>
7. Essex, M., Kanki, P. J., Marlinsk, R., et al. (1990) Antigenic characterization of the human immunodeficiency viruses. *J Am Acad Dermatol* 22 : 1206-1210.
8. Essex, M., McLane, M.F., Lee, T.H. et al. (1983) Antibodies to cell membrane antigens associated with human T-cell leukemia virus in patients with AIDS. *Science* 220 : 859-862.
9. Gallo, R.C., Saluahuiddin, S.Z., Popovic, M., et al. (1984) Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 224 : 500-503.
10. Kenealy, W., Reed, D., Cybulsky, R., et al. (1987) Analysis of human serum antibodies to human immunodeficiency virus (HIV) using recombinant ENV and GAG antigens. *AIDS Res Human Retrovir* 3 : 95-105.
11. Kovacs, A., Xu, J., Rasheed, S., et al. (1995) Comparison of a rapid non-isotopic polymerase chain reaction assay with four commonly used methods for the early diagnosis of human Immunodeficiency Virus type 1 infection in neonates and children. *Pediatr Infect Dis J* 14 : 948-954.
12. Centers for Disease Control (CDC): Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HIV and Recommendations for Postexposure Prophylaxis. *MMWR* 2005; 54(RR-9): 1-17.

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RÉFÉRENCES

1. Essex, M. (1999) Human immunodeficiency viruses in the developing world. *Adv Virus Res* 53 : 71-88.
2. Kanki, P.J., Hopper, J.R. and Essex, M. (1987) The origins of HIV-1 and HTLV-4/HIV-2. *Ann N Y Acad Sci* 511 : 370-375.
3. Nicoli, A., Gill, O.N. (1999) The global impact of HIV infection and disease. *Commun Dis Publ Health* 2 : 85-95.
4. Valdiserri R.O., Holtgrave, D.R., West, G.R. (1999) Promoting early diagnosis and entry into care. *AIDS* 13 : 2317-2330.
5. UNAIDS / WHO Press Release (2005) HIV Infection Rates Decreasing in Several Countries but Globally Number of People Living with HIV Continue to Rise. Disponible à www.unaids.org/epi/2005/doc/docs/PR_EpiUpdateNov05en.pdf : 1-3.
6. Organisation mondiale de la Santé (OMS) : Updated HIV surveillance, estimates, monitoring and evaluation. <http://www.who.int/hiv/topics/me/en/index.html>
7. Essex, M., Kanki, P. J., Marlink, R., et al. (1990) Antigenic characterization of the human immunodeficiency viruses. *J Am Acad Dermatol* 22 : 1206-1210.
8. Essex, M., McLane, M.F., Lee, T.H. et al. (1983) Antibodies to cell membrane antigens associated with human T-cell leukemia virus in patients with AIDS. *Science* 220 : 859-862.
9. Gallo, R.C., Saluahuiddin, S.Z., Popovic, M., et al. (1984) Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 224 : 500-503.
10. Kenealy, W., Reed, D., Cybulsky, R., et.al. (1987) Analysis of human serum antibodies to human immunodeficiency virus (HIV) using recombinant ENV and GAG antigens. *AIDS Res Human Retrovir* 3 : 95-105.
11. Kovacs, A., Xu, J., Rasheed, S., et al. (1995) Comparison of a rapid non-isotopic polymerase chain reaction assay with four commonly used methods for the early diagnosis of human Immunodeficiency Virus type 1 infection in neonates and children. *Pediatr Infect Dis J* 14 : 948-954.
12. Centres de contrôle et de prévention des maladies (CDC) : Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HIV and Recommendations for Postexposure Prophylaxis. *MMWR* 2005; 54(RR-9) : 1-17.

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