## WHO Prequalification of In Vitro Diagnostics Programme PUBLIC REPORT

**Product: Bioline HIV 1/2 3.0<sup>1</sup>**  
**Number: PQDx 0027-012-00**

**Bioline HIV 1/2 3.0** with product codes **03FK10, 03FK16 and 03FK17<sup>2</sup>**, manufactured by **Abbott Diagnostics Korea Inc<sup>3</sup>**, **Rest of the World** regulatory version was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 20 May 2013.

### Summary of prequalification status for Bioline HIV 1/2 3.0

	Date	Outcome
<b>Status on PQ list</b>	20-May-2013	listed
<b>Dossier assessment</b>	11-Aug-2011	MR
<b>Site inspection(s) of quality management system</b>	19-Feb-2013	MR
<b>Product performance evaluation</b>	05-Apr-2013	MR

MR: Meets Requirements

### Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

<sup>1</sup> SD BIOLINE HIV-1/2 3.0 to Bioline HIV 1/2 3.0

<sup>2</sup> Product codes that were initially listed were 03FK10, 03FK16. Product code 03FK17 was added after acceptance of a change request.

<sup>3</sup> change of manufacturer's name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc.

Version	Summary of amendment	Date of report amendment
2.0	Review of text updates in the public report	2013
3.0	Change of inner structure of device to capture a reasonable excess volume of buffer that might be added by the end user, thereby reducing the risk of overflow of buffer.	2013
4.0	Addition of a product code, which corresponds to the IVD supplied with the new safety lancets. Product code is 03FK17, 25 tests per kit, with 25 safety lancets	23-Feb-2017
5.0	Change of product name from SD BIOLINE HIV-1/2 3.0 to Bioline HIV 1/2 3.0 and change of manufacturer's name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc.	20-Aug-2020

#### Intended use<sup>4</sup>:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“The Bioline HIV 1/2 3.0 kit is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 and HIV-2 simultaneously in human serum, plasma or whole blood. The Bioline HIV 1/2 3.0 kit is intended only for professional use and for in vitro diagnostic use. This test may not be suitable for diagnosis of early infection or blood donation screening. Positive samples should be confirmed by a supplemental assay such as ELISA or Western Blot test”*.

#### Assay description:

According to the claim of assay description from Abbott Diagnostics Korea Inc, *“BiolineHIV 1/2 3.0 test contains a membrane strip, which is precoated with recombinant HIV-1 capture antigen (gp41, p24) on test line 1 region and with recombinant HIV-2 capture antigen (gp36) on test line 2 region respectively. The recombinant HIV 1/2 antigen (gp41, p24 and gp36)-colloid gold conjugate and the sample move along the membrane chromatographically to the test region (T) and forms a visible line as the antigen-antibody-antigen gold particle complex forms with high degree of sensitivity and specificity”*.

<sup>4</sup> This product is one that uses Protein A to detect human IgG antibodies. Protein is also able to detect other classes of human antibody (IgA, IgD, IgE and IgM) but not as reliably as it does IgG. This product has been prequalified with respect to its ability to detect human IgG antibodies. Any claim to detect other types of antibodies on this kind of product has not been validated based on WHO prequalification requirements.

**Test kit contents:**

	<b>30T/kit (product code 03FK10)</b>	<b>25T/kit (product code 03FK16)</b>	<b>25T/kit (product code 03FK17)</b>
<b>Test cassettes</b> individually packed in foil pouch with a desiccant	30 test devices	25 test devices	25 test devices
<b>Assay diluent</b> dispensed in plastic bottle	1 x 4ml/bottle	1 x 4ml/bottle	1 x 4ml/bottle
<b>Specimen transfer devices</b> Disposable (20µl)	N/A	25 units of 20 µl	25 units of 20 µl
<b>Lancets</b> Disposable, sterilized	N/A	25 units	25 units <b>(safety lancet)</b>
<b>Alcohol swabs</b> Disposable	N/A	25 units	25 units
<b>Instructions for use</b>	1 copy	1 copy	1 copy

**Storage:**

The test kit should be stored at 1 - 30 °C.

**Shelf-life:**

24 months.

**Warnings/limitations:**

1. The reading time for this product was changed, the revised instructions for use now state: "Time to result is 10 to 20 minutes. After adding the diluent, read the result after 10 minutes but not more than 20 minutes."
2. If the test result is not legible after 10 minutes due to high background colour, read again later but within 20 minutes of adding the diluent. Do not read after 20 minutes.
3. Dual infection of HIV-1 and HIV-2 within one individual is quite rare. Dual reactivity observed in Bioline™ HIV 1/2 3.0, i.e. HIV-1 line and HIV-2 line both reactive, is more likely to be caused by cross-reactivity given certain homology in the amino acid sequences of HIV-1 and HIV-2. To determine the virus type or diagnose a co-infection, confirmatory testing must be performed.

**Prioritization for prequalification**

Based the established eligibility criteria, Bioline HIV 1/2 3.0 was given priority for prequalification assessment.

## Dossier assessment

Abbott Diagnostics Korea Inc submitted a product dossier for Bioline HIV 1/2 3.0 as per the “Instructions for compilation of a product dossier” (PQDx\_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO.

The manufacturer’s response to the nonconformities found during the dossier assessment were accepted on 23 November 2011.

Commitments for prequalification:

1. Analytical performance studies
2. Clinical performance studies
3. Stability studies
4. A new version of the labels and instructions for use.

Commitments for prequalification are under review.

Based on the product dossier assessment findings, the product dossier for Bioline HIV 1/2 3.0 meets WHO prequalification requirements.

## Manufacturing site inspection

A comprehensive second re-inspection was performed at the sites of the legal manufacturer of Bioline HIV 1/2 3.0 at 156-68 Hagal-dong Giheung-gu, Yongin-si, Kyonggi-do 446-930, Republic of Korea and 473-4 Bora-dong Giheung-gu, Yongin-si, Kyonggi-do, 446-904, Republic of Korea in November 2012<sup>5</sup>, as per “Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics.” (PQDx\_014 v1).

Note: the inspection team were not able to review the batch manufacturing records of the lots submitted for the repeat WHO laboratory evaluation for the Bioline HIV 1/2 3.0. The lots for retesting had not been requested at the time of the inspection (November 2012) and were submitted subsequent to the inspection.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted 19 February 2013.

---

<sup>5</sup> A subsequent inspection took place at Production: 65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do 17099 Republic of Korea and 46, Hagalro 15 beon-gil, Giheung-gu, Yongin-si, Gyeonggi-do 17099 Republic of Korea and Warehouse: 19-22, Dongtansandan 3-gil, Dongtan-myeon, Hwaseong-si, Gyeonggi-do 18487, Republic of Korea on 6 to 8 May 2015.

**Commitments for prequalification:**

1. The manufacturer has committed to continuing improvements in the quality management system particularly in the areas of clear lines of authority, identification and traceability, warehousing and clarity of work instructions and batch manufacturing records.
2. The manufacturer has committed to continuing close supervision of the lot release procedures together with ongoing communication over time to finalize any outstanding issues noted in the WHO responses to the inspection findings.

Based on the site inspection and corrective action plan review, the quality management system for BiolineHIV 1/2 3.0 meets WHO prequalification requirements.

**Product performance evaluation**

Bioline HIV 1/2 3.0 was evaluated by WHO at the Institute of Tropical Medicine, Antwerp, Belgium, in the last quarter of 2012 using serum/plasma specimens. From this evaluation, we drew the following conclusions:

Bioline HIV 1/2 3.0 is a lateral flow immunochromatographic assay for the discriminatory detection HIV-1 and HIV-2 antibodies in human serum/plasma and whole blood. A volume of 10µL of serum/plasma or 20µl of whole blood is required to perform the assay. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities and non-laboratory settings. Reading of the results is performed visually i.e. subjectively read.

In this limited evaluation on a panel of 1118 clinically-derived specimens, we found an initial sensitivity (95% CI) of 99.8% (98.8% - 100%) and an initial specificity (95% CI) of 99.9% (99.2% - 100%) compared to the reference assays (Vironostika HIV Ag/ab [bioMérieux] and Enzygnost Anti-HIV 1/2 in parallel; followed by INNO-LIA HIV H/II Score [Innogenetics] ). The final sensitivity (95% CI) was 100% (99.2% - 100%) and the final specificity (95% CI) was 99.9% (99.2% - 100%) compared to the reference assays. Lot to lot variation observed was within the acceptance criteria.

Bioline HIV 1/2 3.0 was unable to discriminate between HIV-1 and HIV-2 for seven HIV-2 specimens, and 22 HIV-1 specimens (6.3% of 460 HIV positive specimens), as two test bands of equal intensity were observed.

For eight seroconversion panels, Bioline HIV 1/2 3.0 detected on average 0.125 specimens later than the benchmark assay (Enzygnost Anti-HIV 1/2 Plus [Siemens Healthcare Diagnostics]).

For the mixed titer panel, Bioline HIV 1/2 3.0 correctly classified all specimens. For the 1<sup>st</sup> International Reference Panel for anti-HIV [NIBSC code 02/210], Bioline HIV 1/2 3.0 detected all HIV-1 subtypes tested (HIV-1 A, HIV-1 B, HIV-C, HIV-1 CRF01\_AE, and HIV-2).

In this study, 0% of the results were recorded as indeterminate. Results were interpreted independently by three technicians; the overall inter-reader variability was 1.9% (0.2% for HIV-1 band, 1.8% for HIV-2 band). The invalid rate was 0%.

Based on these results, the performance evaluation for Bioline HIV 1/2 3.0 meets the WHO prequalification requirements.


## Labelling

1. **Labels**
2. **Instructions for use**

*\*Note: Labelling has been changed as per site name change and product rebranding. But temporarily labelling of legacy brand (SD BIOLINE) will be used in the market according to registration status in each country.*


# 1. Labels

## 1.1 Package box for 03FK10


### 1.2 Package box for 03FK16

### 1.3 Package box for 03FK17

HIV 1/2 3.0

Bioline™

PRODUCT - Bioline™ HCV 25T,  
SAFETY LANCET

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PRINT AREA

**CONTENTS:**

- 25 uses de lancet with desiccant in individual foil pouches
- Asay diluent (1 x 4 ml/1oz)
- 25 capillary pipettes (20 µl)
- 25 safety lancets
- 25 alcohol swabs
- 1 instructions for use

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L19701-00-A0

REF 03FK17

Abbott Diagnostic Reagents Inc.  
6500 Westpark Dr, Irving, TX 75039-2298, USA  
Copyright © 2019, Division of Abbott Diagnostics

© 2019 Abbott. All rights reserved. All trademarks, registered and trademarks of Abbott Diagnostic Group of companies or their respective owners.

## **2. Instructions for use<sup>6</sup>**

---

<sup>6</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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

## PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

- EN Now, open the package and look for the following:

  - Test device with desiccant in individual foil pouch
  - Assay diluent
  - Instructions for use
  - Capillary pipette (20 µl)
  - Safety lancet
  - Alcohol swab

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- EN Carefully read the instructions on how to use the Bioline™ HIV 1/2 3.0 test kit.

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- EN Look at the expiration date at the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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- EN Open the foil pouch and look for the following:

  - Test device
  - Desiccant

Then, label the device with the patient identifier.

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


- EN Dispense 10 µl of plasma or serum specimen or 20 µl of whole blood specimen into the specimen well marked "S"

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
- EN Dispense 4 drops (approximately 120 µl) of assay diluent into the specimen well "S". Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.

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
- EN Interpret the result 10 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 10 min or after 20 min) may provide in false results.

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- EN Clean the area to be lanced with an alcohol swab.

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

- EN Immerse the open end of a new capillary pipette (20 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

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- EN Dispense 20 µl of drawn whole blood specimen in the specimen well marked "S".

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PT Coloque 20 µl da amostra de sangue total colhido no poço da amostra com a marca "S".
- EN Dispense 4 drops (approximately 120 µl) of assay diluent into the specimen well "S". Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.

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