WHO Prequalification of Diagnostics Programme PUBLIC REPORT

Product: Bioline Malaria Ag P.f ¹
PQDx Reference Number: PQDx 0031-012-01

BiolineMalaria Ag P.f with product code **05FK50, 05FK53, 05FK51** and **05FK52**, manufactured by **Abbott Diagnostics Korea Inc²**, **CE-mark regulatory version** was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 06 December 2010.

Summary of prequalification status for Bioline Malaria Ag P.f

	Date	Outcome
Status on PQ list	06-Dec-2010	listed
Dossier assessment	02-Dec-2010	MR
Site inspection(s) of quality management system	Jun-2017 Jul-2015 Nov-2012 Mar-2012 Sep-2010	MR
Product performance evaluation	Round 1 (2009) and Round 5 (2014)	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
1.0-3.0	Editorial amendments and addition of a product code 05FK53	2-Feb-2016
4.0	Changes to the labelling of the cassette	2-May-2016
5.0	Addition of product codes 05FK51 and 05FK52	22-Dec-2016

¹ Product name was changed from SD BIOLINE Malaria Ag P.f and SD BIOLINE Malaria Ag P.f POCT to Bioline Malaria Ag P.f

² Manufacturer's name changed from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc

6.0	Changes to the manufacturer name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc and product name from SD BIOLINE Malaria Ag P.f and SD BIOLINE Malaria Ag P.f POCT to Bioline Malaria P.f	20-Aug-2020
7.0	Addition of affiliate company of Abbott Group (Rapid test manufacturing site in China) for device assembly, pouch sealing and packing process.	15-Jan-2021

Intended use:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of Malaria Plasmodium falciparum in human whole blood. Bioline Malaria Ag P.f test is intended for professional use, only for an initial screening test and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear”*

Test principle:

According to the claim of assay description from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f test device contains a membrane strip, which is precoated with mouse monoclonal antibodies specific to HRP2 of P. falciparum on test line “P.f” region. The mouse monoclonal antibodies specific to HRP2 of P.f - colloid gold conjugate reacts with the Malaria Plasmodium falciparum antigen in the specimen. They move along the membrane chromatographically to the test region “P.f” and form a visible line as the antibody-antigen-antibody gold particle complex with high degree of sensitivity and specificity”*

Test kit contents:

	25T/kit (product code 05FK50)	1T/kit x 25ea (product code 05FK53)	25T/kit (product code 05FK51)	1T/kit x 25ea (product code 05FK52)
Test cassettes individually packed in foil pouch with a desiccant	25 test devices	25x 1 test device	25 test devices	25x 1 test device
Assay diluent dispensed in plastic bottle	1 x 5ml/bottle	25x 180µl/vial	1 x 5ml/bottle	25x 180µl/vial
Specimen transfer devices disposable (5µl)	25 units of 5µl	25 pouches each containing one specimen transfer device (5µl), one lancet, one alcohol swab	25 units of 5µl	25 pouches each containing one specimen transfer device (5µl), one safety lancet , one alcohol swab
Lancets (optional) disposable, sterile	25 units		25 units (safety lancets)	
Alcohol swabs (optional) disposable	25 units		25 units	
Instructions for use	1 copy	1 copy	1 copy	1 copy
Summarized instructions for use	N/A	25 copies	N/A	25 copies

Storage:

The test kit should be stored at 1 – 40 °C.

Shelf-life:

24 months.

Prioritization for prequalification

Based on the WHO product testing results from Round 5, Bioline Malaria Ag P.f was given priority for prequalification assessment.

Dossier assessment

Abbott Diagnostics Korea Inc submitted a product dossier for Bioline Malaria Ag P.f as per the *“Instructions for compilation of a product dossier”* (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2). Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for Bioline Malaria Ag P.f for prequalification.

Commitments for prequalification:

The manufacturer committed to amend and submit additional documentation on the following issues which will be reviewed at the next re-inspection:

1. the risk analysis and control summary reflecting use in resource-limited settings
2. analytical and stability studies
3. revised labels and instructions for use.

The commitments are under review.

Manufacturing site inspection

A comprehensive third re-inspection³ was performed at at the sites of the legal manufacture (Production: 65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, Republic of Korea 446-930 and Warehouse: 19-22, Dongtansandan 3-gil, Dongtan-myeon, Hwaseong-si, Gyeonggi-do, Republic of Korea)in May 2015, as per the *“Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics”* (PQDx_014 v1).

The inspections were based on ‘ISO 13485:2003 Medical devices - Quality management systems - Requirements for regulatory purposes’ and other internationally recognized standards relevant to the manufacture of in vitro diagnostics. In addition, the claims made in the submitted product dossier were verified, particular attention was paid to suitability of product labelling currently in use (including instructions for use and storage requirements), stability testing (in-use, transportation and storage stability), effective mechanisms for customer training, service and feedback, and the adequacy of mechanisms for lot release of the product to customers.

³ Previous site inspections were carried out in September 2010, March 2012, and November 2012.

The inspections found that the manufacturer had an acceptable quality management system and practices in place that ensured the consistent manufacture the above-mentioned products of good quality. For the most recent inspection, the manufacturer's responses to the nonconformities found at the time of the inspection were accepted 24 July 2015.

Product performance evaluation

Bioline Malaria Ag P.f was submitted to Round 1 and Round 5 of WHO Product Testing Programme in 2009 and 2014⁴. The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild type parasite panel, and a mixed panel of *P.falciparum* and non- *P.falciparum* to determine specificity. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

In Round 1, the following results were observed: *P. falciparum* panel detection score (97.5% at 200 parasites/μl), false-positive rates (2.4% for clean negatives, 0.0% for *P. vivax* at 200 parasites/μl, 0.0 % for *P. vivax* at 2000 or 5000 parasites/μl) and invalid rate (0.0%).

In Round 5, the following results were observed: *P. falciparum* panel detection score (95.0% at 200 parasites/μl), false-positive rates (0.0% for clean negatives, 0.0% for *P. vivax* at 200 parasites/μl, 2.9% for *P. vivax* at 2000 or 5000 parasites/μl) and invalid rate (0%), for Bioline Malaria Ag P.f meets the current laboratory evaluation requirements for prequalification.

Summary performance characteristics	Panel detection score	False positive rate (%)		Invalid rate (%)
	200 parasites/μl	200 parasites/μl	Clean negatives	
	Pf	Pv		
Round 1	97.47	0.0	2.38	0.0
Round 5	95.0	0.0	0.0	0.0

Based on these results, the performance evaluation for Bioline Malaria Ag P.f meets the WHO prequalification requirements.

⁴ The same product was also submitted to Round 1 Product Testing.


Labelling

- 1. Labels**
- 2. Instructions for use**


**Note: Labelling has been changed as per site name change and product rebranding. But temporarily labelling of legacy brand (SD BIOLINE) will be used in the market according to registration status in each country.*

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
CONTENTS

EN

- 25 Test strips with desiccant in individual foil pouch
- 250 µgg desiccant
- 25 Diagnostic response pad/membrane (7 µl)
- 25 Safety lancet, 25 µl blood volume
- 25 Control and reference for use
- 1 Instruction for use

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

Malaria Ag P.f. is a registered trademark of Abbott Diagnostic Care Inc. All other trademarks are the property of their respective owners.


1.4 Package box for 05FK53


1.5 Device pouch for 05FK50, 05FK51, 05FK52, 05FK53

1.6 Outer pouch for _05FK52

1.7 Outer pouch for _05FK53

2. Instructions for use⁵

⁵ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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ABOUT THE TEST

[Introduction]

Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is caused by a parasite that is transmitted from one human to another through the bite of infected Anopheles mosquitoes. There are four types of malaria that can infect humans: *Plasmodium falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and are released in another form: merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year. At present, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

[Test principle]

The Bioline[™] Malaria Ag P.f test device contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to HRP2 of *P. falciparum* on test line “P.f” region. The mouse monoclonal antibodies specific to HRP2 of *P.f* -colloid gold conjugate reacts with the Malaria *Plasmodium falciparum* antigen in the specimen. They move along the membrane chromatographically to the test region “P.f” and form a visible line as the antibody-antigen-antibody gold particle complex with high degree of sensitivity and specificity.

25 Test devices with desiccant in individual foil pouch

- Assay diluent (1 x 5 mL/vial),
- 25 Disposable specimen applicators (5 µl), 25 Sterile Lancets, 25 Alcohol swabs
- Instructions for use
- Active ingredients of main component

1 Test strip includes: Gold conjugate : Mouse monoclonal antibodies specific to P.f HRP2 - gold colloid (0.15:0.02 µg), Test line P.f: Mouse monoclonal antibodies specific to P.f HRP2 (0.5:0.1 µg), Control line : Goat anti-mouse IgG (0.1:0.2 µg)

- Assay diluent includes: Bovine serum albumin (q.s.), 200 mM Phosphate Saline Buffer (q.s.), Triton X-100 (q.s.), Sodium azide (q.s.)

[Intended Use]

The Bioline[™] Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of Malaria *Plasmodium falciparum* in human whole blood. Bioline[™] Malaria Ag P.f test is intended for professional use, only for an initial screening test and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear.

Materials provided and active ingredients of main component

- The Bioline[™] Malaria Ag P.f test kit contains the following items to perform the assay:
 - 25 Test devices with desiccant in individual foil pouch
 - Assay diluent (1 x 5 mL/vial),
 - 25 Disposable specimen applicators (5 µl), 25 Sterile Lancets, 25 Alcohol swabs
 - Instructions for use
 - Active ingredients of main component
- 1 test strip includes:
 - Gold conjugate : Mouse monoclonal antibodies specific to P.f HRP2 - gold colloid (0.15:0.02 µg), Test line P.f: Mouse monoclonal antibodies specific to P.f HRP2 (0.5:0.1 µg), Control line : Goat anti-mouse IgG (0.1:0.2 µg)
 - Assay diluent includes: Bovine serum albumin (q.s.), 200 mM Phosphate Saline Buffer (q.s.), Triton X-100 (q.s.), Sodium azide (q.s.)

Materials required but not provided

- Protective gloves, Timer, Biohazard container

Kit storage and stability

- The test kit should be stored at a temperature between 1°C and 40 °C. Do not freeze the kit or its components.
- The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
- Do not use the test kit beyond its expiration date.
- The shelf life of the kit is as indicated on the outer package.
- Do not use the test kit if the package is damaged or the seal is broken.
- When transporting or storing the kit, avoid exposure to high temperature (above 45°C) for a period longer than 1 week.

Warnings

- For in vitro diagnostic use only. Do not reuse the test device and its components.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
- Do not splash or aerosol formation of specimen and assay diluent.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. lancet, specimen applicator, test device) in a biohazard container as if they were infectious waste.
- Do not mix or interchange different specimens.
- Do not mix reagent of different lots or those for other products.
- Do not pipette by mouth.
- Do not eat the desiccant in the foil pouch.
- Do not drink assay diluent.
- The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

Specimen collection and storage

[Collection by venipuncture]

- Use an alcohol swab to collect whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
- If the blood specimen is not immediately tested, it should be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the blood specimen should be tested within 3 days.
- Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.

[Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.

Performance characteristics

- Internal evaluation of Bioline[™] Malaria Ag P.f test
 - A. Bioline[™] Malaria Ag P.f test for the detection of *P. falciparum* malaria according to the level of parasita (Table 1)
 - B. Sensitivity and specificity of Bioline[™] Malaria Ag P.f test. (Table 2)
- Reproducibility of Bioline[™] Malaria Ag P.f has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

Table-1					
Nb. of parasites/µL of blood	Microscopy (No. of positive)	Bioline [™] Malaria Ag P.f (No. of positive)	Sensitivity (%)		
1-50	16	15	93.8%		
51-100	35	35	100.0%		
101-500	96	96	100.0%		
501-1,000	81	81	100.0%		
1,001-5,000	79	79	100.0%		
>5,000	61	61	100.0%		
Total	368	367	99.7%		

Table-2					
Reference method (Microscopic examination)	Result of Bioline [™] Malaria Ag P.f			Sensitivity (%)	
	Positive	Negative	Total		
P.f Positive	367	1	368		
Negative	1	199	200		
Sensitivity (95% CI)				99.7% (98.5 - 100%)	
Specificity (95% CI)				99.5% (97.2 - 99.9%)	

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


- 6 EN** Interpret test results 15 minutes (up to 30 minutes) after adding assay diluent. Reading outside of this time frame (before 15 min or after 30 min) may provide false results.

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


ENGLISH

About the test
[Introduction]
Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is caused by a parasite that is transmitted from one human to another through the bite of infected Anopheles mosquitoes. There are four types of malaria that can infect humans: *Plasmodium falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and are resealed in another form: merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 200 million clinical cases and 0.5 million malaria-caused deaths per year. At present, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

[Test principle]
The Bioline™ Malaria Ag P.F test device contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to HRP2 of *P. falciparum* in the test line region (P.F). The mouse monoclonal antibodies specific to HRP2 of P.F -colloid gold conjugate reacts with the Malaria *Plasmodium falciparum* antigen in the specimen. They move along the membrane chromatographically to the test region "P.F" and form a visible line as the antibody-antigen-antibody gold particle complex with high degree of sensitivity and specificity.

[Precaution]
Anticoagulants including heparin, EDTA and sodium citrate do not affect the test result. Use of other anticoagulants have not been validated. Their use may affect the test result.
Use a new separate disposable capillary pipettes or pipette tips for each specimen in order to avoid cross-contamination of specimens, which could produce erroneous results.
Discard the lancet or alcohol swab if package is pierced or damaged. The use of damaged lancet may cause any infection at the punctured skin area to cease to existing its sterility.
Repeated frozen-thawed cycle for specimen should be avoided.

Test procedure (Refer to figure)
1. Bring all kit components and specimen to room temperature prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. Clean the area to be lanced with an alcohol swab.
4. Twist off or remove the protective cap from the safety lancet to break the seal and discard.
5. Prick the lateral side of the finger with the safety lancet provided. Then, safely dispose of the lancet immediately after use.
6. Using a capillary pipette (5 µl) provided, draw whole blood specimen to black line and then transfer drawn whole blood into the round specimen well. Or, Using a disposable inverted cup (5 µl) provided, pier the circular end of inverted cup into the blood specimen and carefully place the circular end of the inverted cup into the round specimen well.
7. Dispense 4 drops of assay diluent vertically into the square assay diluent well.
8. Interpret test results 15 minutes (up to 30 minutes) after adding assay diluent. Reading outside of this time frame (Before 15 min or after 30 min) may provide false results.

[Intended Use]
The Bioline™ Malaria Ag P.F test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of Malaria *Plasmodium falciparum* in human whole blood. Bioline™ Malaria Ag P.F test is intended for only professional use, only for an initial screening test and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear.

- Materials provided and active ingredients of main component**
- The Bioline™ Malaria Ag P.F test kit contains the following items to perform the assay:
 - 25 Test devices with desiccant in individual foil pouches
 - Assay diluent (1 x 5 mL/vial)
 - 25 Disposable specimen applicators (5 µl), 25 Safety Lancets, 25 Alcohol swabs (option)
 - 1 Instructions for use
 - Active ingredients of main component
 - 1 test strip includes: Gold conjugate; Mouse monoclonal antibodies specific to P.F/HRP2 - gold colloid (0.10±0.02 µg); Test line P.F; Mouse monoclonal antibodies specific to P.F/HRP2 (0.5±0.01 µg); Control line: Goat anti-mouse IgG (1.0±0.2 µg)
 - Assay diluent includes: Bovine serum albumin (g/L), 200 mM Phosphate Saline Buffer (g/L), TRITON™ X-100 (g/L), Sodium azide (g/L)

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- Micropipette, Protective gloves, Tieser, Biohazard container

- Kit storage and stability**
- The test kit should be stored at a temperature between 1°C and 40 °C. Do not freeze the kit or its components. Assay diluent can should be kept firmly sealed between each use.
 - When stored at room temperature, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to the test. Do not open the pouch whilst components come to room temperature.
 - Assay diluent may be opened and resealed for each assay. It is stable until expiration date if kept at 1 - 40 °C.
 - The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from foil pouch.
 - Do not use the test kit beyond its expiration date.
 - The shelf life of the kit is as indicated on the outer package.
 - Do not use the test kit if the pouch is damaged or the seal is broken.
 - When transporting, storing the kit, avoid exposure to high temperature (above 45°C) for a period longer than 1 week.

- Warnings**
- For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
 - The instructions must be followed exactly to achieve accurate results. Any individual performing the test using this product must be trained in its use and must be experienced in laboratory procedures.
 - Do not eat or smoke while handling specimens.
 - Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
 - Dispose of aerosol formation of specimen and assay diluent.
 - Clean up spills thoroughly using an appropriate disinfectant.
 - Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. lancet, specimen applicator, test device) in a biohazard container as if they were infectious waste.
 - Do not mix or interchange different specimens.
 - Do not mix reagent of different lots or those for other products.
 - Do not pipette by mouth.
 - Do not eat the desiccant in the foil pouch.
 - Do not drink assay diluent.
 - The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

- Specimen collection and storage**
[Collection by venipuncture]
- Using antiseptic, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
 - If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
 - If stored at 2 - 8 °C, the blood specimen must be tested within 3 days.
 - Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.
- Performance characteristics**
- Internal evaluation of Bioline™ Malaria Ag P.F test
 - Bioline™ Malaria Ag P.F test for the detection of P. falciparum malaria according to the level of parasite. (Table-1)
 - Sensitivity and specificity of Bioline™ Malaria Ag P.F test. (Table-2)
 - Reproducibility of Bioline™ Malaria Ag P.F has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.
- To evaluate the interference of Bioline™ Malaria Ag P.F test kit with known relevant interfering specimens, the haemolytic specimens, rheumatoid factors-containing specimens and lipaemic, icteric specimens were investigated. In these studies, those specimens did not interfere with the Bioline™ Malaria Ag P.F test.
 - Assay for the detection of Malaria according to the level of parasite.
 - Plasmodium falciparum*: Sensitivity with > 50 parasit/µl of blood is 100%.

Product Disclaimer:
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.

No. of parasites / µL of blood	Microscopy (No. of positive)	Bioline™ Malaria Ag P.F (No. of positive)	Sensitivity (%)
1-50	16	15	93.8 %
51-100	35	35	100.0 %
101-500	96	96	100.0 %
501-1,000	81	81	100.0 %
1,001-5,000	79	79	100.0 %
>5,000	61	61	100.0 %
Total	368	367	99.7%

Reference method (Microscopic examination)	Result of Bioline™ Malaria Ag P.F			
	Positive	Negative	Total	
P.F Positive	367	1	368	
Negative	1	199	200	
Sensitivity (95 % CI)	99.7% (98.5 - 100 %)			
Specificity (95 % CI)	99.5% (97.2 - 99.9 %)			

Método de referência (exames microscópicos)	Resultado de Bioline™ Malaria Ag P.F			
	Positivo	Negativo	Total	
P.F Positivo	367	1	368	
Negativo	1	199	200	
Sensibilidade (IC de 95 %)	99,7% (98,5 - 100,0 %)			
Especificidade (IC de 95 %)	99,5% (97,2 - 99,9 %)			

Bibliography of suggested reading / Bibliographie - Lectures conseillées / Bibliografía de lectura recomendada / Referencias bibliográficas

1. SM Gekunda, JY Carter and Joseph Macharia : FIELD EVALUATION OF Bioline™ -MALARIA ANTIGEN RAPID TEST (2006)

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4. REPORT ON THE EVALUATION OF SD-BIOLINE MALARIA ANTIGEN AT Mt. DARWIN DISTRICT HOSPITAL, ZIMBABWE (2007)

5. http://www.malariaelimination.gov.au/sites/default/files/fileuploads/AuGuideonMalariaEliminationforPolicyMakers.pdf

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ENGLISH

About the test
[Introduction]
Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is caused by a parasite that is transmitted from one human to another through the bite of infected *Anopheles* mosquitoes. There are four types of malaria that can infect humans: *Plasmodium falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and are released in another form: merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year. At present, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

[Text principal]
The Bioline™ Malaria Ag P.f test device contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to HRP2 of *P. falciparum* on test line “P.F” region. The mouse monoclonal antibodies specific to HRP2 of P.f. -colloidal gold conjugate reacts with the Malaria *Plasmodium falciparum* antigen in the specimen. They move along the membrane chromatographically to the test region “P.F” and form a visible line as the antibody-antigen-antibody gold particle complex with high degree of sensitivity and specificity.

[Intended Use]
The Bioline™ Malaria Ag P.f test is a rapid, qualitative test for the detection of histidine-rich protein 2 (HRP2) antigen of *Malaria Plasmodium falciparum* in human whole blood. Bioline™ Malaria Ag P.f test is intended for professional use, only for an initial screening test and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear.

Materials provided and active ingredients of main component
The Bioline™ Malaria Ag P.f test kit contains the following items to perform the assay:
• 25 Test devices with desiccant in individual foil pouch
• 25 Assay diluents
• 25 Disposable specimen applicators (5 µl), 25 Sterile lancets, 25 Alcohol swabs
• 25 Summarized instructions for use
• 1 Instructions for use
3. Active ingredients of main component
• 1 test strip includes : Gold conjugate : Mouse monoclonal antibodies specific to PfHRP2 - gold colloid (0.1:0.02 µg), Test Line P.f. Mouse monoclonal antibodies specific to Pf HRP2 (0.5±0.1 µg), Control Line anti- mouse IgG (0.0:0.2 µg), Control Line
• Assay diluent includes: Bovine serum albumin (g/L), 200 mM Phosphate Saline Buffer (g/L), Triton X-100 (g/L), Sodium azide, (g/L.)

Materials required but not provided
• Micropipette, Protective gloves, Timer, Biohazard container

Kit storage and stability
1. The test kit should be stored at a temperature between 1°C and 40 °C. Do not freeze the kit or its components.
Note: When stored at room temperature, all kit components must be brought to room temperature (15 - 40°C) minimum 30 mins prior to the test.
2. Do not open the pouch whilst components come to room temperature.
3. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is as indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. When transporting or storing the kit, avoid exposure to high temperature (above 45°C) for a period longer than 1 week.

Warnings
1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in the use of the device. Neither the test nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.
Warnings
1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in the use of the device. Neither the test nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Kit storage and stability
1. The test kit should be stored at a temperature between 1°C and 40 °C. Do not freeze the kit or its components.
Note: When stored at room temperature, all kit components must be brought to room temperature (15 - 40°C) minimum 30 mins prior to the test.
2. Do not open the pouch whilst components come to room temperature.
3. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is as indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. When transporting or storing the kit, avoid exposure to high temperature (above 45°C) for a period longer than 1 week.

Warnings
1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in the use of the device. Neither the test nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.
Warnings
1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in the use of the device. Neither the test nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Kit storage and stability
1. The test kit should be stored at a temperature between 1°C and 40 °C. Do not freeze the kit or its components.
Note: When stored at room temperature, all kit components must be brought to room temperature (15 - 40°C) minimum 30 mins prior to the test.
2. Do not open the pouch whilst components come to room temperature.
3. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is as indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. When transporting or storing the kit, avoid exposure to high temperature (above 45°C) for a period longer than 1 week.

Warnings
1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in the use of the device. Neither the test nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.
Warnings
1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in the use of the device. Neither the test nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Limitations and Interferences
1. The test procedure, precautions and interpretation of results for this test must be followed precisely when testing.
2. This test kit detects *Plasmodium falciparum* histidine-rich protein 2 (HRP2) in patient whole blood and is useful as a screening procedure of malaria diagnosis.
3. The test is limited to the detection of antigen to *Malaria Plasmodium falciparum*. Although the test is very accurate in detecting HRP2 antigen of *P. falciparum*, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Interanal quality control
The Bioline™ Malaria Ag P.f test has test line (“P.F”) and control line (“C”) on the surface of the device. Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.
Interanal quality control
The Bioline™ Malaria Ag P.f test has test line (“P.F”) and control line (“C”) on the surface of the device. Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Expected values
The Bioline™ Malaria Ag P.f test has been compared with microscopic examination. The overall accuracy is greater or equal to 99%.

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Table-1				
N° de parasites /µl de sang	Microscopio (No. of positive)	Bioline™ Malaria Ag P.f (No. of positive)	Sensibilité (%)	Spécificité (%)
1-50	16	15	93.8%	
51-100	35	35	100.0%	
101-500	96	96	100.0%	
501-1,000	81	81	100.0%	
1,001-5,000	79	79	100.0%	
>5,000	61	61	100.0%	
Total	368	367	99.7%	
reference panels. All values were identical to reference panel acceptance criteria.				
3. To evaluate the interference of Bioline™ Malaria Ag P.f test kit with known reagents interfering with specimens, the haemolytic specimens, rheumatoid factors-contained specimens and lipemic, icteric specimens were investigated. In these studies, those specimens did not interfere with the Bioline™ Malaria Ag P.f test.				
4. Assay for the detection of Malaria according to the level of parasite.				
• <i>Plasmodium falciparum</i> : Sensitivity with > 50 parasite/µl of blood is 100%.				

Table-1				
N° of parasites /µl of blood	Microscopy (No. of positive)	Bioline™ Malaria Ag P.f (No. of positive)	Sensitivity (%)	Specificity (%)
1-50	16	15	93.8%	
51-100	35	35	100.0%	
101-500	96	96	100.0%	
501-1,000	81	81	100.0%	
1,001-5,000	79	79	100.0%	
>5,000	61	61	100.0%	
Total	368	367	99.7%	
reference panels. All values were identical to reference panel acceptance criteria.				
3. To evaluate the interference of Bioline™ Malaria Ag P.f test kit with known reagents interfering with specimens, the haemolytic specimens, rheumatoid factors-contained specimens and lipemic, icteric specimens were investigated. In these studies, those specimens did not interfere with the Bioline™ Malaria Ag P.f test.				
4. Assay for the detection of Malaria according to the level of parasite.				
• <i>Plasmodium falciparum</i> : Sensitivity with > 50 parasite/µl of blood is 100%.				

Table-2				
Reference method (Microscopic examination)	Result of Bioline™ Malaria Ag P.f	Positive	Negative	Total
P.F Positive		367	1	368
Negative		1	199	200
Sensitivity (95% CI)		99.7% (98.5 - 100%)		
Specificity (95% CI)		99.5% (97.2 - 99.9%)		

2. To evaluate the interference of Bioline™ Malaria Ag P.f test kit with known reagents interfering with specimens, the haemolytic specimens, rheumatoid factors-contained specimens and lipemic, icteric specimens were investigated. In these studies, those specimens did not interfere with the Bioline™ Malaria Ag P.f test.

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