WHO Prequalification of In Vitro Diagnostics Programme PUBLIC REPORT

Product: *CareStart*[™] Malaria HRP2/pLDH (Pf/PAN) COMBO
Number: PQDx 0136-049-00

Abstract

CareStart[™] Malaria HRP2/pLDH (Pf/PAN) COMBO to be labelled as *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT from October 2015, with product codes RMRM-02571 and RMRM-05071 (old product code G0131),¹ manufactured by Access Bio, Inc., at Somerset NJ, USA, rest-of-world regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 28 May 2015.

This public report was amended on 28 June 2018 to reflect the addition of a new product configuration, updates of product codes and inclusion of labelling. The product now has two configurations: **multi kit** (product codes: **RMRM-02571, RMRM-02571CB, RMRM-02571RB, RMRM-02571RI, RMRM-05071, RMRM-05071CB, RMRM-05071RB and RMRM-05071RI**) and **single kit** (product codes: **RMRU-02571, RMRU-02571CB, RMRU-02571RB, RMRU-02571RI, RMRU-04071, RMRU-04071CB, RMRU-04071RB and RMRU-04071RI**). Each multi kit and single kit has 8 different packaging configurations. For further details, please refer to the test kit content table.

Summary of prequalification status for *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

	Date	Outcome
PQ public report amended	28 June 2018 Inclusion of a new configuration (addition of single-use buffer vials and an 'Instruction Card'), updates of product codes and inclusion of labelling.	listed
Status on PQ list	25 May 2015	listed
Dossier assessment	15 April 2015	MR
Inspection status	8-10 February 2017	MR
Laboratory evaluation	2014	MR

MR: Meets Requirements

CareStart[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT was accepted for the WHO list of prequalified in vitro diagnostics on the basis of data submitted and publicly available information.

¹ The old product codes RMRM-02571, RMRM-05071, and G0131 are obsolete.

For the rapid qualitative detection of malaria HRP2 (histidine-rich protein 2) of *Plasmodium falciparum* and pLDH (plasmodium lactate dehydrogenase) of *P.falciparum*, *P.vivax*, *P.ovale* and *P.malariae* in human whole blood as an aid in the diagnosis of malaria infection.

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the Plasmodium species (*P.falciparum*, *P.vivax*, *P.malariae*, and *P.ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P.falciparum*. The conjugate pad is dispensed with antibodies absorbed on gold particles, which are specific to pLDH of PAN and specific to HRP2 of *P. falciparum*.

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is designed for the differentiated diagnosis of *P.falciparum*, *P.vivax*, *P.malariae*, and *P.ovale* infection.

Other clinically available tests are required if the obtained results are questionable. A definitive clinical diagnosis should not be made based on the result of this test, but should only be made by a qualified physician after all clinical and laboratory findings have been evaluated.

Storage:

The test kit should be stored at 1 to 40 °C.

Shelf-life:

30 months.

Reading time:

20 minutes.

Materials needed but not provided:

1. Pair of disposable gloves
2. Timer
3. Sharp box
4. Pencil or pen
5. Sterile gauze or cotton

Consumables: Configuration	Packaging	Item(s)
Multi kit	25 tests/kit	1 × Assay buffer bottle (4 mL)
		25 × Specimen transfer device
		25 × Lancet
		25 × Alcohol swab
		1 × Instructions for Use
	50 tests/kit	1 × Assay buffer bottle (6mL)
		50 × Specimen transfer device
		50 × Lancet
		50 × Alcohol swab
		1 × Instructions for Use
Single kit	25 single kit bags/kit	Each single kit bag contains one of each of the following items:
		- Single use buffer vial
		- Specimen transfer device
		- Lancet
		- Alcohol swab
		- Instruction card
	40 single kit bags /kit	1 × Instructions for Use
		Each single kit bag contains one of each of the following items:
		- Single use buffer vial
		- Specimen transfer device
		- Lancet
		- Alcohol swab
		- Instruction card
		1 × Instructions for use

The test kit contains:

	Packaging	Product code	Item	Instructions for use	Instruction card	RDT Box	Single kit bag (Lab in a pack)	Cassette packaging
Multi kit	25 tests/kit	RMRM-02571	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Inverted cup) 	IFU-RMRM71-EFSP	n/a	LB-KMR2571CIA	n/a	PROL6140
		RMRM-02571CB	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Bulb) 	IFU-RMRM71CB-EFSP		LB-KMR2571CBA		
		RMRM-02571RB	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Bulb) 	IFU-RMRM71RB-EFSP		LB-KMR2571RBA		
		RMRM-02571RI	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Inverted cup) 	IFU-RMRM71RI-EFSP		LB-KMR2571RIA		
	50 tests/kit	RMRM-05071	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Inverted cup) 	IFU-RMRM71-EFSP		LB-KMR5071CIA		
		RMRM-05071CB	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Bulb) 	IFU-RMRM71CB-EFSP		LB-KMR5071CBA		
		RMRM-05071RB	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Bulb) 	IFU-RMRM71RB-EFSP		LB-KMR5071RBA		
		RMRM-05071RI	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Inverted cup) 	IFU-RMRM71RI-EFSP		LB-KMR5071RIA		
Single kit	25 tests/kit	RMRU-02571	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Inverted cup) 	IFU-RMRU71-EFSP	IFU-RMRU71-E	LB-KUR2571CIA	PABO140	
		RMRU-02571CB	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Bulb) 	IFU-RMRU71CB-EFSP	IFU-RMRU71CB-E	LB-KUR2571CBA		
		RMRU-02571RB	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Bulb) 	IFU-RMRU71RB-EFSP	IFU-RMRU71RB-E	LB-KUR2571RBA		
		RMRU-02571RI	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Inverted cup) 	IFU-RMRU71RI-EFSP	IFU-RMRU71RI-E	LB-KUR2571RIA		
	40 tests/kit	RMRU-04071	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Inverted cup) 	IFU-RMRU71-EFSP	IFU-RMRU71-E	LB-KUR4071CIA		
		RMRU-04071CB	<ul style="list-style-type: none"> Lancets (Capped) Specimen transfer devices (Bulb) 	IFU-RMRU71CB-EFSP	IFU-RMRU71CB-E	LB-KUR4071CBA		
		RMRU-04071RB	<ul style="list-style-type: none"> Lancets (Retractable) Specimen transfer devices (Bulb) 	IFU-RMRU71RB-EFSP	IFU-RMRU71RB-E	LB-KUR4071RBA		
		RMRU-04071RI	<ul style="list-style-type: none"> Lancets (Retractable) specimen transfer devices (Inverted cup) 	IFU-RMRU71RI-EFSP	IFU-RMRU71RI-E	LB-KUR4071RIA		

Background information

Access Bio, Inc. submitted an application for prequalification of *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT. Based on the results of the WHO product testing of malaria RDTs Round 5, *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT was given priority for prequalification.

Product dossier assessment

Access Bio, Inc. submitted a product dossier for *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT as per the “Instructions for compilation of a product dossier” (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2). Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT for prequalification.

Manufacturing site inspection

A comprehensive inspection was performed at the site of manufacture (65 Clyde Road Suite A, Somerset NJ, USA), of *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT in February 2017 as per the “Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics” (PQDx_014 v3). The inspection found that the manufacturer had an acceptable quality management system and good manufacturing practices in place that ensured the consistent manufacture of a product of good quality. The manufacturer's responses to the nonconformities found at the time of the inspection were accepted 23 April 2018.

Laboratory evaluation


The fifth round of WHO product testing of RDTs for malaria antigen detection was completed in 2013. The product was evaluated against a *P.falciparum* cultured line panel, *P.falciparum* wild type parasite panel, *P.vivax* wild type parasite panel and a *P.falciparum* negative panel. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

Based on the demonstrated *P. falciparum* panel detection score (90.0% at 200 parasites/μl) and *P. vivax* panel detection score (94.3% at 200 parasites/μl), false-positive rate (0.4% for clean negatives, 1.5% for *P. falciparum* at 200 parasites/μl, 0.7% for *P. vivax* at 200 parasites/μl) and invalid rate (0%), *CareStart*[™] Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT meets the current laboratory evaluation requirements for prequalification.


Labelling

- 1. Labels**
- 2. Instructions for use**
- 3. Instruction Cards**
- 4. RDT box designs**

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LIMITATIONS AND INTERFERENCES:

- The following anticoagulants have been validated for use with this test: heparin, EDTA and citrate.
- This test is designed to detect HRP2 and pLDH antigens of Malaria *Plasmodium* species. Other clinically available tests are required if the obtained results are questionable. A definitive clinical diagnosis should not be made solely based on the result of this test, but should only be made by a qualified physician after all clinical and laboratory findings have been evaluated.
- A faint test line should be read as positive. A false-negative is possible due to a low parasite density.
- This test may still produce a positive result after successful anti-malarial treatment. Therefore, its use is not recommended for monitoring a response to anti-malarial treatment.³⁾
- If specimens cannot be tested immediately, they should be refrigerated at 2 - 8°C for up to 3 days.
- The test may produce a false positive result for a patient with acute schistosomiasis, a high level of rheumatoid factor or presence of human anti-mouse IgG antibody.⁶⁾
- No cross reactivity has been observed when the samples from a patient of chagas, dengue and leishmaniasis.
- Results are stable for 30 minutes after assay time of 20 minutes.
- The prozone effect may cause a false-negative result.

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CareStart™ Malaria Pf /PAN (HRP2/pLDH) Ag Combo RDT

ENGLISH

A RAPID DIAGNOSTIC TEST FOR THE DETECTION OF MALARIA HRP2 (Pf specific) AND pLDH (PAN specific) IN HUMAN BLOOD.

INTENDED USE:

For the rapid qualitative detection of malaria HRP2 (histidine-rich protein 2) of *Plasmodium falciparum* and pLDH (plasmidium lactate dehydrogenase) of *P. falciparum*, *Plasmodium vivax*, *Plasmodium ovale* and *Plasmodium malariae* in human whole blood as an aid in the diagnosis of malaria infection.

SUMMARY:

Malaria is a serious parasitic disease characterized by fever, chills and anemia. It is caused by a bite of infected Anopheles mosquitoes resulting in transmission of protozoan parasites to human. There are four major types of human malaria: *P. falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. When infected, the parasites (called sporozoites) migrated to the liver where they mature and then are released to the blood stream of human infecting the red blood cells. Malaria infection occurs in more than 90 countries worldwide but is mostly prevalent in sub-Saharan Africa. It is estimated that there are over 300 million clinical cases and about 1 million malaria-caused deaths per year.

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the *Plasmodium* species (*P. falciparum*, *P.vivax*, *P. malariae* and/or *P. ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P. falciparum*. The conjugate pad is dispensed with antibodies adsorbed on gold particles.^{(1), (2), (4), (5)}

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is designed for the differentiated diagnosis of *P. falciparum*, *P.vivax*, *P. malariae* and *P. ovale* infection.

WARNINGS AND PRECAUTIONS:

- Kits are for *in vitro* diagnostic and professional use only.
- Read the provided instructions for use before using the kit and follow the provided information when using the kit.
- Wash the hands thoroughly before and after using the kit.
- Wear protective gloves at all time while using the kits and dispose the gloves immediately after each testing.
- Change gloves and wash hands when contacted with potentially infectious materials.
- Do routine clean-up using an appropriate disinfectant.
- Do not eat or smoke while using the kit.
- Use the cassette and accessories (lancet and alcohol swab) immediately after opening the package.
- All provided materials are single-use. Do not re-use any of the materials.
- Do not use the materials from different lots.
- Do not use a cassette, if the packaging of cassette is compromised or the desiccant colour indicator is blue or black.
- Do not use the kit if the expiration date is past.
- The lancet is sterile. If a lancet cap is loose or damaged, do not use the lancet.
- Do not swallow the buffer solution.
- Do not use the alcohol swab if the package of the alcohol swab is damaged.
- Observe the storage condition (1~40°C) indicated on the packaging and RDT box.
- Do not freeze the kit but store the kit in a refrigerator. The refrigerated cassettes need to be warmed to the room temperature prior to use.
- Dispose wastes in accordance with the local regulations.
- Buffer bottle should be tightly closed after each use and stored in a cool area avoiding direct sun light. Discard the buffer solution and the bottle after the expiration date indicated on the bottle.

PERFORMANCE CHARACTERISTICS:

The *CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT has been tested with positive and negative clinical samples confirmed by microscopic examination. The results are shown below:

Specimen	Positive	Negative	Sensitivity
<i>P. falciparum</i>	98	2	98%
<i>P. vivax</i>	96	4	96%

Specimen	Positive	Negative	Specificity
Negative	5	195	97.5%

95% CI, (97.05, 100) for *P. falciparum* positive result, (92.16, 99.84) for *P. vivax* positive result and (94.64, 99.36) for negative result (laboratory study)

PRECISION:

Precision was evaluated by testing 10 replicates of *CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT with three specimens: a negative, a low positive and a strong positive specimen. The test results showed 100% confirmation of the expected results.

DESCRIPTION OF SYMBOL USED / DESCRIPTION DES SYMBOLES UTILISÉS / DESCRIPCIÓN DE LOS SÍMBOLOS UTILIZADOS / DESCRIÇÃO DO S SÍMBOLO USADO

	Single use / Usage unique / Uso único / Uso único		Use by / À utiliser avant / Utilice antes de / Use até
	Keep dry / Garder au sec / Mantêngase seco / Mantenha seco		Manufacturer / Fabricant / Fabricante / Fabricante
	Batch code / Code du lot / Código de lote / Código do lote		Catalogue number / Numéro de catalogue / Número de catálogo / Número de catálogo

CONTENTS OF PRODUCT / CONTENUS DU PRODUIT / CONTENIDOS DEL PRODUCTO / CONTEÚDOS DO PRODUTO

EN *CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains the following items:

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the *Plasmodium* species (*P. falciparum*, *P.vivax*, *P. malariae* and/or *P. ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P. falciparum*. The conjugate pad is dispensed with antibodies adsorbed on gold particles.^{1), 2), 4), 5)}

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WARNINGS AND PRECAUTIONS:

- Kits are for *in vitro* diagnostic and professional use only.
- Read the provided instructions for use before using the kit and follow the provided information when using the kit.
- Wash the hands thoroughly before and after using the kit.
- Wear protective gloves at all time while using the kits and dispose the gloves immediately after each testing.
- Change gloves and wash hands when contacted with potentially infectious materials.
- Do routine clean-up using an appropriate disinfectant.
- Do not eat or smoke while using the kit.
- Use the cassette and accessories (lancet and alcohol swab) immediately after opening the package.
- All provided materials are single-use. Do not re-use any of the materials.
- Do not use the materials from different lots.
- Do not use a cassette, if the packaging of cassette is compromised or the desiccant colour inicator is blue or black.
- Do not use the kit if the expiration date is past.
- The lancet is sterile. If a lancet cap is loose or damaged, do not use the lancet.
- Do not swallow the buffer solution.
- Do not use the alcohol swab if the package of the alcohol swab is damaged.
- Observe the storage condition (1~40°C) indicated on the packaging and RDT box.
- Do not freeze the kit but store the kit in a refrigerator. The refrigerated cassettes need to be warmed to the room temperature prior to use.
- Dispose wastes in accordance with the local regulations.
- Buffer bottle should be tightly closed after each use and stored in a cool area avoiding direct sun light. Discard the buffer solution and the bottle after the expiration date indicated on the bottle.

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CareStart™ Malaria Pf /PAN (HRP2/pLDH) Ag Combo RDT

ENGLISH

A RAPID DIAGNOSTIC TEST FOR THE DETECTION OF MALARIA HRP2 (Pf specific) AND pLDH (PAN specific) IN HUMAN BLOOD

INTENDED USE:

For the rapid qualitative detection of malaria HRP2 (histidine-rich protein 2) of *Plasmodium falciparum* and pLDH (plasmodium lactate dehydrogenase) of *P. falciparum*, *Plasmodium vivax*, *Plasmodium ovale* and *Plasmodium malariae* in human whole blood as an aid in the diagnosis of malaria infection.

SUMMARY:

Malaria is a serious parasitic disease characterized by fever, chills and anemia. It is caused by a bite of infected Anopheles mosquitoes resulting in transmission of protozoan parasites to human. There are four major types of human malaria: *P. falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. When infected, the parasites (called sporozoites) migrated to the liver where they mature and then are released to the blood stream of human infecting the red blood cells. Malaria infection occurs in more than 90 countries worldwide, but is mostly prevalent in sub-Saharan Africa. It is estimated that there are over 300 million clinical cases and about 1 million malaria-caused deaths per year.

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the Plasmodium species (*P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P. falciparum*. The conjugate pad is dispensed with antibodies adsorbed on gold particles.^{10, 21, 41, 51}

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is designed for the differentiated diagnosis of *P. falciparum*, *P. vivax*, *P. malariae* and *P. ovale* infection

WARNINGS AND PRECAUTIONS:

- Kits are for *in vitro* diagnostic and professional use only.
- Read the provided instructions for use before using the kit and follow the provided information when using the kit.
- Wash the hands thoroughly before and after using the kit.
- Wear protective gloves at all time while using the kits and dispose the gloves immediately after each testing.
- Change gloves and wash hands when contacted with potentially infectious materials.
- Do routine clean-up using an appropriate disinfectant.
- Do not eat or smoke while using the kit.
- Use all components immediately after opening the package.
- All provided materials are single-use. Do not re-use any of the materials.
- Do not use the materials from different lots.
- Do not use a cassette, if the packaging of cassette is compromised or the desiccant colour inicator is blue or black.
- Do not use the kit if the expiration date is past.
- The lancet is sterile. If a lancet cap is loose or damaged, do not use the lancet.
- Do not swallow the buffer solution.
- Do not use the alcohol swab if the package of the alcohol swab is damaged.
- Observe the storage condition (1~40°C) indicated on the packaging and RDT box.
- Do not freeze the kit but store the kit in a refrigerator. The refrigerated cassettes need to be warmed to the room temperature prior to use.
- Dispose wastes in accordance with the local regulations.

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the Plasmodium species (*P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P. falciparum*. The conjugate pad is dispensed with antibodies adsorbed on gold particles.^{(1), (2), (4), (5)}

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REFERENCES / RÉFÉRENCES / REFERENCIAS / REFERÊNCIAS

- Valecha N., Eapen A., Usha Devi C., Ravindran J., Aggarwal A., and Subbarao SK. (2002). Field Evaluation of the ICT Malaria P.f./P.v. Immunochromatographic Test in India. *Annals of Tropical Medicine & Parasitology*. 96, 333-336.
- Iqbal J., Hira PR., Sher A., and AL-Enezi AA. (2001). Diagnosis of Imported Malaria by Plasmodium Lactate Dehydrogenase (pLDH) and Histidine-Rich Protein 2 (PfHRP-2)-based Immunocapture Assays. *American Journal of Tropical Medicine and Hygiene*. 64, 20-23.
- Tjitra E., Suprianto S., Dyer ME., Currie BJ. and Anstey NM. (2001). Detection of Histidine-Rich Protein 2 and Panmalarial ICT Malaria P.f./P.v. Test Antigens After Chloroquine Treatment of Uncomplicated falciparum Malaria Does Not Reliably Predict Outcome in Eastern Indonesia. *American Journal of Tropical Medicine and Hygiene*. 65, 593-598.
- Panton LJ., McPhie P., Malloy WL., Welles TE., Taylor DW. and Howard RJ. (1989). Purification and Partial Characterization of an Unusual Protein of Plasmodium falciparum: Histidine-Rich Protein II. *Molecular and Biochemical Parasitology*. 35,149-160.
- Leonard K. Basco, Frederique Marquet, Michael M. Makler, and Jacques Le Bras. (1995). Plasmodium falciparum and Plasmodium vivax: Lactate Dehydrogenase Activity and its Application for in vitro Drug Susceptibility Assay. *Experimental Parasitology*. 80, 260-271.
- Malaria Rapid Diagnostic Test Performance. Results of WHO product testing of malaria RDTs: round 7 (2015-2016). Table A4.7: Phase 2 false-positive rate for *P. falciparum* in samples containing specific non-malarial infectious pathogens. 124.

65 Clyde Road, Suite A Somerset, NJ 08873 U.S.A. Tel: 1-732-873-4040 Fax: 1-732-873-4043 E-mail: info@accessbio.net Web: www.accessbio.net

Instruction Cards

CareStart™ Malaria Pf/PAN(HRP2/pLDH) Ag Combo RDT

REF RMRU-02561CB / RMRU-04061CB / RMRU-02571CB / RMRU-04071CB


TEST PROCEDURE

1. Put on a new pair of gloves.
2. Write the patient's name on the cassette.
3. Clean the area to be pierced using an alcohol swab. Let the alcohol dry completely before proceeding to the next step.

4.  Squeeze the end of a fingertip and pierce the cleaned area of the fingertip using a lancet provided. Discard the lancet in the sharps box.

5. Wipe out the first drop of blood with sterile gauze or cotton.

6.  Collect the blood sample (5µl) using a provided specimen transfer device or a micropipette.

7.  Add 5µl of whole blood into "S" well.

8.  Add 3 drops (60µl) of buffer solution into "A" Well. Start a timer.

9. Read result at 20 minutes.

INTERPRETATION OF THE TEST RESULT

Negative

The presence of a line next to "C" indicates a negative result.

Pf Positive

The presence of two lines (one line in the result window next to "C" and another in the result window next to "1") indicates a positive result for *P. falciparum*.

PAN Positive

The presence of two lines (one line in the result window next to "C" and another line in the result window next to "2") indicates a positive result for *P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*.

Pf Positive or Mixed infection

The presence of three lines (3 lines in the result window next to "C", "2" and "1") indicates a positive result for *P. falciparum* or mixed infection of *P. falciparum* and other malaria (*P. vivax*, *P. ovale* or *P. malariae*).

Invalid

The test is invalid when a line does not appear next to "C". If this occurs, the test should be repeated using a new cassette.

For detailed information about the product and procedure, refer to the full version of instructions for use (Doc. No. IFU-RMRU71CB-EF-SP) included in the RDT box.

TEST PROCEDURE

- Put on a new pair of gloves.
- Write the patient's name on the cassette.
- Clean the area to be pierced using an alcohol swab. Let the alcohol dry completely before proceeding to the next step.

- 

Squeeze the end of a fingertip and pierce the cleaned area of the fingertip using a lancet provided. Discard the lancet in the sharps box.

- Wipe out the first drop of blood with sterile gauze or cotton.

- 

Collect the blood sample (5µl) using a provided specimen transfer device or a micropipette.

- 

Add 5µl of whole blood into "S" well.

- 

Add 3 drops (60µl) of buffer solution into "A" Well. Start a timer.

- Read result at 20 minutes.

INTERPRETATION OF THE TEST RESULT

Negative

The presence of a line next to "C" indicates a negative result.

Pf Positive

The presence of two lines (one line in the result window next to "C" and another in the result window next to "1") indicates a positive result for *P. falciparum*.

PAN Positive

The presence of two lines (one line in the result window next to "C" and another line in the result window next to "2") indicates a positive result for *P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*.

Pf Positive or Mixed infection

The presence of three lines (3 lines in the result window next to "C", "2" and "1") indicates a positive result for *P. falciparum* or mixed infection of *P. falciparum* and other malaria (*P. vivax*, *P. ovale* or *P. malariae*).

Invalid

The test is invalid when a line does not appear next to "C". If this occurs, the test should be repeated using a new cassette.

For detailed information about the product and procedure, refer to the full version of instructions for use (Doc. No. IFU-RMRU71-EFSP) included in the RDT box.

TEST PROCEDURE

[REF] RMRU-02561RB / RMRU-04061RB / RMRU-02571RB / RMRU-04071RB

1. Put on a new pair of gloves.
2. Write the patient's name on the cassette.
3. Clean the area to be pierced using an alcohol swab. Let the alcohol dry completely before proceeding to the next step.

4. Squeeze the end of a fingertip and pierce the cleaned area of the fingertip using a lancet provided. Discard the lancet in the sharps box.

5. Wipe out the first drop of blood with sterile gauze or cotton.

6. Collect the blood sample (5µl) using a provided specimen transfer device or a micropipette.

7. Add 5µl of whole blood into "S" well.

8. Add 3 drops (60µl) of buffer solution into "A" Well. Start a timer.

9. Read result at 20 minutes.

INTERPRETATION OF THE TEST RESULT

Negative

The presence of a line next to "C" indicates a negative result.

Pf Positive

The presence of two lines (one line in the result window next to "C" and another in the result window next to "1") indicates a positive result for *P. falciparum*.

PAN Positive

The presence of two lines (one line in the result window next to "C" and another line in the result window next to "2") indicates a positive result for *P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*.

Pf Positive or Mixed infection

The presence of three lines (3 lines in the result window next to "C", "2" and "1") indicates a positive result for *P. falciparum* or mixed infection of *P. falciparum* and other malaria (*P. vivax*, *P. ovale* or *P. malariae*).

Invalid

The test is invalid when a line does not appear next to "C". If this occurs, the test should be repeated using a new cassette.

For detailed information about the product and procedure, refer to the full version of instructions for use (Doc. No. IFU-RMRU71RB-EF-SP) included in the RDT box.

CareStart™ Malaria Pf/PAN(HRP2/pLDH) Ag Combo RDT

TEST PROCEDURE


1. Put on a new pair of gloves.
2. Write the patient's name on the cassette.
3. Clean the area to be pierced using an alcohol swab. Let the alcohol dry completely before proceeding to the next step.

4. Squeeze the end of a fingertip and pierce the cleaned area of the fingertip using a lancet provided. Discard the lancet in the sharps box.

5. Wipe out the first drop of blood with sterile gauze or cotton.

6. Collect the blood sample (5µl) using a provided specimen transfer device or a micropipette.

7. Add 5µl of whole blood into "S" well.

8. Add 3 drops (60µl) of buffer solution into "A" Well. Start a timer.

9. Read result at 20 minutes.

REF RMRU-02561RI / RMRU-04061RI / RMRU-02571RI / RMRU-04071RI

INTERPRETATION OF THE TEST RESULT

Negative

The presence of a line next to "C" indicates a negative result.

Pf Positive

The presence of two lines (one line in the result window next to "C" and another in the result window next to "1") indicates a positive result for *P. falciparum*.

PAN Positive

The presence of two lines (one line in the result window next to "C" and another line in the result window next to "2") indicates a positive result for *P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*.

Pf Positive or Mixed infection

The presence of three lines (3 lines in the result window next to "C", "2" and "1") indicates a positive result for *P. falciparum* or mixed infection of *P. falciparum* and other malaria (*P. vivax*, *P. ovale* or *P. malariae*).


Invalid

The test is invalid when a line does not appear next to "C". If this occurs, the test should be repeated using a new cassette.

For detailed information about the product and procedure, refer to the full version of instructions for use (Doc. No. IFU-RMRU71RI-EFSP) included in the RDT box.

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

REF RMRM-05071CB IVD ⓘ 50 ⚠️ 1°C 40°C

Contents

- 50 Cassettes, 50 Alcohol swabs,
- 50 Lancets (capped), 50 Specimen transfer devices (bulb type),
- 1 Buffer bottle and 1 Instructions for use

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents
 - 50 Cassettes, 50 Alcohol swabs,
 50 Lancets (capped), 50 Specimen transfer devices (inverted cup type),
 1 Buffer bottle and 1 Instructions for use

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Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents

- 50 Cassettes, 50 Alcohol swabs,
- 50 Lancets (retractable), 50 Specimen transfer devices (bulb type),
- 1 Buffer bottle and 1 Instructions for use

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents

- 50 Cassettes
- 50 Alcohol swabs
- 50 Lancets (retractable)
- 50 Specimen transfer devices (inverted cup type)
- 1 Buffer bottle
- 1 Instructions for use

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents

- 50 Cassettes, 50 Alcohol swabs,
- 50 Lancets (retractable), 50 Specimen transfer devices (inverted cup type),
- 1 Buffer bottle and 1 Instructions for use

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





REF RMRU-04071CB IVD 40°C 1°C

CareStart™
 Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents
 1 Instructions for use and 40 Single kit bags - Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (capped), 1 specimen transfer device (bulb type), 1 buffer vial and 1 instruction card

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Single Kit

Contents - 1 Instructions for use and 40 Single kit bags
 Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (capped), 1 specimen transfer device (bulb type), 1 buffer vial and 1 instruction card

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NO COATING AREA

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents
 1 Instructions for use and 40 Single kit bags - Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (capped), 1 specimen transfer device (inverted cup type), 1 buffer vial and 1 instruction card

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Single Kit

Contents - 1 Instructions for use and 40 Single kit bags
 Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (capped), 1 specimen transfer device (inverted cup type), 1 buffer vial and 1 instruction card

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Single Kit

Contents
 1 Instructions for use and 40 Single kit bags - Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (retractable), 1 specimen transfer device (bulb type), 1 buffer vial and 1 instruction card

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CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Single Kit

Contents - 1 Instructions for use and 40 Single kit bags
 Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (retractable), 1 specimen transfer device (bulb type), 1 buffer vial and 1 instruction card

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REF RMRU-04071RI

CareStart™
Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Contents
 1 Instructions for use and 40 Single kit bags - Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (retractable), 1 specimen transfer device (inverted cup type), 1 buffer vial and 1 instruction card

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CareStart™
Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT

Single Kit

Contents - 1 Instructions for use and 40 Single kit bags
 Each bag contains 1 cassette, 1 alcohol swab, 1 lancet (retractable), 1 specimen transfer device (inverted cup type), 1 buffer vial and 1 instruction card

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