## WHO Prequalification of In Vitro Diagnostics Programme PUBLIC REPORT

**Product: Bioline HIV/Syphilis Duo<sup>1</sup>**  
**Number: PQDx 0179-012-00**

**Bioline HIV/Syphilis Duo** with product codes **06FK30, 06FK35 and 06FK37**, manufactured by **Abbott Diagnostics Korea Inc<sup>2</sup>**, **Rest-of-World regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 28 October 2015.

### Summary of prequalification status for Bioline HIV/Syphilis Duo

	Date	Outcome
<b>Status on PQ list</b>	28 October 2015	listed
<b>Dossier assessment</b>	14 July 2015	MR
<b>Site inspection(s) of quality management system</b>	24 July 2015	MR
<b>Product performance evaluation</b>	Quarter (Q) 3 2014-Q1 2015	MR

MR: Meets Requirements

### Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
4.0	Updated labelling (IFU and labels). the product was rebranded from SD BIOLINE HIV/Syphilis Duo to Alere HIV/Syphilis Duo. Corrections/clarifications on IFU. Alere logo on IFU and boxes.	15 June 2017
5.0	Rebranding of Alere HIV/Syphilis Duo to SD BIOLINE HIV/Syphilis Duo. Standard Diagnostics has been merged to Abbott during the change process of the brand name of this	07 December 2018.

<sup>1</sup> Product name was changed from SD BIOLINE HIV/Syphilis Duo to Bioline™ HIV/Syphilis Duo

<sup>2</sup> Manufacturer's name changed from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc

	product. The rebranding from SD BIOLINE HIV/Syphilis Duo to Alere HIV/Syphilis Duo (as reported in 2 December 2016) did not proceed as planned.	
6.0	Changes to the manufacturer name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc and product name from SD BIOLINE HIV/Syphilis Duo to Bioline HIV/Syphilis Duo.	20 August 2020
7.0	Addition of safety lancet resulting in a new product configuration with product code 06FK37.	22 March 2022

### Intended use:

According to Abbott Diagnostics Korea Inc, *“the Bioline HIV/Syphilis Duo test is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 including subtype-O, HIV-2 and Syphilis (Treponema pallidum) in human serum, plasma or whole blood. The Bioline HIV/Syphilis Duo is intended only for professional use, and for in vitro diagnostic use. This test may not be suitable for diagnosis of early infection or blood donation screening. Reactive specimens should be confirmed by a supplemental assay such as ELISA or Western Blot for antibodies to HIV-1/HIV-2 and the treponemal test (e.g. TPPA, TPHA) and non-treponemal test (e.g. RPR and VDRL) for antibodies to T. pallidum”.*

### Assay description

According to Abbott Diagnostics Korea Inc, *“the Bioline HIV/Syphilis Duo test is a solid phase immunochromatographic assay for the qualitative, simultaneous detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or Treponema pallidum in human serum, plasma or whole blood. The Bioline HIV/Syphilis Duo contains a nitrocellulose membrane strip with two pre-coated test line regions. Test line region HIV is precoated with recombinant HIV-1 capture antigen (gp41), recombinant HIV-2 capture antigen (gp36) and recombinant HIV-sub O antigen; test line region SYP is pre-coated with recombinant Treponema pallidum antigens (17KDa). When specimen and assay diluent is added to the specimen pad, it moves across the antigen-colloid gold conjugate pad. This mixture continues to migrate through the solid phase to the immobilized antigen on the membrane and forms a visible line. This antigen-antibody-antigen gold particle complex forms with high degree of sensitivity and specificity. If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible HIV and C lines in the window. If antibodies to syphilis are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible SYP and C lines in the window. If antibodies to HIV-1 and/or HIV-2 and TP are undetectable, the antigen-colloid gold flows past the window and one control line appears in the window”.*

## Test kit contents

	25 tests/kit Product code 06FK30	25 tests/kit Product code 06FK35	25 tests/kit Product code 06FK37
<b>Test devices</b> Individually foil pouched with a desiccant	25	25	25
<b>Assay diluent</b> Plastic bottle	1 vial of 4mL	1 vial of 4mL	1 vial of 4mL
<b>Specimen transfer devices (capillary pipettes)</b> Disposable, 20µl	N/A	25 units	25 units
<b>Lancets</b>	N/A	25 units of sterile lancets	25 units of sterile safety lancets
<b>Alcohol swabs</b> Disposable	N/A	25 units	25 units
<b>Instructions for Use</b>	1 unit	1 unit	1 unit

## Items required but not provided

- Micropipette,
- Protective gloves
- Timer
- Biohazard container

## Storage

The test kit should be stored at 1 - 30 °C.

## Shelf-life:

24 months.

## Prioritization for prequalification

Abbott Diagnostics Korea Inc submitted an application for prequalification of Bioline HIV/Syphilis Duo. Based on the established eligibility criteria, Bioline HIV/Syphilis Duo was given priority for prequalification assessment.

## Product dossier assessment

Abbott Diagnostics Korea Inc submitted a product dossier for BiolineHIV/Syphilis Duo as per the “Instructions for compilation of a product dossier” (PQDx\_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx\_009 v2). Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for Bioline HIV/Syphilis Duo for prequalification.

## Manufacturing site inspection

A comprehensive inspection was performed at the site of manufacture (Production: 65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, Republic of Korea 446-930 and Warehouse: 19-22, Dongtansandan 3-gil, Dongtan-myeon, Hwaseong-si, Gyeonggi-do, Republic of Korea) of Bioline HIV/Syphilis Duo in May 2015 as per the Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics (PQDx\_014 v1). The inspection found that the manufacturer had an acceptable quality management system and good manufacturing practices in place that ensured the consistent manufacture of a product of good quality. The manufacturer's responses to the nonconformities found at the time of the inspection were accepted 24 July 2015.

## Product performance evaluation

Bioline HIV/Syphilis Duo was evaluated at the Institute of Tropical Medicine, Antwerp, Belgium, on behalf of WHO in the last quarter of 2014/first quarter of 2015 using serum/plasma specimens. From this evaluation, we drew the following conclusions:

Bioline HIV/Syphilis Duo is an immunochromatographic assay for the detection of HIV-1/2 and syphilis antibodies in human serum/plasma, fingerstick/venous whole blood specimens. A volume of 10 µL of serum/plasma is needed to perform the assay. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities and non-laboratory settings. Reading of the results is done visually i.e. subjectively read.

In this limited evaluation a panel of 400 clinically-derived specimens was used. For antibodies to HIV, we observed an initial sensitivity (N=200) of 100% (95% CI: 98.2% - 100%) and an initial specificity (N=200) of 99.5% (95% CI: 97.2% - 100%) compared to the reference assays (Vironostika HIV Ag/Ab [bioMérieux] and Enzygnost anti-HIV 1/2 [Siemens Healthcare Diagnostics] in parallel; followed by INNO-LIA HIV I/II Score [Fujirebio]). The final sensitivity

was 100% (95% CI: 98.2% - 100%) and the final specificity was 99.5% (95% CI: 97.2% - 100%) compared to the reference assays. Lot to lot variation was acceptable.

For the case of antibodies to *Treponema pallidum*, we observed an initial sensitivity\* (N=200) of 86.5% (95% CI: 81.0% - 90.9%) and an initial specificity (N=200) of 99.5% (95% CI: 97.2% - 100%) compared to the reference assays (Vitros Syphilis [Ortho clinical Diagnostics] followed by Serodia TPPA [Fujirebio] to confirm initial reactive results). The final sensitivity was 87.0% (95% CI: 81.5% - 91.3%) and the final specificity was 99.5% (95% CI: 97.2% - 100%) compared to the reference assays. Lot to lot variation was acceptable.

For HIV, eight seroconversion panels were evaluated. Bioline HIV/Syphilis Duo detected on average 0.125 specimens earlier than the benchmark assay (Enzygnost Anti-HIV 1/2 Plus).

For syphilis one seroconversion panel was evaluated. Bioline HIV/Syphilis Duo detected 2 specimens earlier than the benchmark assay (Vitros Syphilis).

For the HIV mixed titer panel, Bioline HIV/Syphilis Duo correctly identified all specimens compared to the reference result.

For the syphilis mixed titer panel, Bioline HIV/Syphilis Duo correctly identified all specimens compared to the reference result.

For the 1st International Reference Panel for anti-HIV [NIBSC code 02/210] Bioline HIV/Syphilis Duo gave an indeterminate result for the HIV-1 group O specimen. Bioline HIV/Syphilis Duo correctly identified all HIV-1 group M subtypes and HIV-2 specimens.

For the 1st International Reference Panel for anti-treponemal [NIBSC code 05/132], Bioline HIV/Syphilis Duo correctly identified all specimens.

In this study, 2.5% of the results for syphilis were recorded as indeterminate, 0% for HIV. Results were interpreted independently by three technicians; the inter-reader variability was 0% for HIV and 4% for syphilis. The invalid rate was 0%.


## **Labelling**

- 1. Labels**
- 2. Instructions for use**


*\*Note: Labelling has been changed as per site name change and product rebranding. But temporarily labelling of legacy brand (SD BIOLINE) will be used in the market according to registration status in each country.*

1. Labels

1.1 Package box for product 06FK30


### 1.2 Package box for product code 06FK35



### 1.3 Package box for product code 06FK37

### 1.4 Device pouch for 06FK30, 06FK35

### 1.5 Bioline HIV/Syphilis Duo test device

### 1.6 Sterile lancet labels ( component in product code 06FK35)


## Sterile Lancet for single use

 LOT NO./Type :20051755-08/28G

 EXP. DATE : 2025-06-21

 QTY :25pcs

[STERILE][R]

Shanghai International Holding Corp. GmbH (Europe)  
Eiffelstrasse 80, 20537 Hamburg, Germany

Shandong Lianfa Medical Plastic Products Co., Ltd.  
No 1 Shuangshan Sanjian Road 250200 Zhangqiu City,  
Jinan, Shandong, PEOPLE'S REPUBLIC OF CHINA

10mm

33.75mm

67.5mm

33.75mm

10mm

10mm

75mm

10mm

75mm

10mm

/28G

LOT NO./Type :  
EXP. DATE :

### Disposable Sterile Lancets

QTY : 25pcs

STERILE R Irradiation Sterilization

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany

**Manufactured by:**  
Tianjin Huahong Technology Co., Ltd.  
**Address:**  
A01, Plant B, No. 278, Hangkong Road, Tianjin Pilot Free Trade  
Zone (Air Port Industrial Park) 300308 Tianjin China

/28G

LOT NO./Type :  
EXP. DATE :

### Disposable Sterile Lancets

QTY : 25pcs

STERILE R Irradiation Sterilization

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany

**Manufactured by:**  
Tianjin Huahong Technology Co., Ltd.  
**Address:**  
A01, Plant B, No. 278, Hangkong Road, Tianjin Pilot Free Trade  
Zone (Air Port Industrial Park) 300308 Tianjin China

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

Twist off and remove the protective cap

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Eiffestraße 80 20537 Hamburg GERMANY


**Manufactured by:**

Shandong Lianfa Medical Plastic Products Co., Ltd.


**Address:**

No 1 Shuangshan Sanjian Road 250200 Zhangqiu City, Jinan,  
Shandong PEOPLE'S REPUBLIC OF CHINA


## Instructions for use :

Twist off the protective cap.


Push the needle into fingertip  
to get blood sample.

Push the cap into the safety lancet back  
and dispose of the used one to a trash bin.

### 1.8 Assay diluent

### 1.9 Alcohol swabs

#### 1.9.1 FA Alcohol Swab-S

Front

Back

### 1.9.2 Alcohol Prep Pad


### 1.9.3 AIM PADS


### 1.10 Capillary tube labels

#### 1.10.1 AIM

#### 1.10.2 Jianguo Kehua

## **2. Instructions for use<sup>3</sup>**

---

<sup>3</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

## **2.1 IFU for product codes 06FK30 and 06FK35.**

**Abbott**

**REF** 06FK30, 06FK35

**Bioline™**

# **HIV/Syphilis Duo**

**Antibodies to HIV-1/2 and *Treponema pallidum* Test**

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**Anticorpos para teste de HIV-1/2 e *Treponema pallidum***

## INTENDED USE

The Bioline™ HIV/Syphilis Duo test is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 including subtype-O, HIV-2 and Syphilis (*Treponema pallidum*) in human serum, plasma or whole blood. The Bioline™ HIV/Syphilis Duo is intended only for professional use, and for *in vitro* diagnostic use. This test may not be suitable for diagnosis of early infection or blood donation screening. Reactive specimens should be confirmed by a supplemental assay such as ELISA or Western Blot for antibodies to HIV-1/HIV-2 and the treponemal test (e.g. TPPA, TPHA) and non-treponemal test (e.g. RPR and VDRL) for antibodies to *T. pallidum*.

## INTRODUCTION

HIV (Human Immunodeficiency Virus) is recognized as the etiologic agent of Acquired Immune Deficiency Syndrome (AIDS). The virus is transmitted by sexual contact, exposure to infected blood, certain body fluids or tissues, and from mother to fetus or child during the perinatal period. HIV-1 has been isolated from patients with AIDS and AIDS related complex, and from healthy persons with high potential risk of developing AIDS. Patients with HIV-2 are found primarily in parts of West Africa. Its course is marked by increasing levels of viral replication and the emergence of more virulent viral strains. HIV-1 and HIV-2 are similar in their morphology, cell tropism, host interaction and generic structure. Serological studies have determined that HIV-1 and HIV-2 have multiple common epitopes in core antigens but much less so in the envelope antigens. This clinical diagnosis of HIV may include the detection of antibodies to HIV-1/2 in human plasma or serum by immunoassay. The presence of HIV can be identified by detection of antibodies to HIV 1/2 in human serum, plasma and whole blood by immunoassay. This advanced assay utilizes recombinant antigens targeted against immunogenic proteins. The major immunoreactive antigens of these proteins are HIV-1 gp41, p24 and HIV-2 gp36.

The spirochetal bacillus *Treponema pallidum* (TP) is the causative agent of the venereal disease syphilis. The presence of TP can be identified by detection of antibodies to TP in human blood by immunoassay. This advanced assay utilizes recombinant antigens targeted against immunogenic TP membrane proteins.<sup>7</sup>

## TEST PRINCIPLE

The Bioline™ HIV/Syphilis Duo test is a solid phase immunochromatographic assay for the qualitative, simultaneous detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or *Treponema pallidum* in human serum, plasma or whole blood. The Bioline™ HIV/Syphilis Duo contains a nitrocellulose membrane strip with two pre-coated test line regions. Test line region HIV is pre-coated with recombinant HIV-1 capture antigen (gp41), recombinant HIV-2 capture antigen (gp36) and recombinant HIV-sub O antigen; test line region SYP is pre-coated with recombinant *Treponema pallidum* antigens (17KDa). When specimen and assay diluent is added to

the specimen pad, it moves across the antigen-colloid gold conjugate pad. This mixture continues to migrate through the solid phase to the immobilized antigen on the membrane and forms a visible line. This antigen-antibody-antigen gold particle complex forms with high degree of sensitivity and specificity. If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible HIV and C lines in the window. If antibodies to syphilis are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible SYP and C lines in the window. If antibodies to HIV-1 and/or HIV-2 and TP are undetectable, the antigen-colloid gold flows past the window and one control line appears in the window.

## MATERIALS PROVIDED AND ACTIVE INGREDIENTS OF MAIN COMPONENTS

1. The Bioline™ HIV/Syphilis Duo test kit contains the following items to perform the assay:

Catalog No.	06FK30	06FK35
<b>Contents</b>	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 4 ml/vial)  3. 1 Instructions for use	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 4 ml/vial)  3. 25 Capillary pipettes (20 µl), 25 Sterile lancets, 25 Alcohol swabs 4. 1 Instructions for use

2. Active ingredients of main components:

1 test strip includes:	Assay diluent includes:
<ul style="list-style-type: none"> <li>- Gold conjugates: Recombinant HIV-1 antigen – gold colloid (0.029±0.006 µg), Recombinant HIV-2 antigen – gold colloid (0.043±0.009 µg), Recombinant Treponema pallidum antigen – gold colloid (0.017±0.003 µg), Chicken IgY– gold colloid (0.038±0.008 µg)</li> <li>- Test line HIV: Recombinant HIV-1 antigen (gp41) (0.384±0.077 µg), Recombinant HIV-2 antigen (gp36) (0.640±0.128 µg), Recombinant HIV - sub O antigen (0.224±0.045 µg)</li> <li>- Test line SYP: Recombinant Treponema pallidum antigen (17KDa) (0.480±0.096 µg)</li> <li>- Control line: Mouse monoclonal anti-chicken IgY (0.640±0.128 µg)</li> </ul>	<ul style="list-style-type: none"> <li>- 20 mM Tris-HCl buffer (100 mM)</li> <li>- Tricine (100 mM)</li> <li>- Casein (0.2 w/v%)</li> <li>- Proclin™ 300 (0.05 w/v%)</li> <li>- Sucrose (5 w/v%)</li> <li>- Sodium azide (0.01 w/v%)</li> </ul>

- Classification according to Regulation (EC) No 1272/2008:

<b>Product identifier</b>	<b>Trade name</b>	Bioline™ HIV/Syphilis Duo Assay diluent		
	<b>Hazardous substance</b>	Proclin™ 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)		
<b>Classification</b>	Skin sensitization Category 1		<b>Signal Word</b>	Warning
<b>Hazard Pictogram</b>			<b>Hazard Statement</b>	H317: May cause an allergic skin reaction
<b>Precautionary statements</b>				
<b>Prevention</b>	P261: Avoid breathing dust/fume/gas/mist/vapours/spray P272: Contaminated work clothing should not be allowed out of the workplace P280: Wear protective gloves/protective clothing/eye protection/face protection			
<b>Response</b>	P302+P352: IF ON SKIN: Wash with plenty of soap and water P321: Specific treatment P333+P313: If skin irritation or rash occurs: Get medical advice/attention P363: Wash contaminated clothing before reuse			
<b>Disposal</b>	P501: Dispose of contents/container in accordance with local/regional/national/international regulation			

## MATERIALS REQUIRED BUT NOT PROVIDED

- Micropipette, Protective gloves, Timer, Biohazard container

## KIT STORAGE AND STABILITY

- The test kit should be stored at a temperature between 1 °C and 30 °C. Do not freeze the kit or its components. Assay diluent cap should be kept firmly sealed between each use.
- Assay diluent may be opened and resealed for each assay. It is stable until expiration date if kept at 1 - 30 °C.
- The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.

4. Do not use the test kit beyond its expiration date. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

## WARNINGS

1. The test devices are for *in vitro* diagnostic use only. Do not reuse the test device.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be proficient.
3. Do not use the pipette by mouth, smoke, drink, eat, apply cosmetics, or handle contact lenses in areas where specimens or kit components are being handled.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
7. Do not mix or interchange different specimens.
8. Do not eat the desiccant from the foil pouch.
9. Avoid splashing or aerosol formation of specimen and assay diluent.
10. Do not mix or interchange components among different lots or those for other products.
11. Do not drink assay diluent.
12. Care should be taken to avoid contamination of the bottle nozzle when dropping assay diluent into the specimen well.
13. The assay diluent contains a proprietary antimicrobial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.
14. The assay diluent contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of these reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
15. Safety data sheet available for professional user on request.

## SPECIMEN COLLECTION AND HANDLING

### 1. Whole blood

#### [Collection by venipuncture]

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
- If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the blood specimen must be tested within 3 days of refrigeration.
- Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.
- Bring blood specimens to room temperature (15 - 30 °C) prior to use.

#### [Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip then prick the lateral side of the finger with the sterile lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new 20 µl capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

### 2. Plasma or serum

- If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8 °C. For storage period longer than 2 weeks, freezing (below -20 °C) is required. Bring plasma or serum specimens to room temperature (15 - 30 °C) prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

#### [Plasma]

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) then centrifuge the tube to generate a plasma specimen.

#### [Serum]

- Using venipuncture, draw whole blood into the collection tube (NOT containing anticoagulants) then leave for 30 minutes to allow blood coagulation to occur. Centrifuge the tube to generate a serum specimen.

\*Anticoagulants including heparin, EDTA and sodium citrate do not affect the test result. Use of other anticoagulants has not been validated. Their use may affect the test result.

## TEST PROCEDURE (REFER TO FIGURE)

1. Allow all kit components and specimens to reach a temperature between 15 °C and 30 °C prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. [Using a micropipette]  
Dispense 10 µl of plasma or serum specimen or dispense 20 µl of whole blood specimen into the specimen well marked “S”.

Or,

[Using a capillary pipette]

Dispense 20 µl of drawn whole blood specimen into the specimen well marked “S”.

**Caution:** Lightly touch the capillary pipette to the specimen pad while dispensing.

4. Dispense 3 drops (approximately 100 µl) of assay diluent into the specimen well marked “S”.

**Caution:** Do not let bottle nozzle touch device in order to avoid cross-contamination. Hold bottle vertically while dispensing. If you do not hold the bottle vertically, it can lead to inaccurate results. Exactly, 3 drops should be added. Adding more than 3 drops may result in reddish color background or an invalid result.

5. As the test begins to work, you will see purple color move across the result window in the center of the test device.
6. Interpret test results 15 - 20 minutes after adding assay diluent. Do not read after 20 minutes.

**Caution:** If the test result is not legible after 15 minutes due to high background color, read again later but within 20 minutes of adding the diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.

## TEST INTERPRETATION (REFER TO FIGURE)

The characters “HIV”, “SYP” and “C” are printed on the test device. They correspond to the positions of the HIV, Syphilis and control lines in the test window.

- **Nonreactive result:**

The presence of only one purple colored line (C) within the result window indicates a nonreactive result.

- **Reactive result:**

**Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

- 1) HIV-1/2 reactive: The presence of both test line HIV and the control line (C) indicates a reactive result for anti-HIV-1/2.
- 2) Syphilis reactive: The presence of both test line SYP and the control line (C) indicates a reactive result for anti-*Treponema pallidum*.

3) HIV-1/2 and Syphilis reactive: The presence of test line HIV, test line SYP and the control line (C) indicates reactive result for both anti-HIV-1/2 and anti-*Treponema pallidum*.

- **Invalid result:**

Absence of the control line (C) and/or presence of a pink/purple smear in the result window indicate an invalid result. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

## TEST LIMITATIONS

1. A reactive result with the Bioline™ HIV/Syphilis Duo test indicates presence of antibodies to *Treponema pallidum*. A reactive result alone is insufficient for diagnosis of Syphilis infection. Results must be interpreted along with other clinical information available to the physician.
2. A reactive result may indicate infection with HIV-1/2. Immunochromatographic testing alone cannot be used to diagnose AIDS. An AIDS diagnosis can only be made in a clinical setting if an individual meets the case definition for AIDS established by the Centers for Disease Control (CDC). Reactive specimens should be confirmed by a supplemental assay such as ELISA or Western Blot test.
3. A nonreactive result does not eliminate the possibility of infection with HIV-1/2 and/or *Treponema pallidum*. The specimen may contain low levels of antibodies that cannot be detected by Bioline™ HIV/Syphilis Duo. If a test result is nonreactive and clinical symptoms persist, additional testing using other clinical methods is recommended.
4. Some known HIV-infected persons taking antiretroviral medication have been shown to produce false negative results when tested by rapid diagnostic tests.<sup>9,10</sup>
5. Where clinical presentation or other data would suggest an inconsistent test result then the individual should be tested by nucleic acid testing (NAT) technologies immediately and/or retested for antibodies to HIV after more than 21 days since the original testing.

## INTERNAL QUALITY CONTROL

The Bioline™ HIV/Syphilis Duo test device has test line HIV, test line SYP and the control line C on the surface of the device. The lines are not visible before applying a specimen. The control line is used for procedural control and shows that the diluent has been applied successfully and that the active ingredients of main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a reactive specimen control.

## PERFORMANCE CHARACTERISTICS

### Performance Summary

Total	HIV	Syphilis
Diagnostic Sensitivity	99.80 %	90.00 %*
Clinical Sensitivity	99.91 %	99.67 %
Diagnostic Specificity	100.00 %	99.90 %
Clinical Specificity	99.67 %	99.72 %

\*Diagnostic Sensitivity on HIV/Syphilis co-infected patient specimens was lower than on clinical specimens not co-infected.

1. Diagnostic sensitivity  
Clinical performance studies were performed by two external evaluation centres, the Institute of Tropical Medicine in Belgium and the German Red Cross in Germany.
- 1) Anti-HIV-1 and anti-HIV-2 reactive specimens  
400 anti-HIV-1 reactive specimens, among which 40 non-B subtype and 100 anti-HIV-2 reactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 500 reactive specimens, is 99.8%.

	Bioline™ HIV/Syphilis Duo		Sensitivity (95 % CI)
	HIV nonreactive	HIV reactive	
Anti-HIV-1 reactive	1	259	
Anti-HIV-1 reactive non-B subtype	0	40	
Anti-HIV-1 and anti-TP reactive	0	100	
Anti-HIV-2 reactive	0	100	
<b>Total</b>	1	499	

- 2) Anti-*Treponema pallidum* reactive specimens  
150 anti-*Treponema pallidum* reactive specimens, among which 100 anti-HIV-1 reactive and 50 anti-HIV-1 nonreactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-*Treponema pallidum* antibody detection, calculated on 150 confirmed reactive specimens is 90.0%.

	Bioline™ HIV/Syphilis Duo		
	Syphilis nonreactive	Syphilis reactive	Sensitivity (95 % CI)
Anti-HIV-1 reactive	13	87	
Anti-HIV-1 nonreactive	2	48	
<b>Total</b>	15	135	90.0 % (84.0 - 94.3 %)

2. Diagnostic specificity

In total 1,000 EDTA plasma and 500 whole blood specimens from blood donors were tested. The specimens originated from 2 collection sites in Germany; Frankfurt and Kassel.

		HIV		Syphilis	
		Nonreactive	Reactive	Nonreactive	Reactive
Blood donors	EDTA Plasma	1,000	0	1,000	0
	Whole blood	500	0	499	1
<b>Specificity (95 % CI)</b>		100 % (99.7 - 100 %)		99.9 % (99.6 - 100 %)	

3. Laboratories in Ghana, Mexico, Laos, Togo, Kenya, and Myanmar participated in the evaluation during 2012–2013. Each site characterized sera using *T. pallidum* particle agglutination assay or *T. pallidum* hemagglutination assay and HIV enzyme immunoassay, Western Blot, and/or HIV antibody rapid tests.

1) Laboratory Performance for Detection of Anti-HIV

Country	Year	n	True Reactive	False Reactive	False Nonreactive	True Nonreactive	Sensitivity Estimate (95 % CI)	Specificity Estimate (95 % CI)
Ghana	2012	400	250	0	0	150	100 % (98.54 - 100 %)	100 % (97.57 - 100 %)
Togo	2013	310	203	0	0	107	100 % (98.20 - 100 %)	100 % (96.61 - 100 %)
Myanmar	2013	245	114	1	0	130	100 % (96.82 - 100 %)	99.24 % (95.82 - 99.98 %)
Kenya	2013	698	345	0	1	352	99.71 % (98.40 - 99.99 %)	100 % (98.96 - 100 %)
Mexico	2013	527	158	0	0	369	100 % (97.69 - 100 %)	100 % (99.01 - 100 %)
Laos	2013	156	53	3	0	100	100 % (93.28 - 100 %)	97.09 % (91.72 - 99.40 %)
<b>Total</b>	-	2,336	1,123	4	1	1,208	99.91 % (99.51 - 100 %)	99.67 % (99.16 - 99.91 %)

## 2) Laboratory performance for Detection of Anti-*Treponema pallidum*

Country	Year	n	True Reactive	False Reactive	False Nonreactive	True Nonreactive	Sensitivity Estimate (95 % CI)	Specificity Estimate (95 % CI)
Ghana	2012	400	250	1	0	149	100 % (98.54 - 100 %)	99.33 % (96.34 - 99.98 %)
Togo	2013	241	88	1	0	152	100 % (95.89 - 100 %)	99.35 % (96.41 - 99.98 %)
Myanmar	2013	200	74	1	1	124	98.67 % (92.79 - 99.97 %)	99.20 % (95.62 - 99.98 %)
Kenya	2013	698	85	0	0	613	100 % (95.75 - 100 %)	100 % (99.40 - 100 %)
Mexico	2013	414	106	1	1	306	99.07 % (94.90 - 99.98 %)	99.67 % (98.20 - 99.99 %)
Laos	2013	106	6	0	0	100	100 % (54.07 - 100 %)	100 % (96.38 - 100 %)
Total	-	2,059	609	4	2	1,444	99.67 % (98.82 - 99.96 %)	99.72 % (99.29 - 99.92 %)

### 4. Sensitivity on seroconversion panels

30 commercially available HIV seroconversion panels were tested showing how early detection can be obtained with the Bioline™ HIV/Syphilis Duo test. 20 HIV Ab panels (pertaining 62 early seroconversion samples) were from the Institute of Tropical Medicine and 10 panels (pertaining 12 early seroconversion samples) were from the Paul-Ehrlich-Institute. The performance of the Bioline™ HIV/Syphilis Duo test on seroconversion panels is comparable to other CE marked Anti-HIV-1/2 screening and rapid tests.

### 5. Analytical specificity

#### • Cross reactivity

The following 16 potential cross-reacting pathogens had no effect on test results of Bioline™ HIV/Syphilis Duo.

Anti-HBs	EBV	HIV (to SYP line)	Syphilis (to HIV line)
<i>Borrelia burgdorferi</i>	HAV	HTLV	Toxoplasma
Chlamydia	HBsAg	Influenza	<i>Trypanosoma cruzi</i> I / II
CMV	HCV	<i>Plasmodium vivax</i>	<i>Trypanosoma gambiense</i>

#### • Interfering substances

The following 9 substances had no effect on test results of Bioline™ HIV/Syphilis Duo.

: Pregnant women, High cholesterol (245 mg/dL), High bilirubin (1.4 mg/dL), Rheumatoid factor (28 IU/ml), Lipemic, Hemolyzed, Autoimmune, Alcoholic cirrhosis, Multiparous pregnancy

#### • Pharmaceutical substances


There was no significant interference with the following 25 drugs on Bioline™ HIV/Syphilis Duo. All drugs were tested at concentrations of 250 µg/ml.

Abacavir	Cholecalciferol	Folic acid	L-ascorbic acid	Pantoprazole
Acetaminophen	Cyclobenzaprine	Hydrochlorothiazide	Magnesium sulfate	Pyrazinamide
Acetylsalicylic acid	Darunavir	Ibuprofen	Metformin	Rifampicin
Amoxicillin	Diclofenac	Iron Chloride	Naproxen	Ritonavir
Aspirin	Ergocalciferol	Isoniazid	Nevirapine	Salicylic acid

6. Reproducibility of the Bioline™ HIV/Syphilis Duo has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

# PREPARATION

Open the foil pouch and look for the following:

1. Test device
2. Desiccant

Then, label the device with the patient identifier.


● SP10 / WB 20  $\mu$ l : Serum 10  $\mu$ l or Plasma 10  $\mu$ l or Whole blood 20  $\mu$ l

□▶●●● : Assay diluent 3 drops


# TEST PROCEDURE

## I. Blood (by venipuncture), Plasma or Serum specimen


- 1 Using a micropipette draw plasma or serum: **10  $\mu$ l** specimen  
or  
whole blood: **20  $\mu$ l** specimen.

- 2** Dispense  
**10  $\mu$ l** of plasma or serum specimen  
or  
**20  $\mu$ l** of whole blood specimen  
into the specimen well marked “S”.

- 3** Dispense **3 drops** (approximately 100  $\mu$ l) of assay diluent into the specimen well marked “S”. Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.


- 4** Interpret test results 15 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.


## II. Blood specimen (with a lancet)


**1** Clean the area to be lanced with an alcohol swab.


**2** Squeeze the fingertip then prick the lateral side of the finger with the sterile lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.


**3** Immerse the open end of a new **20 µl** capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

- 4** Dispense **20  $\mu$ l** of drawn whole blood specimen in the specimen well marked “S”. Lightly touch the capillary pipette to the specimen pad while dispensing.

- 5** Dispense **3 drops** (approximately 100  $\mu$ l) of assay diluent into the specimen well marked “S”. Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.


- 6** Interpret test results 15 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.

## INTERPRETATION

### Nonreactive

The presence of only one purple colored line (C) within the result window indicates a nonreactive result.



### Reactive

 **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.



### HIV-1/2 Reactive

Presence of both “C” and “HIV” lines





## Syphilis Reactive

Presence of both “C” and “SYP” lines



## HIV-1/2 and Syphilis Reactive

Presence of “C”, “HIV” and “SYP” lines

## Invalid

Absence of the control line (C) and/or presence of a pink/purple smear in the result window indicate an invalid result. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

### **Product. Disclaimer:**

While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

### **Warning:**

The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive diagnosis using this product.

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1. Miller JN : Value and limitation of nontreponemal and treponemal tests in the laboratory diagnosis of syphilis. Clin. Obstet Gynecol 18: 191-203, 1975.
2. Lefevre J-C, Bertrand M-A, Baogriaud R : Evaluation of Captia enzyme immunoassays for detection of immunoglobulins G and M to *Treponema pallidum* in syphilis. J Clin Microbiol 28 : 1704-1707, 1990.
3. V.S. Ivanov, Z.K. Suvorova, L.D. Tchikin, A.T. Kozitch and V.T. Ivanov: Effective method for synthetic peptide immobilization that increases the sensitivity and specificity of ELISA procedures Journal of Immunological Methods, 153 (1992) 229-233.
4. Mi Jin Sohn, Young Hae Chong, Ji Eun Chang, Young Ik Lee: Overexpression and simple purification of human immunodeficiency virus-1 gag epitope derived from a recombinant antigen in *E. coli* and its use in ELISA. Journal of Biotechnology 34 (1994) 149-155.
5. Laboratory diagnostics for HIV Infection. Butto S. et al. Centro Nazionale AIDS, Istituto Superiore di Sanita, Rome, Italy. Ann Ist Super Sanita. 2010 Vol. 46 no.1: 24-33.
6. HIV Diagnostic Testing: Evolving Technology and Testing Strategies. Daskalakis D. et al. 2011. New York University Medical Center, USA. Top Antivir Med. 2011 Vol. 19 (1): 18-22.
7. Laboratorial Diagnosis of Syphilis, Syphilis – Recognition, Description and Diagnosis. 2011 Neuzha Satomi Sato et al. (Ed.), ISBN: 978-953-307-554-9, InTech.
8. The laboratory diagnosis of Syphilis. 2005 Ratnam et al. Can J Infect Dis Med Microbiol Vol 16 No 1 January/February 2005.
9. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
10. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.
11. Websites For Peer reviews  
- [http://www.ijgo.org/article/S0020-7292\(15\)00201-5/references](http://www.ijgo.org/article/S0020-7292(15)00201-5/references)  
- [http://sti.bmj.com/content/89/Suppl\\_1/A56.1](http://sti.bmj.com/content/89/Suppl_1/A56.1)  
- <http://www.ncbi.nlm.nih.gov/pubmed/25908677>  
- <http://www.ncbi.nlm.nih.gov/pubmed/25734088>  
- <http://www.ncbi.nlm.nih.gov/pubmed/25339379>  
- <http://www.ncbi.nlm.nih.gov/pubmed/25391806>

**Abbott Diagnostics Korea Inc.**

65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, 17099,  
Republic of Korea

Tel : 82-31-899-2800 | Fax : 82-31-899-2840

[abbott.com/poct](http://abbott.com/poct)

Date issued : 2019. 10  
06FK30/06FK35-04-A0

© 2019 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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## TEST PRINCIPLE

The Bioline™ HIV/Syphilis Duo test is a solid phase immunochromatographic assay for the qualitative, simultaneous detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1/2 and/or *Treponema pallidum* in human serum, plasma or whole blood. The Bioline™ HIV/Syphilis Duo contains a nitrocellulose membrane strip with two pre-coated test line regions. Test line region HIV is pre-coated with recombinant HIV-1 capture antigen (gp41), recombinant HIV-2 capture antigen (gp36) and recombinant HIV-sub O antigen; test line region SYP is pre-coated with recombinant *Treponema pallidum* antigens (17KDa). When specimen and assay diluent is added to the specimen pad, it moves across the antigen-colloid gold conjugate pad. This mixture continues to migrate through the solid phase to the immobilized antigen on the membrane and forms a visible line. This antigen-antibody-antigen gold particle complex forms with high degree

of sensitivity and specificity. If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible HIV and C lines in the window. If antibodies to syphilis are present in the specimen, the antibodies bind to the antigen-colloid gold and to the immobilized recombinant antigens, forming visible SYP and C lines in the window. If antibodies to HIV-1 and/or HIV-2 and TP are undetectable, the antigen-colloid gold flows past the window and one control line appears in the window.

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- The Bioline™ HIV/Syphilis Duo test kit contains the following items to perform the assay:
  - 25 Test devices with desiccant in individual foil pouches
  - Assay diluent (1 x 4 ml/vial) ⚠
  - 25 Capillary pipettes (20 µl), 25 Safety lancets, 25 Alcohol swabs
  - 1 Instructions for use
- Active ingredients of main components:

1 test strip includes:	Assay diluent includes:
<ul style="list-style-type: none"> <li>- Gold conjugates: Recombinant HIV-1 antigen – gold colloid (0.029±0.006 µg), Recombinant HIV-2 antigen – gold colloid (0.043±0.009 µg), Recombinant Treponema pallidum antigen – gold colloid (0.017±0.003 µg), Chicken IgY– gold colloid (0.038±0.008 µg)</li> <li>- Test line HIV: Recombinant HIV-1 antigen (gp41) (0.384±0.077 µg), Recombinant HIV-2 antigen (gp36) (0.640±0.128 µg), Recombinant HIV - sub O antigen (0.224±0.045 µg)</li> <li>- Test line SYP: Recombinant Treponema pallidum antigen (17kDa) (0.480±0.096 µg)</li> <li>- Control line: Mouse monoclonal anti-chicken IgY (0.640±0.128 µg)</li> </ul>	<ul style="list-style-type: none"> <li>- 20 mM Tris-HCl buffer (100 mM)</li> <li>- Tricine (100 mM)</li> <li>- Casein (0.2 w/v%)</li> <li>- Proclin™ 300 (0.05 w/v%)</li> <li>- Sucrose (5 w/v%)</li> <li>- Sodium azide (0.01 w/v%)</li> </ul>

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#### [Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new 20 µl capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

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## TEST PROCEDURE (REFER TO FIGURE)

1. Allow all kit components and specimens to reach a temperature between 15 °C and 30 °C prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. [Using a micropipette]  
Dispense 10 µl of plasma or serum specimen or dispense 20 µl of whole blood specimen into the specimen well marked “S”.  
Or,  
[Using a capillary pipette]  
Dispense 20 µl of drawn whole blood specimen into the specimen well marked “S”.  
⚠ **Caution:** Lightly touch the capillary pipette to the specimen pad while dispensing.
4. Dispense 3 drops (approximately 100 µl) of assay diluent into the specimen well marked “S”.  
⚠ **Caution:** Do not let bottle nozzle touch device in order to avoid cross-contamination. Hold bottle vertically while dispensing. If you do not hold the bottle vertically, it can lead to inaccurate results. Exactly, 3 drops should be added. Adding more than 3 drops may result in reddish color background or an invalid result.
5. As the test begins to work, you will see purple color move across the result window in the center of the test device.
6. Interpret test results 15 - 20 minutes after adding assay diluent. Do not read after 20 minutes.  
⚠ **Caution:** If the test result is not legible after 15 minutes due to high background color, read again later but within 20 minutes of adding the diluent. Reading outside of this time frame (before 15 min or after 20 min) may provide false results.

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- **Reactive result:**

⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

- 1) HIV-1/2 reactive: The presence of both test line HIV and the control line (C) indicates a reactive result for anti-HIV-1/2.
- 2) Syphilis reactive: The presence of both test line SYP and the control line (C) indicates a reactive result for anti-*Treponema pallidum*.
- 3) HIV-1/2 and Syphilis reactive: The presence of test line HIV, test line SYP and the control line (C) indicates reactive result for both anti-HIV-1/2 and anti-*Treponema pallidum*.

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1. Diagnostic sensitivity

Clinical performance studies were performed by two external evaluation centres, the Institute of Tropical Medicine in Belgium and the German Red Cross in Germany.

1) Anti-HIV-1 and anti-HIV-2 reactive specimens

400 anti-HIV-1 reactive specimens, among which 40 non-B subtype and 100 anti-HIV-2 reactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 500 reactive specimens, is 99.8%.

	Bioline™ HIV/Syphilis Duo		Sensitivity (95 % CI)
	HIV nonreactive	HIV reactive	
Anti-HIV-1 reactive	1	259	
Anti-HIV-1 reactive non-B subtype	0	40	
Anti-HIV-1 and anti-TP reactive	0	100	
Anti-HIV-2 reactive	0	100	
<b>Total</b>	1	499	99.8 % (98.9 - 100 %)

2) Anti-*Treponema pallidum* reactive specimens

150 anti-*Treponema pallidum* reactive specimens, among which 100 anti-HIV-1 reactive and 50 anti-HIV-1 nonreactive specimens were tested in the Institute of Tropical Medicine, Belgium. The Diagnostic Sensitivity for anti-*Treponema pallidum* antibody detection, calculated on 150 confirmed reactive specimens is 90.0%.

	Bioline™ HIV/Syphilis Duo		
	Syphilis nonreactive	Syphilis reactive	Sensitivity (95 % CI)
Anti-HIV-1 reactive	13	87	
Anti-HIV-1 nonreactive	2	48	
<b>Total</b>	15	135	90.0 % (84.0 - 94.3 %)

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
- Reproducibility of the Bioline™ HIV/Syphilis Duo has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

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


- 2** Dispense **10  $\mu$ l** of plasma or serum specimen or **20  $\mu$ l** of whole blood specimen into the specimen well marked “S”.

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

**2** Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.


OR


OR


- 3** Immerse the open end of a new **20 µl** capillary pipette in the next blood drop and release the pressure to draw blood into the capillary pipette up to the black line.

- 4** Dispense **20 µl** of drawn whole blood specimen in the specimen well marked "S". Lightly touch the capillary pipette to the specimen pad while dispensing.

- 5** Dispense **3 drops** (approximately 100 µl) of assay diluent into the specimen well marked "S". Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.

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


### Reactive

**⚠ Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

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






C SYP HIV

### **Product. Disclaimer:**

*While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.*

### **Warning:**

*The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive diagnosis using this product.*

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1. Miller JN : Value and limitation of nontreponemal and treponemal tests in the laboratory diagnosis of syphilis. Clin. Obstet Gynecol 18: 191-203, 1975.
2. Lefevre J-C, Bertrand M-A, Baogriaud R : Evaluation of Captia enzyme immunoassays for detection of immunoglobulins G and M to *Treponema pallidum* in syphilis. J Clin Microbiol 28 : 1704-1707, 1990.
3. V.S. Ivanov, Z.K. Suvorova, L.D. Tchikin, A.T. Kozitch and V.T. Ivanov: Effective method for synthetic peptide immobilization that increases the sensitivity and specificity of ELISA procedures Journal of Immunological Methods, 153 (1992) 229-233.
4. Mi Jin Sohn, Young Hae Chong, Ji Eun Chang, Young Ik Lee: Overexpression and simple purification of human immunodeficiency virus-1 gag epitope derived from a recombinant antigen in *E. coli* and its use in ELISA. Journal of Biotechnology 34 (1994) 149-155.
5. Laboratory diagnostics for HIV Infection. Butto S. et al. Centro Nazionale AIDS, Istituto Superiore di Sanita, Rome, Italy. Ann Ist Super Sanita. 2010 Vol. 46 no.1: 24-33.
6. HIV Diagnostic Testing: Evolving Technology and Testing Strategies. Daskalakis D. et al. 2011. New York University Medical Center, USA. Top Antivir Med. 2011 Vol. 19 (1): 18-22.
7. Laboratorial Diagnosis of Syphilis, Syphilis – Recognition, Description and Diagnosis. 2011 Neuza Satomi Sato et al. (Ed.), ISBN: 978-953-307-554-9, InTech.
8. The laboratory diagnosis of Syphilis. 2005 Ratnam et al. Can J Infect Dis Med Microbiol Vol 16 No 1 January/February 2005.
9. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
10. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.
11. Websites For Peer reviews
  - [http://www.ijgo.org/article/S0020-7292\(15\)00201-5/references](http://www.ijgo.org/article/S0020-7292(15)00201-5/references)
  - [http://sti.bmj.com/content/89/Suppl\\_1/A56.1](http://sti.bmj.com/content/89/Suppl_1/A56.1)
  - <http://www.ncbi.nlm.nih.gov/pubmed/25908677>
  - <http://www.ncbi.nlm.nih.gov/pubmed/25734088>
  - <http://www.ncbi.nlm.nih.gov/pubmed/25339379>
  - <http://www.ncbi.nlm.nih.gov/pubmed/25391806>

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© 2021 Abbott. All rights reserved.