WHO Prequalification of In Vitro Diagnostics Programme PUBLIC REPORT

Product: Bioline Malaria Ag P.f (HRP2/pLDH)¹
Number: PQDx 0209-012-00

Bioline Malaria Ag P.f (HRP2/pLDH) with product codes **05FK90,05FK93, 05FK91 and 05FK92** manufactured by **Abbott Diagnostics Korea Inc²**, **CE marked regulatory version**, was accepted for the WHO list of prequalified diagnostics and was listed on 16 October 2015.

Summary of prequalification status for Bioline Malaria Ag P.f (HRP2/pLDH)

	Date	Outcome
Status on PQ list	16-Oct-2015	listed
Dossier assessment	14-Aug-2015	MR
Site inspection(s) of quality management system	24-Jul-2015	MR
Product performance evaluation	Round 3 (2011) and Round 6 (2015)	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
1.0-3.0	Re-introduction of a stable single-use buffer vial (for product code 05FK93).	2-May-2016
4.0	Changes to the labelling of the cassette (for product codes 05FK90, 05FK93).	17-Nov-2016

¹ Product name was changed from SD BIOLINE Malaria Ag P.f(HRP2/pLDH) and SD BIOLINE Malaria Ag P.f (HRP2/pLDH) POCT to Bioline Malaria Ag P.f (HRP2/pLDH)

² Manufacturer's name changed from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc

5.0	Addition of test kit configurations (05FK91, 05FK92) to include safety lancets.	22-Dec-2016
6.0	Changes to the manufacturer name from Standard Diagnostics, Inc to Abbott Diagnostics Korea Inc and product name from SD BIOLINE Malaria Ag P.f (HRP2/pLDH) to Bioline Malaria Ag P.f (HRP2/pLDH).	21-Aug-2020

Intended use:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f (HRP2/pLDH) test is a rapid, qualitative test for the detection of histidine-rich protein II (HRP2) antigen and lactate dehydrogenase from Malaria Plasmodium falciparum in human whole blood. Bioline Malaria Ag P.f (HRP2/pLDH) test is intended for professional use as an initial screening test only. Reactive specimens should be confirmed by a supplemental assay such as microscopic examination of thin blood smear”.*

Test principle:

According to the claim of assay description from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f (HRP2/pLDH) test contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to the histidine-rich protein II (HRP2) and the lactate dehydrogenase of P. falciparum on the test band region. This test device has a letter of "T1", "T2" (test lines) and "C" (control line). on the surface of the device. Neither the test lines nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and should always appear if the test procedure is performed correctly”.*

The test kit contains:

	25T/kit (product code 05FK90)	1T/kit x 25ea (product code 05FK93)	25T/kit (product code 05FK91)	1T/kit x 25ea (product code 05FK92)
Test cassettes individually packed in foil pouch with a desiccant	25 test devices	25x 1 test device	25 test devices	25x 1 test device
Assay diluent dispensed in plastic bottle	1 x 5ml/bottle	25x 180µl/vial	1 x 5ml/bottle	25x 180µl/vial
Specimen transfer devices disposable (5µl)	25 units of 5µl	25 packages each containing one specimen transfer device (5µl), one lancet, one alcohol swab	25 units of 5µl	25 packages each containing one specimen transfer device (5µl), one safety lancet , one alcohol swab
Lancets (optional) disposable, sterile	25 units		25 units (safety lancet)	
Alcohol swabs (optional) disposable	25 units		25 units	
Instructions for use	1 copy	1 copy	1 copy	1 copy
Summarized instructions for use	N/A	25 copies	N/A	25 copies

Storage:

The test kit should be stored at 1 – 40 °C.

Shelf-life:

24 months.

Prioritization for prequalification

Abbott Diagnostics Korea Inc. submitted an application for prequalification of Bioline Malaria Ag P.f (HRP2/pLDH). Based on the results of the WHO product testing of malaria RDTs Round 3, Bioline Malaria Ag P.f (HRP2/pLDH) was given priority for prequalification assessment.

Dossier assessment

Abbott Diagnostics Inc submitted a product dossier for Bioline Malaria Ag P.f (HRP2/pLDH) as per the “Instructions for compilation of a product dossier” (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external

experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2). Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for Bioline Malaria Ag P.f (HRP2/pLDH) for prequalification.

Commitment for prequalification:

Manufacturer will seek specimens from patients infected with *Schistosoma* to test for possible cross-reactions.

Commitment is under review.

Manufacturing site inspection

A comprehensive inspection was performed at the site of manufacture (Production: 65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, Republic of Korea 446-930 and Warehouse: 19-22, Dongtansandan 3-gil, Dongtan-myeon, Hwaseong-si, Gyeonggi-do, Republic of Korea) of Bioline Malaria Ag P.f (HRP2/pLDH) in May 2015³ as per the Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics (PQDx_014 v1). The inspection found that the manufacturer had an acceptable quality management system and good manufacturing practices in place that ensured the consistent manufacture of a product of good quality. The manufacturer's responses to the nonconformities found at the time of the inspection were accepted 24 July 2015.

Product performance evaluation

The product was evaluated in Round 3 and Round 6 of WHO/FIND product testing of RDTs for malaria antigen detection was completed in 2010 and 2015. The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild type parasite panel, and a *P. falciparum* negative panel. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

In Round 3, the following results were observed: *P. falciparum* panel detection score (87.9% at 200 parasites/ μ l), false-positive rates (2.0% for clean negatives, 0% for *P. vivax* at 200 parasites/ μ l, 0% for *P. vivax* at 2000 or 5000 parasites/ μ l) and invalid rate (0%).

In Round 6, the following results were observed: *P. falciparum* panel detection score (overall 88.0% at 200 parasites/ μ l; 87.0% for HRP2 and 52% for pLDH), false-positive rates (0.0% for clean negatives, 0.7% for *P. vivax* at 200 parasites/ μ l, 0% for *P. vivax* at 2000 or 5000 parasites/ μ l) and invalid rate (0%).

³ Previous site inspections were carried out in September 2010, March 2012, and November 2012.

Summary performance characteristics	Panel detection score	False positive rate (%)		Invalid rate (%)
	200 parasites/ μ l	200 parasites/ μ l	Clean negatives	
	Pf	Pv		
Round 3	87.9	0.0	2.0	0.0
Round 6	88.0 (87% for HRP2, 52% for pLDH)	0.7	0.0	0.0

Based on these results, the performance evaluation for Bioline Malaria Ag P.f (HRP2/pLDH) meets the WHO prequalification requirements.


Labelling

- 1. Labels**
- 2. Instructions for use**


**Note: Labelling has been changed as per site name change and product rebranding. But temporarily labelling of legacy brand (SD BIOLINE) will be used in the market according to registration status in each country.*

1. Labels


1.1 Package box for 05FK90



1.3 Package box for 05FK92


1.4 Package box for 05FK93


1.5 Device pouch for 05FK90, 05FK91, 05FK92, 05FK93

1.6 Outer Pouch for_05FK92


1.7 Outer Pouch for_05FK93

2. Instructions for use⁴

⁴ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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ENGLISH

About the test

Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is caused by a parasite that is transmitted from one human to another through an insect. There are four types of malaria that can infect humans: Plasmodium falciparum, P. vivax, P. ovale and P. malariae. In humans, the parasites (called sporozoites) migrate to the liver, where they mature and are released in another form (merozoites). The disease now occurs for most in 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million clinical deaths per year. As a result, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

Test principle

The BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test contains a membrane strip, which is pre-coated with monoclonal antibodies anti-HRP2 and anti-pLDH. The test is performed by the lateral flow methodology of P.falcioparum on the test band region. This test has a letter of "T1", "T2" (test lines) and "C" (control line), on the surface of the device. Neither the test lines nor the control line are visible in the result window. The control line appears only when the procedure is performed correctly and should always appear if the test procedure is performed correctly.

[Intended Use]
The BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test is a rapid, qualitative test for the detection of **histidine-rich protein II (HRP2) antigen** and **lactate dehydrogenase from Malaria Plasmodium falciparum** in human whole blood. **BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test** is intended for professional use only. Results should be confirmed by a supplemental assay such as microscopic examination of thin blood smear.

Materials provided and active ingredients of main component

- The BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test kit contains following items to perform the assay:
 - 25 Disposable swabs with desiccant in individual foil pack
 - 25 Disposable specimen applicator (S.p.)
 - 25 Sterile Lancets, 25 Alcohol swabs
 - 25 Suministrados swabs com dessecante em embalagem individual
 - 25 Aplicadores amostras
 - 25 Agulhas estéreis, 25 Alcoolswabs
 - 25 Suministrados instruções de uso
 - 25 Ingredientes ativos de principal componente
 - 1 Strip isopropil álcool, Gold conjugate: Mouse monoclonal antibodies specific to P.fHRP2 (0.5-10 µg), Gold conjugate: Mouse monoclonal antibodies specific to P.f pLDH (0.5-10 µg), Test line: Mouse monoclonal antibodies specific to P.fHRP2 (0.5-10 µg), Test line: Mouse monoclonal antibodies specific to P.f pLDH (0.5-10 µg), Antigen Diluent: 100 µL, 200 mM Phosphate Saline Buffer (p.s.), Tris, X-100 (g.s.), Sodium azide (g.s.)
- Active ingredients of main component
 - 1 Strip isopropil álcool, Gold conjugate: Mouse monoclonal antibodies specific to P.fHRP2 (0.5-10 µg), Gold conjugate: Mouse monoclonal antibodies specific to P.f pLDH (0.5-10 µg), Test line: Mouse monoclonal antibodies specific to P.fHRP2 (0.5-10 µg), Test line: Mouse monoclonal antibodies specific to P.f pLDH (0.5-10 µg), Antigen Diluent: 100 µL, 200 mM Phosphate Saline Buffer (p.s.), Tris, X-100 (g.s.), Sodium azide (g.s.)

Materials required but not provided

- Microscope, Protective Goggles, Gloves, Biohazard container

Kit storage and stability

- The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the test kit or its components.
- Store in a sealed container at refrigerator, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to use. Do not open the pouch or separate components until ready for use.
- The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
- Do not use test kit beyond its expiration date.
- The shelf life of the kit is indicated on the outer package.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the test kit, avoid exposure to high temperature (above 45 °C) for a period longer than 1 week.

Warnings

- Do not use diagnostic use only. Do not reuse the test device and kit components.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
- Avoid splashing or aerosol formation of specimen and assay diluent.
- Clean up spills thoroughly using appropriate disinfectant.
- Decomantate and dispose of all specimens, reaction kits and potentially contaminated materials (ie. lancet, sample applicator, test device) in a biohazard container as if they were infectious waste.
- Do not mix or interchange different lots or those for other products.
- Do not pipette by mouth.
- Do not eat the desiccant in the foil pouch.
- Do not drink assay diluent.
- The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the eyes and/or skin, washed areas with soap and water immediately. Irritation or signs of toxicity occur, seek medical attention.

Specimen collection and storage

- Collection by venipuncture**
 - Using venipuncture, collect whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
 - Invert the collection tube immediately, then it should be refrigerated at 2 - 8 °C.
 - If stored at 2 - 8 °C, the blood specimen should be tested within 3 days.
 - Using blood specimens that have been stored for more than 3 days can yield a nonspecific reaction.

[Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Wash the hands and the fingers and the lancet with a sterile lancet provided.
- Use the very first drop of blood with sterile gauze or cotton.
- Using a capillary pipette (5 µl) provided, gently squeeze the tube, immerse the capillary end in the blood drop and then gently release the pressure to draw blood into the capillary end.

Performance characteristics

- Internal validity of BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test
 - Sensitivity of BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test for the detection of P.fHRP2

No. of parasitaefol of microscopy (No. of positive)	BioInfect [®] Malaria Ag P.f (HRP2/pLDH) (No. of positive)	Sensitivity
<1-50	14	93.3%
51-100	30	100%
101-500	83	100%
501-1,000	87	100%
1,001-5,000	62	100%
>5,000	35	100%
Total	312	99.7%

- Sensitivity of BioInfect[®] Malaria Ag P.f (HRP2/pLDH) Test for the detection of pLDH specific to P.falcioparum

No. of parasitaefol of microscopy (No. of positive)	BioInfect [®] Malaria Ag P.f (HRP2/pLDH) (No. of positive)	Sensitivity
<1-50	15	8
51-100	30	29
101-500	83	83
501-1,000	87	100%
1,001-5,000	62	100%
>5,000	35	100%
Total	312	97.4%

- Specificity of BioInfect[®] Malaria Ag P.f (HRP2/pLDH) for the detection of antigens of P.falcioparum

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- Reproducibility of BioInfect[®] Malaria Ag P.f (HRP2/pLDH) test has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference plasma specimens criteria.
- Assay for the detection of malaria according to the level of parasite.
 - *Plasmodium falciparum: Sensitivity with >50 parasitaefol of blood is 100%.
 - *Plasmodium falciparum(pLDH): Sensitivity with >100 parasitaefol of blood is 100%.

- Total Sensitivity and Specificity of BioInfect[®] Malaria Ag P.f (HRP2/pLDH) Test

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ENGLISH

About the test

[Intended Use] Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is infected by a parasite that is transmitted from one human to another through the bite of infected *Anopheles* mosquitoes. There are four types of malaria that can infect humans: *Plasmodium falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and are released in another form: merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 200 million clinical cases each year. At present, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

[Test principle]

The BioLine™ Malaria Ag P.f. (HRP2/pLDH) test contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to HRP2 and pLDH. The test is performed by dipping the strip into the blood specimen and placing it on a flat, dry surface. Let the test device with *P. falciparum* on the test band region. This test device has a letter "T", "TZ" (test lines) and "C" (control line), on the surface of the device. Neither the test lines nor the control line is visible in the result window before applying a specimen. The control is used for procedural control and should always appear as the test procedure is performed correctly.

[Intended Use]

The BioLine™ Malaria Ag P.f. (HRP2/pLDH) test is a rapid, qualitative test for the detection of **histidine-rich protein II (HRP2) antigen and lactate dehydrogenase from *Plasmodium falciparum*** in human whole blood. BioLine™ Malaria Ag P.f. (HRP2/pLDH) test is intended for use in a laboratory setting such as microscopic examination of thin blood smears.

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- The test kit: Malaria Ag P.f. (HRP2/pLDH) test kit contains following items to perform the assay:
 - 25 Test devices with desiccant in individual foil pouches
 - Assay diluent (1 x 50 mL)
 - 25 Disposable specimen applicators (5 µl), 25 Safety Lancets, 25 Alcohol swabs
 - Instructions for use

- Active ingredients of main component
 - 1 test strip includes: Gold conjugate; Mouse monoclonal antibodies specific to P.f. HRP2 gold colloid (0.010,0.02 µg), Mouse monoclonal antibodies specific to P.f. pLDH gold colloid (0.010,0.02 µg), Test line 1: Mouse monoclonal antibodies specific to P.f. HRP2 (0.5,0.1 µg), Test line 2: Mouse monoclonal antibodies specific to P.f. pLDH (0.5,0.1 µg), Control line: Goat anti- mouse IgG (0.1,0.2 µg)
 - Assay diluent includes: Bovine serum albumin (q.s.), 200 mM Phosphate Saline Buffer (q.s.), Triton X-100 (q.s.), Sodium azide (q.s.)

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- Microscope, Protective gloves, Timer, Biohazard container

Kit storage and stability

- The test kit should be stored at a temperature between 1°C and 40°C. Do not freeze the kit or its components. Assay diluent cap should be kept firmly sealed between each use.
- When the test kit is opened, the test kit should be stored in the original seal in a temperature (5 - 40°C) minimum 30 mins prior to the test. Do not open the pouch until components come to room temperature.
- When the test kit is opened and resealed for each assay, it is stable until expiration date if kept at 1 - 4°C.
- The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the pouch.
- Do not use test kit beyond its expiration date.
- The shelf life of the kit is indicated on the outer package.
- Do not mix reagent of different lots or those for other products.
- When transporting or storing the kit, avoid exposure to high temperature (above 45°C) for a period longer than 1 week.

Warnings

- For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
- Do not use expired or aerosolized or contaminated reagents and assay kits.
- Clean up spillage thoroughly using an appropriate disinfectant.
- Decomaminate and dispose of all specimens, reaction kits and potentially contaminated materials (e.g. lancets, specimen applicator, test device) in a biohazard container as if they were infectious waste.
- Do not mix or interchange different specimens.
- Do not mix reagent of different lots or those for other products.
- Do not pipette by mouth.
- Do not eat the desiccant in the foil pouch.
- Do not drink assay diluent.
- The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user's normal laboratory safety practices are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

Specimen collection and storage

[Collection by venipuncture] Using venipuncture, collect whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate). Use the product if used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

[Collection using a lancet] Clean the area to be lanced with an alcohol swab. Squeeze the end of the fingertip and place a safety lancet over it. Wipe away the first drop of blood with sterile gauze or cotton. Use the capillary pipette provided, gently squeezing the tube, immerse the open end in the blood drop and then gently release the pressure to draw blood into the pipette.

Performance characteristics

- Intend of evaluation of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test
 - Sensitivity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test for the detection of P.f. HRP2

No. of parasitized of blood (no. of positive)	Microscopy (HRP2/pLDH) (no. of positive)	BioLine™ Malaria Ag P.f. (HRP2/pLDH) (no. of positive)	Sensitivity (%)
1-50	15	14	93.3%
51-100	30	30	100%
101-500	83	83	100%
501-1,000	87	87	100%
1,001-5,000	62	62	100%
>5,000	35	35	100%
Total	312	311	99.7%

- Sensitivity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test for the detection of pLDH specific to P.f. (HRP2/pLDH) test
 - Sensitivity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test for the detection of P.f. HRP2

No. of parasitized of blood (no. of positive)	Microscopy (HRP2/pLDH) (no. of positive)	BioLine™ Malaria Ag P.f. (HRP2/pLDH) (no. of positive)	Sensitivity (%)
1-50	15	8	53.3%
51-100	30	29	96.7%
101-500	83	83	100%
501-1,000	87	87	100%
1,001-5,000	62	62	100%
>5,000	35	35	100%
Total	312	304	97.4%

- Specificity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) Test for the detection of P.f. (HRP2/pLDH) test
 - Specificity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test for the detection of HRP2 of P.f. (HRP2/pLDH) test

No. Type	Potential cross-reacting pathogens	No. Type	Potential cross-reacting pathogens	Specificity (%)
1	Dengue	11	Yellow fever	100%
2	Chikungunya	12	Chikungunya	100%
3	Chikungunya	13	Chikungunya	100%
4	Chikungunya	14	Chikungunya	100%
5	Yersinia	15	Yersinia	100%
6	Hepatitis B	16	Hepatitis B	100%
7	HIV 1 (subtype B)	17	HIV 1 (subtype B)	100%
8	Hepatitis B	18	Hepatitis B	100%
9	Hepatitis C	19	Hepatitis C	100%
10	Hepatitis B	20	Hepatitis B	100%
11	Hepatitis C	21	Hepatitis C	100%
12	Hepatitis B	22	Hepatitis B	100%
13	Hepatitis C	23	Hepatitis C	100%

- Specificity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) Test for the detection of P.f. (HRP2/pLDH) test
 - Specificity of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test for the detection of HRP2 of P.f. (HRP2/pLDH) test

No. Type	Potential cross-reacting pathogens	No. Type	Potential cross-reacting pathogens	Specificity (%)
1	Dengue	11	Yellow fever	100%
2	Chikungunya	12	Chikungunya	100%
3	Chikungunya	13	Chikungunya	100%
4	Chikungunya	14	Chikungunya	100%
5	Yersinia	15	Yersinia	100%
6	Hepatitis B	16	Hepatitis B	100%
7	HIV 1 (subtype B)	17	HIV 1 (subtype B)	100%
8	Hepatitis B	18	Hepatitis B	100%
9	Hepatitis C	19	Hepatitis C	100%
10	Hepatitis B	20	Hepatitis B	100%
11	Hepatitis C	21	Hepatitis C	100%
12	Hepatitis B	22	Hepatitis B	100%
13	Hepatitis C	23	Hepatitis C	100%

- Reproducibility of BioLine™ Malaria Ag P.f. (HRP2/pLDH) test has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panel acceptance criteria.

- Assay for the detection of Malaria according to the label of parasite.
 - *Plasmodium falciparum: Sensitivity with ≥50 parasite/ul of blood is 100%.
 - *Plasmodium falciparum/pLDH: Sensitivity with ≥100 parasite/ul of blood is 100%.

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

5 **EN** Dispense 5 µl of drawn blood into round specimen well.

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
4 **EN** Using a disposable inverted cup provided, dip the circular end of a inverted cup into the blood specimen.

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
5 **EN** Dispense 5 µl of drawn blood into round specimen well touching pad.

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DISPENSE
VERSER
DISPENSE
COLOQUE

EN Use Inverted cup : Let the circular end of the inverted cup touch the pad, then press down lightly.

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7 **EN** Interpret test results 15 minutes (up to 30 minutes) after adding assay diluent. Reading outside of this time frame (before 15 min or after 30 min) may provide false results.

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4 **EN** Open the foil pouch and look for the following:

1. Result window
 2. Specimen well
 3. Assay diluent well
- Then, label the device with the patient identifier.**

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


POSITIVE / POSITIF / POSTIVO

EN **⚠ Caution:** The presence of any line, no matter how faint, the result is considered positive.

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



EN If the control line "C" is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

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