WHO Prequalification of In Vitro Diagnostics

PUBLIC REPORT

Product: Diagnostic kit for HIV (1+2) antibody (colloidal gold) V2

WHO reference number: PQDx 0267-037-00

Diagnostic kit for HIV (1+2) antibody (colloidal gold) V2 with product codes **R-401-50-C-2, KH-R-02, A-GOLD-01, KH-R-03 and R-401-50-C-3** manufactured by **Shanghai Kehua Bio-engineering Co., Ltd, rest of world** regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 21 December 2016.

Summary of WHO prequalification assessment for Diagnostic kit for

HIV (1+2) antibody (colloidal gold) V2

	Date	Outcome
Prequalification listing	21 December 2016	Listed
Dossier assessment	3 June 2016	MR
Site inspection(s) of quality management system	24-26 June 2019	MR
Product performance evaluation and operational characteristics	22 November 2016	MR

MR: Meets requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
1.0-3.0	Changes made on the draft public reports until it was finalized as version 3 that was published December 2016.	21 December 2016
4.0	1. Addition of two package specifications. Each specification contains two bottles of diluent as required by customers. No change with raw materials. No change with test procedure or kit performance. 2. Extension of expiry date/shelf life of reagents from 18 months to 24 months, not including accessories."	29 June 2022

Intended use

According to the claim of intended use from Shanghai Kehua Bio-engineering Co., Ltd, *“Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 is an in vitro, visually read test for the qualitative determination of antibodies against HIV-1 and HIV-2 in human serum, plasma, venous and capillary whole blood. This product is intended for use as an aid to detect antibodies against HIV (1+2) from individuals suspected to be infected by HIV. It is intended to be used by laboratory professional and health care workers at point-of-care.”*

Assay description

According to the claim of assay description from Shanghai Kehua Bio-engineering Co., Ltd, *“Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 adopts the solid phase colloidal gold immunochromatographic technology for the qualitative detection of antibodies against HIV-1/HIV-2. The gold HIV-1 gp160 conjugate and the gold HIV-2 gp36 conjugate are coated to the conjugate pad in advance. The test line (recombinant gp36 and gp41 antigens) and the control line (monoclonal antibody against gp160) are pre-coated on the surface of Nitrocellulose (NC) membrane. When the sample is added to the sample pad, it migrates through the conjugate pad, reconstitutes and mixes with the colloidal gold-antigen complex and moves on till the end of the strip. A reddish-purple test line will be visible in the strip if there are enough antibodies against HIV-1/HIV-2 in the sample. If the specific antibodies are absent, or present at a very low level, no test line appears.”*

Test kit contents:

Component	50 tests (product code R-401-50-C-2)	50 tests (product code KH-R-02)	Accessories kit (product code A-GOLD-01)	50 tests (product code KH-R-03)	50 tests (product code R-401-50-C-3)
HIV (1+2) antibody test cassette	50	50	N/A	50	50
Sample diluent	4ml x 1 bottle	4ml x 1 bottle	N/A	4ml x 2 bottles	4ml x 2 bottles
Safety lancet	50	N/A	50	N/A	50
Alcohol swab	50	N/A	50	N/A	50
Disposable transfer pipette	50	N/A	50	N/A	50
Instructions for use	1	1	1	1	1

N/A: Not Applicable

Items required but not provided

Item
Consumables: Waste container, gloves. If serum/plasma procedure - blood collection tubes, pipette, pipette tips
Equipment: Timer

Storage

The test kit should be stored at 4-30°C.

Shelf-life upon manufacture

24 months

Warnings/limitations

Refer to current version of manufacturer's instructions for use.

Prioritization for prequalification

Based on the established prioritization criteria, Diagnostic kit for HIV (1+2) antibody (colloidal gold) V2 was given priority for WHO prequalification assessment.

Dossier assessment

Shanghai Kehua Bio-engineering Co., Ltd submitted a product dossier for Diagnostic kit for HIV (1+2) antibody (colloidal gold) V2 as per the "Instructions for compilation of a product dossier" (PQDx_018 v1). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 03 June 2016.

Based on the product dossier screening and assessment findings, the product dossier for Diagnostic kit for HIV (1+2) antibody (colloidal gold) V2 meets WHO prequalification requirements.

Manufacturing site inspection

An inspection of Shanghai Kehua Bio-engineering Co., Ltd located at 1189 North Qinzhou Road 200233 Shanghai People Republic of China and at the warehouse located at 388 Minfa Road, Shanghai People Republic of China was conducted from 24 to 26 of June 2019. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality. Routine inspections of the Manufacturing site will be conducted with

copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly.

Note that a WHOPIR reflects the information on the most current inspection performed at a manufacturing site for *in vitro* diagnostic products and gives a summary of the inspection findings.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 15 of November 2019.

Product performance evaluation

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 was evaluated by WHO in the last quarter of 2016 using serum/plasma specimens. From this evaluation, we drew the following conclusions:

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 is an immunochromatographic test for the qualitative detection of HIV-1/2 antibodies in human serum, plasma or whole blood specimens. A volume of 40 µL of specimen is needed to perform the assay. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities and non-laboratory settings. Reading of the results can be done visually i.e. subjectively read.

Performance characteristics in comparison with an agreed reference standard		
	Initial (95% CI)	Final (95% CI)
Sensitivity %	100% (99.2% - 100%)	100% (99.2% - 100%)
Specificity %	100% (99.4% - 100%)	100% (99.4% - 100%)
Invalid rate %	0%	
Inter-reader variability %	0.09%	

Additional performance characteristics	
Sensitivity during seroconversion on eight seroconversion panels in comparison with a benchmark assay; Enzygnost Anti-HIV 1/2 Plus (Siemens Healthcare Diagnostics)	Seroconversion sensitivity index of - 0.25, therefore detection is 0.25 days earlier than the benchmark assay.
Analytical sensitivity on a mixed titer panel in comparison with an agreed reference standard	25 of 25 specimens were correctly classified.
Lot to lot variation on a dilution panel in comparison with an agreed reference standard	Acceptable.

Key operational characteristics	
Validated specimen types	Serum, plasma (EDTA, heparin or sodium citrate), venous whole blood, capillary whole blood.
Number of steps	2 with precision required.
Time to result	16minutes.
Endpoint stability	15-25 minutes.
Internal QC	Yes, the control band appears when only buffer is added to the device (without specimen) and is therefore a reagent addition control.
In-use stability of reagents	The sample diluent should be used within 3 months after opening.

Labelling

- 1. Labels**
- 2. Instructions for use**

1. Labels

Labels of test cassette


Figure 1 Label of aluminum foil pouch for cassette (front side) (120 x 62 mm)


Figure 2 Label of aluminum foil pouch for cassette (back side) (120 x 62 mm)

Label of sample diluent


Figure 3 Label of sample diluent (54 x 22 mm)

Kit box and label for product code R-401-50-C-2

IVD

**Diagnostic Kit for HIV (1+2) Antibody
(Colloidal Gold)**

REF R-401-50-C-2

KHB®

KHB®

REF R-401-50-C-2

**Diagnostic Kit for HIV (1+2) Antibody
(Colloidal Gold)**

IVD

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-64950625
+86-8008203370
Fax: +86-21-64950975
Website: <http://www.skhb.com>

Please refer instructions for use

IVD

KHB®

REF R-401-50-C-2

**Diagnostic Kit for HIV (1+2) Antibody
(Colloidal Gold)**

CONTENT

- 1. HIV (1+2) Antibody test cassette x 50
- 2. One bottle containing 4ml sample diluent
- 3. Safety lancet x 50
- 4. Alcohol pad x 50
- 5. Disposable transfer pipette x 50
- 6. Instructions for use x 1


/ Production date
(YYYY/MM/DD):

/ Lot number:

/ Expiry date
(YYYY/MM/DD):

/ Store at 4°C-30°C

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-64950625
+86-8008203370
Fax: +86-21-64950975
Website: <http://www.skhb.com>

Please refer instructions for use

IVD

Kit box and label for product code KH-R-02


IVD

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

REF KH-R-02

KHB®

KHB®

REF KH-R-02

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

IVD

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-64950625
+86-8008203370
Fax: +86-21-64950975
Website: <http://www.skhb.com>

Please refer instructions for use

IVD

KHB®

REF KH-R-02

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

CONTENT

- 1. HIV (1+2) Antibody test cassettes x 50
- 2. One bottle containing 4ml sample diluent
- 3. Instructions for use x 1


/ Production date (YYYY/MM/DD):

/ Lot number:

/ Expiry date (YYYY/MM/DD):

/ Store at 4°C-30°C

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-64950625
+86-8008203370
Fax: +86-21-64950975
Website: <http://www.skhb.com>

Please refer instructions for use

IVD

Kit box and label for product code A-GOLD-01

CONTENT
 1. Safety lancet x 50
 2. Alcohol pad x 50
 3. Disposable transfer pipette x 50

Accessories kit

REF A-GOLD-01

KHB

KHB

REF A-GOLD-01

Accessories kit

 Shanghai Kehua Bio-engineering Co., Ltd.
 1189 North Qinzhou Road, 200233, Shanghai, P.R. China
 Tel: +86-21-64950625
 +86-8008203370
 Fax: +86-21-64950975
 Website: <http://www.skhb.com>

 **Please refer instructions for use**  

KHB

REF A-GOLD-01

Accessories kit

LOT / Lot number:

 / Expiry date
 (YYYY/MM/DD):

 Shanghai Kehua Bio-engineering Co., Ltd.
 1189 North Qinzhou Road, 200233, Shanghai, P.R. China
 Tel: +86-21-64950625
 +86-8008203370
 Fax: +86-21-64950975
 Website: <http://www.skhb.com>

 **Please refer instructions for use**  

Kit box and label for product code for KH-R-03

IVD

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

REF KH-R-03

KHB

KHB

REF KH-R-03

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

IVD

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-54500744
+86-8008203370
Fax: +86-21-64854051
Website: <http://www.skhb.com>

Please refer instructions for use

IVD

KHB

REF KH-R-03

Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

CONTENT

- 1. HIV (1+2) Antibody test cassette x 50
- 2. Sample diluent (4mL/bottle) x 2
- 3. Instructions for use x 1


/ Production date (YYYY/MM/DD):

/ Lot number:

/ Expiry date (YYYY/MM/DD):

/ Store at 4°C-30°C

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-54500744
+86-8008203370
Fax: +86-21-64854051
Website: <http://www.skhb.com>

Please refer instructions for use

IVD

Kit box and label R-401-50-C-3

IVD

**Diagnostic Kit for HIV (1+2) Antibody
(Colloidal Gold) V2**

REF R-401-50-C-3

KHB

KHB

REF R-401-50-C-3

**Diagnostic Kit for HIV (1+2) Antibody
(Colloidal Gold) V2**

IVD

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-54500744
+86-8008203370
Fax: +86-21-64854051
Website:<http://www.skhb.com>

Please refer instructions for use

IVD

KHB

REF R-401-50-C-3

**Diagnostic Kit for HIV (1+2) Antibody
(Colloidal Gold) V2**

- CONTENT**
1. HIV (1+2) Antibody test cassette x 50
 2. Sample diluent (4mL/bottle) x 2
 3. Safety lancet x 50
 4. Alcohol pad x 50
 5. Disposable transfer pipette x 50
 6. Instructions for use x 1


/ Production date
(YYYY/MM/DD):

/ Lot number:

/ Expiry date
(YYYY/MM/DD):

/ Store at 4°C-30°C

Shanghai Kehua Bio-engineering Co., Ltd.
1189 North Qinzhou Road, 200233, Shanghai, P.R. China
Tel: +86-21-54500744
+86-8008203370
Fax: +86-21-64854051
Website:<http://www.skhb.com>

Please refer instructions for use

IVD

2. Instructions for use¹

¹ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

Accessories kit

Instructions for use

KHB® Shanghai Kehua Bio-engineering Co., Ltd.

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Accessories kit

REF A-GOLD-01 (50 Sets)

△ PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE ALONG WITH KHB DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2 (REF KH-R-02)

Intended Use

The accessories kit is intended to be used for collecting fingerstick whole blood specimen from individuals and transferring the specimen collected for further test.

Materials Provided

- Safety lancet 50 pieces
- Alcohol pad (70% Isopropyl alcohol) 50 pieces
- Disposable transfer pipette 50 pieces
- Instructions for use 1 piece

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- Protective gloves
- Biohazard waste container

△ Warnings and Precautions

The accessories kit is for specimen collection and transfer only. The precautions are included, but not limited to the following:

- Wear gloves during the entire process.
- Do not use the accessories if seal is broken or package is damaged.
- All the accessories are for single use. Do not reuse the accessories.
- The accessories are only used for specimen collection, do not store the specimen in the pipette.
- Clean and disinfect all the areas that may be contaminated by spills of specimens with appropriate disinfectant.
- Decontaminate and dispose of all the accessories and other potentially contaminated materials as if they were infectious wastes, in a biohazard container.
- Application of a sample volume of 40uL has been validated for use with KHB Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 (REF KH-R-02). This volume is attained when the specimen reaches the second graduation mark on the pipette. Use of any other volume has not been validated and may produce incorrect results.

Procedure

1. Choose the ring finger of the patient.
2. Open the alcohol pad and clean the fingertip. Allow the fingertip dry. (Figure 1)
3. Remove the cap of the safety lancet. (Figure 2)
4. Place the lancet on the fingertip, press the lancet against the finger and puncture the skin. (Figure 3)
5. Discard the lancet into a bio-hazard waste container immediately.

Figure 1

Figure 2

Figure 3


Figure 4

6. Squeeze the fingertip to form a bulb and use the disposable transfer pipette to collect the drop of the blood. (Figure 4)
7. Allow the blood to reach the second graduation and the specimen volume will be 40ul. Avoid formation of bubbles during the collection of specimen. (Figure 5)
8. Transfer the specimen to the device immediately after collection. Avoid formation of bubbles.
9. Discard the disposable transfer pipette into a bio-hazard waste container immediately.


Figure 5

KHB® Accessories kit

EN

Key to symbols used

	Caution		Sufficient for "n" use
	Batch code		Manufacturer
	Do not reuse		Consult instruction for use
	Catalogue number or order number		Do not use if package is damaged
	Use by		

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1189 North Qinzhou Road, 200233, Shanghai, P.R. China

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DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2

REF KH-R-02 (50 Tests)

PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE

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Materials Provided

- | | | |
|---|------------------|--|
| • HIV (1+2) antibody test cassette with desiccant | 1 × 50 cassettes | |
| • Sample diluent | 4ml × 1 bottle | |
| • Instruction for use | 1 piece | |

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

1. Take out a test cassette from the foil pouch, and place it on a horizontal surface.
2. Use the specimen of either serum/plasma or whole blood: Add 40 μ l of specimen to the specimen area firstly (Figure 1), then slowly add 1 drop of sample diluent vertically to the same area (Figure 2). Please avoid direct contact of the diluent bottle with the specimen area.
3. Incubate the cassette at room temperature (18-28 $^{\circ}$ C) and read the result after 15 minutes but not more than 25 minutes. Reading the test before 15 minutes or after 25 minutes may cause false result.

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Table 4

HIV Subtype	No. of specimens tested	Positive on Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2
-------------	-------------------------	---

KHB [®]	DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2	For <i>in vitro</i> diagnostic use	EN
A	1	1	
A1	1	1	
A2	1	1	
C	3	3	
CRF01_AE	3	3	
CRF02_AG	3	3	
CRF06_cpx	3	3	
CRF36_cpx	3	3	
D	3	3	
F1	2	2	
F2	1	1	
G	3	3	
H	3	3	
J	3	3	
K	3	3	
Group O	3	3	
Recombinant Form (recombination of subtype B, D, CRF19 and CRF21_A2D)	1	1	
Total	40	40	

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Références

1. **Gallo, R. C., S. Z. Salahuddin, M. Popovic, G. M. Shearer, M. Kaplan, B. F. Haynes, T. J. Palker, R. Redfield, J. Oleske, B. Safai, and et al.** 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* **224**:500-503
2. **Viscidi, R., E. Ellerbeck, L. Garrison, K. Midthun, M. L. Clements, B. Clayman, B. Fergie, and G. Smith.** 1990. Characterization of serum antibody responses to recombinant HIV-1 gp160 vaccine by enzyme immunoassay. NIAID AIDS Vaccine Clinical Trials Network. *AIDS research and human retroviruses* **6**:1251-1256.
3. **Chen, Y. H., Y. Xiao, W. Wu, Q. Wang, G. Luo, and M. P. Dierich.** 2000. HIV-2 transmembrane protein gp36 binds to the putative cellular receptor proteins P45 and P62. *Immunobiology* **201**:317-322..
4. **Brandt, J., L. O. Andersson, and J. Porath.** 1975. Covalent attachment of proteins to polysaccharide carriers by means of benzoquinone. *Biochimica et biophysica acta* **386**:196-202.
5. **Osikowicz, G., M. Beggs, P. Brookhart, D. Caplan, S. Ching, P. Eck, J. Gordon, R. Richerson, S. Sampedro, D. Stimpson, and et al.** 1990. One-step chromatographic immunoassay for qualitative determination of choriogonadotropin in urine. *Clinical chemistry* **36**:1586.
6. **Spielberg, F., C. M. Kabeya, R. W. Ryder, N. K. Kifuani, J. Harris, T. R. Bender, W. L. Heyward, and T. C. Quinn.** 1989. Field testing and comparative evaluation of rapid, visually read screening assays for antibody to human immunodeficiency virus. *Lancet* **1**:580-584.

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DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2

REF KH-R-03 (50 Tests)

 **PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE**

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
Materials Provided

- | | | |
|---|------------------|---|
| • HIV (1+2) antibody test cassette with desiccant | 1 × 50 cassettes |  |
| • Sample diluent | 4mL × 2 bottles |  |
| • Instruction for use | 1 piece | |

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- Specimen transfer device (graduated)
- Tips (Pipette)

- Capillary tubes (graduated)
- Sterile lancet
- Protective gloves
- Specimen and test waste container
- Timer
- Alcohol pad (70% Isopropyl alcohol)
- Collection tube (with or without anticoagulant, see Specimen Collection and Storage Requirements)

Accessories kit (REF: A-GOLD-01) which contains alcohol pad, safety lancet and disposable transfer pipette can be used with the kit.

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[Safety Precautions]

- Wear gloves during the entire testing process.
- Do not pipette by mouth.
- Do not eat or smoke while handling specimens.
- Clean and disinfect all the areas that may be contaminated by spills of specimens or reagents with appropriate disinfectant.
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials as if they were infectious wastes, in a biohazard container.

[Handling Precautions]

- Do not mix and interchange different specimens.
- The sample diluent is specific for each lot of reagent. Do not use diluents from different lots.
- Do not perform the test under environment which leads to rapid evaporation (e.g. high temperature (> 45°C) and low humidity, fast air flow).
- Use of hemolytic specimens, rheumatoid factors-contained specimens and lipidemic, icteric specimens may cause invalid results or false results.
- Pooled specimens or specimens other than specified (e.g. saliva, urine) are strictly forbidden to be tested.
- Ensure that the specified volume of specimen is added by either precision pipette or specimen transfer pipette.

[Warning]

- For *in vitro* diagnosis of the infection of HIV only, this assay will not provide a diagnosis for any other disease.
- The cassette, individually enclosed in a foil pouch and packed 50 tests per kit box, is only for "single use".
- Do not use the test cassette if the pouch is damaged or the seal is broken. Contact the supplier when it happens as this is an important indicator for quality and performance of the IVD.
- The sample diluent in this kit contains Sodium Azide (<0.1%) as the preservative, be cautious while using with the kit. In case of contact with sample diluent to eyes and/or

skin, flush affected area with copious amount of water. Do not dispose the sample diluent directly in drains.

- Care should be taken to avoid the contact between the specimen area of the cassette and the diluent bottle to avoid the contamination of the diluent.

Storage Requirements

- The test cassettes and the sample diluent must be stored at 4-30°C until expiration date.
- Do not use expired devices.
- Test cassette should be used immediately upon opening the pouch.
- Sample diluent should be used within 3 months of opening.

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Whole blood specimens should be brought to room temperature (18-28°C) prior to use. Haemolytic or turbid whole blood specimens should not be used.

- Collection by fingerstick⁽⁶⁾

Squeeze the end of the fingertip and pierce with a sterile lancet.

Draw 40µl of blood with a capillary tube (graduated).

Whole blood specimens collected by fingerstick should be tested immediately.

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- Serum

Draw the whole blood into the collection tube by venipuncture, leave to settle for 2 hours for blood coagulation and then centrifuge blood to get serum supernatant.

Serum and plasma specimens may be stored at 2-8°C for up to 7 days from time of draw, store at -18°C or below for long time storage. Several freeze-thaw cycles should be avoided.

Plasma or Serum specimens should be brought to room temperature prior to use.


Plasma or serum specimens with visible bio-contamination should not be used.

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
1. Take out a test cassette from the foil pouch, and place it on a horizontal surface.
2. Use the specimen of either serum/plasma or whole blood: Add 40 μ l of specimen to the specimen area firstly (Figure 1), then slowly add 1 drop of sample diluent vertically to the same area (Figure 2). Please avoid direct contact of the diluent bottle with the specimen area.
3. Incubate the cassette at room temperature (18-28°C) and read the result after 15 minutes but not more than 25 minutes. Reading the test before 15 minutes or after 25 minutes may cause false result.

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


Limitations of the Procedure

- The kit is designed to detect antibodies against HIV-1 and HIV-2 in human serum, plasma, and whole blood. Specimens other than the specified ones may not supply accurate results.

- The intensity of the T-line does not necessarily correlate to the titer of antibody in the specimen.
- The presence of the control line (C-line) only means the flow of the conjugate.
- When specimens contain high concentration of antibody to HIV-1 are tested on the device, the control line could be absent. In this case, please perform further analysis according to the instructions of Figure 6 in Interpretation of Results.
- As this product is intended to detect antibodies against HIV (1+2) from individuals, clinical diagnosis on AIDS should not be made only based on the test results.
- A non-reactive result should not exclude the possibility of infection caused by HIV-1 or HIV-2. A false negative result can occur in the following circumstances:
 - Negative result may not have detected recently acquired HIV infection.
 - Low levels of antibody (e.g., early seroconversion samples) are below the detection limit of the test.
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration, in exceptional cases this may lead to observation of false negative results.
 - Specimens are not properly stored.
 - False negative result due to very high concentration of antibody(ies) to HIV (which is known as 'hook effect').
- For these reasons, care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.
- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- The product is not validated on specimens from infants, children, patient on anti-retroviral treatment.
- **The product is not specifically designed for detection of antibody to HIV-1 group O.**

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Table 1

Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Specificity
		Negative	Positive	
Blood donor	2500	2500	0	100%
Whole blood donation	500	500	0	100%
Hospitalized patient specimens	200	200	0	100%
Pregnant Women	200	200	0	100%
Total	3400	3400	0	100%

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Table 2

Potentially cross-reacting specimens	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2	
		Negative	Positive
Rheumatoid factor	10	10	0
Cytomegalovirus	5	5	0
Epstein Barr Virus	13	11	2
Malaria	85	85	0
Syphilis	52	51	1
Herpes	5	5	0
HBsAg	124	124	0
Anti-HCV	115	114	1
Anti-HBc	15	15	0
Anti-HTLV1/2	15	15	0
Anti-HEV	5	5	0
Specimens from STI clinics	10	10	0
<i>Trypanosomiasis</i>	24	24	0
<i>E Coli</i> infected	10	10	0
Total	488	484	4

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Table 4

HIV Subtype	No. of specimens tested	Positive on Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2

KHB [®]	DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2	For <i>in vitro</i> diagnostic use	EN
A	1	1	
A1	1	1	
A2	1	1	
C	3	3	
CRF01_AE	3	3	
CRF02_AG	3	3	
CRF06_cpx	3	3	
CRF36_cpx	3	3	
D	3	3	
F1	2	2	
F2	1	1	
G	3	3	
H	3	3	
J	3	3	
K	3	3	
Group O	3	3	
Recombinant Form (recombination of subtype B, D, CRF19 and CRF21_A2D)	1	1	
Total	40	40	

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Shanghai Kehua Bio-engineering Co., Ltd.

1189 North Qinzhou Road, 200233, Shanghai, P.R.China

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Références

1. Gallo, R. C., S. Z. Salahuddin, M. Popovic, G. M. Shearer, M. Kaplan, B. F. Haynes, T. J. Palker, R. Redfield, J. Oleske, B. Safai, and et al. 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* **224**:500-503
2. Viscidi, R., E. Ellerbeck, L. Garrison, K. Midthun, M. L. Clements, B. Clayman, B. Fergie, and G. Smith. 1990. Characterization of serum antibody responses to recombinant HIV-1 gp160 vaccine by enzyme immunoassay. NIAID AIDS Vaccine Clinical Trials Network. *AIDS research and human retroviruses* **6**:1251-1256.
3. Chen, Y. H., Y. Xiao, W. Wu, Q. Wang, G. Luo, and M. P. Dierich. 2000. HIV-2 transmembrane protein gp36 binds to the putative cellular receptor proteins P45 and P62. *Immunobiology* **201**:317-322..
4. Brandt, J., L. O. Andersson, and J. Porath. 1975. Covalent attachment of proteins to polysaccharide carriers by means of benzoquinone. *Biochimica et biophysica acta* **386**:196-202.
5. Osikowicz, G., M. Beggs, P. Brookhart, D. Caplan, S. Ching, P. Eck, J. Gordon, R. Richerson, S. Sampedro, D. Stimpson, and et al. 1990. One-step chromatographic immunoassay for qualitative determination of choriogonadotropin in urine. *Clinical chemistry* **36**:1586.
6. Spielberg, F., C. M. Kabeya, R. W. Ryder, N. K. Kifuani, J. Harris, T. R. Bender, W. L. Heyward, and T. C. Quinn. 1989. Field testing and comparative evaluation of rapid, visually read screening assays for antibody to human immunodeficiency virus. *Lancet* **1**:580-584.

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Principle of the Procedure

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB adopts the solid phase colloidal gold immunochromatographic technology for the qualitative detection of antibodies against HIV-1/HIV-2. The gold HIV-1 gp160⁽²⁾ conjugate and the gold HIV-2 gp36⁽³⁾ conjugate are coated to the conjugate pad in advance. The test line (recombinant gp36 and gp41 antigens) and the control line (monoclonal antibody against gp160) are pre-coated on the surface of Nitrocellulose (NC) membrane. When the specimen is added to the sample pad, it migrates through the conjugate pad, reconstitutes and mixes with the colloidal gold-antigen complex and moves on till the end of the strip. A reddish-purple test line will be visible in the strip if there are enough antibodies against HIV-1/HIV-2 in the specimen. If the specific antibodies are absent, or present at a very low level, no test line appears^(4, 5).

Materials Provided

- | | |
|---|------------------|
| • HIV (1+2) antibody test cassette with desiccant | 1 × 50 cassettes |
| • Sample diluent | 4ml ×1 bottle |
| • Safety lancet | 1 × 50 pieces |
| • Alcohol pad (70% Isopropyl alcohol) | 1 × 50 pieces |
| • Disposable transfer pipette | 1 × 50 pieces |
| • Instruction for use | 1 piece |

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
Figure c


Figure d

1. Choose the ring finger of the patient.
2. Open the alcohol pad and clean the fingertip (Figure a). Allow the fingertip dry.
3. Remove the cap of the safety lancet. (Figure b)
4. Place the lancet on the fingertip, press the lancet against the finger and puncture the skin. (Figure c)
5. Discard the lancet into a bio-hazard waste container immediately.
6. Squeeze the fingertip to form a bulb and use the disposable transfer pipette to collect the drop of the blood. (Figure d)
7. Allow the blood to reach the second graduation and the specimen volume will be 40ul. (Figure e) Avoid formation of bubbles during the collection of specimen.
8. Transfer the specimen to the sample well immediately after collection. Avoid formation of bubbles.
9. Discard the disposable transfer pipette into a bio-hazard waste container immediately.

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

Performance Characteristics

- Specificity

A total of 3400 specimens including serum/plasma and whole blood were tested by KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 1)

Table 1

Population	No. of specimens tested	KHB® Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Specificity
		Negative	Positive	
Blood donor	2500	2500	0	100%
Whole blood donation	500	500	0	100%
Hospitalized patient specimens	200	200	0	100%
Pregnant Women	200	200	0	100%
Total	3400	3400		100%

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DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2

REF R-401-50-C-3 (50 Tests)

⚠ PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE

Intended Use

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB is an *in vitro*, visually read test for the qualitative determination of antibodies against HIV-1 and HIV-2 in human serum, plasma, venous and capillary whole blood. This product is intended for use as an aid to detect antibodies against HIV (1+2) from individuals suspected to be infected by HIV.

It is intended to be used by laboratory professional and health care workers at point-of-care.

Summary and Explanation of the Test

AIDS (Acquired Immunodeficiency Syndrome) is an immunosuppressive disease characterized by depletion of the helper T-cell, which leaves the infected individual susceptible to opportunistic infections and some malignancies. Viruses that have been detected to cause AIDS are HIV-1, first isolated in 1983⁽¹⁾ and HIV-2, a distinctive but related type isolated in 1985. Presence of the virus elicits the production of the antibodies specific to either HIV-1 or HIV-2, and this test is designed to detect the antibodies to indicate the presence of HIV infection.

Principle of the Procedure

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB adopts the solid phase colloidal gold immunochromatographic technology for the qualitative detection of antibodies against HIV-1/HIV-2. The gold HIV-1 gp160⁽²⁾ conjugate and the gold HIV-2 gp36⁽³⁾ conjugate are coated to the conjugate pad in advance. The test line (recombinant gp36 and gp41 antigens) and the control line (monoclonal antibody against gp160) are pre-coated on the surface of Nitrocellulose (NC) membrane. When the specimen is added to the sample pad, it migrates through the conjugate pad, reconstitutes and mixes with the colloidal gold-antigen complex and moves on till the end of the strip. A reddish-purple test line will be visible in the strip if there are enough antibodies against HIV-1/HIV-2 in the specimen. If the specific antibodies are absent, or present at a very low level, no test line appears^(4,5).

Materials Provided

- | | |
|---|------------------|
| • HIV (1+2) antibody test cassette with desiccant | 1 × 50 cassettes |
| • Sample diluent | 4mL × 2 bottles |
| • Safety lancet | 1 × 50 pieces |
| • Alcohol pad (70% Isopropyl alcohol) | 1 × 50 pieces |
| • Disposable transfer pipette | 1 × 50 pieces |
| • Instruction for use | 1 piece |

Cassette
Sample diluent

Materials Required But Not Provided

- Specimen transfer device (graduated)
- Tips (Pipette)
- Protective gloves
- Specimen and test waste container
- Timer

- Collection tube (with or without anticoagulant, see Specimen Collection and Storage Requirements)

⚠ Warnings and Precautions

The test is for *in vitro* diagnostic use only. The test cassette is for single use only. The precautions are included, but not limited to the following:

[Safety Precautions]

- Wear gloves during the entire testing process.
- Do not use the accessories if seal is broken or package is damaged.
- All the accessories are for single use. Do not reuse the accessories.
- Do not pipette by mouth.
- Do not eat or smoke while handling specimens.
- Clean and disinfect all the areas that may be contaminated by spills of specimens or reagents with appropriate disinfectant.
- Decontaminate and dispose of all specimens, reagents, accessories and other potentially contaminated materials as if they were infectious wastes, in a biohazard container.

[Handling Precautions]

- The accessories are only used for specimen collection, do not store the specimen in the pipette.
- The assay is validated for use with 40uL of specimen, obtained when the specimen reaches the second graduation of the pipette. Use of other volumes may not produce the correct result.
- Do not mix and interchange different specimens.
- The sample diluent is specific for each lot of reagent. Do not use diluents from different lots.
- Do not perform the test under environment which leads to rapid evaporation (e.g. high temperature (> 45°C) and low humidity, fast air flow).
- Use of hemolytic specimens, rheumatoid factors-contained specimens and lipidemic, icteric specimens may cause invalid results or false results.
- Pooled specimens or specimens other than specified (e.g. saliva, urine) are strictly forbidden to be tested.
- Ensure that the specified volume of specimen is added by either precision pipette or specimen transfer pipette provided within the test kit.

[Warning]

- For *in vitro* diagnosis of the infection of HIV only, this assay will not provide a diagnosis for any other disease.
- The cassette, individually enclosed in a foil pouch and packed 50 tests per kit box, is only for "single use".
- Do not use the test cassette if the pouch is damaged or the seal is broken. Contact the supplier when it happens as this is an important indicator for quality and performance of the IVD.
- The sample diluent in this kit contains Sodium Azide (<0.1%) as the preservative, be cautious while using with the kit. In case of contact with sample diluent to eyes and/or skin, flush affected area with copious amount of water. Do not dispose the sample diluent directly in drains.
- Care should be taken to avoid the contact between the sample area of the cassette and the diluent bottle to avoid the contamination of the diluent.

Storage Requirements

- The test cassettes and the sample diluent must be stored at 4-30°C until expiration date.

- Do not use expired devices.
- Test cassette should be used immediately upon opening the pouch.
- Sample diluent should be used within 3 months of opening.

Specimen types

Human serum, plasma, and whole blood specimens should be collected through venipuncture (or whole blood through fingerstick).

Plasma specimens and venous whole blood specimens can be prepared with EDTA, heparin or sodium citrate as anticoagulant.

Specimen Collection and Storage Requirements

[Whole Blood]

- Collection by venipuncture

Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture.

Venous whole blood specimens can be stored up to 8 hours at 18-24°C or up to 7 days at 2-8°C if not tested at the same day of specimen collection, and should not be frozen.

Whole blood specimens should be brought to room temperature (18-28°C) prior to use.

Haemolytic or turbid whole blood specimens should not be used.

- Collection by fingerstick⁽⁶⁾

Figure a


Figure b


Figure c


Figure d

1. Choose the ring finger of the patient.
2. Open the alcohol pad and clean the fingertip (Figure a). Allow the fingertip dry.
3. Remove the cap of the safety lancet. (Figure b)
4. Place the lancet on the fingertip, press the lancet against the finger and puncture the skin. (Figure c)
5. Discard the lancet into a bio-hazard waste container immediately.
6. Squeeze the fingertip to form a bulb and use the disposable transfer pipette to collect the drop of the blood. (Figure d)
7. Allow the blood to reach the second graduation and the specimen volume will be 40ul. (Figure e)

Avoid formation of bubbles during the collection of specimen.

8. Transfer the specimen to the sample well immediately after collection. Avoid formation of bubbles.

9. Discard the disposable transfer pipette into a bio-hazard waste container immediately.

[Plasma or Serum]

- Plasma

Draw the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture and lightly mix, then centrifuge blood to get plasma supernatant.

- Serum

Draw the whole blood into the collection tube by venipuncture, leave to settle for 2 hours

Figure e

for blood coagulation and then centrifuge blood to get serum supernatant.


Serum and plasma specimens may be stored at 2-8°C for up to 7 days from time of draw, store at -18°C or below for long time storage. Several freeze-thaw cycles should be avoided. Plasma or Serum specimens should be brought to room temperature prior to use. Plasma or serum specimens with visible bio-contamination should not be used.

Test Procedure


Equilibrate all specimens and the devices to room temperature (18-28°C) before testing.

Sample diluent from different batches cannot be mixed.

1. Take out a test cassette from the foil pouch, and place it on a horizontal surface.
2. Use the specimen of either serum/plasma or whole blood: Add 40µl of specimen to the sample area firstly (Figure 1), then slowly add 1 drop of sample diluent vertically to the same area (Figure 2). Please avoid direct contact of the diluent bottle with the sample area.
3. Incubate the cassette at room temperature (18-28°C) and read the result after 15 minutes but not more than 25 minutes. Reading the test before 15 minutes or after 25 minutes may cause false result.


Interpretation of Results


Positive Result (Figure 3)

Clear visible reddish-purple bands appear both at the control line (C-line) and the test line (T-line) of the cassette. The test result is valid even the band of the control line (C-line) appears much lighter or darker than the band of test line (T-line).


Negative Result (Figure 4)

A reddish-purple band appears only at the control line (C-line) of the cassette (Figure 4).

Invalid Result 1 (Figure 5)

Reddish-purple band appears at neither the control line nor the test line of the cassette (Figure 5). Repeat the test of the specimen. If the result remains the same, then the result is invalid and contact KHB supplier.

**Invalid Result 2 (Figure 6)**

A reddish-purple band appears only at the test line (T-line) of the cassette (Figure 6). Repeat the test and contact KHB supplier if the C-line remains invisible.

Limitations of the Procedure

- The kit is designed to detect antibodies against HIV-1 and HIV-2 in human serum, plasma, and whole blood. Specimens other than the specified ones may not supply accurate results.
- The intensity of the T-line does not necessarily correlate to the titer of antibody in the specimen.
- The presence of the control line (C-line) only means the flow of the conjugate.
- When specimens contain high concentration of antibody to HIV-1 are tested on the device, the control line could be absent. In this case, please perform further analysis according to the instructions of Figure 6 in Interpretation of Results.
- As this product is intended to detect antibodies against HIV (1+2) from individuals, clinical diagnosis on AIDS should not be made only based on the test results.
- A non-reactive result should not exclude the possibility of infection caused by HIV-1 or HIV-2. A false negative result can occur in the following circumstances:
 - Negative result may not have detected recently acquired HIV infection.
 - Low levels of antibody (e.g., early seroconversion samples) are below the detection limit of the test.
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration, in exceptional cases this may lead to observation of false negative results.
 - Specimens are not properly stored.
 - False negative result due to very high concentration of antibody(ies) to HIV (which is known as 'hook effect').

For these reasons, care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- The product is not validated on specimens from infants, children, patient on anti-retroviral treatment.
- The product is not specifically designed for detection of antibody to HIV-1 group O.**

Performance Characteristics

- Specificity

A total of 3400 specimens including serum/plasma and whole blood were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 1)

Table 1

Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Specificity
		Negative	Positive	
Blood donor	2500	2500	0	100%
Whole blood donation	500	500	0	100%
Hospitalized patient specimens	200	200	0	100%
Pregnant Women	200	200	0	100%
Total	3400	3400		100%

A total of 488 potentially cross-reacting specimens were tested, the results are presented in Table 2.

Table 2

Potentially cross-reacting specimens	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2	
		Negative	Positive
Rheumatoid factor	10	10	0
Cytomegalovirus	5	5	0
Epstein Barr Virus	13	11	2
Malaria	85	85	0
Syphilis	52	51	1
Herpes	5	5	0
HBsAg	124	124	0
Anti-HCV	115	114	1
Anti-HBc	15	15	0
Anti-HTLV1/2	15	15	0
Anti-HEV	5	5	0
Specimens from STI clinics	10	10	0
Trypanosomiasis	24	24	0
E Coli infected	10	10	0
Total	488	484	4

- Sensitivity

A total of 750 HIV-1 and HIV-2 antibody positive specimens were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 3)

Table 3

Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Sensitivity
		Positive	Negative	
HIV-1 Ab positive specimens	360	360	0	100%
non-B subtypes HIV-1 Ab positive specimens	40	40	0	100%
HIV-2 Ab positive specimens	100	100	0	100%
HIV Ab positive specimens	250	250	0	100%
Total	750	750	0	100%

A total of 40 HIV-1 Non-B Subtypes positive specimens (in Table 3) were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2, the results are presented in Table 4.

Table 4

HIV Subtype	No. of specimens tested	Positive on KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2	
A	1		1
A1	1		1

KHB [®] DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2		For <i>in vitro</i> diagnostic use	EN
A2	1	1	
C	3	3	
CRF01_AE	3	3	
CRF02_AG	3	3	
CRF06_cpx	3	3	
CRF36_cpx	3	3	
D	3	3	
F1	2	2	
F2	1	1	
G	3	3	
H	3	3	
J	3	3	
K	3	3	
Group O	3	3	
Recombinant Form (recombination of subtype B, D, CRF19 and CRF21_A2D)	1	1	
Total	40	40	

A total of 100 paired HIV-1 antibody positive whole blood and plasma specimens were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 5)

Table 5

Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Sensitivity
		Positive	Negative	
HIV-1 Ab positive whole blood	100	100	0	100%
HIV-1 Ab positive plasma*	100	100	0	100%

*The 100 HIV-1 Ab positive plasma are included in the HIV-1 Ab positive plasma in Table 3.

Seroconversion panels

The sensitivity was evaluated on 30 commercially available seroconversion panels. Overall the Diagnostic kit for HIV (1+2) Antibody (colloidal gold) V2 detect 38 out of the 165 seroconversion panel members, which is 13 specimens more than the least sensitive CE marked antibody test and 15 specimens less than the most sensitive CE marked antibody test.

References

1. Gallo, R. C., S. Z. Salahuddin, M. Popovic, G. M. Shearer, M. Kaplan, B. F. Haynes, T. J. Palker, R. Redfield, J. Oleske, B. Safai, and et al. 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* **224**:500-503
2. Viscidi, R., E. Ellerbeck, L. Garrison, K. Midthun, M. L. Clements, B. Clayman, B. Fennie, and G. Smith. 1990. Characterization of serum antibody responses to recombinant HIV-1 gp160 vaccine by enzyme immunoassay. NIAID AIDS Vaccine Clinical Trials Network. *AIDS research and human retroviruses* **6**:1251-1256.
3. Chen, Y. H., Y. Xiao, W. Wu, Q. Wang, G. Luo, and M. P. Dierich. 2000. HIV-2 transmembrane protein gp36 binds to the putative cellular receptor proteins P45 and P62. *Immunobiology* **201**:317-322..

4. **Brandt, J., L. O. Andersson, and J. Porath.** 1975. Covalent attachment of proteins to polysaccharide carriers by means of benzoquinone. *Biochimica et biophysica acta* **386**:196-202.
5. **Osikowicz, G., M. Beggs, P. Brookhart, D. Caplan, S. Ching, P. Eck, J. Gordon, R. Richerson, S. Sampedro, D. Stimpson, and et al.** 1990. One-step chromatographic immunoassay for qualitative determination of choriogonadotropin in urine. *Clinical chemistry* **36**:1586.
6. **Spielberg, F., C. M. Kabeya, R. W. Ryder, N. K. Kifuani, J. Harris, T. R. Bender, W. L. Heyward, and T. C. Quinn.** 1989. Field testing and comparative evaluation of rapid, visually read screening assays for antibody to human immunodeficiency virus. *Lancet* **1**:580-584.

Key to symbols used

	Caution		Temperature limitation (4-30°C)
	Batch code		Manufacturer
	In vitro Diagnostic use Medical device		Consult instructions for use
	Catalogue number or order number		Do not use if package is damaged
	Use by		Sufficient for "n" use
	Do not reuse		

Shanghai Kehua Bio-engineering Co., Ltd.
 1189 North Qinzhou Road, 200233, Shanghai, P.R.China
 Tel: +86-21- 54500744
 +86-8008203370
 Fax: +86-21- 64854051
 Website: <http://www.skhb.com>

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