WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Bioline Malaria Ag P.f/P.f/P.v
WHO reference number: PQDx 0297-012-00

Bioline Malaria Ag P.f/P.f/P.v¹ with product codes **05FK120** and **05FK123**, manufactured by **Abbott Diagnostics Korea Inc**², CE-marked regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 04 December 2018.

Summary of WHO prequalification assessment for Bioline Malaria Ag P.f/P.f/P.v

	Date	Outcome
Prequalification listing	04-Dec-2018	listed
Dossier assessment	15-May-2018	MR
Site inspection(s) of quality management system	12-Oct-2018	MR
Product performance evaluation	2018	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Updated performance characteristics in the IFU.	29-Jan-2019
3.0	Product name was changed from SD BIOLINE Malaria Ag P.f/P.f/P.v to Bioline Ag P.f/P.f/P.v . Manufacturer's name changed from Standard Diagnostics Inc to Abbott Diagnostics Korea Inc. Manufacturing site inspection commitments for	26-Nov-2020

¹ Formerly called SD BIOLINE Malaria Ag P.f/P.f/P.v to Bioline Ag P.f/P.f/P.v.

² Formerly called Standard Diagnostics Inc.

	prequalification were fulfilled and removed from the public report.	
--	---	--

Intended use:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f/P.f/P.v test kit is an in vitro rapid, qualitative test for the detection of histidine rich protein II (HRP-II) antigen and plasmodium lactate dehydrogenase (pLDH) from Malaria Plasmodium falciparum and pLDH from Malaria Plasmodium vivax in human whole blood. Bioline Malaria Ag P.f/P.f/P.v test is intended for professional use, only for an initial screening test as an aid to diagnosis of clinical malaria disease and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of blood smear.”*

Assay description:

According to the claim of assay description from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f/P.f/P.v test utilizes the principle of immunochromatography. The Bioline Malaria Ag P.f/P.f/P.v test kit contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to histidine-rich protein II (HRP-II) of P. falciparum in test line 1(T1), mouse monoclonal antibodies specific to plasmodium lactate dehydrogenase (pLDH) of P. falciparum in test line 2(T2) and mouse monoclonal antibodies specific to pLDH of P. vivax in test line 3(T3) region. This test device has a letter of T1, T2, T3 and C as “Test Line” and “Control Line” on the surface of the device. All the Test Lines and Control Line in result window are not visible before applying any specimens. The Control Line is used for procedural control. Control Line should always appear if the test procedure is performed properly and the test reagents of control line are working. The principle for reactive results is as follows. As the P.f (HRP-II, pLDH) and P.v reactive specimen flow through the membrane assembly of the cassette after addition of the assay diluent, they form a complex with mouse monoclonal antibodies specific to P.f HRPII- gold conjugates and mouse monoclonal antibodies specific to pan pLDH -gold conjugates. This complex move further on the membrane to the test line region where it is immobilized by mouse monoclonal P.f HRPII antibodies, P.f pLDH antibodies and P.v pLDH antibodies coated on the membrane leading to formation of a red-purple colored band in the respective regions which confirms a reactive test result.*

Absence of a colored band in the test region indicates a non-reactive test result for the corresponding Antigen”.

Test kit contents:

Product code	Tests/kit	Component	Quantity
05FK120	25 T/kit	Test device with desiccant in individual foil pouches	25
		Assay diluent vial	1 × 5 ml/vial
		Disposable inverted cup	25 (5 µL/cup)
		Sterile lancet	25
		Alcohol swabs	25
		Instruction for Use	1
05FK123	1 test/kit (bundle kit)	Each kit contains all necessary accessories (Lancet, alcohol swab, disposable applicator, and assay diluent) for a single test	25 each of 1 test/kit

Items required but not provided:

- Protective gloves
- Timer
- Biohazard container

Storage:

The test kit should be stored at 1°C to 40 °C.

Shelf-life upon manufacture:

24 months.

Warnings/limitations:

- Performance of the Bioline Malaria Ag P.f/P.f/P.v test kit was evaluated during the WHO product testing of malaria RDTs: Round 8, during which it was observed that the *Pf* panel detection score (PDS) based only on the reactivity of the Pf – p-LDH line was 62% at 200 parasites/µl. This is below the WHO recommended performance criterion of greater than or equal to 75% PDS for *P.falciparum* at 200 parasites/µL. Therefore, false negative results can be expected in patients infected with *P.falciparum* parasites, with similar or lower parasitaemias, that do not express HRP2 due to *pfhrp2/3* gene deletions and are therefore reliant only on pf-LDH detection for *P.falciparum* detection.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- The test procedure, precautions and interpretation of results for this test must be followed precisely.

- Refer to current version of manufacturer's instructions for use.

Prioritization for prequalification:

Based on the results of the WHO product testing of malaria RDTs for Round 8, Bioline Malaria Ag P.f/P.f/P.v was given priority for WHO prequalification assessment.

Dossier assessment

Abbott Diagnostics Korea Inc. submitted a product dossier for Bioline Malaria Ag P.f/P.f/P.v as per the "*Instructions for compilation of a product dossier*" (PQDx_018 v1). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 15 May 2018.

Commitments for prequalification:

1. The manufacturer was requested in the stage following the Dossier Review to finalize the translation of the IFU. The manufacturer is requested to demonstrate implementation of the revised IFU within three months of the date of publication of this public report, 04 March 2019.
2. The manufacturer is requested to provide the final report of the evaluation of the real-time stability of the product configuration 05FK123 by September 2020.

The commitments are under review. WHO will verify the implementation of these commitments at the next re-inspection.

Based on the product dossier screening and assessment findings, the product dossier for Bioline Malaria Ag P.f/P.f/P.v meets WHO prequalification requirements.

Manufacturing site inspection

An inspection or a desk assessment of Abbott Diagnostics Korea Inc. located at 46, Hagal-ro 15 beon-gil, Giheung-gu, Yongin-si, Gyeonggi-do 17099, Republic of Korea and 65 Borahagal-ro, Giheung-gu, Yongin-si, Geonggi-do, Republic of Korea was conducted on 26-27 April 2018. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality. Routine inspections of the Manufacturing site will be conducted with copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the

information on the most current inspection or desk assessment performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection or desk assessment can be found at:

https://www.who.int/diagnostics_laboratory/evaluations/PQDxSiteInspection/en/

All published WHOPIRs are with the agreement of the manufacturer.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 12 October 2018.

Product performance evaluation

Bioline Malaria Ag P.f/P.f/P.v was evaluated in the eighth³ round of WHO product testing of RDTs for malaria antigen detection, which was completed in 2018.

Based on the demonstrated overall *P. falciparum* panel detection score against HRP2-expressing parasites (89.0% at 200 parasites/μl), *P. vivax* panel detection score (97.1% at 200 parasites/μl), false-positive rates (0% for clean negatives, 0% for *P. falciparum* at 200 parasites/μl, 0% for *P. vivax* at 200 parasites/μl, 0% for *P. falciparum* at 2000 to 5000 parasites/μl, 0% for *P. vivax* at 2000 to 5000 parasites/μl) and invalid rate (0%), Bioline Malaria Ag P.f/P.f/P.v meets the current laboratory evaluation requirements for prequalification. See Warning/limitations sections regarding use in areas with pfhrp2/3 gene deletions.

Summary performance characteristics	Panel detection score (%)		False positive rate (%)			Invalid rate (%)
	200 parasites/μl		200 parasites/μl		Clean negatives	
	Pf	Pv	Pf	Pv		
Bioline Malaria Ag P.f/P.f/P.v	89 (89/62) ¹	97.1	0	0	0	0

1. Product Panel Detection Score (PDS) of individual HRP2 test line and Pf-pLDH test line, respectively.

³ <https://www.who.int/malaria/publications/atoz/9789241514965/en/>


Labelling

- 1. Labels**
- 2. Instructions for use**


**Note: Labelling has been changed as per site name change and product rebranding. But temporarily labelling of legacy brand (SD BIOLINE) will be used in the market according to registration status in each country.*

1. Labels



1.1 Package box for 05FK120


1.2 Package box for 05FK123

1.3 Device pouch for 05FK120, 05FK123


1.4 Outer Pouch for 05FK123

1.5 Assay diluent label

1.6 Test device image

2. Instructions for use³

³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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

ENGLISH

Intended use

The Bioline™ Malaria Ag P.f/P.f/P.v test kit is a rapid, qualitative test for the detection of **histidine-rich protein II (HRP-II) antigen and plasmodium lactate dehydrogenase (pLDH) from Malaria Plasmodium falciparum and pLDH from Malaria Plasmodium vivax** in human whole blood.

REF 05FK123

Materials provided

- Test device
- Assay diluent
- Inverted cup (5 µl)
- Sterile lancet
- Alcohol swab

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the kit or its components.
Note : When stored at refrigerator, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to the test. Do not open the pouch whilst components come to room temperature.
2. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
3. Do not use the test kit beyond its expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.
6. When transporting or storing the kit, avoid exposure to high temperature (above 45 °C) for a period longer than 1 week.

Specimen collection using a lancet

1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the end of the fingertip and pierce with a sterile lancet provided.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. With a disposable inverted cup (5 µl) provided, dip the circular end of an inverted cup into the blood specimen and carefully place the circular end of the inverted cup into the round specimen well.

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© 2019 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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