## WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

**Product: STANDARD Q Malaria P.f Ag Test**  
**WHO reference number: PQDx 0346-117-00**

**STANDARD Q Malaria P.f Ag Test** with product codes **09MAL10D**, manufactured by **SD Biosensor, Inc, CE-Mark regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 6 March 2020.

### Summary of WHO prequalification assessment for STANDARD Q Malaria P.f Ag Test

	Date	Outcome
<b>Prequalification listing</b>	06-Mar-2020	listed
<b>Dossier assessment</b>	29-Oct-2019	MR
<b>Site inspection(s) of quality management system</b>	17-19-Apr-2019	MR
<b>Product performance evaluation</b>	2018	MR

MR: Meets Requirements

### Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendment	Date of report amendment
2.0	1. Change of IFU from IFU in English to multilanguage IFU (EN, ES, FR, PT), 2. Change of outer package, 3. Add additional supplier for sterile lancet.	18-May-2021

### Intended use

According to the claim of intended use from SD Biosensor, Inc, "*STANDARD Q Malaria P.f Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of Plasmodium falciparum (P. falciparum) specific Histidine Rich Protein 2 (HRP-2) in human*

*capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f Ag Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing."*

### Assay description

According to the claim of assay description from SD Biosensor, Inc, "STANDARD Q Malaria P.f Ag Test contains two pre-coated lines, "P.f" (*P. falciparum*) as test line and "C" as control line on the surface of the nitrocellulose membrane. Both test line and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-P. falciparum HRP-2 is coated on the test line region and monoclonal anti-chicken IgY is coated on the control line region. During the test, the *P. falciparum* specific HRP-2 in the specimen reacts to the gold-conjugated monoclonal anti-Malaria HRP-2 and binds to it. Any *P. falciparum* specific HRP-2 antigen-antibody gold particle complex also migrates with the buffer and is immobilized by monoclonal anti-P. falciparum HRP-2 at the test line to formation of a violet test colored band which confirms a positive result. Absence of this violet colored band indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working."

### Test kit contents

Component	25 tests (product code 09MAL10D)
Test device (individually in a foil pouch with desiccant)	25
Buffer Bottle	1 x 4 mL
Inverted cup (5µl)	Pack of 25
Sterile lancet	Pack of 25
Alcohol swab	Pack of 25
Instructions for use	1

### Items required but not provided

- Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood
- Micropipette and tip
- Timer
- PPE (Personal Protective Equipment)
- Pen/pencil
- Extra lancets and alcohol swabs
- Biosafety sharps container
- Biohazard container
- Sterile gauze

**Storage**

The test kit should be stored at 2-40°C.

**Shelf-life upon manufacture**

24 months.

**Warnings/limitations**

Refer to current version of manufacturer's instructions for use.

**Prioritization for prequalification**

Based on the established eligibility criteria, STANDARD Q Malaria P.f Ag Test was given priority for WHO prequalification assessment.

**Dossier assessment**

SD Biosensor, Inc submitted a product dossier for STANDARD Q Malaria P.f Ag Test as per the "*Instructions for compilation of a product dossier*" (PQDx\_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 29 October 2019.

**Commitment for prequalification**

Manufacturer commits to submit the final report of the ongoing shipping stability study by May 2021. WHO will follow-up on implementation of this commitment at the next re-inspection, or at the date assigned in the commitment review letter.

Based on the product dossier screening and assessment findings, the product dossier for STANDARD Q Malaria P.f Ag Test meets WHO prequalification requirements.

**Manufacturing site inspection**

An inspection of **SD Biosensor Inc.**, located at 74, *Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea* was conducted between 17-19 April 2019. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality

management system and manufacturing practices in place that would support the manufacture of the product of consistent quality. Routine inspections of the Manufacturer will be conducted with copies of these WHO Public Inspection Reports (WHOPIRs) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. To note that a WHOPIR reflects the information on the most current inspection performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection can be found at:

[https://www.who.int/diagnostics\\_laboratory/evaluations/PQDxSiteInspection/en/](https://www.who.int/diagnostics_laboratory/evaluations/PQDxSiteInspection/en/)

All published WHOPIRs are with the agreement of the manufacturer.

## Product performance evaluation

**STANDARD Q Malaria P.f Ag Test**<sup>1</sup> was included in the eighth round of WHO product testing of RDTs for malaria antigen detection, which was completed in 2018. The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild-type parasite panel, *P. vivax* wild-type parasite panel and a negative panel. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

Based on the demonstrated *P. falciparum* panel detection score (87.0% at 200 parasites/ $\mu$ l), false-positive rates (0% for clean negatives, 0% for *P. vivax* at 200 parasites/ $\mu$ l, 0% for *P. vivax* at 2000 / $\mu$ l) and invalid rate (0%), STANDARD Q Malaria P.f Ag Test meets the current laboratory evaluation requirements for prequalification.

Summary performance characteristics	Panel detection score		False positive rate (%)			Invalid rate (%)
	200 parasites/ $\mu$ l		200 parasites/ $\mu$ l		Clean negatives	
	Pf	Pv	Pf	Pv		
STANDARD Q Malaria P.f Ag Test	87.0	NA	NA	0	0	0

Based on these results, **STANDARD Q Malaria P.f Ag Test** meets the current performance evaluation requirements for prequalification.

<sup>1</sup> The product code of the product included in the WHO product testing of malaria RDTs was 09MAL10B, which corresponded to the same product, with a different kit configuration (i.e. number of tests) <https://www.who.int/malaria/publications/atoz/9789241514965/en/>

## **Labelling**

- 1. Labels**
- 2. Instructions for use**


## **1. Labels**

### **1.1 Device Package**

3

SD BIOSENSOR

www.sdbiosensor.com

CONTENTS

1) Test device (individually in a foil pouch with desiccant) x 25  
 2) Inverted cup (5µl) x 25  
 3) Buffer bottle x 1  
 4) Sterile lancet x 25  
 5) Alcohol swab x 25  
 6) Instructions for use x 1

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

Manufacturing site  
74, Daepyeong-ro, Daepyeong-eup, Daepyeong-gun, Gyeonggi-do, 10770, REPUBLIC OF KOREA  
+82-31-300-0400 www.sdbiosensor.com

EC REP

Authorized Representative  
**MIT Promed Consulting GmbH**  
Altenhofstrasse 69 66388 St. Ingbert Germany  
Phone: +49 6894 581020 Fax: +49 6894 581021

STANDARD™

### 1.2 Foil pouch

1.3 Buffer label

**STANDARD Q Malaria P.f Ag**

Product Code: **GMLD1-3**

---

**Buffer bottle**

LOT XXXXXX Vol : 4ml

YYYY.MM.DD. YYYY.MM.DD.

IVD

2°C (36°F) 40°C (104°F)

L36RT1ENR1 Issue date : 2017.10

 **STANDARD**<sup>TM</sup>

1.4 Inverted cup(5 µl) label

**Inverted cup (5µl)**

LOT No. : XXXXXX

EXP : YYYY.MM.DD.

Quantity : 25PCS

IVD

Σ<sub>25</sub>

L15IC1ENR2 Issue date : 2020.02

 **SD BIOSENSOR**

1.5 Sterile Lancet labels

# Disposable Sterile Lancets

 LOT No. : **XXXXXX**

 MFG Date : **YYYY.MM.DD.**

 EXP Date : **YYYY.MM.DD.**

**INTENDED USE**  
To obtain a capillary blood specimen from the fingertip.

**INSTRUCTIONS FOR USE**  
To use, twist-off the protective cap.

**CAUTION**  
The lancet is guaranteed sterile while protective cap is sealed to the base.  
Do not use if the seal has been damaged or broken.

 **Manufactured by**  
Beijing Ruicheng Medical Supplies Co., Ltd.  
No. 558 Zhangzikou, Yangsong Town,  
Huairou District, 101400 Beijing, China

Product code : 01GL25

## 25PCS

## 28G

B06LCT2ENR2  
Issue Date : 2020.02

 **Authorized Representative**  
Lotus NL B.V.  
Koningin Julianaplein 10, 1e Verd,  
2595AA, The Hague, Netherlands.  
Tel: +31645171879, +31626669008

# Disposable Sterile Lancets

 LOT No. : **XXXXXX**

 MFG Date : **YYYY.MM.DD.**

 EXP Date : **YYYY.MM.DD.**

**INTENDED USE**  
To obtain a capillary blood specimen from the fingertip.

**INSTRUCTIONS FOR USE**  
To use, twist-off the protective cap.

**CAUTION**  
The lancet is guaranteed sterile while protective cap is sealed to the base.  
Do not use if the seal has been damaged or broken.

 **Manufactured by**  
Tianjin Huahong Technology Co., Ltd.  
A01, Plant B No. 278, Hangkong Road,  
Tianjin Pilot Free Trade Zone  
(Air Port Industrial Park),  
300308 Tianjin, China

Product code : 01GL27

## 25PCS

## 28G

B05SL2ENR0  
Issue Date : 2020.07

 **Authorized Representative**  
Shanghai International Holding Corp.  
GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany  
Tel: +49-40-2513175

### 1.6 Alcohol swab label

Front

Back

20

1.7 Cassette image

## **2. Instructions for use<sup>2</sup>**

---

<sup>2</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

**STANDARD Q Malaria P.f Ag**  
STANDARD™ Q Malaria P.f Ag Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

**KIT CONTENTS**

**MATERIALS REQUIRED BUT NOT PROVIDED**

- Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood
- Micropipette and tip
- Timer
- PPE (Personal Protective Equipment)
- Pen/pencil
- Extra lancets and alcohol swabs
- Biosafety sharps container
- Biohazard container
- Sterile gauze

**SPECIMEN COLLECTION AND PREPARATION**

**Capillary whole blood**

- Capillary whole blood should be collected aseptically by fingerstick.
- Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
- Squeeze the end of the fingertip and pierce with a sterile lancet.
- Wipe the first drop with a sterile gauze.
- Collect the capillary whole blood and it must be tested immediately after collection.

**Venous whole blood**

- Collect the venous whole blood into the commercially available anticoagulant tube such as heparin, EDTA, or sodium citrate by venipuncture.
- If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/ 36-46°F, the specimen can be used for testing within 72 hours after collection.
- Do not use hemolyzed blood specimen.


Anticoagulants such as heparin, EDTA, or sodium citrate do not affect the test result. As known relevant interference, hemolytic specimen, rheumatoid factors-contained specimen and lipemic, icteric specimen may lead to inaccurate results. Use separate inverted cup for each specimen in order to avoid cross-contamination of either specimens which can cause erroneous results.

**PREPARATION AND TEST PROCEDURE**

**Preparation**

- Carefully read the instructions for using the STANDARD Q Malaria P.f Ag Test.
- Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
- Open the foil pouch, and check the test device and the desiccant inside the foil pouch.


**Test Procedure**


Capillary and venous whole blood can be used for the STANDARD Q Malaria P.f Ag Test.

**1. Test with a sterile lancet**


- Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.



- Dry and pierce the wiped fingertip with a sterile lancet to bleed.

- Wipe away the first drop of blood with sterile gauze.

- Take the inverted cup(5µl) provided and collect 5µl of the capillary whole blood specimen by dipping the circular end of the inverted cup into the specimen.
- To apply the collected capillary whole blood specimen to the specimen well(circle) of the test device, place the circular end of the inverted cup on the specimen well pad and then press down lightly.


**1-1. Test by venipuncture**


- Gently invert 4-5 times for proper mixing.

- Using a micropipette, draw 5µl of venous whole blood specimen.

- Apply the collected venous whole blood specimen to the specimen well(circle) of the test device.

**2. Dropping of Buffer**

- Hold the buffer bottle at 90° (perpendicular) to the test device and not touching the specimen well to avoid contamination. Add 3 drops (90µl) of the buffer into the buffer well (square) of the test device.


**3. Reading Time**

- Read the test result in 15 minutes. The test can be read up to 30 minutes.

Do not read test results after 30 minutes. It may give false results.


**INTERPRETATION OF TEST RESULT**

The presence of only one colored band ("C" Control line) within the result window indicates a negative result.

**P.f. (Plasmodium falciparum) Positive**  
The presence of two colored bands ("C" Control line and "P.f." Plasmodium falciparum line) within the result window, no matter which band appears first, indicates a P.f. positive result.

The presence of any line, no matter how faint the result is considered positive.

If the control band ("C" Control line) is not visible within the result window, the result is considered invalid. The directions may not have been followed correctly or the test device may have deteriorated. Re-test with a new specimen and a new test device.

- A colored band will appear in the top section of the result window to show that the test is working properly. This band is the control line (C).
- A colored band will appear in the lower section of the result window. This band is the Plasmodium falciparum line (P.f.).

**\*This test only indicates the presence of Histidine-Rich Protein 2 (HRP2) of Plasmodium falciparum in human capillary and venous whole blood specimen and should not be used as the sole criteria for the diagnosis of malaria. As other diagnostic tests, all test results should be considered with other clinical history available to the physician.**

**EXPLANATION AND SUMMARY**

**Introduction**

Malaria remains an important cause of illness and death in children and adults in countries in which it is endemic. Globally, an estimated 3.2 billion people in 97 countries and territories are at risk of being infected with malaria and developing disease, and 1.2 billion are at high risk (>1 in 1,000 cases of getting malaria in a year). According to the World Malaria Report 2015, there were 210 million cases of malaria globally in 2015 (uncertainty range 149-303 million) and 438,000 malaria deaths (range 236,000-635,000), representing a decrease in malaria cases and deaths of 37% and 60% since 2000, respectively. The protozoal parasites that cause malaria are from Plasmodium falciparum, vivax, ovale and malariae with the first two species causing the most infections worldwide. Classic symptoms of malaria include fever, headache, chills, vomiting, shivering and convulsions. In some rare forms of P. falciparum, the patient may present with delirium or coma. Several anemia is often attributed to the cause of death from malaria. Accurate and prompt diagnosis of malaria is of utmost importance due to the morbidity associated with the other malaria forms. Rapid diagnostic test is an ideal diagnostic tool for malaria diagnosis in that it can provide a rapid determination if the patient is infected with malaria allowing for accurate treatment and improved outcomes. STANDARD Q Malaria P.f Ag Test, a reliable and sensitive screening test, would enhance the accuracy of the diagnosis of malaria infection and thus make clinical treatment decision effectively.

P.f.: Plasmodium falciparum  
P.v.: Plasmodium vivax  
P.o.: Plasmodium ovale  
P.m.: Plasmodium malariae

**Intended use**

STANDARD Q Malaria P.f Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of Plasmodium falciparum (P. falciparum) specific Histidine Rich Protein 2 (HRP-2) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f Ag Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.

**Test principle**

STANDARD Q Malaria P.f Ag Test contains two pre-coated lines, "P.f." (P. falciparum) as test line and "C" as control line on the surface of the nitrocellulose membrane. Both test line and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-P. falciparum HRP-2 is coated on the test line region and monoclonal anti-chicken IgY is coated on the control line region. During the test, the P. falciparum specific HRP-2 in the specimen reacts to the gold-conjugated monoclonal anti-Malaria HRP-2 and binds to it. Any P. falciparum specific HRP-2 antigen-antibody gold particle complex also migrates with the buffer and is immobilized by monoclonal anti-P. falciparum HRP-2 at the test line to formation of a violet colored band which confirms a positive result. Absence of this violet colored band indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

**KIT STORAGE AND STABILITY**

- Store the sealed pouch and the buffer provided in the kit at 2-40°C/36-104°F out of the direct sunlight for the duration of its shelf life.
- Do not open the foil pouch until you are ready to perform a test.
- Close the buffer cap tightly after using, and then store it at 2-40°C/36-104°F out of the direct sunlight. It is stable until the expiry date of the kit and the buffer label after opening its cap, if it is tightly closed.

**WARNINGS AND PRECAUTIONS**

- Do not freeze.
- Do not use beyond the expiration date.
- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use buffer of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- In the pouch, there is a desiccant containing a humidity indicator, the humidity indicator should be yellow. If the humidity indicator is green, throw away the test device and take another test device. Throw away the desiccant in the nonsharps (non-infectious) disposal.
- Before testing, check that the lot number of buffer on the bottle matches the lot number of buffer on the package and the pouch you are using.
- The buffer contains 0.095% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.
- For in vitro diagnostic use only.

**LIMITATION OF TEST**

- The test procedure, precautions, and interpretation of results for this test must be followed strictly when testing.
- STANDARD Q Malaria P.f Ag Test is designed for use only on human capillary and venous whole blood specimens.
- Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a

**QUALITY CONTROL**

- A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent is reactive.
- Control materials are not supplied with this test kit. However, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
- If there is a problem with the result such as invalid result, retest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

**Potential Cross-reacting Substances**

STANDARD Q Malaria P.f Ag Test had no cross-reaction with the following microorganisms such as other hepatitis (Hepatitis A virus, Hepatitis B virus, Hepatitis C virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, E. coli (K99), Salmonella Typhi), Streptococcus pneumoniae, Streptococcus pyogenes, Legionella pneumophila, Treponema pallidum. Also, this product had no cross-reaction with HIV Serconversion panel, HCV positive plasma, HCV core Ag positive plasma, Hbs Ag positive plasma, HAV IgM positive plasma, Dengue IgM positive plasma, Tick borne encephalitis IgM positive plasma, Chikungunya positive plasma, Zika virus positive plasma, Salmonella typhi IgM positive plasma, Leishmania positive plasma, Japanese Encephalitis positive plasma, Brucella IgM positive plasma and West Nile Virus positive plasma.

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