WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: NxTek Eliminate Malaria Pf and NxTek Malaria Pf¹
WHO reference number: PQDx 0349-012-00

NxTek Eliminate Malaria Pf with product codes **05FK140, 05FK141, 05FK142, and 05FK143, and NxTek Malaria Pf** with product code **05FK144**, manufactured by **Abbott Diagnostics Korea Inc²**, **CE-marked regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 12 April 2019.

Summary of WHO prequalification assessment NxTek Eliminate Pf and NxTek Malaria Pf

	Date	Outcome
Prequalification listing	12 April 2019	listed
Dossier assessment	12 January 2019	MR
Site inspection(s) of the quality management system	26 April 2018	MR
Product performance evaluation	2016	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarised in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	The product name was changed from Alere Malaria Ag P.f to NxTek Eliminate Pf. The manufacturer's name was changed from Standard Diagnostics Inc. to Abbott Diagnostics Korea Inc. Manufacturing site inspection commitments for prequalification were fulfilled and removed from the public report.	26 November 2020

¹ The new name and product code for the product for the Ghanaian and Nigerian markets with product name NxTek Malaria Pf and product code 05FK144. The product is the same as the prequalified product, the NxTek Eliminate Pf.

² Manufacturer's name changed from Standard Diagnostics Inc to Abbott Diagnostics Korea Inc.

3.0	Labelling change that included a 2D Barcode System for traceability. It is used for UDI (Unique Device identification), serialisation and aggregation for each component and product to improve patient safety and healthcare business processes.	7 March 2022
4.0	Introduction of a new name and product code for the product for the Ghanaian and Nigerian markets with product name NxTek Malaria Pf and product code 05FK144. Changes to alcohol, sterile lancet, and inverted cup labels.	11 September 2023

Intended use

According to the claim of intended use from Abbott Diagnostics Korea Inc, “*NxTek Eliminate Malaria Pf is an ultra sensitive, qualitative in vitro test for the detection of histidine-rich protein 2 (HRP2) antigen of Plasmodium falciparum malaria in human whole blood. The test is intended for professional use and the point-of-care testing to aid in the diagnosis of malaria infection.*”

Assay description

According to the claim of assay description from Abbott Diagnostics Korea Inc, “*NxTek Eliminate Malaria Pf test contains a membrane strip, which is pre-coated with immobilised poly streptavidin on the test line and mouse monoclonal anti-chicken IgY on the control line. Two types of conjugates (mouse rFab fragment specific to HRP2 of P.f carboxyl modified latex conjugate and chicken IgY gold conjugate) move upward on the membrane chromatographically and react with biotinylated sheep rFab specific to HRP2 of P. falciparum and pre-coated mouse monoclonal anti-chicken IgY respectively. For a positive result, biotinylated sheep recombinant fragment antigen-binding (rFab fragment) specific to HRP2 of P. falciparum and mouse rFab fragment specific to HRP2 of P.f – carboxyl modified latex (CML) conjugate will form in the test line on the P.f region of the result window. Neither the test or control lines are visible in the result window prior to applying the patient specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional.*”

Test kit contents:**NxTek Eliminate Malaria Pf**

Product code	Test/kit	Component	Quantity
05FK140	25T/Kit	Test device	25
		Assay diluent	1
		Disposable inverted cup	25
		Alcohol swab	25
		Sterile lancet	25
		Instructions for Use (IFU)	1
05FK141	25T/kit	Test device	25
		Assay diluent	1
		Disposable inverted cup	25
		Alcohol swab	25
		Safety lancet	25
		Instructions for Use	1
05FK142	1T/kit x 25 each	Each kit contains a test device, an assay diluent, a disposable inverted cup, an alcohol swab, a safety lancet and a summarised IFU for a single test.	25 each of 1 test/kit
	In addition to the summarised IFU in each kit, a bundle pack IFU is included in the package box.		
05FK143	1T/kit x 25 each	Each kit contains a test device, an assay diluent, a disposable inverted cup, an alcohol swab, a sterile lancet and a summarised IFU for a single test.	25 each of 1 test/kit
	In addition to the summarised IFU in each kit, a bundle pack IFU is included in the package box.		

NxTek Malaria Pf

Product code	Test/kit	Component	Quantity
05FK144	25T/Kit	Test device	25
		Assay diluent	1
		Disposable inverted cup	25
		Alcohol swab	25
		Sterile lancet	25
		Instructions for Use (IFU)	1

Items required but not provided

- Protective gloves
- Timer
- Biohazard container
- Sterile gauze or cotton
- If collection by venipuncture: Venipuncture blood collection tube, Micropipette, Pipette tips

Storage

The test kit should be stored at 1 °C and 30 °C.

Shelf-life upon manufacture

12 months.

Warning/limitation

Please refer to the instructions for use attached to this public report.

Prioritisation for prequalification

Based on the WHO product testing of malaria RDTs results for Round 7 and eligibility criteria, NxTek Eliminate Malaria Pf was prioritised for the WHO prequalification assessment.

Dossier assessment

Abbott Diagnostics Korea Inc. submitted a product dossier for NxTek Eliminate Malaria Pf as per the “*Instructions for compilation of a product dossier*” (PQDx_018). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 14 January 2019.

Based on the product dossier screening and assessment findings, the product dossier for NxTek Eliminate Malaria Pf meets WHO prequalification requirements.

Manufacturing site inspection

An inspection or a desk assessment of Abbott Diagnostics Korea Inc. located at 46, Hagal-ro 15 beon-gil, Giheung-gu, Yongin-si, Gyeonggi-do 17099, Republic of Korea and 65 Borahagal-ro, Giheung-gu, Yongin-si, Geonggi-do, Republic of Korea was conducted on 26-27 April 2018. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality.

Routine inspections of the Manufacturing site will be conducted with copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current inspection or desk assessment performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection or desk assessment can be found at:

<https://extranet.who.int/pgweb/inspection-services/prequalification-reports/whopirs-vitro-diagnostics>

All published WHOPIRs are with the agreement of the manufacturer.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 12 October 2018.

Product performance evaluation

The seventh round of WHO product testing of RDTs for malaria antigen detection was completed in 2016. The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild-type parasite panel and *Plasmodium spp.* negative panel. Thermal stability was assessed after two months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

Based on the demonstrated *P. falciparum* panel detection score (98% at 200 parasites/μl), false-positive rates (0.9% for clean negatives, 0% for *P. vivax* at 200 parasites/μl, 0% for *P. vivax* at 2000 parasites/μl) and invalid rate (0.1%), **NxTek Eliminate Malaria Pf** meets the current laboratory evaluation requirements for prequalification.

Summary performance characteristics	Panel detection score	False positive rate (%)		Invalid rate (%)
	200 parasites/μl	200 parasites/μl	Clean negatives	
	<i>Pf</i>	<i>Pv</i>		
NxTek Eliminate Malaria Pf	98%	0%	0.9%	0.1%


Labelling

1. Labels
2. Instructions for use


**Note: Labelling has been changed per site name change and product rebranding. However, temporary labelling of the legacy brand (Alere) will be used in the market according to registration status in each country.*

1. NxTek Eliminate Malaria Pf Labels


1.1.1 Package box for 05FK140


1.1.2 Package box for 05FK141

1.1.3 Package box for 05FK142


1.1.4 Package box for 05FK143


1.1.5 Device pouch for 05FK140, 05FK141, 05FK142, 05FK143


Front

Back



1.1.6 Assay diluent label

1.1.7 Test device image


Product codes 05FK140 and 05FK141

Product codes 05FK142 and 05FK143

1.1.8 Alcohol swab labels

1.1.9 Lancet labels

(01)16957138602815
(1)YYMMDD
(17)YYMMDD
(10)XXXXXXXX

LOT
[Symbol]

EC REP
Shanghai International Holding Corp. GmbH (Europe)
Elffestrasse 80, 20637 Hamburg, Germany
Tel: 0049-40-2513175
E-mail: shholding@hotmail.com

Tianjin Huahong Technology Co., Ltd.
A01, Plant B, No.278, Hangtong Road,
Tianjin Pilot Free Trade Zone (A1 Port
Industrial Park), 300205 Tianjin, China
Tel: 0086 22 58775333 400 8587753
Fax: 0086 22 58775353
E-mail: huahongsales@163.com

Twist-off lancet (use with lancing device)
25 pcs (28G I)

[Instructions]
1. Unscrew the lancing device cap from the body of the lancing device.
2. Insert a sterile lancet into the lancing device firmly.
3. Twist off the protective cap of the lancet.
4. Carefully screw the lancing device cap back to the lancing device.
5. Adjust the puncture depth.
6. Put the lancing device against the side of the finger and press the release button to prick your fingertip. If blood can't be collected, please adjust the puncture depth according to the instruction of the lancing device, and repeat this step.
7. Apply a drop of blood to test strip.
8. Recap lancet and discard it into an appropriate container.
9. Place the lancing device cap back to the lancing device.

[Structure] The twist-off lancets consist of needle, main body and protective cap.
[Performance] The device is sterilized by irradiation.
[Intended use] The twist-off lancet (use with lancing device) is a single use, sterile medical device, which is designed for use of micro blood sampling puncture to obtain capillary blood samples from the fingertip. The twist-off lancet needs using with lancing device.
[Intended users] Healthcare personnel, only when used with single lancing device to obtain blood sample from single patient.
[Intended patient population] Patient who needs to obtain capillary blood samples from fingertip, with no restriction to gender, weight or race.
[Indication] The lancet is intended for people who need capillary blood sampling.
[Storage] The products should be stored in a well-ventilated room and kept away from light, moisture and corrosive gases.
[Contraindications] None.


 The lancet and lancing device must be compatible with each other. When you use Huahong's twist-off lancet with other manufacturers' lancing devices, please be sure that these products are compatible. We recommend that you shall consult a health care professional before use.

[Cautions]
1. Read the instructions carefully before use.
2. Use the lancets within the expiry date.
3. Do not use if protective cap is removed or broken previously.
4. Do not store lancet in lancing device, discard lancet after use, do not re-use.
5. Lancet is only for single-use, and re-use of the device could lead to cross-infection and damage of the needle.
6. If the used needle is not recapped, it might cause sharp injury of other person and might cause infection or cross-infection.
7. If the needle is broken and embedded into the finger, please consult a healthcare professional immediately.
8. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

[LOT] Batch code  Date of manufacture 


 Use-by date  Manufacturer

 Do not re-use  Do not resterilize

 Medical device  Caution

 Sterilized using irradiation

 Consult instructions for use

 Do not use if package is damaged and consult instructions for use

EC REP Authorized representative in the European community

CE CE Marking **CE 0197**

Revision: 2022-09-27

1.1.10 Disposable inverted cup

Inverted cup 80mm (5µl)

Code No.: P65-003

LOT No.:

Quantity: 25PCS

MFG Date:

AIM, Inc.
101-908, Digital Empire II, 88, Sinwon-ro, Yeongtong-gu
Suwon-si, Gyeonggi-do, 16581, South Korea

EC REP Cmc Medical Devices & Drugs S.L.
C/ Horacio Lergam 18 CP 29006 Málaga-Spain

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1°C | 30°C
 34°F | 86°F

©2021 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

1.2.2 Foil pouch label


Front

Back

1.2.3 Nx Tek Malaria Pf Device

1.2.4 Nx Tek Malaria Pf Buffer label

Note: Disposable inverted cup sampling device, alcohol swab and sterile lancet. Please refer to the labels above.

2. Instructions for use³

² ~~English version of the IFU~~ was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

2.1 IFU for 05FK140

Type of specimen	No.	Potential interfering substance	Test concentration	
Emergent substance	11	Amodiaquine	393 ng/ml	
	12	Desethylamodiaquine	751 ng/ml	
	13	Arimether	540 ng/ml	
	14	Lumefantrine	28300 ng/ml	
	15	Artesunate	29544 ng/ml	
	16	Dihydroartemisinin	2043 ng/ml	
	17	Atovaquone	13270 ng/ml	
	18	Primaquine	730 ng/ml	
	19	Clindamycin	14 ng/ml	
	20	Dihydroartemisinin	698 ng/ml	
	21	Piperaquine	730 ng/ml	
	22	Doxycycline	0.0069 ng/ml	
	23	Mefloquine	3279 ng/ml	
	24	Quinine	295 ng/ml	
	25	Sulfadoxine	19000 ng/ml	
			Pyrimethamine	860 ng/ml

*There was no interference until biotin concentration 375 ng/ml, but if higher concentrations are present in the blood, it can be affected to the performance of the NxTek™ Eliminate Malaria Pf. Cross-reactivity testing was performed with normal whole blood spiked with serum, plasma or culture fluid specimens which contain potential cross-reacting pathogens. The following 22 pathogens had no effect on test results of NxTek™ Eliminate Malaria Pf. Although no cross reactivity occurred in this internal study, the possibility cannot be excluded completely.

3. Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptospira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity. Reproducibility of NxTek™ Eliminate Malaria Pf test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), day to day (five days), between run (three different test operators), between sites (three sites), and lot-to-lot (3 different lots).

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
Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptospira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity.

Reproducibility of NxTek™ Eliminate Malaria Pf test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), day to day (five days), between run (three different test operators), between sites (three sites), and lot-to-lot (3 different lots).

4. Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptospira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity. Reproducibility of NxTek™ Eliminate Malaria Pf test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), day to day (five days), between run (three different test operators), between sites (three sites), and lot-to-lot (3 different lots).

Type d'échantillon	N°	Substance potentiellement interférente	Concentration du test
Substance étrangère	11	Amodiaquine	393,2 ng/ml
	12	Diéthylamodiaquine	751 ng/ml
	13	Artemether	540 ng/ml
	14	Luméfantérine	28 300 ng/ml
	15	Artésunate	29 544 ng/ml
	16	Dihydroartémisine	2 043 ng/ml
	17	Atovaquone	13 270 ng/ml
	8	Rubella	751 ng/ml
	9	West Nile Virus (WNV)	14 ng/ml
	18	Chloroquine	1 430 ng/ml
	19	Clindamycine	14 ng/ml
	20	Dihydroartémisine	698 ng/ml
	21	Pipérasquine	730 ng/ml
	22	Doxycycline	0,0069 ng/ml
	23	Méfloquine	3 279 ng/ml
	24	Quinine	295 ng/ml
	25	Sulfadoxine	19 000 ng/ml
		Pyriméthamine	860 ng/ml
		Chagas	
		HAMA (anticorps humain anti-souris)	
		Syphilis	

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*There was no interference until boitin concentration 375 ng/ml, but if higher concentrations are present in the blood, it can be affected to the performance of the NxTek™ Eliminate Malaria Pf. Cross-reactivity testing was performed with normal whole blood spiked with serum, plasma or culture fluid specimens which contain potential cross-reacting pathogens. The following 22 pathogens had no effect on test results of NxTek™ Eliminate Malaria Pf. Although no cross reactivity occurred in this internal study, the possibility cannot be excluded completely.

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PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

1 **EN** Open the package and look for the following:

1. Test device with desiccant in individual foil pouch
2. Assay diluent
3. Disposable inverted cup (5 µl)
4. Safety lancet
5. Alcohol swab
6. Instructions for use

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
3 **EN** Next, look at the expiration date on the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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4 **EN** Open the foil pouch and look for the following:
1. Result window
2. Specimen well
3. Assay diluent well
Then, label the device with the patient identifier.

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2 **EN** Dispense 4 drops of assay diluent into the round specimen well.

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

3 **EN** Dispense 4 drops of assay diluent into the square assay diluent well.

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1 **EN** Clean the area to be lanced with an alcohol swab.

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2 **EN** Twist off or remove the protective cap from the safety lancet to break the seal and discard.

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3 **EN** Prick the lateral side of the finger with the safety lancet provided. Then, safely dispose of the lancet immediately after.

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Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptospirosis, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity. Reproducibility of NxTek™ Eliminate Malaria Pf test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), day to day (five days), between run (three different test operators), between sites (three sites), and lot-to-lot (3 different lots).

INTERPRETATION / INTERPRÉTATION / INTERPRETACIÓN / INTERPRETAÇÃO

NEGATIVE / NÉGATIF / NEGATIVO

5 **EN** The presence within the result window of one colored line at "C" Control line (violet) indicates a negative result.

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POSITIVE / POSITIF / POSITIVO

5 **EN** The presence within the result window of two colored lines; one each at "C" Control line (violet) and "P.F" Test line (any blue) indicates a positive result.

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5 **EN** If the "C" Control line (violet) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

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
*There was no interference until boatin concentration 375 ng/ml, but if higher concentrations are present in the blood, it can be affected to the performance of the NxTek™ Eliminate Malaria Pf. Cross-reactivity testing was performed with normal whole blood spiked with serum, plasma or culture fluid specimens which contain potential cross-reacting pathogens. The following 22 pathogens had no effect on test results of NxTek™ Eliminate Malaria Pf. Although no cross reactivity occurred in this internal study, the possibility cannot be excluded completely.

Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptospira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity. Reproducibility of NxTek™ Eliminate Malaria Pf test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), day to day (five days), between run (three different test operators), between sites (three sites), and lot-to-lot (3 different lots).

PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

- 1** **EN** Open the package and look for the following:
1. Test device with desiccant in individual foil pouch
 2. Assay diluent
 3. Disposable inverted cup (5 µl)
 4. Safety lancet
 5. Alcohol swab
 6. Instructions for use
 7. Summarized instructions for use

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- 3** **EN** Next, look at the expiration date on the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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- 4** **EN** Open the foil pouch and look for the following:
1. Result window
 2. Specimen well
 3. Assay diluent well
- Then, label the device with the patient identifier.

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ENGLISH

INTENDED USE

NxTek™ Eliminate Malaria Pf is an ultra sensitive, qualitative *in vitro* test for the detection of histidine-rich protein 2 (HRP2) antigen of *Plasmodium falciparum* malaria in human whole blood. The test is intended for professional use and the point-of care testing to aid in the diagnosis of malaria infection.

REF 05FK142

MATERIALS PROVIDED

- Test device
- Assay diluent
- Inverted cup (5 µl)
- Safety lancet
- Alcohol swab

KIT STORAGE AND STABILITY

1. The test kit should be stored at a temperature between 1 °C and 30 °C. Do not expose the kit or its components to freezing or high temperatures.
Note : When stored at refrigerator, all kit components must be brought to room temperature (15 - 30 °C) minimum 30 minutes prior to performing the test. Do not open the pouch whilst components come to room temperature.
2. The test device is sensitive to both heat and humidity. Perform the test immediately after opening the foil pouch, to avoid exposure to humidity and moisture.
3. Do not use test kit beyond its expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

SPECIMEN COLLECTION USING A LANCET

1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the fingertip then prick the lateral side of the finger with the safety lancet provided. Immediately, safely dispose of the lancet.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. Using a new disposable inverted cup (5 µl) provided, dip the circular end of the inverted cup into the blood specimen.
5. The specimen collected must be used immediately. The specimen collected cannot be stored.

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© 2021 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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

REF 05FK142

TEST PROCEDURE / PROCÉDURE DE TEST

- EN** Dispense 5 µl of drawn whole blood specimen into the round specimen well.

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
- EN** Dispense all of the assay diluent from the diluent tube into the square well of test device.

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- EN** Read the test result after 20 minutes.

⚠ Caution: Do not read test results after 30 minutes.; reading after 30 minutes can yield erroneous results.

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

3. Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptosira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity. Reproducibility of NxTek™ Eliminate Malaria Pf test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), day to day (five days), between run (three different test operators), between sites (three sites), and lot-to-lot (3 different lots).

PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

1. Open the package and look for the following:
 1. Test device with desiccant in individual foil pouch
 2. Assay diluent
 3. Disposable inverted cup (5 µl)
 4. Sterile lancet
 5. Alcohol swab
 6. Instructions for use
 7. Summarized instructions for use

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3. First, carefully read the instructions for using NxTek™ Eliminate Malaria Pf test kit.

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
5. Next, look at the expiration date on the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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7. Open the foil pouch and look for the following:
 1. Result window
 2. Specimen well
 3. Assay diluent well
 Then, label the device with the patient identifier.

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Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptosira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity.

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4. Hepatitis A, measles, Influenza A & B, Tick borne encephalitis, Leishmania, Leptosira, M. tuberculosis, Schistosoma, Toxoplasma, Brucella or Human African Trypanosomiasis (HAT) specimens are not investigated for cross-reactivity.

Type d'échantillon	N°	Substance potentiellement interférente	Concentration du test
Substance antigénique	11	Armodiaquina	39,3 ng/ml
	12	Desethyl amodiaquina	751 ng/ml
	13	Artemether	540 ng/ml
	14	Lumefantrina	28 300 ng/ml
	15	Artesunato	29 644 ng/ml
	16	Dihydroartémisine	2 043 ng/ml
	17	Atovaquone	13 270 ng/ml
	8	Rubella	751 ng/ml
	9	West Nile Virus (WNV)	1 430 ng/ml
	18	Chloroquine	1 430 ng/ml
	19	Cinamylquina	751 ng/ml
	20	Dihydroartémisine	698 ng/ml
	21	Piperazina	730 ng/ml
	22	Dioxyprina	0,0069 ng/ml
	23	Maltoquina	3 279 ng/ml
	24	Primaquina	295 ng/ml
	25	Sulfadoxina	217 800 ng/ml
		Pyriméthamine	860 ng/ml

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SPECIMEN COLLECTION USING A LANCET

1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the fingertip then prick the lateral side of the finger with the sterile lancet provided. Immediately, safely dispose of the lancet.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. Using a new disposable inverted cup (5 µl) provided, dip the circular end of the inverted cup into the blood specimen.
5. The specimen collected must be used immediately. The specimen collected cannot be stored.

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