WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Wondfo HIV Self-Test WHO reference number: PQDx 0357-004-01

Wondfo HIV Self-Test with product codes W006P0058, W006P0059, and W006P0060, manufactured by Guangzhou Wondfo Biotech Co., Ltd, Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 13 July 2022.

Summary of WHO prequalification assessment for Wondfo HIV Self-Test¹

	Date	Outcome
Prequalification listing for the Wondfo HIV Self-Test	13 July 2022	Listed
Prequalification listing for the One Step HIV1/2 Whole Blood/Serum/Plasma Test	29 November 2018	Listed
Dossier assessment	26 April 2018	MR
Site inspection(s) of the quality management system	12 October 2018	MR
Product performance evaluation	First-quarter of 2018	MR

MR: Meets Requirements

*Change notification

In 2021, Guangzhou Wondfo Biotech Co., Ltd submitted a change notification for their prequalified product One Step HIV1/2 Whole Blood/Serum/Plasma Test to introduce a new configuration with an intended use specific for HIV self-testing(Wondfo HIV Self-Test). The new format was adapted from the corresponding professional use product (One Step HIV1/2 Whole Blood/Serum/Plasma Test), for which a WHO prequalification assessment has already taken place. Guangzhou Wondfo Biotech Co., Ltd generated additional data to meet requirements for self-testing as set out in the WHO Technical Specifications Series document TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and self-testing².

¹ Dossier assessment and product performance evaluation for the Wondfo HIV Self-Test were considered from the previous assessment of professional use product, One Step HIV1/2 Whole Blood/Serum/Plasma Test, which was prequalified in 2018. Based on the product dossier assessment and product performance evaluation, Wondfo HIV Self-Test meets WHO prequalification requirements. Please refer to the WHO Prequalification of Diagnostics Programme Public Report for One Step HIV1/2 Whole Blood/Serum/Plasma Test.

<https://extranet.who.int/pqweb/content/public-report-one-step-hiv12-whole-bloodserumplasma-test-pqdx-0357-004-00>

² <https://apps.who.int/iris/bitstream/handle/10665/251857/9789241511742-eng.pdf;jsessionid=153ABC9D88E7623A1AD1DF946A22B4C8?sequence=1>

Intended use

According to the intended use claim from Guangzhou Wondfo Biotech Co., Ltd, “*the Wondfo HIV Self-Test is a single-use in vitro diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as a self-test and/or by medical professionals. For in vitro diagnostic use only.*”

Assay description

According to the claim of assay description from Guangzhou Wondfo Biotech Co., Ltd, “*Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into their respective wells, they migrate along the device by capillary action. When the levels of HIV-1 or HIV-2 antibodies are at or above the target cut-off level, HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen(gp36/41). The complex is captured by HIV recombinant antigen (gp41and gp36) immobilized in the test region (T). This produces a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV 1/2 antibody are zero or below the target cut-off level, there is no visible colored band in the test region (T). This indicates a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C), if the test has been performed properly.*”

Test kit contents

Components		Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)		1	1x20	1x100
	Desiccant (pcs)		1	1x20	1x100
Accessories	Dropper (pcs)		1	1x20	1x100
	Buffer (vial)		1	1x20	1x100
	IFU (pcs)		1	1x20	1x100
	Blood Lancet for Single Use (pcs)		1	1x20	1x100
	Alcohol prep pad (pcs)		1	1x20	1x100
	Cotton swab (pcs)		1	1x20	1x100
	Disposal bag (pcs)		1	/	/

Materials required but not provided

- Timer
- Tissue
- Bandage.

Storage

The test kit must be stored at 2 -30 °C.

Shelf-life upon manufacture

24 months (test kit and buffer).

Warnings/limitations

Refer to manufacturer's instructions for use.

Labelling


- 1. Labels**
- 2. Instructions for use**

1. Labels


1.1 Labels for product code W006P0058

1. Test Cassette Foil Pouch for W006P0058

2. Accessories Pouch for W006P0058

3. Disposal Bag for W006P0058

Wondfo HIV Self-Test

LOT XXXXXXXXX
REF W006P0058

XXXX-XX-XX
XXXX-XX-XX

Rev. A2 Rel.: 2021/10/09

Contents:

- | | |
|--------------------|-------------------------------|
| 1 Test cassette | 1 Desiccant pouch |
| 1 Dropper | 1 Instruction for use (IFU) |
| 1 Buffer | 1 Cotton swab |
| 1 Alcohol prep pad | 1 Blood Lancet for Single Use |

 **Guangzhou Wondfo Biotech Co., Ltd.**
No. 8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: www.wondfo.com.cn

4. Foil Pouch for W006P0058

Rev. A3 Rel.: 2022/01/06

Wondfo

Leading POCT Manufacturer

Contents:

- | | |
|--------------------|-------------------------------|
| 1 Test cassette | 1 Desiccant pouch |
| 1 Dropper | 1 Instruction for use (IFU) |
| 1 Buffer | 1 Cotton swab |
| 1 Alcohol prep pad | 1 Blood Lancet for Single Use |
| 1 Disposal bag | |

REF W006P0058
LOT XXXXXXXXX
XXXX-XX-XX
XXXX-XX-XX

Wondfo

Leading POCT Manufacturer

Wondfo HIV Self-Test

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR SELF-TESTING USE


Wondfo
Leading POCT Manufacturer


Guangzhou Wondfo Biotech Co., Ltd.
Add: No. 8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063 Website: www.wondfo.com.cn

5. The Box for W006P0058


6. Round Sticker for W006P0058

1.2 Labels for product code W006P0059

1. Test Cassette Foil Pouch for W006P0059

2. Accessories Pouch for W006P0059

Wondfo HIV Self-Test

LOT XXXXXXXXX
REF W006P0059

XXXX-XX-XX
XXXX-XX-XX

Rev. A2 Rel.: 2021/10/09

Contents:

- | | |
|--------------------|-------------------------------|
| 1 Test cassette | 1 Desiccant pouch |
| 1 Dropper | 1 Instruction for use (IFU) |
| 1 Buffer | 1 Cotton swab |
| 1 Alcohol prep pad | 1 Blood Lancet for Single Use |

 **Guangzhou Wondfo Biotech Co., Ltd.**
No. 8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: www.wondfo.com.cn

3. Foil Pouch for W006P0059

Wondfo HIV Self-Test

20
Tests

CONTENTS:


- 20 Individual pouches, each containing:
- 1 Test cassette
- 1 Desiccant pouch
- 1 Dropper
- 1 Instruction for use (IFU)
- 1 Blood Lancet for Single Use
- 1 Alcohol prep pad
- 1 Buffer
- 1 Cotton Swab


IVD

Version: A3
Rel.: 2022/07/11

4. Secondary Packaging Label of W006P0059


5. The Box for W006P0059


6. Round Sticker for W006P0059

1.3 Labels for product code W006P0060

1. Test Cassette Foil Pouch for W006P0060

2. Accessories Pouch for W006P0060

Wondfo HIV Self-Test

LOT XXXXXXXXX
REF W006P0060

XXXX-XX-XX
XXXX-XX-XX

Rev. A2 Rel.: 2021/10/09





Contents:

- | | |
|--------------------|-------------------------------|
| 1 Test cassette | 1 Desiccant pouch |
| 1 Dropper | 1 Instruction for use (IFU) |
| 1 Buffer | 1 Cotton swab |
| 1 Alcohol prep pad | 1 Blood Lancet for Single Use |

 **Guangzhou Wondfo Biotech Co., Ltd.**
No. 8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: www.wondfo.com.cn

3. Foil Pouch for W006P0060

	    <p> ART NO.: C/NO.: QTY: PCS N.W: KGS G.W: KGS MEAS: 36×28×21 CM </p>		    <p> ART NO.: C/NO.: QTY: PCS N.W: KGS G.W: KGS MEAS: 36×28×21 CM </p>

4. The Carton (inner packaging container) for W006P0060

Wondfo HIV Self-Test

Wondfo
Leading POCT Manufacturer

IVD

Version: A2
Rel.: 2021/10/09

CONTENTS:

100 Individual pouches, each containing:

1 Test cassette	1 Desiccant pouch
1 Cotton swab	1 Alcohol prep pad
1 Blood Lancet for Single Use	1 Buffer
1 Instruction for use (IFU)	1 Dropper

REF W006P0060

LOT XXXXXXXXXX

 XXXX-XX-XX

 XXXX-XX-XX

 100

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: www.wondfo.com.cn

5. Secondary Packaging Label for W006P0060

2.0 Instructions for use³


³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

Wondfo HIV Self-Test

Instruction for use

- You must follow the test procedure carefully to get an accurate result.
- Sit in a clean, well-lit area and ensure you have all contents before beginning the test.
- For single-use only. Do not open foil pouch containing cassette until ready to test.

Contents

Watch the operation video

Materials may be required but not provided.

Preparation

Questions and answers

1. What is HIV?

The human immunodeficiency virus (HIV) targets cells of the immune system, called CD4 cells, which help the body respond to infection. Within the CD4 cell, HIV replicates and in turn, damages and destroys the cell. Without effective treatment of a combination of antiretroviral (ARV) drugs, the immune system will become weakened to the point that it can no longer fight infection and disease.

2. How is HIV transmitted?

HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:

- Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
- Blood transfusion of contaminated blood.
- Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments.
- From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.

If a person living with HIV is on ART (Antiretroviral therapy), which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

3. When do I need to test myself?

People who have been exposed to HIV or are at risk of HIV should seek testing. Most widely-used HIV diagnostic tests detect antibodies produced by the person as part of their immune response to fight HIV. In most cases, people develop antibodies to HIV within 28 days of infection. During this time, people may receive a false-negative result using an antibody test. While people at risk of HIV should test as soon as possible, negative test results within 28 days of exposure should be confirmed with additional testing no more than three months later.

4. What should I do if I get a positive result?

You need to follow up with a health care worker to get additional testing to confirm the result. At that time your local clinic and the provider will suggest the next steps that need to be taken.

Product information

CATALOG NO.

W006P0058 W006P0059 W006P0060

INTENDED USE

The Wondfo HIV Self-Test is a single-use *in vitro* diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as self-test and/or by medical professionals.

For *in vitro* diagnostic use only.

PRINCIPLES OF THE PROCEDURE

Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into the respective wells, they will migrate along the device by capillary action. The HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen (gp36/41), and the complex is then captured by the HIV recombinant antigen (gp41 and gp36) immobilized in the test region (T). When the levels of HIV-1 or HIV-2 antibodies are at or above the limit of detection (LOD) of the assay, it will produce a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV-1 or HIV-2 antibodies are zero or below the LOD, there is no visible colored band in the test region (T) indicating a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C).

WARNINGS AND CAUTIONS

- Do not** use if the test kit beyond expiration date.
 - Do not** use if the pouch is punctured or improperly sealed.
 - Do not** use for self-testing if you are under 12 years old.
 - Do not** use for self-testing if you have a bleeding disorder.
 - Do not** use for self-testing if you are already diagnosed as HIV positive.
 - Do not** open the pouch until you are ready to perform the test.
- Wash your hands and ensure that they are clean and dry before starting test. Adequate lighting is required to read the test results.

KIT CONTENT

There are 3 configurations of the test kits, 1 test/kit, 20 tests/kit and 100 tests/kit. The kit components are provided as below:

Components		Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)		1	1x20	1x100
	Desiccant (pcs)		1	1x20	1x100
	Dropper (pcs)		1	1x20	1x100
	Buffer (vial)		1	1x20	1x100
Accessories	IFU (pcs)		1	1x20	1x100
	Blood Lancet for Single Use (pcs)		1	1x20	1x100
	Alcohol prep pad (pcs)		1	1x20	1x100
	Cotton swab (pcs)		1	1x20	1x100
	Disposal bag (pcs)		1	/	/

One test strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

STORAGE AND STABILITY

- The test kit can be stored at 2-30 °C for 24 months.
- Use the test cassette within 1 hour after opening the pouch.
- Keep away from sunlight, moisture and heat.
- Use the kit at 10-30 °C.

LIMITATIONS OF THE PROCEDURE

- The test is designed for detecting human fingerstick whole blood.
- The test is limited to the qualitative detection of HIV-1 and HIV-2 antibodies.
- The assay procedure and result interpretation must be followed closely when testing. For optimal test performance, proper specimen collection is critical. Failure to follow the procedure may lead to inaccurate test results.
- False-negative results can occur in the following conditions:
 - Patients exposed to HIV less than 3 months.
 - Patients under HIV treatment (Antiretroviral therapy).
 - If the quantity of antibodies for HIV present in the specimen is below the detection limit of the assay.
- False-positive results can occur in the following conditions:
 - Patients have participated in a HIV vaccine clinical trial.
 - The presence of the control line only means that migration of added liquid occurred. It does not guarantee that:
 - The correct specimen has been used.
 - The specimen has been applied correctly.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but instead should be determined by a healthcare provider in conjunction with clinical findings and the results from other laboratory tests and evaluations. Results from the Wondfo HIV Self-Test should not be used as the sole basis for diagnosis.

PERFORMANCE CHARACTERISTICS

In the clinical study, 900 participants whose HIV status were unknown were given the Wondfo HIV Self-Test to test. The results were compared to the 4th generation laboratory test. The laboratory results shown that a total of 77 participants were HIV positive, 822 participants were HIV negative and 1 undetermined. A total of 43 participants (5 HIV positive, 37 HIV negative and 1 undetermined) were excluded from the performance analysis. The comparison of results was as follows:

- 95.8% of participants (69/72) correctly reported the result as positive. This means that 3 participants infected with HIV reported negative result. This is called a false negative.
- 99.6% of participants (782/785) correctly reported the result as negative. This means that 3 participants not infected with HIV reported positive result. This is called a false positive.
- 4.7% of participants (42/900) failed to obtain a result, 1 participant's HIV infection status was not confirmed during the clinical study so it was excluded from the analysis.


REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or selftesting, Geneva: World Health Organization; 2016.


SYMBOLS KEY

	In vitro diagnostic medical device		Consult instructions for use		Expiry date YYYY-MM-DD
	Contains sufficient for <n> tests		Date of manufacture YYYY-MM-DD		Keep dry
	Batch code		Temperature limit 2°C		Keep away from sunlight
	Manufacturer		Do not re-use		Product code/ Catalogue number
	Caution				

How to use the test kit (for fingerstick whole blood use)

How to read the test

DISPOSE

Throw away all contents in the normal trash or dispose according to local infectious waste disposal laws and regulations.

Manufacturer information

Guangzhou Wondfo Biotech Co., Ltd.

Address: No.8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R. China

Tel: +86-20-32296083 400-888-5268 (Toll Free)

Fax: +86-20-32296063

E-mail: global@wondfo.com.cn

Website: www.wondfo.com.cn

Please contact the manufacturer or your local distributor if you have any questions related to the product.

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RÉFÉRENCES

1. WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization ; 2019.
2. Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
3. WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or self-testing, Geneva: World Health Organization; 2016.

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