## WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

### Product: Determine HBsAg 2 WHO reference number: PQDx 0451-013-00

**Determine HBsAg 2** with product codes **7D2942, 7D2943 and 7D2943 SET**, manufactured by **Abbott Diagnostics Medical Co.Ltd<sup>1</sup>**, **Rest-of-World regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 2 September 2019.

### Summary of WHO prequalification assessment for Determine HBsAg

	Date	Outcome
<b>Prequalification listing</b>	2 September 2019	listed
<b>Dossier review</b>	N/A	N/A
<b>Site inspection(s) of the quality management system</b>	24 May 2019	MR
<b>Product performance evaluation</b>	Quarter 2 of 2019	MR

MR: Meet Requirements

N/A: Not Applicable

### Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Change of the name of the manufacturer from Alere Medical Co., Ltd. to Abbott Diagnostics Medical Co., Ltd. Products themselves are completely the same. Correction on regulatory version from CE-Marked to Rest-of-World.	23 January 2020
3.0	Updated the labelling of the Chase Buffer (7D2243) included in each SET product.	14 July 2023

<sup>1</sup> The company name changed from Alere Medical Co.,Ltd. to Abbott Diagnostics Medical Co., Ltd.

**Intended use:**

According to the claim of intended use by the manufacturer, "*Determine HBsAg 2 is an in vitro, visually read, qualitative immunoassay for the detection of Hepatitis B Surface Antigen (HBsAg) in human capillary and venous whole blood, plasma or serum. The test is intended as an aid in the diagnosis of HBV infection through detection of HBsAg from infected individuals. The test is not intended for screening blood donors. The test is for professional use only.*"

**Assay description:**

According to the manufacturer's claim, "*Determine HBsAg 2 is an immunochromatographic test for the qualitative detection of Hepatitis B Surface Antigen (HBsAg). A specimen is added to the sample pad. The specimen mixes with biotinylated anti-HBsAg mouse monoclonal antibodies and black particles coated with anti-HBsAg mouse monoclonal antibodies. This mixture migrates along the solid phase to the immobilized avidin at the patient bar. If HBsAg is present in the specimen, the antigen binds to the biotinylated anti-HBsAg mouse monoclonal antibodies and the black particles coated with anti-HBsAg mouse monoclonal antibodies. This complex binds to the immobilized avidin forming a black bar on the test strip. If HBsAg is not present, black particles flow past the patient bar and no black bar is formed on the test strip. To ensure assay validity, a procedural control bar is incorporated in the assay device on the test strip. The red control bar will appear.*"

**Test kit contents:**

Component	20 tests (product code 7D2942)	100 tests (product code 7D2943)	100 tests (product code 7D2943SET)
Determine HBsAg 2 test cards	2 cards of 10 tests/card	10 cards of 10 tests/card	10 cards of 10 testes/card
Chase buffer (7D2243)	NA	NA	1 bottle x 2.5ml
EDTA capillary tube (7D2222)	NA	NA	1 pack for 100 tests
Sterile lancets	NA	NA	1 box for 100 tests

**Items required but not provided:**

- Disposable gloves
- Timing device
- Micropipette capable of delivering 50 µL (not required for fingerstick method)

**Storage:**

The test kit should be stored at 2-30 °C.

**Shelf-life upon manufacture:**

18 months.

**Warnings/limitations:**

Refer to the instructions for use.

**Prioritization for prequalification**

Based on the established eligibility criteria, **Determine HBsAg 2** was given priority for the WHO prequalification assessment.

**Product dossier assessment**

In accordance with the WHO procedure for abridged prequalification assessment, Abbott Diagnostics Medical Co. Ltd (formerly called Alere Medical Co. Ltd) was not required to submit a product dossier for **Determine HBsAg 2** as per the *"Instructions for compilation of a product dossier"* (PQDx\_018 version 3). Notwithstanding, certain aspects of the product dossier previously submitted for stringent regulatory review were reviewed by an assessor during the site inspection.

**Manufacturing site inspection**

In accordance with the WHO procedure for abridged prequalification assessment, a shortened inspection with fewer inspectors was conducted at the site of manufacture ( 357 Matsuhidai, Matsudo-shi, 270-2214, Chiba-ken Japan) for **Determine HBsAg 2** from 26-27 November 2018 as per the *"Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics"* (PQDx\_014 version 4).

The inspection found that the manufacturer had an acceptable quality management system and good manufacturing practices in place that ensured the consistent manufacture of a product of good quality.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 24 May 2019.

Based on the site inspection and corrective action plan review, the quality management system for Determine HBsAg 2 meets WHO prequalification requirements.

## Product performance evaluation

Determine HBsAg 2 (Abbott Diagnostics Medical Co. Ltd, formerly Alere Medical Co.Ltd) was evaluated by WHO at the Virus Reference Department, Public Health England, London, UK, in the 2nd quarter of 2019 using plasma specimens. From this evaluation, we drew the following conclusions:

Determine HBsAg 2 is a rapid immunochromatographic assay for the detection of HBsAg in human serum, plasma or whole blood. A volume of 50 µL of specimen is needed to perform the assay. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities. Reading the results can be done visually.

In this limited evaluation on a panel of 514 clinically-derived serum/plasma specimens of European, African, Latin American and Asian origin, compared to the reference algorithm (Monolisa Ag HBs Plus [Bio-Rad Laboratories] and Hepanostika HBsAg Uniform II [bioMérieux]; followed by Hepanostika HBsAg Uniform II Confirmatory Assay [bioMérieux]), the following performance characteristics were estimated:

Performance characteristics in comparison with an agreed reference standard		
	Initial (95% CI)	Final (95% CI)
Sensitivity % (N=201)	100 (98.2-100)	100 (98.2-100)
Specificity % (N=313)	99.7 (98.2-100)	100 (98.8-100)
Invalid rate % (N=869)	0.12	
Inter-reader variability %	0	

In addition, analytical performance characteristics were assessed using commercially available and locally-made panels and the following results were obtained:

Additional performance characteristics	
Sensitivity during seroconversion on 5 seroconversion panels in comparison with a benchmark assay (Monolisa Ag HBs Plus)	Seroconversion sensitivity index of 0. Therefore Determine HBsAg 2 detected HBsAg on average 0 specimens earlier than the benchmark assay.
Analytical sensitivity on WHO Biological Preparation panel 03/262 and dilutions series of 3rd WHO International standard 12/226	Detected a concentration of 0.52 IU/ml (panel 03/262) and 0.125 IU/ml (dilutions of 12/226)

Analytical sensitivity on a low titer panel (PHA105, SeraCare)	10 of 14 specimens with low titer (< 1 IU/mL) were detected
Analytical sensitivity on the 1st WHO International Reference Panel for HBV Genotypes for HBsAg assays (PEI code number 6100/09)	All 15 HBsAg genotypes in the panel were detected.
Analytical sensitivity on HBV mutant panel	11 of 12 HBsAg mutants were detected to a concentration of 1 IU/mL
Lot to lot variation on a dilution panel in comparison with an agreed reference standard	Acceptable

<b>Key operational characteristics</b>	
Validated specimen types	Serum, plasma (EDTA), venous and capillary whole blood
Number of steps	1 with precision (serum/plasma) or 2 with/without precision (venous/capillary whole blood)
Time to result	15 minutes
Endpoint stability	15 minutes (read within a maximum of 30 minutes)
Internal QC	Yes, reagent addition control
In-use stability of reagents	The test should be used within 2 hours after removing the protective foil cover from the test strip.

## **Labelling**

- 1. Labels**
- 2. Instructions for use**

## **1. Labels**

**1.1 Determine HBsAg 2 Determine HBsAg (product code 7D2942) label.**

**EN**

Determine™ HBsAg 2 is an *in vitro*, visually read, qualitative immunoassay for the detection of Hepatitis B Surface Antigen (HBsAg) in human capillary and venous whole blood, plasma or serum. For professional use only.

7D2243 Chase Buffer is required for whole blood testing.

**Kit contains:**

Determine™ HBsAg 2 Test Cards: 2 cards (containing 10 tests/card)

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1/4

EN

Determine™ HBsAg 2 is an *in vitro*, visually read, qualitative immunoassay for the detection of Hepatitis B Surface Antigen (HBsAg) in human capillary and venous whole blood, plasma or serum.  
For professional use only.

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© 2019 Abbott. All rights reserved.  
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. July, 2019


4 571226 474570

241787/R2

DETERMINE™  
HBsAg 2 SET


DETERMINE™  
HBsAg 2 SET

DETERMINE™  
HBsAg 2 SET

### 1.3.2 Outside package images

#### Determine HBsAg 2 100T Box

### 1.3.4 Sterile lancet label (7D2233)

**Product Code:**  
7D2233

**medicore™**  
Nashville, TN  
1-800-327-8894

**LOT :** \_\_\_\_\_

**⏰ :** \_\_\_\_\_

**CE 0197** **STERILE R**

Single Use Only 

**Torjin Huorong Technology Co., Ltd.**  
Air Port Industrial Park No. 1 Floor, Building C,  
Jing'Er Road No. 226 Airport Economic Area  
300308 Torjin China  
Tel: 0086-22-5875333

**EC REP** **Yiheng Xu**  
34 Toronho House, Surrey  
Grove Road, London, UK SE16 7AJ  
Tel: 44+(0)2072370150 0044+70 2110813

Made in China

**100 Sterile Lancets**

- Fine gauge, tri-bevel tip for virtually painless sampling
- Consistent depth penetration
- Universal design fits almost all lancing devices

**28 gauge**

**Instructions for Use:**


**Do not use if safety seal has been broken. Instructions for use:**

1. Wash hands with soap and warm water. Dry thoroughly.
2. Select fingertip site, slightly off-center.
3. Remove protective lancet cover above (steps 1 & 2). If seal is already broken discard lancet.
4. Insert lancet into lancing device, use to puncture skin.
5. Recap lancet, discard in appropriate container.

### **1.3.5 Chase buffer labels (7D2243)**

**Chase buffer box label**

<p>Chase Buffer  Black</p> <p>Box</p> <p>Size: (W)22mm x (L)44mm x (H)56mm</p>	<p>PN: 1135945401</p>
---	-----------------------

REF 7D2243

1156253101  
April 2022Abbott Diagnostics Medical Co., Ltd.  
357 Matsuhidai, Matsudo-shi, Chiba,  
270 -2214, Japan  
Tel +81 47 311 5750Abbott Rapid Dx International Limited  
Parkmore East Business Park, Ballybritt,  
Galway, Ireland, H91 VK7E

DETERMINE™

CHASE BUFFER

Abbott

EN

## Intended Purpose

Chase Buffer is intended for use with Determine™ in vitro diagnostic rapid test devices. Chase Buffer enables the migration of the specimen and conjugate through the test device when testing whole blood specimens. It is specially formulated and optimised to ensure the Determine™ test performs correctly and controls assay flow. It is for professional use and near-patient testing, and not for self-testing.

## For full procedure refer to the diagnostic assay IFUs:

Determine™ HIV-1/2 (7D23)  
Determine™ HIV Early Detect (7D28)  
Determine™ HBsAg 2 (7D29)  
Determine™ Syphilis TP (7D24)  
Determine™ HBsAg (7D25)

When adding Chase Buffer to the sample pad, hold the bottle vertically. One drop to be used per test (one bottle of Chase Buffer can be used for 100 tests). Avoid drop with bubble, do not add buffer to specimen.

## USER TRAINING, QUALIFICATION AND EXPERIENCE REQUIREMENTS

This Chase Buffer is designed for use by professionals in both laboratory and near-patient settings including laboratory technologists. Suitable personnel shall have a consummate level of training and experience in specimen preparation, the running and interpretation of specimens, management of all supplies and waste disposal as per local regulations and generation of reports.

## Caution: Do not eat, drink, smoke or apply make-up while using this product.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

## Contents

1 Bottle (2.5 mL) Chase Buffer contains:  
**CHASE BUFFER** phosphate buffer solution (20mM pH7.4) containing preservatives Nipasept™ and A56620 (sarafloxacin).

## Storage Instruction

Recap the Chase Buffer to avoid evaporation or spillage and store at 2-30 °C until expiration date.

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## Advice Line/Infotelefon/Asistencia/Conseil/Assistenza/Rådgivning

For further information, please contact your distributor, or call Abbott Technical Specialists: Europe: Tel: +44 161 483 9032. Email: EME.TechSupport@abbott.com

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IVD

241609/R10

**EN** EDTA Capillary Tubes are to be used in the fingerstick collection of whole blood specimens for use in **Determine™** products. Refer to the package insert of the diagnostic assay for full procedure. Bring to room temperature before opening the cap. For professional use only.  
**Caution: Risk of glass breakage and injury. Handle with care. Store in a dry location.**

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© 2019 Abbott. All rights reserved.  
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Produced by:  
Drummond Scientific Co.  
Broomall, PA USA 19008

**Abbott Diagnostics  
Medical Co., Ltd.**  
357 Matsuhidai, Matsudo-shi,  
Chiba, 270-2214, Japan

Exp. Lot

4 571226 470176

## **2. Instructions for use<sup>2</sup>**

---

<sup>2</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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### Invalid/No válido Non valide/Inválido

 Control bar (red) Patient bar (black)

This instruction for use must be read carefully prior to use. All procedures must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this IFU.

### NAME AND INTENDED USE

Determine™ HBsAg 2 is an *in vitro*, visually read, qualitative immunoassay for the detection of Hepatitis B Surface Antigen (HBsAg) in human capillary and venous whole blood, plasma or serum. The test is intended as an aid in the diagnosis of HBV infection through detection of HBsAg from infected individuals. The test is not intended for screening blood donors. The test is for professional use only.

### SUMMARY AND EXPLANATION OF THE TEST

Hepatitis B virus (HBV), is a DNA virus transmitted percutaneously, sexually and perinatally. Worldwide it causes deaths from cirrhosis, liver failure and hepatocellular carcinoma<sup>1</sup>. HBsAg seropositivity is the first serological marker to appear in acute HBV infection<sup>2</sup>.

### BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Determine™ HBsAg 2 is an immunochromatographic test for the qualitative detection of Hepatitis B Surface Antigen (HBsAg).

A specimen is added to the sample pad. The specimen mixes with biotinylated anti-HBsAg mouse monoclonal antibodies and black particles coated with anti-HBsAg mouse monoclonal antibodies. This mixture migrates along the solid phase to the immobilized avidin at the patient bar.

If HBsAg is present in the specimen, the antigen binds to the biotinylated anti-HBsAg mouse monoclonal antibodies and the black particles coated with anti-HBsAg mouse monoclonal antibodies. This complex binds to the immobilized avidin forming a black bar on the test strip. If HBsAg is not present, black particles flow past the patient bar and no black bar is formed on the test strip.

To ensure assay validity, a procedural control bar is incorporated in the assay device on the test strip. The red control bar will appear.

### CONTENTS

- Determine™ HBsAg 2 20 Test (7D2942): 2 cards (10 tests/card)
- Determine™ HBsAg 2 100 Test (7D2943): 10 cards (10 tests/card)
- Determine™ HBsAg 2 SET (7D2943SET) 100 Test for testing whole blood samples
- Determine™ HBsAg 2 (7D2943) 10 cards (10 tests/card)
- Chase Buffer (7D2243) 1 bottle of 2.5 ml
- EDTA Capillary Tubes (7D2222) 1 pack for 100 tests
- Blood Lancet (sterilized) (7D2233) 1 box for 100 tests

### ACCESSORIES (required)

#### For testing Whole Blood Samples

Chase Buffer (7D2243) 1 Bottle (2.5 mL) containing phosphate buffered saline, preservative and antimicrobial agent

#### For testing Whole Blood Samples (fingerstick assay)

Single use sterile lancet  
EDTA capillary tube (7D2222)

### Materials Required But Not Provided

- Disposable gloves
- Timing device
- Micropipette capable of delivering 50 µL (not required for fingerstick method)
- Alcohol swab, gauze pad

### WARNING AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

Use the Chase Buffer bottle with attention to avoid contamination of the nozzle. Don't touch nozzle to the sample or sample pad.

Safety data sheet available for professional users on request.

### CAUTION:

When handling specimens and reagents, use appropriate biosafety practices<sup>4,5</sup>. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using suitable disinfectant, such as 0.5% sodium hypochlorite<sup>3,4</sup>.
- Decontaminate and dispose of all specimens, used test strips, and other potentially contaminated materials in accordance with local regulations<sup>3,4</sup>.

### STORAGE

Store Determine™ HBsAg 2 Test Cards and Chase Buffer at 2-30°C until expiration date. After open the Chase buffer, re-cap and store at 2-30°C to avoid evaporation or spillage.

- When handled and stored as directed, kit components are stable until expiration date. Do not use kit components beyond expiration date.
- Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.
- Do not use wet devices or damaged packages.

### SPECIMEN COLLECTION

#### Serum, Plasma, and Whole Blood Collection by Venipuncture

Use EDTA collection tubes for whole blood and plasma specimens.

- Collect human whole blood by aseptic venipuncture.
- To obtain serum, separate from the clot. To obtain plasma, separate from the packed cells. Separate specimens as soon as possible to avoid any hemolysis.

#### Whole Blood Collection by fingerstick<sup>4</sup> (See Fig.1)

#### Use EDTA capillary tube (7D2222)

**CAUTION:** Glass capillaries may be damaged during transportation or when in use. Handle with care in order to avoid injury when removing from the package as well as during use and during disposal.

### EDTA CAPILLARY TUBES

Before collecting a fingerstick specimen, place a capillary tube on a clean dry surface.

1. Choose the fingertip of the middle, ring, or index finger (whichever is the least callused).

Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.

2. Clean fingertip with alcohol; allow to air dry.
3. Position the hand palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose the lancet in an appropriate biohazard sharps container.
4. Wipe away the first drop of blood with a sterile gauze pad.
5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times.
6. Touch the tip of the capillary tube to the drop of blood. Avoid air bubbles. Fill the tube with whole blood up to between the two marked lines (50 µL).

### SPECIMEN STORAGE

- Store serum and plasma specimens at 2-8°C and run the test within 7 days of collection. If testing is delayed more than 7 days, the specimens should be frozen (-20°C or colder).
- Avoid repeated freeze/thaw cycles.
- If serum or plasma specimens show particulate matter or turbidity, centrifuge at 10,000g for 5 min at room temperature before sampling. Carefully take the 50 µL test sample from the supernatant. If a lipid layer is formed on the surface of the liquid, ensure that the sample is taken from the clear liquid below that layer.
- For whole blood collected by venipuncture, store at 2-8°C. Do not freeze whole blood specimens. Run the test within 2 days (48 hours) of collection. Mix the specimen well by gentle inversion of the tube immediately before testing.
- For whole blood collected by fingerstick, test immediately.

### TEST PROCEDURE

This test should be performed at 18 to 40°C.

1. Remove the desired number of test strips from the 10-test card by bending and tearing at the perforation.
  - To preserve the lot number which appears on the left side of the test card, remove individual test strips starting from the right side of the test card. The lot number and expiry date are not printed on the individual test strips.
2. Label or write the patient identification on the top white area of the device.
3. Remove the protective foil cover from each test.
  - After removing the protective foil cover from each test strip, start the assay within 2 hours.
4. For serum or plasma samples:

- a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
- b. Wait a minimum of 15 minutes (30 minutes maximum) from addition of the sample and read result.

For whole blood (venipuncture) samples:

- a. Apply 50 µL of sample (precision pipette) to the middle of the sample pad (marked by the arrow symbol).
- b. Wait one minute to allow the sample to be absorbed, and then apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.
- c. Wait a minimum of 15 minutes (30 minutes maximum) from addition of the sample and read result.

For whole blood (fingerstick) samples using EDTA capillary tubes:

- a. Place the capillary tube containing the blood sample to the middle of the sample pad (marked by the arrow symbol) at an upright (vertical) position.
- b. Wait until all the blood is transferred from the capillary tube to the sample pad. Then immediately apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.
  - Caution:** do not lift the capillary tube from the sample pad before all the blood has been transferred – a bubble may form which will prevent the complete transfer of sample and invalidate the test. It may take more than one minute for full transfer of the sample.
- c. Dispose the used capillary tube as biohazardous material according to local regulations.
- d. Wait a minimum of 15 minutes (30 minutes maximum) from addition of the sample and read result.

### NOTE:

- Running the test in high temperature/low humidity may affect the appearance of the patient bar. If the test strip is partially dried and it is difficult to read the result at 15 to 30 minutes, the test should be repeated using a new test strip and result read at 15 minutes. When the test strip is partially dried, it appears as mixed white spot and grayish area. At 40°C/ 40% RH or low humidity conditions, the test strip dries at 30 minutes.
- If serum or plasma sample does not flow or shows abnormal flow, such as stopping in the middle of the window, centrifuge the specimen and repeat the test with a new test strip.

### QUALITY CONTROL

To ensure assay validity, a procedural control system is incorporated in the device. This red control bar will appear in the window. If the control bar does not appear by assay completion, the test result is invalid. Repeat the test using a new test strip.


### INTERPRETATION OF RESULTS (See pictures)

#### NOTE:

- Interpret any visible red bar (even very faint) in the window as a valid result.
- The test result is valid even if the black patient bar appears lighter or darker than the red control bar. Any black patient bar, no matter how faint, is interpreted as reactive.
- Even if the color of control bar is dark red, it is valid.
- A test which gives very high background should be considered invalid.
- If an invalid result occurs repeatedly, or for technical assistance, contact your local distributor or call the Advice Line.
- When a whole blood venipuncture or whole blood fingerstick result is reactive with Determine™ HBsAg 2, it is recommended that the result is followed up with a plasma or serum Determine™ HBsAg 2 test or a laboratory reference test.

**Reactive (Two bars: Red bar and Black bar)**

Two bars: one red bar and one black bar appear in the window. The red bar corresponds to control bar and the black bar corresponds to patient bar, respectively.

**Non-reactive (One bar: Red bar)**

One red bar appears and no black bar appears in the window.

**INVALID (No bar or One bar: Black bar)**

If there is no red control bar in the window, and even if a black patient bar appears in the window, the result is invalid. Repeat the test using a new test strip.

**LIMITATIONS OF THE PROCEDURES**

- The Determine™ HBsAg 2 test is designed to detect Hepatitis B Surface Antigen (HBsAg) in human serum, plasma and whole blood. Other body fluids or pooled specimens may not give accurate results and should not be used.
- The intensity of the patient bar does not necessarily correlate with the titer of antigen in the specimen.
- No test provides absolute assurance that a sample does not contain low levels of HBsAg such as those present at a very early stage of infection.
- If the result is non-reactive at 15 minutes and reactive at 30 minutes, the patient has very low HBsAg below the detection limit.
- A non-reactive result at any time does not preclude the possibility of exposure to or infection with hepatitis B virus.
- The use of anticoagulants other than EDTA for whole blood or plasma specimens has not been validated for use with the Determine™ HBsAg 2 and may give incorrect results.
- For diagnostic purposes and to differentiate acute HBV infection from chronic HBV infection, the detection of HBsAg must be correlated with patient symptoms and other hepatitis B viral serological markers and nucleic acid testing (NAT).
- Biotin treatment higher than 5mg per day may lead to decreased patient bar intensity. Biotin concentrations up to 50 ng/mL in serum or plasma did not impact the sensitivity.

**PERFORMANCE CHARACTERISTICS**

**Sensitivity**

Sensitivity was evaluated by testing confirmed HBsAg positive specimens, commercial seroconversion panels, HBV genotype and serotype panels.

**1. HBsAg Positive Specimens**

A total of 437 confirmed HBsAg positive specimens were tested with Determine™ HBsAg 2 and a commercially available rapid HBsAg test (Table I). The sensitivity (95% CI) of Determine™ HBsAg 2 on this population of specimens was calculated to be 98.4% (96.7-99.4%). The non-reactive specimens with Determine™ HBsAg 2 were 0.06, 0.06, 0.09, 0.10, 0.13, 0.33 and 1.29 IU/mL by a quantitative HBsAg test kit.

**Table I: HBsAg Positive Specimens**

Types	Number of Specimens tested	Reactive by Determine™ HBsAg 2	Reactive by alternative rapid HBsAg test
HBsAg positive specimens	369	362	294
HBV genotype determined specimens*	51	51	43
HBsAg serotype determined specimens**	17	17	17
Plasma sensitivity (95% CI)	369	98.6% (364/369) (96.9-99.6%)	81.3% (300/369) (76.9-85.2%)
Serum sensitivity (95% CI)	68	97.1% (66/68) (89.8-99.6%)	79.4%(54/68) (67.9-88.3%)
Total Sensitivity (95% CI)	437	98.4% (430/437) (96.7-99.4%)	81.0% (354/437) (77.0-84.6%)

\* Genotypes: A, B, C, D, E, F, G, H, A/E, A/F, D/E, D/F and D/G  
 \*\* Serotypes: adw2, adr, ayw2, ayw3 and ayw4 including 15 samples of # 6100/09 PEI Hepatitis B virus genotype panel

**2. HBV seroconversion panels**

The sensitivity of Determine™ HBsAg 2 was evaluated using 32 sets of seroconversion panels; each including early seroconversion panel members. The results were compared with the results of a commercially available rapid immunochromatographic HBsAg test kit and a quantitative HBsAg test kit (Chemiluminescent micro-particle immunoassay (CMIA)). On all sets, Determine™ HBsAg 2 detected HBsAg earlier than a rapid immunochromatographic test kit. On 4 of 32 sets, Determine™ HBsAg 2 detected HBsAg earlier than a CMIA kit. Detection of 4 of 32 sets by Determine™ HBsAg 2 was delayed by 1 bleed date when compared to CMIA kit.

The first reactive date followed by continuously reactive results was compared with the first reactive date with the Determine™ HBsAg 2.

**3. Analytical sensitivity of HBsAg**

The analytical sensitivity of Determine™ HBsAg 2 was evaluated by testing WHO International Standard for HBsAg (NIBSC code 12/226). The results demonstrated that the test could detect a concentration of 0.1 IU/mL HBsAg.

**4. HBsAg Mutant Detection**

HBsAg mutant susceptibility was evaluated with Determine™ HBsAg 2. A panel consisting of 14 different recombinant HBsAg mutant panels was tested. The mutant panel consisted of following mutations, P120Q, T123A, T126N, T126S, Q129R, Q129H, Q129L, M133H, M133L, K141E, P142S, T143K, D144A and G145R. All 14 samples were detected with Determine™ HBsAg 2.

**Specificity**

A total of 1650 confirmed negative serum or plasma specimens were tested with the Determine™ HBsAg 2 and the specificity (95% CI) was determined (Table II). The specificity was 99.6% (99.1-99.8%).

**Table II: HBsAg Negative Specimens**

Population	Number of tested specimens	Determine™ HBsAg 2	Rapid HBsAg Test
		Non-reactive	Non-reactive
Seronegative specimens*	1027	1027	1027
Clinical specimens	213	209	210
Pregnant women	206	205	205
Disease states other than HBV and potentially interfering substances**	204	202	202
Plasma specificity (95% CI)	530	100% (530/530) (99.3-100%)	100%(530/530) (99.3-100%)
Serum specificity (95% CI)	1120	99.4% (1113/1120) (98.7-99.8%)	99.5%(1114/1120) (98.8-99.8%)
Total Specificity (95% CI)	1650	99.6%(1643/1650) (99.1-99.8%)	99.6%(1644/1650) (99.2-99.9%)

\* Including specimens collected in USA (1240) and Africa (20)  
 \*\*Rheumatoid factor (15/15), antinuclear antibody (10/10), systemic lupus erythematosus (10/10), high cholesterol (9/9), high total protein (3/3), high IgM (10/10), high IgG (4/4), human anti-mouse IgG (15/15), other infections - Hepatitis A virus (6/6), Hepatitis C virus (15/15), HTLV I (9/9), Cytomegalovirus (10/10), Toxoplasma IgG (10/10), Syphilis (8/8), Herpes simplex virus 1/2 (10/10) and Epstein-Barr virus (9/10), Flu vaccinated patients (10/10), non-virus liver disease and dialysis patient specimens (19/20)

**Sample type**

All specimen matrices (serum, plasma, whole blood venipuncture and whole blood fingerstick) were tested.

**Table III: HBsAg sensitivity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens**

Population	No. of matched specimens tested	Type of Specimens and No. of reactive or non-reactive by Determine™ HBsAg 2				
		Serum	Plasma	Whole blood venipuncture	Whole blood fingerstick**	
HBsAg positive	145	Reactive	142	143	141	141
		Non-reactive	3	2	4	4
		Sensitivity (95% CI)	97.9% (94.1-99.6%)	98.6% (95.1-99.8%)	97.2%*** (93.1-99.2%)	97.2%*** (93.1-99.2%)
HBsAg negative	203	Reactive	1	0	0 <3**>	0 <3**>
		Non-reactive	202	202*	203 <200***>	203 <200***>
		Specificity (95% CI)	99.5% (97.3-100%)	100% (98.2-100%)	100% <98.5%***> (98.2-100%)	100% <98.5%***> (98.2-100%)

\* There was no plasma result for 1 study subject.  
 \*\* EDTA capillary tubes were used for 225 specimens. MICROSAFE® Tubes were used for 123 specimens. Prospective multiple (matched) specimens of whole blood, serum and plasma from 348 individuals (145 HBsAg positive) in European countries were tested with Determine™ HBsAg 2 (Table III). Whole blood initially reactive results were adjudicated utilizing the serum or plasma results. The true negative but initially reactive whole blood results are listed in parentheses <> in Table III.

The results obtained from all specimen matrices showed correlation, except for the following. One positive specimen (0.27 IU/mL) was reactive with plasma but non-reactive with serum, venipuncture whole blood and fingerstick.

One positive specimen (0.14 IU/mL) was reactive with serum and plasma but non-reactive with venipuncture whole blood and fingerstick. One negative specimen was non-reactive with plasma, venipuncture whole blood and fingerstick but reactive with serum. Three negative venipuncture whole blood specimens were reactive. Three negative fingerstick specimens were reactive.

All whole blood true negative but reactive specimens were adjudicated utilizing the serum or plasma results; all were non-reactive on serum and plasma testing; the final assessment is that the samples were non-reactive.

Two seropositive specimens (0.08 IU/ml and 0.14 IU/ml) were non-reactive with all specimen matrices at 15 minutes. At 30 minutes, results for serum and plasma were reactive.

\*\*\* When using the cut-off 0.13 IU/mL, the sensitivity at 15 and at 30 minutes for fingerstick and venous whole blood samples was 97.9%.

**Advice Line (See back page)**

**BIBLIOGRAPHY (See back page)**

**The manufacturing process produces different lot numbers for the kit and test cards; these lot numbers are traceable.**

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**BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/BIBLIOGRAFIA**

1. Yun-Fan Liaw, Chia-Ming Chu (2009) Hepatitis B virus infection. Lancet 373:582-592
2. Jules L. Dienstag (2008) Hepatitis B Virus Infection. The New England Journal of Medicine 259:1486-1500
3. Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-6th GP42-A6 September 2008
4. Clinical and Laboratory Standards Institute, Clinical Laboratory Waste Management: Approved Guideline-Third Edition, GP05-A3. January 2011
5. EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1-5, R1-R3, A1-A24.

**Advice Line**

For further information, please contact your distributor, or call to one of the following Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+ 44 (0) 161 483 9032	EMEproductsupport@abbott.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@abbott.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCISproductsupport@abbott.com
Latin America	+ (57) 2 661 8797	LAPproductsupport@abbott.com

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