

WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Abbott RealTime High Risk HPV
WHO reference number: PQDx 0455-180-00

Abbott RealTime High Risk HPV with product codes **02N09-092** and **02N09-080**, manufactured by **Abbott GmbH (formerly called Abbott GmbH & Co.KG)**, CE Mark **regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 10 October 2019.

Summary of WHO prequalification assessment for Abbott RealTime High Risk HPV

	Date	Outcome
Prequalification listing	10 October 2019	listed
Dossier review	N/A	N/A
Site inspection(s) of quality management system	23 September 2019	MR
Product performance evaluation	N/A	N/A

MR: Meet Requirements

N/A: Not Applicable

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Updated the labelling to refer to its new legal entity Name, Abbott GmbH.	20 October 2021

Intended use:

According to the claim of intended use from Abbott GmbH & Co.KG *“The Abbott RealTime High Risk HPV is a qualitative in vitro test for the detection of DNA from 14 high risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in clinical specimens. The assay specifically identifies HPV genotypes 16 and 18 while concurrently detecting the other high risk genotypes at clinically relevant infection levels.*

The Abbott RealTime High Risk HPV is indicated:

a) To screen patients with ASC-US (atypical squamous cells of undetermined significance) cervical cytology results to determine the need for referral to colposcopy. The results of this test are not intended to prevent women from proceeding to colposcopy.

b) To be used with cervical cytology to adjunctively screen to assess the presence or absence of high risk HPV genotypes.

c) To be used as a first-line primary screening test to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease.

d) To assess the presence or absence of HPV genotypes 16 and 18 to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease with or without cervical cytology.

The results from the Abbott RealTime High Risk HPV, together with the physician’s assessment of cytology history, other risk factors, and professional guidelines, may be used to guide patient management.”

Assay description:

According to the claim of assay description from Abbott GmbH & Co.KG *“The Abbott RealTime HR HPV assay uses the Abbott m2000sp instrument, the Abbott m24sp instrument, or the manual sample preparation method for processing samples and the Abbott m2000rt instrument for amplification and detection. A primer mix consisting of 3 forward primers and 2 reverse primers targeting a conserved L1 region is used to amplify HPV targets. Signal for 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) is generated with the use of fluorescent labeled probes. Internal Control (IC) amplicons are generated with a primer set targeting an endogenous human beta globin sequence and are detected with the IC specific probe. The Abbott RealTime HR HPV assay detects the endogenous human beta globin sequence as sample validity control for cell adequacy, sample extraction and amplification efficiency. Probes for HPV 16, HPV 18, non-HPV 16/18 genotypes (Other HR HPV) and IC are labeled with different fluorophores allowing their signals to be distinguishable in a single reaction.*

Test kit contents:

Component	96 tests (product code 2N09-092)
Abbott RealTime High Risk HPV Amplification Reagent Kit	
Amplification reagent pack	4 packs x 24 tests/ pack
DNA polymerase in a Buffered solution with stabilizers	0.070 mL x 1 bottle
HPV oligonucleotide reagent	0.502 mL x 1 bottle
Activation reagent	0.778 mL x 1 bottle

Items required but not provided:

Component	Description
Abbott RealTime High Risk HPV Control Kit	Product code 2N09-80
Abbott RealTime High Risk HPV Negative Control	0.5 mL x 12 vials
Abbott RealTime High Risk HPV Positive Control	0.5 mL x 12 vials
Materials for manual sample preparation (Assay Protocol I)	
Abbott <i>m</i> Sample Preparation System _{DNA} for RealTime High Risk HPV	3N92
Abbott Optical Adhesive Covers	04J71-75
Abbott Adhesive Cover Applicator	9K32-01
Abbott 96-Well Optical Reaction Plate	04J71-70
Abbott Splash-Free Support Base	09K31-01
Calibrated precision pipettes	10 µL to 1 000 µL
Aerosol barrier pipette tips	20 µL to 1 000µL
Single-use DNase-free tube or container	General laboratory material
Materials for Abbott m24sp (Assay Protocol II)	

Abbott m24sp instrument containing the scripts necessary to run the Abbott RealTime HR HPV assay (m24sp Database v 3.0 or higher)	50-148470 or higher
Abbott <i>m</i> Sample Preparation System _{DNA}	06K12-24
Calibrated precision pipettes	10 µL to 1 000 µL
Aerosol barrier pipette tips	20 µL to 1 000µL
Sample input tubes	General laboratory equipment
1000 µL disposable tips	04J71-10
200 µL disposable tips	04J71-17
Vortex mixer	General laboratory equipment
USP grade 190 to 200 proof ethanol (95 to 100% ethanol: Do not use ethanol that contains denaturants.)	General laboratory material
Abbott Optical Adhesive Covers	04J71-75
Abbott Adhesive Cover Applicator	9K32-01
Abbott 96-Deep-Well Plate	04J71-30
Abbott Splash-Free Support Base	09K31-01
13 mm Sample Racks	04J72-82
1.5 mL Reaction Vessels and Output Tubes (1.5 mL screw top microfuge tubes and caps)	4J71-50
Abbott 96-Well Optical Reaction Plate	04J71-70
Abbott Splash-Free Support Base	09K31-01
Calibrated precision pipettes	10 µL to 1 000 µL
Aerosol barrier pipette tips	20 µL to 1 000µL
Single-use DNase-free tube or container	
Materials for Abbott m2000sp (Assay Protocol III)	
Abbott m2000sp instrument with Software Version 3.0 or higher	50-148393 or higher
Abbott <i>m</i> Sample Preparation System _{DNA}	06K12-24
5 mL Reaction Vessels	4J71-20
Calibrated precision pipettes	10 µL to 1 000 µL
Aerosol barrier pipette tips	20 µL to 1 000µL
Sample input tubes	General laboratory material
1000 µL disposable tips	04J71-10
200 µL disposable tips	04J71-17
Vortex mixer	

USP grade 190 to 200 proof ethanol (95 to 100% ethanol): Do not use ethanol that contains denaturants.	General laboratory material
Abbott Optical Adhesive Covers	04J71-75
Abbott Adhesive Cover Applicator	9K32-01
Abbott Splash-Free Support Base	09K31-01
Master Mix Tube	04J71-80
200 mL Reagent Vessels	4J71-60
Abbott 96-Deep-Well Plate	04J71-30
Abbott 96-Well Optical Reaction Plate	04J71-70
Materials for Abbott m2000rt	
Abbott m2000rt instrument with Software Version 3.0 or higher	50-148392 or higher
Abbott m2000rt Optical Calibration Kit	4J71-93

Other materials

- Biological safety cabinet approved for working with infectious materials
- Sealable plastic bags
- DNase-free water†
- Microcentrifuge Tubes†
- Cotton Tip Applicators (Puritan or equivalent) †

NOTE: † These three items are used in the procedure for Monitoring the Laboratory for the presence of Contamination. Refer to the QUALITY CONTROL PROCEDURES section of the package insert.

Storage:

- Abbott RealTime High Risk HPV Amplification Reagent Kit (product code 02N09-092) must be stored at -25 to -15°C when not in use.
- Abbott RealTime High Risk HPV Control Kit (product code 2N09-80) must be stored at -10°C or colder
- Reagents and controls are shipped on dry ice.

Shelf-life upon manufacture:

18 months.

Warnings/limitations:

Refer to the latest version of instructions for use.

Prioritization for prequalification

Based on the established prioritization criteria, Abbott RealTime High Risk HPV was given priority for WHO prequalification assessment.

Product dossier assessment

In accordance with the WHO procedure for abridged prequalification assessment, Abbott GmbH & Co.KG was not required to submit a product dossier for the Abbott RealTime High Risk HPV as per the *“Instructions for compilation of a product dossier”* (PQDx_018 version 3). Notwithstanding, certain aspects of the product dossier previously submitted for stringent regulatory review were reviewed by an assessor during the site inspection.

Commitment for prequalification: To provide a revised IFU by 31 March 2022 with amended intended use to state the population the product is intended to be used for.

Manufacturing site inspection

In accordance with the WHO procedure, an inspection of a manufacturing site(s) may be waived by WHO in writing under defined circumstances such as; a recent inspection with appropriate scope by a WHO-recognized national regulatory authority or by a Medical Devices Single Audit Program (MDSAP) participating Auditing Organisation as per the *“Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics”* (PQDx_014 version 4)..

The Abbott GmbH & Co. KG site located at Max-Planck-Ring 2, Wiesbaden, 65205, Germany was inspected by the Medical Device Single Audit Program (MDSAP), audit report (dated 4-7 June 2018). The site was found compliant and meet all requirements of ISO 13485: 2016.

An inspection of the site was waived in light of the MDSAP inspection.

Based on the MDSAP report, the quality management system for Abbott RealTime High Risk HPV meets WHO prequalification requirements.

Product performance evaluation

In accordance with the WHO procedure for prequalification assessment at the date of prioritization and given the fact that the Abbott RealTime High Risk HPV assay is used as the benchmark assay in the WHO evaluation protocol for HPV core antigen assays, it was decided that WHO will not conduct the performance evaluation of this assay as part of the prequalification assessment process.

Consequently, the laboratory evaluation of Abbott RealTime High Risk HPV was waived.


Labelling

- 1. Labels**
- 2. Instructions for use**

1. Labels

1.1 Abbott RealTime High Risk HPV Amplification Reagent Kit Label

List Number: 02N09-092



Amplification Reagent Kit

Abbott RealTime High Risk HPV

52-602801/R3

REF



LOT

GTIN



IVD

REF 02N09-092



Amplification Reagent Kit

Abbott RealTime High Risk HPV

(de) In-vitro-Test.

Inhalt:

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent Pack (4 Packungen, 24 Tests/Packung)

Jede Reagenzpackung enthält:

- 1 Fläschchen (0,070 ml) AmpliTaq Gold Enzym (5,4 bis 5,9 Einheiten/µl) in einer gepufferten Lösung mit Stabilisatoren.
- 1 Fläschchen (0,502 ml) HPV Oligonucleotidreagenz. < 0,1 % synthetische Oligonucleotide und < 1 % dNTPs, in einer gepufferten Lösung mit einem Referenzfarbstoff. Konservierungsmittel: Natriumazid und 0,16 % ProClin 950.
- 1 Fläschchen (0,778 ml) Aktivierungsreagenz. 38 mmol/l Magnesiumchlorid in einer gepufferten Lösung. Konservierungsmittel: Natriumazid und 0,15 % ProClin 950.

ProClin und AmpliTaq Gold sind Eigentum der Rechteinhaber.

(fr) Test *in vitro*.

Composition :

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent Pack (4 coffrets, 24 tests/coffret)

Chaque coffret-réactifs contient :

- 1 flacon (0,070 ml) d'enzyme AmpliTaq Gold (5,4 à 5,9 unités/µl) dans une solution tampon avec des stabilisants.
- 1 flacon (0,502 ml) de réactif d'oligonucléotides HPV. < 0,1 % d'oligonucléotides synthétiques et < 1 % de dNTPs, dans une solution tampon avec un fluorochrome de référence. Conservateurs : azide de sodium et ProClin 950 à 0,16 %.
- 1 flacon (0,778 ml) de réactif d'activation. 38 mmol/l de chlorure de magnésium dans une solution tampon. Conservateurs : azide de sodium et ProClin 950 à 0,15 %.

ProClin et AmpliTaq Gold sont la propriété de leurs détenteurs respectifs.

(es) Análisis *in vitro*.

Contenido:

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent Pack (4 envases, 24 pruebas/envase)

Cada equipo de reactivos contiene:

- 1 frasco (0,070 ml) AmpliTaq Gold Enzyme (5,4 unidades/µl a 5,9 unidades/µl) en solución tamponada con estabilizantes.
- 1 frasco (0,502 ml) HPV Oligonucleotide Reagent (reactivo de oligonucleótidos del VPH). < 0,1% de oligonucleótidos sintéticos y < 1% de dNTPs, en solución tamponada con fluoróforo de referencia. Conservantes: azida sódica y ProClin 950 al 0,16%.
- 1 frasco (0,778 ml) Activation Reagent (reactivo de activación). 38 mmol/l de cloruro de magnesio en solución tamponada. Conservantes: azida sódica y ProClin 950 al 0,15%.

ProClin y AmpliTaq Gold están a nombre de su propietario.

(it) Test *in vitro*.

Contenuto:

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent Pack (4 confezioni, 24 test/confezione)

Ciascuna confezione del reagente contiene:

- 1 fialone (0,070 ml) di AmpliTaq Gold Enzyme (enzima AmpliTaq Gold, da 5,4 a 5,9 unità/µl) in una soluzione tamponata con stabilizzanti.
- 1 fialone (0,502 ml) di HPV Oligonucleotide Reagent (reagente di oligonucleotidi dell'HPV). Oligonucleotidi sintetici < 0,1% e dNTP < 1% in una soluzione tamponata con un colorante di riferimento. Conservanti: sodio azoturo e ProClin 950 allo 0,16%.
- 1 fialone (0,778 ml) di Activation Reagent (reagente di attivazione). 38 mM di cloruro di magnesio in soluzione tamponata. Conservanti: sodio azoturo e ProClin 950 allo 0,15%.

ProClin e AmpliTaq Gold sono proprietà dei relativi titolari.

Abbott RealTime
High Risk HPV

REF 2N09

IVD

24 Tests

AMPLIFICATION REAGENT PACK



LOT

 Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

51-602328/R3

Top Edge

Abbott RealTime High Risk HPV


REF 02N09

IVD


AMPLIFICATION REAGENT PACK



24 Tests



LOT

 Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

51-602806/R3

1.2 Abbott RealTime High Risk HPV Amplification Reagent Pack Label

List Number: 02N09-90

ProCin et AmpliTaQ Gold sont la propriété de leurs détenteurs respectifs.

• 1 flacon (0,778 ml) de réactif d'activation. Conservateur: azide de sodium et ProCin 950 à 0,15%.

• 1 flacon (0,502 ml) de réactif d'oligonucléotides synthétiques. Conservateur: azide de sodium et ProCin 950 à 0,16%.

• 1 flacon (0,502 ml) d'HPV Oligonucleotide Reagent (réactif de oligonucleotides synthétiques). Conservateur: azide de sodium et ProCin 950 à 0,16%.

• 1 flacon (0,070 ml) AmpliTaQ Gold Enzyme (enzyma AmpliTaQ Gold de 5,4 à 5,9 unités/µl) dans une solution tampon avec des stabilisants.

Chaque coffret de réactifs contient :

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent (4 emballages, 24 tests/emballage)

Contenu:

(pt) Para uso em diagnóstico *in vitro*. Abbott RealTime High Risk HPV é um ensaio qualitativo *in vitro* para a detecção de ADN de 14 genótipos de papilomavírus humano (HPV) de alto risco 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 e 68 em amostras clínicas.

ProCin y AmpliTaQ Gold están a nombre de sus propietarios.

(fr) Pour diagnostic *in vitro*, Abbott RealTime High Risk HPV est un dosage qualitatif *in vitro* pour la détection de l'ADN des 14 génotypes de papillomavirus humain (HPV) à haut risque dans les échantillons cliniques.

Composition :

• 1 flacon (0,778 ml) de réactif d'activation. Conservateur: azide de sodium et ProCin 950.

• 1 flacon (0,502 ml) de réactif d'oligonucléotides synthétiques. Conservateur: azide de sodium et ProCin 950 à 0,16%.

• 1 flacon (0,502 ml) d'HPV Oligonucleotide Reagent (réactif de oligonucleotides synthétiques). Conservateur: azide de sodium et ProCin 950 à 0,16%.

• 1 flacon (0,070 ml) AmpliTaQ Gold Enzyme (5,4 à 5,9 unités/µl) en solution tamponnée avec des stabilisants.

Chaque coffret de réactifs contient :

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent (4 emballages, 24 tests/emballage)

Contenu:

(en) Para uso em diagnóstico *in vitro*. Abbott RealTime High Risk HPV é um ensaio qualitativo *in vitro* para a detecção de DNA de 14 genótipos de papilomavírus humano (HPV) de alto risco 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 e 68 em amostras clínicas.

ProCin and AmpliTaQ Gold are property of their respective owners.

(de) *In-vitro*-Diagnostik. Abbott RealTime High Risk HPV ist ein qualitativer *In-vitro*-Test zum Nachweis von DNA der 14 Genotypen (HPV) in Patientenproben.

Inhalt:

• 1 Flasche (0,778 ml) Aktivierungsmittel. Natriumazid und 0,15% ProCin 950.

• 1 Flasche (0,502 ml) HPV Oligonucleotide Reagent. < 0,16% synthetische Oligonucleotide und < 0,16% dNTPs, in einer gepufferten Lösung mit einem Referenzfarbstoff.

• 1 Flasche (0,502 ml) HPV Oligonucleotide Reagent. < 0,16% synthetische Oligonucleotide und < 0,16% dNTPs, in einer gepufferten Lösung mit einem Referenzfarbstoff.

• 1 Flasche (0,070 ml) AmpliTaQ Gold Enzym (5,4 bis 5,9 Einheiten/µl) in einer gepufferten Lösung mit Stabilisatoren.

Jedes Reagenzpaket enthält:

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent (4 Packungen, 24 Tests/Packung)

Inhalt:

(fr) Pour diagnostic *in vitro*, Abbott RealTime High Risk HPV est un test qualitatif *in vitro* pour la détection de l'ADN des 14 génotypes de papillomavirus humain (HPV) à haut risque dans les échantillons cliniques.

Composition :

• 1 flacon (0,778 ml) de réactif d'activation. Conservateur: azide de sodium et ProCin 950 à 0,15%.

• 1 flacon (0,502 ml) de réactif d'oligonucléotides synthétiques. Conservateur: azide de sodium et ProCin 950 à 0,16%.

• 1 flacon (0,502 ml) d'HPV Oligonucleotide Reagent (réactif de oligonucleotides synthétiques). Conservateur: azide de sodium et ProCin 950 à 0,16%.

• 1 flacon (0,070 ml) AmpliTaQ Gold Enzyme (5,4 à 5,9 unités/µl) en solution tamponnée avec des stabilisants.

Chaque coffret de réactifs contient :

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent (4 emballages, 24 tests/emballage)

Contenu:

(en) For *In Vitro* Diagnostic Use. The Abbott RealTime High Risk HPV is a qualitative *in vitro* test for the detection of DNA from 14 high risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in clinical specimens.

Contents:

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent Pack (4 packs, 24 tests/pack)

Each Reagent Pack contains:

- 1 Bottle (0.070 mL) AmpliTaQ Gold Enzyme (5.4 to 5.9 Units/µL) in a buffered solution with stabilizers.
- 1 Bottle (0.502 mL) HPV Oligonucleotide Reagent. <0.1% synthetic oligonucleotides and <1% dNTPs, in a buffered solution with a reference dye. Preservatives: sodium azide and 0.16% ProCin 950.
- 1 Bottle (0.778 mL) Activation Reagent. 38mM magnesium chloride in a buffered solution. Preservatives: sodium azide and 0.15% ProCin 950.

ProCin and AmpliTaQ Gold are property of their respective owners.

(pt) Para utilização em diagnóstico *in vitro*. O Abbott RealTime High Risk HPV é um ensaio qualitativo *in vitro* para a detecção de ADN de 14 genótipos de papilomavírus humano (HPV) em amostras clínicas.

Conteúdo:

AMPLIFICATION REAGENT PACK Abbott RealTime High Risk HPV Amplification Reagent Pack (4 embalagens, 24 testes/embalagem)

Cada embalagem de reagentes contém:

- 1 Frasco (0,070 ml) de enzima AmpliTaQ Gold (5,4 a 5,9 unidades/µl) numa solução tamponada com estabilizadores.
- 1 Frasco (0,502 ml) de reagente de oligonucleotídios do HPV. <0,1% de oligonucleotídios sintéticos e <1% de dNTPs, numa solução tamponada com um corante de referência. Conservantes: azida sódica e 0,16% de ProCin 950.
- 1 Frasco (0,778 ml) de reagente de activação. 38mM de cloreto de magnésio numa solução tamponada. Conservantes: azida sódica e 0,15% de ProCin 950.

ProCin e AmpliTaQ Gold são propriedade dos respectivos titulares.

Amplification Reagent Kit



REF 2N09-90

IVD



51-6-02327/R4

GTIN

LOT

REF

Abbott RealTime High Risk HPV

Amplification Reagent Kit



H317, EUH032
P261, P272, P280, P302 + P352,
P333 + P313, P362 + P364, P501

1.3 Abbott RealTime High Risk HPV Control Kit Label

List Number: 2N09-80

ProClin est la propriété de son détenteur.

tampón conteniendo de l'ADN enterrador. Conservantes: azida de sodium et ProClin 950 a 0,15 % < 0,01 % d'ADN non infectueux avec du VPH et des séquences beta-globine dans une solution.

2. **CONTROL +** Abbott RealTime High Risk HPV Positive Control (12 flacons de 0,5 ml chacun). Tampón conteniendo de l'ADN enterrador. Conservantes: azida de sodium et ProClin 950 a 0,15 % (chacun). < 0,01 % d'ADN non infectueux avec une séquence beta-globine dans une solution.

1. **CONTROL -** Abbott RealTime High Risk HPV Negative Control (12 flacons de 0,5 ml chacun). DNA non infectueux con DNA carrier. Conservantes: azida de sodium et ProClin 950 a 0,15 %.

2. **CONTROL +** Abbott RealTime High Risk HPV Positive Control (12 flacons de 0,5 ml per flacon). DNA non infectueux con DNA carrier. Conservantes: azida de sodium et ProClin 950 a 0,15 %.

1. **CONTROL -** Abbott RealTime High Risk HPV Negative Control (12 flacons, 0,5 ml per flacon). DNA non infectuos con secuencia de beta-globina en solución tamporada con DNA portador. Conservantes: azida sódica y ProClin 950 al 0,15 %.

2. **CONTROL +** Control positivo Abbott RealTime High Risk HPV (12 frascos, 0,5 ml por frasco). < 0,01% de DNA no infeccioso con secuencia de beta-globina y ProClin 950 al 0,15%.

1. **CONTROL -** Control negativo Abbott RealTime High Risk HPV (12 frascos, 0,5 ml por frasco). < 0,01% de DNA no infeccioso con secuencia de beta-globina en solución tamporada con DNA portador. Conservantes: azida sódica y ProClin 950 al 0,15%.

Contenido:

(fr) Per uso diagnostic *in vitro*. Les Abbott RealTime High Risk HPV Controls sont utilisés pour établir la validité du test Abbott RealTime High Risk HPV lors de la détection de l'ADN du virus du papillome humain (VPH) à haut risque dans les échantillons cliniques.

(it) Per uso diagnostico *in vitro*. Gli Abbott RealTime High Risk HPV Controls vengono utilizzati per stabilire la validità della seduta analitica del dosaggio Abbott RealTime High Risk HPV nella rilevazione del DNA di papillomavirus umano (HPV) ad alto rischio in campioni clinici.

ProClin está a nombre de su propietario.

ProClin ist Eigentum des Rechteinhabers.

Träger-DNA, Konservierungsmittel: Natriumazid und 0,15 % ProClin 950.

2. **CONTROL +** Abbott RealTime High Risk HPV Positive Control (12 Flaschen, je 0,5 ml). < 0,01 % nicht infektiöse DNA mit HPV und Beta-Globin-Sequenzen in einer gepufferten Lösung mit Träger-DNA, Konservierungsmittel: Natriumazid und 0,15 % ProClin 950.

1. **CONTROL -** Abbott RealTime High Risk HPV Negative Control (12 Flaschen, je 0,5 ml). < 0,01 % nicht infektiöse DNA mit Beta-Globin-Sequenz in einer gepufferten Lösung mit Träger-DNA, Konservierungsmittel: Natriumazid und 0,15 % ProClin 950.

Inhalt:

(de) In-vitro-Diagnosikum. Die Abbott RealTime High Risk HPV Kontrollen dienen zur Sicherstellung der Testgültigkeit des Abbott RealTime High Risk HPV Assays beim Nachweis von DNA des humanen Hochrisiko-Papillomavirus (HPV) in Patientenproben.

(es) Para uso en diagnóstico *in vitro*. Abbott RealTime High Risk HPV Controls se utilizan para establecer la validez del procesamiento de los ensayos Abbott RealTime High Risk HPV en la detección del DNA del papilomavirus humano (VPH) de alto riesgo, en muestras clínicas.

ProClin is property of its owner.


Abbott RealTime High Risk HPV

Control Kit

(en) For *In Vitro* Diagnostic Use. The Abbott RealTime High Risk HPV Controls are used to establish run validity of the Abbott RealTime High Risk HPV assay when used for the detection of high risk human papillomavirus (HPV) DNA in clinical specimens.

Contents:

- CONTROL -** Abbott RealTime High Risk HPV Negative Control (12 vials, 0.5 mL per vial). <0.01% noninfectious DNA with Beta Globin sequence in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.
- CONTROL +** Abbott RealTime High Risk HPV Positive Control (12 vials, 0.5 mL per vial). <0.01% noninfectious DNA with HPV and Beta Globin sequences in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.




(pt) Para utilização *in vitro*. Os Abbott RealTime High Risk HPV Controls destinam-se a estabelecer a validade do ensaio Abbott RealTime High Risk HPV quando utilizado para a detecção de ADN do papilomavirus humano de alto risco (HPV) em amostras clínicas.

Conteúdo:

- CONTROL -** Abbott RealTime High Risk HPV Negative Control (12 frascos, 0,5 ml por frasco). <0,01% de ADN não-infeccioso com sequência de betaglobulina numa solução tamporada com ADN portador. Conservantes: azida sódica e 0,15% de ProClin 950.
- CONTROL +** Abbott RealTime High Risk HPV Positive Control (12 frascos, 0,5 ml por frasco). <0,01% de ADN não-infeccioso com sequências de HPV e de betaglobulina numa solução tamporada com ADN portador. Conservantes: azida sódica e 0,15% de ProClin 950.

ProClin é propriedade do respectivo titular.



Abbott RealTime High Risk HPV

Control Kit

H317, EUH032
P261, P272, P280, P302 + P352,
P333 + P313, P362 + P364, P501


Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany

www.abbottmolecular.com



1.4 Abbott RealTime High Risk HPV Negative Control Vial Label

List Number: 2N09Z



REF 2N09Z
IVD

0.5mL

Abbott RealTime
High Risk HPV



CONTROL -



LOT


51-602332/R3

Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden, Germany



1.5 Abbott RealTime High Risk HPV Positive Control Vial Label

List Number: 2N09A



REF 2N09A
IVD

0.5mL

Abbott RealTime
High Risk HPV



CONTROL +



LOT

51-602333/R3

 Abbott GmbH
Max-Eyth-Ring 2
65205 Wiesbaden, Germany

 **Abbott**

2. Instructions for use¹

¹ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages

2.1 Abbott RealTime High Risk HPV Control Kit Package Insert

List Number: 2N09-80

FINALIDAD DE USO

Para uso en diagnóstico *in vitro*. Abbott RealTime High Risk HPV Controls (controles) se utilizan para establecer la validez del procesamiento del ensayo Abbott RealTime High Risk HPV para la detección del DNA del papilomavirus humano (VPH) de alto riesgo en muestras clínicas.

CONTENUTO

- CONTROL-** Abbott RealTime High Risk HPV Negative Control [control negativo] (n° de ref.: 2N09Z) (12 frascos, 0,5 ml cada uno). < 0,01% de DNA no infeccioso con secuencia de beta-globina en solución tamponada con DNA portador. Conservantes: azida sódica y ProClin 950 al 0,15%.
- CONTROL+** Abbott RealTime High Risk HPV Positive Control [control positivo] (n° de ref.: 2N09A) (12 frascos, 0,5 ml cada uno). < 0,01% de DNA no infeccioso con secuencias de beta-globina y VPH en solución tamponada con DNA portador. Conservantes: azida sódica y ProClin 950 al 0,15%.

- Abbott RealTime High Risk HPV Control Kit (equipo de controles) sólo puede utilizarse con el ensayo Abbott RealTime High Risk HPV (n° de ref.: 2N09).


PRECAUCIONES

- IVD**
- No utilizar una vez transcurrida la fecha de caducidad.

Los componentes de Abbott RealTime High Risk HPV Control Kit (equipo de controles, n° de ref.: 2N09-80) contienen lo siguiente:

- 2-metil-2H-isotiazol-3-ona
- Azida sódica

Se aplican las siguientes advertencias:

	Atención	
EUH032	En contacto con ácidos libera gases muy tóxicos.	
H317	Puede provocar una reacción alérgica en la piel.	
P261	Evitar respirar la niebla/los vapores/el aerosol.	
P280	Llevar guantes/prendas/gafas de protección.	
P272	Las prendas de trabajo contaminadas no podrán sacarse del lugar de trabajo.	
P302+P352	EN CASO DE CONTACTO CON LA PIEL: lavar con agua abundante.	
P333+P313	En caso de irritación o erupción cutánea: consultar a un médico.	
P362+P364	Quitar las prendas contaminadas y lavarlas antes de volver a usarlas.	
P501	Eliminar el contenido/el recipiente conforme a las normativas locales.	



CONDICIONES PARA EL TRANSPORTE

Transportar con nieve carbónica.

ProClin está a nombre de su propietario.



Mayo 2020
© 2008, 2020 Abbott Molecular Inc.

www.abbottmolecular.com

FINALITÀ D'USO

Per uso diagnostico *in vitro*. Gli Abbott RealTime High Risk HPV Controls vengono utilizzati per stabilire la validità della seduta analitica del dosaggio Abbott RealTime High Risk HPV nella rilevazione del DNA di papilomavirus umano (HPV) ad alto rischio in campioni clinici.

CONTENUTO

- CONTROL-** Abbott RealTime High Risk HPV Negative Control (n. di listino 2N09Z) (12 flaconi, 0,5 ml per flacone). DNA non infettivo <0,01% con sequenza di beta globina in una soluzione tamponata con DNA carrier. Conservanti: sodio azoturo e ProClin 950 allo 0,15%.
- CONTROL+** Abbott RealTime High Risk HPV Positive Control (n. di listino 2N09A) (12 flaconi, 0,5 ml per flacone). DNA non infettivo <0,01% con sequenze di HPV e di beta globina in una soluzione tamponata con DNA carrier. Conservanti: sodio azoturo e ProClin 950 allo 0,15%.

- L'Abbott RealTime High Risk HPV Control Kit deve essere utilizzato solamente con il dosaggio Abbott RealTime High Risk HPV (n. di listino 2N09).


PRECAUZIONI

- IVD**
- Non usare oltre la data di scadenza.

I componenti di Abbott RealTime High Risk HPV Control Kit (n. di listino 2N09-80) contengono i seguenti componenti:

- 2-metil-2H-isotiazol-3-one
- sodio azoturo

Si applicano le seguenti avvertenze:

	Attenzione	
EUH032	A contatto con acidi libera gas molto tossici.	
H317	Può provocare una reazione allergica della pelle.	
P261	Evitare di respirare la nebbia/i vapori/aerosol.	
P280	Indossare guanti/indumenti protettivi/Proteggere gli occhi.	
P272	Gli indumenti da lavoro contaminati non devono essere portati fuori dal luogo di lavoro.	
P302+P352	IN CASO DI CONTATTO CON LA PELLE: lavare abbondantemente con acqua.	
P333+P313	In caso di irritazione o eruzione della pelle, consultare un medico.	
P362+P364	Togliere gli indumenti contaminati e lavarli prima di indossarli nuovamente.	
P501	Smaltire il contenuto/recipiente in conformità alla regolamentazione locale.	



CONDIZIONI DI SPEDIZIONE

Spedire su ghiaccio secco.

ProClin è proprietà del suo titolare.



Maggio 2020
© 2008, 2020 Abbott Molecular Inc.

www.abbottmolecular.com

FINALIDADE DE USO

Para utilização em diagnóstico *in vitro*. Os controlos Abbott RealTime High Risk HPV são utilizados para estabelecer a validade do ensaio Abbott RealTime High Risk HPV quando utilizado para a deteção de ADN do papilomavirus humano de alto risco (HPV) em amostras clínicas.

CONTEÚDO

- CONTROL-** Abbott RealTime High Risk HPV Negative Control (Nº de Lista 2N09Z) (12 frascos, 0,5 ml por frasco). <0,01% de ADN não infeccioso com sequência de betaglobina numa solução tamponada com ADN transportador. Conservantes: azida sódica e 0,15% de ProClin 950.
- CONTROL+** Abbott RealTime High Risk HPV Positive Control (Nº de Lista 2N09A) (12 frascos, 0,5 ml por frasco). <0,01% de ADN não infeccioso com sequências de HPV e de betaglobina numa solução tamponada com ADN transportador. Conservantes: azida sódica e 0,15% de ProClin 950.

- O Abbott RealTime High Risk HPV Control Kit só pode ser utilizado com o ensaio Abbott RealTime High Risk HPV (Nº de Lista 2N09).


PRECAUÇÕES

- IVD**
- Não utilizar após o final do prazo de validade.

Os componentes do Abbott RealTime High Risk HPV Control Kit (Nº de lista 2N09-80) contêm os seguintes componentes:

- 2-metil-2H-isotiazol-3-ona
- Azida sódica

Aplicam-se os seguintes avisos:


	Atenção	
EUH032	Em contacto com ácidos liberta gases muito tóxicos.	
H317	Pode provocar uma reação alérgica cutânea.	
P261	Evitar respirar as névoas/vapores/aerossóis.	
P280	Usar luvas de proteção/vestuário de proteção/proteção ocular.	
P272	A roupa de trabalho contaminada não pode sair do local de trabalho.	
P302+P352	SE ENTRAR EM CONTACTO COM A PELE: lavar abundantemente com água.	
P333+P313	Em caso de irritação ou erupção cutânea: consulte um médico.	
P362+P364	Retirar a roupa contaminada e lavá-la antes de a voltar a usar.	
P501	Eliminar o conteúdo/recipiente em conformidade com os regulamentos locais.	



CONDIÇÕES DE TRANSPORTE

Transportar em gelo seco.

ProClin é propriedade do respetivo titular.



Maio 2020
© 2008, 2020 Abbott Molecular Inc.

www.abbottmolecular.com

 **Abbott RealTime
High Risk HPV**



REF 2N09-80






G59198R04

C2N099



Read Highlighted Changes: Revised May 2020.
Bitte Änderungen beachten: Überarbeitet im Mai 2020.
Faire attention aux modifications: Révision de mai 2020.
Consulte las modificaciones marcadas: Revisado en mayo de 2020.
Fare attenzione alle modifiche: Revisione di maggio 2020.
Consultar as alterações assinaladas: Revisto em maio de 2020.

Controls

Key to symbols used Erläuterung der verwendeten Symbole / Légende des symboles utilisés / Clave de los símbolos utilizados / Legenda dei simboli utilizzati / Legenda dos símbolos utilizados			
REF	Reference Number / Bestellnummer / Référence / Número de referencia / Numero di listino / Número de referência	LOT	Lot Number / Chargenbezeichnung / Número de lot / Número de lote / Numero di lotto / Número de lote
IVD	<i>In Vitro</i> Diagnostic Medical Device / In-vitro-Diagnostikum / Dispositif médical de diagnostic <i>in vitro</i> / Producto sanitario para diagnóstico <i>in vitro</i> / Dispositivo medico-diagnostico <i>in vitro</i> / Dispositivo médico para diagnóstico <i>in vitro</i>	CONTROL-	Negative Control / Negative Kontrolle / Contrôle négatif / Control negativo / Controllo negativo / Controllo negativo
	Consult instructions for use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Consultare le istruzioni per l'uso / Consultar as instruções de utilização	CONTROL+	Positive Control / Positive Kontrolle / Contrôle positif / Control positivo / Controllo positivo / Controllo positivo
	Upper limit of temperature / Temperaturobergrenze / Conserver jusqu'à / Limite superior de temperatura / Limite di temperatura superiore / Limite superior de temperatura		Warning / Achtung / Mise en garde / Atención / Attenzione / Atenção
			Manufacturer / Hersteller / Fabricant / Fabricante / Fabricante / Fabricante

 **Abbott**

INTENDED USE

For *In Vitro* Diagnostic Use. The Abbott RealTime High Risk HPV Controls are used to establish run validity of the Abbott RealTime High Risk HPV assay when used for the detection of high risk human papillomavirus (HPV) DNA in clinical specimens.

CONTENTS

- CONTROL-** Abbott RealTime High Risk HPV Negative Control (List No. 2N09Z) (12 vials, 0.5 mL per vial). <0.01% noninfectious DNA with Beta Globin sequence in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.
 - CONTROL+** Abbott RealTime High Risk HPV Positive Control (List No. 2N09A) (12 vials, 0.5 mL per vial). <0.01% noninfectious DNA with HPV and Beta Globin sequences in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.
- The Abbott RealTime High Risk HPV Control Kit must only be used with the Abbott RealTime High Risk HPV assay (List No. 2N09).


PRECAUTIONS

- IVD**
- Do not use beyond expiration date.

Components of the Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80) contain the following components:

- 2-Methyl-2H-isothiazol-3-one
- Sodium azide

The following warnings apply:

	Warning	
EUH032	Contact with acids liberates very toxic gas.	
H317	May cause an allergic skin reaction.	
P261	Avoid breathing mist / vapours / spray.	
P280	Wear protective gloves / protective clothing / eye protection.	
P272	Contaminated work clothing should not be allowed out of the workplace.	
P302+P352	IF ON SKIN: Wash with plenty of water.	
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.	
P362+P364	Take off contaminated clothing and wash it before reuse.	
P501	Dispose of contents / container in accordance with local regulations.	



SHIPPING CONDITIONS

Ship on dry ice.

ProClin is property of its owner.



Mai 2020
© 2008, 2020 Abbott Molecular Inc.
www.abbottmolecular.com

VERWENDUNGSZWECK

Zur Verwendung als In-vitro-Diagnostikum. Die Abbott RealTime High Risk HPV Controls dienen zur Sicherstellung der Testgültigkeit des Abbott RealTime High Risk HPV Assays beim Nachweis von DNA des humanen Hochrisiko-Papillomavirus (HPV) in Patientenproben.

INHALT

- CONTROL-** Abbott RealTime High Risk HPV Negative Control (Best.-Nr. 2N09Z) (12 Fläschchen, je 0,5 ml). <0,01 % nicht infektiöse DNA mit Beta-Globin-Sequenz in einer gepufferten Lösung mit Träger-DNA. Konservierungsmittel: Natriumazid und 0,15 % ProClin 950.
 - CONTROL+** Abbott RealTime High Risk HPV Positive Control (Best.-Nr. 2N09A) (12 Fläschchen, je 0,5 ml). <0,01 % nicht infektiöse DNA mit HPV und Beta-Globin-Sequenzen in einer gepufferten Lösung mit Träger-DNA. Konservierungsmittel: Natriumazid und 0,15 % ProClin 950.
- Der Abbott RealTime High Risk HPV Control Kit darf nur mit dem Abbott RealTime High Risk HPV Assay (Best.-Nr. 2N09) verwendet werden.


VORSICHTSMASSNAHMEN

- IVD**
- Nicht über das Verfallsdatum hinaus verwenden.

Komponenten des Abbott RealTime High Risk HPV Control Kit (Best.-Nr. 2N09-80) enthalten die folgenden Bestandteile:

- 2-Methyl-2H-isothiazol-3-on
- Natriumazid

Es gelten die folgenden Gefahrenhinweise:


	Achtung	
EUH032	Entwickelt bei Berührung mit Säure sehr giftige Gase.	
H317	Kann allergische Hautreaktionen verursachen.	
P261	Einatmen von Nebel / Dampf / Aerosol vermeiden.	
P280	Schutzhandschuhe / Schutzkleidung / Augenschutz tragen.	
P272	Kontaminierte Arbeitskleidung nicht außerhalb des Arbeitsplatzes tragen.	
P302+P352	BEI BERÜHRUNG MIT DER HAUT: Mit viel Wasser waschen.	
P333+P313	Bei Hautreizung oder -ausschlag: Ärztlichen Rat einholen / ärztliche Hilfe hinzuziehen.	
P362+P364	Kontaminierte Kleidung ausziehen und vor erneutem Tragen waschen.	
P501	Inhalt / Behälter gemäß den geltenden gesetzlichen Vorschriften entsorgen.	



TRANSPORTBEDINGUNGEN

Auf Trockeneis versenden.

ProClin ist Eigentum des Rechteinhabers.



Mai 2020
© 2008, 2020 Abbott Molecular Inc.
www.abbottmolecular.com

DOMAINE D'APPLICATION

Pour diagnostic *in vitro*. Les Abbott RealTime High Risk HPV Controls sont utilisés pour établir la validité du test Abbott RealTime High Risk HPV lors de la détection de l'ADN du papillomavirus humain (HPV) à haut risque dans les échantillons cliniques.

COMPOSITION

- CONTROL-** Abbott RealTime High Risk HPV Negative Control (Réf. 2N09Z) (12 flacons de 0,5 ml chacun). < 0,01 % d'ADN non infectieux avec une séquence Beta-Globine dans une solution tampon contenant de l'ADN entraîneur. Conservateurs : azide de sodium et ProClin 950 à 0,15 %.
 - CONTROL+** Abbott RealTime High Risk HPV Positive Control (Réf. 2N09A) (12 flacons de 0,5 ml chacun). < 0,01 % d'ADN non infectieux avec de l'HPV et des séquences Beta-Globine dans une solution tampon contenant de l'ADN entraîneur. Conservateurs : azide de sodium et ProClin 950 à 0,15 %.
- L'Abbott RealTime High Risk HPV Control Kit ne doit être utilisé qu'avec le test Abbott RealTime High Risk HPV (Réf. 2N09).


PRECAUTIONS

- IVD**
- Ne pas les utiliser au-delà de leur date de péremption.

L'Abbott RealTime High Risk HPV Control Kit (Réf. 2N09-80) contient les composants suivants :

- 2-méthyl-2H-isothiazole-3-one
- Azide de sodium

Les mises en garde suivantes s'appliquent :

	Mise en garde	
EUH032	Au contact d'un acide, dégage un gaz très toxique.	
H317	Peut provoquer une allergie cutanée.	
P261	Eviter de respirer les brouillards / vapeurs / aérosols.	
P280	Porter des gants de protection / des vêtements de protection / un équipement de protection des yeux.	
P272	Les vêtements de travail contaminés ne devraient pas sortir du lieu de travail.	
P302+P352	EN CAS DE CONTACT AVEC LA PEAU : Laver abondamment à l'eau.	
P333+P313	En cas d'irritation ou d'éruption cutanée : Consulter un médecin.	
P362+P364	Enlever les vêtements contaminés et les laver avant réutilisation.	
P501	Eliminer le contenu / récipient conformément aux réglementations locales.	



CONDITIONS D'EXPEDITION

Expédier sur de la carboglace.

ProClin est la propriété de son détenteur.



Mai 2020
© 2008, 2020 Abbott Molecular Inc.
www.abbottmolecular.com



Abbott RealTime High Risk HPV

CE 
en

REF 2N09
G59272R05
B2N090

Read Highlighted Changes: Revised May 2020.

Customer Service

International: Call your Abbott Representative

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

 **Abbott**

Key to symbols used

REF

List Number

CONTROL -

Negative Control

IVD

In Vitro Diagnostic Medical Device

CONTROL +

Positive Control

LOT

Lot Number

AMPLIFICATION REAGENT PACK

Amplification Reagent Pack



Expiration Date



Store at -10°C or colder



Consult instructions for use



Manufacturer



Warning

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

CONTENTS

NAME	4
INTENDED USE	4
SUMMARY AND EXPLANATION OF THE TEST.....	4
BIOLOGICAL PRINCIPLES OF THE PROCEDURE.....	5
REAGENTS	8
WARNINGS AND PRECAUTIONS	9
REAGENT STORAGE AND HANDLING INSTRUCTIONS.....	13
INSTRUMENTS/METHODS.....	13
SPECIMEN COLLECTION AND HANDLING INSTRUCTIONS	14
ASSAY PROCEDURE.....	14
ASSAY PROTOCOL I: MANUAL SAMPLE PREPARATION METHOD AND <i>m2000rt</i> INSTRUMENT	19
ASSAY PROTOCOL II: <i>m24sp</i> AND <i>m2000rt</i> INSTRUMENTS	23
ASSAY PROTOCOL III: <i>m2000sp</i> AND <i>m2000rt</i> INSTRUMENTS	28
QUALITY CONTROL PROCEDURES.....	31
RESULTS.....	34
LIMITATIONS OF THE PROCEDURE.....	36
SPECIFIC PERFORMANCE CHARACTERISTICS.....	36
BIBLIOGRAPHY	49

NAME

Abbott RealTime High Risk HPV

INTENDED USE

The Abbott RealTime High Risk HPV is a qualitative *in vitro* test for the detection of DNA from 14 high risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in clinical specimens.

SUMMARY AND EXPLANATION OF THE TEST

HPV is a small, non-enveloped, double-stranded DNA virus (approximately 8,000 base pairs) that replicates in the nucleus of squamous epithelial cells and induces hyperproliferative lesions.¹ HPV infections are among the most common sexually transmitted infections.² Most HPV infections have a benign clinical consequence and are cleared spontaneously.³ However, persistent HPV infection may result in progression to cervical cancer.⁴⁻⁷ More than one hundred different HPV genotypes have been identified, among which over forty infect mucosal and genital epithelia.⁸ Genital HPV genotypes are generally classified into high risk (HR) and low risk (LR) groups based on their carcinogenic potential. HR HPV genotypes are associated with invasive cervical cancer or its immediate precursor (high-grade squamous intraepithelial lesion, cervical intraepithelial neoplasia or carcinoma *in situ*), whereas LR HPV genotypes induce benign lesion and are not associated with cervical cancer.⁹⁻¹² Approximately 70% of invasive cervical cancer cases worldwide are caused by HPV 16 and HPV 18.¹³ Infection by HPV 16 or HPV 18 is associated with higher risk of disease progression compared to other HR HPV genotypes.¹⁴ Compared with cervical screening methods identifying cytological abnormalities, molecular tests that specifically detect the presence of HR HPV DNA in cervical cells can potentially increase sensitivity and cost-effectiveness of cervical cancer screening programs.¹⁵⁻²⁰ Furthermore, HPV DNA tests can be effectively used in triaging patients with equivocal cytology, in post-therapeutic follow-up and in monitoring vaccine efficacy.²¹⁻²³

The Abbott RealTime HR HPV assay is a qualitative *in vitro* test that amplifies and detects HR HPV DNA in cervical cells collected in liquid media. The detection of fourteen HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) is achieved through a primer mix targeting a conserved region of HPV genomes and single-stranded DNA probes. The assay can differentiate between HPV 16, HPV 18 and non-HPV 16/18 genotypes (Other HR HPV).

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The Abbott RealTime HR HPV assay uses the Abbott *m2000sp* instrument, the Abbott *m24sp* instrument or the manual sample preparation method for processing samples and the Abbott *m2000rt* instrument for amplification and detection. A primer mix consisting of three forward primers and two reverse primers targeting a conserved L1 region is used to amplify HPV targets. Signal for fourteen HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) is generated with the use of fluorescent labeled probes. Internal Control (IC) amplicons are generated with a primer set targeting an endogenous human beta globin sequence and are detected with the IC specific probe. The Abbott RealTime HR HPV assay detects the endogenous human beta globin sequence as sample validity control for cell adequacy, sample extraction and amplification efficiency. Probes for HPV 16, HPV 18, non-HPV 16/18 genotypes (Other HR HPV) and IC are labeled with different fluorophores allowing their signals to be distinguishable in a single reaction.

Sample Preparation

The purpose of sample preparation is to extract, concentrate, and purify the target DNA molecules for amplification. The Abbott *mSample* Preparation System_{DNA} uses magnetic particle technology to capture nucleic acids and washes the particles to remove unbound sample components. The bound nucleic acids are eluted and are then ready for amplification.

NOTE: One *mSample* Preparation System_{DNA} kit is sufficient to complete 4 x 48 (192) HPV sample preparations.

Two automated instrument systems, the *m2000sp* or the *m24sp*, can be used to prepare samples for the Abbott RealTime HR HPV assay. The *m2000sp* provides automated sample eluate transfer and reaction assembly in the Abbott 96-Well Optical Reaction Plate, while the *m24sp* requires manual sample eluate transfer and reaction assembly.

Alternatively, samples can be prepared manually following the instructions in “Manual Sample Preparation Using the ABBOTT *m*Sample Preparation System^{DNA} for RealTime High Risk HPV” (List No. 3N92). The manual sample preparation method requires manual transfer of the eluted samples to a 96-Well Optical Reaction Plate and manual reaction assembly before amplification.

Reagent Preparation and Reaction Plate Assembly

The *m2000sp* combines the Abbott RealTime HR HPV Amplification Reagent components (HPV Oligonucleotide Reagent, AmpliTaq Gold Enzyme, and Activation Reagent). The *m2000sp* dispenses the resulting master mix to the 96-Well Optical Reaction Plate along with aliquots of the nucleic acid samples prepared by the *m2000sp*. The plate is ready, after manual application of the optical seal, for transfer to the *m2000rt*.

The *m24sp* users and manual sample preparation method users manually combine the Abbott RealTime HR HPV Amplification Reagent components to create the amplification master mix and transfer aliquots of the master mix and sample eluates to the 96-Well Optical Reaction Plate. The plate is ready, after manual application of the optical seal, for transfer to the *m2000rt*.

Amplification

During the amplification reaction on the *m2000rt*, the target DNA is amplified by AmpliTaq Gold polymerase enzyme in the presence of dNTPs and magnesium. The AmpliTaq Gold polymerase enzyme is a thermophilic enzyme that has been modified in its active site by a molecule that renders it inactive. When the enzyme is heated prior to the initiation of PCR, the inhibitory molecule is cleaved from the enzyme allowing it to regain its activity. In this way, the enzyme is only active at temperatures where specific DNA-DNA interactions occur. This greatly reduces non-specific PCR artifacts such as primer dimers. In the Abbott RealTime HR HPV assay, the AmpliTaq Gold enzyme is first activated at 92°C for 10 minutes. During each subsequent round of thermal cycling, a high temperature is used to melt double-stranded DNA strands apart, followed by a low temperature where primers anneal to their respective targets and are extended to generate double-stranded DNA products. Exponential amplification of the products is achieved through repeated cycling between high and low temperatures, resulting in a billion-fold or greater amplification of target sequences. Amplification of both targets (HPV and IC) takes place simultaneously in the same reaction.

The target sequence for the Abbott RealTime HR HPV assay is in the conserved L1 region of the HPV genomes. A primer mix consisting of three forward primers and two reverse primers is designed to hybridize to the consensus regions among HPV genotypes of approximately 150 bases. The IC target sequence is a region of 136 bases in the endogenous human beta globin gene.

Detection

During the last 38 cycles of amplification, in an additional reading step, the temperature is lowered further to allow fluorescence detection of amplification products as the HPV and IC probes anneal to their targets (referred to as real-time fluorescence detection). The HPV and IC probes are single-stranded DNA oligonucleotides modified with a fluorescent moiety covalently linked to one end of the probe and a quenching moiety to the other end. In the absence of HPV or IC target sequences, the probes adopt a series of random conformations, some of which bring the quencher close enough to the excited fluorophore to absorb its energy before it can be fluorescently emitted. When a probe binds to its complementary sequence in the target, the fluorophore and the quencher are held apart, allowing fluorescent emission and detection by the *m2000rt*.

Signal for fourteen HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) is generated with the use of fluorescent labeled probes. IC signal is generated with an IC specific probe. Probes for HPV 16, HPV 18, Other HR HPV and IC are labeled with different fluorophores allowing their distinct signals to be simultaneously detected and distinguishable in a single reaction. Signals for HPV 16, HPV 18, Other HR HPV, and IC are detected in VIC, NED, FAM, and Cy5 channels, respectively.

Assay Results

The Abbott RealTime HR HPV assay is a qualitative assay. Results are reported as detected or not detected. In addition, each detected signal (HPV 16, HPV 18, or Other HR HPV) is also listed in the reported result. Refer to the “**RESULTS**” section of the package insert for further details.

REAGENTS

The Abbott RealTime HR HPV assay consists of two kits:

- Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 2N09-90)
- Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)

Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 2N09-90)

AMPLIFICATION REAGENT PACK (4 packs, 24 tests/pack)

Each Reagent Pack contains:

- 1 Bottle (0.070 mL) AmpliTaq Gold Enzyme (5.4 to 5.9 Units/ μ L) in a buffered solution with stabilizers.
- 1 Bottle (0.502 mL) HPV Oligonucleotide Reagent. < 0.1% synthetic oligonucleotides and < 1% dNTPs, in a buffered solution with a reference dye. Preservatives: sodium azide and 0.16% ProClin 950.
- 1 Bottle (0.778 mL) Activation Reagent. 38 mM magnesium chloride in a buffered solution. Preservatives: sodium azide and 0.15% ProClin 950.

NOTE: The Abbott RealTime Reagent components (enzyme, oligonucleotide reagent, activation reagent) are intended for single-use only and unused reagents should be discarded.

Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)

Abbott RealTime High Risk HPV Negative Control

- **CONTROL -** (12 vials, 0.5 mL per vial)
< 0.01% noninfectious DNA with Beta Globin sequence in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.

Abbott RealTime High Risk HPV Positive Control

- **CONTROL +** (12 vials, 0.5 mL per vial)
< 0.01% noninfectious DNA with HPV and Beta Globin sequences in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.

NOTE: The Negative and Positive Controls are intended for single-use only and unused reagents should be discarded.

WARNINGS AND PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

Refer to the *m2000sp* (List No. 9K20), *m24sp* (List No. 3N09) and *m2000rt* (List No. 9K25) Operations Manuals, Hazards Section, and “Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime HR HPV” (List No. 3N92) for instructions on safety precautions.

- There are no human sourced materials in any of the Abbott RealTime HR HPV Amplification Reagents or Controls.
- This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled with appropriate biosafety practices. Wear disposable gloves while handling specimens and wash hands thoroughly afterwards. Use of protective eyewear is recommended.

The following warnings apply to the HPV Oligonucleotide Reagent, Activation Reagent and the controls.



Warning

Hazard-determining components of labeling:

2-Methyl-2H-isothiazol-3-one

Sodium Azide

H317

May cause an allergic skin reaction.

EUH032

Contact with acids liberates very toxic gas.

P261

Avoid breathing mist / vapours / spray.

P280

Wear protective gloves / protective clothing / eye protection / face protection.

P272

Contaminated work clothing should not be allowed out of the workplace.

P302+P352

IF ON SKIN: Wash with plenty of water.

P333+P313

If skin irritation or rash occurs: Get medical advice / attention.

P362+P364

Take off contaminated clothing and wash it before reuse.

P501

Dispose of contents / container in accordance with local regulations.

Specimen Collection and Handling Precautions

- Specimens collected in PreservCyt Solution (Cytoc Corporation) can be used with the Abbott RealTime HR HPV assay. Users must follow the manufacturer's instructions for collecting and handling cervical specimens in PreservCyt Solution.
- Specimens collected in SurePath Preservative Fluid (TriPath Imaging, Inc.) can be used with the Abbott RealTime HR HPV assay. Either the sample from the original SurePath collection vial or the remaining cell pellet sample obtained after slide preparation with TriPath Imaging PrepStain Slide Processor can be used for testing. Users must follow the manufacturer's instructions for collecting, handling and processing cervical specimens in SurePath Preservative Fluid.
- Specimens collected with the Abbott Cervi-Collect Specimen Collection Kit can be used with the Abbott RealTime HR HPV assay. Users must follow the instructions in the Abbott Cervi-Collect Specimen Collection Kit Package Insert (List No. 4N73) for collecting and handling cervical specimens.

Laboratory Precautions

- During preparation of samples, compliance with good laboratory practices is essential to minimize the risk of cross-contamination between samples as well as the inadvertent introduction of nucleases into samples during and after the extraction procedure. Proper aseptic technique should always be used when working with DNA.
- Work area and instrument platforms must be considered potential sources of contamination. Change gloves after having contact with potential contaminants (such as DNases, specimens, eluates, and/or amplified product) before handling unopened reagents, negative control, positive control, or specimens. Refer to the *m24sp*, *m2000sp* and *m2000rt* Operations Manuals for instrument cleaning procedures.
- Wear appropriate personal protective equipment at all times.
- Use powder-free gloves.
- To reduce the risk of nucleic acid contamination due to aerosols formed during pipetting, pipettes with aerosol barrier tips must be used for all pipetting. The length of the tip should be sufficient to prevent contamination of the pipette barrel. While pipetting, care should be taken to avoid touching the pipette barrel to the inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.

- Change aerosol barrier pipette tips between ALL manual liquid transfers.
- Clean and disinfect spills of specimens and reagents as stated in the following manuals: the *m24sp* Operations Manual, the *m2000sp* Operations Manual, the *m2000rt* Operations Manual, and “Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime High Risk HPV”.

Contamination Precautions

- Amplification reactions such as PCR are sensitive to accidental introduction of product from previous amplification reactions. Incorrect results could occur if either the clinical specimen or the reagents used become contaminated by accidental introduction of even a few molecules of amplification product. Measures to reduce the risk of contamination in the laboratory include physically separating the activities involved in performing PCR in compliance with good laboratory practices.
- The use of three dedicated areas within the laboratory is recommended for performing the Abbott RealTime HR HPV assay with the *m24sp* or manual sample preparation using the *mSample* Preparation System_{DNA} and the *m2000rt*:
 - The Reagent Preparation Area is dedicated to combining the Abbott RealTime HR HPV Amplification Reagent components to create the amplification master mix and transferring aliquots of the master mix to the 96-Well Optical Reaction Plate. Laboratory coats, pipettes, and pipette tips used in the Reagent Preparation Area must remain in this area and not be moved to either the Sample Preparation Area or the Amplification Area. Do not bring target or amplification product into the Reagent Preparation Area.
 - The Sample Preparation Area is dedicated to processing samples (specimens, Abbott RealTime HR HPV Controls) and to adding processed samples and controls to the 96-Well Optical Reaction Plate. All reagents used in the Sample Preparation Area should remain in this dedicated area at all times. Laboratory coats, pipettes, pipette tips, and vortex mixers used in the Sample Preparation Area must remain in this area and not be moved to either the Reagent Preparation Area or the Amplification Area. Do not bring amplification product into the Sample Preparation Area.
 - The Amplification Area is dedicated to the amplification and detection of amplified product. Laboratory coats and equipment used in the Amplification Area must remain in this area and not be moved to either the Reagent Preparation Area or the Sample Preparation Area.

- Only two dedicated areas, Sample Preparation Area and Amplification Area, are recommended when the *m2000sp* and *m2000rt* are used.
- If the *m2000sp* run is aborted, dispose of all commodities and reagents according to the *m2000sp* Operations Manual. If the *m24sp* run is aborted, dispose of all commodities and reagents (if not being reused) according to the *m24sp* Operations Manual. If the manual sample preparation procedure is incorrectly performed or is interrupted at any point so that the timing of the steps exceeds the recommended timing per the manual instructions, dispose of all commodities and reagents (if not being reused) according to the instructions in "Manual Sample Preparation Using the ABBOTT *m*Sample Preparation System_{DNA} for RealTime High Risk HPV".
- If the *m2000sp* master mix addition protocol is aborted after amplification reagents are added to the 96-Well Optical Reaction Plate, seal the 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the *m2000sp* Operations Manual, Hazards section, along with the gloves used to handle the plate. Do not import the test order onto the *m2000rt*. If manual preparation of the PCR reaction mix is aborted after amplification reagents are added to the 96-Well Optical Reaction Plate, seal the 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to laboratory guidelines, along with the gloves used to handle the plate.
- For all completed, interrupted or aborted *m2000rt* runs, dispose of the 96-Well Optical Reaction Plate in a sealable plastic bag according to the *m2000rt* Operations Manual along with the gloves used to handle the plate.
- **Autoclaving the sealed Reaction Plate will not degrade the amplified product and may contribute to the release of the amplified product by opening the sealed plate. The laboratory area can become contaminated with amplified product if the waste materials are not carefully handled and contained.**
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local, state, and federal regulations.^{24,25} All materials should be handled in a manner that minimizes the chance of potential contamination of the work area.


Contamination from External Deoxy-Uracil (dU)-Containing Amplified Product

- HPV amplification assays containing dU may cause contamination and inaccurate results in the Abbott RealTime HR HPV. When negative controls are persistently reactive or where contamination with dU-containing HPV amplified product is likely to have occurred, it is recommended that the laboratory uses a contamination control procedure. This procedure (List No. 2N09-66) is available through your Abbott representative.


REAGENT STORAGE AND HANDLING INSTRUCTIONS

NOTE: Care must be taken to separate the Abbott RealTime High Risk HPV Amplification Reagent Kit that is in use from direct contact with specimens and Abbott RealTime High Risk HPV Control Kit reagents.

Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 2N09-90)

- 
- 10°C
- The Abbott RealTime High Risk HPV Amplification Reagent Pack must be stored at -10°C or colder when not in use.
 - Reagents are shipped on dry ice.

Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)

- 
- 10°C
- The Abbott RealTime High Risk HPV Negative and Positive Controls must be stored at -10°C or colder.
 - Reagents are shipped on dry ice.

INSTRUMENTS/METHODS

The Abbott RealTime HR HPV assay is performed with manual sample preparation method or on the *m24sp* or the *m2000sp* for sample extraction and the *m2000rt* for amplification and detection. Refer to "Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime High Risk HPV" or the *m24sp*, the *m2000sp* or the *m2000rt* Operations Manuals for detailed operating procedures.

The appropriate database containing sample preparation protocols must be installed on the *m24sp* prior to performing the assay. For detailed information on database installation, refer to the *m24sp* Operations Manual.

The Abbott RealTime HR HPV application files must be installed on the *m2000rt* and/or *m2000sp* from the Abbott RealTime High Risk HPV *m2000* System ROW Combined Application CD-ROM (List No. 4N05) prior to performing the assay. For detailed information on application file installation, refer to the *m2000sp* and the *m2000rt* Operations Manuals, Operating Instructions section.

SPECIMEN COLLECTION AND HANDLING INSTRUCTIONS

Specimen Collection

Specimens collected in PreservCyt Solution (Cytoc Corporation) or SurePath Preservative Fluid (TriPath Imaging, Inc.), or collected with Abbott Cervi-Collect Specimen Collection Kit (Abbott List No. 4N73) can be used with the Abbott RealTime HR HPV assay. For SurePath specimens, either the sample from the original SurePath collection vial or the remaining cell pellet sample obtained after cytological processing can be used. Users must follow the respective manufacturer's instructions for collecting cervical specimens in PreservCyt Solution or SurePath Preservative Fluid. Users must follow the instructions in the Abbott Cervi-Collect Specimen Collection Kit Package Insert (List No. 4N73) for collecting cervical specimens with the Abbott Cervi-Collect Specimen Collection Kit.

Specimen Transport and Storage

Cervical specimens collected in PreservCyt Solution can be transported at 15-30°C or 2-8°C and may be stored for up to 4 months at 15-30°C or up to 6 months at 2-8°C and -10°C or colder following collection.

Cervical specimens collected in SurePath Preservative Fluid (the sample from the original SurePath collection vial or the remaining cell pellet sample obtained after cytological processing) can be transported at 15-30°C or 2-8°C and may be stored for up to 2 months at 15-30°C or up to 6 months at 2-8°C and -10°C or colder following collection.

Cervical specimens collected with the Abbott Cervi-Collect Specimen Collection Kit can be transported at 2°C to 30°C and may be stored for up to 14 days at 2°C to 30°C or up to 90 days at -10°C or colder. Thaw specimens at 2°C to 30°C. Specimens should not undergo more than four freeze/thaw cycles.

Time and temperature conditions for storage must be adhered to during transport. For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical, diagnostic, or biological specimens.

ASSAY PROCEDURE

This Abbott RealTime HR HPV package insert contains three assay protocols:

- Samples prepared for amplification using the manual sample preparation method following **ASSAY PROTOCOL I**.
- Samples prepared for amplification using the *m24sp* instrument following **ASSAY PROTOCOL II**.
- Samples prepared for amplification using the *m2000sp* instrument following **ASSAY PROTOCOL III**.

Materials Provided

- Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 2N09-90)

Materials Required But Not Provided

- Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)
- Abbott RealTime High Risk HPV *m2000* System ROW Combined Application CD-ROM (List No. 4N05)
- **Materials for Manual Sample Preparation (Assay Protocol I)**

Sample Preparation Area

- Refer to the Materials and Equipment Required Section of "Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime High Risk HPV" (List No. 3N92).
- Abbott Optical Adhesive Covers (List No. 4J71-75)
- Abbott Adhesive Cover Applicator (List No. 9K32-01)

Reagent Preparation Area

- Abbott 96-Well Optical Reaction Plate (List No. 4J71-70)
- Abbott Splash Free Support Base (List No. 9K31-01)
- Calibrated Pipettes capable of delivering 20-1000 μ L
- Aerosol Barrier Pipette Tips for 20-1000 μ L Pipettes
- Single-use DNase-free tube or container
- **Materials for *m24sp* (Assay Protocol II)**

Sample Preparation Area

- Abbott *m24sp* instrument containing the scripts necessary to run the Abbott RealTime HR HPV assay (*m24sp* Database v 3.0 or higher)
- Abbott *mSample* Preparation System_{DNA} (List No. 6K12)

NOTE: One kit is sufficient to complete 192 HPV sample preparations.

- Calibrated Pipettes capable of delivering 20-1000 μ L
- Aerosol Barrier Pipette Tips for 20-1000 μ L Pipettes
- Sample input tubes (refer to “**ASSAY PROTOCOL II**” section for details)
- 1000 μ L Disposable Aerosol Barrier Pipette Tips (List No. 4J71-10)
- 200 μ L Disposable Aerosol Barrier Pipette Tips (List No. 4J71-15)
- Vortex mixer
- USP Grade 190-200 Proof Ethanol (95%-100% Ethanol). **Do not use ethanol that contains denaturants.**
- Abbott Optical Adhesive Covers (List No. 4J71-75)
- Abbott Adhesive Cover Applicator (List No. 9K32-01)
- Abbott 96 Deep-Well Plate (List No. 4J71-30)
- Abbott Splash Free Support Base (List No. 9K31-01)
- 13 mm Sample Racks
- 1.5 mL Reaction Vessels and Output Tubes (1.5 mL screw top microfuge tubes and caps, List No. 4J71-50 or equivalent)

Reagent Preparation Area

- Abbott 96-Well Optical Reaction Plate (List No. 4J71-70)
 - Abbott Splash Free Support Base (List No. 9K31-01)
 - Calibrated Pipettes capable of delivering 20-1000 μ L
 - Aerosol Barrier Pipette Tips for 20-1000 μ L Pipettes
 - Single-use DNase-free tube or container
- **Materials for *m2000sp* (Assay Protocol III)**

Sample Preparation Area

- Abbott *m2000sp* instrument with Software Version 3.0 or higher
- Abbott *mSample Preparation System*_{DNA} (List No. 6K12)

NOTE: One kit is sufficient to complete 192 HPV sample preparations.

- 5 mL Reaction Vessels (List No. 4J71-20)
- Calibrated Pipettes capable of delivering 20-1000 μ L
- Aerosol Barrier Pipette Tips for 20-1000 μ L Pipettes
- Sample input tubes (refer to “**ASSAY PROTOCOL III**” section for details)
- 1000 μ L Disposable Aerosol Barrier Pipette Tips (List No. 4J71-10)
- 200 μ L Disposable Aerosol Barrier Pipette Tips (List No. 4J71-15)
- Vortex mixer
- USP Grade 190-200 Proof Ethanol (95%-100% Ethanol). **Do not use ethanol that contains denaturants.**
- Abbott Optical Adhesive Covers (List No. 4J71-75)
- Abbott Adhesive Cover Applicator (List No. 9K32-01)
- Abbott Splash Free Support Base (List No. 9K31-01)
- Master Mix Tube (List No. 4J71-80)
- 200 mL Reagent Vessels (List No. 4J71-60)
- Abbott 96 Deep-Well Plate (List No. 4J71-30)
- Abbott 96-Well Optical Reaction Plate (List No. 4J71-70)
- 13 mm Sample Racks
- **Materials for *m2000rt***
 - Abbott *m2000rt* instrument with Software Version 3.0 or higher
 - Abbott *m2000rt* Optical Calibration Kit (List No. 4J71-93)

Other Materials

- Biological safety cabinet approved for working with infectious materials
- Sealable plastic bags
- DNase-free water[†]
- Microcentrifuge Tubes[†]
- Cotton Tip Applicators (Puritan or equivalent)[†]

[†] **NOTE:** These three items are used in the procedure for **Monitoring the Laboratory for the Presence of Contamination**. Refer to the **“QUALITY CONTROL PROCEDURES”** section of the package insert.

Procedural Precautions

- Read the instructions in the package insert carefully before processing samples.
- Do not use kits or reagents beyond expiration date.
- Control kit lots and amplification reagent kit lots can be used interchangeably. Components contained within a kit are intended to be used together. For example, do not use the negative control from control kit lot X with the positive control from control kit lot Y.
- Amplification Reagent components (enzyme, oligonucleotide reagent and activation reagent) and Controls are for single-use only and should be discarded after use. Use new reagent vessels and new reaction vessels, for every new Abbott RealTime HR HPV assay run. At the end of each run, discard all these remaining reagents as stated in the following manuals: the *m24sp* Operations Manual, the *m2000sp* Operations Manual, and “Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime High Risk HPV”.
- The Abbott RealTime HR HPV Controls must be processed with the specimens to be tested. The use of the Abbott RealTime HR HPV Controls is integral to the performance of the Abbott RealTime HR HPV assay. Refer to the **“QUALITY CONTROL PROCEDURES”** section in the package insert for details.
- Use only USP Grade 190-200 Proof Ethanol (95%-100% Ethanol) to prepare the *mWash* 2_{DNA} sample preparation reagent. **Do not use ethanol that contains denaturants.**
- Use aerosol barrier pipette tips or disposable pipettes only one time when pipetting. To prevent contamination to the pipette barrel while pipetting, care should be taken to avoid touching the pipette barrel to inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.

- Replace any empty or partially used 200 μ L and 1000 μ L disposable tips on the *m2000sp* or *m24sp* with full trays before every run. Refer to the *m2000sp* and *m24sp* Operations Manuals, Operating Instructions section.
- Monitoring procedures for the presence of amplification product can be found in the “**QUALITY CONTROL PROCEDURES**” section in the package insert.
- To reduce the risk of nucleic acid contamination, clean and disinfect spills of specimens, reagents and controls by using a detergent solution followed by a tuberculocidal disinfectant such as 1.0% (v/v) sodium hypochlorite or other suitable disinfectant.

ASSAY PROTOCOL I: MANUAL SAMPLE PREPARATION METHOD AND *m2000rt* INSTRUMENT

Refer to the “**WARNINGS AND PRECAUTIONS**” section of the package insert before preparing samples.

1. Vortex each specimen for 15-20 seconds. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial. Immediately transfer 400 μ L of each specimen to a reaction tube.

NOTE: Ensure that the volume of the cell pellet SurePath samples after cytological processing is approximately 2.8 mL. The sample volume must be adjusted to 6 mL using SurePath Preservative Fluid prior to vortexing and transferring.

2. Thaw control reagents at 15-30°C or at 2-8°C; see “**QUALITY CONTROL PROCEDURES**” section of the package insert.
 - Vortex each assay control for 15-20 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
 - Once thawed, assay controls can be stored at 2-8°C for up to 24 hours before use.
3. Thaw amplification reagents at 15-30°C or at 2-8°C and store at 2-8°C until required for the amplification master mix procedure.
 - Once thawed, the amplification reagents can be stored at 2-8°C for up to 24 hours if not used immediately.

NOTE: A maximum of 96 reactions can be performed per run.

For up to 24 reactions use: one tube of Positive Control, one tube of Negative Control, one Amplification Reagent Pack, and one set of *mSample Preparation System_{DNA}* reagents.

For 25 to 48 reactions use: one tube of Positive Control, one tube of Negative Control, two Amplification Reagent Packs, and one set of *mSample Preparation System_{DNA}* reagents.

For 49 to 72 reactions use: one tube of Positive Control, one tube of Negative Control, three Amplification Reagent Packs, one bottle of *mMicroparticle_{DNA}* and *mLysis_{DNA}* Buffer, and two bottles of *mWash1_{DNA}* Buffer, *mWash2_{DNA}* Buffer and *mElution_{DNA}* Buffer.

For 73 to 96 reactions use: one tube of Positive Control, one tube of Negative Control, four Amplification Reagent Packs, one bottle of *mMicroparticle_{DNA}* and *mLysis_{DNA}* Buffer, and two bottles of *mWash1_{DNA}* Buffer, *mWash2_{DNA}* Buffer and *mElution_{DNA}* Buffer.

Sample Preparation Area

4. Refer to the Extraction Protocol section of “Manual Sample Preparation Using the ABBOTT *mSample Preparation System_{DNA}* for RealTime High Risk HPV” for sample preparation procedures.

NOTE: *mSample Preparation System_{DNA}* reagents can be used up to three times within 14 days for a total of 48 samples when stored tightly capped at 15°C to 30°C. If reusing the *mSample Preparation System_{DNA}* reagents, mark the *mWash2_{DNA}* bottle to indicate that ethanol has already been added. Once prepared, do not add more ethanol to the *mWash2_{DNA}* bottle at any time. If reusing the *mSample Preparation System_{DNA}* reagents, after removing the caps from all the *mSample Preparation System_{DNA}* reagents, store the caps on a clean, absorbent surface for recapping after the run.

NOTE: The assembly of the amplification master mix and sample eluates into the 96-Well Optical Reaction Plate (step 12) must be initiated within one hour after completion of Sample Preparation.

Amplification Area

5. Switch on and initialize the *m2000rt*. The *m2000rt* requires a 15-minute warm-up prior to starting a run. Refer to the *m2000rt* Operations Manual, Operating Instructions section.
6. Create the *m2000rt* test order. Refer to the Operating Instructions section of the *m2000rt* Operations Manual. From the Protocol screen, select the appropriate application file.

Reagent Preparation Area

NOTE: All reagent preparation must take place in the dedicated Reagent Preparation Area. Refer to the “Contamination Precautions” section of the package insert before preparing reagents. Change gloves before handling the amplification reagents.

7. Prepare the amplification master mix.
 - Each Amplification Reagent Pack supports up to 24 reactions.
 - Prior to opening the amplification reagents, ensure that the contents of the Amplification Reagent Pack are at the bottom by tapping the Amplification Reagent Pack in an upright position on the bench to bring the liquid to the bottom of the vials.
 - Identify the amplification reagents as follows:
 - Activation Reagent (Reagent 1): clear bottle, teal cap
 - Oligonucleotide Reagent (Reagent 2): black bottle, white cap
 - AmpliTaq Gold Enzyme (Reagent 3): clear bottle, white cap
 - Remove and discard caps.
 - Prepare the master mix by using a **PIPETTE DEDICATED FOR REAGENT USE ONLY** to add 278 μL of the HPV Activation Reagent (Reagent 1) and 402 μL of the HPV Oligonucleotide Reagent (Reagent 2) together in the AmpliTaq Gold Enzyme bottle (Reagent 3). Mix the Enzyme vial containing the reaction mixture (master mix) by gently pipetting up and down six times. Avoid creating foam.
 - If performing 25 to 48 reactions, prepare the amplification master mix from two Amplification Reagent Packs. If performing 49 to 72 reactions, prepare the amplification master mix from three Amplification Reagent Packs. If performing 73 to 96 reactions, prepare the amplification master mix from four Amplification Reagent Packs.
- NOTE: The *m2000rt* protocol (step 14) must be initiated within one hour of the addition of the Activation Reagent into the AmpliTaq Gold Enzyme bottle (step 7).**
8. Pipette the contents of the master mix from the Enzyme bottle(s) into a single-use DNase-free tube. Mix by gently pipetting up and down six times. Avoid creating foam.

9. Prior to addition of master mix and sample, insert a 96-Well Optical Reaction Plate onto a Splash Free Support Base to prevent contamination.
 - Contamination of the bottom of the 96-Well Optical Reaction Plate with fluorescent materials could potentially interfere with the HPV assay. The 96-Well Optical Reaction Plate should be held and transported with the Splash Free Support Base to minimize contamination.
10. Using a **DEDICATED PIPETTE**, dispense 25 μ L aliquots of the amplification master mix into each well of the 96-Well Optical Reaction Plate that will be used depending on the number of samples to be run, including controls. A calibrated repeat pipettor may be used. Visually verify that 25 μ L has been dispensed into each well.

NOTE: Remaining activated master mix can be recapped and stored at -10°C or colder for up to 14 days and reused at a later time if the volume is sufficient. The activated master mix should not undergo more than three freeze/thaw cycles. The frozen master mix can be thawed at room temperature for up to one hour prior to the initiation of amplification and detection on the *m2000rt*.
11. Transfer the 96-Well Optical Reaction Plate on the Splash Free Support Base to the Sample Preparation Area.

Sample Preparation Area

12. In the Sample Preparation Area, transfer 25 μ L of sample eluate to the 96-Well Optical Reaction Plate on the Splash Free Support Base. **Use a separate pipette tip for each sample eluate transfer.** Visually verify that a total of 50 μ L has been dispensed into each well.
13. Seal the 96-Well Optical Reaction Plate according to the instructions in the *m2000rt* Operations Manual.

Amplification Area

14. Place the 96-Well Optical Reaction Plate in the *m2000rt* and initiate the RealTime HR HPV protocol as described in the *m2000rt* Operations Manual, Operating Instructions section. At the completion of the run, assay results are reported on the *m2000rt*. Refer to the “**RESULTS**” section of the package insert for further details.
15. After the *m2000rt* has completed the amplification and detection protocol, remove the 96-Well Optical Reaction Plate and dispose of according to the instructions in the “**Contamination Precautions**” section of the package insert. Place the 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the *m2000rt* Operations Manual along with the gloves used to handle the plate.

Post Processing Procedures

1. Refer to the Clean Up section of “Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime High Risk HPV”.
2. Clean the Splash Free Support Base before next use, according to the *m2000rt* Operations Manual.

ASSAY PROTOCOL II: *m24sp* AND *m2000rt* INSTRUMENTS

Refer to the “**WARNINGS AND PRECAUTIONS**” section of the package insert before preparing samples.

1. Vortex each specimen for 15-20 seconds. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial. Immediately transfer the specimens to the sample input tubes.

NOTE: Ensure that the volume of the cell pellet SurePath samples after cytological processing is approximately 2.8 mL. The sample volume must be adjusted to 6 mL using SurePath Preservative Fluid prior to vortexing and transferring.

- For specimens collected in PreservCyt Solution or SurePath Preservative Fluid, to ensure that 400 μ L of each specimen is transferred by the *m24sp* to the reaction vessel:
 - transfer a minimum of 500 μ L of each specimen if using Master Mix Tubes or Abbott Transport Tubes as sample input tubes.
 - transfer a minimum of 700 μ L of each specimen if using 5 mL Reaction Vessels or any other 13 mm round bottom non-skirted tubes as sample input tubes.
 - For specimens collected with the Abbott Cervi-Collect Specimen Collection Kit, load the tubes without cap directly on the *m24sp* (these specimens do not require a transfer).
2. Thaw control reagents at 15-30°C or at 2-8°C; see “**QUALITY CONTROL PROCEDURES**” section of the package insert.
 - Vortex each assay control for 15-20 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
 - Once thawed, assay controls can be stored at 2-8°C for up to 24 hours before use.

3. Thaw amplification reagents at 15-30°C or at 2-8°C and store at 2-8°C until required for the amplification master mix procedure.

- Once thawed, the amplification reagents can be stored at 2-8°C for up to 24 hours if not used immediately.

NOTE: A maximum of 24 reactions can be performed per run. For up to 24 reactions, use one tube of Positive Control, one tube of Negative Control, one Amplification Reagent Pack, and one set of *mSample Preparation System_{DNA}* reagents.

Sample Preparation Area

4. Place the controls and patient specimens into the *m24sp* sample rack, as described in the *m24sp* Operations Manual, Operating Instructions section.

CAUTION: Use only 13 mm sample racks. Do NOT skip any positions in a sample rack. Load specimens and controls into the 13 mm sample racks in consecutive positions, starting with the third position in the first sample rack. Fill all positions in each sample rack without skipping any positions before loading specimens into the next sample rack.

Insert specimen and control tubes into sample racks carefully to avoid splashing. Ensure that each tube is placed securely in the sample rack so that the bottom of the tube reaches the inside bottom of the rack.

Load filled sample racks onto the *m24sp* in consecutive sample rack positions, with the first rack farthest to the right on the worktable, and if needed, the second rack to the left of the first rack.

5. Open the *mSample Preparation System_{DNA}* reagent pack. Prepare the *mWash_{2DNA}* by adding 70 mL of USP Grade 190-200 Proof Ethanol (95%-100% Ethanol) to the *mWash_{2DNA}* bottle as described in the *mSample Preparation System_{DNA}* product information. **Do not use ethanol that contains denaturants.** Gently invert each reagent bottle to ensure a homogenous solution. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.

NOTE: *mSample Preparation System_{DNA}* reagents can be used up to three times within 14 days for a total of 48 samples when stored tightly capped at 15°C to 30°C. If reusing the *mSample Preparation System_{DNA}* reagents, mark the *mWash_{2DNA}* bottle to indicate that ethanol has already been added. Once prepared, do not add more ethanol to the *mWash_{2DNA}* bottle at any time.

6. Initiate the *m24sp* protocol as described in the *m24sp* Operations Manual, Operating Instruction section. From the Protocol screen, select the appropriate script to run the HPV assay depending on desired output vessels (*m24sp_HP_V_DNA_Tube* for 1.5 mL tubes or *m24sp_HP_V_DNA_DWP* for 96 Deep-Well Plate). When prompted by the instrument, vigorously mix or vortex the *mMicroparticle_{DNA}* bottle until the *mMicroparticles_{DNA}* are fully resuspended. Put the *mMicroparticle_{DNA}* bottle on the deck of the instrument in the designated position.

NOTE: If reusing the *mSample Preparation System_{DNA}* reagents, after removing the caps from all the *mSample Preparation System_{DNA}* reagents, store the caps on a clean, absorbent surface for recapping after the run.

NOTE: The assembly of the amplification master mix and sample eluates into the 96-Well Optical Reaction Plate (step 14) must be initiated within one hour after completion of Sample Preparation.

Amplification Area

7. Switch on and initialize the *m2000rt*. The *m2000rt* requires a 15-minute warm-up prior to starting a run. Refer to the *m2000rt* Operations Manual, Operating Instructions section.
8. Create the *m2000rt* test order. Refer to the Operating Instructions section of the *m2000rt* Operations Manual. From the Protocol screen, select the appropriate application file.

Reagent Preparation Area

NOTE: All reagent preparation must take place in the dedicated Reagent Preparation Area. Refer to the “Contamination Precautions” section of the package insert before preparing reagents. Change gloves before handling the amplification reagents.

9. Prepare the amplification master mix.
 - Each Amplification Reagent Pack supports up to 24 reactions.
 - Prior to opening the amplification reagents, ensure that the contents of the Amplification Reagent Pack are at the bottom by tapping the Amplification Reagent Pack in an upright position on the bench to bring the liquid to the bottom of the vials.

- Identify the amplification reagents as follows:
 - Activation Reagent (Reagent 1): clear bottle, teal cap
 - Oligonucleotide Reagent (Reagent 2): black bottle, white cap
 - AmpliTaq Gold Enzyme (Reagent 3): clear bottle, white cap
- Remove and discard caps.
- Prepare the master mix by using a **PIPETTE DEDICATED FOR REAGENT USE ONLY** to add 278 μL of the HPV Activation Reagent (Reagent 1) and 402 μL of the HPV Oligonucleotide Reagent (Reagent 2) together in the AmpliTaq Gold Enzyme bottle (Reagent 3). Mix the Enzyme vial containing the reaction mixture (master mix) by gently pipetting up and down six times. Avoid creating foam.

NOTE: The *m2000rt* protocol (step 16) must be initiated within one hour of the addition of the Activation Reagent into the AmpliTaq Gold Enzyme Reagent bottle (step 9).

10. Pipette the contents of the master mix from the Enzyme bottle(s) into a single-use DNase-free tube. Mix by gently pipetting up and down six times. Avoid creating foam.
11. Prior to addition of master mix and sample, insert a 96-Well Optical Reaction Plate onto a Splash Free Support Base to prevent contamination.
 - Contamination of the bottom of the 96-Well Optical Reaction Plate with fluorescent materials could potentially interfere with the HPV assay. The 96-Well Optical Reaction Plate should be held and transported with the Splash Free Support Base to minimize contamination.
12. Using a **DEDICATED PIPETTE**, dispense 25 μL aliquots of the amplification master mix into each well of the 96-Well Optical Reaction Plate that will be used depending on the number of samples to be run, including controls. A calibrated repeat pipettor may be used. Visually verify that 25 μL has been dispensed into each well.

NOTE: Remaining activated master mix can be recapped and stored at -10°C or colder for up to 14 days and reused at a later time if the volume is sufficient. The activated master mix should not undergo more than three freeze/thaw cycles. The frozen master mix can be thawed at room temperature for up to one hour prior to the initiation of amplification and detection on the *m2000rt*.
13. Transfer the 96-Well Optical Reaction Plate on the Splash Free Support Base to the Sample Preparation Area.

Sample Preparation Area

14. In the Sample Preparation Area, transfer 25 μL of sample eluate to the 96-Well Optical Reaction Plate on the Splash Free Support Base. **Use a separate pipette tip for each sample eluate transfer.** Visually verify that a total of 50 μL has been dispensed into each well.
15. Seal the 96-Well Optical Reaction Plate according to the instructions in the *m2000rt* Operations Manual.

Amplification Area

16. Place the 96-Well Optical Reaction Plate in the *m2000rt* and initiate the RealTime HR HPV protocol as described in the *m2000rt* Operations Manual, Operating Instructions section. At the completion of the run, assay results are reported on the *m2000rt*. Refer to the **“RESULTS”** section of the package insert for further details.
17. After the *m2000rt* has completed the amplification and detection protocol, remove the 96-Well Optical Reaction Plate and dispose of according to the instructions in the **“Contamination Precautions”** section of the package insert. Place the 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the *m2000rt* Operations Manual along with the gloves used to handle the plate.

Post Processing Procedures

1. At the end of each run, remove and discard all remaining reagents from the *m24sp* worktable as stated in the *m24sp* Operations Manual.
2. Decontaminate and dispose of all specimens, reagents (except for amplification master mix when applicable), and other potentially contaminated materials in accordance with local, state, and federal regulations.
3. Clean the Splash Free Support Base before next use, according to the *m2000rt* Operations Manual.

ASSAY PROTOCOL III: *m2000sp* AND *m2000rt* INSTRUMENTS

Refer to the “**WARNINGS AND PRECAUTIONS**” section of the package insert before preparing samples.

1. Vortex each specimen for 15-20 seconds. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial. Immediately transfer the specimens to the sample input tubes.

NOTE: Ensure that the volume of the cell pellet SurePath samples after cytological processing is approximately 2.8 mL. The sample volume must be adjusted to 6 mL using SurePath Preservative Fluid prior to vortexing and transferring.

- For specimens collected in PreservCyt Solution or SurePath Preservative Fluid, to ensure that 400 μ L of each specimen is transferred by the *m2000sp* to the reaction vessel:
 - transfer a minimum of 500 μ L of each specimen if using Master Mix Tubes or Abbott Transport Tubes as sample input tubes.
 - transfer a minimum of 700 μ L of each specimen if using 5 mL Reaction Vessels or any other 13 mm round bottom non-skirted tubes as sample input tubes.
 - For specimens collected with the Abbott Cervi-Collect Specimen Collection Kit, load the tubes without cap directly on the *m2000sp* (these specimens do not require a transfer).
2. Thaw control reagents at 15-30°C or at 2-8°C; see “**QUALITY CONTROL PROCEDURES**” section of the package insert.
 - Vortex each assay control for 15-20 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
 - Once thawed, assay controls can be stored at 2-8°C for up to 24 hours before use.
 3. Thaw amplification reagents at 15-30°C or at 2-8°C and store at 2-8°C until required for the amplification master mix procedure.
 - Once thawed, the amplification reagents can be stored at 2-8°C for up to 24 hours if not used immediately.

NOTE: A maximum of 96 reactions can be performed per run.

For up to 24 reactions use: one tube of Positive Control, one tube of Negative Control, one Amplification Reagent Pack, and one set of *m*Sample Preparation System_{DNA} reagents.

For 25 to 48 reactions use: one tube of Positive Control, one tube of Negative Control, two Amplification Reagent Packs, and one set of *mSample Preparation System*_{DNA} reagents.

For 49 to 72 reactions use: one tube of Positive Control, one tube of Negative Control, three Amplification Reagent Packs, one bottle of *mMicroparticle*_{DNA} and *mLysis*_{DNA} Buffer, and two bottles of *mWash1*_{DNA} Buffer, *mWash2*_{DNA} Buffer and *mElution*_{DNA} Buffer.

For 73 to 96 reactions use: one tube of Positive Control, one tube of Negative Control, four Amplification Reagent Packs, one bottle of *mMicroparticle*_{DNA} and *mLysis*_{DNA} Buffer, and two bottles of *mWash1*_{DNA} Buffer, *mWash2*_{DNA} Buffer and *mElution*_{DNA} Buffer.

NOTE: *mSample Preparation System*_{DNA} is for single-use only and should be discarded after use. Use newly opened reagents for every new Abbott RealTime HR HPV assay run.

4. Place the controls and the patient specimens into the *m2000sp* sample rack.

CAUTION: Use only 13 mm sample racks. Do NOT skip any positions in a sample rack. Load specimens and controls into the 13 mm sample racks in consecutive positions, starting with the first position in the first sample rack. Fill all positions in each sample rack without skipping any positions before loading specimens into the next sample rack.

Insert specimen and control tubes into sample racks carefully to avoid splashing. If used, bar codes on tube labels must face right for scanning. Ensure that each tube is placed securely in the sample rack so that the bottom of the tube reaches the inside bottom of the rack.

Load filled sample racks onto the *m2000sp* in consecutive sample rack positions, with the first rack farthest to the right on the worktable, and any additional rack progressively to the left of the first rack.

5. Open the *mSample Preparation System*_{DNA} reagent pack(s). Prepare the *mWash2*_{DNA} by adding 70 mL of USP Grade 190-200 Proof Ethanol (95%-100% Ethanol) to the *mWash2*_{DNA} bottle as described in the *mSample Preparation System*_{DNA} product information. **Do not use ethanol that contains denaturants.** Gently invert each reagent bottle to ensure a homogenous solution and pour the contents into the appropriate reagent vessels per the *m2000sp* Operations Manual, Operating Instructions section. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.

NOTE: Before pouring the *m*Microparticles_{DNA} into the 200 mL reagent vessels, vigorously mix or vortex until the *m*Microparticles_{DNA} are fully resuspended.

6. Initiate the sample extraction protocol as described in the *m2000sp* Operations Manual, Operating Instructions section.
7. While the *m2000sp* is performing sample preparation, switch on and initialize the *m2000rt*. The *m2000rt* requires a 15-minute warm-up prior to starting a run. Refer to the *m2000rt* Operations Manual, Operating Instructions section.

NOTE: Once sample preparation is completed, the master mix protocol should be started within one hour.

8. Load the amplification reagents and the master mix tube on the *m2000sp* worktable.
 - Prior to opening the amplification reagents, ensure that the contents of the Amplification Reagent Pack(s) are at the bottom by tapping the Amplification Reagent Pack(s) in an upright position on the bench to bring the liquid to the bottom of the vials.
 - Remove and discard vial caps.

NOTE: Change gloves before handling the amplification reagents.

9. Initiate the *m2000sp* Master Mix Addition protocol as described in the *m2000sp* Operations Manual, Operating Instructions section.
10. After the *m2000sp* has completed the addition of samples and amplification reagents, seal the 96-Well Optical Reaction Plate according to the instructions in the *m2000sp* Operations Manual.
 - Contamination of the bottom of the 96-Well Optical Reaction Plate with fluorescent materials could potentially interfere with the HPV assay. The 96-Well Optical Reaction Plate should be held and transported with the Splash Free Support Base to minimize contamination.

NOTE: Within one hour of starting the master mix protocol, the sealed 96-Well Optical Reaction Plate should be transferred to the *m2000rt* to begin amplification/detection.

11. Place the 96-Well Optical Reaction Plate in the *m2000rt* and initiate the RealTime HR HPV assay protocol as described in the *m2000rt* Operations Manual, Operating Instructions section. At the completion of the run, assay results are reported on the *m2000rt*. Refer to the “**RESULTS**” section of the package insert for further details.

12. After the *m2000rt* has completed the amplification and detection protocol, remove the 96-Well Optical Reaction Plate and dispose of according to the instructions in the “**Contamination Precautions**” section of the package insert. Place the 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the *m2000rt* Operations Manual along with the gloves used to handle the plate.

Post Processing Procedures

1. At the end of each run, remove and discard all remaining reagents from the *m2000sp* worktable as stated in the *m2000sp* Operations Manual.
2. Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local, state, and federal regulations.
3. Clean the Splash Free Support Base before next use, according to the *m2000rt* Operations Manual.

QUALITY CONTROL PROCEDURES

***m2000rt* Optical Calibration**

Optical calibration of the *m2000rt* is required for the accurate measurement and discrimination of dye fluorescence during the Abbott RealTime HR HPV assay.

The following *m2000rt* Optical Calibration Plates are used to calibrate the *m2000rt* for the Abbott RealTime HR HPV assay:

- FAM Plate (Carboxyfluorescein)
- Cy5 Plate (Cyanine)
- NED Plate (ABI proprietary dye)
- ROX Plate (Carboxy-X-rhodamine)
- VIC Plate (Proprietary dye)

Refer to the Calibration Procedures section in the *m2000rt* Operations Manual for a detailed description of how to perform an *m2000rt* Optical Calibration.

Detection of Inhibition and/or Cell Inadequacy

The Abbott RealTime HR HPV assay detects the endogenous human beta globin sequence as Internal Control (IC) signal to evaluate cell adequacy, sample extraction and amplification efficiency. A flag or an error code is displayed when IC cycle number (CN) value of a sample or control exceeds the established range.

Negative and Positive Controls

A Negative Control and a Positive Control are required for every run to verify that the sample processing, the amplification, and the detection steps are performed correctly. The Abbott RealTime HR HPV controls need to be processed together with the samples prior to running the amplification portion of the assay.

The Negative Control is formulated with DNA containing IC sequence. The only signal detected for Negative Control should be the IC signal in the Cy5 channel. The Positive Control is formulated with DNA containing HPV 16, HPV 18, HPV 58 and IC sequences. All four signals (VIC signal for HPV 16, NED signal for HPV 18, FAM signal for HPV 58, and Cy5 signal for IC) should be detected in the Positive Control. A flag is displayed when a control result is out of range. If Negative or Positive Controls are out of range, all of the samples and controls from that run must be reprocessed, beginning with sample preparation.

HR HPV must not be detected in the Negative Control. HR HPV detected in the Negative Control is indicative of contamination from other samples or amplified product introduced during sample preparation or during preparation of the 96-Well Optical Reaction Plate. To remove contamination, clean the *m24sp* or *m2000sp* and the *m2000rt* according to the *m24sp*, the *m2000sp* and the *m2000rt* Operations Manuals. For manual sample preparation, clean the equipment according to the instructions in "Manual Sample Preparation Using the ABBOTT *mSample* Preparation System_{DNA} for RealTime High Risk HPV". Following cleaning, repeat sample processing for controls and samples following the appropriate sample preparation protocol outlined in the package insert.

IC results for the Negative Control and Positive Control that are outside the validity limit indicate the occurrence of inhibition during sample preparation or during the amplification reaction steps of the assay. Repeat the processing for controls and samples following the appropriate sample preparation protocol outlined in the package insert.

Monitoring the Laboratory for the Presence of Contamination

It is recommended that the following procedure be done at least once a month to monitor laboratory surfaces and equipment for contamination. It is very important to test all areas that may have been exposed to processed samples and controls and/or amplification product. This includes routinely handled objects such as pipettes, *m24sp*, *m2000sp* and *m2000rt* function keys, bench surfaces and other equipment that may be present in the work areas.

1. Add 0.8 mL DNase-free water to a new Master Mix Tube.
2. Saturate the cotton tip of an applicator (Puritan or equivalent) in the DNase-free water from the Master Mix Tube.
3. Using the saturated cotton tip of the applicator, wipe the area to be monitored using a sweeping motion. Place the applicator into the Master Mix Tube.
4. Swirl the cotton tip in DNase-free water 10 times, then press the applicator along the inside of the tube so that the liquid drains back into the solution at the bottom of the Master Mix Tube. Discard the applicator.
5. Cap the Master Mix Tube and vortex.
6. Remove the caps from the Master Mix Tubes and test the sample according to the appropriate assay procedure section of the package insert.
7. Contamination is indicated by the detection of HR HPV in the swab samples.
 - If contamination is present, the instrument will report "HR HPV Detected" (disregard IC flag if present).
 - If there is no contamination, the instrument will report "Not Detected" or no result will be displayed (disregard error codes 4951 or 4952 if present).
8. If contamination is detected on the equipment, follow the cleaning and decontaminating guidelines given in that equipment's operations manual. If HR HPV is detected on surfaces, clean the contaminated areas with 1.0% (v/v) sodium hypochlorite solution, followed by 70% ethanol or water.

NOTE: Chlorine solutions may pit equipment and metal. Use sufficient amounts or repeated applications of 70% ethanol until chlorine residue is no longer visible.
9. Repeat testing of the contaminated area by following steps 1 through 6.
10. If the presence of contamination is detected again, repeat steps 8 and 9 until no HR HPV amplification is detected.

RESULTS

The Abbott RealTime HR HPV assay is a qualitative assay. A minimum of one Negative Control and one Positive Control are required with each run. The Negative Control serves to verify that HR HPV DNA contamination of the Negative Control did not occur during the sample preparation and set-up of the amplification reaction. If HR HPV signal is detected for the Negative Control, the -QC flag is displayed next to all sample results for the run. Samples with the -QC flag may have been similarly contaminated with analyte during processing. If the Negative Control is not processed, the -QC flag is indicated next to all sample results for that run.

The IC signal in samples serves to confirm that each sample had sufficient cell input for accurate HR HPV detection and was processed correctly and to indicate whether inhibitors of amplification are present. If the IC is out of range (i.e. IC CN not generated or greater than or equal to a fixed cutoff cycle) and HR HPV is detected, the sample will have an interpretation of "HR HPV Detected". An IC flag will be reported next to the result. If the IC is out of range and HR HPV is not detected, no result will be reported and an error code will be generated. The sample with the error code must be retested starting with sample preparation.

For more information about error codes and flags, refer to the *m2000rt* Operations Manual Version 3.0 and Operations Manual Addendum Version 3.0.

Results Reporting

Three HPV signals corresponding to HPV 16, HPV 18 and Other HR HPV are evaluated for each sample. Each signal is either determined as "Detected" if the CN is less than a fixed assay cutoff cycle or is determined as "Not Detected" if the CN is not generated or the CN is greater than or equal to the assay cutoff cycle. All the detected signals (HPV 16, HPV 18 or Other HR HPV) are reported in the sample result with the respective CN values (in parenthesis after the target result). Samples with any of the three HR HPV signals detected will have an interpretation of "HR HPV Detected". Samples with all three HR HPV signals not detected will have an interpretation of "Not Detected".

Assay results and interpretations will look similar to the following examples:

Sample ID	Results	Interpretation	Explanation
1	HPV 16 (20.76)	HR HPV Detected	HPV 16 is detected with a CN of 20.76 HPV 18 and Other HR HPV are not detected
2	HPV 18 (21.20)	HR HPV Detected	HPV 18 is detected with a CN of 21.20 HPV 16 and Other HR HPV are not detected
3	Other HR HPV (14.48)	HR HPV Detected	Other HR HPV is detected with a CN of 14.48 HPV 16 and HPV 18 are not detected
4	HPV 16 (22.20); Other HR HPV (17.21)	HR HPV Detected	HPV 16 and Other HR HPV are detected with CN of 22.20 and 17.21, respectively HPV 18 is not detected
5	HPV 18 (18.67); Other HR HPV (15.88)	HR HPV Detected	HPV 18 and Other HR HPV are detected with CN of 18.67 and 15.88, respectively HPV 16 is not detected
6	HPV 16 (24.51); HPV 18 (23.11)	HR HPV Detected	HPV 16 and HPV 18 are detected with CN of 24.51 and 23.11, respectively Other HR HPV is not detected
7	HPV 16 (21.35); HPV 18 (22.60); Other HR HPV (19.45)	HR HPV Detected	HPV 16 and HPV 18 and Other HR HPV are detected with CN of 21.35, 22.60, and 19.45, respectively
8	Not Detected	Not Detected	HR HPV is not detected

LIMITATIONS OF THE PROCEDURE

- For *In Vitro* Diagnostic Use Only.
- This method has been tested using clinically-collected PreservCyt and SurePath liquid pap and Abbott Cervi-Collect specimens. Performance with other specimen types has not been evaluated.
- Optimal performance of this test requires appropriate specimen collection, handling, and storage (refer to the “**SPECIMEN COLLECTION AND HANDLING INSTRUCTIONS**” section of the package insert).
- Use of the Abbott RealTime HR HPV assay is limited to personnel who have been trained on the use of the *m24sp* or *m2000sp* or manual sample preparation method for sample extraction and *m2000rt* for amplification and detection.
- The instruments and assay procedures reduce the risk of contamination by amplification product. However, nucleic acid contamination from the controls, specimens, and amplification product must be controlled by good laboratory practice and careful adherence to the procedures specified in the package insert.
- A negative result does not preclude the possibility of infection because results are dependent on appropriate specimen collection. Test results may be affected by improper specimen collection, technical error, or specimen mix-up.
- As with any diagnostic test, results from the Abbott RealTime HR HPV assay should be interpreted in conjunction with other clinical and laboratory findings.

SPECIFIC PERFORMANCE CHARACTERISTICS

Genotype Inclusivity and Partial Genotyping

The ability of the Abbott RealTime HR HPV assay to detect 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) and to distinguish HPV 16 and HPV 18 from the other 12 HR HPV genotypes was evaluated. Fifty-one samples containing HPV DNA targets from each of the 14 genotypes individually and in combinations were tested as listed in Table 1. Results from 51 samples that included 14 samples with single genotype, 25 samples with two genotypes and 12 samples with three genotypes were reported accurately; the presence or absence of HPV 16 and HPV 18 DNA was accurately determined in each case.

Table 1: Genotype Detection and Partial Genotyping Capability

Sample No.	HPV Genotype	Reported Result
1	HPV16	HPV 16
2	HPV18	HPV 18
3	HPV31	Other HR HPV
4	HPV33	Other HR HPV
5	HPV35	Other HR HPV
6	HPV39	Other HR HPV
7	HPV45	Other HR HPV
8	HPV51	Other HR HPV
9	HPV52	Other HR HPV
10	HPV56	Other HR HPV
11	HPV58	Other HR HPV
12	HPV59	Other HR HPV
13	HPV66	Other HR HPV
14	HPV68	Other HR HPV
15	HPV16 and HPV18	HPV 16; HPV 18
16	HPV16 and HPV31	HPV 16; Other HR HPV
17	HPV16 and HPV33	HPV 16; Other HR HPV
18	HPV16 and HPV35	HPV 16; Other HR HPV
19	HPV16 and HPV39	HPV 16; Other HR HPV
20	HPV16 and HPV45	HPV 16; Other HR HPV
21	HPV16 and HPV51	HPV 16; Other HR HPV
22	HPV16 and HPV52	HPV 16; Other HR HPV

Sample No.	HPV Genotype	Reported Result
23	HPV16 and HPV56	HPV 16; Other HR HPV
24	HPV16 and HPV58	HPV 16; Other HR HPV
25	HPV16 and HPV59	HPV 16; Other HR HPV
26	HPV16 and HPV66	HPV 16; Other HR HPV
27	HPV16 and HPV68	HPV 16; Other HR HPV
28	HPV18 and HPV31	HPV 18; Other HR HPV
29	HPV18 and HPV33	HPV 18; Other HR HPV
30	HPV18 and HPV35	HPV 18; Other HR HPV
31	HPV18 and HPV39	HPV 18; Other HR HPV
32	HPV18 and HPV45	HPV 18; Other HR HPV
33	HPV18 and HPV51	HPV 18; Other HR HPV
34	HPV18 and HPV52	HPV 18; Other HR HPV
35	HPV18 and HPV56	HPV 18; Other HR HPV
36	HPV18 and HPV58	HPV 18; Other HR HPV
37	HPV18 and HPV59	HPV 18; Other HR HPV
38	HPV18 and HPV66	HPV 18; Other HR HPV
39	HPV18 and HPV68	HPV 18; Other HR HPV
40	HPV16 and HPV18 and HPV31	HPV 16; HPV 18; Other HR HPV
41	HPV16 and HPV18 and HPV33	HPV 16; HPV 18; Other HR HPV
42	HPV16 and HPV18 and HPV35	HPV 16; HPV 18; Other HR HPV
43	HPV16 and HPV18 and HPV39	HPV 16; HPV 18; Other HR HPV
44	HPV16 and HPV18 and HPV45	HPV 16; HPV 18; Other HR HPV
45	HPV16 and HPV18 and HPV51	HPV 16; HPV 18; Other HR HPV

Sample No.	HPV Genotype	Reported Result
46	HPV16 and HPV18 and HPV52	HPV 16; HPV 18; Other HR HPV
47	HPV16 and HPV18 and HPV56	HPV 16; HPV 18; Other HR HPV
48	HPV16 and HPV18 and HPV58	HPV 16; HPV 18; Other HR HPV
49	HPV16 and HPV18 and HPV59	HPV 16; HPV 18; Other HR HPV
50	HPV16 and HPV18 and HPV66	HPV 16; HPV 18; Other HR HPV
51	HPV16 and HPV18 and HPV68	HPV 16; HPV 18; Other HR HPV

Clinical Sensitivity and Specificity in Referral Population: Disease Detection

A total of 512 PreservCyt liquid pap specimens from a referral population were tested with the Abbott RealTime HR HPV assay and the hc2 High-Risk HPV DNA Test (HC2). The clinical sensitivity and specificity for detection of disease were determined for both assays. The presence of disease was defined by a histology result of Cervical Intraepithelial Neoplasia (CIN) 2 or greater. For subjects lacking a histology evaluation, disease status was defined by a liquid based cytology (LBC) result of high grade squamous intraepithelial lesion (HSIL) or greater at enrollment in colposcopy clinics. Disease prevalence in this population was 24.6%. Of 126 disease positive subjects, 121 were detected by the Abbott RealTime HR HPV assay and 119 were detected by HC2. Of 386 disease negative specimens, 154 were not detected by the Abbott RealTime HR HPV assay and 147 were not detected by HC2. The sensitivity of the Abbott RealTime HR HPV assay for detection of disease was 96.0% and of HC2 was 94.4%. The specificity of the Abbott RealTime HR HPV assay in this referral population was 39.9% and of HC2 was 38.1% (Table 2).

Table 2: Clinical Performance for Detection of Disease in Referral Population

Test	Sensitivity (95% CI)	Specificity (95% CI)	Positive Predictive Value	Negative Predictive Value
Abbott RealTime HR HPV	96.0% (91.0-98.7%)	39.9% (35.0-45.0%)	34.3%	96.9%
HC2	94.4% (88.9-97.7%)	38.1% (33.2-43.1%)	33.2%	95.5%

A total of 128 subjects from this population had a cytology result of atypical squamous cells of undermined significance (ASCUS). In this ASCUS population, the clinical sensitivity was 100% for both the Abbott RealTime HR HPV assay and HC2. The clinical specificity of the Abbott RealTime HR HPV assay was 46.2% and of HC2 was 45.2%.

Clinical Sensitivity and Specificity in Referral Population: High Risk HPV Detection

The sensitivity and specificity of the Abbott RealTime HR HPV assay for detection of HR HPV were evaluated by testing 517 PreservCyt liquid pap specimens collected from a referral population. The high risk HPV status of cervical specimens was determined by the concordance between the Abbott RealTime HR HPV assay and HC2 tests and by further analysis of the specimens with discordant results using LINEAR ARRAY HPV Genotyping Test (Linear Array). A total of 337 specimens were detected by both assays and 136 specimens were not detected by either assay. The results of 44 discordant specimens were resolved by Linear Array. Of the 363 HR HPV positive specimens, 354 were detected by the Abbott RealTime HR HPV assay and 346 were detected by HC2. Of the 154 HR HPV negative specimens, 153 were not detected by the Abbott RealTime HR HPV assay and 137 were not detected by HC2. The sensitivity of the Abbott RealTime HR HPV assay for detection of HR HPV was 97.5% and of HC2 assay was 95.3%. The specificity of the Abbott RealTime HR HPV assay was 99.4% and of HC2 assay was 89.0% (Table 3).

Table 3: Sensitivity and Specificity for Detection of HR HPV

Test	Sensitivity (95% CI)	Specificity (95% CI)
Abbott RealTime HR HPV	97.5% (95.3-98.9%)	99.4% (96.4-100%)
HC2	95.3% (92.6-97.2%)	89.0% (82.9-93.4%)

Clinical Specificity in General Screening Population with Normal Cytology (Age ≥ 30 years)

A total of 362 cytologically normal PreservCyt liquid pap specimens collected in a general screening population from women 30 years of age or older were tested with the Abbott RealTime HR HPV assay and HC2. The clinical specificity was determined for both assays. The Abbott RealTime HR HPV assay detected 4.1% of the specimens and HC2 detected 3.0% of the specimens. In this population, the clinical specificity of the Abbott RealTime HR HPV assay and HC2 were 95.9% and 97.0% respectively.

Table 4: Clinical Specificity in General Screening Population with Normal Cytology (Age ≥ 30 years)

Test	Number detected/tested	Detection Rate	Specificity (95% CI)
Abbott RealTime HR HPV	15/362	4.1%	95.9% (93.3-97.7%)
HC2	11/362	3.0%	97.0% (94.6-98.5%)

Analytical Sensitivity for High Risk HPV Genotypes

Analytical sensitivity of the Abbott RealTime HR HPV assay was determined by testing HPV DNA from each of 14 HR HPV genotypes in the presence of human cellular DNA in PreservCyt Solution. Four hundred microliters of sample is used per assay. For each genotype, a minimum of 4 levels, with 9 replicates at each level were tested. Testing was performed with three lots of amplification reagents on three *m2000* Systems.

Probit analysis determined that with a probability of greater than 95%, HPV 16, 18, 35, 39, 45, 51, 59, 66, and 68 can be detected at 500 copies per assay, HPV 31, 33, 52, and 56 can be detected at 2,000 copies per assay and HPV 58 at 5,000 copies per assay.

Analytical Specificity (Cross-reactivity)

A panel of bacteria, viruses and fungi were evaluated for potential cross-reactivity in the Abbott RealTime HR HPV assay (Table 5). The panel included 15 Low Risk HPV genotypes and other organisms that can be found in the female anogenital tract. Human cellular DNA was also evaluated for potential cross-reactivity. Each potential cross-reactant was spiked into HPV negative samples at concentrations (per 0.4 mL sample input) shown in Table 5. Purified nucleic acids were used except where noted. Cross-reactivity was not observed with any of the organisms tested.

Table 5: Cross-reactivity Panel

Organisms	Concentration	Organisms	Concentration
<i>Bacteroides fragilis</i>	10 ⁷ genomic copies	HPV 6	10 ⁷ genomic copies
<i>Candida albicans</i> ¹	10 ⁷ CFU	HPV 11	10 ⁷ genomic copies
<i>Chlamydia trachomatis</i> ¹	10 ⁷ EBs	HPV 13	10 ⁷ genomic copies
<i>Corynebacterium genitalium</i>	10 ⁷ genomic copies	HPV 26	10 ⁷ genomic copies
<i>Enterobacter cloacae</i>	10 ⁷ genomic copies	HPV 30	10 ⁷ genomic copies
<i>Enterococcus faecalis</i>	10 ⁷ genomic copies	HPV 32	10 ⁷ genomic copies
<i>Escherichia coli</i>	10 ⁷ genomic copies	HPV 40	10 ⁷ genomic copies
<i>Gardnerella vaginalis</i>	10 ⁷ genomic copies	HPV 42	10 ⁷ genomic copies
<i>Haemophilis ducreyi</i>	10 ⁷ genomic copies	HPV 43	10 ⁷ genomic copies
<i>Lactobacillus acidophilus</i>	10 ⁷ genomic copies	HPV 44	10 ⁷ genomic copies
<i>Mycoplasma genitalium</i>	10 ⁷ genomic copies	HPV 53	10 ⁷ genomic copies
<i>Mycoplasma hominis</i>	10 ⁷ genomic copies	HPV 54	10 ⁷ genomic copies
<i>Neisseria gonorrhoeae</i>	10 ⁷ genomic copies	HPV 55	10 ⁷ genomic copies
<i>Neisseria meningitides</i>	10 ⁷ genomic copies	HPV 57	10 ⁷ genomic copies
<i>Proteus mirabilis</i>	10 ⁷ genomic copies	HPV 61	10 ⁷ genomic copies
<i>Staphylococcus aureus</i>	10 ⁷ genomic copies	HSV-I	10 ⁷ genomic copies
<i>Staphylococcus epidermidis</i>	10 ⁷ genomic copies	HSV-II	10 ⁷ genomic copies
<i>Streptococcus pneumoniae</i>	10 ⁷ genomic copies	HBV	10 ⁷ genomic copies
<i>Trichomonas vaginalis</i>	10 ⁶ genomic copies	HCV ²	10 ⁶ viral RNA copies
<i>Ureaplasma urealyticum</i>	10 ⁷ genomic copies	HIV-1	10 ⁶ viral RNA copies
Human Cellular DNA	10 ⁷ genomic copies		

¹Cultured microorganisms. ²Clinical specimen

Reproducibility

The reproducibility of the Abbott RealTime HR HPV assay was evaluated by testing a panel of 20 well-characterized clinical specimen pools (10 HR HPV positive and 10 HR HPV negative). The twenty panel members were tested by two operators. Each operator, using a unique combination of reagent lot and instrument pair, tested two replicates of each panel member per day for four days for a total of eight replicates. Percent (%) Agreement results, based on comparison of the Abbott RealTime HR HPV results to expected results, for each panel member individually and for overall negative and positive panels are shown in Table 6. For positive samples, results for each HPV signal (HPV 16, HPV 18, and Other HR HPV) were accurately reported for all replicates. The overall agreement for 319 results compared with expected results was 100%. The agreement for 159 comparisons between the two operators using two different reagent lots and two instruments was 100%.

Table 6: Reproducibility

Panel No.	Expected Result	N	% Detected	% Agreement
1	Not Detected	16	0	100
2	Not Detected	16	0	100
3	Not Detected	16	0	100
4	Not Detected	16	0	100
5	Not Detected	16	0	100
6	Not Detected	16	0	100
7	Not Detected	16	0	100
8	Not Detected	16	0	100
9	Not Detected	16	0	100
10	Not Detected	16	0	100
11	HR HPV Detected (Other HR HPV)	16	100	100
12	HR HPV Detected (HPV 16; HPV 18)	16	100	100
13	HR HPV Detected (HPV 16)	16	100	100

Panel No.	Expected Result	N	% Detected	% Agreement
14	HR HPV Detected (HPV 16; Other HR HPV)	16	100	100
15	HR HPV Detected (Other HR HPV)	16	100	100
16	HR HPV Detected (Other HR HPV)	16	100	100
17	HR HPV Detected (Other HR HPV)	15 [^]	100	100
18	HR HPV Detected (Other HR HPV)	16	100	100
19	HR HPV Detected (Other HR HPV)	16	100	100
20	HR HPV Detected (Other HR HPV)	16	100	100
Negative Samples (panels 1-10)		160	0	100
Positive Samples (panels 11-20)		159	100	100

[^] Invalid reaction was excluded from the analysis.

Reproducibility Between Manual, *m24sp* and *m2000sp* Sample Preparation Methods

Three different sample processing options are available for the Abbott RealTime HR HPV assay: manual, *m24sp*, and *m2000sp*. The reproducibility between *m2000sp* and manual sample preparation methods and between *m2000sp* and *m24sp* was determined by testing separate aliquots of the same cervical specimens using these different sample preparation methods. For each comparison 110 PreservCyt liquid pap specimens were tested. Agreements between *m2000sp* and manual sample preparation methods (Table 7) and between *m2000sp* and *m24sp* (Table 8) were both 100%.

Table 7: Agreement Between *m2000sp* and Manual Sample Preparation

	Manual Sample Preparation	
	Detected	Not Detected
<i>m2000sp</i>	55	0
	0	55

		<i>m24sp</i>	
		Detected	Not Detected
<i>m2000sp</i>	Detected	55	0
	Not Detected	0	55

Potentially Interfering Substances

The potential for interference in the Abbott RealTime HR HPV assay was assessed with substances that may be present in cervical specimens. HR HPV negative samples and HR HPV positive samples were tested in the presence or absence of each of the substances listed in Table 9. Blood and mucus were spiked into PreservCyt solution at a concentration of 5%, all other substances at a concentration of 0.5%. Interference was not observed with any of the substances tested.

Table 9: Potentially Interfering Substances Tested

Blood	Monistat-1 Day or Night Treatment
Mucus	Norforms Deodorant Suppositories
CLOTRIMAZOLE Vaginal Cream (2%)	Terazol-3 Vaginal Cream
Delfen Vaginal Contraceptive Foam	Vagi-gard Povidone Iodine Medicated Douche
Gynecort 1% Hydrocortisone Anti-itch Creme	Vagisil Anti-Itch Creme
K-Y Jelly	Vagisil Intimate Lubricant
Lubrin	Yeast Gard Homeopathic Vaginal Suppositories
MetroGel-Vaginal	Zovirax Cream (Acyclovir) 5%
Miconazole Nitrate Suppository	

Performance for Detection of High Risk HPV with Specimens Collected Using Abbott Cervi-Collect Specimen Collection Kit

Specimens collected with the Abbott Cervi-Collect Specimen Collection Kit were tested with the Abbott RealTime HR HPV assay. Specimens collected in PreservCyt Solution from the same subjects were tested with the Abbott RealTime HR HPV assay and with HC2. A total of 153 paired specimens that had sufficient volume for all three tests were included in the analysis. The high risk HPV status of cervical specimens was determined by the concordance of the Abbott RealTime HR HPV and HC2 results, and by further analysis of the specimens with discordant results using Linear Array. Among the 70 HR HPV positive specimens, the detection rates were 92.9%, 98.6% and 84.3% for Abbott RealTime HR HPV using Cervi-Collect specimens, Abbott RealTime HR HPV using PreservCyt liquid pap specimens and HC2, respectively (Table 10). Among the 83 HR HPV negative specimens, the detection rates were 3.6%, 2.4% and 3.6% for Abbott RealTime HR HPV using Cervi-Collect specimens, Abbott RealTime HR HPV using PreservCyt liquid pap specimens and HC2, respectively (Table 10).

Table 10: HR HPV Detection

Test	HR HPV Positive (N=70)		HR HPV Negative (N=83)	
	Number detected	% Detected (95% CI)	Number detected	% Detected (95% CI)
Abbott RealTime HR HPV with Cervi-Collect	65	92.9 (84.1-97.6)	3	3.6 (0.8-10.2)
Abbott RealTime HR HPV with PreservCyt Liquid Pap	69	98.6 (92.3-100)	2	2.4 (0.3-8.4)
HC2 with PreservCyt Liquid Pap	59	84.3 (73.6-91.9)	3	3.6 (0.8-10.2)

The agreement in Abbott RealTime HR HPV results with specimens collected in Cervi-Collect versus specimens collected in PreservCyt Solution from the same patients was 94.4% (Table 11).

Table 11: Agreement between Cervi-Collect and PreservCyt Liquid Pap Specimens

		Abbott RealTime HR HPV Cervi-Collect	
		Detected	Not Detected
Abbott RealTime HR HPV PreservCyt Liquid Pap	Detected	69	6
	Not Detected	3	83

Agreement = 94.4% (152/161)

Performance for Detection of High Risk HPV with Specimens Collected in SurePath Preservative Fluid

To assess the performance of the Abbott RealTime HR HPV assay with specimens collected in SurePath Preservative Fluid, a total of 263 specimens collected in SurePath Preservative Fluid were tested with the Abbott RealTime HR HPV assay and the HC2 test. Both the sample from the original SurePath collection vial and the remaining cell pellet sample obtained after cytological processing were tested with the Abbott RealTime HR HPV assay. The cell pellet sample was tested with HC2 per manufacturer's instructions. The high risk HPV status of cervical specimens was determined by the concordance of the Abbott RealTime HR HPV and HC2 results, and by further analysis of the specimens with discordant results using Linear Array. Among the 138 HR HPV positive specimens, the detection rates were 98.6%, 97.1% and 99.3% for Abbott RealTime HR HPV using the sample from the original collection vial, Abbott RealTime HR HPV using the cell pellet sample and HC2, respectively (Table 12). Among the 125 HR HPV negative specimens, the detection rates were 0.0%, 0.0% and 13.6% for Abbott RealTime HR HPV using the sample from the original collection vial, Abbott RealTime HR HPV using the cell pellet sample and HC2, respectively (Table 12).

Table 12: HR HPV Detection

Test	HR HPV Positive (N=138)		HR HPV Negative (N=125)	
	Number detected	% Detected (95% CI)	Number detected	% Detected (95% CI)
Abbott RealTime HR HPV with SurePath samples from Original Collection Vial	136	98.6 (94.9-99.8)	0	0 (0.0-2.9)
Abbott RealTime HR HPV with SurePath samples from Cell Pellet	134	97.1 (92.7-99.2)	0	0 (0.0-2.9)
HC2 with SurePath samples from Cell Pellet	137	99.3 (96.0-100)	17	13.6 (8.1-20.9)

The agreement in Abbott RealTime HR HPV results with the SurePath sample from the original collection vial versus the cell pellet sample was 99.2% (Table 13).

Table 13: Agreement between SurePath Samples from Original Collection Vial and Cell Pellet

Abbott RealTime HR HPV Original Collection Vial		Abbott RealTime HR HPV Cell Pellet	
		Detected	Not Detected
Detected	Detected	134	2
	Not Detected	0	129

Agreement = 99.2% (263/265)

BIBLIOGRAPHY

1. Howley PM. Papillomaviridae: the viruses and their replication. In: Fields BN, Knipe DM, Howley PM, eds. *Virology*, 3rd ed. Philadelphia, Lippincott-Raven Publishers 1996:947-78.
2. CDC. Genital HPV Infection - CDC Fact Sheet. 2008; <http://www.cdc.gov/std/HPV/STDFact-HPV.htm>.
3. zur Hausen H. Papillomaviruses and cancer: from basic studies to clinical application. *Nat Rev Cancer*. 2002;2:342-50.
4. Walboomers JMM, Jacobs MV, Manos MM, *et al*. Human papillomavirus is a necessary cause of invasive cervical cancer worldwide. *J Pathol* 1999;189:12–19.
5. Snijders PJ, Steenbergen RD, Heideman DA, *et al*. HPV-mediated cervical carcinogenesis: concepts and clinical implications. *J Pathol*. 2006;208:152-64.
6. Kjaer SK, van den Brule AJC, Paull G, *et al*. Type specific persistence of high risk human papillomavirus (HPV) as indicator of high grade cervical squamous intraepithelial lesions in young women: population based prospective follow up study. *BMJ* 2002;325:572–578.
7. Cuschieri KS, Cubie HA, Whitley MW, *et al*. Persistent high risk HPV infection associated with development of cervical neoplasia in a prospective population study. *J Clin Pathol*. 2005;58:946-50.
8. de Villiers EM, Fauquet C, Broker TR, *et al*. Classification of papillomaviruses, *Virology* 2004;324:17-27.
9. IARC Monographs on the evaluation of carcinogenic risks to humans. Human Papillomaviruses. Lyon: *International Agency for Research on Cancer* 2007; Volume 90.
10. Muñoz N, Bosch FX, de Sanjosé S, *et al*. Epidemiologic classification of human papillomavirus types associated with cervical cancer. *N Engl J Med*. 2003;348:518-27.
11. Clifford GM, Smith JS, Plummer M, *et al*. Human papillomavirus types in invasive cervical cancer worldwide: a meta-analysis. *Br J Cancer*. 2003;88:63-73.
12. Muñoz N, Castellsagué X, de González AB, *et al*. Chapter 1: HPV in the etiology of human cancer. *Vaccine*. 2006;24 Suppl 3:S3/1-10.

13. Smith JS, Lindsay L, Hoots B, *et al.* Human papillomavirus type distribution in invasive cervical cancer and high-grade cervical lesions: a meta-analysis update. *Int J Cancer.* 2007;121:621-32.
14. Khan MJ, Castle PE, Lorincz AT, *et al.* The elevated 10-year risk of cervical precancer and cancer in women with human papillomavirus (HPV) type 16 or 18 and the possible utility of type-specific HPV testing in clinical practice. *J Natl Cancer Inst.* 2005;97:1072-9.
15. Davies P, Arbyn M, Dillner J, *et al.* A report on the current status of European research on the use of human papillomavirus testing for primary cervical cancer screening. *Int J Cancer.* 2006;118:791-6.
16. Cuzick J, Clavel C, Petry KU, *et al.* Overview of the European and North American studies on HPV testing in primary cervical cancer screening. *Int J Cancer.* 2006;119:1095-101.
17. Mayrand MH, Duarte-Franco E, Rodrigues I, *et al.* Human papillomavirus DNA versus Papanicolaou screening tests for cervical cancer. *N Engl J Med.* 2007;357:1579-88.
18. Goldie SJ, Gaffikin L, Goldhaber-Fiebert JD, *et al.* Cost-effectiveness of cervical-cancer screening in five developing countries. *N Engl J Med.* 2005;353:2158-68.
19. Kim JJ, Wright TC, Goldie SJ. Cost-effectiveness of human papillomavirus DNA testing in the United Kingdom, The Netherlands, France, and Italy. *J Natl Cancer Inst.* 2005;97:888-95.
20. Goldie SJ, Kim JJ, Wright TC. Cost-effectiveness of human papillomavirus DNA testing for cervical cancer screening in women aged 30 years or more. *Obstet Gynecol.* 2004;103:619-31.
21. Cuschieri KS, Cubie HA. The role of human papillomavirus testing in cervical screening. *J Clin Virol.* 2005;32 Suppl 1:S34-42.
22. Franco EL, Cuzick J. Cervical cancer screening following prophylactic human papillomavirus vaccination. *Vaccine.* 2008;26 Suppl 1:A16-23.
23. Stanley M, Villa LL. Monitoring HPV vaccination. *Vaccine.* 2008;26 Suppl 1:A24-7.

24. National Committee for Clinical Laboratory Standards. Clinical Laboratory Waste Management: Approved Guideline - Second Edition. NCCLS Document GP5-A2. Wayne, PA: NCCLS, 2002;22(3):1-23, 32-44.
25. US Environmental Protection Agency. EPA Guide for Infectious Waste Management Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1-5-5, R1-R3, A1-A24.

Abbott *m*, *m2000*, *m2000rt*, *m2000sp*, and Cervi-Collect are trademarks of Abbott Laboratories in various jurisdictions. AmpliTaq Gold, ProClin, FAM, ROX, NED, VIC, Cy5, PreservCyt, the Spirit design, Celera, PrepStain, SurePath, TriPath Imaging, hc2 High-Risk HPV DNA Test, Linear Array, Delfen, Gynecort, K-Y Jelly, Lubrin, MetroGel-Vaginal, Monistat, Norforms, Terazol, Vagi-gard, Vagisil, Yeast Gard and Zovirax are property of their respective owners.

www.abbottmolecular.com



May 2020

© 2008, 2020 Abbott Laboratories

2.2 2.2 Abbott RealTime High Risk HPV Amplification Reagent Kit IFU

List Number: 02N09-092



Read Highlighted Changes:
Revised May 2020

Key to Symbols Used	
	Reference Number
	Lot Number
	In Vitro Diagnostic Medical Device
	Use By
	Negative Control
	Positive Control
AMPLIFICATION REAGENT PACK	
	Amplification Reagent Pack
	Upper Limit of Temperature
	Temperature Limit
	WARNING
	Contains sufficient for <n> tests
	Consult instructions for use
	Manufacturer

See REAGENTS section for a full explanation of symbols used in reagent component naming.

Customer Service

International: Call your Abbott Representative

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

Abbott RealTime High Risk HPV

INTENDED USE

The Abbott RealTime High Risk HPV is a qualitative in vitro test for the detection of DNA from 14 high risk human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 in clinical specimens. The assay specifically identifies HPV genotypes 16 and 18 while concurrently detecting the other high risk genotypes at clinically relevant infection levels.

The Abbott RealTime High Risk HPV is indicated:

- To screen patients with ASC-US (atypical squamous cells of undetermined significance) cervical cytology results to determine the need for referral to colposcopy. The results of this test are not intended to prevent women from proceeding to colposcopy.

- To be used with cervical cytology to adjunctively screen to assess the presence or absence of high risk HPV genotypes.
- To be used as a first-line primary screening test to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease.
- To assess the presence or absence of HPV genotypes 16 and 18 to identify women at increased risk for the development of cervical cancer or the presence of high-grade disease with or without cervical cytology.

The results from the Abbott RealTime High Risk HPV, together with the physician's assessment of cytology history, other risk factors, and professional guidelines, may be used to guide patient management.

SUMMARY AND EXPLANATION OF THE TEST

HPV is a small, non-enveloped, double-stranded DNA virus (approximately 8,000 base pairs) that replicates in the nucleus of squamous epithelial cells and induces hyperproliferative lesions.¹ HPV infections are among the most common sexually transmitted infections.² Most HPV infections have a benign clinical consequence and are cleared spontaneously.³ However, persistent HPV infection may result in progression to cervical cancer.⁴⁻⁷ More than one hundred different HPV genotypes have been identified, among which over forty infect mucosal and genital epithelia.⁸ Genital HPV genotypes are generally classified into high risk (HR) and low risk (LR) groups based on their carcinogenic potential. HR HPV genotypes are associated with invasive cervical cancer or its immediate precursor (high-grade squamous intraepithelial lesion, cervical intraepithelial neoplasia or carcinoma *in situ*), whereas LR HPV genotypes induce benign lesion and are not associated with cervical cancer.⁹⁻¹² Approximately 70% of invasive cervical cancer cases worldwide are caused by HPV 16 and HPV 18.¹³ Infection by HPV 16 or HPV 18 is associated with higher risk of disease progression compared to other HR HPV genotypes.¹⁴ Compared with cervical screening methods identifying cytological abnormalities, molecular tests that specifically detect the presence of HR HPV DNA in cervical cells can potentially increase sensitivity and cost-effectiveness of cervical cancer screening programs.¹⁵⁻²⁰ Furthermore, HPV DNA tests can be effectively used in triaging patients with equivocal cytology, in post-therapeutic follow-up and in monitoring vaccine efficacy.²¹⁻²³

The Abbott RealTime HR HPV assay is a qualitative in vitro test that amplifies and detects HR HPV DNA in cervical cells collected in liquid media. The detection of 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) is achieved through a primer mix targeting a conserved region of HPV genomes and single-stranded DNA probes. The assay can differentiate between HPV 16, HPV 18 and non-HPV 16/18 genotypes (Other HR HPV).

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The Abbott RealTime HR HPV assay uses the Abbott *m2000sp* instrument, the Abbott *m24sp* instrument, or the manual sample preparation method for processing samples and the Abbott *m2000rt* instrument for amplification and detection. A primer mix consisting of 3 forward primers and 2 reverse primers targeting a conserved L1 region is used to amplify HPV targets. Signal for 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) is generated with the use of fluorescent labeled probes. Internal Control (IC) amplicons are generated with a primer set targeting an endogenous human beta globin sequence and are detected with the IC specific probe. The Abbott RealTime HR HPV assay detects the endogenous human beta globin sequence as sample validity control for cell adequacy, sample extraction and amplification efficiency. Probes for HPV 16, HPV 18, non-HPV 16/18 genotypes (Other HR HPV) and IC are labeled with different fluorophores allowing their signals to be distinguishable in a single reaction.

Sample Preparation

The purpose of sample preparation is to extract, concentrate, and purify the target DNA molecules for amplification.

The Abbott *mSample Preparation System_{DNA}* uses magnetic particle technology to capture nucleic acids and washes the particles to remove unbound sample components. The bound nucleic acids are eluted and are then ready for amplification.

NOTE: One Abbott *mSample Preparation System_{DNA}* kit is sufficient to complete 4 x 48 (192) HPV sample preparations.

Two automated instrument systems, the Abbott *m2000sp* or the Abbott *m24sp*, can be used to prepare samples for the Abbott RealTime HR HPV assay. The Abbott *m2000sp* provides automated sample eluate transfer and reaction assembly in the Abbott 96-Well Optical Reaction Plate, while the Abbott *m24sp* requires manual sample eluate transfer and reaction assembly.

Alternatively, samples can be prepared manually following the instructions in **Manual Sample Preparation Using the Abbott *mSample Preparation System_{DNA}* for RealTime High Risk HPV** (List No. 3N92). The manual sample preparation method requires manual transfer of the eluted samples to a Abbott 96-Well Optical Reaction Plate and manual reaction assembly before amplification.

Reagent Preparation and Reaction Plate Assembly

The Abbott *m2000sp* combines the Abbott RealTime HR HPV Amplification Reagent components (HPV Oligonucleotide Reagent, DNA Polymerase, and Activation Reagent). The Abbott *m2000sp* dispenses the resulting master mix to the Abbott 96-Well Optical Reaction Plate along with aliquots of the nucleic acid samples prepared by the Abbott *m2000sp*. The plate is ready, after manual application of the optical seal, for transfer to the Abbott *m2000rt*.

The Abbott *m24sp* users and manual sample preparation method users manually combine the Abbott RealTime HR HPV Amplification Reagent components to create the amplification master mix and transfer aliquots of the master mix and sample eluates to the Abbott 96-Well Optical Reaction Plate. The plate is ready, after manual application of the optical seal, for transfer to the Abbott *m2000rt*.

Amplification

During the amplification reaction on the Abbott *m2000rt*, the target DNA is amplified by DNA Polymerase in the presence of dNTPs and magnesium. The DNA Polymerase is a thermophilic enzyme that has been modified in its active site by a molecule that renders it inactive. When the enzyme is heated prior to the initiation of PCR, the inhibitory molecule is cleaved from the enzyme allowing it to regain its activity. In this way, the enzyme is only active at temperatures where specific DNA-DNA interactions occur. This greatly reduces non-specific PCR artifacts such as primer dimers. In the Abbott RealTime HR HPV assay, the DNA Polymerase is first activated at 92°C for 10 minutes. During each subsequent round of thermal cycling, a high temperature is used to melt double-stranded DNA strands apart, followed by a low temperature where primers anneal to their respective targets and are extended to generate double-stranded DNA products. Exponential amplification of the products is achieved through repeated cycling between high and low temperatures, resulting in a billion-fold or greater amplification of target sequences. Amplification of both targets (HPV and IC) takes place simultaneously in the same reaction.

The target sequence for the Abbott RealTime HR HPV assay is in the conserved L1 region of the HPV genomes. A primer mix consisting of three forward primers and two reverse primers is designed to hybridize to the consensus regions among HPV genotypes of approximately 150 bases. The IC target sequence is a region of 136 bases in the endogenous human beta globin gene.

Detection

During the last 38 cycles of amplification, in an additional reading step, the temperature is lowered further to allow fluorescence detection of amplification products as the HPV and IC probes anneal to their targets (referred to as real-time fluorescence detection). The HPV and IC probes are single-stranded DNA oligonucleotides modified with a fluorescent moiety covalently linked to one end of the probe and a quenching moiety to the other end. In the absence of HPV or IC target sequences, the probes adopt a series of random conformations, some of which bring the quencher close enough to the excited fluorophore to absorb its energy before it can be fluorescently emitted. When a probe binds to its complementary sequence in the target, the fluorophore and the quencher are held apart, allowing fluorescent emission and detection by the Abbott *m2000rt*.

Signal for 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) is generated with the use of fluorescent labeled probes. IC signal is generated with an IC specific probe. Probes for HPV 16, HPV 18, Other HR HPV and IC are labeled with different fluorophores allowing their distinct signals to be simultaneously detected and distinguishable in a single reaction. Signals for HPV 16, HPV 18, Other HR HPV, and IC are detected in VIC, NED, FAM, and Cy5 channels, respectively.

Assay Results

The Abbott RealTime HR HPV assay is a qualitative assay. Results are reported as detected or not detected. In addition, each detected signal (HPV 16, HPV 18, or Other HR HPV) is also listed in the reported result. Refer to the **RESULTS** section of the package insert for further details.

REAGENTS

The Abbott RealTime HR HPV assay consists of 2 kits:

- Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 02N09-092)
- Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)

Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 02N09-092)

AMPLIFICATION REAGENT PACK (4 packs, 24 tests/pack)

Each Reagent Pack contains:

- 1 Bottle (0.070 mL) DNA Polymerase (5.4 to 5.9 Units/ μ L) in a buffered solution with stabilizers.
- 1 Bottle (0.502 mL) HPV Oligonucleotide Reagent. < 0.1% synthetic oligonucleotides and < 1% dNTPs, in a buffered solution with a reference dye. Preservatives: sodium azide and 0.16% ProClin 950.
- 1 Bottle (0.778 mL) Activation Reagent. 38 mM magnesium chloride in a buffered solution. Preservatives: sodium azide and 0.15% ProClin 950.

NOTE: The Abbott RealTime Reagent components (enzyme, oligonucleotide reagent, activation reagent) are intended for single-use only and unused reagents should be discarded.

Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)

Abbott RealTime High Risk HPV Negative Control

- **CONTROL₋** (12 vials, 0.5 mL per vial)
< 0.01% noninfectious DNA with Beta Globin sequence in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.

Abbott RealTime High Risk HPV Positive Control

- **CONTROL₊** (12 vials, 0.5 mL per vial)
< 0.01% noninfectious DNA with HPV and Beta Globin sequences in a buffered solution with carrier DNA. Preservatives: sodium azide and 0.15% ProClin 950.

NOTE: The Negative and Positive Controls are intended for single-use only and unused reagents should be discarded.

WARNINGS AND PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use
- In Vitro Test

Safety Precautions

Refer to the Abbott *m2000sp* (List No. 9K20), Abbott *m24sp* (List No. 3N09), and Abbott *m2000rt* (List No. 9K25) Operations Manuals, Hazards Section, and **Manual Sample Preparation Using the Abbott *mSample Preparation System_{DNA}* for RealTime HR HPV** (List No. 3N92) for instructions on safety precautions.

- There are no human sourced materials in any of the Abbott RealTime HR HPV Amplification Reagents or Abbott RealTime HR HPV Controls.
- This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled with appropriate biosafety practices. Wear disposable gloves while handling specimens and wash hands thoroughly afterwards. Use of protective eyewear is recommended.

The controls, HPV Oligonucleotide Reagent and Activation Reagent contain a methylisothiazolone (which is a component of ProClin). The following warnings apply to the HPV Oligonucleotide Reagent, Activation Reagent and the controls.



Warning

Hazard-determining components of labeling:

2-Methyl-2H-isothiazol-3-one	
Sodium Azide	
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
P261	Avoid breathing mist / vapours / spray.
P280	Wear protective gloves / protective clothing / eye protection / face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P501	Dispose of contents / container in accordance with local regulations.

Specimen Collection and Handling Precautions

- Specimens collected in PreservCyt Solution (Hologic, Inc.) can be used with the Abbott RealTime HR HPV assay. Users must follow the manufacturer's instructions for collecting and handling cervical specimens in PreservCyt Solution.
- Specimens collected in SurePath Preservative Fluid (TriPath Imaging, Inc.) can be used with the Abbott RealTime HR HPV assay. Either the sample from the original SurePath collection vial or the remaining cell pellet sample obtained after slide preparation with TriPath Imaging PrepStain Slide Processor can be used for testing. Users must follow the manufacturer's instructions for collecting, handling and processing cervical specimens in SurePath Preservative Fluid.
- Specimens collected with the Abbott Cervi-Collect Specimen Collection Kit can be used with the Abbott RealTime HR HPV assay. Users must follow the instructions in the Abbott Cervi-Collect Specimen Collection Kit Package Insert (List No. 4N73) for collecting and handling cervical specimens.

Laboratory Precautions

- During preparation of samples, compliance with good laboratory practices is essential to minimize the risk of cross-contamination between samples as well as the inadvertent introduction of nucleases into samples during and after the extraction procedure. Proper aseptic technique should always be used when working with DNA.
- Work area and instrument platforms must be considered potential sources of contamination.
Change gloves after having contact with potential contaminants (such as DNases, specimens, eluates, and/or amplified product) before handling unopened reagents, negative control, positive control, or specimens. Refer to the Abbott *m24sp*, Abbott *m2000sp*, and Abbott *m2000rt* Operations Manuals for instrument cleaning procedures.
- Wear appropriate personal protective equipment at all times.
- Use powder-free gloves.

- To reduce the risk of nucleic acid contamination due to aerosols formed during pipetting, pipettes with aerosol barrier tips must be used for all pipetting. The length of the tip should be sufficient to prevent contamination of the pipette barrel. While pipetting, care should be taken to avoid touching the pipette barrel to the inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.
- Change aerosol barrier pipette tips between ALL manual liquid transfers.
- Clean and disinfect spills of specimens and reagents as stated in the following manuals: the Abbott *m24sp* Operations Manual, the Abbott *m2000sp* Operations Manual, the Abbott *m2000rt* Operations Manual, and **Manual Sample Preparation Using the Abbott *mSample* Preparation System_{DNA} for RealTime High Risk HPV.**

Contamination Precautions

- Amplification reactions such as PCR are sensitive to accidental introduction of product from previous amplification reactions. Incorrect results could occur if either the clinical specimen or the reagents used become contaminated by accidental introduction of even a few molecules of amplification product. Measures to reduce the risk of contamination in the laboratory include physically separating the activities involved in performing PCR in compliance with good laboratory practices.
- The use of 3 dedicated areas within the laboratory is recommended for performing the Abbott RealTime HR HPV assay with the Abbott *m24sp* or manual sample preparation using the Abbott *mSample* Preparation System_{DNA} and the Abbott *m2000rt*:
 - The Reagent Preparation Area is dedicated to combining the Abbott RealTime HR HPV Amplification Reagent components to create the amplification master mix and transferring aliquots of the master mix to the Abbott 96-Well Optical Reaction Plate. Laboratory coats, pipettes, and pipette tips used in the Reagent Preparation Area must remain in this area and not be moved to either the Sample Preparation Area or the Amplification Area. Do not bring target or amplification product into the Reagent Preparation Area.
 - The Sample Preparation Area is dedicated to processing samples (specimens, Abbott RealTime HR HPV Controls) and to adding processed samples and controls to the Abbott 96-Well Optical Reaction Plate. All reagents used in the Sample Preparation Area should remain in this dedicated area at all times. Laboratory coats, pipettes, pipette tips, and vortex mixers used in the Sample Preparation Area must remain in this area and not be moved to either the Reagent Preparation Area or the Amplification Area. Do not bring amplification product into the Sample Preparation Area.
 - The Amplification Area is dedicated to the amplification and detection of amplified product. Laboratory coats and equipment used in the Amplification Area must remain in this area and not be moved to either the Reagent Preparation Area or the Sample Preparation Area.
- Only 2 dedicated areas, Sample Preparation Area and Amplification Area, are recommended when the Abbott *m2000sp* and Abbott *m2000rt* are used.
- If the Abbott *m2000sp* run is aborted, dispose of all commodities and reagents according to the Abbott *m2000sp* Operations Manual. If the Abbott *m24sp* run is aborted, dispose of all commodities and reagents (if not being reused) according to the Abbott *m24sp* Operations Manual. If the manual sample preparation procedure is incorrectly performed or is interrupted at any point so that the timing of the steps exceeds the recommended timing per the manual instructions, dispose of all commodities and reagents (if not being reused) according to the instructions in **Manual Sample Preparation Using the Abbott *mSample* Preparation System_{DNA} for RealTime High Risk HPV.**

- If the Abbott *m2000sp* master mix addition protocol is aborted after amplification reagents are added to the Abbott 96-Well Optical Reaction Plate, seal the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the Abbott *m2000sp* Operations Manual, Hazards section, along with the gloves used to handle the plate. Do not import the test order onto the Abbott *m2000rt*. If manual preparation of the PCR reaction mix is aborted after amplification reagents are added to the Abbott 96-Well Optical Reaction Plate, seal the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to laboratory guidelines, along with the gloves used to handle the plate.
- For all completed, interrupted or aborted Abbott *m2000rt* runs, dispose of the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.
- **Autoclaving the sealed reaction plate will not degrade the amplified product and may contribute to the release of the amplified product by opening the sealed plate. The laboratory area can become contaminated with amplified product if the waste materials are not carefully handled and contained.**
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local, state, and federal regulations.^{24,25} All materials should be handled in a manner that minimizes the chance of potential contamination of the work area.


Contamination from External Deoxy-Uracil (dU)-Containing Amplified Product

- HPV amplification assays containing dU may cause contamination and inaccurate results in the Abbott RealTime HR HPV. When negative controls are persistently reactive or where contamination with dU-containing HPV amplified product is likely to have occurred, it is recommended that the laboratory uses a contamination control procedure. This procedure (List No. 2N09-66) is available through your Abbott representative.


REAGENT STORAGE AND HANDLING INSTRUCTIONS

NOTE: Care must be taken to separate the Abbott RealTime High Risk HPV Amplification Reagent Kit that is in use from direct contact with specimens and Abbott RealTime High Risk HPV Control Kit reagents.

Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 02N09-092)

- 
- The Abbott RealTime High Risk HPV Amplification Reagent Pack must be stored at -25 to -15°C when not in use.
 - Reagents are shipped on dry ice.

Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)

- 
- The Abbott RealTime High Risk HPV Negative and Positive Controls must be stored at -10°C or colder.
 - Reagents are shipped on dry ice.

INSTRUMENTS/METHODS

The Abbott RealTime HR HPV assay is performed with manual sample preparation method or on the Abbott *m24sp* or the Abbott *m2000sp* for sample extraction and the Abbott *m2000rt* for amplification and detection. Refer to **Manual Sample Preparation Using the Abbott *mSample Preparation System_{DNA}* for RealTime High Risk HPV** or the Abbott *m24sp*, the Abbott *m2000sp* or the Abbott *m2000rt* Operations Manuals for detailed operating procedures.

The appropriate database containing sample preparation protocols must be installed on the Abbott *m24sp* prior to performing the assay. For detailed information on database installation, refer to the Abbott *m24sp* Operations Manual.

The Abbott RealTime HR HPV application files must be installed on the Abbott *m2000rt* and/or Abbott *m2000sp* from the Abbott RealTime High Risk HPV Abbott *m2000* System ROW Combined Application CD-ROM (List No. 4N05) prior to performing the assay. For detailed information on application file installation, refer to the Abbott *m2000sp* and the Abbott *m2000rt* Operations Manuals, Operating Instructions section.

SPECIMEN COLLECTION AND HANDLING INSTRUCTIONS

Specimen Collection

Specimens collected in PreservCyt Solution (Hologic, Inc.) or SurePath Preservative Fluid (TriPath Imaging, Inc.), or collected with Abbott Cervi-Collect Specimen Collection Kit (Abbott List No. 4N73) can be used with the Abbott RealTime HR HPV assay. For SurePath specimens, either the sample from the original SurePath collection vial or the remaining cell pellet sample obtained after cytological processing can be used. Users must follow the respective manufacturer's instructions for collecting cervical specimens in PreservCyt Solution or SurePath Preservative Fluid. Users must follow the instructions in the Abbott Cervi-Collect Specimen Collection Kit Package Insert (List No. 4N73) for collecting cervical specimens with the Abbott Cervi-Collect Specimen Collection Kit.

Specimen Transport and Storage

Cervical specimens collected in PreservCyt Solution can be transported at 15 to 30°C or 2 to 8°C and may be stored for up to 4 months at 15 to 30°C or up to 6 months at 2 to 8°C and -10°C or colder following collection.

Cervical specimens collected in SurePath Preservative Fluid (the sample from the original SurePath collection vial or the remaining cell pellet sample obtained after cytological processing) can be transported at 15 to 30°C or 2 to 8°C and may be stored for up to 2 months at 15 to 30°C or up to 6 months at 2 to 8°C and -10°C or colder following collection.

Cervical specimens collected with the Abbott Cervi-Collect Specimen Collection Kit can be transported at 2 to 30°C and may be stored for up to 14 days at 2 to 30°C or up to 90 days at -10°C or colder. Thaw specimens at 2 to 30°C. Specimens should not undergo more than four freeze/thaw cycles.

Time and temperature conditions for storage must be adhered to during transport. For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical, diagnostic, or biological specimens.

ASSAY PROCEDURE

This Abbott RealTime HR HPV package insert contains 3 assay protocols:

- Samples prepared for amplification using the manual sample preparation method following **ASSAY PROTOCOL I**.
- Samples prepared for amplification using the Abbott *m24sp* instrument following **ASSAY PROTOCOL II**.
- Samples prepared for amplification using the Abbott *m2000sp* instrument following **ASSAY PROTOCOL III**.

Materials Provided

- Abbott RealTime High Risk HPV Amplification Reagent Kit (List No. 02N09-092)

Materials Required But Not Provided

- Abbott RealTime High Risk HPV Control Kit (List No. 2N09-80)
- Abbott RealTime High Risk HPV *m2000* System ROW Combined Application CD-ROM (List No. 4N05)
- **Materials for Manual Sample Preparation (Assay Protocol I)**

Sample Preparation Area

- Refer to the Materials and Equipment Required Section of **Manual Sample Preparation Using the Abbott *mSample Preparation System_{DNA}* for RealTime High Risk HPV** (List No. 3N92).
- Abbott Optical Adhesive Covers (List No. 04J71-75)
- Abbott Adhesive Cover Applicator (List No. 9K32-01)

Reagent Preparation Area

- Abbott 96-Well Optical Reaction Plate (List No. 04J71-70)
- Abbott Splash-Free Support Base (List No. 09K31-01)
- Calibrated precision pipettes capable of delivering 10 µL to 1000 µL
- 20 µL to 1000 µL aerosol barrier pipette tips for precision pipettes
- Single-use DNase-free tube or container

Materials for Abbott *m24sp* (Assay Protocol II)

Sample Preparation Area

- Abbott *m24sp* instrument containing the scripts necessary to run the Abbott RealTime HR HPV assay (*m24sp* Database v 3.0 or higher)
- Abbott *mSample Preparation System_{DNA}* (List No. 06K12-24)

NOTE: One kit is sufficient to complete 192 HPV sample preparations.

- Calibrated precision pipettes capable of delivering 10 µL to 1000 µL
- 20 µL to 1000 µL aerosol barrier pipette tips for precision pipettes
- Sample input tubes (refer to **ASSAY PROTOCOL II** section for details)
- 1000 µL disposable tips (List No. 04J71-10)
- 200 µL disposable tips (List No. 04J71-17)
- Vortex mixer
- USP grade 190 to 200 proof ethanol (95 to 100% ethanol). **Do not use ethanol that contains denaturants.**
- Abbott Optical Adhesive Covers (List No. 04J71-75)
- Abbott Adhesive Cover Applicator (List No. 9K32-01)
- Abbott 96-Deep-Well Plate (List No. 04J71-30)
- Abbott Splash-Free Support Base (List No. 09K31-01)
- 13 mm Sample Racks
- 1.5 mL Reaction Vessels and Output Tubes (1.5 mL screw top microfuge tubes and caps, List No. 4J71-50 or equivalent)

Reagent Preparation Area

- Abbott 96-Well Optical Reaction Plate (List No. 04J71-70)
- Abbott Splash-Free Support Base (List No. 09K31-01)
- Calibrated precision pipettes capable of delivering 10 µL to 1000 µL
- 20 µL to 1000 µL aerosol barrier pipette tips for precision pipettes
- Single-use DNase-free tube or container

Materials for Abbott m2000sp (Assay Protocol III)

Sample Preparation Area

- Abbott m2000sp instrument with Software Version 3.0 or higher
- Abbott mSample Preparation System_{DNA} (List No. 06K12-24)

NOTE: One kit is sufficient to complete 192 HPV sample preparations.

- 5 mL Reaction Vessels (List No. 4J71-20)
 - Calibrated precision pipettes capable of delivering 10 µL to 1000 µL
 - 20 µL to 1000 µL aerosol barrier pipette tips for precision pipettes
 - Sample input tubes (refer to **ASSAY PROTOCOL III** section for details)
 - 1000 µL disposable tips (List No. 04J71-10)
 - 200 µL disposable tips (List No. 04J71-17)
 - Vortex mixer
 - USP grade 190 to 200 proof ethanol (95 to 100% ethanol). **Do not use ethanol that contains denaturants.**
 - Abbott Optical Adhesive Covers (List No. 04J71-75)
 - Abbott Adhesive Cover Applicator (List No. 9K32-01)
 - Abbott Splash-Free Support Base (List No. 09K31-01)
 - Master Mix Tube (List No. 04J71-80)
 - 200 mL Reagent Vessels (List No. 4J71-60)
 - Abbott 96-Deep-Well Plate (List No. 04J71-30)
 - Abbott 96-Well Optical Reaction Plate (List No. 04J71-70)
 - 13 mm Sample Racks
- **Materials for Abbott m2000rt**
- Abbott m2000rt instrument with Software Version 3.0 or higher
 - Abbott m2000rt Optical Calibration Kit (List No. 4J71-93)

Other Materials

- Biological safety cabinet approved for working with infectious materials
- Sealable plastic bags
- DNase-free water[†]
- Microcentrifuge Tubes[†]
- Cotton Tip Applicators (Puritan or equivalent)[†]

[†] **NOTE: These three items are used in the procedure for Monitoring the Laboratory for the Presence of Contamination. Refer to the QUALITY CONTROL PROCEDURES section of the package insert.**

Procedural Precautions

- Read the instructions in the package insert carefully before processing samples.
- Do not use kits or reagents beyond expiration date.
- Control kit lots and amplification reagent kit lots can be used interchangeably. Components contained within a kit are intended to be used together. For example, do not use the negative control from control kit lot X with the positive control from control kit lot Y.

- Amplification Reagent components (enzyme, oligonucleotide reagent and activation reagent) and Controls are for single-use only and should be discarded after use. Use new reagent vessels and new reaction vessels, for every new Abbott RealTime HR HPV assay run. At the end of each run, discard all these remaining reagents as stated in the following manuals: the Abbott m24sp Operations Manual, the Abbott m2000sp Operations Manual, and **Manual Sample Preparation Using the Abbott mSample Preparation System_{DNA} for RealTime High Risk HPV.**
- The Abbott RealTime HR HPV Controls must be processed with the specimens to be tested. The use of the Abbott RealTime HR HPV Controls is integral to the performance of the Abbott RealTime HR HPV assay. Refer to the **QUALITY CONTROL PROCEDURES** section in the package insert for details.
- Use only USP grade 190 to 200 proof ethanol (95 to 100% ethanol) to prepare the mWash 2_{DNA} sample preparation reagent. **Do not use ethanol that contains denaturants.**
- Use aerosol barrier pipette tips or disposable pipettes only one time when pipetting. To prevent contamination to the pipette barrel while pipetting, care should be taken to avoid touching the pipette barrel to inside of the sample tube or container. The use of extended aerosol barrier pipette tips is recommended.
- Replace any empty or partially used 200 µL and 1000 µL disposable tips on the Abbott m2000sp or Abbott m24sp with full trays before every run. Refer to the Abbott m2000sp and Abbott m24sp Operations Manuals, Operating Instructions section.
- Monitoring procedures for the presence of amplification product can be found in the **QUALITY CONTROL PROCEDURES** section in the package insert.
- To reduce the risk of nucleic acid contamination, clean and disinfect spills of specimens, reagents and controls by using a detergent solution followed by a tuberculocidal disinfectant such as 1.0% (v/v) sodium hypochlorite or other suitable disinfectant.

ASSAY PROTOCOL I: MANUAL SAMPLE PREPARATION METHOD AND ABBOTT m2000rt INSTRUMENT

Refer to the **WARNINGS AND PRECAUTIONS** section of the package insert before preparing samples.

1. Vortex each specimen for 15 to 20 seconds. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial. Immediately transfer 400 µL of each specimen to a reaction tube.
NOTE: Ensure that the volume of the cell pellet SurePath samples after cytological processing is approximately 2.8 mL. The sample volume must be adjusted to 6 mL using SurePath Preservative Fluid prior to vortexing and transferring.
2. Thaw control reagents at 15 to 30°C or at 2 to 8°C; see **QUALITY CONTROL PROCEDURES** section of the package insert.
 - Vortex each assay control for 15 to 20 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
 - Once thawed, assay controls can be stored at 2 to 8°C for up to 24 hours before use.
3. Thaw amplification reagents at 15 to 30°C or at 2 to 8°C and store at 2 to 8°C until required for the amplification master mix procedure.
 - Once thawed, the amplification reagents can be stored at 2 to 8°C for up to 24 hours if not used immediately.

NOTE: A maximum of 96 reactions can be performed per run.

For up to 24 reactions use: 1 tube of positive control, 1 tube of negative control, 1 amplification reagent pack, and 1 set of Abbott mSample Preparation System_{DNA} reagents.

For 25 to 48 reactions use: 1 tube of positive control, 1 tube of negative control, 2 amplification reagent packs, and 1 set of Abbott mSample Preparation System_{DNA} reagents.

For 49 to 72 reactions use: 1 tube of positive control, 1 tube of negative control, 3 amplification reagent packs, 1 bottle of mMicroparticle_{DNA} and mLysis_{DNA} Buffer, and 2 bottles of mWash 1_{DNA} Buffer, mWash 2_{DNA} Buffer and mElution Buffer_{DNA}.

For 73 to 96 reactions use: 1 tube of positive control, 1 tube of negative control, 4 amplification reagent packs, 1 bottle of mMicroparticle_{DNA} and mLysis_{DNA} Buffer, and 2 bottles of mWash 1_{DNA} Buffer, mWash 2_{DNA} Buffer and mElution Buffer_{DNA}.

Sample Preparation Area

- Refer to the Extraction Protocol section of **Manual Sample Preparation Using the Abbott mSample Preparation System_{DNA} for RealTime High Risk HPV** for sample preparation procedures.

NOTE: Abbott mSample Preparation System_{DNA} reagents can be used up to 3 times within 14 days for a total of 48 samples when stored tightly capped at 15 to 30°C. If reusing the Abbott mSample Preparation System_{DNA} reagents, mark the mWash 2_{DNA} bottle to indicate that ethanol has already been added. Once prepared, do not add more ethanol to the mWash 2_{DNA} bottle at any time. If reusing the Abbott mSample Preparation System_{DNA} reagents, after removing the caps from all the Abbott mSample Preparation System_{DNA} reagents, store the caps on a clean, absorbent surface for recapping after the run. **NOTE:** The assembly of the amplification master mix and sample eluates into the Abbott 96-Well Optical Reaction Plate (step 12) must be initiated within 1 hour after completion of Sample Preparation.

Amplification Area

- Switch on and initialize the Abbott m2000rt. The Abbott m2000rt requires a 15-minute warm-up prior to starting a run. Refer to the Abbott m2000rt Operations Manual, Operating Instructions section.
- Create the Abbott m2000rt test order. Refer to the Operating Instructions section of the Abbott m2000rt Operations Manual. From the Protocol screen, select the appropriate application file.

Reagent Preparation Area

NOTE: All reagent preparation must take place in the dedicated Reagent Preparation Area. Refer to the Contamination Precautions section of the package insert before preparing reagents. Change gloves before handling the amplification reagents.

- Prepare the amplification master mix.
 - Each amplification reagent pack supports up to 24 reactions.
 - Prior to opening the amplification reagents, ensure that the contents of the amplification reagent pack are at the bottom by tapping the amplification reagent pack in an upright position on the bench to bring the liquid to the bottom of the vials.
 - Identify the amplification reagents as follows:
 - Activation Reagent (Reagent 1): clear bottle, teal cap
 - Oligonucleotide Reagent (Reagent 2): black bottle, white cap
 - DNA Polymerase (Reagent 3): clear bottle, white cap
 - Remove and discard caps.
 - Prepare the master mix by using a **PIPETTE DEDICATED FOR REAGENT USE ONLY** to add 278 µL of the HPV Activation Reagent (Reagent 1) and 402 µL of the HPV Oligonucleotide Reagent (Reagent 2) together in the DNA Polymerase bottle (Reagent 3). Mix the Enzyme vial containing the reaction mixture (master mix) by gently pipetting up and down 6 times. Avoid creating foam.
 - If performing 25 to 48 reactions, prepare the amplification master mix from 2 amplification reagent packs. If performing 49 to 72 reactions, prepare the amplification master mix from 3 amplification reagent packs. If performing 73 to 96 reactions, prepare the amplification master mix from 4 amplification reagent packs.

NOTE: The Abbott m2000rt protocol (step 14) must be initiated within 1 hour of the addition of the activation reagent into the DNA Polymerase bottle (step 7).

- Pipette the contents of the master mix from the enzyme bottle(s) into a single-use DNase-free tube. Mix by gently pipetting up and down 6 times. Avoid creating foam.
- Prior to addition of master mix and sample, insert an Abbott 96-Well Optical Reaction Plate onto an Abbott Splash-Free Support Base to prevent contamination.
 - Contamination of the bottom of the Abbott 96-Well Optical Reaction Plate with fluorescent materials could potentially interfere with the HPV assay. The Abbott 96-Well Optical Reaction Plate should be held and transported with the Abbott Splash-Free Support Base to minimize contamination.
- Using a **DEDICATED PIPETTE**, dispense 25 µL aliquots of the amplification master mix into each well of the Abbott 96-Well Optical Reaction Plate that will be used depending on the number of samples to be run, including controls. A calibrated repeat pipettor may be used. Visually verify that 25 µL has been dispensed into each well.

NOTE: Remaining activated master mix can be recapped and stored at –10°C or colder for up to 14 days and reused at a later time if the volume is sufficient. The activated master mix should not undergo more than 3 freeze/thaw cycles. The frozen master mix can be thawed at room temperature for up to 1 hour prior to the initiation of amplification and detection on the Abbott m2000rt.

- Transfer the Abbott 96-Well Optical Reaction Plate on the Abbott Splash-Free Support Base to the Sample Preparation Area.

Sample Preparation Area

- In the Sample Preparation Area, transfer 25 µL of sample eluate to the Abbott 96-Well Optical Reaction Plate on the Abbott Splash-Free Support Base. Use a separate pipette tip for each sample eluate transfer. Visually verify that a total of 50 µL has been dispensed into each well.
- Seal the Abbott 96-Well Optical Reaction Plate according to the instructions in the Abbott m2000rt Operations Manual.

Amplification Area

- Place the Abbott 96-Well Optical Reaction Plate in the Abbott m2000rt and initiate the Abbott RealTime HR HPV protocol as described in the Abbott m2000rt Operations Manual, Operating Instructions section. At the completion of the run, assay results are reported on the Abbott m2000rt. Refer to the **RESULTS** section of the package insert for further details.
- After the Abbott m2000rt has completed the amplification and detection protocol, remove the Abbott 96-Well Optical Reaction Plate and dispose of according to the instructions in the **Contamination Precautions** section of the package insert. Place the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the Abbott m2000rt Operations Manual along with the gloves used to handle the plate.

Post Processing Procedures

- Refer to the Clean Up section of **Manual Sample Preparation Using the Abbott mSample Preparation System_{DNA} for RealTime High Risk HPV**.
- Clean the Abbott Splash-Free Support Base before next use, according to the Abbott m2000rt Operations Manual.

ASSAY PROTOCOL II: ABBOTT m24sp AND ABBOTT m2000rt INSTRUMENTS

Refer to the **WARNINGS AND PRECAUTIONS** section of the package insert before preparing samples.

- Vortex each specimen for 15 to 20 seconds. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial. Immediately transfer the specimens to the sample input tubes.

NOTE: Ensure that the volume of the cell pellet SurePath samples after cytological processing is approximately 2.8 mL. The sample volume must be adjusted to 6 mL using SurePath Preservative Fluid prior to vortexing and transferring.

- For specimens collected in PreservCyt Solution or SurePath Preservative Fluid, to ensure that 400 µL of each specimen is transferred by the Abbott m24sp to the reaction vessel:
 - transfer a minimum of 500 µL of each specimen if using Master Mix Tubes or Abbott Transport Tubes as sample input tubes.
 - transfer a minimum of 700 µL of each specimen if using 5 mL Reaction Vessels or any other 13 mm round bottom non-skirted tubes as sample input tubes.
 - For specimens collected with the Abbott Cervi-Collect Specimen Collection Kit, load the tubes without cap directly on the Abbott m24sp (these specimens do not require a transfer).
- Thaw control reagents at 15 to 30°C or at 2 to 8°C; see **QUALITY CONTROL PROCEDURES** section of the package insert.
 - Vortex each assay control for 15 to 20 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
 - Once thawed, assay controls can be stored at 2 to 8°C for up to 24 hours before use.

- Thaw amplification reagents at 15 to 30°C or at 2 to 8°C and store at 2 to 8°C until required for the amplification master mix procedure.
 - Once thawed, the amplification reagents can be stored at 2 to 8°C for up to 24 hours if not used immediately.

NOTE: A maximum of 24 reactions can be performed per run. For up to 24 reactions, use 1 tube of positive control, 1 tube of negative control, 1 amplification reagent pack, and 1 set of Abbott *mSample Preparation System*_{DNA} reagents.

Sample Preparation Area

- Place the controls and patient specimens into the Abbott *m24sp* sample rack, as described in the Abbott *m24sp* Operations Manual, Operating Instructions section.

CAUTION: Use only 13 mm sample racks. Do NOT skip any positions in a sample rack. Load specimens and controls into the 13 mm sample racks in consecutive positions, starting with the third position in the first sample rack. Fill all positions in each sample rack without skipping any positions before loading specimens into the next sample rack.

Insert specimen and control tubes into sample racks carefully to avoid splashing. Ensure that each tube is placed securely in the sample rack so that the bottom of the tube reaches the inside bottom of the rack.

Load filled sample racks onto the Abbott *m24sp* in consecutive sample rack positions, with the first rack farthest to the right on the worktable, and if needed, the second rack to the left of the first rack.

- Open the Abbott *mSample Preparation System*_{DNA} reagent pack. Prepare the *mWash 2*_{DNA} by adding 70 mL of USP grade 190 to 200 proof ethanol (95 to 100% ethanol) to the *mWash 2*_{DNA} bottle as described in the Abbott *mSample Preparation System*_{DNA} product information. Do not use ethanol that contains denaturants. Gently invert each reagent bottle to ensure a homogenous solution. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.

NOTE: Abbott *mSample Preparation System*_{DNA} reagents can be used up to 3 times within 14 days for a total of 48 samples when stored tightly capped at 15 to 30°C. If reusing the Abbott *mSample Preparation System*_{DNA} reagents, mark the *mWash 2*_{DNA} bottle to indicate that ethanol has already been added. Once prepared, do not add more ethanol to the *mWash 2*_{DNA} bottle at any time.

- Initiate the Abbott *m24sp* protocol as described in the Abbott *m24sp* Operations Manual, Operating Instruction section. From the Protocol screen, select the appropriate script to run the HPV assay depending on desired output vessels (*m24sp_HP*_{DNA}_Tube for 1.5 mL tubes or *m24sp_HP*_{DNA}_DWP for 96-Deep-Well Plate). When prompted by the instrument, vigorously mix or vortex the *mMicroparticle*_{DNA} bottle until the *mMicroparticles*_{DNA} are fully resuspended. Put the *mMicroparticle*_{DNA} bottle on the deck of the instrument in the designated position.

NOTE: If reusing the Abbott *mSample Preparation System*_{DNA} reagents, after removing the caps from all the Abbott *mSample Preparation System*_{DNA} reagents, store the caps on a clean, absorbent surface for recapping after the run.

NOTE: The assembly of the amplification master mix and sample eluates into the Abbott 96-Well Optical Reaction Plate (step 14) must be initiated within 1 hour after completion of Sample Preparation.

Amplification Area

- Switch on and initialize the Abbott *m2000rt*. The Abbott *m2000rt* requires a 15-minute warm-up prior to starting a run. Refer to the Abbott *m2000rt* Operations Manual, Operating Instructions section.
- Create the Abbott *m2000rt* test order. Refer to the Operating Instructions section of the Abbott *m2000rt* Operations Manual. From the Protocol screen, select the appropriate application file.

Reagent Preparation Area

NOTE: All reagent preparation must take place in the dedicated Reagent Preparation Area. Refer to the Contamination Precautions section of the package insert before preparing reagents. Change gloves before handling the amplification reagents.

- Prepare the amplification master mix.
 - Each amplification reagent pack supports up to 24 reactions.
 - Prior to opening the amplification reagents, ensure that the contents of the amplification reagent pack are at the bottom by tapping the amplification reagent pack in an upright position on the bench to bring the liquid to the bottom of the vials.

- Identify the amplification reagents as follows:
 - Activation Reagent (Reagent 1): clear bottle, teal cap
 - Oligonucleotide Reagent (Reagent 2): black bottle, white cap
 - DNA Polymerase (Reagent 3): clear bottle, white cap
- Remove and discard caps.
- Prepare the master mix by using a **PIPETTE DEDICATED FOR REAGENT USE ONLY** to add 278 µL of the HPV Activation Reagent (Reagent 1) and 402 µL of the HPV Oligonucleotide Reagent (Reagent 2) together in the DNA Polymerase bottle (Reagent 3). Mix the enzyme vial containing the reaction mixture (master mix) by gently pipetting up and down 6 times. Avoid creating foam.

NOTE: The Abbott *m2000rt* protocol (step 16) must be initiated within 1 hour of the addition of the activation reagent into the DNA Polymerase reagent bottle (step 9).

- Pipette the contents of the master mix from the enzyme bottle(s) into a single-use DNase-free tube. Mix by gently pipetting up and down 6 times. Avoid creating foam.
- Prior to addition of master mix and sample, insert an Abbott 96-Well Optical Reaction Plate onto an Abbott Splash-Free Support Base to prevent contamination.
 - Contamination of the bottom of the Abbott 96-Well Optical Reaction Plate with fluorescent materials could potentially interfere with the HPV assay. The Abbott 96-Well Optical Reaction Plate should be held and transported with the Abbott-Splash Free Support Base to minimize contamination.

- Using a **DEDICATED PIPETTE**, dispense 25 µL aliquots of the amplification master mix into each well of the Abbott 96-Well Optical Reaction Plate that will be used depending on the number of samples to be run, including controls. A calibrated repeat pipettor may be used. Visually verify that 25 µL has been dispensed into each well.

NOTE: Remaining activated master mix can be recapped and stored at -10°C or colder for up to 14 days and reused at a later time if the volume is sufficient. The activated master mix should not undergo more than 3 freeze/thaw cycles. The frozen master mix can be thawed at room temperature for up to 1 hour prior to the initiation of amplification and detection on the Abbott *m2000rt*.

- Transfer the Abbott 96-Well Optical Reaction Plate on the Abbott Splash-Free Support Base to the Sample Preparation Area.

Sample Preparation Area

- In the Sample Preparation Area, transfer 25 µL of sample eluate to the Abbott 96-Well Optical Reaction Plate on the Abbott Splash-Free Support Base. Use a separate pipette tip for each sample eluate transfer. Visually verify that a total of 50 µL has been dispensed into each well.
- Seal the Abbott 96-Well Optical Reaction Plate according to the instructions in the Abbott *m2000rt* Operations Manual.

Amplification Area

- Place the Abbott 96-Well Optical Reaction Plate in the Abbott *m2000rt* and initiate the Abbott RealTime HR HPV protocol as described in the Abbott *m2000rt* Operations Manual, Operating Instructions section. At the completion of the run, assay results are reported on the Abbott *m2000rt*. Refer to the RESULTS section of the package insert for further details.
- After the Abbott *m2000rt* has completed the amplification and detection protocol, remove the Abbott 96-Well Optical Reaction Plate and dispose of according to the instructions in the Contamination Precautions section of the package insert. Place the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.

Post Processing Procedures

- At the end of each run, remove and discard all remaining reagents from the Abbott *m24sp* worktable as stated in the Abbott *m24sp* Operations Manual.
- Decontaminate and dispose of all specimens, reagents (except for amplification master mix when applicable), and other potentially contaminated materials in accordance with local, state, and federal regulations.
- Clean the Abbott Splash-Free Support Base before next use, according to the Abbott *m2000rt* Operations Manual.

ASSAY PROTOCOL III: ABBOTT *m2000sp* AND ABBOTT *m2000rt* INSTRUMENTS

Refer to the **WARNINGS AND PRECAUTIONS** section of the package insert before preparing samples.

1. Vortex each specimen for 15 to 20 seconds. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial. Immediately transfer the specimens to the sample input tubes.

NOTE: Ensure that the volume of the cell pellet SurePath samples after cytological processing is approximately 2.8 mL. The sample volume must be adjusted to 6 mL using SurePath Preservative Fluid prior to vortexing and transferring.

- For specimens collected in PreservCyt Solution or SurePath Preservative Fluid, to ensure that 400 µL of each specimen is transferred by the Abbott *m2000sp* to the reaction vessel:
 - transfer a minimum of 500 µL of each specimen if using Master Mix Tubes or Abbott Transport Tubes as sample input tubes.
 - transfer a minimum of 700 µL of each specimen if using 5 mL Reaction Vessels or any other 13 mm round bottom non-skirted tubes as sample input tubes.
 - For specimens collected with the Abbott Cervi-Collect Specimen Collection Kit, load the tubes without cap directly on the Abbott *m2000sp* (these specimens do not require a transfer).
2. Thaw control reagents at 15 to 30°C or at 2 to 8°C; see **QUALITY CONTROL PROCEDURES** section of the package insert.
 - Vortex each assay control for 15 to 20 seconds before use. Ensure that the contents of each vial are at the bottom after vortexing by tapping the vials on the bench to bring liquid to the bottom of the vial.
 - Once thawed, assay controls can be stored at 2 to 8°C for up to 24 hours before use.
 3. Thaw amplification reagents at 15 to 30°C or at 2 to 8°C and store at 2 to 8°C until required for the amplification master mix procedure.
 - Once thawed, the amplification reagents can be stored at 2 to 8°C for up to 24 hours if not used immediately.

NOTE: A maximum of 96 reactions can be performed per run.

For up to 24 reactions use: 1 tube of positive control, 1 tube of negative control, 1 amplification reagent pack, and 1 set of Abbott *mSample Preparation System_{DNA}* reagents.

For 25 to 48 reactions use: 1 tube of positive control, 1 tube of negative control, 2 amplification reagent packs, and 1 set of Abbott *mSample Preparation System_{DNA}* reagents.

For 49 to 72 reactions use: 1 tube of positive control, 1 tube of negative control, 3 amplification reagent packs, 1 bottle of *mMicroparticle_{DNA}* and *mLysis_{DNA}* Buffer, and two bottles of *mWash 1_{DNA}* Buffer, *mWash 2_{DNA}* Buffer and *mElution Buffer_{DNA}*.

For 73 to 96 reactions use: 1 tube of positive control, 1 tube of negative control, 4 amplification reagent packs, 1 bottle of *mMicroparticle_{DNA}* and *mLysis_{DNA}* Buffer, and 2 bottles of *mWash 1_{DNA}* Buffer, *mWash 2_{DNA}* Buffer and *mElution Buffer_{DNA}*.

NOTE: Abbott *mSample Preparation System_{DNA}* is for single-use only and should be discarded after use. Use newly opened reagents for every new Abbott RealTime HR HPV assay run.

4. Place the controls and the patient specimens into the Abbott *m2000sp* sample rack.

CAUTION: Use only 13 mm sample racks. Do NOT skip any positions in a sample rack. Load specimens and controls into the 13 mm sample racks in consecutive positions, starting with the first position in the first sample rack. Fill all positions in each sample rack without skipping any positions before loading specimens into the next sample rack.

Insert specimen and control tubes into sample racks carefully to avoid splashing. If used, bar codes on tube labels must face right for scanning. Ensure that each tube is placed securely in the sample rack so that the bottom of the tube reaches the inside bottom of the rack.

Load filled sample racks onto the Abbott *m2000sp* in consecutive sample rack positions, with the first rack farthest to the right on the worktable, and any additional rack progressively to the left of the first rack.

5. Open the Abbott *mSample Preparation System_{DNA}* reagent pack(s). Prepare the *mWash 2_{DNA}* by adding 70 mL of USP grade 190 to 200 proof ethanol (95 to 100% ethanol) to the *mWash 2_{DNA}* bottle as described in the Abbott *mSample Preparation System_{DNA}* product information. **Do not use ethanol that contains denaturants.** Gently invert each reagent bottle to ensure a homogenous solution and pour the contents into the appropriate reagent vessels per the Abbott *m2000sp* Operations Manual, Operating Instructions section. If crystals are observed in any of the reagent bottles upon opening, allow the reagent to equilibrate at room temperature until the crystals disappear. Do not use the reagents until the crystals have dissolved.

NOTE: Before pouring the *mMicroparticles_{DNA}* into the 200 mL reagent vessels, vigorously mix or vortex until the *mMicroparticles_{DNA}* are fully resuspended.

6. Initiate the sample extraction protocol as described in the Abbott *m2000sp* Operations Manual, Operating Instructions section.
7. While the Abbott *m2000sp* is performing sample preparation, switch on and initialize the Abbott *m2000rt*. The Abbott *m2000rt* requires a 15-minute warm-up prior to starting a run. Refer to the Abbott *m2000rt* Operations Manual, Operating Instructions section.

NOTE: Once sample preparation is completed, the master mix protocol should be started within 1 hour.

8. Load the amplification reagents and the master mix tube on the Abbott *m2000sp* worktable.
 - Prior to opening the amplification reagents, ensure that the contents of the amplification reagent pack(s) are at the bottom by tapping the amplification reagent pack(s) in an upright position on the bench to bring the liquid to the bottom of the vials.
 - Remove and discard vial caps.

NOTE: Change gloves before handling the amplification reagents.

9. Initiate the Abbott *m2000sp* Master Mix Addition protocol as described in the Abbott *m2000sp* Operations Manual, Operating Instructions section.
10. After the Abbott *m2000sp* has completed the addition of samples and amplification reagents, seal the Abbott 96-Well Optical Reaction Plate according to the instructions in the Abbott *m2000sp* Operations Manual.
 - Contamination of the bottom of the Abbott 96-Well Optical Reaction Plate with fluorescent materials could potentially interfere with the HPV assay. The Abbott 96-Well Optical Reaction Plate should be held and transported with the Abbott Splash-Free Support Base to minimize contamination.

NOTE: Within 1 hour of starting the master mix protocol, the sealed Abbott 96-Well Optical Reaction Plate should be transferred to the Abbott *m2000rt* to begin amplification/detection.

11. Place the Abbott 96-Well Optical Reaction Plate in the Abbott *m2000rt* and initiate the Abbott RealTime HR HPV assay protocol as described in the Abbott *m2000rt* Operations Manual, Operating Instructions section. At the completion of the run, assay results are reported on the Abbott *m2000rt*. Refer to the **RESULTS** section of the package insert for further details.
12. After the Abbott *m2000rt* has completed the amplification and detection protocol, remove the Abbott 96-Well Optical Reaction Plate and dispose of according to the instructions in the **Contamination Precautions** section of the package insert. Place the Abbott 96-Well Optical Reaction Plate in a sealable plastic bag and dispose of according to the Abbott *m2000rt* Operations Manual along with the gloves used to handle the plate.

Post Processing Procedures

1. At the end of each run, remove and discard all remaining reagents from the Abbott *m2000sp* worktable as stated in the Abbott *m2000sp* Operations Manual.
2. Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local, state, and federal regulations.
3. Clean the Abbott Splash-Free Support Base before next use, according to the Abbott *m2000rt* Operations Manual.

QUALITY CONTROL PROCEDURES

Abbott m2000rt Optical Calibration

Optical calibration of the Abbott m2000rt is required for the accurate measurement and discrimination of dye fluorescence during the Abbott RealTime HR HPV assay.

The following Abbott m2000rt Optical Calibration Plates are used to calibrate the Abbott m2000rt for the Abbott RealTime HR HPV assay:

- FAM Plate (Carboxyfluorescein)
- Cy5 Plate (Cyanine)
- NED Plate (ABI proprietary dye)
- ROX Plate (Carboxy-X-rhodamine)
- VIC Plate (Proprietary dye)

Refer to the Calibration Procedures section in the Abbott m2000rt Operations Manual for a detailed description of how to perform an Abbott m2000rt Optical Calibration.

Detection of Inhibition and/or Cell Inadequacy

The Abbott RealTime HR HPV assay detects the endogenous human beta globin sequence as Internal Control (IC) signal to evaluate cell adequacy, sample extraction and amplification efficiency. A flag or an error code is displayed when IC cycle number (CN) value of a sample or control exceeds the established range.

Negative and Positive Controls

A negative control and a positive control are required for every run to verify that the sample processing, the amplification, and the detection steps are performed correctly. The Abbott RealTime HR HPV controls need to be processed together with the samples prior to running the amplification portion of the assay.

The negative control is formulated with DNA containing IC sequence. The only signal detected for negative control should be the IC signal in the Cy5 channel. The positive control is formulated with DNA containing HPV 16, HPV 18, HPV 58 and IC sequences. All 4 signals (VIC signal for HPV 16, NED signal for HPV 18, FAM signal for HPV 58, and Cy5 signal for IC) should be detected in the positive control. A flag is displayed when a control result is out of range. If negative or positive controls are out of range, all of the samples and controls from that run must be reprocessed, beginning with sample preparation.

HR HPV must not be detected in the negative control. HR HPV detected in the negative control is indicative of contamination from other samples or amplified product introduced during sample preparation or during preparation of the Abbott 96-Well Optical Reaction Plate. To remove contamination, clean the Abbott m24sp or Abbott m2000sp and the Abbott m2000rt according to the Abbott m24sp, the Abbott m2000sp, and the Abbott m2000rt Operations Manuals. For manual sample preparation, clean the equipment according to the instructions in **Manual Sample Preparation Using the Abbott mSample Preparation System_{DNA} for RealTime High Risk HPV**. Following cleaning, repeat sample processing for controls and samples following the appropriate sample preparation protocol outlined in the package insert.

IC results for the negative control and positive control that are outside the validity limit indicate the occurrence of inhibition during sample preparation or during the amplification reaction steps of the assay. Repeat the processing for controls and samples following the appropriate sample preparation protocol outlined in the package insert.

Monitoring the Laboratory for the Presence of Contamination

It is recommended that the following procedure be done at least once a month to monitor laboratory surfaces and equipment for contamination. It is very important to test all areas that may have been exposed to processed samples and controls and/or amplification product. This includes routinely handled objects such as pipettes, Abbott m24sp, Abbott m2000sp, and Abbott m2000rt function keys, bench surfaces and other equipment that may be present in the work areas.

1. Add 0.8 mL DNase-free water to a new Master Mix Tube.
2. Saturate the cotton tip of an applicator (Puritan or equivalent) in the DNase-free water from the Master Mix Tube.
3. Using the saturated cotton tip of the applicator, wipe the area to be monitored using a sweeping motion. Place the applicator into the Master Mix Tube.
4. Swirl the cotton tip in DNase-free water 10 times, then press the applicator along the inside of the tube so that the liquid drains back into the solution at the bottom of the Master Mix Tube. Discard the applicator.

5. Cap the Master Mix Tube and vortex.
6. Remove the caps from the Master Mix Tubes and test the sample according to the appropriate assay procedure section of the package insert.
7. Contamination is indicated by the detection of HR HPV in the swab samples.
 - If contamination is present, the instrument will report "HR HPV Detected" (disregard IC flag if present).
 - If there is no contamination, the instrument will report "Not Detected" or no result will be displayed (disregard error codes 4951 or 4952 if present).
8. If contamination is detected on the equipment, follow the cleaning and decontaminating guidelines given in that equipment's operations manual. If HR HPV is detected on surfaces, clean the contaminated areas with 1.0% (v/v) sodium hypochlorite solution, followed by 70% ethanol or water.

NOTE: Chlorine solutions may pit equipment and metal. Use sufficient amounts or repeated applications of 70% ethanol until chlorine residue is no longer visible.
9. Repeat testing of the contaminated area by following steps 1 through 6.
10. If the presence of contamination is detected again, repeat steps 8 and 9 until no HR HPV amplification is detected.

RESULTS

The Abbott RealTime HR HPV assay is a qualitative assay. A minimum of 1 negative control and 1 positive control are required with each run. The negative control serves to verify that HR HPV DNA contamination of the negative control did not occur during the sample preparation and set-up of the amplification reaction. If HR HPV signal is detected for the negative control, the -QC flag is displayed next to all sample results for the run. Samples with the -QC flag may have been similarly contaminated with analyte during processing. If the negative control is not processed, the -QC flag is indicated next to all sample results for that run.

The IC signal in samples serves to confirm that each sample had sufficient cell input for accurate HR HPV detection and was processed correctly and to indicate whether inhibitors of amplification are present. If the IC is out of range (i.e. IC CN not generated or greater than or equal to a fixed cutoff cycle) and HR HPV is detected, the sample will have an interpretation of "HR HPV Detected". An IC flag will be reported next to the result. If the IC is out of range and HR HPV is not detected, no result will be reported and an error code will be generated. The sample with the error code must be retested starting with sample preparation.

For more information about error codes and flags, refer to the Abbott m2000rt Operations Manual Version 3.0 and Operations Manual Addendum Version 3.0.

Results Reporting

Three HPV signals corresponding to HPV 16, HPV 18 and Other HR HPV are evaluated for each sample. Each signal is either determined as "Detected" if the CN is less than a fixed assay cutoff cycle or is determined as "Not Detected" if the CN is not generated or the CN is greater than or equal to the assay cutoff cycle. All the detected signals (HPV 16, HPV 18 or Other HR HPV) are reported in the sample result with the respective CN values (in parenthesis after the target result). Samples with any of the 3 HR HPV signals detected will have an interpretation of "HR HPV Detected". Samples with all 3 HR HPV signals not detected will have an interpretation of "Not Detected".

Assay results and interpretations will look similar to the following examples:

Sample ID	Results	Interpretation	Explanation
1	HPV 16 (20.76)	HR HPV Detected	HPV 16 is detected with a CN of 20.76 HPV 18 and Other HR HPV are not detected
2	HPV 18 (21.20)	HR HPV Detected	HPV 18 is detected with a CN of 21.20 HPV 16 and Other HR HPV are not detected
3	Other HR HPV (14.48)	HR HPV Detected	Other HR HPV is detected with a CN of 14.48 HPV 16 and HPV 18 are not detected

Sample ID	Results	Interpretation	Explanation
4	HPV 16 (22.20); Other HR HPV (17.21)	HR HPV Detected	HPV 16 and Other HR HPV are detected with CN of 22.20 and 17.21, respectively HPV 18 is not detected
5	HPV 18 (18.67); Other HR HPV (15.88)	HR HPV Detected	HPV 18 and Other HR HPV are detected with CN of 18.67 and 15.88, respectively HPV 16 is not detected
6	HPV 16 (24.51); HPV 18 (23.11)	HR HPV Detected	HPV 16 and HPV 18 are detected with CN of 24.51 and 23.11, respectively Other HR HPV is not detected
7	HPV 16 (21.35); HPV 18 (22.60); Other HR HPV (19.45)	HR HPV Detected	HPV 16 and HPV 18 and Other HR HPV are detected with CN of 21.35, 22.60, and 19.45, respectively
8	Not Detected	Not Detected	HR HPV is not detected

LIMITATIONS OF THE PROCEDURE

- For In Vitro Diagnostic Use Only.
- This method has been tested using clinically-collected PreservCyt and SurePath liquid pap and Abbott Cervi-Collect specimens. Performance with other specimen types has not been evaluated.
- Optimal performance of this test requires appropriate specimen collection, handling, and storage (refer to the **SPECIMEN COLLECTION AND HANDLING INSTRUCTIONS** section of the package insert).
- Use of the Abbott RealTime HR HPV assay is limited to personnel who have been trained on the use of the Abbott *m24sp* or Abbott *m2000sp* or manual sample preparation method for sample extraction and Abbott *m2000rt* for amplification and detection.
- The instruments and assay procedures reduce the risk of contamination by amplification product. However, nucleic acid contamination from the controls, specimens, and amplification product must be controlled by good laboratory practice and careful adherence to the procedures specified in the package insert.
- A negative result does not preclude the possibility of infection because results are dependent on appropriate specimen collection. Test results may be affected by improper specimen collection, technical error, or specimen mix-up.
- As with any diagnostic test, results from the Abbott RealTime HR HPV assay should be interpreted in conjunction with other clinical and laboratory findings.

SPECIFIC PERFORMANCE CHARACTERISTICS

Genotype Inclusivity and Partial Genotyping

The ability of the Abbott RealTime HR HPV assay to detect 14 HR HPV genotypes (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) and to distinguish HPV 16 and HPV 18 from the other 12 HR HPV genotypes was evaluated. Fifty-one samples containing HPV DNA targets from each of the 14 genotypes individually and in combinations were tested as listed in Table 1. Results from 51 samples that included 14 samples with single genotype, 25 samples with 2 genotypes and 12 samples with 3 genotypes were reported accurately; the presence or absence of HPV 16 and HPV 18 DNA was accurately determined in each case.

Table 1: Genotype Detection and Partial Genotyping Capability

Sample No.	HPV Genotype	Reported Result
1	HPV 16	HPV 16
2	HPV 18	HPV 18
3	HPV 31	Other HR HPV
4	HPV 33	Other HR HPV
5	HPV 35	Other HR HPV
6	HPV 39	Other HR HPV
7	HPV 45	Other HR HPV
8	HPV 51	Other HR HPV
9	HPV 52	Other HR HPV
10	HPV 56	Other HR HPV
11	HPV 58	Other HR HPV
12	HPV 59	Other HR HPV
13	HPV 66	Other HR HPV
14	HPV 68	Other HR HPV
15	HPV 16 and HPV 18	HPV 16; HPV 18
16	HPV 16 and HPV 31	HPV 16; Other HR HPV
17	HPV 16 and HPV 33	HPV 16; Other HR HPV
18	HPV 16 and HPV 35	HPV 16; Other HR HPV
19	HPV 16 and HPV 39	HPV 16; Other HR HPV
20	HPV 16 and HPV 45	HPV 16; Other HR HPV
21	HPV 16 and HPV 51	HPV 16; Other HR HPV
22	HPV 16 and HPV 52	HPV 16; Other HR HPV
23	HPV 16 and HPV 56	HPV 16; Other HR HPV
24	HPV 16 and HPV 58	HPV 16; Other HR HPV
25	HPV 16 and HPV 59	HPV 16; Other HR HPV
26	HPV 16 and HPV 66	HPV 16; Other HR HPV
27	HPV 16 and HPV 68	HPV 16; Other HR HPV
28	HPV 18 and HPV 31	HPV 18; Other HR HPV
29	HPV 18 and HPV 33	HPV 18; Other HR HPV
30	HPV 18 and HPV 35	HPV 18; Other HR HPV
31	HPV 18 and HPV 39	HPV 18; Other HR HPV
32	HPV 18 and HPV 45	HPV 18; Other HR HPV
33	HPV 18 and HPV 51	HPV 18; Other HR HPV
34	HPV 18 and HPV 52	HPV 18; Other HR HPV
35	HPV 18 and HPV 56	HPV 18; Other HR HPV
36	HPV 18 and HPV 58	HPV 18; Other HR HPV
37	HPV 18 and HPV 59	HPV 18; Other HR HPV
38	HPV 18 and HPV 66	HPV 18; Other HR HPV
39	HPV 18 and HPV 68	HPV 18; Other HR HPV
40	HPV 16 and HPV 18 and HPV 31	HPV 16; HPV 18; Other HR HPV
41	HPV 16 and HPV 18 and HPV 33	HPV 16; HPV 18; Other HR HPV
42	HPV 16 and HPV 18 and HPV 35	HPV 16; HPV 18; Other HR HPV
43	HPV 16 and HPV 18 and HPV 39	HPV 16; HPV 18; Other HR HPV
44	HPV16 and HPV18 and HPV 45	HPV 16; HPV 18; Other HR HPV
45	HPV 16 and HPV 18 and HPV 51	HPV 16; HPV 18; Other HR HPV
46	HPV 16 and HPV 18 and HPV 52	HPV 16; HPV 18; Other HR HPV
47	HPV 16 and HPV 18 and HPV 56	HPV 16; HPV 18; Other HR HPV
48	HPV 16 and HPV 18 and HPV 58	HPV 16; HPV 18; Other HR HPV
49	HPV 16 and HPV 18 and HPV 59	HPV 16; HPV 18; Other HR HPV
50	HPV 16 and HPV 18 and HPV 66	HPV 16; HPV 18; Other HR HPV
51	HPV 16 and HPV 18 and HPV 68	HPV 16; HPV 18; Other HR HPV

Clinical Sensitivity and Specificity in Referral Population: Disease Detection

A total of 512 PreservCyt liquid pap specimens from a referral population were tested with the Abbott RealTime HR HPV assay and the hc2 High-Risk HPV DNA Test (HC2). The clinical sensitivity and specificity for detection of disease were determined for both assays. The presence of disease was defined by a histology result of Cervical Intraepithelial Neoplasia (CIN) 2 or greater. For subjects lacking a histology evaluation, disease status was defined by a liquid based cytology (LBC) result of high grade squamous intraepithelial lesion (HSIL) or greater at enrollment in colposcopy clinics. Disease prevalence in this population was 24.6%. Of 126 disease positive subjects, 121 were detected by the Abbott RealTime HR HPV assay and 119 were detected by HC2. Of 386 disease negative specimens, 154 were not detected by the Abbott RealTime HR HPV assay and 147 were not detected by HC2. The sensitivity of the Abbott RealTime HR HPV assay for detection of disease was 96.0% and of HC2 was 94.4%. The specificity of the Abbott RealTime HR HPV assay in this referral population was 39.9% and of HC2 was 38.1% (Table 2).

Table 2: Clinical Performance for Detection of Disease in Referral Population

Test	Sensitivity (95% CI)	Specificity (95% CI)	Positive Predictive Value	Negative Predictive Value
Abbott RealTime HR HPV	96.0% (91.0-98.7%)	39.9% (35.0-45.0%)	34.3%	96.9%
HC2	94.4% (88.9-97.7%)	38.1% (33.2-43.1%)	33.2%	95.5%

A total of 128 subjects from this population had a cytology result of ASC-US. In this ASC-US population, the clinical sensitivity was 100% for both the Abbott RealTime HR HPV assay and HC2. The clinical specificity of the Abbott RealTime HR HPV assay was 46.2% and of HC2 was 45.2%.

Clinical Sensitivity and Specificity in Referral Population: High Risk HPV Detection

The sensitivity and specificity of the Abbott RealTime HR HPV assay for detection of HR HPV were evaluated by testing 517 PreservCyt liquid pap specimens collected from a referral population. The high risk HPV status of cervical specimens was determined by the concordance between the Abbott RealTime HR HPV assay and HC2 tests and by further analysis of the specimens with discordant results using LINEAR ARRAY HPV Genotyping Test (Linear Array). A total of 337 specimens were detected by both assays and 136 specimens were not detected by either assay. The results of 44 discordant specimens were resolved by Linear Array. Of the 363 HR HPV positive specimens, 354 were detected by the Abbott RealTime HR HPV assay and 346 were detected by HC2. Of the 154 HR HPV negative specimens, 153 were not detected by the Abbott RealTime HR HPV assay and 137 were not detected by HC2. The sensitivity of the Abbott RealTime HR HPV assay for detection of HR HPV was 97.5% and of HC2 assay was 95.3%. The specificity of the Abbott RealTime HR HPV assay was 99.4% and of HC2 assay was 89.0% (Table 3).

Table 3: Sensitivity and Specificity for Detection of HR HPV

Test	Sensitivity (95% CI)	Specificity (95% CI)
Abbott RealTime HR HPV	97.5% (95.3-98.9%)	99.4% (96.4-100%)
HC2	95.3% (92.6-97.2%)	89.0% (82.9-93.4%)

Clinical Specificity in General Screening Population with Normal Cytology (Age ≥ 30 years)

A total of 362 cytologically normal PreservCyt liquid pap specimens collected in a general screening population from women 30 years of age or older were tested with the Abbott RealTime HR HPV assay and HC2. The clinical specificity was determined for both assays. The Abbott RealTime HR HPV assay detected 4.1% of the specimens and HC2 detected 3.0% of the specimens. In this population, the clinical specificity of the Abbott RealTime HR HPV assay and HC2 were 95.9% and 97.0% respectively.

Table 4: Clinical Specificity in General Screening Population with Normal Cytology (Age ≥ 30 years)

Test	Number detected/ tested	Detection Rate	Specificity (95% CI)
Abbott RealTime HR HPV	15/362	4.1%	95.9% (93.3-97.7%)
HC2	11/362	3.0%	97.0% (94.6-98.5%)

Clinical Sensitivity and Specificity from Additional Studies on Referral Populations: Disease Detection

The clinical sensitivity and specificity for detection of disease among referral populations were determined for Abbott RealTime HR HPV assay in 4 studies²⁶⁻²⁹ in comparison with HC2. All specimens were collected in PreservCyt Solution. The results from peer-reviewed literature are summarized in Table 5.

Table 5: Clinical Performance for Detection of Disease in Referral Populations

Study	Number Positive ^a	Number Disease Negative ^a	Sensitivity (95% CI)		Specificity (95% CI)	
			Abbott RealTime HR HPV	HC2	Abbott RealTime HR HPV	HC2
1 ²⁶	229	473	97.8% (95.0-99.3%)	95.6% (92.1-97.9%)	32.8% (28.6-37.2%)	35.7% (31.4-40.2%)
2 ²⁷	39	76	90.0% (85.0-95.0%)	95.0% (91.0-99.0%)	50.0% (41.0-59.0%)	50.0% (41.0-59.0%)
3 ²⁸	359	740	93.3% (90.1-95.6%)	96.3% (93.8-98.0%)	27.3% (24.1-30.7%)	19.5% (16.7-22.6%)
4 ²⁹	156	163	92.4% (87.0-96.0%)	91.7% (86.3-95.5%)	61.7% (53.8-69.2%)	58.6% (50.6-66.3%)

^a Disease positive specimens were generally defined as having a histology result of CIN2 or greater. Disease negative specimens were generally defined as having a histology result of less than CIN2.

Clinical Sensitivity and Specificity from Additional Studies on ASC-US Populations: Disease Detection

The clinical sensitivity and specificity for detection of disease among patients with ASC-US or equivalent cytology results were determined for Abbott RealTime HR HPV assay in 2 studies^{26,30} in comparison with HC2. All specimens were collected in PreservCyt Solution. The results from peer-reviewed literature are summarized in Table 6.

Table 6: Clinical Performance for Detection of Disease in ASC-US Populations

Study	Number Positive	Number Disease Negative	Sensitivity (95% CI)		Specificity (95% CI)	
			Abbott RealTime HR HPV	HC2	Abbott RealTime HR HPV	HC2
1 ²⁶	52	141	96.2% (86.8-99.5%)	94.2% (84.1-98.8%)	33.3% (25.6-41.8%)	39.0% (30.9-47.6%)
2 ^{a,30}	37	240	97.3% ^b	97.4% ^b	39.6% ^b	33.6% ^b

^a Study subjects included in this data set had a cytology result of borderline dyskaryosis, which correlates with atypical squamous cells.³¹
^b The 95% CI range is not reported by the publication.³⁰

Clinical Sensitivity and Specificity in Screening Populations: Disease Detection

The clinical sensitivity and specificity for detection of disease among screening populations were determined for Abbott RealTime HR HPV assay in 3 studies³²⁻³⁴ in comparison with benchmark tests. All specimens were collected in PreservCyt Solution. The results from peer-reviewed literature are summarized in Table 7.

Table 7: Clinical Performance for Detection of Disease in Screening Populations (Age ≥ 30 years)

Study	Number Positive	Number Disease Negative	Sensitivity (95% CI)		Specificity (95% CI)	
			Abbott RealTime HR HPV	Benchmark Test ^a	Abbott RealTime HR HPV	Benchmark Test ^a
1 ³²	38	3,091	100% (86.5-100%)	97.4% (86.2-99.9%)	93.3% (92.4-94.2%)	91.8% (90.8-92.7%)
2 ³³	68	858	95.6% (87.2-98.6%)	98.5% (90.3-99.8%)	92.0% (90.0-93.5%)	91.8% (89.9-93.4%)
3 ³⁴	16	4,629 ^b	100% (79.4-100%)	100% (79.4-100%)	90.3% (89.4-91.1%)	88.8% (87.9-89.7%)

^a Benchmark test for Studies 1 and 3 was HC2. Benchmark test for Study 2 was GP5+/6+ PCR.

^b Based on the specimens tested with Abbott RealTime HR HPV.

Accuracy in Identification of HPV 16 and/or HPV 18 in Women with Cervical Disease

The performance of the Abbott RealTime HR HPV in identification of HPV 16 and/or HPV 18 in cervical disease (CIN2 or greater) is evaluated based on the results from a referral population.²⁶ Out of 229 specimens with cervical disease, 210 had a valid Abbott RealTime HR HPV result with an interpretation of "HR HPV Detected" and a valid Linear Array result that reported one or more of the high risk HPV genotypes targeted by Abbott RealTime HR HPV. The overall agreement for detection of HPV 16 and/or HPV 18 between the Abbott RealTime HR HPV and Linear Array tests was 100% (210/210).

Table 8: Genotyping Accuracy for HPV 16 and/or HPV 18

Linear Array	HPV 16 and/or HPV 18 Reported ^c	Abbott RealTime HR HPV	
		HPV 16 and/or HPV 18 Detected ^a	Other HR HPV Detected ^b
	Non-HPV 16/18 High Risk Genotype(s) Reported ^d	0	57

^a These specimens were detected for HPV 16 and/or HPV 18 signal(s) with or without Other HR HPV signal detected.

^b These specimens were not detected for HPV 16 or HPV 18 signal and detected for Other HR HPV signal.

^c These specimens were reported with HPV 16 and/or HPV 18 genotype(s) with or without non-HPV 16/18 high risk HPV genotype(s) reported.

^d These specimens were reported with 1 or more of the non-HPV 16/18 high risk HPV genotypes that are targeted by Abbott RealTime HR HPV. HPV 16 or HPV 18 was not reported.

Estimate of Relative Disease Risk Associated with Different Genotype Results

The relative risks of having cervical disease (CIN2 or greater) were estimated for HPV 16 and/or HPV 18 Detected vs. Other HR HPV Detected results based on data obtained in a referral population,²⁶ an ASC-US Population,²⁶ and a screening population (women 30 years of age or older).³²

Table 9: Relative Risk of Cervical Disease Associated with Different Genotype Results (HPV 16 and/or HPV 18 Detected vs Other HR HPV Detected)

Study	Relative Risk	95% CI
Referral Population	2.1	(1.7, 2.7)
ASC-US Population	2.6	(1.5, 4.6)
Screening Population (Age ≥ 30 years)	2.5	(1.4, 4.4)

Analytical Sensitivity for High Risk HPV Genotypes

Analytical sensitivity of the Abbott RealTime HR HPV assay was determined by testing HPV DNA from each of 14 HR HPV genotypes in the presence of human cellular DNA in PreservCyt Solution. Four hundred microliters of sample is used per assay. For each genotype, a minimum of 4 levels, with 9 replicates at each level were tested. Testing was performed with 3 lots of amplification reagents on 3 Abbott m2000 RealTime Systems.

Probit analysis determined that with a probability of greater than 95%, HPV 16, 18, 35, 39, 45, 51, 59, 66, and 68 can be detected at 500 copies per assay, HPV 31, 33, 52, and 56 can be detected at 2,000 copies per assay and HPV 58 at 5,000 copies per assay.

Analytical Specificity (Cross-reactivity)

A panel of bacteria, viruses and fungi were evaluated for potential cross-reactivity in the Abbott RealTime HR HPV assay (Table 10). The panel included 15 Low Risk HPV genotypes and other organisms that can be found in the female anogenital tract. Human cellular DNA was also evaluated for potential cross-reactivity. Each potential cross-reactant was spiked into HPV negative samples at concentrations (per 0.4 mL sample input) shown in Table 10. Purified nucleic acids were used except where noted. Cross-reactivity was not observed with any of the organisms tested.

Table 10: Cross-reactivity Panel

Organisms	Concentration	Organisms	Concentration
<i>Bacteroides fragilis</i>	10 ⁷ genomic copies	HPV 6	10 ⁷ genomic copies
<i>Candida albicans</i> ^a	10 ⁷ CFU	HPV 11	10 ⁷ genomic copies
<i>Chlamydia trachomatis</i> ^a	10 ⁷ EBs	HPV 13	10 ⁷ genomic copies
<i>Corynebacterium genitalium</i>	10 ⁷ genomic copies	HPV 26	10 ⁷ genomic copies
<i>Enterobacter cloacae</i>	10 ⁷ genomic copies	HPV 30	10 ⁷ genomic copies
<i>Enterococcus faecalis</i>	10 ⁷ genomic copies	HPV 32	10 ⁷ genomic copies
<i>Escherichia coli</i>	10 ⁷ genomic copies	HPV 40	10 ⁷ genomic copies
<i>Gardnerella vaginalis</i>	10 ⁷ genomic copies	HPV 42	10 ⁷ genomic copies
<i>Haemophilis ducreyi</i>	10 ⁷ genomic copies	HPV 43	10 ⁷ genomic copies
<i>Lactobacillus acidophilus</i>	10 ⁷ genomic copies	HPV 44	10 ⁷ genomic copies
<i>Mycoplasma genitalium</i>	10 ⁷ genomic copies	HPV 53	10 ⁷ genomic copies
<i>Mycoplasma hominis</i>	10 ⁷ genomic copies	HPV 54	10 ⁷ genomic copies
<i>Neisseria gonorrhoeae</i>	10 ⁷ genomic copies	HPV 55	10 ⁷ genomic copies
<i>Neisseria meningitides</i>	10 ⁷ genomic copies	HPV 57	10 ⁷ genomic copies
<i>Proteus mirabilis</i>	10 ⁷ genomic copies	HPV 61	10 ⁷ genomic copies
<i>Staphylococcus aureus</i>	10 ⁷ genomic copies	HSV-I	10 ⁷ genomic copies
<i>Staphylococcus epidermidis</i>	10 ⁷ genomic copies	HSV-II	10 ⁷ genomic copies
<i>Streptococcus pneumoniae</i>	10 ⁷ genomic copies	HBV	10 ⁷ genomic copies
<i>Trichomonas vaginalis</i>	10 ⁶ genomic copies	HCV ^b	10 ⁶ viral RNA copies
<i>Ureaplasma urealyticum</i>	10 ⁷ genomic copies	HIV-1	10 ⁶ viral RNA copies
Human Cellular DNA	10 ⁷ genomic copies		

^a Cultured microorganisms.

^b Clinical specimen

Reproducibility

The reproducibility of the Abbott RealTime HR HPV assay was evaluated by testing a panel of 20 well-characterized clinical specimen pools (10 HR HPV positive and 10 HR HPV negative). The 20 panel members were tested by 2 operators. Each operator, using a unique combination of reagent lot and instrument pair, tested 2 replicates of each panel member per day for 4 days for a total of 8 replicates. Percent (%) Agreement results, based on comparison of the Abbott RealTime HR HPV results to expected results, for each panel member individually and for overall negative and positive panels are shown in Table 11. For positive samples, results for each HPV signal (HPV 16, HPV 18, and Other HR HPV) were accurately reported for all replicates. The overall agreement for 319 results compared with expected results was 100%. The agreement for 159 comparisons between the 2 operators using 2 different reagent lots and 2 instruments was 100%.

Table 11: Reproducibility

Panel No.	Expected Result	N	% Detected	% Agreement
1	Not Detected	16	0	100
2	Not Detected	16	0	100
3	Not Detected	16	0	100
4	Not Detected	16	0	100
5	Not Detected	16	0	100
6	Not Detected	16	0	100
7	Not Detected	16	0	100
8	Not Detected	16	0	100
9	Not Detected	16	0	100
10	Not Detected	16	0	100
11	HR HPV Detected (Other HR HPV)	16	100	100
12	HR HPV Detected (HPV 16; HPV 18)	16	100	100
13	HR HPV Detected (HPV 16)	16	100	100
14	HR HPV Detected (HPV 16; Other HR HPV)	16	100	100
15	HR HPV Detected (Other HR HPV)	16	100	100
16	HR HPV Detected (Other HR HPV)	16	100	100
17	HR HPV Detected (Other HR HPV)	15 ^a	100	100
18	HR HPV Detected (Other HR HPV)	16	100	100
19	HR HPV Detected (Other HR HPV)	16	100	100
20	HR HPV Detected (Other HR HPV)	16	100	100
Negative Samples (panels 1-10)		160	0	100
Positive Samples (panels 11-20)		159	100	100

^a Invalid reaction was excluded from the analysis.

Reproducibility Between Manual, Abbott *m24sp*, and Abbott *m2000sp* Sample Preparation Methods

Three different sample processing options are available for the Abbott RealTime HR HPV assay: manual, Abbott *m24sp*, and Abbott *m2000sp*. The reproducibility between Abbott *m2000sp* and manual sample preparation methods and between Abbott *m2000sp* and Abbott *m24sp* was determined by testing separate aliquots of the same cervical specimens using these different sample preparation methods. For each comparison 110 PreservCyt liquid pap specimens were tested. Agreements between Abbott *m2000sp* and manual sample preparation methods (Table 12) and between Abbott *m2000sp* and Abbott *m24sp* (Table 13) were both 100%.

Table 12: Agreement Between Abbott *m2000sp* and Manual Sample Preparation

		Manual Sample Preparation	
		Detected	Not Detected
Abbott <i>m2000sp</i>	Detected	55	0
	Not Detected	0	55

Table 13: Agreement Between Abbott *m2000sp* and Abbott *m24sp*

		Abbott <i>m24sp</i>	
		Detected	Not Detected
Abbott <i>m2000sp</i>	Detected	55	0
	Not Detected	0	55

Potentially Interfering Substances

The potential for interference in the Abbott RealTime HR HPV assay was assessed with substances that may be present in cervical specimens. HR HPV negative samples and HR HPV positive samples were tested in the presence or absence of each of the substances listed in Table 14. Blood and mucus were spiked into PreservCyt solution at a concentration of 5%, all other substances at a concentration of 0.5%. Interference was not observed with any of the substances tested.

Table 14: Potentially Interfering Substances Tested

Blood
Mucus
CLOTRIMAZOLE Vaginal Cream (2%)
Deifen Vaginal Contraceptive Foam
Gynecort 1% Hydrocortisone Anti-itch Creme
K-Y Jelly
Lubrin
MetroGel-Vaginal
Miconazole Nitrate Suppository
Monistat-1 Day or Night Treatment
Norforms Deodorant Suppositories
Terazol-3 Vaginal Cream
Vagi-gard Povidone Iodine Medicated Douche
Vagisil Anti-Itch Creme
Vagisil Intimate Lubricant
Yeast Gard Homeopathic Vaginal Suppositories
Zovirax Cream (Acyclovir) 5%

Performance for Detection of High Risk HPV with Specimens Collected Using Abbott Cervi-Collect Specimen Collection Kit

Specimens collected with the Abbott Cervi-Collect Specimen Collection Kit were tested with the Abbott RealTime HR HPV assay. Specimens collected in PreservCyt Solution from the same subjects were tested with the Abbott RealTime HR HPV assay and with HC2. A total of 153 paired specimens that had sufficient volume for all 3 tests were included in the analysis. The high risk HPV status of cervical specimens was determined by the concordance of the Abbott RealTime HR HPV and HC2 results, and by further analysis of the specimens with discordant results using Linear Array. Among the 70 HR HPV positive specimens, the detection rates were 92.9%, 98.6% and 84.3% for Abbott RealTime HR HPV using Cervi-Collect specimens, Abbott RealTime HR HPV using PreservCyt liquid pap specimens and HC2, respectively (Table 15). Among the 83 HR HPV negative specimens, the detection rates were 3.6%, 2.4% and 3.6% for Abbott RealTime HR HPV using Cervi-Collect specimens, Abbott RealTime HR HPV using PreservCyt liquid pap specimens and HC2, respectively (Table 15).

Table 15: HR HPV Detection

Test	HR HPV Positive (N=70)		HR HPV Negative (N=83)	
	Number detected	% Detected (95% CI)	Number detected	% Detected (95% CI)
Abbott RealTime HR HPV with Cervi-Collect	65	92.9 (84.1-97.6)	3	3.6 (0.8-10.2)
Abbott RealTime HR HPV with PreservCyt Liquid Pap	69	98.6 (92.3-100)	2	2.4 (0.3-8.4)
HC2 with PreservCyt Liquid Pap	59	84.3 (73.6-91.9)	3	3.6 (0.8-10.2)

The agreement in Abbott RealTime HR HPV results with specimens collected in Cervi-Collect versus specimens collected in PreservCyt Solution from the same patients was 94.4% (Table 16).

Table 16: Agreement between Cervi-Collect and PreservCyt Liquid Pap Specimens

	Abbott RealTime HR HPV Cervi-Collect	
	Detected	Not Detected
Abbott RealTime HR HPV	69	6
PreservCyt Liquid Pap	3	83

Agreement = 94.4% (152/161)

Performance for Detection of High Risk HPV with Specimens Collected in SurePath Preservative Fluid

To assess the performance of the Abbott RealTime HR HPV assay with specimens collected in SurePath Preservative Fluid, a total of 263 specimens collected in SurePath Preservative Fluid were tested with the Abbott RealTime HR HPV assay and the HC2 test. Both the sample from the original SurePath collection vial and the remaining cell pellet sample obtained after cytological processing were tested with the Abbott RealTime HR HPV assay. The cell pellet sample was tested with HC2 per manufacturer's instructions. The high risk HPV status of cervical specimens was determined by the concordance of the Abbott RealTime HR HPV and HC2 results, and by further analysis of the specimens with discordant results using Linear Array. Among the 138 HR HPV positive specimens, the detection rates were 98.6%, 97.1% and 99.3% for Abbott RealTime HR HPV using the sample from the original collection vial, Abbott RealTime HR HPV using the cell pellet sample and HC2, respectively (Table 17). Among the 125 HR HPV negative specimens, the detection rates were 0.0%, 0.0% and 13.6% for Abbott RealTime HR HPV using the sample from the original collection vial, Abbott RealTime HR HPV using the cell pellet sample and HC2, respectively (Table 17).

Table 17: HR HPV Detection

Test	HR HPV Positive (N=138)		HR HPV Negative (N=125)	
	Number detected	% Detected (95% CI)	Number detected	% Detected (95% CI)
Abbott RealTime HR HPV with SurePath samples from Original Collection Vial	136	98.6 (94.9-99.8)	0	0 (0.0-2.9)
Abbott RealTime HR HPV with SurePath samples from Cell Pellet	134	97.1 (92.7-99.2)	0	0 (0.0-2.9)
HC2 with SurePath samples from Cell Pellet	137	99.3 (96.0-100)	17	13.6 (8.1-20.9)

The agreement in Abbott RealTime HR HPV results with the SurePath sample from the original collection vial versus the cell pellet sample was 99.2% (Table 18).

Table 18: Agreement between SurePath Samples from Original Collection Vial and Cell Pellet

	Abbott RealTime HR HPV Cell Pellet	
	Detected	Not Detected
Abbott RealTime HR HPV Original Collection Vial	134	2
	0	129

Agreement = 99.2% (263/265)

BIBLIOGRAPHY

- Howley PM. Papillomaviridae: the viruses and their replication. In: Fields BN, Knipe DM, Howley PM, eds. *Virology*, 3rd ed. Philadelphia, Lippincott-Raven Publishers 1996:947-78.
- CDC. Genital HPV Infection - CDC Fact Sheet. 2008; <http://www.cdc.gov/std/HPV/STDFact-HPV.htm>.
- zur Hausen H. Papillomaviruses and cancer: from basic studies to clinical application. *Nat Rev Cancer*. 2002;2:342-50.
- Walboomers JMM, Jacobs MV, Manos MM, et al. Human papillomavirus is a necessary cause of invasive cervical cancer worldwide. *J Pathol* 1999;189:12-19.
- Snijders PJ, Steenbergen RD, Heideman DA, et al. HPV-mediated cervical carcinogenesis: concepts and clinical implications. *J Pathol*. 2006;208:152-64.
- Kjaer SK, van den Brule AJC, Paull G, et al. Type specific persistence of high risk human papillomavirus (HPV) as indicator of high grade cervical squamous intraepithelial lesions in young women: population based prospective follow up study. *BMJ* 2002;325:572-578.
- Cuschieri KS, Cubie HA, Whitley MW, et al. Persistent high risk HPV infection associated with development of cervical neoplasia in a prospective population study. *J Clin Pathol*. 2005;58:946-50.
- de Villiers EM, Fauquet C, Broker TR, et al. Classification of papillomaviruses. *Virology* 2004;324:17-27.
- IARC Monographs on the evaluation of carcinogenic risks to humans. Human Papillomaviruses. Lyon: *International Agency for Research on Cancer* 2007; Volume 90.
- Muñoz N, Bosch FX, de Sanjosé S, et al. Epidemiologic classification of human papillomavirus types associated with cervical cancer. *N Engl J Med*. 2003;348:518-27.
- Clifford GM, Smith JS, Plummer M, et al. Human papillomavirus types in invasive cervical cancer worldwide: a meta-analysis. *Br J Cancer*. 2003;88:63-73.
- Muñoz N, Castellsagué X, de González AB, et al. Chapter 1: HPV in the etiology of human cancer. *Vaccine*. 2006;24 Suppl 3:S3/1-10.
- Smith JS, Lindsay L, Hoots B, et al. Human papillomavirus type distribution in invasive cervical cancer and high-grade cervical lesions: a meta-analysis update. *Int J Cancer*. 2007;121:621-32.
- Khan MJ, Castle PE, Lorincz AT, et al. The elevated 10-year risk of cervical precancer and cancer in women with human papillomavirus (HPV) type 16 or 18 and the possible utility of type-specific HPV testing in clinical practice. *J Natl Cancer Inst*. 2005;97:1072-9.
- Davies P, Arbyn M, Dillner J, et al. A report on the current status of European research on the use of human papillomavirus testing for primary cervical cancer screening. *Int J Cancer*. 2006;118:791-6.
- Cuzick J, Clavel C, Petry KU, et al. Overview of the European and North American studies on HPV testing in primary cervical cancer screening. *Int J Cancer*. 2006;119:1095-101.
- Mayrand MH, Duarte-Franco E, Rodrigues I, et al. Human papillomavirus DNA versus Papanicolaou screening tests for cervical cancer. *N Engl J Med*. 2007;357:1579-88.
- Goldie SJ, Gaffikin L, Goldhaber-Fiebert JD, et al. Cost-effectiveness of cervical-cancer screening in five developing countries. *N Engl J Med*. 2005;353:2158-68.
- Kim JJ, Wright TC, Goldie SJ. Cost-effectiveness of human papillomavirus DNA testing in the United Kingdom, The Netherlands, France, and Italy. *J Natl Cancer Inst*. 2005;97:888-95.
- Goldie SJ, Kim JJ, Wright TC. Cost-effectiveness of human papillomavirus DNA testing for cervical cancer screening in women aged 30 years or more. *Obstet Gynecol*. 2004;103:619-31.
- Cuschieri KS, Cubie HA. The role of human papillomavirus testing in cervical screening. *J Clin Virol*. 2005;32 Suppl 1:S34-42.

22. Franco EL, Cuzick J. Cervical cancer screening following prophylactic human papillomavirus vaccination. *Vaccine*. 2008;26 Suppl 1:A16-23.
23. Stanley M, Villa LL. Monitoring HPV vaccination. *Vaccine*. 2008;26 Suppl 1:A24-7.
24. CLSI. Clinical Laboratory Waste Management; Approved Guideline - Third Edition. CLSI document GP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2011.
25. US Environmental Protection Agency. EPA Guide for Infectious Waste Management Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1-5-5, R1-R3, A1-A24.
26. Huang S, Erickson B, Tang N, *et al*. Clinical performance of Abbott RealTime High Risk HPV test for detection of high-grade cervical intraepithelial neoplasia in women with abnormal cytology. *J Clin Virol*. 2009;45(suppl 1):S19-23.
27. Halfon P, Benmoura D, Agostini A, *et al*. Evaluation of the clinical performance of the Abbott RealTime High-Risk HPV for carcinogenic HPV detection. *J Clin Virol*. 2010;48(4):246-50.
28. Szarewski A, Mesher D, Cadman L, *et al*. Comparison of seven tests for high-grade cervical intraepithelial neoplasia in women with abnormal smears: the Predictors 2 study. *J Clin Microbiol*. 2012;50(6):1867-73.
29. Jentschke M, Soergel P, Lange V, *et al*. Evaluation of a new multiplex real-time polymerase chain reaction assay for the detection of human papillomavirus infections in a referral population. *Int J Gynecol Cancer*. 2012;22(6):1050-6.
30. Mesher D, Szarewski A, Cadman L, *et al*. Comparison of human papillomavirus testing strategies for triage of women referred with low-grade cytological abnormalities. *Eur J Cancer*. 2013;49(9):2179-86. Supplementary table S1.
31. IARC Handbooks of Cancer Prevention. Volume 10. Cervix Cancer Screening. Lyon, France: International Agency for Research on Cancer; 2005.
32. Poljak M, Ostrbenk A, Seme K, Učakar V, Hillemanns P, Bokal EV, Jancar N, Klavs I. Comparison of clinical and analytical performance of the Abbott Realtime High Risk HPV test to the performance of hybrid capture 2 in population-based cervical cancer screening. *J Clin Microbiol*. 2011;49(5):1721-9.
33. Hesselink AT, Meijer CJLM, Poljak M, Berkhof J, van Kemenade FJ, van der Salm ML, Bogaarts M, Snijders PJF, Heideman DAM. Clinical validation of the Abbott RealTime High Risk (HR) HPV assay according to the guidelines for human papillomavirus DNA test requirements for cervical screening. *J Clin Microbiol*. 2013;51(7):2409-10.
34. Cuzick J, Cadman L, Mesher D, Austin J, Ashdown-Barr L, Ho L, Terry G, Liddle S, Wright C, Lyons D, Szarewski A. Comparing the performance of six human papillomavirus tests in a screening population. *Br J Cancer*. 2013;108(4):908-13. Supplementary table A1.

Abbott *m*, *m2000*, *m2000rt*, *m2000sp*, and Cervi-Collect are trademarks of Abbott Laboratories in various jurisdictions. ProClin, FAM, ROX, NED, VIC, Cy5, PreservCyt, PrepStain, SurePath, TriPath Imaging, hc2 High-Risk HPV DNA Test, Linear Array, Delfen, Gynecort, K-Y Jelly, Lubrin, MetroGel-Vaginal, Monistat, Norforms, Terazol, Vagi-gard, Vagisil, Yeast Gard and Zovirax are property of their respective owners.

www.abbottmolecular.com

 Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580



May 2020
© 2014, 2020 Abbott Laboratories

