

WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

Product: parascreen - Rapid test for Malaria Pan/Pf
WHO reference number: PQDx 0291-025-00

The parascreen - Rapid test for Malaria Pan/Pf with product codes 503030025(1T), 503030010, 503030025, 503030050, and 503030100, manufactured by Zephyr Biomedicals – A Division of Tulip Diagnostics (P) Ltd., Rest-of-World version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 03 June 2020.

Summary of WHO's prequalification assessment for the parascreen - Rapid test for Malaria Pan/Pf

	Date	Outcome
Prequalification listing	3 June 2020	listed
Dossier assessment	21 January 2020	MR
Product performance evaluation	2018	MR

MR: Meets Requirements

Report amendments or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	<ol style="list-style-type: none"> 1. Change of storage conditions of the product from 4 - 30 °C to 4 - 40 °C. 2. Change in use of specimen applicator from sample loops to Inverted cups. 3. Change in use of lancets from stainless steel sterile blood lancets to blood lancet plastic (twist off). 4. Updates on the status of commitments to prequalification. 	11 July 2022
3.0	Addition of a pack size, with single tests in a pack of 25 tests (product code 503030025(1T)).	3 October 2023
4.0	1. Removal of the CE mark from the WHO PQ products.	13 April 2026

	2. Additional approved suppliers of sterile blood lancets for the WHO PQ products (PQC-IVD-2025-0008).	
--	--	--

Intended use

According to the claim of intended use from Zephyr Biomedicals – A Division of Tulip Diagnostics (P) Ltd., “*parascreen [- Rapid test for Malaria Pan/Pf] is a rapid, qualitative, two site sandwich immunoassay utilizing capillary and venous whole blood specimens of symptomatic patients for the detection of P.falciparum specific histidine rich protein-2 (Pf. HRP-2) and Plasmodium Lactate Dehydrogenase (pLDH) antigens produced by P.falciparum, P.vivax, P.ovale, and P.malariae species. It is used in the diagnosis of malaria for differentiation of P.falciparum and other malaria species. parascreen [- Rapid test for Malaria Pan/Pf] is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating at point-of-care in resource-limited settings. This product is not intended for self-testing, and it is not for blood donor screening. The test is not automated; it needs to be performed and interpreted manually by the user*”.

Test kit contents

Component	25 tests (product code 503030025 (1T))	10 tests (product code 503030010)	25 tests (product code 503030025)	50 tests (product code 503030050)	100 tests (product code 503030100)
Pouch-sealed test with desiccant and specimen transfer device	25 x 1T (25 single kit test/kit)	10	25	50	100
Assay buffer bottle(s)	25 x 1.0 mL	1 bottle (total volume 3.0 mL)	1 bottle (total volume 4.0 mL)	2 bottles (total volume 8.0 mL)	4 bottles (total volume 16.0 mL)
Alcohol swabs	25	10	25	50	100
Sterile lancets	25	10	25	50	100
Instructions for use	25	1	1	1	1

Pictorial instructions for use	25	1	1	1	1
---------------------------------------	----	---	---	---	---

Items required but not provided

Item	Description
Consumables: Disposable micropipette tips Venipuncture blood collection kit Additional alcohol swabs Additional sterile lancets	NA
Durables: Biohazard waste container	NA
Equipment: Calibrated micropipette Timer	Micropipette should be capable of delivering 5µL of specimen.

Storage

The test kit must be stored between 4 and 40 °C.

Shelf-life upon manufacture¹

24 months.

Dossier assessment

Zephyr Biomedicals – A Division of Tulip Diagnostics (P) Ltd. submitted a product dossier for parascreen - Rapid test for Malaria Pan/Pf as per the “*Instructions for compilation of a product dossier*” (PQDx_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

Based on the product dossier screening and assessment findings, the product dossier for parascreen - Rapid test for Malaria Pan/Pf meets WHO prequalification requirements.

Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pgweb/vitro-diagnostics/who-public-inspection-reports>

Product performance evaluation

parascreen - Rapid test for Malaria Pan/Pf was evaluated in the eighth² round of WHO product testing of RDTs for malaria antigen detection, which was completed in 2018.

parascreen - Rapid test for Malaria Pan/Pf was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild-type parasite panel, *P. vivax* wild-type parasite panel, and a *P. falciparum* and *P. vivax* negative panel.

Performance Characteristics		
	<i>P. falciparum</i>	<i>P. vivax</i>
Panel detection score at 200 parasites/ μ L (N=100)	91	94.3
False positive results % (N= 208)	0.5	
Invalid rate % (N= 1210)	0	
Inter-reader variability %*	Not applicable	
Lowest concentration of HRP2/pLDH was detected using	Not applicable	

² <https://www.who.int/malaria/publications/atoz/9789241514965/en/>

the 1 st WHO International standard for Pf antigens (NIBSC code: 16/376)*	
--	--

** Not applicable for assays evaluated in WHO product testing of RDTs for malaria antigen detection*

Operational characteristics and ease of use

This assay does not require laboratory equipment and can be performed in laboratories with limited facilities or non-laboratory settings.

Key operational characteristics	
Number of steps*	2 steps in total
Time to result	20 minutes
Endpoint stability (interval)	10 minutes (the test can be read between 20 and 30 minutes after the addition of diluent)
Internal QC	The test has an internal control line. The presence of the control line indicates that migration of liquid has occurred; however, it does not guarantee that the correct specimen type or volume was added or that the test procedure was followed correctly.

** Definition: each action required to obtain a result (excluding specimen collection, device preparation – opening the pouch), e.g., for RDTs: add specimen, add buffer (2 steps).*

Labelling review

The labelling submitted for the parascreen - Rapid test for Malaria Pan/Pf assay was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.


The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

Controlled Labelling References

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
Outer box artwork	10 Test Carton	VER-03	21/02/2022	0320/VER-03
	25 Test Carton	VER-03	14/09/2020	0320/VER-03
	50 Test Carton	VER-03	10/04/2023	0320/VER-03
	100 Test Carton	VER-03	10/04/2023	0320/VER-03
	25x1 Test Carton	VER-02	09/08/2023	0723/VER-02
	1 Test Carton	VER-02	07/03/2026	0326/VER-02
Pouch / Device label	Pouch	VER-07	07/03/2026	0326/VER-07
Reagent bottle labels	Buffer label	VER-07	07/03/2026	VER-07
	Buffer label – 1ml	VER-02	07/03/2026	VER-02
Lancets	Blood lancets (Shangdong)	REV.20250621	N/A	N/A
	Blood lancets (Promised)	REV 01.	N/A	N/A
Instructions for Use (IFU)	IFU	VER-09	07/03/2026	PI/50303/0326/VER-09
Pictorial IFU	Pictorial Representation	VER-08	07/03/2026	PR/50303/0326/VER-08
Carton Label	10 Test Kit label	VER-07	07/03/2026	0326/VER-07
	25 Test Kit label	VER-08	07/03/2026	0326/VER-08
	50 Test Kit label	VER-07	07/03/2026	0326/VER-07
	100 Test Kit label	VER-07	07/03/2026	0326/VER-07
	25x1 Tests Kit label	VER-02	07/03/2026	0326/VER-02

Labels

1 Tests Carton
Size: 70 mm x 22 mm X 130 mm



Lot No. :
Mfg. Date :
Expiry Date :
M.R.P. (₹) :
(Inclusive of all taxes)



IN VITRO DIAGNOSTIC MEDICAL DEVICE
NOT FOR MEDICINAL USE
REFER PACKAGE INSERT
FOR INSTRUCTIONS
BEFORE USE



Manufactured by:
Zephyr Biomedicals
A Division of Tulip Diagnostics (P) Ltd.
M 46-47, Phase III B, Verna Industrial Estate,
Verna, Goa - 403 722, INDIA.
Regd. Office: Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O.,
Goa - 403 202, INDIA.

An ISO 13485
Certified Company

0326/VER-02



25 x 1 Test Carton
Size: 310 mm x 135 mm X 145 mm



Manufactured by:
Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.
M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.
Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O., Goa - 403 202, INDIA.
Email address: sales@tulipgroup.com
Tel. : (0832) 2458546, (0832) 2458547




This side up

Lot No. :
 Mfg. Date :
 Expiry Date :
M.R.P. (₹) :
(Inclusive of all taxes)


An ISO 13485
Certified Company




IVD IN VITRO DIAGNOSTIC MEDICAL DEVICE
NOT FOR MEDICINAL USE

REFER PACKAGE INSERT FOR
INSTRUCTIONS BEFORE USE



25 Tests Carton
Size: 200 mm x 132 mm X 70 mm



Manufactured by:
Zephyr Biomedicals
A Division of Tulip Diagnostics (P) Ltd.
M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.
Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O., Goa - 403 202, INDIA.
Email address: sales@tulipgroup.com
Tel. : (0832) 2458546, (0832) 2458547


 **Lot No.** :
 **Mfg. Date** :
 **Expiry Date** :
M.R.P. (₹) :
(Inclusive of all taxes)



This side up

An ISO 13485
Certified Company

 **IN VITRO DIAGNOSTIC MEDICAL DEVICE**
NOT FOR MEDICINAL USE




 **REFER PACKAGE INSERT FOR**
INSTRUCTIONS BEFORE USE




100 Tests Carton
Size: 240 mm x 128 mm X 190 mm




This side up


 **Lot No.** :
 **Mfg. Date** :
 **Expiry Date** :
M.R.P. (₹) :
(Inclusive of all taxes)


Manufactured by:
Zephyr Biomedicals
A Division of Tulip Diagnostics (P) Ltd.
M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.
Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O., Goa - 403 202, INDIA.
Email address: sales@tulipgroup.com
Tel. : (0832) 2458546, (0832) 2458547

An ISO 13485
Certified Company




 **IN VITRO DIAGNOSTIC MEDICAL DEVICE
NOT FOR MEDICINAL USE**

 **REFER PACKAGE INSERT FOR
INSTRUCTIONS BEFORE USE**




10 Tests Carton (Device)
Size: 132 mm x 90 mm X 70 mm



Manufactured by:
Zephyr Biomedicals
A Division of Tulip Diagnostics (P) Ltd.
M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.
Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O., Goa - 403 202, INDIA.
Email address: sales@tulipgroup.com
Tel. : (0832) 2458546, (0832) 2458547


 **Lot No.** :
 **Mfg. Date** :
 **Expiry Date** :
M.R.P. (₹) :
(Inclusive of all taxes)



An ISO 13485
Certified Company

 **IN VITRO DIAGNOSTIC MEDICAL DEVICE
NOT FOR MEDICINAL USE**




 **REFER PACKAGE INSERT FOR
INSTRUCTIONS BEFORE USE**



50 Tests Carton
Size: 200 mm x 132 mm X 132 mm



This side up

 **Lot No.** :
 **Mfg. Date** :
 **Expiry Date** :
M.R.P. (₹) :
(Inclusive of all taxes)



Manufactured by:

Zephyr Biomedicals


A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O., Goa - 403 202, INDIA.

Email address: sales@tulipgroup.com
Tel. : (0832) 2458546, (0832) 2458547

An ISO 13485
Certified Company

 **IN VITRO DIAGNOSTIC MEDICAL DEVICE
NOT FOR MEDICINAL USE**

 **REFER PACKAGE INSERT FOR
INSTRUCTIONS BEFORE USE**

Size : 100 x 78 mm

Contents:

- 10 x Test devices.
- 1 x Clearing buffer bottle.
- 10 x Disposable specimen applicators.
- 1 x Instructions for use.
- 1 x Pictorial Instructions for use.
- 10 x Alcohol swabs.
- 10 x Sterile lancets.

Mfg. Lic. No.: MFG/IVD/2020/000043

REF Cat No.: 503030010



4°C -40°C Store at 4°C to 40°C
Do Not Freeze

(DO NOT ACCEPT IF THE SEAL IS BROKEN)

10 TESTS



parascreen®

DEVICE

RAPID TEST FOR MALARIA
Pan / Pf

0326/VER-07

Size : 100 x 78 mm

Contents:

Mfg. Lic. No.: MFG/IVD/2020/000043



Cat No.: 503030025

- 25 x Test devices.
- 1 x Clearing buffer bottle.
- 25 x Disposable specimen applicators.
- 1 x Instructions for use.
- 1 x Pictorial Instructions for use.
- 25 x Alcohol swabs.
- 25 x Sterile lancets.

4°C Store at 4°C to 40°C
Do Not Freeze



(DO NOT ACCEPT IF THE SEAL IS BROKEN)

25 TESTS







parascreen®

DEVICE

RAPID TEST FOR MALARIA
Pan / Pf

0326/VER.08

70 x 18 mm

IVD	For in vitro diagnostic use only Not for medicinal use Read package insert for instructions before use	parascreen [®] CLEARING BUFFER	4°C-40°C Store at 4°C to 40°C Do Not Freeze
			
Manufactured by Zephyr Biomedicals A Division of Taty Diagnostics (P) Ltd. 18 66-67, Phase II B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA. Regd. Office: Changan, Tulp Block, D, Antonio De Rego Bldg, Alto Bertemuz, Bertemuz, Coimbra P.C., 304-403 202, PORTUGAL.		BUF	Mfg. Lic. No.: MFG/IVD/2020/000043
			Lot No. : Mfg. Date : Expiry Date : Vol. :
			
			VER-07

Size : 125 x 90 mm

Contents:

- 50 x Test devices.
- 2 x Clearing buffer bottles.
- 50 x Disposable specimen applicators.
- 1 x Instructions for use.
- 1 x Pictorial Instructions for use.
- 50 x Alcohol swabs.
- 50 x Sterile lancets.

Mfg. Lic. No.: MFG/IVD/2020/000043

REF Cat No.: 503030050



4°C - 40°C Store at 4°C to 40°C
Do Not Freeze



(DO NOT ACCEPT IF THE SEAL IS BROKEN)



50 TESTS



parascreen®

DEVICE

RAPID TEST FOR MALARIA
Pan / Pf

0326/VER-07

Size : 150 x 120 mm


Contents:

- 100 x Test devices.
- 4 x Clearing buffer bottles.
- 100 x Disposable specimen applicators.
- 1 x Instructions for use.
- 1 x Pictorial Instructions for use.
- 100 x Alcohol swabs.
- 100 x Sterile lancets.

Mfg. Lic. No.: MFG/IVD/2020/000043

REF Cat No.: 503030100




 Store at 4°C to 40°C
 Do Not Freeze



(DO NOT ACCEPT IF THE SEAL IS BROKEN)



100 TESTS



parascreen®

DEVICE

RAPID TEST FOR MALARIA
Pan / Pf

0326/VER-07

Size : 70 x 130 mm

parascreen®

DEVICE

RAPID TEST FOR MALARIA
Pan / Pf

CONTAINS ONE MEMBRANE TEST ASSEMBLY
WITH DESICCANT

IVD IN VITRO DIAGNOSTIC MEDICAL DEVICE
NOT FOR MEDICINAL USE

 FOR SINGLE USE ONLY. DO NOT REUSE

 STORE BETWEEN 4°C TO 40°C

 REFER INSTRUCTIONS FOR USE



Manufactured by:
Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitarjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex P.O., Goa - 403 202, INDIA.

Email address: sales@tulipgroup.com Tel. : (0832) 2458546, (0832) 2458547

Mfg. Lic. No. : MFG/IVD/2020/000043




 :
 Lot No. :
 Mfg. Date :
 Expiry Date :



0326/VER-07


PANTONE Reflex Blue C
Size:50*42mm


Blood Lancets

Model / Specification : I / 28G

	LOT NO.	:	XXXXXXXX
	MFG. DATE	:	XXXX-XX-XX
	EXP. DATE	:	XXXX-XX-XX
	QTY	:	25pcs

  0197


 **Linkfar Healthcare GmbH**
Niederheinstraße 71, 40474 Düsseldorf, Germany
TEL: +49-21138332888

 **Shandong Lianfa Medical Plastic Products Co., Ltd.**
No. 3, Shuangshan Sanjiao Road, Zhangqiu, Jinan City, 250000, Shandong P. R. China

Exclusively Imported & Marketed in India By: RadSun Biotech Corporation
Pn D-31, Kh No. 391 Block-D Indrapuri Loni Ghaziabad, Uttar Pradesh-201102

REV:20250621

Size 200 mm x 75 mm



Instructions for use and pictorial package insert³

³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

RAPID TEST FOR MALARIA

Pan/Pf

DEVICE

INTENDED USE

parascreen® is a rapid, qualitative, two site sandwich immunoassay utilizing capillary and venous whole blood specimens of symptomatic patients for the detection of *P.falciparum* specific histidine rich protein-2 (Pf. HRP-2) and Plasmodium Lactate Dehydrogenase (pLDH) antigens produced by *P.falciparum*, *P.vivax*, *P.ovale* and *P.malariae* species. It is used in the diagnosis of malaria for differentiation of *P.falciparum* and other malaria species.

parascreen® is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating at point-of-care in resource-limited settings. This product is not intended for self-testing and it is not for blood donor screening. The test is not automated; it needs to be performed and interpreted manually by the user.

SUMMARY

Four species of the Plasmodium parasites are responsible for malaria infections in human viz. *P. falciparum*, *P.vivax*, *P.ovale* and *P.malariae*. Of these, *P. falciparum* and *P.vivax* are the most prevalent. Early detection and differentiation of malaria is of utmost importance due to incidence of cerebral malaria and drug resistance associated with falciparum malaria and due to the morbidity associated with the other malarial forms.

parascreen® detects the presence of Pan malaria specific pLDH released from parasitised blood cells, for the detection of all malarial parasites. Whereas, for the detection of *P. falciparum* malaria, **parascreen®** utilises the detection of *P. falciparum* specific histidine rich protein-2 (Pf. HRP-2) which is a water soluble protein that is released from parasitised red blood cells of infected individuals. In the absence of *P.falciparum* specific Pf. HRP-2, the presence of Pan malaria specific band points to the presence of other malarial species such as *P.vivax*, *P.ovale* or *P.malariae*. Speciation is done and results inferred in the context of prevalence rates of the malarial species prevalent in the particular region.

PRINCIPLE

parascreen® utilizes the principle of agglutination of antibodies/ antisera with respective antigen in immuno-chromatography format along with use of nano gold particles as agglutination revealing agent. As the test specimen flows through the membrane assembly of the device after addition of the clearing buffer, the colored colloidal gold conjugates of the Agglutinating sera for HRP-2 / Agglutinating sera for Pan malaria specific pLDH - colloidal gold conjugate complexes the proteins in the lysed specimen. This complex moves further on the membrane to the test region where it is immobilised by the Agglutinating sera for HRP-2 / Agglutinating sera for Pan malaria specific pLDH coated on the membrane leading to formation of pink-purple colored band/s which confirms a positive test result. Absence of this colored band/s in the test region indicates a negative test result.

The unreacted conjugate along with the rabbit globulin-colloidal gold conjugate and unbound complex if any, move further on the membrane and are subsequently immobilised by Agglutinating sera for rabbit globulin coated on the membrane at the control region, forming a pink-purple band. The control band formation is based on the 'Rabbit / Agglutinating sera for Rabbit globulin' system. Since it is completely independent of the analyte detection system, it facilitates formation of consistent control band signal independent of the analyte concentration. This control band serves to validate the test performance.

REAGENTS AND MATERIALS SUPPLIED

parascreen® kit contains:

- A. Individual pouches, each containing:
 1. **DEVICE** Test Device: Membrane assembly pre-dispensed with Agglutinating sera for HRP-2 - colloidal gold conjugate, Agglutinating sera for Pan malaria specific pLDH - colloidal gold conjugate, rabbit globulin-colloidal gold conjugate, Agglutinating sera for HRP-2, Agglutinating sera for Pan malaria specific pLDH and Agglutinating sera for rabbit globulin at the respective regions.
 2. Desiccant pouch.
- B. **PIPETTE** Disposable Plastic Specimen Applicator.
- C. **BUF** Clearing Buffer in a dropper bottle.
- D. Instructions for use.
- E. Pictorial instructions for use.
- F. Alcohol swabs – 70% Isopropyl alcohol
- G. Sterile lancets.

Product codes	REF	503030025(1T)	503030010	503030025	503030050	503030100
Pouch Sealed tests	▽	25 x 1T (25 single kit test/kit)	10	25	50	100
Disposable Plastic Specimen Applicator		25	10	25	50	100
Clearing buffer bottles		25 x 1.0ml	01 x 3.0ml	01 x 4.0ml	02 x 4.0ml	04 x 4.0ml
Alcohol swabs		25	10	25	50	100
Sterile lancets		25	10	25	50	100
Instructions for use		25	01	01	01	01
Pictorial instructions for use		25	01	01	01	01

MATERIALS REQUIRED BUT NOT PROVIDED

Calibrated micropipette capable of delivering 5µl specimen accurately, disposable micropipette tips.

Permanent marker Pen/pencil, disposable gloves, timer.

Biosafety sharps container and Biohazard waste container (for potentially infectious waste).

Venipuncture blood collection kit (if whole blood is collected by venepuncture).

Additional alcohol swabs (if any included in the kit are found dry) and additional sterile lancets (if any included in the kit have the sterility seal broken).

STORAGE AND STABILITY

The sealed pouches in the test kit & the kit components may be stored between 4°C to 40°C till the duration of the shelf life as indicated on the pouch/ carton. DO NOT FREEZE. After first opening of the clearing buffer bottle, it can be stored between 4°C to 40°C for the remaining duration of its shelf life.

WARNINGS

Read the instructions carefully before performing the test.

For in vitro diagnostic use only. NOT FOR MEDICINAL USE. For professional use.

The test is for aiding in diagnosis of malaria infection and not for screening which requires confirmation.

Do not use beyond expiry date.

Do not use components from different lots of the product.

The device, specimen applicator, alcohol swab and blood lancet are for single use only.

Contact with the contents of desiccant pouch containing, among other substances, cobalt chloride (CAS# 7646-79-9) should be kept to a minimum. Inhalation / swallowing may cause harm.

Handle all specimens as potentially infectious.

Follow standard biosafety guidelines for handling and disposal of potentially infectious material.

Clearing buffer contains Sodium Azide(0.1%), avoid skin contact with this reagent. Azide may react with lead and copper in the plumbing system and form highly explosive metal oxides. Flush with large volumes of water to prevent azide build up in the plumbing.

SPECIMEN COLLECTION AND PREPARATION


For specimen collection, refer to pictorial instructions for use.

Fresh capillary/venous wholeblood from finger prick / puncture should be used as a test specimen. However, fresh anti-coagulated venous whole blood may also be used as a test specimen. Using standard blood collection practices, collect venous whole blood into the commercially available anti-coagulant tube such as EDTA or CPDA or Heparin or Oxalate or Tri-sodium Citrate. If immediate testing is not possible then the specimen may be stored at 2°C to 8°C for upto 72 hours before testing and should be brought to room temperature (20°C to 30°C) before use on the test. Clotted,

hemolysed or lipaemic whole blood specimens should not be used for performing the test.

TESTING PROCEDURE AND INTERPRETATION OF RESULTS

- Bring the **parascreen**® kit components to room temperature (20°C to 30°C) before testing.
- Open the cap to retrieve the device and desiccant pouch. Check the color of the desiccant. It should be blue, if it has turned colorless or pink, discard the device and use another device. **Once opened, the device must be used immediately.**
- Label the test device with patient identifier.
- Place the testing device on a flat horizontal surface.
- Tighten the cap of the clearing buffer bottle provided with the kit in the clockwise direction to pierce the buffer bottle nozzle.
- Specimen application**
 - 6.1 Venous whole blood:** Evenly mix the anti-coagulated whole blood by gentle swirling. Dip the specimen applicator provided in the kit into the whole blood. Ensuring that an applicator full of blood is retrieved, immediately blot the blood so collected in the specimen port 'A' (This delivers approximately 5µl of the whole blood specimen). Alternatively, 5µl of the anti-coagulated venous whole blood specimen may be delivered in the specimen port 'A' using a micropipette.
 - 6.2 Capillary whole blood:** Touch the specimen applicator provided in the kit into the whole blood on the finger prick. Ensuring that an applicator full of blood is retrieved, immediately blot the blood so collected in the specimen port 'A' (Care should be taken that whole blood specimen is not clotted and transfer to the specimen port is immediate). Alternatively, 5µl of the capillary finger-prick whole blood specimen may be delivered in the specimen port 'A' using a micropipette.
- Note:** Ensure that the whole blood from the specimen applicator has been completely taken up at the specimen port 'A'.
- Immediately dispense **two drops** of clearing buffer into buffer port 'B' holding the buffer bottle vertically and switch on the timer. To avoid contamination of clearing buffer bottle, do not touch the buffer port 'B' with the tip of clearing buffer bottle.
- Read the results at the end of **20 minutes** as follows :



NEGATIVE for malaria: Only one pink-purple band appears in the control window 'C'.



POSITIVE for P. falciparum: In addition to the control band, two pink-purple bands appear at regions 'Pf' and 'Pan' in the test window 'T'.

OR

In addition to the control band, a pink-purple band appear only at region 'Pf' in the test window 'T'.

Appearance of coloured bands of any intensity (faint to dark) at 'Pf' and/or 'Pan' should be considered as positive result for *P. falciparum*.



POSITIVE for Mixed infection (P. falciparum and P. vivax or P. malariae or P. ovale) : In addition to the control band, two pink-purple bands appear at regions 'Pf' and 'Pan' in the test window 'T'. Appearance of coloured bands of any intensity (faint to dark) at 'Pf' and 'Pan' should be considered as positive result for Mixed infection.



POSITIVE for Other species (non falciparum): In addition to the control band, one pink-purple band appears only at region 'Pan' in the test window 'T'. Appearance of a coloured band of any intensity (faint to dark) at 'Pan' should be considered as positive result *P. vivax* or *P. malariae* or *P. ovale* malaria.



INVALID RESULT: The test should be considered invalid if no bands appear on the device. The test should also be considered invalid if only test bands (Pan and/or Pf) appear and no control band appears. Repeat the test with a new device ensuring that the test procedure has been followed accurately.

CAUTION : Do not read results after 30 minutes as it may show erroneous results.

PERFORMANCE CHARACTERISTICS

A. Analytical Performance Study

A1. Potentially interfering exogenous and endogenous substances :

The following Potentially interfering substances have no impact on test results of **parascreen**® :

Type of Specimen		Sr. No.	Potential Interfering substances
Endogenous substance		1	Total Protein
		2	Bilirubin, conjugated
		3	Cholesterol
		4	Triglycerides
		5	Haemoglobin
Common Drugs	Antibiotic	1	Amoxicillin
		2	Ciprofloxacin
	Anti-inflammatory	1	Aspirin
		2	Ibuprofen
Exogenous Substance	Anti-Malaria Drugs	1	Chloroquine
		2	Doxycycline
		3	ACT
		4	Primaquine
		5	Mefloquine
		6	Sulfadoxine
		7	Pyrimethamine
	Anti-TB Drugs	1	Ethambutol
		2	Isoniazide
		3	Rifampin
	Anti-Retroviral Drugs	1	Lamivudine
		2	Efavirenz
		3	Emtricitabine
		4	Tenofovir
		5	Atazanavir

A2. Cross Reacting infections, disease and medical conditions:

The following 17 potential cross reacting infections/diseases/conditions did not affect the performance of **parascreen**®.

Potential Cross reacting infections/diseases/conditions			
1	<i>T. cruzi</i>	10	<i>Toxoplasma gondii</i>
2	Dengue virus	11	Influenza A/B
3	<i>Leishmania spp</i>	12	Yellow fever virus
4	<i>Brucella spp</i>	13	<i>Leptospira spp</i>
5	Measles virus (Rubeola virus)	14	<i>Treponema pallidum</i>
6	HAV	15	HAMA
7	HBV	16	ANA
8	HCV	17	Rheumatoid factor
9	HIV-1/HIV-2		

A3. Precision (Repeatability)

Within run, precision was determined using 10 replicates of 5 different venous whole blood specimens in 03 different lots of **parascreen**[®] which is summarized below:

*Quality control Panel	Accuracy (%)
Malaria Negative	100%
<i>P.falciparum</i> Positive (Moderate Positive)	100%
<i>P.falciparum</i> Positive (Weak Positive)	100%
<i>P.vivax</i> Positive (Moderate Positive)	100%
<i>P.vivax</i> Positive (Weak Positive)	100%

A4. Precision (Reproducibility)

Between run, precision was determined using 5 different blinded venous whole blood specimens in 3 different lots of **parascreen**[®] X 3 different operators X 3 different sites X 5 different days which is summarized below:

*Quality control Panel	Accuracy (%)			
	Between Day	Between Operator	Between lot	Between site
Malaria Negative	100%	100%	100%	100%
<i>P.falciparum</i> Positive (Moderate Positive)	100%	100%	100%	100%
<i>P.falciparum</i> Positive (Weak Positive)	100%	97.7%	100%	100%
<i>P.vivax</i> Positive (Moderate Positive)	100%	100%	100%	100%
<i>P.vivax</i> Positive (Weak Positive)	100%	100%	100%	100%

*Quality control panel specimens have been confirmed by microscopy as malaria negative and malaria positive. Malaria positive specimens were classified as moderate or weak positive based on respective parasite counts as determined by microscopy.

A5. Analytical Sensitivity

The sensitivity of **parascreen**[®] for *P.falciparum* is 100 parasites/μl and for *P.vivax* is 200 parasite/μl based on microscopy results.

B. Clinical Performance study: Diagnostic Specificity and Diagnostic Sensitivity

B1. In an in-house study, a panel of 251 venous whole blood specimens whose results were earlier confirmed with microscopy were tested with **parascreen**[®]. The results obtained are as follows:

Specimens	Total no. of specimens tested	parascreen [®]		Sensitivity (95% CI)	Specificity (95% CI)
		Positive	Negative		
<i>P.falciparum</i>	16	16	0	100% (79.41% to 100.00%)	-
<i>P.vivax</i>	25	25	0	100% (86.28% to 100.00%)	-
Malaria Negative	210	0	210	-	100% (98.26% to 100.00%)

B2. External evaluation studies:

Table 1

Study Site	Total Number of Malaria Negative specimens Tested	Specimen Type		Number of specimens Negative by Microscopy	Number of specimens Negative in parascreen [®]	Number of specimens falsely Positive in parascreen [®]
		Population type	Mode of Collection			
Jharkhand, India	985	Hospitalized Patients	Finger prick/ venous phlebotomy	985	985	0
Maharashtra, India.	1000	Blood Donors	Venous whole blood	1000	1000	0
Goa, India.	39	Symptomatic/ Asymptomatic Individuals	Capillary Whole Blood	39	39	0
			Venous Whole Blood	39	39	0
Odisha, India	545	Pregnant Women	Venous Whole Blood	545	545	0
Odisha, India	497	Neonates	Heel Prick	497	497	0
Based on above data:						
Total Nos. tested		Overall Specificity		95% Confidence Interval		
3105		100%		99.88% to 100.00%		

Table 2

Study Site	Total Number of Malaria Positive specimens Tested	Population type	Specimen Type		Number of specimens Positive by Microscopy	Number of specimens Positive in parascreen [®]	Number of specimens falsely Negative in parascreen [®]
			Mode of Collection	Species Type			
India	403	Hospitalized Patients	Finger prick/ venous phlebotomy	<i>P.falciparum</i>	403	403	0
	312			<i>P.vivax</i>	312	312	0
	26	Symptomatic/ Asymptomatic Individuals	Capillary and Venous Whole Blood	<i>P.falciparum</i>	26	26	0
	29			<i>P.vivax</i>	29	29	0
	06			<i>P.falciparum</i> + <i>P.vivax</i>	06	06	0
	01	Pregnant Women	Venous Whole Blood	<i>P.falciparum</i>	01	01	0
	04			<i>P.vivax</i>	04	04	0
	01	Neonates	Heel Prick	<i>P.falciparum</i>	01	01	0
02	<i>P.vivax</i>			02	02	0	
Malawi	25	Hospitalized Patients	Venous whole blood	<i>P.ovale</i>	25	25	0
	25			<i>P.malariae</i>	25	25	0
Malaysia	90	Hospitalized Patients	Venous whole blood	<i>P. knowlesi</i>	90	74	16
Based on above data:							
Plasmodium species		Total Nos. tested		Overall Sensitivity		95% Confidence Interval	
<i>P.falciparum</i>		431		100%		99.15% to 100.00%	
<i>P.vivax</i>		347		100%		98.94% to 100.00%	
<i>P.ovale</i>		25		100%		84.56% to 100.00%	
<i>P.malariae</i>		25		100%		84.56% to 100.00%	
<i>P. knowlesi</i>		90		82.22%		72.74% to 89.48%	
<i>P.falciparum</i> + <i>P.vivax</i>		06		100%		54.07% to 100.00%	

LIMITATIONS OF THE TEST

1. As with all diagnostic tests, the test result must always be correlated with clinical findings.
2. The results of test are to be interpreted within the epidemiological, clinical and therapeutic context.
3. Any modification to the above procedure and / or use of other reagents will invalidate the test procedure.
4. Hook effect may be observed at parasite density $\geq 3,00,000$ parasite/ μ l. In such cases, repeat the test by using different dilutions of same specimen. Other clinical data (e.g symptoms, travel history, risky factors) should be used in conjunction with the test results.
5. Interference due to presence of heterophile antibodies in patient's specimen can lead to erroneous analyte detection in immunoassay, has been reported in various studies. **parascreen**[®] uses HETEROPHILIC BLOCKING REAGENT (HBR) to inhibit majority of these interferences.
6. Potential cross-reacting diseases such as HAT, Tick-borne Encephalitis and those caused by *Schistosoma spp* have not been tested in this product, and their associated interference in **parascreen**[®] is not known.
7. Due to limited evidence, the manufacturer does not claim the limit of detection of **parascreen**[®] for *P.knowlesi* though it is found to detect these Plasmodium species as low as 226 -300 parasites/ μ l.
8. In case of mixed infection (*P.falciparum*, with other malarial species), both, 'Pf' and 'Pan' malaria bands will be positive. Hence, differentiation of infection due to *P.vivax*, *P.ovale* or *P.malariae* cannot be done.
9. In *P.falciparum* malaria infection, HRP-2 is not secreted in gametogony stage. Hence, in "Carriers", the HRP-2 band may be absent.











WARRANTY

This product is designed to perform as described on the label and Instructions for use. The manufacturer disclaims any implied warranty of use and sale for any other purpose. In the event of performance changes or product malfunction, please contact manufacturer.

BIBLIOGRAPHY

1. Howard, R.J., et al., 1986: Secretion of a Malarial Histidine-rich Protein (Pf. HRP II) from *Plasmodium falciparum*-infected Erythrocytes. J. Cell Biol., 103, 1269-1277.
2. Parra, M.E., et al., 1991: Identification of *Plasmodium falciparum* Histidine-Rich Protein 2 in the Plasma of Humans with Malaria. J. Clin. Microbiol., 29, 1629-1634.
3. Rodriguez-Del Valle, M., et al., 1991: Detection of Antigens and Antibodies in the Urine of Humans with *Plasmodium falciparum* Malaria. J. Clin. Microbiol., 29, 1236-1242.
4. Piper, R. C., et al., (1999) Immuno-capture diagnostic assays for malaria utilizing *Plasmodium* Lactate Dehydrogenase (pLDH) Am. J. Trop. Med. Hyg. 60(1) 109-118.
5. Hunte-Cooke A., et al., (1999) Comparison of a Parasite Lactate Dehydrogenase-based Immunochromatographic Antigen Detection assay (OptiMAL[®]) with Microscopy for the Detection of Malaria Parasites in Human Blood Samples. Am J. Trop Med 60(2). 173-176.
6. Quintana M., et al., (1998) Malaria diagnosis by dipstick assay in a Honduran Population with coendemic *Plasmodium falciparum* and *Plasmodium vivax*. Am. J. Trop. Med. Hyg. 59(6) 868-871.
7. Palmer, C. J., (1998) Evaluation of OptiMal test for rapid diagnosis of *Plasmodium vivax* and *Plasmodium falciparum*. J. Clin Microbiol. 36(1) 203-206.
8. Moody A., et al., (2000) Performance of the OptiMAL[®] malaria antigen capture dipstick for malaria diagnosis and treatment monitoring. British Journal of Hematology, 109, 1-5.
9. Data on file: Zephyr Biomedicals.

SYMBOL KEYS

 Temperature Limitation	 Manufacturer	DEVICE Device	 Harmful if swallowed. Do not breathe vapour, if swallowed, seek medical advice immediately and show this container or label. Avoid release to the environment. Refer to special instructions.
 Use by	 Consult Instructions for use	PIPETTE Disposable Plastic Specimen Applicator	
 Date of Manufacture	REF Catalogue Number	BUF Clearing Buffer	
LOT Batch Number / Lot Number	IVD <i>In vitro</i> Diagnostic Medical Device	 This side up	
 Contains sufficient for <n> tests	 Do not reuse	 Do not use if package is damaged	



Manufactured by:

Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.


M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.


Email address: sales@tulipgroup.com

Tel. : (0832) 2458546, (0832) 2458547


Interpretation of Results




Band on 'C' Area
NEGATIVE
for Malaria




Band on 'C' & 'T'-Pan Area
POSITIVE
for Non *Pf.* Malaria




Band on 'C' & 'T'-Pf Area
POSITIVE
for *Pf.* Malaria




Band on 'C' & 'T'-Pan, Pf Area
POSITIVE
for *Pf.* or Mixed infection
(*Pf.* or *Pv.* or *Pm.* or *Po.*)




No band on 'C' & 'T' Area
INVALID
Repeat the test



No band on 'C'.
Band on 'T'-Pan Area
INVALID
Repeat the test




No band on 'C'.
Band on 'T'-Pf Area
INVALID
Repeat the test




No band on 'C'.
Band on 'T'-Pan & Pf Area
INVALID
Repeat the test

Note : Coloured bands of any intensity (faint to dark) at 'Pan' and/or 'Pf' in addition to Control band should be considered as positive result.


Test Procedure using Capillary whole blood specimen




Open pouch at cut mark and remove the device & desiccant pouch. Once opened, the device must be used immediately.
Note: Before opening the pouch, check the date of expiration on the back side of the pouch. Do not use beyond expiry date.




For blood collection, hold the specimen applicator vertically as shown above Fig A. With the pricked finger facing downwards towards the specimen applicator, allow the blood drop to touch the cup side of specimen applicator, Fig.B. Ensure that the cup full of blood is retrieved Fig.C.




Component details are shown above. Check the colour of the desiccant. It should be blue, if it has turned colourless or pink, discard the device & use a new device.




Blot the collected blood (5µl) in the specimen port 'A' by touching applicator vertically straight onto the specimen pad. Ensure that the blood from the applicator has been completely taken up at the specimen pad.



Clean finger to be pricked with an alcohol swab. Allow to dry. "Avoid finger with callus".



Hold the clearing buffer bottle vertically straight over the buffer port 'B' without having contact with the pad to avoid contamination & add exactly 2 drops of buffer onto the buffer port. Start the stop watch.



Take a sterile lancet. Open the lancet by twisting the cap. Squeeze the finger tip then prick the lateral side of the tip with sterile lancet. Safely dispose off the used lancet as biohazardous sharps waste.



Read the test Results at the end of 20 minutes. Do not read test results after 30 minutes.

For further information, contact :

Zephyr Biomedicals


A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.


Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

E-mail address: sales@tulipgroup.com Tel. : (0832) 2458546, (0832) 2458547


Test Procedure for Venous Whole Blood Specimen




Open pouch at cut mark and remove the device & desiccant pouch. Once opened, the device must be used immediately.
Note: Before opening the pouch, check the date of expiration on the back side of the pouch. Do not use beyond expiry date.




Component details are shown above. Check the colour of the desiccant. It should be blue, if it has turned colourless or pink, discard the device & use a new device.




Assemble equipment & include needle & syringe.



Apply a tourniquet, about 4 to 5 finger widths above the selected venepuncture site.




Ask the patient to form a fist.




Put on well fitting non-sterile gloves.




Disinfect the site using 70% Isopropyl Alcohol and let it dry completely.




Enter the syringe in vein swiftly at a 30° angle.




Release the tourniquet once sufficient blood has been collected.




Withdraw the needle gently and then give the patient a clean cotton.




Label the tube with respective patient's name.




Transfer blood to the labeled tube.




Place the cap in the tube.




Invert specimen tube back and forth 5-10 times to mix the additives with the whole blood.



Dip the cup side of the specimen applicator in the tube. Ensure that an applicator full of blood is retrieved as shown in the above figure.



Blot the collected blood (5µl) in the specimen port 'A' by touching applicator vertically straight onto specimen pad. Ensure that the blood from the applicator has been completely taken up at the specimen pad.



Hold the clearing buffer bottle vertically straight over the buffer port 'B' without having contact with the pad to avoid contamination & add exactly 2 drops of buffer onto the buffer port. Start the stop watch.



Read the test result at the end of 20 minutes. Do not read test result after 30 minutes.