## WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

**Product: STANDARD Q HCV Ab Test**  
**WHO reference number: PQDx 0360-117-00**

STANDARD Q HCV Ab Test with product code 09HCV10D, manufactured by SD Biosensor, Inc., Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 5 March 2020.

### Summary of WHO Prequalification Assessment for the STANDARD Q HCV Ab Test

	Date	Outcome
<b>Prequalification listing</b>	5 March 2020	listed
<b>Dossier assessment</b>	30 May 2019	MR
<b>Product performance evaluation</b>	Quarter 2-2018	MR

MR: Meets Requirements

#### Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendments	Date of report amendment
2.0	Design change and addition of multi-language translation (EN, ES, FR, PT) to IFU and change of Test kit package.	2 June 2023
3.0	1. Addition of a manufacturer for Cutting, Assembly, Buffer Preparation & dispensing, sealing, printing, packing, and shipping process. 2. Addition of suppliers for raw materials: <ul style="list-style-type: none"> <li>• Upper device,</li> <li>• Lower device,</li> <li>• Desiccant</li> <li>• Aluminum Pouch,</li> <li>• Buffer bottle</li> <li>• Capillary tube</li> <li>• Outer package/IFU/Label.</li> </ul>	6 March 2024

	3. Changes to the Outer package, Instructions for Use (IFU) and Labels.	
4.0	The addition of the manufacturing site is located in the public report of Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India.	23 April 2024
5.0	Updated labels of the capillary tube.	5 August 2025
6.0	Established the PKG label, including the print area for the procurement (PQC-IVD-2026-0031).	14 May 2026

### Intended use

According to the claim of intended use from SD Biosensor Inc, “*STANDARD Q HCV Ab Test is a rapid chromatographic immunoassay for the qualitative detection of specific antibodies to HCV present in human serum, plasma or whole blood. This test is for in vitro professional diagnostic use and intended as an aid to early diagnosis of HCV infection in patients with clinical symptoms of HCV infection. It provides only an initial screening test result. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HCV infection. STANDARD Q HCV Ab Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing*”.

### Test kit contents

Component	25 tests (product code 09HCV10D)
Test device (individually in a foil pouch with desiccant)	25
Buffer Bottle	1 x 4 mL
Capillary tube (20µL)	Pack of 25
Instructions for use	1

### Items required but not provided

- Micropipette and tip
- Blood collection tube
- Alcohol swab
- Lancet
- PPE (Personal Protective Equipment)
- Biohazard container
- Timing device

### Storage

The test kit must be stored between 2 and 40°C.

**Shelf-life upon manufacture<sup>1</sup>**

24 months.

**Warnings/limitations**

During the performance evaluation, operators noted that reading at 5 minutes (i.e. minimum recommended time) was difficult due to the strong red background and thin control line observed at that time and that 8 minutes allowed for sufficient clearing of the background colour. For more warnings, please refer to the instructions for use attached to this public report.

**Dossier assessment**

SD Biosensor Inc. submitted a product dossier for the STANDARD Q HCV Ab Test as per the "Instructions for compilation of a product dossier" (PQDx\_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 30 May 2019.

Based on the product dossier screening and assessment findings, the product dossier for the STANDARD Q HCV Ab Test meets WHO prequalification requirements.

**Manufacturing site inspection**

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pgweb/vitro-diagnostics/who-public-inspection-reports>

---

<sup>1</sup> The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

Based on the site inspection and corrective action plan review, the quality management system for the STANDARD Q HCV Ab Test meets WHO prequalification requirements.

### Product performance evaluation

STANDARD Q HCV Ab Test was evaluated by the National Reference Laboratory, Melbourne, Australia, on behalf of WHO in the 2nd quarter of 2018, according to protocol PQDx\_040, version 6.0.

### Clinical performance evaluation

In this limited laboratory-based evaluation of clinical performance characteristics, a panel of 483 serum/plasma specimens was used. The specimens were characterized using the following reference algorithm: Murex Anti-HCV version 4.0 (DiaSorin) and Monalisa anti-HCV Plus version 2.0 (Bio-Rad) in parallel, followed by CHIRON RIBA HCV 3.0 Strip Immunoblot Assay or MP diagnostics HCV Blot 3.0 on initially reactive specimens.

Clinical performance characteristics in comparison with an agreed reference standard		
	Initial (95% CI)	Final (95% CI)
Sensitivity % (N=163)	99.4 (96.6-100)	100 (97.8-100)
Specificity % (N= 320)	99.7 (98.3-100)	100 (98.9-100)
Invalid rate % (N= 483)	0.2	
Inter-reader variability % (N= 483)	0.8	

### Analytical performance evaluation

Analytical performance characteristics	
Sensitivity during seroconversion on 4 seroconversion panels (PHV913, PHV919, PHV920, PHV922, SeraCare) in comparison with a benchmark assay (Murex Anti-HCV version 4.0)	Of a total of 26 specimens, 20 were detected by the assay under evaluation; versus 21 specimens detected by the benchmark assay. Seroconversion sensitivity index of 0.25; therefore, detection is 0.25 specimens later than the benchmark assay.
Analytical sensitivity on a mixed titer panel (0810-0175, SeraCare)	All 16 specimens were correctly classified.

Analytical sensitivity on a low titer panel (0810-0192, SeraCare)	10 of 11 specimens were correctly classified. One low-titer specimen was not detected.
Lot to lot variation on a dilution panel	Lot to lot variation was within +/- 1 two-fold dilutions for 8 of 10 dilution series. For the 2 other dilution series, there was a 2 two-fold dilution difference.

### Operational characteristics and ease of use

This assay does not require laboratory equipment and can be performed in laboratories with limited facilities or non-laboratory settings.

The assay was found easy to use by the operators performing the evaluation.

Key operational characteristics	
Number of steps*	2 steps in total 1 step with precision pipetting (for serum/plasma/venous whole blood)
Time to result	10 minutes (if the test read at 8 minutes)
Endpoint stability (interval)	15 minutes (the test can be read between 5 and 20 minutes after the addition of diluent)
Internal QC	Yes, the control line on the test device (reagent control)

\* Definition: each action required to obtain a result (excluding specimen collection, device preparation – opening the pouch), e.g. for RDTs: add specimen, add buffer (2 steps).

Based on these results, the performance evaluation for the STANDARD Q HCV Ab Test meets the WHO prequalification requirements.

### Labelling review

The labelling submitted for the STANDARD Q HCV Ab Test was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

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
# STANDARD Q HCV Ab\_25T\_09HCV10D\_인도제조소 O

Unit : mm


<b>Classification</b>	Package	<b>Color</b>	3 Color (K, Pantone 2415C, 2945C)
<b>Document number</b>	PINOB25HCV1ML4R0-WHO	<b>Post-processing</b>	UV Coating + Uncoated Area / 3-side sealed
<b>Size</b>	W155 * D124 * H71	<b>Issue date</b>	2025.10.27
<b>Material/gram</b>	GC2 / 270g	<b>Department</b>	Product Planning Division

# STANDARD Q HCV Ab\_25T\_09HCV10D\_인도제조소 X

Unit : mm

<b>Classification</b>	Package	<b>Color</b>	3 Color (K, Pantone 2415C, 2945C)
<b>Document number</b>	PINXB25HCV1ML4R0-WHO	<b>Post-processing</b>	UV Coating + Uncoated Area / 3-side sealed
<b>Size</b>	W155 * D124 * H71	<b>Issue date</b>	2025.10.27
<b>Material/gram</b>	GC2 / 270g	<b>Department</b>	Product Planning Division

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


# STANDARD Q HCV Ab\_25T\_09HCV10D\_한국생산용

Unit : mm

<b>Classification</b>	Package	<b>Color</b>	3 Color (K, Pantone 2415C, 2945C)
<b>Document number</b>	PB25HCV1ML4R0-WHO	<b>Post-processing</b>	Gloss Lamination / 3-side sealed
<b>Size</b>	W155 * D124 * H71	<b>Issue date</b>	2025.10.27
<b>Material/gram</b>	GC2 / 270g	<b>Department</b>	Product Planning Division

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



EC REP Cmc Medical Devices & Drugs S.L.  
C/ Horacio Lengo n18 C.P. 29008 Málaga-Spain

(01)08809899620106

P000302502REV.0

### **Instructions for use<sup>2</sup>**

---

<sup>2</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

EN REF: QHCV018  
Cat. No.: 09HCV10D


## STANDARD Q HCV Ab

STANDARD™ Q HCV Ab Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

SD BIOSENSORS

### KIT CONTENTS

### MATERIALS REQUIRED BUT NOT PROVIDED

- Micropipette and tip
- Blood collection tube
- Alcohol swab
- Lancet
- PPE (Personal Protective Equipment)
- Biohazard container
- Timing device

### SPECIMEN COLLECTION AND STORAGE

#### Serum

- Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- Store the serum in a refrigerator(2-8°C/36-46°F) or room temperature(15-25°C/59-77°F) until 72 hours.
- Do not freeze-thawing Serum more than 3 times.

#### Plasma

- Collect the venous whole blood into the commercially available anticoagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
- Store the plasma in a refrigerator(2-8°C/36-46°F) or room temperature(15-25°C/59-77°F) until 72 hours.
- Do not freeze-thawing Plasma more than 3 times.

#### Whole blood

- Capillary whole blood should be collected aseptically by fingerstick.
- Select the finger that is free from callus. Gently rub the finger to warm it to stimulate blood circulation.
- Clean the area to be lanced with an alcohol swab and dry it completely.
- Pierce with a sterile lancet. The finger should NOT be squeezed.
- Collect the capillary whole blood to the black line of the capillary tube for the testing.
- The capillary whole blood must be tested immediately after collection.

#### Venous whole blood

- Collect the venous whole blood into the commercially available anticoagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
- Store the venous whole blood in a refrigerator(2-8°C/36-46°F) or room temperature(15-25°C/59-77°F) until 60 hours.
- Do not use hemolyzed blood specimen.


**CAUTION**

- Anticoagulants such as heparin, EDTA or sodium citrate do not affect the test result.
- Use separate disposable materials for each specimen in order to avoid cross-contamination which can cause erroneous results.
- Collect the specimen followed by CLSI recommended collection method.

### PREPARATION AND TEST PROCEDURE

#### Preparation

- Carefully read the instructions for using the STANDARD Q HCV Ab Test.
- Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
- Open the foil pouch, and check the test device and the desiccant inside the foil pouch.

#### Test Procedure

##### [For using a micropipette]

- Collecting of Specimen  
Collect the 10µl of serum/plasma or 20µl of venous whole blood using a micropipette.


##### Adding of specimen

- Hold the micropipette at 90° angle to the testing device and add collected specimen into the specimen well.

### 3. Dropping of Buffer

- Hold the buffer bottle at 90° angle to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.

### 4. Reading Time

- Read the test results after 5 minutes. The test can be read up to 20 minutes.


### CAUTION

- Do not read test results after 20 minutes. It may give false results.
- If the T line does not appear until 5 minutes, make sure to observe the appearance of the T line until 20 minutes.

### [For using a capillary tube (20µl)]

#### 1. Collecting of Specimen

- Collect the 20µl of capillary whole blood to the black line of a capillary tube.

#### 2. Adding of specimen

- Place the capillary tube containing the blood specimen to the middle of the specimen pad at an upright (vertical) position.
- Wait until all the blood is transferred from the capillary tube to the specimen pad.

### 3. Dropping of Buffer

- Hold the buffer bottle at 90° to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.

### 4. Reading Time

- Read the test results after 5 minutes. The test can be read up to 20 minutes.

**CAUTION**

- Do not read test results after 20 minutes. It may give false results.
- If the T line does not appear until 5 minutes, make sure to observe the appearance of the T line until 20 minutes.

### INTERPRETATION OF TEST RESULT

\*"C" Control Line | "T" Test Line

#### Negative

The presence of only one colored band ("C" Control line) within the result window indicates a negative result.

#### Positive

The presence of two colored bands ("C" Control line and "T" Test line) within the result window, no matter which band appears first, indicates a HCV Ab positive result. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a HCV Ab positive result.

#### Invalid

If the control band ("C" Control line) is not visible within the result window, the result is considered invalid. The directions may not have been followed correctly or the test device may have deteriorated. Re-test with a new specimen and a new test device.

**\*This test only indicates the presence of HCV antibodies in human serum and plasma specimen and should not be used as the sole criteria for the diagnosis of HCV. As with other diagnostic tests, all test results should be considered with other clinical history available.**

### EXPLANATION AND SUMMARY

#### Introduction

Hepatitis C virus (HCV) is one of several hepatitis viruses that can cause inflammation of the liver. It is a bloodborne virus and is most commonly transmitted through unsafe injection practices, inadequate sterilization of medical equipment and the transfusion of unscreened blood and blood products. HCV can cause both acute and chronic hepatitis infection. Acute HCV infection is a short-term viral infection, and is usually asymptomatic. About 15-45% of infected persons spontaneously improve or resolve the infection within just several months without treatment. However, the remaining 55-85% of infected persons will develop chronic HCV infection. Chronic HCV infection is a serious disease that can result in long-term problems in the liver, including liver damage, cancer, even death. According to the World Health Organization (WHO), about 130-150 million people globally have chronic HCV infection, with more than 350,000 people dying from Hepatitis C-related liver diseases each year. Antiviral medicines can cure approximately 90% of persons with HCV, thereby reducing the risk of death, but access to diagnosis is difficult. Establishing best practices for early diagnosis for HCV infection can prevent health problems that may result from infection and prevent transmission of the virus. STANDARD Q HCV Ab Test provides significantly fast, easy and accurate system to detect the specific antibodies to HCV in human serum, plasma or whole blood.

### Intended use

STANDARD Q HCV Ab Test is a rapid chromatographic immunoassay for the qualitative detection of specific antibodies to HCV present in human serum, plasma or whole blood. This test is for *in vitro* professional diagnostic use and intended as an aid to early diagnosis of HCV infection in patients with clinical symptoms of HCV infection. It provides only an initial screening test result. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HCV infection. STANDARD Q HCV Ab Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.

### Test principle

STANDARD Q HCV Ab Test has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimen. Monoclonal anti-NS3 and monoclonal anti-Core are coated on the control line region and monoclonal anti-human IgG is coated on the test line region. Four recombinant HCV antigens from the core, NS3, NS4 and NS5 regions conjugated with colloidal gold particles are used as detectors for HCV antibodies. During the test, HCV antibodies in the specimen interact with recombinant HCV antigens conjugated with colloidal gold particles making antibody-antigen gold particle complex. This complex migrates to the membrane via capillary action until the test line, where it will be captured by the monoclonal anti-human IgG. A violet test line will be visible in the result window if HCV antibodies are present in the specimen. The intensity of violet test line will vary depending upon the concentration of HCV antibodies present in the specimen. If HCV antibodies are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

### KIT STORAGE AND STABILITY

- Store the kit unopened at 2-40°C/36-104°F, out of direct sunlight.
- Do not open the aluminum pouch until you are ready to perform a test.
- Discard the testing device if it has been exposed for more than 30 minutes after opening the pouch and not used.
- Do not use the kit after 24 months from the manufacturing date.
- The buffer provided in the kit is stable until the expiry date of the kit after opening its cap, if it is tightly closed.

### WARNINGS AND PRECAUTIONS

- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use buffer of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- A moisture indicator with desiccant within foil pouch is to absorb moisture and keep humidity from affecting products. If the desiccant is green, discard the test device.

### TEST PROCEDURE

#### Preparation

- Carefully read the instructions for using the STANDARD Q HCV Ab Test.
- Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
- Allow the STANDARD Q HCV Ab Test components and specimen to come to room temperature(15-25°C/59-77°F) for 30min prior to testing.
- Check that the test device packaging is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from orange to green), throw away the test device and take another test device packaging. If the color of the buffer does not show a change, you can use the test. Throw away the buffer bottle in the non-sharps (non-infectious) disposal container.
- Procedure method should be followed for the specific specimen type being tested.

#### Test procedure

- For serum/plasma/venous whole blood specimen**
  - Collect the 10µl of serum/plasma or 20µl of venous whole blood specimen using a micropipette.
  - Add the collected specimen to the specimen well of the test device.
  - Add 3 drops of buffer into the specimen well of the test device.
  - Read the test results after 5 minutes. Test can be read up to 20 minutes.

- For capillary whole blood specimen**

- Clean a fingertip by wiping with an alcohol swab.
- Dry and pierce the wiped fingertip with a lancet to bleed.
- Collect the 20µl of capillary whole blood to the black line of the capillary tube (20µl).
- Add the collected whole blood to the specimen well of the test device.
- Add 3 drops of buffer into the specimen well of the test device.
- Read the test results after 5 minutes. Test can be read up to 20 minutes.

**CAUTION**

- Do not read test results after 20 minutes. It may give false results.
- Place the test device on the flat surface after applying the specimen.
- Use correct amount of specimen and buffer.

### LIMITATION OF TEST

- The test should be used for the detection of HCV antibodies in human serum, plasma or whole blood specimens.
- Neither the quantitative value nor the ratio of HCV antigen to antibody can be determined by this qualitative test.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- Although the STANDARD Q HCV Ab Test is very accurate for detecting HCV antibodies, specimens with antibody titres below the device's limit of detection may produce false negative results.
- This test only indicates the presence of HCV antibodies in human serum, plasma or whole blood specimen and should not be used as the sole criteria for the diagnosis of HCV. As with other diagnostic tests, all test results should be considered with other clinical history available to the physician.
- Specimen from patients with elevated IgG (from approximately 140%) may show a lower positive signal.

### Diagnostic Sensitivity

The Diagnostic Sensitivity of STANDARD Q HCV Ab Test, calculated on 413 positive specimens, was 100.00%(413/413). The specimens used are confirmed by CE marked comparator devices.

	STANDARD Q HCV Ab Test	
	Negative	Positive
Anti-HCV positive serum/plasma	0	210
Anti-HCV positive whole blood / plasma couples	0	101
Anti-HCV with known genotype	0	102
TOTAL	0	413
Sensitivity	413/413 = 100.00%	[95% CI: 99.08% - 100.00%]

### QUALITY CONTROL

- A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent is reactive.
- Control materials are not supplied with this test kit. However, it is recommended that the positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
- If there is a problem with the result such as invalid result, retest with a new kit and/or new specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

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### Diagnostic Specificity

The Diagnostic Specificity of STANDARD Q HCV Ab Test, calculated on 1500 negative specimens, was 97.67%(1465/1500). The specificity on plasma and whole blood was 97.20%(972/1000) and 98.60%(493/500) respectively. We claim that the performance of STANDARD Q HCV Ab Test for plasma and whole blood is equivalent because the test result for both specimen type are within the range of each 95% confidence interval. The donor specimens used are confirmed by CE marked comparator devices.

	STANDARD Q HCV Ab Test	
	Negative	Positive
Blood donor EDTA plasma specimen	972	28
Blood donor whole blood specimen	493	7
Specificity of plasma	972/1000 = 97.20%	[95% CI: 95.98% - 98.06%]
Specificity of whole blood	493/500 = 98.60%	[95% CI: 97.14% - 99.32%]
Total Specificity	1465/1500 = 97.67%	[95% CI: 96.77% - 98.32%]

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
### KIT CONTENTS

#### ■ For Cat. No.: 09HCV10D

1. Test device (individually in a foil pouch with desiccant) x 25
2. Capillary tube (20 µl) x 25
3. Buffer bottle x 1
4. Instructions for use x 1

#### ■ For Cat. No.: 09HCV11D

1. Test device (individually in a foil pouch with desiccant) x 25
2. Capillary tube (20 µl) x 25
3. Buffer bottle x 1
4. Sterile lancet x 25
5. Alcohol swab x 25
6. Instructions for use x 1

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### SPECIMEN COLLECTION AND STORAGE

#### ■ Serum

1. Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
2. Store the serum in a refrigerator (2 – 8°C / 36 – 46°F) or room temperature (15 – 25°C / 59 – 77°F) until 72 hours.
3. Do not freeze-thawing Serum more than 3 times.

#### ■ Plasma

1. Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
2. Store the plasma in a refrigerator (2 – 8°C / 36 – 46°F) or room temperature (15 – 25°C / 59 – 77°F) until 72 hours.
3. Do not freeze-thawing Serum more than 3 times.

#### ■ Whole blood

##### ■ Capillary whole blood

1. Capillary whole blood should be collected aseptically by fingertip.
2. Select the finger that is free from callus. Gently rub the finger to warm it to stimulate blood circulation.
3. Clean the area to be lanced with an alcohol swab and dry it completely.
4. Pierce with a sterile lancet. The finger should NOT be squeezed.
5. Collect the capillary whole blood to the black line of the capillary tube for the testing.
6. The capillary whole blood must be tested immediately after collection.

##### ■ Venous whole blood


1. Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
2. Store the venous whole blood in a refrigerator (2 – 8°C / 36 – 46°F) or room temperature (15 – 25°C / 59 – 77°F) until 60 hours.
3. Do not use hemolyzed blood specimen.

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### PREPARATION AND TEST PROCEDURE

#### ■ Preparation

1. Carefully read the instructions for using the STANDARD Q HCV Ab Test.
2. Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed. Open the foil pouch, and check the test device and the desiccant inside the foil pouch.
3. Allow the STANDARD Q HCV Ab Test components and specimen to come to room temperature (15 – 25°C / 59 – 77°F) for 30 min prior to testing.
4. Check that the test device packaging is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from yellow to green), throw away the test device and take another test device packaging. If the color of the buffer does not show a change, you can use the test. Throw away the buffer bottle in the non-sharps (non-infectious) disposal container.
5. Procedure method should be followed for the specific specimen type being tested.

### PERFORMANCE CHARACTERISTICS

#### ■ Potential Interfering Substances

STANDARD Q HCV Ab Test was evaluated with the following interfering substances present in specimen in order to assess their potential effect on the assay performance as per CLSI guideline EP7-A2. It was not affected by interfering materials such as artemether-lumefantrine (Malaria drug), doxycycline hyclate (Malaria drug), quinine (Malaria drug), lamivudine (Retroviral medication), ribavirin (HCV drug), daclatasvir (HCV drug), acetaminophen, acetylsalicylic acid, Ibuprofen, erythromycin (antibiotic), ciprofloxacin (antibiotic), caffeine, ethanol, bilirubin (Unconjugate). There was no interference with human anti-mouse antibody, whole blood of pregnant women and whole blood having elevated levels of haemoglobin and elevated levels of C-reactive protein.

#### ■ Potential Cross-reacting Substances

STANDARD Q HCV Ab Test had no cross-reaction with the following microorganism such as other hepatitis (Hepatitis A virus, Hepatitis B virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue fever virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, Escherichia coli, Salmonella), HbsAg positive plasma, Rubella, CMV, Tick borne encephalitis, West Nile Virus, Treponema pallidum, HAV, Leishmania, Brucella positive plasma.

STANDARD Q HCV Ab Test had no cross-reaction with the following microorganism such as other hepatitis (Hepatitis A virus, Hepatitis B virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue fever virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, Escherichia coli, Salmonella), HbsAg positive plasma, Rubella, CMV, Tick borne encephalitis, West Nile Virus, Treponema pallidum, HAV, Leishmania, Brucella positive plasma.

STANDARD Q HCV Ab Test had no cross-reaction with the following microorganism such as other hepatitis (Hepatitis A virus, Hepatitis B virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue fever virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, Escherichia coli, Salmonella), HbsAg positive plasma, Rubella, CMV, Tick borne encephalitis, West Nile Virus, Treponema pallidum, HAV, Leishmania, Brucella positive plasma.


STANDARD Q HCV Ab Test had no cross-reaction with the following microorganism such as other hepatitis (Hepatitis A virus, Hepatitis B virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue fever virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, Escherichia coli, Salmonella), HbsAg positive plasma, Rubella, CMV, Tick borne encephalitis, West Nile Virus, Treponema pallidum, HAV, Leishmania, Brucella positive plasma.

### ■ Test Procedure

#### 1. Collecting of Specimen


##### ■ For serum/plasma/venous whole blood specimen

- ① Collect the 10 µl of serum/plasma or 20 µl of venous whole blood specimen using a micropipette.

##### ■ For capillary whole blood specimen


- ① Clean a fingertip by wiping with an alcohol swab.
- ② Dry and pierce the wiped fingertip with a lancet to bleed.
- ③ Collect the 20 µl of capillary whole blood to the black line of the capillary tube (20 µl).

#### 2. Adding of Specimen


##### ■ For serum/plasma/venous whole blood specimen

Add the collected specimen to the specimen well of the test device. Hold the micropipette at 90° angle to the testing device and add collected specimen into the specimen well.


##### ■ For capillary whole blood specimen

Add the collected whole blood to the specimen well of the test device. Place the capillary tube containing the blood specimen to the middle of the specimen pad at an upright (vertical) position. Wait until all the blood is transferred from the capillary tube to the specimen pad.

#### 3. Dropping of Buffer

Hold the buffer bottle at 90° to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.

#### 4. Read the test results after 5 minutes. Test can be read up to 20 minutes.

- Do not read test results after 20 minutes. It may give false results.
- Place the test device on the flat surface after applying the specimen.
- Use correct amount of specimen and buffer.
- If the T line does not appear until 5 minutes, make sure to observe the appearance of the T line until 20 minutes.

### INTERPRETATION OF TEST RESULT

Test result	Example	Description
Negative		The presence of only one colored band ("C" Control line) within the result window indicates a negative result.
Positive		The presence of two colored bands ("C" Control line and "T" Test line) within the result window, no matter which band appears first, indicates a positive result. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a HCV Ab positive result.
Invalid		If the control band ("C" Control line) is not visible within the result window, the result is considered invalid. The directions may not have been followed correctly or the test may have deteriorated. Re-test with a new patient specimen and a new test device.

- This test only indicates the presence of HCV antibodies in human serum and plasma specimen and should not be used as the sole criteria for the diagnosis of HCV. As with other diagnostic tests, all test results should be considered with other clinical history available.
- The presence of any line no matter how faint it is, should be considered as a line formed.
- This test is for screening purposes. Confirmatory testing according to national guidelines is recommended to confirm the infection status.

### QUALITY CONTROL

1. A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent are reactive.
2. Control materials are not supplied with this test kit. However, it is recommended that the positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
3. If there is a problem with the result such as invalid result, retest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

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#### ■ Test Principle

STANDARD Q HCV Ab Test has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimen. Monoclonal anti-NS3 and monoclonal anti-Core are coated on the control line region and monoclonal anti-human IgG is coated on the test line region. Four recombinant HCV antigens from the core, NS3, NS4 and NS5 regions conjugated with colloidal gold particles are used as detectors for HCV antibodies. During the test, HCV antibodies in the specimen interact with recombinant HCV antigens conjugated with colloidal gold particles making antibody-antigen gold particle complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the monoclonal anti-human IgG. A violet test line would be visible in the result window if HCV antibodies are present in the specimen. The intensity of violet test line will vary depending upon the concentration of HCV antibodies present in the specimen. If HCV antibodies are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

### KIT STORAGE AND STABILITY

1. Store the kit unopened at 2 – 40°C / 36 – 104°F, out of direct sunlight.
2. Do not open the aluminum pouch until you are ready to perform a test.
3. Discard the testing device if it has been exposed for more than 30 minutes after opening the pouch and not used.
4. Do not use the kit after 24 months from the manufacturing date.
5. The buffer provided in the kit is stable until the expiry date of the kit after opening its cap, if it is tightly closed.

### WARNINGS AND PRECAUTIONS

1. Do not re-use the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use buffer of another lot.
4. Do not smoke, drink or eat while handling specimen.
5. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Handle all specimens as if they contain infectious agents.
8. Observe established precautions against microbiological hazards throughout testing procedures.
9. Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
10. A moisture indicator with desiccant within foil pouch is to absorb moisture and keep humidity from affecting products. If the desiccant is green, discard the test device.

### LIMITATION OF TEST

1. The test should be used for the detection of HCV antibodies in human serum, plasma or whole blood specimens.
2. Neither the quantitative value nor the ratio of HCV antigen to antibody can be determined by this qualitative test.
3. Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
4. Although the STANDARD Q HCV Ab Test is very accurate for detecting HCV antibodies, specimens with antibody titres below the device's limit of detection may produce false negative results.
5. This test only indicates the presence of HCV antibodies in human serum, plasma or whole blood specimen and should not be used as the sole criteria for the diagnosis of HCV. As with other diagnostic tests, all test results should be considered with other clinical history available to the physician.
6. Specimen from patients with elevated IgG (from approximately 140%) may show a lower positive signal.

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### SYMBOLE

REF	Reference number	Caution	Use by	LOT	Batch code	Consult Instructions for Use	Do not re-use
IVD	<i>In vitro</i> Diagnostics	Note	Manufacturer		Date of manufacture	Contains Sufficient for <=3> Tests	Keep away from sunlight
	Indicate that you should keep the product dry		To indicate the temperature limitations in which the transport package has to be kept and handled.		Do not use if packaging is damaged		Authorized representative in the european community

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