## WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

**Product: STANDARD Q HIV 1/2 Ab 3-Line Test**  
**WHO reference number: PQDx 0383-117-00**

The STANDARD Q HIV 1/2 Ab 3-Line Test with product codes 09HIV30D and 09HIV30DM, manufactured by SD Biosensor, Inc., Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 10 June 2020.

### Summary of WHO prequalification assessment for STANDARD Q HIV 1/2 Ab 3-Line Test

	Date	Outcome
<b>Prequalification listing</b>	10 June 2020	listed
<b>Dossier assessment</b>	25 May 2020	MR
<b>Product performance evaluation</b>	Quarter 4-2019	MR

MR: Meets Requirements

### Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendment	Date of report amendment
2.0	1. Addition of a manufacturer for Cutting, Assembly, Buffer Preparation & dispensing, sealing, printing, packing, and shipping process. 2. Addition of suppliers for raw materials: <ul style="list-style-type: none"> <li>• Upper device,</li> <li>• Lower device,</li> <li>• Desiccant,</li> <li>• Aluminum Pouch,</li> <li>• Buffer bottle,</li> <li>• Capillary tube,</li> <li>• Lancet,</li> <li>• Alcohol swab,</li> <li>• Outer package/IFU/Label.</li> </ul>	7 March 2024

	3. Changes to the Outer package, Instructions for Use (IFU) and labels.	
3.0	The addition of the manufacturing site is located in the public report of Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India.	7 May 2024
4.0	Addition of a manufacturing facility for the packing and shipping process at Colexa Biosensor Limited in Nigeria.	17 January 2025
5.0	Updated labels of capillary tube, sterile lancet, and alcohol swab.	5 August 2025
6.0	<p>1. New Package labels for Print area addition Established the PKG label, adding the print area to deal with the requirements for the procurement.</p> <p>2. Test device pouch revision Established the test device pouch label for the exclusive product for better identification.</p> <p>3. RMS Lancet pouch revision Update of Manufacturer's Address, EU Representative symbol. Addition of the Korean language and the KGMP symbol to fulfil the requirements for domestic distribution (PQC-IVD-2025-0136).</p>	13 May 2026

### Intended use

According to the claim of intended use from SD Biosensor Inc., *"STANDARD Q HIV 1/2 Ab 3-Line Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O and HIV 2 in human serum, plasma or venous and capillary whole blood. The test is for in vitro diagnostic use and intended as an aid to diagnosis of HIV in symptomatic, asymptomatic populations and persons at risk of HIV infection. The test is not intended for use on neonates or infants below 2 years. The test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource limited countries. The product is not for blood donor screening. The test kit is not automated and does not require any additional instruments. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HIV infection."*

**Test kit contents**

Component	25 tests (product code 09HIV30D )	25 tests (product code 09HIV30DM)
Test device (individually in a foil pouch with desiccant)	25	25
Buffer Bottle	1 x 4 mL	1 x 4 mL
Capillary tube (20 µL)	Pack of 25	
Instructions for use	1	1
Sterile lancet	25	
Alcohol swab	25	

**Items required but not provided****For Cat. No.: 09HIV30D**

- Anti-coagulant tube containing heparin, EDTA or sodium citrate
- Micropipette and tip
- PPE (Personal Protective Equipment)
- Biohazard container

**For Cat. No.: 09HIV30DM**

- Anti-coagulant tube containing heparin, EDTA or sodium citrate
- Micropipette and tip
- Sterile lancet
- Alcohol swab
- PPE (Personal Protective Equipment)
- Biohazard container
- Capillary whole blood collection tool

**Storage**

The test kit must be stored between 2 and 40 °C.

**Shelf-life upon manufacture<sup>1</sup>**

24 months.

<sup>1</sup> The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

## Dossier assessment

SD Biosensor Inc. submitted a product dossier for STANDARD Q HIV 1/2 Ab 3-Line Test as per the "Instructions for compilation of a product dossier" (PQDx\_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the discrepancies found during the dossier review were accepted on 25 May 2020.

Based on the product dossier screening and assessment findings, the product dossier for the STANDARD Q HIV 1/2 Ab 3-Line Test meets WHO prequalification requirements.

## Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pgweb/vitro-diagnostics/who-public-inspection-reports>

Based on the site inspection and corrective action plan review, the quality management system for the STANDARD Q HIV 1/2 Ab 3-Line Test meets WHO prequalification requirements.

## Product performance evaluation

STANDARD Q HIV 1/2 Ab 3-Line Test (SD Biosensor) was evaluated by the Institute of Tropical Medicine, Belgium, on behalf of WHO in the 4th quarter of 2019, according to protocol PQDx\_030, version 10.

### Clinical performance evaluation

In this limited laboratory-based evaluation of clinical performance characteristics, a panel of 1200 serum/plasma specimens was used. The specimens were characterized using the following reference algorithm: Vironostika HIV Ag/Ab (bioMérieux) and Enzygnost Anti-HIV 1/2 Plus (Siemens Healthcare Diagnostics), followed by INNO-LIA HIV I/II Score (Fujirebio) for initially reactive results.

Clinical performance characteristics in comparison with an agreed reference standard		
	Initial (95% CI)	Final (95% CI)
Sensitivity % (N=470)	100 (99.2-100)	100 (99.2-100)
Specificity % (N=730)	99.2 (98.2-99.7)	99.3 (98.4-99.8)
Invalid rate % (N=1200)	0.08	
Inter-reader variability % (N=1200)	2.1	

Among 449 HIV-1 positive specimens, STANDARD Q HIV 1/2 Ab 3-Line Test showed the presence of both HIV 1 and HIV-2 lines in 63 (14 %) specimens, although in most cases (n=60), the HIV-2 line was weaker than the HIV-1 line, which is interpreted as HIV-1 positive result according to the IFU of the assay.

Among 21 HIV-2 positive specimens, STANDARD Q HIV 1/2 Ab 3-Line Test showed the presence of both HIV-1 and HIV-2 lines in 18 (86 %) specimens, with medium to strong reactivity in most cases (n=14).

For inter-reader variability, disagreements (n=25) mostly concerned the HIV-2 line in HIV-1 positive specimens (n=21) or the HIV-1 line in HIV-2 positive specimens (n=2) or in negative specimens (n=2).

### Analytical performance evaluation

Analytical performance characteristics	
Sensitivity during seroconversion on 8 seroconversion panels in comparison with a benchmark assay (Enzygnost Anti-HIV 1/2 Plus)	Of a total of 52 specimens, 23 were detected by the assay under evaluation versus 21 specimens detected by the benchmark assay. Seroconversion sensitivity index of -0.25. Therefore, detection is 0.25 specimens earlier than the benchmark assay.
Analytical sensitivity on a mixed titer panel (PRB205, SeraCare)	All 25 specimens were correctly classified.

Analytical sensitivity on the 1 <sup>st</sup> International Reference Panel for anti-HIV (NIBSC code 02/210)	All 6 HIV subtypes/groups in the panel were detected.
Lot to lot variation on a dilution panel	Lot to lot variation was within +/- 1 two-fold dilutions for 8 dilution series. For the other 2 dilution series, the last dilution tested was reactive on both lots, so the end-point could not be determined.

### Operational characteristics and ease of use

This assay does not require laboratory equipment and can be performed in laboratories with limited facilities or in non-laboratory settings.

The assay was found easy to use by the operators performing the evaluation.

Key operational characteristics	
Number of steps*	2 steps in total 1 step with precision pipetting (for serum/plasma/venous whole blood)
Time to result	12 minutes
Endpoint stability (interval)	10 minutes (the test can be read between 10 and 20 minutes after the addition of diluent)
Internal QC	Yes, reagent addition control

\* Definition: each action required to obtain a result (excluding specimen collection, device preparation – opening the pouch), e.g. for RDTs: add specimen, add buffer (2 steps).

Based on these results, the performance evaluation for the STANDARD Q HIV 1/2 Ab 3-Line Test meets the WHO prequalification requirements.

### Labelling review

The labelling submitted for the STANDARD Q HIV 1/2 Ab 3-Line Test was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries. The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

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




# STANDARD Q HIV 1/2 Ab 3-Line 25T\_09HIV30D\_인도제조소 O

Unit : mm


<b>Classification</b>	Package	<b>Color</b>	3 Color (K, Pantone 2415C, 1797C)
<b>Document number</b>	PINOB25HIV3ML4R0-WHO	<b>Post-processing</b>	UV Coating + Uncoated Area / 3-side sealed
<b>Size</b>	W155 * D124 * H71	<b>Issue date</b>	2025.10.27
<b>Material/gram</b>	GC2 / 270g	<b>Department</b>	Product Planning Division

# STANDARD Q HIV 1/2 Ab 3-Line 25T\_09HIV30D\_인도제조소 X

Unit : mm

<b>Classification</b>	Package	<b>Color</b>	3 Color (K, Pantone 2415C, 1797C)
<b>Document number</b>	PINXB25HIV3ML4R0-WHO	<b>Post-processing</b>	UV Coating + Uncoated Area / 3-side sealed
<b>Size</b>	W155 * D124 * H71	<b>Issue date</b>	2025.10.24
<b>Material/gram</b>	GC2 / 270g	<b>Department</b>	Product Planning Division

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

# STANDARD Q HIV 1/2 Ab 3-Line 25T\_09HIV30D\_한국생산용

Unit : mm


<b>Classification</b>	Package	<b>Color</b>	3 Color (K, Pantone 2415C, 1797C)
<b>Document number</b>	PB25HIV3ML4R0-WHO	<b>Post-processing</b>	Gloss Lamination / 3-side sealed
<b>Size</b>	W155 * D124 * H71	<b>Issue date</b>	2025.10.27
<b>Material/gram</b>	GC2 / 270g	<b>Department</b>	Product Planning Division

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
SD BIOSENSOR

www.sdbiosensor.com

- CONTENTS
- 1) Test device (individually in a foil pouch with desiccant) x 25
  - 2) Capillary tube (20µl) x 25
  - 3) Buffer bottle x 1
  - 4) Sterile Lancet x 25
  - 5) Alcohol swab x 25
  - 6) Instructions for use x 1

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**25 TESTS KIT**

**STANDARD Q**  
**HIV 1/2 Ab 3-Line Test**

Manufactured by  
**SD Biosensor, Inc.**  
Head office  
C-40b5th, 16, Deogyong-daero 155beon-gil, Yeongtong-gu, Seongju, Gyeongju-do, 6690, REPUBLIC OF KOREA  
Manufacturing site  
74, Osonggaengmyeong 4-ro, Osong-eup, Heungdeok-si, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA  
Locally Packaged in Nigeria by  
**Colexa Biosensor Limited**  
31/33 Industrial Avenue, Ikeja, Lagos, Nigeria.

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EC REP Cmc Medical Devices & Drugs S.L.  
C/ Horacio Lengo n18 C.P. 29008 Málaga-Spain

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


# Disposable Sterile Lancets\_Pouch\_28G\_25pcs

Unit : mm

Classification	Label	Color	1 color (K)
Document number	B06LCT2ENR4	Post-processing	
Size	W90 * H110	Issue date	2025.06.24
Material/gram	Printing on the pouch	Department	Product Planning Division

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To use, twist-off the protective cap.

### CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the based.  
Do not use if the seal has been damaged or broken.

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 Manufactured by  
**Beijing Ruicheng Medical Supplies Co., Ltd.**  
Building 5, No.8 Yanqi West Road, Yanqi  
Economic Development Zone, Huairou  
District 101407 Beijing P.R. China

 Authorized Representative  
**Lotus Nl B.V.**  
Koningin Julianaplein 10, 1e Verd,  
2595AK The Hague, Netherlands.  
Tel: +31644168999  
Email: info@lotusnl.com

B06LCT2ENR4  
Issue date: 2025.06

# Disposable Sterile Lancets

 LOT No. : **XXXXXX**

Product code : 01GL27

 MFG Date : **YYYY.MM.DD.**

 EXP Date : **YYYY.MM.DD.**

**25 PCS**  
**28G**

## INTENDED USE

To obtain a capillary blood specimen from the fingertip.

## INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

## CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base.  
Do not use if the seal has been damaged or broken.

B05SL2ENR0

Issue Date : 2020.07

### Manufactured by

Tianjin Huahong Technology Co., Ltd.  
A01, Plant B No. 278, Hangkong Road,  
Tianjin Pilot Free Trade Zone  
(Air Port Industrial Park),  
300308 Tianjin, China

### Authorized Representative

Shanghai International Holding Corp.  
GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany  
Tel:+49-40-2513175

STERILE R

**SURFACE JOB**

110MM

220MM

110MM

# Disposable Sterile Lancets

Product code :SM.01

**25PCS**

**28G**

 LOT No. : XXXXXXXXX

 MFG Date : YYYY-MM

 EXP Date : YYYY-MM

## INTENDED USE

To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

## CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base. Do not use if the seal has been damaged or broken.

Issue Date : 2023.09

 Manufactured by

Phoenix Innovative Healthcare Manufacturers Pvt. Ltd.  
E1-209, Shil Mahape Road, Electronic Zone, MIDC,  
TTC Industrial Area, Mahape, Navi Mumbai - 400 710,  
Maharashtra, India

 Authorized Representative

Med Path GmbH  
Miles-van-der-Rohe-Strasse 8 80907  
Munich, Germany

CE 0197

 STERILE

25 Pcs.

100MM

size: 95mm x 100mm No.256

95mm

## Blood Lancet

 LOT No. : 240311

Product code : BL-28G

 MFG Date : 2024-03-11

**25 PCS**

 EXP Date : 2029-02-11

**28G**

### INTENDED USE

To obtain the capillary blood specimen for test.

IFU-02-BL-01-NO.256

### INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

Issue Date: 2023-09-15

### CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base.

Do not use if the seal has been damaged or broken.

### Manufactured by

Guangzhou iCare Medical Technology Co., Ltd.  
First floor A No.8, Lianhua Port Industrial Zone,  
Lotus Mountain Bonded Area, Shilou Town,  
Panyu District , Guangzhou,  
People's Republic of China.

Tel: +86-020-31169913 Fax: +86-020-84881433

### Authorized Representative

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175 Fax: +49-40-2513175

100mm

0123

**STERILE R**

### **Instructions for use<sup>2</sup>**

---

<sup>2</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.


# STANDARD Q HIV 1/2 Ab 3-Line Test

STANDARD™ Q HIV 1/2 Ab 3-Line Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST.

**SD BIOSENSOR**

## KIT CONTENTS

### ■ For Cat. No.: 09HIV30D

- Test device (individually in a foil pouch with desiccant) x 25
- Capillary tube (20µl) x 25
- Buffer bottle x 1
- Sterile lancet x 25
- Alcohol swab x 25
- Instructions for use x 1

### ■ For Cat. No.: 09HIV30DM

- Test device (individually in a foil pouch with desiccant) x 25
- Buffer bottle x 1
- Instructions for use x 1

## MATERIALS REQUIRED BUT NOT PROVIDED

### ■ For Cat. No.: 09HIV30D

- Anti-coagulant tube containing heparin, EDTA or sodium citrate
- Micropipette and tip
- PPE (Personal Protective Equipment)
- Biohazard container

### ■ For Cat. No.: 09HIV30DM

- Anti-coagulant tube containing heparin, EDTA or sodium citrate
- Micropipette and tip
- Sterile lancet
- Alcohol swab
- PPE (Personal Protective Equipment)
- Biohazard container
- Capillary whole blood collection tool

## SPECIMEN COLLECTION AND PREPARATION

### ■ Serum

- Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant. If serum in the plain tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 3 days after collection. For prolonged storage, it should be at below -40°C / -40°F up to 3 freeze-thaw cycles.
- It should be brought to room temperature prior to use.

### ■ Plasma

- Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
- If plasma in an anti-coagulant tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 3 days after collection. For prolonged storage, it should be at below -40°C / -40°F up to 3 freeze-thaw cycles.
- It should be brought to room temperature prior to use.

### ■ Whole blood

#### [Capillary whole blood]

- Capillary whole blood should be collected aseptically by fingertip.
- Select the finger that is not calloused and gently rub the finger to warm it to stimulate blood circulation. And then clean the fingertip by wiping with an alcohol swab.
- Wait until the fingertip is dried completely and squeeze the end of the fingertip and pierce the wiped fingertip with a sterile lancet to bleed.
- Gently squeeze capillary tube and immerse open end in the center of a blood drop and release the capillary tube slowly to draw up the blood up to the 20µl black marking line on the capillary tube.
- The capillary whole blood must be tested immediately after collection.


#### [Venous whole blood]

- Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
- If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 1 – 2 days after collection.
- Do not freeze venous whole blood specimen.
- Do not use hemolyzed blood specimen.

- Anticoagulants such as heparin, EDTA or sodium citrate do not affect the test result.
- As known relevant interference, haemolytic specimen, rheumatoid factors-contained specimen and lipaemic, icteric specimen can lead to impair the test results.

## PREPARATION AND TEST PROCEDURE

### ■ Preparation

- Carefully read the instructions for using the STANDARD Q HIV 1/2 Ab 3-Line Test.
- Check that the expiry date on the back of the foil pouch has not passed.
- Allow the STANDARD Q HIV 1/2 Ab 3-Line Test components and specimen to come to room temperature(15 – 30°C / 59 – 86°F) for 30min prior to testing.
- Check that the test device package is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from yellow to green), throw away the test device and take another test device packaging.
- Procedure method should be followed for the specific specimen type being tested.

## PERFORMANCE CHARACTERISTICS

### ■ Diagnostic Sensitivity

[Plasma] The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 457 positive plasma specimens, is 99.78% (456 / 457) with a Wilson 95% confidence interval of [98.77% - 99.96%].

	STANDARD Q HIV 1/2 Ab 3-Line Test	
	HIV Negative	HIV Positive
Anti-HIV-1 positive	1	294
Anti-HIV-1 positive non-B subtype	0	40
Anti-HIV-2 positive	0	122
<b>Sensitivity</b>	<b>456/457= 99.78%</b>	

[Serum] The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 623 serum specimens, is 99.84% (622 / 623) with a Wilson 95% confidence interval of [99.10% - 99.97%].

	STANDARD Q HIV 1/2 Ab 3-Line Test	
	HIV Negative	HIV Positive
Anti-HIV-1 positive	1	589
Anti-HIV-1 positive non-B subtype	0	10
Anti-HIV-2 positive	0	23
<b>Sensitivity</b>	<b>622/623= 99.84%</b>	

### ■ Test Procedure

#### [For serum/plasma/venous whole blood specimen]

- Collect the 10 µl of serum/plasma or 20 µl of venous whole blood specimen using a micropipette.

#### [For capillary whole blood specimen]

- Select the finger that is not calloused and gently rub the finger to warm it to stimulate blood circulation. And then clean the fingertip by wiping with an alcohol swab.
- Wait until the fingertip is dried completely and pierce the wiped fingertip with a sterile lancet to bleed.
- Gently squeeze capillary tube and immerse open end in the center of a blood drop and release the capillary tube slowly to draw up the blood up to the 20 µl black marking line on the capillary tube.

### 2. Adding of Specimen

#### [For serum/plasma/venous whole blood specimen]

Add the collected specimen to the specimen well of the test device using the micropipette.

#### [ For capillary whole blood specimen]

Add the collected specimen to the specimen well of the test device. Let the end of the capillary tube to touch the pad, then press down lightly.

### 3. Dropping of Buffer

Hold the buffer bottle at 90° angle to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.

### 4. Reading Time

Read the test results between 10 to 20 minutes after adding buffer

- Do not read test results after 20 minutes. It may give false results.
- Place the test device on the flat surface before applying the specimen.
- Use correct amount of specimen and buffer.
- If the T line does not appear until 10 minutes, make sure to observe the appearance of the T line until 20 minutes.

## INTERPRETATION OF TEST RESULT

Test result	Example	Description
<b>Negative</b>		The presence of only "C" line indicates a negative result.
<b>HIV-1</b>		The presence of two lines as "C" and "T1" line indicates a positive result for HIV-1.
		In case of the presence of three lines as "C", "T1" and "T2", if the intensity of the "T1" line is stronger than "T2" line, it should be interpreted as HIV-1 positive.
		The presence of two lines as "C" and "T2" line indicates a positive result for HIV-2.
<b>Positive</b>		In case of the presence of three lines as "C", "T2" and "T1", if the intensity of the "T2" line is stronger than "T1" line, it should be interpreted as HIV-2 positive.
		The presence of three lines as "C" and equivalent intensity of "T2" and "T1", it should be interpreted as HIV-1 and HIV-2 both positive.
<b>Invalid</b>		No presence of "C" line indicates an invalid result. It is recommended that the specimen be retested using a new device.

- Interpret faint lines of "T1" and/or "T2" as the positive for HIV-1 and HIV-2.
- In case the intensity of "T1" and "T2" line is similar, confirmatory testing is recommended to confirm the virus type or a co-infection.
- The presence of any line no matter how faint it is, should be considered as a line formed.
- Positive results should be considered in conjunction with the clinical history and other data available.
- This test is for screening purposes. Confirmatory testing according to national guidelines is recommended to confirm the infection status.

## QUALITY CONTROL

- A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent are reactive.
- Control materials are not supplied with this test kit. However, it is recommended that the positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
- If there is a problem with the result such as invalid result, retest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

## EXPLANATION AND SUMMARY

### ■ Introduction

AIDS is caused by two known types of HIV (human immunodeficiency virus), HIV type 1 and HIV type 2. HIV type 1 (HIV-1) is found in patients with AIDS, AIDS-related complex (ARC), and asymptomatic infected individuals at high risk for AIDS. The virus is transmitted by sexual contact, by exposure to infected blood or blood products, or from an infected mother to her fetus or infant. The infection of HIV type 2 (HIV-2) is endemic only in West Africa, and it has been identified in individuals who had sexual relations with individuals from that geographic region. HIV-2 is similar to, but distinct from, HIV-1. Both viruses have similar morphology and lymphotropism, and the modes of transmission appear to be identical. Within the two major HIV types, there is significant variation. HIV-1 has been divided into four groups: group M (for major) including at least ten subtypes; group O (for outlier), group P, and group N (for non-M, non-O). Similarly, the HIV-2 strains have been classified into at least five subtypes (A through E). STANDARD Q HIV 1/2 Ab 3-Line Test can detect antibodies to both HIV-1 and HIV-2, which first appear significantly later, 20-45 days after infection. Detecting HIV earlier with STANDARD Q HIV 1/2 Ab 3-Line Test is helping to prevent future transmission during extremely infectious stage.

### ■ Intended use

STANDARD™ Q HIV 1/2 Ab 3-Line Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O and HIV-2 in human serum, plasma or venous and capillary whole blood. The test is for *in vitro* diagnostic use and intended as an aid to diagnosis of HIV in symptomatic, asymptomatic populations and persons at risk of HIV infection. The test is not intended for use on neonates or infants below 2 years. The test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited countries. The product is not for blood donor screening. The test kit is not automated and does not require any additional instruments. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HIV Virus infection.

### ■ Test Principle

STANDARD Q HIV 1/2 Ab 3-Line Test has "T1", "T2" and "C" line pre-coated with recombinant HIV-1 GP41 protein / recombinant HIV-1 subtype O GP41, recombinant HIV-2 GP36 protein and monoclonal anti-chicken IgY respectively. The anti-HIV-1/anti-HIV-1 subtype O in patient specimens interacts with the recombinant HIV-1 GP41-gold / recombinant HIV-1 subtype O GP41-gold and the anti-HIV-2 in patient specimens interacts with the recombinant HIV-2 GP36-gold in the conjugation pad. The complex of gold conjugated antigens and antibodies moves along the membrane chromatographically to the membrane with buffer and is captured by the recombinant HIV-1 and HIV-2 antigens on the test regions (T1 and T2). If the antibodies against HIV-1/2 are in the patient specimen, visible lines are formed in the test region. The chicken IgY antibody conjugated with colloidal gold particles are used as detectors for "C" control line. The control line should always appear if the test procedure is performed properly.

## KIT STORAGE AND STABILITY

- Store the sealed pouch and the buffer provided in the kit at 2 – 40°C / 36 – 104°F out of the direct sunlight for the duration of its shelf life.
- Do not open the foil pouch until you are ready to perform a test. Test device can be used immediately after being opened.
- Close the buffer cap tightly after using, and then store it at 2 – 40°C / 36 – 104°F out of the direct sunlight. It is stable until the expiry date of the kit and the buffer label after opening its cap, if it is tightly closed.

## WARNINGS AND PRECAUTIONS

- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use after the expiration date.
- Do not use the buffer bottle of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Allow Kit components to reach room temperature (15 – 30°C / 59 – 86°F) before performing a test.
- Follow necessary precautions when handling specimens with this test. Use personal protective equipment (PPE) consistent with current guidelines.
- The buffer contains 0.01% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.

## LIMITATION OF TEST

- All three test lines ("T1", "T2" and "C") may develop when tested with specimen containing high titers of HIV-1 antibodies. Hence, reactive test bands for both HIV-1 and HIV-2 may not indicate mixed infection but may result from the cross reactivity of HIV-1 and HIV-2 because of the similarity of their genomic structure.
- A positive result indicates presence of antibodies to HIV-1 and/or HIV-2. However, a positive result does not indicate a conclusive HIV infection diagnosis. A positive result should be confirmed by a supplemental test.
- A negative result does not eliminate the possibility of infection with HIV-1/2. The STANDARD Q HIV 1/2 Ab 3-Line Test may not detect extremely low concentration of the antibody to HIV 1/2. For negative result with this kit, additional test using other clinical method is necessary.
- The test results alone should not be used in diagnosis of infection with HIV-1/2. For overall clinical diagnosis, results must be interpreted with the patient's clinical symptoms, and clinical history.

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### PERFORMANCE CHARACTERISTICS

- Owen SM et al. Alternative Algorithms for Human Immunodeficiency Virus Infection Diagnosis Using Tests That Are Licensed in the United States. J Clin Microbiol 46:1588-1595, 2008.
- Barre-Sinoussi F, Chermann JC, Rey F, et al: Isolation of a T-lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). Science 220:868-871, 1983.
- Centers for Disease Control: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Rep 37:371-388, 1988.
- Hoff R, Weblein BJ, Schwenzler M, et al: Specific antibodies to HIV-2 detected in an anonymous newborn blood specimen from Massachusetts. Fourth Consensus Conference on Testing for Human Retroviruses, March 1989.
- Charneau P, Borman AM, Quillan C, et al: Isolation and envelope sequence of a highly divergent.
- Perspectives in Disease Prevention and Health Promotion Update: Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and Other Bloodborne Pathogens in Health-Care Settings. MMWR, June 24, 1988;37(24):377-388.

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
# STANDARD Q HIV 1/2 Ab 3-Line Test

STANDARD™ Q HIV 1/2 Ab 3-Line Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

## KIT CONTENTS

- **For Cat. No.: 09HIV30D**
  - Test device (individually in a foil pouch with desiccant) x 25
  - Capillary tube (20µl) x 25
  - Buffer bottle x 1
  - Sterile lancets x 25
  - Alcohol swab x 25
  - Instructions for use x 1
- **For Cat. No.: 09HIV30DM**
  - Test device (individually in a foil pouch with desiccant) x 25
  - Buffer bottle x 1
  - Instructions for use x 1

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- **For Cat. No.: 09HIV30D**
  - Anti-coagulant tube containing heparin, EDTA or sodium citrate
  - Micropipette and tip
  - PPE (Personal Protective Equipment)
  - Biohazard container
- **For Cat. No.: 09HIV30DM**
  - Anti-coagulant tube containing heparin, EDTA or sodium citrate
  - Micropipette and tip
  - Sterile lancet
  - Alcohol swab
  - PPE (Personal Protective Equipment)
  - Biohazard container
  - Capillary whole blood collection tool

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- **Serum**
  - Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
  - If serum in the plain tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 3 days after collection. For prolonged storage, it should be at below -40°C / -40°F up to 3 freeze-thaw cycles.
  - It should be brought to room temperature prior to use.
- **Plasma**
  - Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
  - If plasma in an anti-coagulant tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 3 days after collection. For prolonged storage, it should be at below -40°C / -40°F up to 3 freeze-thaw cycles.
  - It should be brought to room temperature prior to use.

- **Whole blood (Capillary whole blood)**
  - Capillary whole blood should be collected aseptically by fingertip.
  - Select the finger that is not calloused and gently rub the finger to warm it to stimulate blood circulation. And then clean the fingertip by wiping with an alcohol swab.
  - Wait until the fingertip is dried completely and squeeze the end of the fingertip and pierce the wiped fingertip with a sterile lancet to bleed.
  - Gently squeeze capillary tube and immerse open end in the center of a blood drop and release the capillary tube slowly to draw up the blood to the 20µl black marking line on the capillary tube.
  - The capillary whole blood must be tested immediately after collection.
- **Venous whole blood)**
  - Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
  - If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 1 – 2 days after collection.
  - Do not freeze venous whole blood specimen.
  - Do not use hemolyzed blood specimen.

- ⚠️ **Anticoagulants** such as heparin, EDTA or sodium citrate do not affect the test result.
- ⚠️ **As known relevant interference**, haemolytic specimen, rheumatoid factors-contained specimen and lipaemic, icteric specimen can lead to impair the test results.

## PREPARATION AND TEST PROCEDURE

- **Preparation**
    - Result window
    - Specimen well
- 

- Carefully read the instructions for using the STANDARD Q HIV 1/2 Ab 3-Line Test.
- Check that the expiry date on the back of the foil pouch has not passed.
- Allow the STANDARD Q HIV 1/2 Ab 3-Line Test components and specimen to come to room temperature (15 – 30°C / 59 – 86°F) for 30min prior to testing.
- Check that the test device package is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from yellow to green), throw away the test device and take another test device packaging.
- Procedure method should be followed for the specific specimen type being tested.

## PERFORMANCE CHARACTERISTICS

### Diagnostic Sensitivity

[Plasma]  
The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 457 positive plasma specimens, is 99.78% (456 / 457) with a Wilson 95% confidence interval of [98.77% - 99.96%].

	STANDARD Q HIV 1/2 Ab 3-Line Test	
	HIV Negative	HIV Positive
Anti-HIV-1 positive	1	254
Anti-HIV-1 positive non-B subtype	0	40
Anti-HIV-2 positive	0	122
<b>Sensitivity</b>	<b>456/457= 99.78%</b>	


[Serum]  
The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 623 serum specimens, is 99.84% (622 / 623) with a Wilson 95% confidence interval of [99.10% - 99.97%].

	STANDARD Q HIV 1/2 Ab 3-Line Test	
	HIV Negative	HIV Positive
Anti-HIV-1 positive	1	589
Anti-HIV-1 positive non-B subtype	0	10
Anti-HIV-2 positive	0	23
<b>Sensitivity</b>	<b>622/623= 99.84%</b>	

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


#### 2. Adding of Specimen

- [For serum/plasma/venous whole blood specimen]
- Add the collected specimen to the specimen well of the test device using the micropipette.

#### [For capillary whole blood specimen]

- Add the collected specimen to the specimen well of the test device. Let the end of the capillary tube to touch the pad, then press down lightly.

#### 3. Dropping of Buffer

- Hold the buffer bottle at 90° angle to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.

#### 4. Reading Time

- Read the test results between 10 to 20 minutes after adding buffer

- ⚠️ Do not read test results after 20 minutes. It may give false results.
- ⚠️ Place the test device on the flat surface before applying the specimen.
- ⚠️ Use correct amount of specimen and buffer.
- ⚠️ If the line does not appear until 10 minutes, make sure to observe the appearance of the line until 20 minutes.

## INTERPRETATION OF TEST RESULT

Test result	Example	Description
Negative		The presence of only "C" line indicates a negative result.
Positive	HIV-1 	The presence of two lines as "C" and "T1" line indicates a positive result for HIV-1. In case of the presence of three lines as "C", "T1" and "T2", if the intensity of the "T1" line is stronger than "T2" line, it should be interpreted as HIV-1 positive.
	HIV-2 	The presence of two lines as "C" and "T2" line indicates a positive result for HIV-2. In case of the presence of three lines as "C", "T2" and "T1", if the intensity of the "T2" line is stronger than "T1" line, it should be interpreted as HIV-2 positive.
HIV-1 & HIV-2 	The presence of three lines as "C" and equivalent intensity of "T2" and "T1", it should be interpreted as HIV-1 and HIV-2 both positive.	
Invalid		No presence of "C" line indicates an invalid result. It is recommended that the specimen be retested using a new device.

- ⚠️ Interpret faint lines of "T1" and/or "T2" as the positive for HIV-1 and HIV-2.
- ⚠️ In case the intensity of "T1" and "T2" line is similar, confirmatory testing is recommended to confirm the virus type or a co-infection.
- ⚠️ The presence of any line no matter how faint it is, should be considered as a line formed.
- ⚠️ Positive results should be considered in conjunction with the clinical history and other data available.
- ⚠️ This test is for screening purposes. Confirmatory testing according to national guidelines is recommended to confirm the infection status.

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## EXPLANATION AND SUMMARY

### Introduction

AIDS is caused by two known types of HIV (human immunodeficiency virus). HIV type 1 and HIV type 2. HIV type 1 (HIV-1) is found in patients with AIDS, AIDS-related complex (ARC), and asymptomatic infected individuals at high risk for AIDS. The virus is transmitted by sexual contact, by exposure to infected blood or blood products, or from an infected mother to her fetus or infant. The infection of HIV type 2 (HIV-2) is endemic only in West Africa, and it has been identified in individuals who had sexual relations with individuals from that geographic region. HIV-2 is similar to, but distinct from, HIV-1. Both viruses have similar morphology and lymphotropism, and the modes of transmission appear to be identical. Within the two major HIV types, there is significant variation. HIV-1 has been divided into four groups: group M (for major)-including at least ten subtypes, group O (for outlier), group P, and group N (for non-M, non-O). Similarly, the HIV-2 strains have been classified into at least five subtypes (A through E). STANDARD Q HIV 1/2 Ab 3-Line Test can detect antibodies to both HIV-1 and HIV-2, which first appear significantly later, 20-45 days after infection. Detecting HIV earlier with STANDARD Q HIV 1/2 Ab 3-Line Test is helping to prevent virus transmission during extremely infectious stage.

### Intended use

STANDARD™ Q HIV 1/2 Ab 3-Line Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O and HIV 2 in human serum, plasma or venous and capillary whole blood. The test is for *in vitro* diagnostic use and intended as an aid to diagnosis of HIV in symptomatic, asymptomatic populations and persons at risk of HIV infection. The test is not intended for use on neonates or infants below 2 years. The test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited countries. The product is not for blood donor screening. The test kit is not automated and does not require any additional instruments. More specific alternative diagnosis methods should be performed in order to ensure the confirmation of HIV Virus infection.

### Test Principle

STANDARD Q HIV 1/2 Ab 3-Line Test has "T1", "T2" and "C" line pre-coated with recombinant HIV-1 GP41 Ig1 / recombinant HIV-1 subtype O GP41, recombinant HIV-2 GP36 protein and monoclonal anti-chicken IgY respectively. The anti-HIV-1/anti-HIV-1 subtype O in patient specimens interacts with the recombinant HIV-1 GP41-gold / recombinant HIV-1 subtype O GP41-gold and the anti-HIV-2 in patient specimens interacts with the recombinant HIV-2 GP36-gold in the conjugation pad. The complex of gold conjugated antigens and antibodies moves along the membrane chromatographically to the membrane with buffer and is captured by the recombinant HIV-1 and HIV-2 antigens on the test regions (T1 and T2). If the antibodies against HIV-1/2 are in the patient specimen, visible lines are formed in the test region. The chicken IgY antibody conjugated with colloidal gold particles are used as detectors for "C" control line. The control line should always appear if the test procedure is performed properly.

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## WARNINGS AND PRECAUTIONS

- Do not reuse the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use after the expiration date.
- Do not use the buffer bottle of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Allow Kit components to reach room temperature (15 – 30°C / 59 – 86°F) before performing a test.
- Follow necessary precautions when handling specimens with this test. Use personal protective equipment (PPE) consistent with current guidelines.
- The buffer contains 0.01% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.

## LIMITATION OF TEST

- All three test lines ("T1", "T2" and "C") may be obtained when tested with specimen containing high titers of HIV-1 antibodies. Hence, reactive test bands for both HIV-1 and HIV-2 may not indicate mixed infection but may result from the cross reactivity of HIV-1 and HIV-2 because of the similarity of their genomic structure.
- A positive result indicates presence of antibodies to HIV-1 and/or HIV-2. However, a positive result does not indicate a conclusive HIV infection diagnosis. A positive result should be confirmed by a supplemental test.
- A negative result does not eliminate the possibility of infection with HIV-1/2. The STANDARD Q HIV 1/2 Ab 3-Line Test may not detect extremely low concentration of the antibody to HIV 1/2. For negative result with this kit, additional test using other clinical method is necessary.
- The test results alone should not be used in diagnosis of infection with HIV-1/2. For overall clinical diagnosis, results must be interpreted with the patient's clinical symptoms, and clinical history.

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
## STANDARD Q HIV 1/2 Ab 3-Line

STANDARD™ Q HIV 1/2 Ab 3-Line Test  
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

### KIT CONTENTS

- **For Cat. No.: 09HIV30D**
- 1 Test device (individually in a foil pouch with desiccant) x 25
  - 2 Capillary tube (20µl) x 25
  - 3 Buffer bottle x 1
  - 4 Sterile lancet x 25
  - 5 Alcohol swab x 25
  - 6 Instructions for use x 1

- **For Cat. No.: 09HIV30DM**
- 1 Test device (individually in a foil pouch with desiccant) x 25
  - 2 Buffer bottle x 1
  - 3 Instructions for use x 1

### MATERIALS REQUIRED BUT NOT PROVIDED

- **For Cat. No.: 09HIV30D**
- 1 Anti-coagulant tube containing heparin, EDTA or sodium citrate
  - 2 Micropipette and tip
  - 3 PPE (Personal Protective Equipment)
  - 4 Biohazard container

- **For Cat. No.: 09HIV30DM**
- 1 Anti-coagulant tube containing heparin, EDTA or sodium citrate
  - 2 Micropipette and tip
  - 3 Sterile lancet
  - 4 Alcohol swab
  - 5 PPE (Personal Protective Equipment)
  - 6 Biohazard container
  - 7 Capillary whole blood collection tool

### SPECIMEN COLLECTION AND PREPARATION

- **Serum**
- 1 Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
  - 2 If serum in the plain tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 3 days after collection. For prolonged storage, it should be at below -40°C / -40°F up to 3 months. It should be brought to room temperature prior to use.
- **Plasma**
- 1 Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
  - 2 If plasma in an anti-coagulant tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 3 days after collection. For prolonged storage, it should be at below -40°C / -40°F up to 3 freeze-thaw cycles.
  - 3 It should be brought to room temperature prior to use.


- **Whole blood**
- **Capillary whole blood**
- 1 Capillary whole blood should be collected aseptically by fingertip.
  - 2 Select the finger that is not calloused and gently rub the finger to warm it to stimulate blood circulation. And then clean the fingertip by wiping with an alcohol swab.
  - 3 Wait until the fingertip is dried completely and squeeze the end of the fingertip and pierce the wiped fingertip with a sterile lancet to bleed.
  - 4 Gently squeeze capillary tube and immerse open end in the center of a blood drop and release the capillary tube slowly to draw up the blood up to the 20µl black marking line on the capillary tube.
  - 5 The capillary whole blood must be tested immediately after collection.

- **Venous whole blood**
- 1 Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
  - 2 If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2 – 8°C / 36 – 46°F, the specimen can be used for testing within 1 – 2 days after collection.
  - 3 Do not freeze venous whole blood specimen.
  - 4 Do not use hemolyzed blood specimen.

- Anticoagulants such as heparin, EDTA or sodium citrate do not affect the test result.
- As known relevant interference, haemolytic specimen, rheumatoid factors-contained specimen and lipemic, icteric specimen can lead to impair the test results.

### PREPARATION AND TEST PROCEDURE

#### ■ Preparation

- 1 Carefully read the instructions for using the STANDARD Q HIV 1/2 Ab 3-Line Test.
- 2 Check that the expiry date on the back of the foil pouch has not passed.
- 3 Allow the STANDARD Q HIV 1/2 Ab 3-Line Test components and specimen to come to room temperature (15 – 30°C / 59 – 86°F) for 30min prior to testing.
- 4 Check that the test device package is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from yellow to green), throw away the test device and take another test device packaging.
- 5 Procedure method should be followed for the specific specimen type being tested.

### PERFORMANCE CHARACTERISTICS

#### ■ Diagnostic Sensitivity

■ **Plasma**  
The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 457 positive plasma specimens, is 99.78% (456 / 457) with a Wilson 95% confidence interval of [98.77% - 99.96%].

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■ **Serum**  
The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 623 serum specimens, is 99.84% (622 / 623) with a Wilson 95% confidence interval of [99.10% - 99.97%].


	STANDARD Q HIV 1/2 Ab 3-Line Test	
	HIV Negative	HIV Positive
Anti-HIV-1 positive	1	583
Anti-HIV-1 positive non-B subtype	0	10
Anti-HIV-2 positive	0	23
<b>Sensitivity</b>	<b>622/623= 99.84%</b>	

### ■ Test Procedure

#### 1. Collecting of Specimen


[For serum/plasma/venous whole blood specimen]

- 1 Collect the 10 µl of serum/plasma or 20 µl of venous whole blood specimen using a micropipette.


#### [For capillary whole blood specimen]


- 1 Select the finger that is not calloused and gently rub the finger to warm it to stimulate blood circulation. And then clean the fingertip by wiping with an alcohol swab.

- 2 Wait until the fingertip is dried completely and pierce the wiped fingertip with a sterile lancet to bleed.


- 3 Gently squeeze capillary tube and immerse open end in the center of a blood drop and release the capillary tube slowly to draw up the blood up to the 20 µl black marking line on the capillary tube.

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
### 4. Reading Time

Read the test results between 10 to 20 minutes after adding buffer.

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### INTERPRETATION OF TEST RESULT

If the T line does not appear until 10 minutes, make sure to observe the appearance of the T line until 20 minutes.

Test result	Example	Description
<b>Negative</b>		The presence of only "C" line indicates a negative result.
<b>HIV-1</b>		- The presence of two lines as "C" and "T1" line indicates a positive result for HIV-1. - In case of the presence of three lines as "C", "T1" and "T2", if the intensity of the "T1" line is stronger than "T2" line, it should be interpreted as HIV-1 positive.
		- The presence of two lines as "C" and "T2" line indicates a positive result for HIV-2. - In case of the presence of three lines as "C", "T2" and "T1", if the intensity of the "T2" line is stronger than "T1" line, it should be interpreted as HIV-2 positive.
<b>Positive</b>		The presence of three lines as "C" and equivalent intensity of "T2" and "T1", it should be interpreted as HIV-1 and HIV-2 both positive.
<b>Invalid</b>		No presence of "C" line indicates an invalid result. It is recommended that the specimen be retested using a new device.

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### QUALITY CONTROL

1. A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent are reactive.
2. Control materials are not supplied with this test kit. However, it is recommended that the positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
3. If there is a problem with the result such as invalid result, retest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

### EXPLANATION AND SUMMARY

#### ■ Introduction

AIDS is caused by two known types of HIV (human immunodeficiency virus), HIV type 1 and HIV type 2. HIV type 1 (HIV-1) is found in patients with AIDS, AIDS-related complex (ARC), and asymptomatic infected individuals at high risk for AIDS. The virus is transmitted by sexual contact, by exposure to infected blood or blood products, or from an infected mother to her fetus or infant. The infection of HIV type 2 (HIV-2) is endemic only in West Africa, and it has been identified in individuals who had sexual relations with individuals from that geographic region. HIV-2 is similar to, but distinct from, HIV-1. Both viruses have similar morphology and lymphotropism, and the modes of transmission appear to be identical. Within the two major HIV types, there is significant variation. HIV-1 has been divided into four groups: group M (for major)-including at least ten subtypes, group O (for outlier), group P, and group N (for non-M, non-O). Similarly, the HIV-2 strains have been classified into at least five subtypes (A through E). STANDARD Q HIV 1/2 Ab 3-Line Test can detect antibodies to both HIV-1 and HIV-2, which first appear significantly later, 20-45 days after infection. Detecting HIV earlier with STANDARD Q HIV 1/2 Ab 3-Line Test is helping to prevent future transmission during extremely infectious stage.

#### ■ Intended use

STANDARD™ Q HIV 1/2 Ab 3-Line Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O and HIV 2 in human serum, plasma or venous and capillary whole blood. The test is for *in vitro* diagnostic use and intended as an aid to diagnosis of HIV in symptomatic, asymptomatic populations and persons at risk of HIV infection. The test is not intended for use on neonates or infants below 2 years. The test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited countries. The product is not for blood donor screening. The test kit is not automated and does not require any additional instruments. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HIV Virus infection.

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### KIT STORAGE AND STABILITY

1. Store the sealed pouch and the buffer provided in the kit at 2 – 40°C / 36 – 104°F out of the direct sunlight for the duration of its shelf life.
2. Do not open the foil pouch until you are ready to perform a test. Test device can be used immediately after being opened.
3. Close the buffer cap tightly after using, and then store it at 2 – 40°C / 36 – 104°F out of the direct sunlight. It is stable until the expiry date of the kit and the buffer label after opening its cap, if it is tightly closed.

### WARNINGS AND PRECAUTIONS

1. Do not re-use the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use after the expiration date.
4. Do not use the buffer bottle of another lot.
5. Do not smoke, drink or eat while handling specimen.
6. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
7. Clean up spills thoroughly using an appropriate disinfectant.
8. Handle all specimens as if they contain infectious agents.
9. Observe established precautions against microbiological hazards throughout testing procedures.
10. Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
11. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
12. Allow kit components to reach room temperature (15 – 30°C / 59 – 86°F) before performing a test.
13. Follow necessary precautions when handling specimens with this test. Use personal protective equipment (PPE) consistent with current guideline.
14. The buffer contains 0.1% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.

### LIMITATION OF TEST

1. All three test lines ("T1", "T2" and "C") may develop when tested with specimen containing high titers of HIV-1 antibodies. Hence, reactive test bands for both HIV-1 and HIV-2 may not indicate mixed infection but may result from the cross reactivity of HIV-1 and HIV-2 because of the similarity of their genomic structure.
2. A positive result indicates presence of antibodies to HIV-1 and/or HIV-2. However, a positive result does not indicate a conclusive HIV infection diagnosis. A positive result should be confirmed by supplemental test.
3. A negative result does not eliminate the possibility of infection with HIV-1/2. The STANDARD Q HIV 1/2 Ab 3-Line Test may not detect extremely low concentration of the antibody to HIV 1/2. For negative result with this kit, additional test using other clinical method is necessary.
4. The test results alone should not be used in diagnosis of infection with HIV-1/2. For overall clinical diagnosis, results must be interpreted with the patient's clinical symptoms, and clinical history.

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### PERFORMANCE CHARACTERISTICS

- Owen SM et al. Alternative Algorithms for Human Immunodeficiency Virus Infection Diagnosis Using Tests That Are Licensed in the United States. J Clin Microbiol 46:1588-1595, 2008.
- Barre-Sinoussi F, Chermann JC, Rey F, et al. Isolation of a T-lymptropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). Science 220:868-871, 1983.
- Centers for Disease Control: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Rep 37:377-388, 1988.
- Hoff R, Weiblen BJ, Schwenzer M, et al. Specific antibodies to HIV-2 detected in an anonymous newborn blood specimen from Massachusetts. Fourth Consensus Conference on Testing for Human Retroviruses, March 1989.
- Charneau P, Borman AM, Quilliant C, et al: Isolation and envelope sequence of a highly divergent.
- Perspectives in Disease Prevention and Health Promotion Update: Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and Other Bloodborne Pathogens in Health-Care Settings. MMWR, June 24, 1988;37(24):377-388.

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