WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: STANDARD Q HIV/Syphilis Combo Test
WHO reference number: PQDx 0382-117-00

STANDARD Q HIV/Syphilis Combo Test with product code **09HIV20D**, manufactured by **SD Biosensor, Inc., Rest-of-World regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 26 May 2020.

Summary of WHO prequalification assessment for STANDARD Q HIV/Syphilis Combo Test

	Date	Outcome
Prequalification listing	26 May 2020	listed
Dossier assessment	1 May 2020	MR
Site inspection(s) of the quality management system	24-26 May 2023	MR
Product performance evaluation	Quarter 4-2019	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendment	Date of report amendment
2.0	1. Addition of a manufacturer for Cutting, Assembly, Buffer Preparation & dispensing, sealing, printing, packing, and shipping process. 2. Addition of suppliers for raw materials: <ul style="list-style-type: none"> • Upper device, • Lower device, • Desiccant, • Aluminum Pouch, • Buffer bottle, • Capillary tube, • Lancet, 	7 March 2024

	<ul style="list-style-type: none"> • Alcohol swab, • Outer package/IFU/Label. <p>3. Changes to the Outer package, Instructions for Use (IFU) and labels.</p>	
--	--	--

Intended use

According to the claim of intended use from SD Biosensor Inc., *“STANDARD Q HIV/Syphilis Combo Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O, HIV 2 and Syphilis (Treponema pallidum) in human serum, plasma or whole blood. The test is for in vitro diagnostic use and intended as an aid to early diagnosis of HIV and Syphilis infection for HIV or Syphilis infected patients, patients with signs and symptoms for HIV and Syphilis and persons at risk. It provides only an initial screening test result. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HIV Virus and Syphilis infection.”*

Assay description

According to the claim of assay description from SD Biosensor Inc, *“STANDARD Q HIV/Syphilis Combo Test has “H1”, “H2”, “SYP” and “C” line region pre-coated with recombinant HIV-1 GP41 protein / recombinant HIV-1 subtype O GP41, recombinant HIV-2 GP36 protein, recombinant p17 Treponema pallidum protein (recombinant TPP 17 protein) and monoclonal anti-HIV-1 / monoclonal anti-syphilis respectively. The anti-HIV-1/anti-HIV-1 subtype O in patient sample interacts with the recombinant HIV-1 GP41-gold / recombinant HIV-1 subtype O GP41-gold and the anti-HIV-2 in patient sample interacts with the recombinant HIV-2 GP36-gold in the conjugation pad. The anti-syphilis in patient sample interacts with the recombinant TPP 17 protein-gold. The complex moves along the membrane chromatographically with assay diluent and is captured by the recombinant HIV-1 and HIV-2 antigens and/or recombinant TPP 17 antigen on the each test line (H1, H2, SYP). If the antibodies against HIV 1/ 2 and/or syphilis are in the patient sample, visible lines are formed in the each test line. The control line should always appear if the test procedure is performed properly.”*

Test kit contents

Component	25 tests (product code 09HIV20D)
Test device (individually in a foil pouch with desiccant)	25
Buffer Bottle	1 x 4 mL
Capillary tube (20µl)	Pack of 25
Instructions for use	1
Sterile lancet	25
Alcohol swabs	25

Items required but not provided

- Micropipette and tip
- Blood collection tube
- PPE (Personal Protective Equipment)
- Biohazard container
- Timing device

Storage

The test kit should be stored at 2-40°C.

Shelf-life upon manufacture

24 months.

Warnings/limitations

Please refer to the IFU attached to this public report.

Prioritization for Prequalification

Based on the established eligibility criteria, STANDARD Q HIV/Syphilis Combo Test was given priority for the WHO prequalification assessment.

Dossier assessment

SD Biosensor Inc. submitted a product dossier for STANDARD Q HIV/Syphilis Combo Test as per the “Instructions for compilation of a product dossier” (PQDx_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The Manufacturer's responses to the non-conformities found during the dossier review were accepted on 1 May 2020.

Commitment for Prequalification

SD Biosensor, Inc. committed to providing the interim study report and raw data for device stability studies on 24 November 2020 and the final report and raw data on 23 March 2022.

Based on the product dossier screening and assessment findings, the product dossier for STANDARD Q HIV/Syphilis Combo Test meets WHO prequalification requirements.

Manufacturing site inspection

An inspection of SD Biosensor Inc. located at 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea was conducted from 24-26 May 2023. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of the product of consistent quality. Routine inspections of the Manufacturer will be conducted with copies of these WHO Public Inspection Reports (WHOPIRs) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current inspection performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection can be found at:

<https://extranet.who.int/prequal/inspection-services/prequalification-reports/whopirs-vitro-diagnostics>

All published WHOPIRs are with the manufacturer's agreement.

Product performance evaluation

STANDARD Q HIV/Syphilis Combo Test was evaluated at the Institute of Tropical Medicine, Belgium, on behalf of WHO in the 4th quarter of 2019, according to protocol PQDx_150, version 4.1.

Clinical performance evaluation

In this limited laboratory-based evaluation of clinical performance characteristics, a panel of 400 serum/plasma specimens was used. The specimens were characterized using the following reference algorithms. For HIV: Vironostika HIV Ag/Ab (bioMérieux) and Enzygnost Anti-HIV 1/2 Plus (Siemens Healthcare Diagnostics) or Genscreen HIV-1/2 Version 2 (Bio-Rad) in parallel, followed by INNO-LIA HIV I/II Score (Fujirebio Inc.) on initially reactive specimens. For *Treponema pallidum*: Vitros Syphilis TPA Assay (Ortho Clinical Diagnostics), followed by SERODIA-TP.PA (Fujirebio Inc.).

Clinical performance characteristics in comparison with an agreed reference standard				
	HIV-1/2		Syphilis	
	Initial (95% CI)	Final (95% CI)	Initial (95% CI)	Final (95% CI)
Sensitivity % (N=200)	100 (98.2-100)	100 (98.2-100)	95.0 (91.0-97.6)	95.5 (91.6-97.9)
Specificity % (N= 200)	99.0 (96.4-99.9)	99.5 (97.2-100)	99.5 (97.2-100)	99.5 (97.2-100)
Invalid rate % (N= 400)	0%			
Inter-reader variability % (N= 400)	3.25%*		1.0%	

* All 13 disagreements on HIV-1/2 results were on the HIV-2 line in HIV-1 positive specimens.

Out of 200 HIV-1 positive specimens, STANDARD Q HIV/Syphilis Combo Test showed the presence of the HIV-2 line in 30 (15 %) specimens, although in most cases (n=28), the HIV-2 line was weaker than the HIV-1 line, which is interpreted as HIV-1 positive result according to the IFU of the assay.

Analytical performance evaluation

Analytical performance characteristics		
	HIV-1/2	Syphilis
Sensitivity during seroconversion in comparison with a benchmark assay (Enzygnost Anti-HIV 1/2 Plus)	Of a total of 52 specimens in 8 panel, 23 were detected by the assay under evaluation, versus 21 specimens detected by the benchmark assay (Enzygnost Anti-HIV 1/2 Plus). Seroconversion sensitivity index of -0.25. Therefore, detection is 0.25 specimens earlier than the benchmark assay.	Of a total of 9 specimens in 1 panel, 5 were detected by the assay under evaluation versus 2 specimens detected by the benchmark assay (Vitro Syphilis TPA Assay).
Analytical sensitivity on mixed titer panels	All 25 specimens of panel PRB-205 (SeraCare) were correctly classified.	All 17 specimens of panel PSS-202 (SeraCare) were correctly classified.
Analytical sensitivity on WHO reference preparation panels	All 6 HIV subtypes/groups in the 1 st International Reference Panel for anti-HIV (NIBSC code 02/210) were detected.	The 1 st International Standard for human syphilitic plasma IgG (NIBSC code 05/122) was detected.

Lot to lot variation on a dilution panel	Lot to lot variation was within +/- 1 two-fold dilutions for all 10 dilution series.	Lot to lot variation was within +/- 1 two-fold dilutions for all 10 dilution series.
--	--	--

Operational characteristics and ease of use

This assay does not require laboratory equipment and can be performed in laboratories with limited facilities or in non-laboratory settings.

The assay was found easy to use by the operators performing the evaluation.

Key operational characteristics	
Number of steps*	2 steps in total 1 step with precision pipetting (only for serum/plasma)
Time to result	15 minutes
Endpoint stability (interval)	5 minutes (the test can be read between 15 and 20 minutes after the addition of diluent)
Internal QC	Yes, reagent addition control

* Definition: each action required to obtain a result (excluding specimen collection, device preparation – opening the pouch), e.g. for RDTs: add specimen, add buffer (2 steps).

Based on these results, the performance evaluation for the STANDARD Q HIV/Syphilis Combo Test meets the WHO prequalification requirements.

Labelling


- 1. Labels**
- 2. Instructions for use**

1. Labels


1.1 Device Package Artwork

STANDARD Q

HIV/Syphilis Combo 25T


STANDARD Q HIV/Syphilis Combo 25T

SD BIOSENSOR

www.sdbiosensor.com

- CONTENTS**
- 1) Test device (individually in a foil pouch with desiccant) x 25
 - 2) Capillary tube (20µl) x 25
 - 3) Buffer bottle x 1
 - 4) Sterile lancet x 25
 - 5) Alcohol swab x 25
 - 6) Instructions for use x 1

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
STANDARD Q HIV/Syphilis Combo Test

STANDARD™

SD Biosensor, Inc.
C-4th&5th, 16, Deogyang-dae-ro 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16696, REPUBLIC OF KOREA

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INTENDED USE
To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE
To use, twist-off the protective cap.

CAUTION
The lancet is guaranteed sterile while protective cap is sealed to the base.
Do not use if the seal has been damaged or broken.

 **Manufactured by**
Tianjin Huahong Technology Co., Ltd.
A01, Plant B No. 278, Hangkong Road,
Tianjin Pilot Free Trade Zone
(Air Port Industrial Park),
300308 Tianjin, China

 **Authorized Representative**
Shanghai International Holding Corp.
GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
Tel:+49-40-2513175

 **0197**    

Disposable Sterile Lancets

Product code : SM.01

LOT No. : 21080495

MFG Date : 2023-09

EXP Date : 2028-08

25PCS

28G

INTENDED USE
To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE
To use, twist-off the protective cap.

CAUTION
The lancet is guaranteed sterile while protective cap is sealed to the base. Do not use if the seal has been damaged or broken.

Issue Date : 2023.09

Manufactured by
Phoenix Innovative Healthcare Manufacturers Pvt. Ltd.
EL-209, Shil Mahape Road, Electronic Zone, MIDC,
TTC Industrial Area, Mahape, Navi Mumbai - 400 710,
Maharashtra, India

Authorized Representative
Med Path Gm6H
Miles-van-der-Rohe-Strasse 8,80807
Munich, Germany

CE 0197 STERILE R

Disposable Sterile Lancets

Product code : 01GL25

LOT No. :

MFG Date :

EXP Date :

25PCS

28G

INTENDED USE
To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE
To use, twist-off the protective cap.

CAUTION
The lancet is guaranteed sterile while protective cap is sealed to the base. Do not use if the seal has been damaged or broken.

B06LCT2ENR2
Issue Date : 2020.02

Manufactured by
Beijing Ruicheng Medical Supplies Co., Ltd.
No. 558 Zhangzikou, Yangsong Town,
Huairou District, 101400 Beijing, China

Authorized Representative
Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague, Netherlands.
Tel: +31645171879, +31626669008

CE 0197 STERILE R

1.6 Alcohol swab labels

Front

Back

PANTONE Reflex Blue U

PANTONE Orange 021 U


2 3 4 5 6 7 8

2 3 4 5 6 7 8

F /
F

A →

2. Instructions for use¹

¹ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.


STANDARD Q HIV/Syphilis Combo

STANDARD™ Q HIV/Syphilis Combo Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

- 1) Micropipette and tip
- 2) Blood collection tube
- 3) PPE (Personal Protective Equipment)
- 4) Biohazard container

SPECIMEN COLLECTION AND STORAGE

Serum

1. Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
2. If serum in the plain tube is stored in a refrigerator at 2–8°C / 36–46°F, the specimen can be used for testing within 4 days after collection. For prolonged storage, it should be at below -40°C / -40°F.
3. It should be brought to room temperature prior to use.

Plasma

1. Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
2. If plasma in an anti-coagulant tube is stored in a refrigerator at 2–8°C / 36–46°F, the specimen can be used for testing within 4 days after collection. For prolonged storage, it should be at below -40°C / -40°F.
3. It should be brought to room temperature prior to use.

Whole blood

(Capillary whole blood)

1. Capillary whole blood should be collected aseptically by fingertip.
2. Select the finger that is free from callus. Gently rub the finger to warm it to stimulate blood circulation.
3. Squeeze the end of the fingertip and pierce with a sterile lancet.
4. Collect the capillary whole blood to the black line of the capillary tube for the testing.
5. The capillary whole blood must be tested immediately after collection.

(Venous whole blood)


1. Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2–8°C / 36–46°F, the specimen can be used for testing within 1–2 days after collection.
3. Do not use hemolyzed blood specimen.

- Anticoagulants such as heparin, EDTA or sodium citrate do not affect the test result.
- As known relevant interference, haemolytic specimens, rheumatoid factors-contained specimens and lipaemic, icteric specimens can lead to impair the test results.

PREPARATION AND TEST PROCEDURE

Preparation


1. Carefully read the instructions for using the STANDARD Q HIV/Syphilis Combo Test.
2. Check the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
3. Open the foil pouch, and check the test device and the color indicator desiccant pack in foil pouch.
4. Allow the STANDARD Q HIV/Syphilis Combo Test components and specimen to come to room temperature (15–20°C / 59–68°F) for 30 min prior to testing.
5. Check that the test device packaging is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from orange to green), throw away the test device and take another test device packaging. If the color of the buffer bottle does not show a change, you can use the test. Throw away the buffer bottle in the non-sharps (non-infectious) disposal container.
6. Procedure method should be followed for the specific specimen type being tested.

Test Procedure


1. Collecting of Specimen

- 1-1. For serum/plasma/venous whole blood specimen
Collect the 10 µl of serum/plasma and 20µl of venous whole blood specimen using a micropipette.

- 1-2. For capillary whole blood specimen

- 1) Clean a fingertip by wiping with an alcohol swab.


PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity


1. HIV Ab detection: The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 637 positive specimens, is 100% (637/637) with a Wilson 95% confidence interval of [99.4% - 100.0%].

	STANDARD Q HIV/Syphilis Combo Test	
	HIV negative	HIV positive
Anti-HIV-1 positive/Anti-Tp negative	247	
Anti-HIV-1/Anti-Tp positive	0	250
Anti-HIV-1 positive non-B subtype	0	40
Anti-HIV-2 positive	0	100
Sensitivity	637/637 = 100%	

- 2) Dry and pierce the wiped fingertip with a lancet to bleed.


- 3) Collect the 20 µl of capillary whole blood to the black line of the capillary tube (20 µl).

2. Adding of Specimen


2-1. For serum/plasma/venous whole blood specimen

Add the collected specimen to the specimen well of the test device.


2-2. For capillary whole blood specimen

Add the collected specimen to the specimen well of the test device.


3. Dropping of Buffer

Hold the buffer bottle at 90° angle to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.

4. Reading Time

Read the test results between 15 to 20 minutes after adding Buffer.

- Do not read test results after 20 minutes. It may give false results.
- Place the test device on the flat surface after applying the specimen.
- Use correct amount of specimen and buffer.
- If the T line does not appear until 15 minutes, make sure to observe the appearance of the T line until 20 minutes.

INTERPRETATION OF TEST RESULT

Test result	Example	Description
Negative		The presence of only "C" line indicates a negative result.
HIV-1 Positive		The presence of two lines as "C" and "H1" line indicates a positive result for HIV-1.
		In case of the presence of three lines as "C", "H1" and "H2", if the intensity of the "H1" line is stronger than "H2" line, it should be interpreted as HIV-1 positive.
HIV-2 Positive		The presence of two lines as "C" and "H2" line indicates a positive result for HIV-2.
		In case of the presence of three lines as "C", "H2" and "H1", if the intensity of the "H2" line is stronger than "H1" line, it should be interpreted as HIV-2 positive.
Positive		The presence of two lines as "C" and "SYP" line indicates a positive result for Syphilis.
		Interpret faint lines of "SYP" as the positive for Syphilis.
HIV-1 & Syphilis Positive		The presence of three lines as "C", "H1" and "SYP" line indicates a positive result for HIV-1 and Syphilis.
HIV-2 & Syphilis Positive		The presence of three lines as "C", "H2" and "SYP" line indicates a positive result for HIV-2 and Syphilis.
Invalid		No presence of "C" line indicates an invalid result even in the presence of other lines. It is recommended that the specimen be retested using a new device.

- In case the intensity of 'H1' and 'H2' line is similar, confirmatory testing is recommended to confirm the virus type or a co-infection.
- Positive results should be considered in conjunction with the clinical history and other data available.
- This test is for screening purposes. Confirmatory testing according to national guidelines is recommended to confirm the infection status.

EXPLANATION AND SUMMARY

Introduction

It is common that co-infection of HIV and syphilis occur in similar patient group since they are both are sexually transmitted infections. Furthermore, they have an influence on each other in many ways. Many studies show that syphilis can increase the transmission of HIV; the presence of genital ulcers can increase shedding of HIV; syphilis also increases HIV-1 viral load and decreases CD4 T-cell count. It is known that the manifestation and symptoms of syphilis are different in HIV-infected patient group, so therapeutics of syphilis should be changed in HIV co-infected patients. For these reasons, there are no doubts about the usefulness of management of HIV and syphilis in tandem. Especially, it is essential to check whether pregnant women have the HIV/syphilis infection due to the risk of mother-to-child transmission (MCT) of HIV and syphilis infection. STANDARD Q HIV/Syphilis Combo Test is able to quickly screen the HIV/syphilis infection and gives the advantage of managing these two diseases at the same time.

Intended use

STANDARD Q HIV/Syphilis Combo Test is a rapid chromatographic immunoassay for the

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BIBLIOGRAPHY

1. Owen SM et al. Alternative Algorithms for Human Immunodeficiency Virus Infection Diagnosis Using Tests That Are Licensed in the United States. J Clin Microbiol 46:1588-1595, 2008.
2. Barre-Sinoussi F, Chermann JC, Rey F, et al. Isolation and envelope sequence of a highly divergent.
3. Centers for Disease Control. Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Rep 37:377-388, 1988.
4. Hoff R, Weiblen BJ, Schwertzer M, et al: Specific antibodies to HIV-2 detected in an anonymous newborn blood specimen from Massachusetts. Fourth Consensus Conference on Testing for Human Retroviruses, March 1989.
5. Chansue J, Borran AM, Quillian C, et al. Isolation and envelope sequence of a highly divergent.
6. Fathi, D, Dupin, N (September-October 2010). "Origins of syphilis and management in the immunocompetent patient: facts and controversies". Clinics in Dermatology 28 (5): 533-8.
7. Miller JN : value and limitation of nontreponemal and treponemal tests in the laboratory diagnosis of syphilis. Clin. Obstet Gynecol 18: 191-203, 1975.
8. Syphilis - CDC Fact Sheet (Detailed). CDC. November 2, 2015. Retrieved 3 February 2016.
9. Alexander, JM; Sheffield, JS; Sanchez, PJ; Mayfield, J; Wendel GD, Jr (January 1999). "Efficacy of treatment for syphilis in pregnancy". Obstetrics and gynecology 93 (1): 5-8.
10. A Gerber et al., Recombinant Treponema pallidum antigens in syphilis serology. Immunobiology, 196(C):535-49, 1996-1997.
11. Perspectives in Disease Prevention and Health Promotion Update: Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and Other Bloodborne Pathogens in Health-Care Settings. MMWR, June 24,1988/37(24):377-388.

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
SYMBOL

REF	Reference number		Caution		Use by	LOT	Batch code		Consult Instructions for Use		Do not re-use
IVD	In vitro Diagnostics		Note		Manufacturer		Date of manufacture		Contains Sufficient for >10 Tests		Do not use if packaging is damaged
	Indicate that you should keep the product dry		To indicate the temperature limitations in which the transport package has to be kept and handled.		Keep away from sunlight		Do not use if packaging is damaged				

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- Micropipette and tip
- Blood collection tube
- PPE (Personal Protective Equipment)
- Biohazard container

SPECIMEN COLLECTION AND STORAGE


- ### ■ Serum
- Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
 - If serum in the plain tube is stored in a refrigerator at 2–8°C / 36–46°F, the specimen can be used for testing within 4 days after collection. For prolonged storage, it should be at below -40°C / -40°F.
 - It should be brought to room temperature prior to use.
- ### ■ Plasma
- Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
 - If plasma in an anti-coagulant tube is stored in a refrigerator at 2–8°C / 36–46°F, the specimen can be used for testing within 4 days after collection. For prolonged storage, it should be at below -40°C / -40°F.
 - It should be brought to room temperature prior to use.
- ### ■ Whole blood (Capillary whole blood)
- Capillary whole blood should be collected aseptically by fingertip.
 - Select the finger that is free from callus. Gently rub the finger to warm it to stimulate blood circulation.
 - Squeeze the end of the fingertip and pierce with a sterile lancet.
 - Collect the capillary whole blood to the black line of the capillary tube for the testing.
 - The capillary whole blood must be tested immediately after collection.

- ### ■ Venous whole blood
- Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture.
 - If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2–8°C / 36–46°F, the specimen can be used for testing within 1–2 days after collection.
 - Do not use hemolyzed blood specimen.

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PREPARATION AND TEST PROCEDURE


- ### ■ Preparation
- Carefully read the instructions for using the STANDARD Q HIV/Syphilis Combo Test.
 - Check the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
 - Open the foil pouch, and check the test device and the color indicator desiccant pack in foil pouch.
 - Allow the STANDARD Q HIV/Syphilis Combo Test components and specimen to come to room temperature (15–30°C / 59–86°F) for 30 min prior to testing.
 - Check that the test device packaging is not damaged. If damaged, discard the test and use another test. If a humidity indicator inside shows saturation (color changed from orange to green), throw away the test device and take another test device packaging. If the color of the buffer bottle does not show a change, you can use the test. Throw away the buffer bottle in the non-sharps (non-infectious) disposal container.
 - Procedure method should be followed for the specific specimen type being tested.

■ Test Procedure


1. Collecting of Specimen

- 1-1. For serum/plasma/venous whole blood specimen
Collect the 10 µl of serum/plasma and 20µl of venous whole blood specimen using a micropipette.

1-2. For capillary whole blood specimen


- Clean a fingertip by wiping with an alcohol swab.


PERFORMANCE CHARACTERISTICS

- ### ■ Diagnostic sensitivity
- HIV Ab detection: The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 637 positive specimens, is 100% (637/637) with a Wilson 95% confidence interval of [99.4% - 100.0%].
 - Syphilis Ab detection: The Diagnostic Sensitivity for anti-Treponema pallidum antibody detection, calculated on 400 positive specimens, is 98.8% (395/400) with a Wilson 95% confidence interval of [97.1% - 99.5%].

	STANDARD Q HIV/Syphilis Combo Test	
	HIV negative	HIV positive
Anti-HIV-1 positive/Anti-Tp negative	247	0
Anti-HIV-1/Anti-Tp positive	0	250
Anti-HIV-1 positive non-B subtype	0	40
Anti-HIV-2 positive	0	100
Sensitivity	637/637 = 100%	


- Dry and pierce the wiped fingertip with a lancet to bleed.


2. Adding of Specimen

- 2-1. For serum/plasma/venous whole blood specimen
Add the collected specimen to the specimen well of the test device.

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
- Add the collected specimen to the specimen well of the test device.


3. Dropping of Buffer

- Hold the buffer bottle at 90° angle to the test device without touching the specimen well to avoid contamination of the buffer.
Add 3 drops of the buffer into the specimen well of the test device.

4. Reading Time

- Read the test results between 15 to 20 minutes after adding Buffer.

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INTERPRETATION OF TEST RESULT

Test result	Example	Description
Negative		The presence of only "C" line indicates a negative result.
HIV-1 Positive		The presence of two lines as "C" and "H1" line indicates a positive result for HIV-1. In case of the presence of three lines as "C", "H1" and "H2", if the intensity of the "H1" line is stronger than "H2" line, it should be interpreted as HIV-1 positive.
		The presence of two lines as "C" and "H2" line indicates a positive result for HIV-2. In case of the presence of three lines as "C", "H1" and "H2", if the intensity of the "H2" line is stronger than "H1" line, it should be interpreted as HIV-2 positive. Interpret faint lines of "H1" and/or "H2" as the positive for HIV-1 and/or HIV-2.
Positive		The presence of two lines as "C" and "SYP" line indicates a positive result for Syphilis. Interpret faint lines of "SYP" as the positive for Syphilis.
		The presence of three lines as "C", "H1" and "SYP" line indicates a positive result for HIV-1 and Syphilis.
	The presence of three lines as "C", "H2" and "SYP" line indicates a positive result for HIV-2 and Syphilis.	
Invalid		No presence of "C" line indicates an invalid result even in the presence of other lines. It is recommended that the specimen be retested using a new device.
		No presence of "C" line indicates an invalid result even in the presence of other lines. It is recommended that the specimen be retested using a new device.

- In case the intensity of "H1" and "H2" line is similar, confirmatory testing is recommended to confirm the virus type or a co-infection.
- Positive results should be considered in conjunction with the clinical history and other data available.
- This test is for screening purposes. Confirmatory testing according to national guidelines is recommended to confirm the infection status.

EXPLANATION AND SUMMARY

■ Introduction

It is common that co-infection of HIV and syphilis occur in similar patient group since they are both are sexually transmitted infections. Furthermore, they have an influence on each other in many ways. Many studies show that syphilis can increase the transmission of HIV, the presence of genital ulcers can increase shedding of HIV; syphilis also increases HIV-1 viral load and decreases CD4⁺ cell count. It is known that the manifestation and symptoms of syphilis are different in HIV-infected patient group, so therapeutics of syphilis should be changed in HIV co-infected patients. For these reasons, there are no doubts about the usefulness of management of HIV and syphilis in tandem. Especially, it is essential to check whether pregnant women have the HIV/syphilis infection due to the risk of mother-to-child transmission (MCT) of HIV and syphilis infection. STANDARD Q HIV/Syphilis Combo Test is able to quickly screen the HIV/syphilis infection and gives the advantage of managing these two diseases at the same time.

■ Intended use

STANDARD Q HIV/Syphilis Combo Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O, HIV-2 and Syphilis (*Treponema pallidum*) in human serum, plasma or whole blood. The test is for in vitro diagnostic use and intended as an aid to diagnosis of HIV and Syphilis infection for HIV or Syphilis infected patients, patients with signs and symptoms for HIV and Syphilis and persons at risk including pregnant women. The test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited countries. It provides only an initial screening test result. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of HIV Virus and Syphilis infection.

■ Test principle

STANDARD Q HIV/Syphilis Combo Test has "H1", "H2", "SYP" and "C" line region pre-coated with recombinant HIV-1 GP41 protein / recombinant HIV-1 subtype O GP41, recombinant HIV-2 GP36 protein, recombinant p17 *Treponema pallidum* protein (recombinant TPP 17 protein) and monoclonal anti-HIV-1 / monoclonal anti-syphilis respectively. The anti-HIV-1/anti-HIV-1 subtype O in patient specimen interacts with the recombinant HIV-1 GP41-gold / recombinant HIV-1 subtype O GP41-gold and the anti-HIV-2 in patient specimen interacts with the recombinant HIV-2 GP36-gold in the conjugation pad. The anti-syphilis in patient specimen interacts with the recombinant TPP 17 protein-gold. The complex moves along the membrane chromatographically with buffer and is captured by the recombinant HIV-1 and HIV-2 antigens and/or recombinant TPP 17 antigen on the each test line (H1, H2, SYP). If the antibodies against HIV 1/ 2 and/or syphilis are in the patient specimen, visible lines are formed in the each test line. The control line should always appear if the test procedure is performed properly.

KIT STORAGE AND STABILITY

- Store the kit unopened at 2–40°C / 36–104°F, out of direct sunlight.
- Do not open the aluminum pouch until you are ready to perform a test.
- The buffer bottle provided in the kit is stable until the expiry date of the kit after opening its cap, if it is tightly closed.

WARNINGS AND PRECAUTIONS

- Do not reuse the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use after the expiration date.
- Do not use buffer bottle of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly when afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Allow Kit components to reach room temperature (15–30°C) before performing a test.
- Follow necessary precautions when handling specimens with this test. Use personal protective equipment (PPE) consistent with guideline.
- The buffer contains 0.01% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.
- For in vitro diagnostic use only.

LIMITATION OF TEST

- All three test lines ("H1", "H2" and "C") may develop when tested with specimen containing high titers of HIV-1 antibodies. Hence, reactive test bands for both HIV-1 and HIV-2 may not indicate mixed infection but may result from the cross reactivity of HIV-1 and HIV-2 because of the similarity of their genomic structure.
- A positive result indicates presence of antibodies to HIV-1/2 and/or Syphilis. However, a positive result does not indicate a conclusive HIV and Syphilis infection diagnosis. A positive result should be confirmed by a supplemental test.
- A negative result does not eliminate the possibility of infection with HIV-1/2 and Syphilis. The STANDARD Q HIV/Syphilis combo Test may not detect extremely low concentration of the antibody to HIV-1/2 and Syphilis. For negative result with this kit additional test using other clinical method is necessary.
- The test results alone should not be used in diagnosis of infection with HIV-1/2 and Syphilis. For overall clinical diagnosis, results must be used in conjunction with the patient's clinical symptoms, and clinical history.

QUALITY CONTROL

- A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent are reactive.
- Control materials are not supplied with this test kit. However, it is recommended that the positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

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BIBLIOGRAPHY

- Dwen SM et al. Alternative Algorithms for Human Immunodeficiency Virus Infection Diagnosis Using Tests That Are Licensed in the United States. J Clin Microbiol 46:1588-1595, 2008.
- Barre-Sinoussi F, Chermann JC, Rey F, et al. Isolation of a T-lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). Science 220:868-871, 1983.
- Centers for Disease Control. Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Rep 37:377-388, 1988.
- Hoff R, Weiblen BJ, Schwertler M, et al: Specific antibodies to HIV-2 detected in an anonymous newborn blood specimen from Massachusetts. Fourth Consensus Conference on Testing for Human Retroviruses, March 1989.
- Chamone D, Borman AM, Quilley C, et al: Isolation and envelope sequence of a highly divergent.
- Farhi, D, Dupin, N (September-October 2010). "Origins of syphilis and management in the immunocompetent patient: facts and controversies". Clinics in Dermatology 28 (5): 533-8.
- Miller JN : value and limitation of nontreponemal and treponemal tests in the laboratory diagnosis of syphilis. Clin. Obstet Gynecol 18: 191-203, 1975.
- Syphilis - CDC Fact Sheet (Detailed). CDC. November 2, 2015. Retrieved 3 February 2016.
- Alexander, JM; Sheffield, JS; Sanchez, PJ; Mayfield, J; Wendel GD, Jr (January 1999). "Efficacy of treatment for syphilis in pregnancy". Obstetrics and gynecology 93 (1): 5-8.
- A Gerber et al., Recombinant Treponema pallidum antigens in syphilis serology. Immunobiology, 196(5):535-49, 1996-1997.
- Perspectives in Disease Prevention and Health Promotion Update: Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and Other Bloodborne Pathogens in Health-Care Settings. MMWR, June 24, 1988/37(24):377-388.

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SYMBOL

REF	Reference number	Caution	Use by	LOT	Batch code	Consult Instructions for Use	Do not re-use
IVD	In vitro Diagnostics	Note	Manufacturer		Date of manufacture	Contains Sufficient for >10 Tests	
	Indicate that you should keep the product dry	To indicate the temperature limitations in which the transport package has to be kept and handled.	Keep away from sunlight		Do not use if packaging is damaged		

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

SPECIMEN COLLECTION AND STORAGE

- Serum: Collect whole blood by venipuncture into commercially available tubes WITHOUT anti-coagulant, and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
Plasma: Collect the venous whole blood into the commercially available anti-coagulant tube such as heparin, EDTA or sodium citrate by venipuncture and centrifuge blood to get plasma specimen.
Whole blood (Capillary whole blood): Capillary whole blood should be collected aseptically by fingerstick.


2. Dry and pierce the wiped fingertip with a lancet to bleed.


3. Collect the 20 µl of capillary whole blood to the black line of the capillary tube (20 µl).

2. Adding of Specimen
2-1. For serum/plasma/venous whole blood specimen
Add the collected specimen to the specimen well of the test device.


2-2. For capillary whole blood specimen
Add the collected specimen to the specimen well of the test device.

3. Dropping of Buffer
Hold the buffer bottle at 90° angle to the test device without touching the specimen well to avoid contamination of the buffer. Add 3 drops of the buffer into the specimen well of the test device.


4. Reading Time
Read the test results between 15 to 20 minutes after adding Buffer.

INTERPRETATION OF TEST RESULT

Table with 3 columns: Test result, Example, and Description. It details various outcomes such as HIV-1 Positive, HIV-2 Positive, Syphilis Positive, and Invalid results based on the presence of lines C, H1, H2, and SYP.

1-2. For capillary whole blood specimen
1. Clean a fingertip by wiping with an alcohol swab.

EXPLANATION AND SUMMARY

Introduction
It is common that co-infection of HIV and syphilis occur in similar patient group since they are both sexually transmitted infections. Furthermore, they have an influence on each other in many ways. Many studies show that syphilis can increase the transmission of HIV: the presence of genital ulcers can increase shedding of HIV; syphilis also increases HIV-1 viral load and decreases CD4 T-cell count. It is known that the manifestation and symptoms of syphilis are different in HIV-infected patient group, so therapeutics of syphilis should be changed in HIV co-infected patients.

Intended use
STANDARD Q HIV/Syphilis Combo Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O, HIV-2 and Syphilis (Treponema pallidum) in human serum, plasma or whole blood. The test is for in vitro diagnostic use and intended as an aid to diagnosis of HIV and Syphilis infection for HIV or Syphilis infected patients, patients with signs and symptoms for HIV and Syphilis and persons at risk including pregnant women.

Test principle
STANDARD Q HIV/Syphilis Combo Test has "H1", "H2", "SYP" and "C" line region pre-coated with recombinant HIV-1 GP41 protein / recombinant HIV-1 subtype O GP41, recombinant HIV-2 GP36 protein, recombinant p17 Treponema pallidum protein (recombinant TPP 17 protein) and monoclonal anti-HIV-1 / monoclonal anti-syphilis respectively. The anti-HIV-1 (HIV-1 subtype O) in patient specimen interacts with the recombinant HIV-1 GP41-gold / recombinant HIV-1 subtype O GP41-gold and the anti-HIV-2 in patient specimen interacts with the recombinant HIV-2 GP36-gold in the conjugation pad.

KIT STORAGE AND STABILITY

- 1. Store the kit unopened at 2 - 40°C / 36 - 104°F, out of direct sunlight. Do not reuse the test kit.
2. Do not open the aluminum pouch until you are ready to perform a test.
3. The buffer bottle provided in the kit is stable until the expiry date of the kit after opening its cap, if it is tightly closed.

WARNINGS AND PRECAUTIONS

- 1. Do not reuse the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use after the expiration date.
4. Do not use buffer bottle of another lot.
5. Do not smoke, drink or eat while handling specimen.
6. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly when afterwashes.
7. Clean up spills thoroughly using an appropriate disinfectant.
8. Handle all specimens as if they contain infectious agents.
9. Observe established precautions against microbiological hazards throughout testing procedures.
10. Dispose of all specimens and materials used to perform the test as bio-hazard waste.
11. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products.
12. Allow Kit components to reach room temperature (15 - 30°C) before performing a test.
13. Follow necessary precautions when handling specimens with this test. Use personal protective equipment (PPE) consistent with guideline.

LIMITATION OF TEST

- 1. All three test lines ("H1", "H2" and "C") may develop when tested with specimen containing high titers of HIV-1 antibodies. Hence, reactive test bands for both HIV-1 and HIV-2 may not indicate mixed infection but may result from the cross reactivity of HIV-1 and HIV-2 because of the similarity of their genomic structure.
2. A positive result indicates presence of antibodies to HIV-1/2 and/or Syphilis. However, a positive result does not indicate a conclusive HIV and Syphilis infection diagnosis. A positive result should be confirmed by a supplemental test.
3. A negative result does not eliminate the possibility of infection with HIV-1/2 and Syphilis. The STANDARD Q HIV/Syphilis combo Test may not detect extremely low concentration of the antibody to HIV-1/2 and Syphilis. For negative result with this kit additional test using other clinical method is necessary.
4. The test results alone should not be used in diagnosis of infection with HIV-1/2 and Syphilis. For overall clinical diagnosis, results must be used in conjunction with the patient's clinical symptoms, and general history.

QUALITY CONTROL

- 1. A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent are reactive.
2. Control materials are not supplied with this test kit. However, it is recommended that the positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

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
Table with 3 columns: Test result, Example, and Description. It details various outcomes such as HIV-1 Positivo, HIV-2 Positivo, Syphilis Positivo, and Invalid results based on the presence of lines C, H1, H2, and SYP.

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PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity
1. HIV Ab detection: The Diagnostic Sensitivity for anti-HIV antibody detection, calculated on 637 positive specimens, is 100% (637/637) with a Wilson 95% confidence interval of [99.4% - 100.0%].

Table showing diagnostic sensitivity for HIV negative and HIV positive specimens. Results show 0% false positive and 100% sensitivity.

Table showing diagnostic sensitivity for Syphilis negative and Syphilis positive specimens. Results show 100% sensitivity.

Diagnostic specificity:
The Diagnostic Specificity for anti-HIV Ab detection, calculated on 1500 negative specimens, is 99.9% (1499/1500) with a Wilson 95% confidence interval of [99.6% - 100.0%].

Table showing diagnostic specificity for Blood donors EDTA-K3 plasma and EDTA-K3 whole blood. Results show 99.9% specificity.

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BIBLIOGRAPHY

1. Owen SM et al. Alternative Algorithms for Human Immunodeficiency Virus Infection Diagnosis Using Tests That Are Licensed in the United States. J Clin Microbiol 46:1588-1595, 2008.
2. Barre-Sinoussi F, Chermann JC, Rey F, et al. Isolation of a T lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). Science 220:868-871, 1983.
3. Centers for Disease Control. Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Rep 37:377-388, 1988.
4. Hoff R, Weiblen BJ, Schwertler M, et al. Specific antibodies to HIV-2 detected in an anonymous newborn blood specimen from Massachusetts. Fourth Consensus Conference on Testing for Human Retroviruses, March 1989.
5. Charnous P, Berman AM, Quillan C, et al. Isolation and envelope sequence of a highly divergent.
6. Farhi, D, Dupin, N (September-October 2010). "Origins of syphilis and management in the immunocompetent patient: facts and controversies". Clinics in Dermatology 28 (5): 533-8.
7. Miller JN : value and limitation of nontreponemal and treponemal tests in the laboratory diagnosis of syphilis. Clin. Obstet Gynecol 18: 191-203, 1975.
8. Syphilis - CDC Fact Sheet (Detailed)". CDC. November 2, 2015. Retrieved 3 February 2016.
9. Alexander, JM; Sheffield, JS; Sanchez, PJ; Mayfield, J; Wendel GD, Jr (January 1999). "Efficacy of treatment for syphilis in pregnancy". Obstetrics and gynecology 93 (1): 5-8.
10. A Cerber et al., Recombinant Treponema pallidum antigens in syphilis serology. Immunology, 1963:525-69, 1996-1991.
11. Perspectives in Disease Prevention and Health Promotion Update: Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and Other Bloodborne Pathogens in Health-Care Settings. MMWR, June 24,1988(37)(4):377-388.

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