WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

Product: STANDARD Q Malaria P.f Ag Test
WHO reference number: PQDx 0346-117-00

STANDARD Q Malaria P.f Ag Test with product code 09MAL10D, manufactured by SD Biosensor, Inc, CE-Mark regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 6 March 2020.

Summary of the WHO Prequalification Assessment for the STANDARD Q Malaria P.f Ag Test

	Date	Outcome
Prequalification listing	6 March 2020	listed
Dossier assessment	29 October 2019	MR
Product performance evaluation	2018	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendment	Date of report amendment
2.0	<ol style="list-style-type: none"> 1. Change of IFU from IFU in English to multilanguage IFU (EN, ES, FR, PT), 2. Change of outer package, 3. Added supplier for sterile lancet. 	17 May 2021
3.0	<ol style="list-style-type: none"> 1. Addition of Manufacturer for cutting, assembly, sealing, buffer preparation & dispensing, packing, and shipping process. 2. Addition of suppliers for raw materials: <ul style="list-style-type: none"> • Upper device, • Lower device, • Desiccant, • Aluminum Pouch, 	7 March 2024

	<ul style="list-style-type: none"> • Buffer bottle, • Inverted cup, • Lancet, • Alcohol swab, • Outer package/IFU/Label. <p>3. Changes to the Outer package, Instructions for Use (IFU) and Labels.</p>	
4.0	The addition of the manufacturing site is located in the public report of Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India.	19 April 2024.
5.0	Updated labels of inverted cup, sterile lancet and alcohol swab (PQC-IVD-2024-0028).	4 August 2025
6.0	<p>1. Package revision for M mark addition Updated the PKG label, adding the “M” mark for better identification.</p> <p>2. New Package labels for Print area addition Established the PKG label, adding the print area to deal with the requirements for the procurement.</p> <p>3. Test device pouch revision Established the test device pouch label to the exclusive product for the better identification.</p> <p>4. RMS Lancet pouch revision Update of Manufacturer’s Address, EU Representative symbol. Addition of Korean language and KGMP symbol to fulfill the requirements for domestic distribution (PQC-IVD-2025-0135).</p>	9 April 2026

Intended use

According to the claim of intended use from SD Biosensor, Inc, “*STANDARD Q Malaria P.f Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of Plasmodium falciparum (P. falciparum) specific Histidine Rich Protein 2 (HRP-2) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f Ag Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.*”

Test kit contents

Component	25 tests (product code 09MAL10D)
Test device (individually in a foil pouch with desiccant)	25
Buffer Bottle	1 x 4 mL
Inverted cup (5µL)	Pack of 25
Sterile lancet	Pack of 25
Alcohol swab	Pack of 25
Instructions for use	1

Items required but not provided

- Anti-coagulant tube containing heparin, EDTA or sodium citrate for the collection of venous whole blood
- Micropipette and tip
- Timer
- PPE (Personal Protective Equipment)
- Pen/pencil
- Extra lancets and alcohol swabs
- Biosafety sharps container
- Biohazard container
- Sterile gauze

Storage

The test kit must be stored between 2 and 40 °C.

Shelf-life upon manufacture¹

24 months.

Dossier assessment

SD Biosensor, Inc submitted a product dossier for STANDARD Q Malaria P.f Ag Test as per the “*Instructions for compilation of a product dossier*” (PQDx_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 29 October 2019.

Based on the product dossier screening and assessment findings, the product dossier for STANDARD Q Malaria P.f Ag Test meets WHO prequalification requirements.

Manufacturing site inspection

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pgweb/vitro-diagnostics/who-public-inspection-reports>

Based on the site inspection and corrective action plan review, the quality management system for the STANDARD Q Malaria P.f Ag Test meets WHO prequalification requirements.

Product performance evaluation

STANDARD Q Malaria P.f Ag Test² was included in the eighth round of WHO product testing of RDTs for malaria antigen detection, which was completed in 2018. The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild-type parasite panel, *P. vivax* wild-type parasite panel and a negative panel. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

Based on the demonstrated *P. falciparum* panel detection score (87.0% at 200 parasites/ μ L), false-positive rates (0% for clean negatives, 0% for *P. vivax* at 200 parasites/ μ L, 0% for *P. vivax* at 2000 / μ L) and invalid rate (0%), STANDARD Q Malaria P.f Ag Test meets the current laboratory evaluation requirements for Prequalification.

² The product code of the product included in the WHO product testing of malaria RDTs was 09MAL10B, which corresponded to the same product, with a different kit configuration (i.e. number of tests) <https://www.who.int/malaria/publications/atoz/9789241514965/en/>

Summary performance characteristics	Panel detection score		False positive rate (%)			Invalid rate (%)
	200 parasites/ μ L		200 parasites/ μ L		Clean negatives	
	Pf	Pv	Pf	Pv		
STANDARD Q Malaria P.f Ag Test	87.0	NA	NA	0	0	0

Based on these results, STANDARD Q Malaria P.f Ag Test meets the current performance evaluation requirements for Prequalification.

Labelling review

The labelling submitted for the STANDARD Q Malaria P.f Ag Test was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

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

SD BIOSENSOR

STANDARD Q
Malaria P.f Ag Test


Manufactured by
SD Biosensor, Inc.
Head office
C-Winkoh, 16, Deogyong-daru 1550boon-gil, Yeongtong-gu, Seongsi, Gyeonggi-do, 16696, REPUBLIC OF KOREA
Manufacturing site
74, Osongsangmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA
Authorized Representative
MT Promed Consulting GmbH
Ernst-Heckel-Strasse 7 66386 St. Ingbert Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

STANDARD

STANDARD Q Malaria P.f Ag 25T_09MAL10D_인도제조소 O

Unit : mm


Classification	Package	Color	3 Color (K, Pantone 2415C, 1495C)
Document number	INOB25MAL1ML4R1-WHO	Post-processing	UV Coating + Uncoated Area / 3-side sealed
Size	W155 * D124 * H71	Issue date	2025.10.27
Material/gram	GC2 / 270g	Department	Product Planning Division

STANDARD Q Malaria P.f Ag 25T_09MAL10D_인도제조소 X


Unit : mm

Classification	Package	Color	3 Color (K, Pantone 2415C, 1495C)
Document number	INXB25MAL1ML4R1-WHO	Post-processing	UV Coating + Uncoated Area / 3-side sealed
Size	W155 * D124 * H71	Issue date	2025.10.27
Material/gram	GC2 / 270g	Department	Product Planning Division

SD BIOSENSOR

www.sdbiosensor.com

- CONTENTS
- 1) Test device (individually in a foil pouch with desiccant) x 25
 - 2) Inverted cup (5 µl) x 25
 - 3) Buffer bottle x 1
 - 4) Sterile lancet x 25
 - 5) Alcohol swab x 25
 - 6) Instructions for use x 1

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SD BIOSENSOR

STANDARD Q
Malaria P.f Ag Test


SD Biosensor, Inc.
C-4thSt, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Authorized Representative
MT Promed Consulting GmbH
Ernst-Heckel-Strasse 7 66386 St. Ingbert Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

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
STANDARD Q Malaria P.f Ag 25T_09MAL10D_인도제조소 O


Unit : mm

Classification	Package	Color	3 Color (K, Pantone 2415C, 1495C)
Document number	PINOB25MAL1ML4R0-WHO	Post-processing	UV Coating + Uncoated Area / 3-side sealed
Size	W155 * D124 * H71	Issue date	2025.10.27
Material/gram	GC2 / 270g	Department	Product Planning Division

SD BIOSENSOR

www.sdbiosensor.com

- CONTENTS**
 1) Test device (individually in a foil pouch with desiccant) x 25
 2) Inverted cup (5 µl) x 25
 3) Buffer bottle x 1
 4) Sterile lancet x 25
 5) Alcohol swab x 25
 6) Instructions for use x 1

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STANDARD Q
Malaria P.f Ag Test

STANDARD


SD Biosensor, Inc.
 16, Deogyong-daero 150beon-gil, Yeongtong-gu, Seowon-si, Gyeonggi-do, 16694, REPUBLIC OF KOREA
Manufacturing site
 Plot No.38, Sector 4, IMT Manesar, Gurgaon, Haryana 122052, INDIA
 Manufacturing License No.: MFG/MD/2018/00005
Authorized Representative
MT Promed Consulting GmbH
 Ernst-Heckel-Strasse 7 66386 St. Ingbert Germany
 Phone : +49 6894 581020, Fax : +49 6894 581021

Uncoated Area

STANDARD Q Malaria P.f Ag 25T_09MAL10D_인도제조소 X

Unit : mm



Classification	Package	Color	3 Color (K, Pantone 2415C, 1495C)
Document number	PINXB25MAL1ML4R0-WHO	Post-processing	UV Coating + Uncoated Area / 3-side sealed
Size	W155 * D124 * H71	Issue date	2025.10.27
Material/gram	GC2 / 270g	Department	Product Planning Division

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


25 PCS
28G

INTENDED USE

To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base.

Do not use if the seal has been damaged or broken.

B05SL2ENR0

Issue Date : 2020.07

Manufactured by

Tianjin Huahong Technology Co., Ltd.
A01, Plant B No. 278, Hangkong Road,
Tianjin Pilot Free Trade Zone
(Air Port Industrial Park),
300308 Tianjin, China

Authorized Representative

Shanghai International Holding Corp.
GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
Tel:+49-40-2513175

STERILE R

Disposable Sterile Lancets_Pouch_28G_25pcs

Unit : mm

Classification	Label	Color	1 color (K)
Document number	B06LCT2ENR4	Post-processing	
Size	W90 * H110	Issue date	2025.06.24
Material/gram	Printing on the pouch	Department	Product Planning Division

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To use, twist-off the protective cap.

CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the based.
Do not use if the seal has been damaged or broken.

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 Manufactured by
Beijing Ruicheng Medical Supplies Co., Ltd.
Building 5, No.8 Yanqi West Road, Yanqi
Economic Development Zone, Huairou
District 101407 Beijing P.R. China

 Authorized Representative
Lotus Nl B.V.
Koningin Julianaplein 10, 1e Verd,
2595AK The Hague, Netherlands.
Tel: +31644168999
Email: info@lotusnl.com

B06LCT2ENR4
Issue date: 2025.06

SURFACE JOB

110MM

220MM

110MM

Disposable Sterile Lancets

Product code :SM.01

25PCS

28G

 LOT No. : XXXXXXXXX

 MFG Date : YYYY-MM

 EXP Date : YYYY-MM

INTENDED USE

To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base. Do not use if the seal has been damaged or broken.

Issue Date : 2023.09

 Manufactured by

Phoenix Innovative Healthcare Manufacturers Pvt. Ltd.
E1-209, Shil Mahape Road, Electronic Zone, MIDC,
TTC Industrial Area, Mahape, Navi Mumbai - 400 710,
Maharashtra, India

 Authorized Representative

Med Path GmbH
Miles-van-der-Rohe-Strasse 8 80907
Munich, Germany

CE 0197

 STERILE

25 Pcs.

100MM

size: 95mm x 100mm No.256

95mm

Blood Lancet

 LOT No. : 240311

Product code : BL-28G

 MFG Date : 2024-03-11

25 PCS

 EXP Date : 2029-02-11

28G

INTENDED USE

To obtain the capillary blood specimen for test.

INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

IFU-02-BL-01-NO.256

Issue Date: 2023-09-15

CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base.

Do not use if the seal has been damaged or broken.

Manufactured by

Guangzhou iCare Medical Technology Co., Ltd.
First floor A No.8, Lianhua Port Industrial Zone,
Lotus Mountain Bonded Area, Shilou Town,
Panyu District , Guangzhou,
People's Republic of China.

Tel: +86-020-31169913 Fax: +86-020-84881433

Authorized Representative

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175 Fax: +49-40-2513175

100mm

0123

STERILE

Instructions for use³

³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

STANDARD Q Malaria P.f Ag

STANDARD™ Q Malaria Pf Ag Test


DE, IT

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

- 1) Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood
- 2) Micropipette and tip
- 3) Timer
- 4) PPE (Personal Protective Equipment)
- 5) Pen/pencil
- 6) Extra lancets and alcohol swabs
- 7) Biosafety sharps container
- 8) Biohazard container
- 9) Sterile gauze

SPECIMEN COLLECTION AND PREPARATION

Capillary whole blood

1. Capillary whole blood should be collected aseptically by fingertick.
2. Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
3. Squeeze the end of the fingertip and pierce with a sterile lancet.
4. Wipe the first drop with a sterile gauze.
5. Collect the capillary whole blood and it must be tested immediately after collection.

Venous whole blood

1. Collect the venous whole blood into the commercially available anticoagulant tube such as heparin, EDTA, or sodium citrate by venipuncture.
2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/ 36-46°F, the specimen can be used for testing within 72 hours after collection.
3. Do not use hemolyzed blood specimen.


CAUTION

Anticoagulants such as heparin, EDTA, or sodium citrate do not affect the test result. As known relevant interference, hemolytic specimen, rheumatoid factors-contained specimen and lipaemic, icteric specimen may lead to inaccurate results. Use separate inverted cup for each specimen in order to avoid cross-contamination of either specimens which can cause erroneous results.

PREPARATION AND TEST PROCEDURE

Preparation

1. Carefully read the instructions for using the STANDARD Q Malaria Pf Ag Test.
2. Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
3. Open the foil pouch, and check the test device and the desiccant inside the foil pouch.

CAUTION

Bring the kit contents and the specimens to room temperature before testing.

Test Procedure

Capillary and venous whole blood can be used for the STANDARD Q Malaria Pf Ag Test.

1. Test with a sterile lancet

1. Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
-
2. Dry and pierce the wiped fingertip with a sterile lancet to bleed.
-

PERFORMANCE CHARACTERISTICS

Limit of Detection (LOD)

	P.f Line
	<i>Plasmodium falciparum</i>
Limit of detection (parasites/µL)	183

Potential Interfering Substances

STANDARD Q Malaria Pf Ag Test was evaluated with the following interfering substances present in specimen in order to assess their potential effect on the assay performance as per CLSI guideline EP7-A2. It was not possible to evaluate the effect of the following substances: artesunate-lumefantrine (malaria drug), doxycycline hydrate (malaria drug), lamivudine (retroviral medication), acetaminophen, ibuprofen, acetylsalicylic acid, bilirubin (uncjugated), caffeine, ethanol. There was no interference with Human anti-mouse antibody, whole blood of pregnant women, infant and neonate specimen and whole blood having elevated levels of hemoglobin, elevated levels of C-reactive protein and lipidemic specimen.

EXPLANATION AND SUMMARY

Introduction

Malaria remains an important cause of illness and death in children and adults in countries in which it is endemic. Globally, an estimated 3.2 billion people in 97 countries and territories are at risk of being infected with malaria and developing disease, and 1.2 billion are at high risk (>1 in 1,000 chance of getting malaria in a year). According to the World Malaria Report 2015, there were 210 million cases of malaria globally in 2015 (uncertainty range 149-303 million) and 438,000 malaria deaths (range 236,000-635,000), representing a decrease in malaria cases and deaths of 37% and 60% since 2000, respectively. The protozoal parasites that cause malaria are from *Plasmodium falciparum*, *vivax*, *ovale* and *malariae* with the first two species causing the most infections worldwide. Classic symptoms of malaria include fever, headache, chills, vomiting, shivering and convulsions. In some rare forms of *P. falciparum*, the patient may present with delirium or coma. Several anemia is often attributed to the cause of death from malaria. Accurate and prompt diagnosis of malaria is of utmost importance due to the morbidity associated with the other malaria forms. Rapid diagnostic test is an ideal diagnostic tool for malaria diagnosis in that it can provide a rapid determination if the patient is infected with malaria allowing for accurate treatment and improved outcomes. STANDARD Q Malaria Pf Ag Test, a reliable and sensitive screening test, would enhance the accuracy of the diagnosis of malaria infection and thus make clinical treatment decision effectively.

NOTE

<i>Pf.: Plasmodium falciparum</i>
<i>Pv.: Plasmodium vivax</i>
<i>Po.: Plasmodium ovale</i>
<i>Pm.: Plasmodium malariae</i>

Intended use

STANDARD Q Malaria Pf Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of *Plasmodium falciparum* (*P. falciparum*) specific Histidine Rich Protein Z (HRP-2) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria Pf Ag Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.

Test principle

STANDARD Q Malaria Pf Ag Test contains two pre-coated lines, "P.f" (*P. falciparum*) as test line and "C" as control line on the surface of the nitrocellulose membrane. Both test line and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-*P. falciparum* HRP-2 is coated on the test line region and monoclonal anti-chicken IgY is coated on the control line region. During the test, the *P. falciparum* specific HRP-2 in the specimen reacts to the gold-conjugated monoclonal anti-Malaria HRP-2 and binds to it. Any *P. falciparum* specific HRP-2 antigen-antibody gold particle complex also migrates with the buffer and is immobilized by monoclonal anti-*P. falciparum* HRP-2 at the test line to formation of a violet test colored band which confirms a positive result. Absence of this violet colored band indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

KIT STORAGE AND STABILITY

1. Store the sealed pouch and the buffer provided in the kit at 2-40°C/36-104°F out of the direct sunlight for the duration of its shelf life.
2. Do not open the foil pouch until you are ready to perform a test.
3. Close the buffer cap tightly after using, and then store it at 2-40°C/36-104°F out of the direct sunlight. It is stable until the expiry date of the kit and the buffer label after opening its cap, if it is tightly closed.

WARNINGS AND PRECAUTIONS

1. Do not freeze.
2. Do not use beyond the expiration date.
3. Do not re-use the test kit.
4. Do not use the test kit if the pouch is damaged or the seal is broken.
5. Do not use buffer of another lot.
6. Do not smoke, drink or eat while handling specimen.
7. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.
8. Clean up spills thoroughly using an appropriate disinfectant.
9. Handle all specimens as if they contain infectious agents.
10. Observe established precautions against microbiological hazards throughout testing procedures.
11. Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
12. In the pouch, there is a desiccant containing a humidity indicator, the humidity indicator should be yellow. If the humidity indicator is green, throw away the test device and take another test device. Throw away the desiccant in the nonsharps (non-infectious) disposal.
13. Before testing, check that the lot number of buffer on the bottle matches the lot number of buffer on the package and the pouch you are using.
14. The buffer contains 0.095% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.
15. For *in vitro* diagnostic use only.

LIMITATION OF TEST

1. The test procedure, precautions, and interpretation of results for this test must be followed strictly when testing.
2. STANDARD Q Malaria Pf Ag Test is designed for use only on human capillary and venous whole blood specimens.
3. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

QUALITY CONTROL

1. A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent is reactive.
2. Control materials are not supplied with this test kit. However, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
3. If there is a problem with the result such as invalid result, retest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

***This test only indicates the presence of Histidine-Rich protein 2 (HRP2) of Plasmodium falciparum in human capillary and venous whole blood specimen and should not be used as the sole criteria for the diagnosis of malaria. As other diagnostic tests, all test results should be considered with other clinical history available.**

Potential Cross-reacting Substances

STANDARD Q Malaria Pf Ag Test had no cross-reaction with the following microorganisms such as other hepatitis (Hepatitis A virus, Hepatitis B virus, Hepatitis C virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, E. coli (K99), *Salmonella* Typhi), *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Legionella pneumophila*, *Treponema pallidum*. Also, this product had no cross-reaction with HIV Seroconversion panel, HCV positive plasma, HCV core Ag positive plasma, Hbs Ag positive plasma, HAV IgM positive plasma, Dengue IgM positive plasma, Tick borne encephalitis IgM positive plasma, Chikungunya positive plasma, Zika virus positive plasma, *Salmonella* typhi IgM positive plasma, *Leishmania* positive plasma, Japanese Encephalitis positive plasma, *Brucella* IgM positive plasma and West Nile Virus positive plasma.

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