WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: STANDARD Q Malaria P.f/P.v Ag Test
WHO reference number: PQDx 0348-117-00

STANDARD Q Malaria P.f/P.v Ag Test with product code **09MAL20D**, manufactured by **SD Biosensor, Inc, CE-Mark regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 6 March 2020.

Summary of WHO prequalification assessment for STANDARD Q Malaria P.f/P.v Ag Test

	Date	Outcome
Prequalification listing	6 March 2020	listed
Dossier assessment	29 October 2018	MR
Site inspection(s) of the quality management system	24-26 May 2023 14-16 June 2023	MR
Product performance evaluation	2018	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendments	Date of report amendment
2.0	1. Change of IFU from IFU in English to multilingual IFU (EN, ES, FR, PT), 2. Change of outer package.	18 May 2021.
3.0	1. Addition of manufacturer site for Cutting, Assembly, Buffer Preparation & dispensing, sealing, printing, packing and shipping process. 2. Addition of suppliers for raw materials. (1) Upper device, (2) Lower device, (3) Aluminium Pouch, (4) Inverted cup, (5) Lancet, (6) Alcohol swab". 3. Update lancet label.	19 February 2024.

4.0	The addition of the manufacturing site is located in the public report of Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India.	23 April 2024
5.0	Updated labels of inverted cup, sterile lancet and alcohol swab.	5 August 2025

Intended use

According to the claim of intended use from SD Biosensor, Inc., "*STANDARD Q Malaria P.f/P.v Ag Test is a rapid and membrane-based immunochromatography for the qualitative detection of Plasmodium falciparum (P. falciparum) specific Histidine Rich Protein 2 (HRP-2) and Plasmodium vivax specific Plasmodium lactate dehydrogenase (pLDH) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f/P.v Ag Test is intended to be used by trained healthcare or laboratory professionals or other healthcare workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.*"

Assay description

According to the claim of assay description from SD Biosensor, Inc., "*STANDARD Q Malaria P.f/P.v Ag Test contains two-pre-coated lines, "P.f" (P. falciparum line), "P.v" (P. vivax) as test lines and "C" as control line on the surface of the nitrocellulose membrane. The test lines and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-P. falciparum HRP-2 is coated on the P.f test line region, monoclonal anti-Malaria P. vivax LDH is coated on the P.v test line region and monoclonal anti-chicken IgY is coated on the control line region. During the test, the P. falciparum specific HRP-2 and/or P. vivax specific pLDH in the specimen react to the gold-conjugated monoclonal anti-Malaria HRP-2 and/or gold-conjugated monoclonal anti-Malaria P. vivax LDH and then bind to them respectively. Any P. falciparum specific HRP-2 antigen-antibody gold particle complex and/or P. vivax specific pLDH antigen-antibody gold particle complex also migrate with the buffer and are immobilized by monoclonal anti-P. falciparum HRP-2 and/or monoclonal anti-Malaria P. vivax LDH at the two individual test lines to formation of violet test colored band(s) which confirms a positive result. Absence of this violet colored bands indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.*"

Test kit contents

Component	25 tests (product code 09MAL20D)
Test device (individually in a foil pouch with desiccant)	25
Buffer Bottle	1 x 4 mL
Inverted cup (5µl)	Pack of 25
Sterile lancet	Pack of 25
Alcohol swab	Pack of 25
Instructions for use	1

Items required but not provided

- Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood,
- Micropipette and tip,
- Timer,
- PPE (Personal Protective Equipment),
- Pen/pencil,
- Extra lancets and alcohol swabs,
- Biosafety sharps container,
- Biohazard container,
- Sterile gauze.

Storage

The test kit must be stored at 2-40°C.

Shelf-life upon manufacture

24 months.

Warnings/limitations

Refer to the current version of the Manufacturer's instructions for use.

Prioritization for Prequalification

Based on the established eligibility criteria, the STANDARD Q Malaria P.f/P.v Ag Test was given priority for the WHO prequalification assessment.

Dossier assessment

SD Biosensor, Inc. submitted a product dossier for STANDARD Q Malaria P.f/P.v Ag Test as per the "*Instructions for compilation of a product dossier*" (PQDx_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the discrepancies found during dossier screening and assessment findings were accepted on 29 October 2018.

Based on the product dossier screening and assessment findings, the product dossier for STANDARD Q Malaria P.f/P.v Ag product name meets WHO prequalification requirements.

Manufacturing site inspection

An inspection of SD Biosensor Inc., located at 74, *Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea* was conducted between 24-26 May 2023 and Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India between 14-16 June 2023. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of the product of consistent quality. Routine inspections of the Manufacturer will be conducted with copies of these WHO Public Inspection Reports (WHOPIRs) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current inspection performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection can be found at:

<https://extranet.who.int/prequal/inspection-services/prequalification-reports/whopirs-vitro-diagnostics>

All published WHOPIRs are with the manufacturer's agreement.

Product performance evaluation

STANDARD Q Malaria P.f/P.v Ag Test¹ was included in the eighth round of WHO product testing of RDTs for malaria antigen detection, which was completed in 2018. The product

¹ The product code of the product included in the WHO product testing of malaria RDTs was 09MAL20B, which corresponded to the same product, with a different kit configuration (i.e. number of tests) <https://www.who.int/malaria/publications/atoz/9789241514965/en/>

was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild-type parasite panel, *P. vivax* wild-type parasite panel and a negative panel. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

Based on the demonstrated *P. falciparum* panel detection score (85.0% at 200 parasites/ μ l), *P. vivax* panel detection score (100% at 200 parasites/ μ l), false-positive rates (0% for clean negatives, 0.5% for *P. falciparum* at 200 parasites/ μ l, 0% for *P. vivax* at 200 parasites/ μ l, 0.5% for *P. falciparum* at 2000 parasites/ μ l, 0% for *P. vivax* at 2000 / μ l) and invalid rate (0%), STANDARD Q Malaria P.f/P.v Ag Test meet the current laboratory evaluation requirements for prequalification.

Summary performance characteristics	Panel detection score		False positive rate (%)			Invalid rate (%)
	200 parasites/ μ l		200 parasites/ μ l		Clean negatives	
	Pf	Pv	Pf	Pv		
STANDARD Q Malaria P.f/P.v Ag Test	85.0	100	0.5	0	0	0

Based on these results, the STANDARD Q Malaria P.f/P.v Ag Test meets the current performance evaluation requirements for prequalification.

Labelling

- 1. Labels**
- 2. Instructions for use**

1. Labels

1.1 Device Package

Note: Three versions (Figures 1.1a-1.1c) are valid and may be found in distribution.

STANDARD Q Malaria P.f/P.v Ag 25T


Figure 1.1a


Figure 1.1b

STANDARD Q Malaria P.f/P.v Ag 25T


Figure 1.1c

1.2 Foil pouch

Note: Both versions(Figure 1.2a and Figure 1.2b) are valid and may be found in distribution.


Figure 1.2a


Figure 1.2b

1.3 Buffer labels

Note: Both versions(Figure 1.3a and Figure 1.3b) are valid and may be found in distribution.


Figure 1.3a


Figure 1.3b

1.4 Inverted cup(5 μ l) label

Note: Previous version (Figure 1.4a) and updated version (Figure 1.4b). Both label versions are valid and may be found in distribution.

Figure 1.4a

Figure 1.4b

1.5 Sterile Lancet label

Note: Previous versions (Figure 1.5a and 1.5b) and updated versions (Figure 1.5c and 1.5d). Both label versions are valid and may be found in distribution.

Figure 1.5a

Disposable Sterile Lancets

 LOT No. : XXXXXX Product code : 01GL27
 MFG Date : YYYY.MM.DD. **25PCS**
 EXP Date : YYYY.MM.DD. **28G**

INTENDED USE
To obtain a capillary blood specimen from the fingertip.

INSTRUCTIONS FOR USE
To use, twist-off the protective cap.

CAUTION
The lancet is guaranteed sterile while protective cap is sealed to the base. Do not use if the seal has been damaged or broken.

 **Manufactured by**
 Tianjin Huahong Technology Co., Ltd.
 A01, Plant B No. 278, Hangkong Road,
 Tianjin Pilot Free Trade Zone
 (Air Port Industrial Park),
 300308 Tianjin, China

 **Authorized Representative**
 Shanghai International Holding Corp.
 GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel:+49-40-2513175

 0197
 
 
 
 

Figure 1.5b

Blood Lancet

 LOT No. : 240311 Product code : BL-28G
 MFG Date : 2024-03-11 **25 PCS**
 EXP Date : 2029-02-11 **28G**

INTENDED USE
To obtain the capillary blood specimen for test.

INSTRUCTIONS FOR USE
To use, twist-off the protective cap.

CAUTION
The lancet is guaranteed sterile while protective cap is sealed to the base. Do not use if the seal has been damaged or broken.

 **Manufactured by**
 Guangzhou iCare Medical Technology Co., Ltd.
 First floor A No.8, Lianhua Port Industrial Zone,
 Lotus Mountain Bonded Area, Shilou Town,
 Panyu District, Guangzhou,
 People's Republic of China.
 Tel: +86-020-31169913 Fax: +86-020-84881433

 **Authorized Representative**
 Shanghai International Holding Corp. GmbH (Europe)
 Eiffestrasse 80, 20537 Hamburg, Germany
 Tel: +49-40-2513175 Fax: +49-40-2513175

 0123
 
 
 
 

Figure 1.5c


Figure 1.5d

1.6 Alcohol swab label

Note: Note: Updated version (Figure 1.6a) and previous version (Figure 1.6b). Both label versions are valid and may be found in distribution.

Figure 1.6a


Figure 1.6b

1.7 Cassette image

2. Instructions for use²

² English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.


STANDARD Q
Malaria P.f/P.v Ag

STANDARD™ Q Malaria P.f/P.v Ag Test

DE, IT

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

**KIT CONTENTS****MATERIALS REQUIRED BUT NOT PROVIDED**

- Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood
- Micropipette and tip
- PPE (Personal Protective Equipment)
- Pen/pencil
- Extra lancets and alcohol swabs
- Biosafety sharps container
- Biohazard container
- Sterile gauze

SPECIMEN COLLECTION AND PREPARATION**Capillary whole blood**

- Capillary whole blood should be collected aseptically by fingerstick.
- Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
- Squeeze the end of the fingertip and pierce with a sterile lancet.
- Wipe the first drop with a sterile gauze.
- Collect capillary whole blood and it must be tested immediately after collection.

Venous whole blood

- Collect venous whole blood into a commercially available anti-coagulant tube such as heparin, EDTA, or sodium citrate by venipuncture.
- If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 72 hours after collection.
- Do not use hemolyzed blood specimen.

CAUTION


Anticoagulants such as heparin, EDTA, or sodium citrate do not affect the test result.

As known relevant interference, hemolytic specimen, rheumatoid factors-containing specimen and lipaemic, icteric specimen may lead to inaccurate results.

Use separate inverted cup for each specimen in order to avoid cross-contamination of either specimens which can cause erroneous results.

PREPARATION AND TEST PROCEDURE**Preparation**

- Carefully read the instructions for using the STANDARD Q Malaria P.f/P.v Ag Test.
- Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
- Open the foil pouch, and check the test device and the desiccant inside the foil pouch.

**Test Procedure**


Capillary and venous whole blood can be used for the STANDARD Q Malaria P.f/P.v Ag Test.

1. Test with a sterile lancet


- Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.


- Dry and pierce the wiped fingertip with a sterile lancet to bleed.


- Wipe away the first drop of blood with sterile gauze.


- Take the inverted cup (5µl) provided and collect 5µl of the capillary whole blood specimen by dipping the circular end of the inverted cup into the specimen.


- To apply the collected capillary whole blood specimen to the specimen well(circle) of the test device, place the circular end of the inverted cup on the specimen well pad and then press down lightly.

**1-1. Test by venipuncture**


- Gently invert 4-5 times for proper mixing.


- Using a micropipette, draw 5µl of venous whole blood specimen.

- Apply the collected venous whole blood specimen to the specimen well(circle) of the test device.

**2. Dropping of Buffer**

- Hold the buffer bottle at 90° (perpendicular) to the test device and not touching the specimen well to avoid contamination. Add 3 drops (90µl) of the buffer into the buffer well(square) of the test device.

**3. Reading Time**


- Read the test result in 15 minutes. The test can be read up to 30 minutes.


- Do not read test results after 30 minutes. It may give false results.
- If the T line does not appear until 15 minutes, make sure to observe the appearance of the T line until 30 minutes.

INTERPRETATION OF TEST RESULT

* "C" Control Line | "P.v" P.v Line | "P.f" P.f Line

The presence of only one colored band ("C" Control line) within the result window indicates a negative result.

P.f. Positive	P.v. Positive	Mixed P.f. and P.v. Positive
The presence of two colored bands ("C" Control line and "P.f" <i>Plasmodium falciparum</i> line) within the result window, no matter which band appears first, indicates a P.f. positive result.	The presence of two colored bands ("C" Control line and "P.v" <i>Plasmodium vivax</i> line) within the result window, no matter which band appears first, indicates a P.v. positive result.	The presence of three colored bands ("C" Control line, "P.f" <i>Plasmodium falciparum</i> line and "P.v" <i>Plasmodium vivax</i> line) within the result window, no matter which band appears first, indicates a mixed P.f. and P.v. positive result.

CAUTION

The presence of any line, no matter how faint the result is considered positive.

Invalid

If the control band ("C" Control line) is not visible within the result window, the result is considered invalid. The directions may not have been followed correctly or the test device may have deteriorated. Re-test with a new specimen and a new test device.


***NOTE • P.f.: Plasmodium falciparum, • P.v.: Plasmodium vivax**

- A colored band will appear in the top section of the result window to show that the test is working properly. This band is the control line (C).
- Colored bands will appear in the middle and lower section of the result window. These bands are the *Plasmodium vivax* line (P.v) and the *Plasmodium falciparum* line (P.f) respectively.

***This test only indicates the presence of Histidine-Rich Protein-2 (HRP2) of Plasmodium falciparum and Plasmodium lactate dehydrogenase (pLDH) of Plasmodium vivax in human capillary and venous whole blood specimen and should not be used as the sole criteria for the diagnosis of malaria. As other diagnostic tests, all test results should be considered with other clinical history available.**

EXPLANATION AND SUMMARY**Introduction**

Malaria remains an important cause of illness and death in children and adults in countries in which it is endemic. Globally, an estimated 3.2 billion people in 97 countries and territories are at risk of being infected with malaria and dengue disease, and 1.2 billion are at high risk (>1 in 1,000 chance of getting malaria in a year). According to the World Malaria Report 2015, there were 210 million cases of malaria globally in 2015 (uncertainty range 149-303 million) and 438,000 malaria deaths (range 236,000-635,000), representing a decrease in malaria cases and deaths of 37% and 60% since 2,000, respectively. The prototypical parasite *Plasmodium falciparum* is from *Plasmodium falciparum*, *vivax*, *ovale* and *malariae* with the first two species causing the most infections worldwide. Classic symptoms of malaria include fever, headache, chills, vomiting, shivering and convulsions. In some rare forms of *P. falciparum*, the patient may present with delirium or coma. Severe anemia is often attributed to the cause of death from malaria. Accurate and prompt diagnosis of malaria is of utmost importance due to the morbidity associated with the other malarial forms. Rapid diagnostic test is an ideal diagnostic tool for malaria diagnosis in that it can provide a rapid determination if the patient is infected with malaria allowing for accurate treatment and improved outcomes. STANDARD Q Malaria P.f/P.v Ag Test, a reliable and sensitive screening test, would enhance the accuracy of the diagnosis of malaria infection and thus make clinical treatment decision effectively.

**Intended use**

STANDARD Q Malaria P.f/P.v Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of *Plasmodium falciparum* (*P. falciparum*) specific Histidine Rich Protein 2 (HRP-2) and *Plasmodium vivax* specific Plasmodium lactate dehydrogenase (pLDH) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f/P.v Ag Test is intended to be used by trained healthcare or laboratory professionals or other healthcare workers who have received appropriate training. This product can be used by trained by providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.

Test principle

STANDARD Q Malaria P.f/P.v Ag Test contains two pre-coated lines, "P.f" (*P. falciparum*), "P.v" (*P. vivax*) as test lines and "C" as control line on the surface of the nitrocellulose membrane. The test lines and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-*P. falciparum* HRP-2 is coated on the P.f test line region, monoclonal anti-Malaria *P. vivax* LDH is coated on the P.v test line region and monoclonal anti-*IgG* is coated on the control line region. During the test, the *P. falciparum* specific HRP-2 and/or *P. vivax* specific pLDH in the specimen react to the gold-conjugated monoclonal anti-Malaria HRP-2 and/or gold-conjugated monoclonal anti-Malaria *P. vivax* LDH and then bind to them respectively. Any *P. falciparum* specific HRP-2 antigen-antibody gold particle complex and/or *P. vivax* specific pLDH antigen-antibody gold particle complex also migrate with the buffer and are immobilized by monoclonal anti-*P. falciparum* HRP-2 and/or monoclonal anti-Malaria *P. vivax* LDH at the two individual test lines to form of violet test colored band(s) which confirms a positive result. Absence of this violet colored bands indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

KIT STORAGE AND STABILITY

- Store the sealed pouch and the buffer provided in the kit at 2-40°C/36-104°F out of the direct sunlight for the duration of its shelf life.
- Do not open the foil pouch until you are ready to perform a test.
- Close the buffer cap tightly after using, and then store it at 2-40°C/36-104°F out of the direct sunlight. It is stable until the expiry date of the kit and the buffer label after opening its cap, if it is tightly closed.

WARNINGS AND PRECAUTIONS

- Do not freeze.
- Do not use beyond the expiration date.
- Do not reuse the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use buffer of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- In the pouch, there is a desiccant containing a humidity indicator; the humidity indicator should be yellow. If the humidity indicator is green, throw away the test device and take another test device. Throw away the desiccant in the nonsharps (non-infectious) disposal.
- Before testing, check that the lot number of buffer on the bottle matches the lot number of buffer on the package and the pouch you are using.
- The buffer contains 0.095% sodium azide as a preservative which may

Potential Cross-reacting Substances

STANDARD Q Malaria P.f/P.v Ag Test had no cross-reaction with the following microorganisms such as other hepatitis (Hepatitis A virus, Hepatitis B virus, Hepatitis C virus), Respiratory infectious disease (Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, *E. coli* (k99), *Salmonella* Typhi), *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Legionella pneumophila*, *Treponema pallidum*. Also, this product had no cross-reaction with HIV Seroconversion panel, HCV positive plasma, HCV core Ag positive plasma, Hbs Ag positive plasma, HAV IgM positive plasma, Dengue IgM positive plasma, Tick borne encephalitis IgM positive plasma, Chikungunya positive plasma, Zika virus positive plasma, *Salmonella* typhi IgM positive plasma, *Leishmania* positive plasma, Japanese Encephalitis positive plasma, *Brucella* IgM positive plasma and West Nile Virus positive plasma.

Diagnostic Sensitivity & Specificity**1. Venous whole blood**

1618 venous whole blood specimen were tested using the STANDARD Q Malaria P.f/P.v Ag Test to determine its clinical sensitivity and specificity.

be toxic if ingested. When disposed of through a sink, flush with large quantities of water.

15. In vitro diagnostic use only.**LIMITATION OF TEST**

- The test procedure, precautions, and interpretation of results for this test must be followed strictly when testing.
- STANDARD Q Malaria P.f/P.v Ag Test is designed for use only on human capillary and venous whole blood specimens.
- Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

QUALITY CONTROL

- A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent is reactive.
- Control materials are not supplied with this test kit. However, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
- If there is a problem with the result such as invalid result, retest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

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BIBLIOGRAPHY

- World Malaria Report 2015, World Health Organization, p.8-12, 2015.
- Confronting *Plasmodium vivax* malaria, World Health Organization, p.1-3, 2015.
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition, U.S. Department of Health and human Services, CDC, NIH, Washington DC, P.182-184, 2007
- Guidelines for the treatment of malaria, World Health Organization, p.5-6, 2006.
- Recommended selection criteria for procurement of malaria rapid diagnostic tests, World Health Organization, p.1-4, 2015.
- C. Beaudic, G.W. Long, M. C. Eby, C. Hoffman, G.W. Long, W.R. Weiss, S.M. Maret, A.J. Oloo, Diagnosis of malaria by detection of *Plasmodium falciparum* HRP-2 antigen with a rapid diagnostic antigen-capture assay, The Lancet. Vol. 343, Issue: 8897, 1994, p.564-568.
- A. Moody, Rapid diagnostic tests for malaria parasites, Clin Microbiol Rev. Vol. 15, Issue 1, 2002, p.66-78.
- LB Ochola, P. Vounatsou, T. Smith, MLH Mabaso, CJ Newton, and The reliability of diagnostic techniques in the diagnosis and

management of malaria in the absence of a gold standard, Lancet Infect Dis. Vol. 6, Issue 9, 2006, p.582-588.

- K. Abba, M. A. M. S. Olliaro, C.M. Neang, S.M. Jackson et al. Rapid diagnostic tests for diagnosing uncomplicated *P. falciparum* malaria in endemic countries. Cochrane Database Syst Rev. Vol. 18, Issue 12, 2011.
- MU Ahmed, M. A. M. S. Olliaro, AKM Shah Suzzaman, MU Alam, AH Khan et al. Rapid Diagnosis of malaria by Antigen Detection. Bangladesh Journal of Medical Microbiology, Bangladesh J Med Microbiol. Vol. 3, Issue 1, 2009, p. 14-16.

SYMBOL

REF	Reference number	Caution	Use by	LOT	Batch code	Consult Instructions for Use	Do not re-use
IVD	In vitro Diagnostics	Note	Manufacturer		Date of manufacture	Contains sufficient for <n> Tests	Keep away from sunlight
	Indicate that you should keep the product dry	To indicate the temperature limitations in which the transport package has to be kept and handled.			Full fill the requirements of Directive 98/79/EC on in vitro diagnostic medical devices	Do not use if packaging is damaged	

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Manufacturing site

Plot no. 38, Sector 4, IIMT Manesar, Gurgaon, Haryana 122052, INDIA
Please contact us for any complaints/inquiries/suggestions via email (ts@sdbiosensor.com), phone (+82-80-970-9700) or website (www.sdbiosensor.com).

L23MMALZINRO-WHO

Issue date: 2023.11


STANDARD Q Malaria P.f/P.v Ag
STANDARD™ Q Malaria P.f/P.v Ag Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST **SD BIOSENSOR**

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
Anticoagulants such as heparin, EDTA, or sodium citrate do not affect the test result. As known relevant interference, hemolytic specimen, rheumatoid factors-contained specimen and lipemic, icteric specimen may lead to inaccurate results. Use separate inverted cup for each specimen in order to avoid cross-contamination of either specimens which can cause erroneous results.

PREPARATION AND TEST PROCEDURE

- Preparation**
- Carefully read the instructions for using the STANDARD Q Malaria P.f/P.v Ag Test.
 - Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
 - Open the foil pouch, and check the test device and the desiccant inside the foil pouch.

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



3. Wipe away the first drop of blood with sterile gauze.

4. Take the inverted cup (5µl) provided and collect 5µl of the capillary whole blood specimen by dipping the circular end of the inverted cup into the specimen.


5. To apply the collected capillary whole blood specimen to the specimen well(circle) of the test device, place the circular end of the inverted cup on the specimen well pad and then press down lightly.


1-1. Test by venipuncture

- Gently invert 4-5 times for proper mixing.
- Using a micropipette, draw 5µl of venous whole blood specimen.

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

3. Apply the collected venous whole blood specimen to the specimen well(circle) of the test device.

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INTERPRETATION OF TEST RESULT

*"C" Control Line "Pv" P.v Line "Pf" P.f Line
Negative
Positive
Pf. Positive
Pv. Positive
Mixed Pf. and Pv. Positive
CAUTION

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- *NOTE • *Pf.*: *Plasmodium falciparum*. • *Pv.*: *Plasmodium vivax*
- A colored band will appear in the top section of the result window to show that the test is working properly. This band is the control line (C).
 - Colored bands will appear in the middle and lower section of the result window. These bands are the *Plasmodium vivax* line (Pv) and the *Plasmodium falciparum* line (Pf) respectively.

*This test only indicates the presence of Histidine-Rich Portein-2 (HRP2) of *Plasmodium falciparum* and *Plasmodium lactate dehydrogenase* (pLDH) of *Plasmodium vivax* in human capillary and venous whole blood specimen and should not be used as the sole criteria for the diagnosis of malaria. As other diagnostic tests, all test results should be considered with other clinical history available.

EXPLANATION AND SUMMARY

Introduction

Malaria remains an important cause of illness and death in children and adults in countries in which it is endemic. Globally, an estimated 3.2 billion people in 97 countries and territories are at risk of being infected with malaria and developing disease, and 1.2 billion are at high risk (v) in 1,000 chance of getting malaria in a year). According to the World Malaria Report 2015, there were 210 million cases of malaria globally in 2015 (uncertainty range 149-303 million) and 438,000 malaria deaths (range 236,000-635,000), representing a decrease in malaria cases and deaths of 27% and 60% since 2,000, respectively. The protozoal parasites that cause malaria are from *Plasmodium falciparum*, *vivax*, *ovale* and *malariae* with the first two species causing the most infections worldwide. Classic symptoms of malaria include fever, headache, chills, vomiting, shivering and convulsions. In some rare forms of *P. falciparum*, the patient may present with delirium or coma. Severe anemia is often attributed to the cause of death from malaria. Accurate and prompt diagnosis of malaria is of utmost importance due to the morbidity associated with the other malarial forms. Rapid diagnostic test is an ideal diagnostic tool for malaria diagnosis in that it can provide a rapid determination if the patient is infected with malaria although not accurate outcomes and improved outcomes. STANDARD Q Malaria P.f/P.v Ag Test, a reliable and sensitive screening test, would enhance the accuracy of the diagnosis of malaria infection and thus make clinical treatment decision effectively.

Pf.: *Plasmodium falciparum*
Pv.: *Plasmodium vivax*
Po.: *Plasmodium ovale*
Pa.: *Plasmodium malariae*

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WARNINGS AND PRECAUTIONS

- Do not freeze.
- Do not use beyond the expiration date.
- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use buffer of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- In the pouch, there is a desiccant containing a humidity indicator; the humidity indicator should be yellow. If the humidity indicator is green, throw away the test device and take another test device. Throw away the desiccant in the nonsharps (non-infectious) disposal.
- Before testing, check that the lot number of buffer on the bottle matches the lot number of buffer on the package and the pouch you are using.
- The buffer contains 0.095% sodium azide as a preservative which may

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STANDARD Q Malaria P.f/P.v Ag Test had no cross-reaction with the following microorganisms such as other hepatitis(Hepatitis A virus, Hepatitis B virus, Hepatitis C virus), Respiratory infectious disease(Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, *E. coli* (k99), *Salmonella* Typhi), *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Legionella pneumophila*, *Treponema pallidum*. Also, this product had no cross-reaction with HIV Seroconversion panel, HCV positive plasma, HCV core Ag positive plasma, Hbs Ag positive plasma, HAV IgM positive plasma, Dengue IgM positive plasma, Tick borne encephalitis IgM positive plasma, Chikungunya positive plasma, Zika virus positive plasma, *Salmonella* typhi IgM positive plasma, *Leishmania* positive plasma, Japanese Encephalitis positive plasma, *Bruceella* IgM positive plasma and West Nile Virus positive plasma.

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- For *in vitro* diagnostic use only.

LIMITATION OF TEST

- The test procedure, precautions, and interpretation of results for this test must be followed strictly when testing.
- STANDARD Q Malaria P.f/P.v Ag Test is designed for use only on human capillary and venous whole blood specimens.
- Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of as a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

QUALITY CONTROL

- A colored line appearing in the control line is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagent is reactive.
- Control materials are not supplied with this test kit. However, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
- If there is a problem with the result such as invalid result, rerest with a new kit and/or specimen. If the problem is repeated, contact SD Biosensor through your local distributor.

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Potential Interfering Substances

STANDARD Q Malaria P.f/P.v Ag Test was evaluated with the following interfering substances present in specimen in order to assess their potential effect on the assay performance as per CLSI guideline EP7-A2. It was not affected by interfering materials such as artemether-lumefantrine (malaria drug), doxycycline hyclate (malaria drug), lamivudine (retroviral medication), acetaminophen, erythromycin, Ibuprofen, acetylsalicylic acid, bilirubin (unconjugate), caffeine, ethanol. There was no interference with Human anti-mouse antibody, whole blood of pregnant women, infant and neonate specimen and whole blood having elevated levels of hemoglobin, elevated levels of C-reactive protein and lipidemic specimen.

[Clinical sensitivity]

-Clinical sensitivity for Malaria *Pf.*: 99.59% (487/489, 95% CI: 98.53-99.95%)
-Clinical sensitivity for Malaria *P.v.*: 100% (123/123, 95% CI: 97.05-100.00%)
[Clinical specificity] 100% (1006/1006, 95% CI: 99.63-100.00%)

STANDARD Q Malaria P.f/P.v Ag Test	Reference Test (Microscopy)			Sum
	P.f.	P.v.	Negative	
P.f line Positive	487	0	0	487
P.v line Positive	0	123	0	123
Negative	2	0	1006	1008
Total	489	123	1006	1618

2. Capillary whole blood

605 capillary whole blood specimen were tested using the STANDARD Q Malaria P.f/P.v Ag Test to determine its clinical sensitivity and specificity.

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
KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

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SPECIMEN COLLECTION AND PREPARATION

- Capillary whole blood**
- Capillary whole blood should be collected aseptically by fingerstick.
 - Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
 - Squeeze the end of the fingertip and pierce with a sterile lancet.
 - Wipe the first drop with a sterile gauze.
 - Collect capillary whole blood and it must be tested immediately after collection.

Venous whole blood

- Collect venous whole blood into a commercially available anti-coagulant tube such as heparin, EDTA, or sodium citrate by venipuncture.
- If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/ 36-46°F, the specimen can be used for testing within 72 hours after collection.
- Do not use hemolyzed blood specimen.

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
As known relevant interference, hemolytic specimen, rheumatoid factors-contained specimen and lipemic, icteric specimen may lead to inaccurate results. Use separate inverted cup for each specimen in order to avoid cross-contamination of either specimens which can cause erroneous results.

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CAUTION

Bring the kit contents and the specimens to room temperature before testing.

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1. Test with a sterile lancet

- Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
- Dry and pierce the wiped fingertip with a sterile lancet to bleed.

- Wipe away the first drop of blood with sterile gauze.
- Take the inverted cup (5µl) provided and collect 5µl of the capillary whole blood specimen by dipping the circular end of the inverted cup into the specimen.
- To apply the collected capillary whole blood specimen to the specimen well(circle) of the test device, place the circular end of the inverted cup on the specimen well pad and then press down lightly.

1-1. Test by venipuncture

- Gently invert 4-5 times for proper mixing.
- Using a micropipette, draw 5µl of venous whole blood specimen.
- Apply the collected venous whole blood specimen to the specimen well(circle) of the test device.

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3. Reading Time


- Read the test result in 15 minutes. The test can be read up to 30 minutes.
- Do not read test results after 30 minutes. It may give false results.
 - If the T line does not appear until 15 minutes, make sure to observe the appearance of the T line until 30 minutes.

INTERPRETATION OF TEST RESULT

*"C" Control Line	"Pv" Pv Line	"Pf" Pf Line	Result
+	-	-	Negative
+	+	-	Positive
+	-	+	Mixed Pf, and Pv. Positive
+	+	+	Mixed Pf, and Pv. Positive

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- *NOTE • Pf: *Plasmodium falciparum*, • Pv: *Plasmodium vivax*
- A colored band will appear in the top section of the result window to show that the test is working properly. This band is the control line (C).
 - These bands are the *Plasmodium vivax* line (P.v) and the *Plasmodium falciparum* line (P.f) respectively.

*This test only indicates the presence of Histidine-Rich Protein-2 (HRP2) of *Plasmodium falciparum* and *Plasmodium lactate dehydrogenase (pLDH) of Plasmodium vivax* in human capillary and venous whole blood specimens and should not be used as the sole criteria for the diagnosis of malaria. As other diagnostic tests, all test results should be considered with other clinical history available.

EXPLANATION AND SUMMARY

Introduction

Malaria remains an important cause of illness and death in children and adults in countries in which it is endemic. Globally, an estimated 3.2 billion people in 97 countries and territories are at risk of being infected with malaria and developing disease, and 1.2 billion are at high risk (>1 in 1,000 chance of getting malaria in a year). According to the World Malaria Report 2015, there were 210 million cases of malaria globally in 2015 (uncertainty range 149-303 million) and 438,000 malaria deaths (range 236,000-635,000), representing a decrease in malaria cases and deaths of 37% and 60% since 2,000, respectively. The protocol parasites that cause malaria are from *Plasmodium falciparum*, *vivax*, *ovale* and *malariae* with the first two species causing the most infectious worldwide. Classic symptoms of malaria include fever, headache, chills, vomiting, shivering and convulsions. In some rare forms of *P. falciparum*, the patient may present with delirium or coma. Severe anemia is often attributed to the cause of death from malaria. Accurate and prompt diagnosis of malaria is of utmost importance due to the morbidity associated with the other malarial forms. Rapid diagnostic test is an ideal diagnostic tool for malaria diagnosis in that it can provide a rapid determination if the patient is infected with malaria allowing for accurate treatment and improved outcomes. STANDARD Q Malaria P.f/P.v Ag Test, a reliable and sensitive screening test, would enhance the accuracy of the diagnosis of malaria infection and thus make clinical treatment decision effectively.

NOTE

Pf: *Plasmodium falciparum*
Pv: *Plasmodium vivax*
Po: *Plasmodium ovale*
Pm: *Plasmodium malariae*

Intended use

STANDARD Q Malaria P.f/P.v Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of *Plasmodium falciparum* (*P. falciparum*) specific Histidine Rich Protein 2 (HRP-2) and *Plasmodium vivax* specific Plasmodium lactate dehydrogenase (pLDH) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f/P.v Ag Test is intended to be used by trained healthcare or laboratory professionals or other healthcare workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.

Test principle

STANDARD Q Malaria P.f/P.v Ag Test contains two pre-coated lines, "Pf" (*P. falciparum*), "Pv" (*P. vivax*) as test lines and "C" as control line on the surface of the nitrocellulose membrane. The test lines and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-*P. falciparum* HRP-2 is coated on the Pf test line region, monoclonal anti-Malaria *P. vivax* LDH is coated on the Pv test line region and monoclonal anti-chicken IgY is coated on the control line region. During the test, the *P. falciparum* specific HRP-2 and/or *P. vivax* specific pLDH in the specimen react to the gold-conjugated monoclonal anti-Malaria HRP-2 and/or gold-conjugated monoclonal anti-Malaria *P. vivax* LDH and then bind to them respectively. Any *P. falciparum* specific HRP-2 antigen-antibody gold particle complex and/or *P. vivax* specific pLDH antigen-antibody gold particle complex also migrate with the buffer and are immobilized by monoclonal anti-*P. falciparum* HRP-2 and/or monoclonal anti-Malaria *P. vivax* LDH at the two individual test lines to formation of violet test colored bands) which confirms a positive result. Absence of this violet colored bands indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

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WARNINGS AND PRECAUTIONS

- Do not freeze.
- Do not use beyond the expiration date.
- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use buffer of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- In the pouch, there is a desiccant containing a humidity indicator. The humidity indicator should be yellow. If the humidity indicator is green, throw away the test device and take another test device. Throw away the desiccant in the nonsharps (non-infectious) disposal.
- Before testing, check that the lot number of buffer on the bottle matches the lot number of buffer on the package and the pouch you are using.
- The buffer contains 0.095% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.
- For *in vitro* diagnostic use only.

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STANDARD Q Malaria P.f/P.v Ag Test had no cross-reaction with the following microorganisms such as other hepatitis(Hepatitis A virus, Hepatitis B virus, Hepatitis C virus), Respiratory infectious disease(Influenza A, B), Vector-borne infectious disease (Chikungunya virus, Zika virus, Dengue virus, Yellow fever virus), Intestinal infectious disease (Rotavirus, *E. coli* (k99), *Salmonella* Typhi), *Streptococcus pyogenes*, *Legionella pneumophila*, *Treponema pallidum*. Also, this product had no cross-reaction with HIV Seroconversion panel, HCV positive plasma, HCV core Ag positive plasma, Hbs Ag positive plasma, HAV IgM positive plasma, Dengue IgM positive plasma, Tick borne encephalitis IgM positive plasma, Chikungunya positive plasma, Zika virus positive plasma, *Salmonella* typhi IgM positive plasma, *Leishmania* positive plasma, Japanese Encephalitis positive plasma, *Bruceella* IgM positive plasma and West Nile Virus positive plasma.

Diagnostic Sensitivity & Specificity

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STANDARD Q Malaria P.f/P.v Ag Test was evaluated with the following interfering substances present in specimen in order to assess their potential effect on the assay performance as per CLSI guideline EP7-A2. It was not affected by interfering materials such as artemether-lumefantrine (malaria drug), doxycycline hyclate (malaria drug), lamivudine (retroviral medication), acetaminophen, erythromycin, Ibuprofen, acetylsalicylic acid, bilirubin (unconjugate), caffeine, ethanol. There was no interference with Human anti-mouse antibody, whole blood of pregnant women, infant and neonate specimen and whole blood having elevated

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- Recommended selection criteria for procurement of malaria rapid diagnostic tests, World Health Organization, p.1-4, 2015.
- Beattie, G.W Long, Mc Elroy, Hoffman, G.W Long, W.R Weiss, S.M Maret, A.J Oloo, Diagnosis of malaria by detection of *Plasmodium falciparum* HRP-2 antigen with a rapid dipstick antigen-capture assay. The Lancet, Vol. 343, Issue. 8897, 1994, p.564-568.
- A. Moody, Rapid diagnostic tests for malaria parasites. Clin Microbiol Rev, Vol. 15, Issue 1, 2002, p.66-78.
- LB Ochoha, P Youmansou, T Smith, MLM Mabaso, CRJC Newton. The reliability of diagnostic techniques in the diagnosis and

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SYMBOL

REF	Reference number		Caution		Use by	LOT	Batch code		Consult Instructions for Use		Do not re-use
IVD	<i>In vitro</i> Diagnostics		Note		Manufacturer		Date of manufacture		Contains Sufficient for n-Tests		Keep away from sunlight
	Indicate that you should keep the product dry		To indicate the temperature limitations in which the transport package has to be kept and handled.		Fullfill the requirements of Directive 98/79/EC on <i>in vitro</i> diagnostic medical devices		Do not use if packaging is damaged				

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