WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

Product: Wondfo HIV Self-Test
WHO reference number: PQDx 0357-004-01

Wondfo HIV Self-Test with product codes W006P0058, W006P0059, W006P0060, and W006P0079, manufactured by Guangzhou Wondfo Biotech Co., Ltd, Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 13 July 2022.

Summary of WHO's Prequalification Assessment for the Wondfo HIV Self-Test¹

	Date	Outcome
Prequalification listing for the Wondfo HIV Self-Test	13 July 2022	Listed
Prequalification listing for the One Step HIV1/2 Whole Blood/Serum/Plasma Test	29 November 2018	Listed
Dossier assessment	26 April 2018	MR
Product performance evaluation	First-quarter of 2018	MR

MR: Meets Requirements

*Change notification

In 2021, Guangzhou Wondfo Biotech Co., Ltd submitted a change notification for their prequalified product, One Step HIV1/2 Whole Blood/Serum/Plasma Test, to introduce a new configuration with an intended use specific for HIV self-testing (Wondfo HIV Self-Test). The new format was adapted from the corresponding professional-use product (One Step HIV1/2 Whole Blood/Serum/Plasma Test), for which a WHO prequalification assessment has already been completed. Guangzhou Wondfo Biotech Co., Ltd generated additional data to meet requirements for self-testing as set out in the WHO Technical Specifications Series document TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and self-testing².

¹ Dossier assessment and product performance evaluation for the Wondfo HIV Self-Test were considered from the previous assessment of professional use product, One Step HIV1/2 Whole Blood/Serum/Plasma Test, which was prequalified in 2018. Based on the product dossier assessment and product performance evaluation, Wondfo HIV Self-Test meets WHO prequalification requirements. Please refer to <https://extranet.who.int/prequal/WHOPR/public-report-one-step-hiv12-whole-bloodserumplasma-test-pqdx-0357-004-00>

² <https://apps.who.int/iris/bitstream/handle/10665/251857/9789241511742-eng.pdf;jsessionid=153ABC9D88E7623A1AD1DF946A22B4C8?sequence=1>

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendment and change request reference, where applicable.	Date of report amendment
2.0	Added local languages to IFU. The changes included: 1. The length of the box for W006P0059(20T/kit) increased by 2 cm, from 230X185X105mm to 250X185X105mm; another artwork of the box has not changed. 2. Two optional sizes of IFU, the prequalified 420*285 mm and the additional 570*420 mm.	19 July 2024
3.0	Change the section 'PERFORMANCE CHARACTERISTICS' of the current Instructions for Use by supplementing clinical evaluation data obtained from healthcare professionals. Addition of Portuguese to the multilingual IFU (PQC-IVD-2025-0002). New configuration containing a disposable bag and plaster. Addition of French and Portuguese to the labels (PQC-IVD-2025-0063).	15 January 2026

Intended use:

According to the intended use claim from Guangzhou Wondfo Biotech Co., Ltd, *“the Wondfo HIV Self-Test is a single-use in vitro diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as a self-test and/or by medical professionals. For in vitro diagnostic use only.”*

Test kit contents

Catalog No.		W006P0058	W006P0059	W006P0060	W006P0079
Components					
Test cassette pouch	Test cassette (pcs)	1	1x20	1x100	1x20
	Desiccant (pcs)	1	1x20	1x100	1x20
Accessories	Dropper (pcs)	1	1x20	1x100	1x20
	Buffer (vial)	1	1x20	1x100	1x20
	IFU (pcs)	1	1x20	1x100	1x20
	Blood Lancet for Single Use (pcs)	1	1x20	1x100	1x20
	Alcohol prep pad (pcs)	1	1x20	1x100	1x20
	Cotton swab (pcs)	1	1x20	1x100	1x20
	Disposal bag (pcs)	1	/	/	1x20
	Plaster (pcs)	/	/	/	1x20

Materials required but not provided:

- Timer
- Tissue

Storage

The test kit must be stored between 2 and 30 °C.

Shelf-life upon manufacture³

24 months (test kit and buffer).

Labelling review

The labelling submitted for the Wondfo HIV Self-Test was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

³ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

Controlled Labelling References

Catalog No.	Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.
W006P0058	Outer box artwork	Box artwork	A5	2025/03/25	DLA-W006(3)-01-01
	Pouch / Device label	Foil pouch	A4	2025/03/25	DAT-W006(3)-01-01
	Test cassette foil pouch	Test cassette foil pouch	A3	2023/03/31	DLC-W006(3)-01-01
	Accessory labeling (e.g., pipettes, buffer caps)	Accessories pouch	A2	2025/03/25	DAV-W006(3)-01-01
W006P0059	Outer box artwork	Box artwork	230*185*105mm: A2; 250*185*105mm: A1	2023/03/31 ; 2025/03/25	DLA-W006(3)-01-02
	Outer box label	Box label	A5	2025/03/25	DLB-W006(3)-01-02
	Pouch / Device label	Foil pouch	A4	2025/03/25	DAT-W006(3)-01-02
	Test cassette foil pouch	Test cassette foil pouch	A3	2023/03/31	DLC-W006(3)-01-02
	Accessory labeling (e.g., pipettes, buffer caps)	Same as W006P0058			
W006P0060	Outer box artwork	Box artwork	/	/	DLA-W006(3)-01-03
	Outer box label	Box label	A4	2025/03/25	DLB-W006(3)-01-03
	Pouch / Device label	Foil pouch	A4	2025/03/25	DAT-W006(3)-01-03
	Test cassette foil pouch	Test cassette foil pouch	A3	2023/03/31	DAV-W006(3)-01-01
	Accessory labeling (e.g., pipettes, buffer caps)	Same as W006P0058			

W006P0079	Outer box artwork	Box artwork	A1	2023/03/31	DLA-W006(3)-01-04
	Outer box label	Box label	A5	2025/03/25	DLB-W006(3)-01-04
	Pouch / Device label	Foil pouch	A4	2025/03/25	DAT-W006(3)-01-04
	Test cassette foil pouch	Test cassette foil pouch	A3	2023/03/31	DLC-W006(3)-01-04
	Accessory labeling (e.g., pipettes, buffer caps)	Same as W006P0058			
/	Instructions for Use	IFU	A8	2025/05/28	IFU-W006(3)-01-01


Labels

1.1 Labels for product code W006P0058


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXXXX which is nine characters.

Front

Back

Foil Pouch for W006P0058


Note: The lot number have three formats. Format 1 is WXXXXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXXX which is nine characters.


Test Cassette Foil Pouch for W006P0058


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front


Back

Disposal Bag for W006P0058


Round Sticker

1.2 Labels for product code W006P0059

Note: The lot number have three formats. Format 1 is WXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

The image displays a comprehensive vector illustration of a Wondfo IVD (In Vitro Diagnostic) box layout, designed for a 250 x 185 x 105mm product. The layout is presented as a flat, unfolded template, showing the front, back, and side panels. The design is clean and professional, featuring a white background with grey and red accents. The front panel prominently displays the Wondfo logo in blue, with the tagline "Leading POCT Manufacturer" below it. A large, stylized "IVD" label is positioned in the upper right corner. The back panel features the same Wondfo logo and tagline, but with a large, empty rectangular area for a barcode or additional labeling. The side panels are decorated with a series of parallel grey and red diagonal stripes. The top and bottom flaps of the box are also shown, with the top flap featuring a red and grey striped pattern. The layout includes various graphical elements such as a red outline of a box on the front panel, a yellow outline of a box on the back panel, and a red outline of a box on the side panel. The overall design is modern and functional, suitable for a medical or scientific product.

The image displays a comprehensive vector illustration of a Wondfo IVD (In Vitro Diagnostic) box layout, designed for a 250 x 185 x 105mm product. The layout is presented as a flat, unfolded template, showing the front, back, and side panels. The design is clean and professional, featuring a white background with grey and red accents. The front panel prominently displays the Wondfo logo in blue, with the tagline "Leading POCT Manufacturer" below it. A large, stylized "IVD" label is positioned in the upper right corner. The back panel features the same Wondfo logo and tagline, but with a large, empty rectangular area for a barcode or additional labeling. The side panels are decorated with a series of parallel grey and red diagonal stripes. The top and bottom flaps of the box are also shown, with the top flap featuring a red and grey striped pattern. The layout includes various graphical elements such as a red outline of a box on the front panel, a yellow outline of a box on the back panel, and a red outline of a box on the side panel. The overall design is modern and functional, suitable for a medical or scientific product.

Wondfo HIV Self-Test

EN

CONTENTS:

20 Individual pouches

each containing:

1 Test cassette

1 Desiccant pouch

1 Instruction for use

1 Buffer

1 Dropper

1 Alcohol prep pad

1 Blood Lancet for Single Use

1 Cotton swab

1 Disposal bag

(W006P0058, W006P0079 ONLY)

1 Plaster

(W006P0079 ONLY)

FR

CONTENUS:

20 Pochettes unitaires

contenant chacune :

1 Cassette de test

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1 Mode d'emploi

1 Solution tampon

1 Compte-gouttes

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1 Cotton-tige

1 Sac poubelle

(W006P0058, W006P0079 UNIQUEMENT)

1 Pansement

(W006P0079 UNIQUEMENT)

PT

CONTEÚDO:

20 bolsas individuais

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1 Cassete de teste

1 Bolsa dessecante

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1 Tampão

1 Conta-gotas

1 Toallete desinfetante

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1 Cotonete de algodão

1 Bolsa de descarte

(W006P0058, W006P0079 APENAS)

1 Pansamento


(W006P0079 APENAS)


Version: A5
Rel.: 2025/03/25

Secondary Packaging Label of W006P0059

Front

Back

Foil Pouch for W006P0059


Note: The lot number have three formats. Format 1 is WXXXXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXXX which is nine characters.


Test Cassette Foil Pouch for W006P0059


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front

Back

Round Sticker

1.3 Labels for product code W006P0060

	ART NO.: C/NO.: QTY: PCS N.W: KGS G.W: KGS MEAS: 36×28×21 CM		ART NO.: C/NO.: QTY: PCS N.W: KGS G.W: KGS MEAS: 36×28×21 CM

Secondary packaging of W006P0060

Wondfo HIV Self-Test

Wondfo
Leading POCT Manufacturer

IVD

Version: A4
Rel.: 2025/03/25

EN
CONTENTS:
100 Individual pouches
each containing:
1 Test cassette
1 Desiccant pouch
1 Instruction for use
1 Buffer
1 Dropper
1 Alcohol prep pad
1 Blood Lancet for Single Use
1 Cotton swab

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

Secondary Packaging Label for W006P0060

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front


Back

Foil Pouch for W006P0060


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.


Test Cassette Foil Pouch for W006P0060

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front

Back

1.4 Labels for product code W006P0079

250 x 185 x 105mm

Note: The lot number have three formats. Format 1 is WXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXXX which is nine characters.

Version: A1
Rel.: 2023/03/31

Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

IVD

IVD

IVD

IVD

REF W006P0079
LOT WXXXXXXXXX
YYYY-MM-DD
YYYY-MM-DD

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296063 Fax: 0086-20-32296063
Website: en.wondfo.com


Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296063 Fax: 0086-20-32296063
Website: en.wondfo.com

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296063 Fax: 0086-20-32296063
Website: en.wondfo.com

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296063 Fax: 0086-20-32296063
Website: en.wondfo.com

Weave Working For Your Health

The image displays a comprehensive vector illustration of a Wondfo IVD (In Vitro Diagnostic) box layout, designed for a 250 x 185 x 105mm product. The layout is presented as a flat, unfolded template, showing the front, back, and side panels. The design is clean and professional, featuring a white background with grey and red accents. The front panel prominently displays the Wondfo logo, the text 'Leading POCT Manufacturer', and a large 'IVD' label. It also includes a barcode area and a date field. The back panel features the Wondfo logo, the text 'Leading POCT Manufacturer', and a large 'IVD' label. It also includes a barcode area and a date field. The side panels are decorated with a repeating pattern of grey and red diagonal stripes. The top and bottom flaps of the box are also shown, with the top flap featuring a red and grey striped pattern. The overall design is modern and functional, suitable for a medical device packaging.

Wondfo HIV Self-Test

EN

CONTENTS:

20 Individual pouches

each containing:

1 Test cassette

1 Desiccant pouch

1 Instruction for use

1 Buffer

1 Dropper

1 Alcohol prep pad

1 Blood Lancet for Single Use

1 Cotton swab

1 Disposal bag

(W006P0058, W006P0079 ONLY)

1 Plaster

(W006P0079 ONLY)

FR

CONTENUS:

20 Pochettes unitaires

contenant chacune :

1 Cassette de test

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1 Mode d'emploi

1 Solution tampon

1 Compte-gouttes

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1 Cotton-tige

1 Sac poubelle

(W006P0058, W006P0079 UNIQUEMENT)

1 Pansement

(W006P0079 UNIQUEMENT)

PT

CONTEÚDO:

20 bolsas individuais

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1 Cassete de teste

1 Bolsa dessecante

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1 Tampão

1 Conta-gotas

1 Toalhete desinfetante

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1 Cotonete de algodão

1 Bolsa de descarte

(W006P0058, W006P0079 APENAS)

1 Pansamento


(W006P0079 APENAS)


Version: A5
Rel.: 2025/03/25

Secondary Packaging Label of W006P0079

Front

Back

Foil Pouch for W006P0079


Note: The lot number have three formats. Format 1 is WXXXXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXXX which is nine characters.


Test Cassette Foil Pouch for W006P0079

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front

Back

Sterile Wound Plaster



LOT


CE 0123

STERILE EO

Disposable Bag for W006P0079

Round Sticker

Instructions for use⁴

⁴ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.


Wondfo

HIV Self-Test

Instruction for use

- You must follow the test procedure carefully to get an accurate result.
- Sit in a clean, well-lit area and ensure you have all contents before beginning the test.
- For single-use only. Do not open foil pouch containing cassette until ready to test.

Contents

Watch the operation video

Materials may be required but not provided.



Timer


Tissue


Preparation

How to use the test kit (for fingerstick whole blood use)

How to read the test

Questions and answers

1. What is HIV?

The human immunodeficiency virus (HIV) targets cells of the immune system, called CD4 cells, which help the body respond to infection. Within the CD4 cell, HIV replicates and in turn, damages and destroys the cell. Without effective treatment of a combination of antiretroviral (ARV) drugs, the immune system will become weakened to the point that it can no longer fight infection and disease.

2. How is HIV transmitted?

HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:

- Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
- Blood transfusion of contaminated blood.
- Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments.
- From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.

If a person living with HIV is on ART (Antiretroviral therapy), which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

3. When do I need to test myself?

People who have been exposed to HIV or are at risk of HIV should seek testing. Most widely-used HIV diagnostic tests detect antibodies produced by the person as part of their immune response to fight HIV. In most cases, people develop antibodies to HIV within 28 days of infection. During this time, people may receive a false-negative result using an antibody test. While people at risk of HIV should test as soon as possible, negative test results within 28 days of exposure should be confirmed with additional testing no more than three months later.

4. What should I do if I get a positive result?

You need to follow up with a health care worker to get additional testing to confirm the result. At that time your local clinic and the provider will suggest the next steps that need to be taken.

Product information

CATALOG NO.

W006P0058 W006P0059 W006P0060 W006P0079

INTENDED USE

The Wondfo HIV Self-Test is a single-use *in vitro* diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as self-test and/or by medical professionals.

For *in vitro* diagnostic use only.

PRINCIPLES OF THE PROCEDURE

Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into the respective wells, they will migrate along the device by capillary action. The HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen (gp36/41), and the complex is then captured by the HIV recombinant antigen (gp41 and gp36) immobilized in the test region (T). When the levels of HIV-1 or HIV-2 antibodies are at or above the limit of detection (LOD) of the assay, it will produce a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV-1 or HIV-2 antibodies are zero or below the LOD, there is no visible colored band in the test region (T) indicating a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C).

WARNINGS AND CAUTIONS

Do not use if the test kit beyond expiration date.
Do not use if the pouch is punctured or improperly sealed.
Do not use for self-testing if you are under 12 years old.
Do not use for self-testing if you have a bleeding disorder.
Do not use for self-testing if you are already diagnosed as HIV positive.
Do not open the pouch until you are ready to perform the test.
Wash your hands and ensure that they are clean and dry before starting test.
Adequate lighting is required to read the test results.

KIT CONTENT

There are 4 configurations of the test kits, 1 test/kit, 20 tests/kit and 100 tests/kit. The kit components are provided as below:

Catalog No.		W006P0058	W006P0059	W006P0060	W006P0079
Components					
Accessories	Test cassette (pcs)	1	1x20	1x100	1x20
	Desiccant (pcs)	1	1x20	1x100	1x20
	Dropper (pcs)	1	1x20	1x100	1x20
	Buffer (vial)	1	1x20	1x100	1x20
	IFU (pcs)	1	1x20	1x100	1x20
	Blood Lancet for Single Use (pcs)	1	1x20	1x100	1x20
	Alcohol prep pad (pcs)	1	1x20	1x100	1x20
	Cotton swab (pcs)	1	1x20	1x100	1x20
	Disposal bag (pcs)	1	/	/	1x20
	Plaster (pcs)	/	/	/	1x20

One test strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

STORAGE AND STABILITY

- The test kit can be stored at 2-30 °C until the stated Expiry Date on the kit.
- Use the test cassette within 1 hour after opening the pouch.
- Keep away from sunlight, moisture and heat.
- Use the kit at 10-30 °C.

LIMITATIONS OF THE PROCEDURE

- The test is designed for testing human fingerstick whole blood.
- The test is limited to the qualitative detection of HIV-1 and HIV-2 antibodies.
- The assay procedure and result interpretation must be followed closely when testing. For optimal test performance, proper specimen collection is critical. Failure to follow the procedure may lead to inaccurate test results.
- False-negative results can occur in the following conditions:
 - Patients exposed to HIV less than 3 months.
 - Patients under HIV treatment (Antiretroviral therapy).
 - If the quantity of antibodies for HIV present in the specimen is below the detection limit of the assay.
- False-positive results can occur in the following conditions:
 - Patients have participated in a HIV vaccine clinical trial.
- The presence of the control line only means that migration of added liquid occurred. It does not guarantee that:
 - The correct specimen has been used.
 - The specimen has been applied correctly.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but instead should be determined by a healthcare provider in conjunction with clinical findings and the results from other laboratory tests and evaluations. Results from the Wondfo HIV Self-Test should not be used as the sole basis for diagnosis.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity:

A clinical study of the Wondfo HIV Self-Test was conducted in 5 different test institutions. Healthcare professionals tested 3911 samples with the Wondfo HIV Self-Test. Laboratory testing confirmed that 1172 samples were HIV positive and 2739 were HIV negative.

- Sensitivity: 100% (95%CI: 99.67% - 100%) - the Wondfo HIV Self-Test correctly identified all 1172 HIV positive samples.
- Specificity: 99.96% (95%CI: 99.79% - 99.99%) - the Wondfo HIV Self-Test correctly identified 2738 of the 2739 HIV negative samples.

Usability Study:

A usability study was conducted in which 900 untrained lay participants whose HIV status was previously unknown tested themselves using the Wondfo HIV Self-Test. The participants' HIV true status was determined by testing with a 4th generation (HIV Ag/Ab) laboratory test.

- 42 (4.7%) of the 900 participating lay-users, failed to obtain a Self-Test result, including 5 whose laboratory test was HIV positive. An additional participant's 4th generation HIV test status was not determined before the study closed.
- 69 (95.8%) of the remaining 72 laboratory confirmed HIV positive lay-users correctly reported a positive Self-Test.
- 782 (99.6%) of the remaining 785 laboratory confirmed HIV negative lay-users correctly reported a negative Self-Test.

REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or selftesting, Geneva: World Health Organization; 2016.

SYMBOLS KEY

	<i>In vitro</i> diagnostic medical device		Consult instructions for use		Use-by date
	Contains sufficient for <n> tests		Date of manufacture		Keep dry
	Batch code		Temperature limit		Keep away from sunlight
	Manufacturer		Do not re-use		Catalogue number
	Caution		Do not use if package is damaged and consult instructions for use		

Manufacturer information

Guangzhou Wondfo Biotech Co., Ltd.

Add: No.8 Lizhishan Road, Science City, Huangpu District, 510663,Guangzhou, P.R. China

Tel: +86-20-32296083 400-888-5268 (Toll Free)

Fax: +86-20-32296063

E-mail: global@wondfo.com.cn

Website: en.wondfo.com

Please contact the manufacturer or your local distributor if you have any questions related to the product.

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RÉFÉRENCES

- WHO. TGS-5 Designing instruction for use for in vitro diagnostic medical devices, Geneva: World Health Organization ; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or selftesting. Geneva: World Health Organization; 2016.

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