WHO Prequalification of In Vitro Diagnostics PUBLIC ASSESSMENT REPORT

Product: Wondfo Malaria P.f (HRP2/pLDH) Test WHO reference number: PQDx 0651-004-00

Wondfo Malaria P.f (HRP2/pLDH) Test with product codes W037P0041, W037P0022, W037P0028, W037P0030, W037P0024, W037P0033, W037P0035, and W037P0037, manufactured by Guangzhou Wondfo Biotech Co., Ltd, Rest-of-World version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 30 July 2024.

Summary of WHO prequalification assessment for the Wondfo Malaria P.f (HRP2/pLDH) Test

	Date	Outcome
Prequalification listing	30 July 2024	listed
Dossier assessment	29 July-2024	MR
Site inspection(s) of quality management system	13 February 2024	MR
Product performance evaluation	Quarter 4 2022	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Reason for amendment and change request reference, where applicable.	Date of report amendment
2.0	<p>The manufacturer has proposed the following changes:</p> <ol style="list-style-type: none"> 1. Changes to the accessories, including: <ul style="list-style-type: none"> • replacement of lancet supplier; • change to the manufacturing address of the inverted cup supplier; 2. Addition of languages to IFU and labels. 3. Change to the pack size of the test kit (addition of 	19 March 2026

	configuration 10 tests/kit)(PQC-IVD-2025-0067).	
--	---	--

Intended use

According to the intended use claim from Guangzhou Wondfo Biotech Co., Ltd, “*The Wondfo Malaria P.f (HRP2/pLDH) Test is an immunochromatographic, rapid diagnostic assay intended for the qualitative detection of P. falciparum specific histidine rich protein-2 (P.f HRP2) and P. falciparum specific pLDH (P.f pLDH) in human fingerstick whole blood and venous whole blood. It is intended to be used by laboratory professionals and at point-of-care by healthcare workers or trained lay providers. It is used as an aid for diagnosis of P. falciparum Malaria infection for symptomatic patients including paediatric testing. Results from the Wondfo Malaria P.f (HRP2/pLDH) Test should not be used as the sole basis for diagnosis.*

For in vitro diagnostic use only. For professional use only.”

Test kit contents

Component	10 Test (T) (W037P0041)	25 T (W037P0022)	40 T (W037P0024)	25 T (W037P0028)	25 T (W037P0030)	40 T (W037P0033)	40 T (W037P0035)	100 T (W037P0037)
Cassette (pieces (pcs))	10	25	40	25	25	40	40	100
Instructions	1	1	1	1	1	1	1	1
Buffer	1x5mL/bottle	25x0.3 mL/vial	40x0.3 mL/vial	2x5 mL/bottle	2x5 mL/bottle	4x5 mL/bottle	4x5 mL/bottle	8x5 mL/bottle
Quick guide (pcs)	/	25 pouches, each containing one quick guide:	40 pouches, each containing one quick guide:	/	/	/	/	/
Inverted cup (pcs)	10	one inverted cup,	one inverted cup,	25	25	40	40	100
Alcohol prep pad (pcs)	10	one alcohol prep pad,	one alcohol prep pad,	25	25	40	40	100
Sterile lancet (pcs)	10	one sterile lancet.	one sterile lancet.	25	/	40	/	100
Blood lancet for single use (pcs)	/	/	/	/	25	/	40	/

Items required but not provided

- Timer
- Protective gloves
- Specimen and test waste container

- Sterile cotton swab/ball

Note: If whole blood is collected by venipuncture, venipuncture blood collection materials are required in addition to a precision pipette and pipette tips.

Storage

The test kit must be stored between 1 and 40 °C.

Shelf-life upon manufacture¹

24 months.

Warnings/limitations

Refer to the current version of the manufacturer's instructions for use attached to this public assessment report.

The WHO Global Malaria Programme recommends the use of tests detecting pLDH and a panel detection of at least 75% on the panel of non-HRP2 expressing specimens in areas with high prevalence of *pfhrp2/3* gene deletions². This test does not meet these criteria and is not recommended for use in such settings.

Dossier assessment

Guangzhou Wondfo Biotech Co., Ltd submitted a product dossier for the Wondfo Malaria P.f (HRP2/pLDH) Test as per the "Instructions for compilation of a product dossier" (PQDx_018). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO. The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 29 July 2024.

Based on the product dossier screening and assessment findings, the product dossier for the Wondfo Malaria P.f (HRP2/pLDH) Test meets WHO prequalification requirements.

Manufacturing site inspection

¹ The assigned device shelf-life is based on stability data generated from the date of manufacture. The finished goods shelf-life, calculated from the date of packaging completion, may be shorter depending on the time elapsed between manufacture and final packaging of the device.

² <https://iris.who.int/bitstream/handle/10665/258972/WHO-HTM-GMP-2017.18-eng.pdf?sequence=1>

The inspection of the manufacturing site(s) was conducted to assess whether the manufacturer's quality management system (QMS) and manufacturing practices are in alignment with:

- (i) applicable international standards, such as ISO 13485 (Medical devices – Quality management systems – Requirements for regulatory purposes);
- (ii) the manufacturer's own documented procedures and quality requirements; and
- (iii) other relevant international standards and guidelines applicable to in vitro diagnostic (IVD) medical devices. The WHO's Public Inspection Reports are accessible at:

<https://extranet.who.int/pqweb/vitro-diagnostics/who-public-inspection-reports>

Based on the site inspection and corrective action plan review, the quality management system for the Wondfo Malaria P.f (HRP2/pLDH) Test meets WHO prequalification requirements.

Product performance evaluation

Wondfo Malaria P.f (HRP2/pLDH) Test was evaluated by the Centers for Disease Control and Prevention on behalf of WHO in the fourth quarter of 2022, according to protocol PQDx_317, version 3.1.

Clinical performance evaluation

Wondfo Malaria P.f (HRP2/pLDH) Test was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild type parasite panel, *P. vivax* wild type parasite panel, *hrp2/3* deletion panel, a parasite-negative panel (including clean negatives from afebrile donors and specimens with possible cross-reactants), and the WHO International standard for *P. falciparum* antigens.

Performance characteristics			
	<i>P. falciparum</i>	<i>P. vivax</i>	Pf – <i>hrp2</i> deletion panel
Panel detection score [#] at 200 parasites/μL	93/100, 93.0 %	NA	21/60, 35.0%
Proportion detected by line at 200 parasites/μL	HRP2 line: 98.0% Pf-pLDH line: 55.0% (N=400)	NA	HRP2 line: 16.7% Pf-pLDH line: 52.5% (N=240)
False positive results %	All negative specimens: 0/200, 0.0% Of which, clean negative specimens: 0/104, 0% No false Pf result on the <i>P. vivax</i> specimens at 200 and 2000 parasites/μL		
Invalid rate %	0/1322, 0.0 %		

(N= 1322)	
Inter-reader variability on the wild-type and negative panels % (N= 1010)	HRP2 test line: 1/1010, 0.1 % Pf-pLDH test line: 166/1010, 16.4%
Lowest concentration of HRP2/pLDH detected using the 1 st WHO International standard for Pf antigens (NIBSC code: 16/376)*	HRP2 test line: Lot A: 7.8 IU/mL Lot B: 3.9 IU/mL Pf-pLDH test line: Lot A: 15.6 IU/mL Lot B: 15.6 IU/mL

The panel detection score is the proportion of specimens at a specified concentration that gave positive results in all 4 tests performed with this specimen (2 tests on each of 2 lots).

* Not applicable for assays evaluated in WHO product testing of RDTs for malaria antigen detection.

Operational characteristics and ease of use

This assay does not require laboratory equipment and can be performed in laboratories with limited facilities or in non-laboratory settings.

The assay was found easy to use by the operators performing the evaluation.

However, the operators noted that 383/1010 (37.9 %) of the tests performed on the wild-type and negative panels showed anomalies (234 with incomplete clearing, 38 with red background, 51 with red background obscuring test line, 58 with patchy line, 1 with strip misplaced, and 1 with failed migration).

Key operational characteristics	
Specimen type and volume	5 µL of fingerstick whole blood or venous whole blood collected into a blood collection tube with anticoagulant (EDTA, sodium citrate, or heparin sodium)
Number of steps*	2 steps in total 1 step with specimen transfer device (precision pipette was used during the evaluation)
Time to result	15 minutes
Endpoint stability (interval)	15 minutes (the test can be read between 15 and 30 minutes after addition of diluent)
Internal QC	Yes, reagent addition control

* Definition: each action required to obtain a result (excluding specimen collection, device preparation – opening the pouch), e.g. for RDTs: add specimen, add buffer (2 steps).

Based on these results, the performance evaluation for the Wondfo Malaria P.f (HRP2/pLDH) Test meets the WHO prequalification requirements.

Labelling review

The labelling submitted for the Wondfo Malaria P.f (HRP2/pLDH) Test was reviewed by WHO staff and external technical experts appointed by WHO. The review evaluated the labelling for clarity and consistency with the information submitted in the product dossier, alignment with international guidance and standards, and suitability for the intended users and settings in WHO Member States, including low- and middle-income countries.

The table below provides traceability of the labelling documents reviewed during the assessment, including document titles, version numbers, approval dates, and control identifiers.

Controlled Labelling References

Document Type	Document Title	Version / Revision	Date Approved	Controlled Document No.	Applicable Catalog No.	
Outer box artwork	Box artwork	A3	2024/06/25	DLA-W037(2)-01-01	W037P0022	
				DLA-W037(2)-01-02	W037P0024	
				DLA-W037(2)-01-03	W037P0028	
				DLA-W037(2)-01-04	W037P0030	
				DLA-W037(2)-01-05	W037P0033	
				DLA-W037(2)-01-06	W037P0035	
				DLA-W037(2)-01-07	W037P0037	
		A1	2025/06/06	DLA-W037(2)-01-08	W037P0041	
Outer box label	Box label	A3	2025/06/06	DLB-W037(2)-01-01	W037P0022	
				DLB-W037(2)-01-02	W037P0024	
				DLB-W037(2)-01-03	W037P0028	
				DLB-W037(2)-01-05	W037P0033	
				DLB-W037(2)-01-07	W037P0037	
				A4	DLB-W037(2)-01-04	W037P0030
				DLB-W037(2)-01-06	W037P0035	
		A1	DLB-W037(2)-01-08	W037P0041		
Pouch / Device label	Pouch	A3	2024/06/25	DLC-W037(2)-01-01	All	
Accessory labeling	Accessory pouch	A3	2025/06/06	DAV-W037(2)-01-01	W037P0022	
				DAV-W037(2)-01-02	W037P0024	
Accessory labeling (buffer)	Buffer label	A2	2024/06/25	DLD-W037(2)-01-01	W037P0022	
				DLD-W037(2)-01-02	W037P0024	
				DLD-W037(2)-01-03	W037P0028	
					W037P0030	
					W037P0033	
					W037P0035	
					W037P0037	

					W037P0041
Accessory labeling (inverted cup)	Dropper label	A3	2025/06/06	DLX-W037(2)-01-01	W037P0028
		A1		DLX-W037(2)-01-02	W037P0037
		A3		DLX-W037(2)-01-03	W037P0041
				DLX-W037(2)-01-04	W037P0030
				DLX-W037(2)-01-05	W037P0033
Accessory labeling (lancet labels)	Lancet	A2	2024/12/11	DLV-W037(2)-01-01	W037P0035
				DLV-W037(2)-01-02	W037P0028
		A0	2025/07/01	DLV-W037(2)-01-03	W037P0037
				DLV-W037(2)-01-04	W037P0033
				DLV-W037(2)-01-05	W037P0030
A1	2025/06/06	DLV-W037(2)-01-05	W037P0035		
Quick Guide	Quick Guide	A2	2024/9/26	DLI-W037(2)-01-01	W037P0041
					W037P0022
Instructions for Use (IFU)	IFU	A7	2025/6/6	IFU-W037(2)-01-01	W037P0024
					All

Labels

Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 25 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 25 Vials of P.f buffer (0.3mL/vial)
3. 1 Instruction for use
4. 25 Alcohol prep pads
5. 25 Sterile lancets
6. 25 Inverted cups
7. 25 Quick guides

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



Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 40 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 40 Vials of P.f buffer (0.3mL/vial)
3. 1 Instruction for use
4. 40 Alcohol prep pads
5. 40 Sterile lancets
6. 40 Inverted cups
7. 40 Quick guides

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
Version: A3
Rel.:2025/06/06

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8, Lixin Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32290663
Website: en.wondfo.com

Wondfo
Leading POCT Manufacturer

FOR IN VITRO DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8, Lixin Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32290663
Website: en.wondfo.com

Version: A3
Rel.: 2024/06/25

IVD

Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

FOR IN VITRO DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

REF: W037P0024
LOT: XXXXXXXX
YYMM-DD
YYMM-DD

We Are Working For Your Health

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Wondfo
Leading POCT Manufacturer

Wondfo Malaria P.f (HRP2/pLDH) Test

LOT XXXXXXXX

YYYY-MM-DD

REF W037P0024

YYYY-MM-DD

Version: A3

FOR IN VITRO DIAGNOSTIC USE ONLY

Rel.: 2024/06/25

Wondfo

IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

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


Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 25 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 2 Bottles of buffer (5mL/bottle)
3. 1 Instruction for use
4. 25 Alcohol prep pads
5. 25 Sterile lancets
6. 25 Inverted cups

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Wondfo
Leading POCT Manufacturer

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32293962 Fax: +86-20-32296663
Website: en.wondfo.com

Version: A3
Rel.: 2024/06/25

Wondfo
Leading POCT Manufacturer

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

REF W037P0028
LOT XXXXXXXX
YYYY-MM-DD
YYYY-MM-DD

We Are Working For Our Health

IVD

Wondfo
Leading POCT Manufacturer

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Wondfo
Leading POCT Manufacturer

Wondfo Malaria P.f (HRP2/pLDH) Test

LOT XXXXXXXX

⌚ YYYY-MM-DD

REF W037P0028

📄 YYYY-MM-DD

Version: A3

Rel.: 2024/06/25

FOR IN VITRO DIAGNOSTIC USE ONLY

Wondfo

IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: +86-20-32053962 Fax: +86-20-32296063

Website: en.wondfo.com

Sterile Lancet

25
pcs

Rev. A2
Rel.: 2024/12/11

Beijing Ruicheng Medical Supplies Co., Ltd.
Building 5, No. 8, Yanqi West Road, Yanqi
Economic Development Zone, Huairou District
101407 Beijing P.R. China

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd, 2595AA, The
Hague, Netherlands.

STERILE R

LOT

Manufactured for Wondfo Biotech by Beijing Ruicheng Medical Supplies Co., Ltd.

Inverted cup

 Foshan Yuyang Medical Instrument Co.,Ltd.
Room 501,5th Floor, Building B, No. 1 Kesheng Road,Lishui
Town,Nanhai District, Foshan City,Guangdong, China

LOT

Manufactured for Wondfo Biotech by
Foshan Yuyang Medical Instrument Co.,Ltd.

Language: English
Rev.A3 Ref.20250608

size:38*20mm

**Buffer for Wondfo Malaria
P.f (HRP2/pLDH) Test**

Lot: XXXXXXXX

Ref: W037P0028

Exp: YYYY-MM-DD

Vol:5ml

Store at:1-40°C

Version:A2

Rel.:2024/06/25

size:38*20mm

Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 25 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 2 Bottles of buffer (5mL/bottle)
3. 1 Instruction for use
4. 25 Alcohol prep pads
5. 25 Blood lancets for single use
6. 25 Inverted cups

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Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

REF W037P0030
LOT XXXXXXXX
YYYY-MM-DD
YYYY-MM-DD

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32293962 Fax: +86-20-32296683
Website: en.wondfo.com

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32293962 Fax: +86-20-32296683
Website: en.wondfo.com

Version: A3
Rel.: 2024/06/25

WeAreWorkingForOur Health

IVD

Wondfo
Leading POCT Manufacturer

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Wondfo
Leading POCT Manufacturer

Wondfo Malaria P.f (HRP2/pLDH) Test

LOT XXXXXXXX

 YYYY-MM-DD

REF W037P0030

 YYYY-MM-DD

Version: A3

Rel.: 2024/06/25

FOR IN VITRO DIAGNOSTIC USE ONLY

Wondfo

IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: +86-20-32053962 Fax: +86-20-32296063

Website: en.wondfo.com

Lancet (21G)

Blood Lancet for Single Use

25
pcs

Hebei Xinle Sci&Tech CO.,LTD
No. 2, Xingye Street, Xinle City, 050700 Shijiazhuang City,
Hebei province, PEOPLE'S REPUBLIC OF CHINA
Tel:0086-311-88655818 Fax:0086-311-88655819

Rev.A0
Rel.:2025/07/01

MD STERILE R

EU REP

Shanghai International Holding Corp.GmbH(Europe)
Eiffestrasse 80, 20537 Hamburg,Germany
Tel:0049-40-2513175 Fax:0049-40-255726
E-mail:shholding@hotmail.com

⊗ ⊗ CE 0123

Model: XLII-2118

LOT

Manufactured for Wondfo Biotech by Hebei Xinle Sci&Tech CO.,LTD

www.hbxinle.com

Inverted cup

 Foshan Yuyang Medical Instrument Co.,Ltd.
Room 501,5th Floor, Building B, No. 1 Kesheng Road,Lishui
Town,Nanhai District, Foshan City,Guangdong, China

LOT

Manufactured for Wondfo Biotech by
Foshan Yuyang Medical Instrument Co.,Ltd.

Language: English
Rev.A3 Ref.20250608

size:38*20mm

**Buffer for Wondfo Malaria
P.f (HRP2/pLDH) Test**

Lot: XXXXXXXX

Ref: W037P0030

Exp: YYYY-MM-DD

Vol:5ml

Store at:1-40°C

Version:A2

Rel.:2024/06/25

size:38*20mm

Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 40 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 4 Bottles of buffer (5mL/bottle)
3. 1 Instruction for use
4. 40 Alcohol prep pads
5. 40 Sterile lancets
6. 40 Inverted cups

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
Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishiwan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32295063
Website: en.wondfo.com

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishiwan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32295063
Website: en.wondfo.com

Version: A3
Rel.: 2024/06/25

Wondfo
Leading POCT Manufacturer

We're Working For Your Health

IVD

REF	W037P0033
LOT	XXXXXXXX
	YYYY-MM-DD
	YYYY-MM-DD

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Wondfo
Leading POCT Manufacturer

Wondfo Malaria P.f (HRP2/pLDH) Test

LOT XXXXXXXX

YYYY-MM-DD

REF W037P0033

YYYY-MM-DD

Version: A3

Rel.: 2024/06/25

FOR IN VITRO DIAGNOSTIC USE ONLY

Wondfo

IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: +86-20-32053962 Fax: +86-20-32296063

Website: en.wondfo.com

Sterile Lancet

40
pcs

Rev. A2
Rel.: 2024/12/11

Beijing Ruicheng Medical Supplies Co., Ltd.
Building 5, No. 8, Yanqi West Road, Yanqi
Economic Development Zone, Huairou District
101407 Beijing P.R. China

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd, 2595AA, The
Hague, Netherlands.

CE

0197

STERILE R

LOT

Manufactured for Wondfo Biotech by Beijing Ruicheng Medical Supplies Co., Ltd.

Inverted cup

 Foshan Yuyang Medical Instrument Co., Ltd.
Room 501, 5th Floor, Building B, No. 1 Kesheng Road, Lishui
Town, Nanhai District, Foshan City, Guangdong, China

LOT

Manufactured for Wondfo Biotech by
Foshan Yuyang Medical Instrument Co., Ltd.

Language: English
Rev A3 (Rev. 2025/06/06)

size:38*20mm

**Buffer for Wondfo Malaria
P.f (HRP2/pLDH) Test**

Lot: XXXXXXXX

Ref: W037P0033

Exp: YYYY-MM-DD

Vol:5ml

Store at:1-40°C

Version:A2

Rel.:2024/06/25

size:38*20mm

Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 40 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 4 Bottles of buffer (5mL/bottle)
3. 1 Instruction for use
4. 40 Alcohol prep pads
5. 40 Blood lancets for single use
6. 40 Inverted cups

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
Wondfo
Leading POCT Manufacturer

Wondfo
Leading POCT Manufacturer

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

FOR *IN VITRO* DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishiwan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32296063
Website: en.wondfo.com

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishiwan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32296063
Website: en.wondfo.com

Version: A3
Rel.: 2024/06/25

We're Working For You
Health

IVD

Leading POCT Manufacturer
Wondfo

REF	W037P0035
LOT	XXXXXXXX
	YYYY-MM-DD
	YYYY-MM-DD

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Wondfo
Leading POCT Manufacturer

Wondfo Malaria P.f (HRP2/pLDH) Test

LOT XXXXXXXX

XXXX-XX-XX

REF W037P0035

XXXX-XX-XX

Version: A3

Rel.: 2024/06/25

FOR IN VITRO DIAGNOSTIC USE ONLY

Wondfo

IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: +86-20-32053962 Fax: +86-20-32296063

Website: en.wondfo.com

Lancet (21G)

Blood Lancet for Single Use

40
pcs

Hebei Xinle Sci&Tech CO.,LTD
No. 2, Xingye Street, Xinle City, 050700 Shijiazhuang City,
Hebei province, PEOPLE'S REPUBLIC OF CHINA
Tel:0086-311-88655818 Fax:0086-311-88655819

Rev.A0
Rel.:2025/07/01

MD STERILE R

EU REP

Shanghai International Holding Corp.GmbH(Europe)
Eiffestrasse 80, 20537 Hamburg,Germany
Tel:0049-40-2513175 Fax:0049-40-255726
E-mail:shholding@hotmail.com

⊗ ⊗ CE 0123

Model: XLII-2118

LOT

Manufactured for Wondfo Biotech by Hebei Xinle Sci&Tech CO.,LTD

www.hbxinle.com

Inverted cup

 Foshan Yuyang Medical Instrument Co., Ltd.
Room 501, 5th Floor, Building B, No. 1 Kesheng Road, Lishui
Town, Nanhai District, Foshan City, Guangdong, China

LOT

Manufactured for Wondfo Biotech by
Foshan Yuyang Medical Instrument Co., Ltd.

Language: English
Rev A3 Ref: 2025/06/06

size:38*20mm

**Buffer for Wondfo Malaria
P.f (HRP2/pLDH) Test**

Lot: XXXXXXXX

Ref: W037P0035

Exp: YYYY-MM-DD

Vol:5ml

Store at:1-40°C

Version:A2

Rel.:2024/06/25

size:38*20mm

Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 100 Individual pouches, each containing:
 - Test cassette
 - Desiccant pouch
2. 8 Bottles of buffer (5mL/bottle)
3. 1 Instruction for use
4. 100 Alcohol prep pads
5. 100 Sterile lancets
6. 100 Inverted cups

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


Version: A3
Rel.: 2025/06/06

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32296063
Website: en.wondfo.com

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China
Tel: +86-20-32053962 Fax: +86-20-32296063
Website: en.wondfo.com

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Wondfo
Leading POCT Manufacturer

Wondfo Malaria P.f (HRP2/pLDH) Test

LOT XXXXXXXX
REF W037P0037

YYYY-MM-DD
YYYY-MM-DD

Version: A3

Rel.: 2024/06/25

FOR IN VITRO DIAGNOSTIC USE ONLY

Wondfo

IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: +86-20-32053962 Fax: +86-20-32296063

Website: en.wondfo.com

Sterile Lancet

25
PCS

Rev. A2
Rel.: 2024/12/11

Beijing Ruicheng Medical Supplies Co.,Ltd.
Building 5, No. 8, Yanqi West Road, Yanqi
Economic Development Zone, Huairou District
101407 Beijing P.R. China

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd, 2595AA, The
Hague, Netherlands.

0197
STERILE R

Manufactured for Wondfo Biotech by Beijing Ruicheng Medical Supplies Co., Ltd.

Inverted cup

Foshan Yuyang Medical Instrument Co.,Ltd.
Room 501,5th Floor, Building B, No. 1 Kesheng Road,Lishui
Town,Nanhai District, Foshan City,Guangdong, China

LOT

Manufactured for Wondfo Biotech by
Foshan Yuyang Medical Instrument Co.,Ltd.

Language: English
Rev.A3 Ref.20250608

size:38*20mm

**Buffer for Wondfo Malaria
P.f (HRP2/pLDH) Test**

Lot: XXXXXXXX

Ref: W037P0037

Exp: YYYY-MM-DD

Vol:5ml

Store at:1-40°C

Version:A2

Rel.:2024/06/25

size:38*20mm

Wondfo Malaria P.f (HRP2/pLDH) Test

CONTENTS:

1. 10 Individual pouches, each containing:
Test cassette
Desiccant pouch
2. 1 Bottles of buffer(5mL/bottle)
3. 1 Instruction for use
4. 10 Alcohol prep pads
5. 10 Sterile lancets
6. 10 Inverted cups

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
IVD

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: +86-20-32053962 Fax: +86-20-32296063

Website: en.wondfo.com

75mm

10mm

95mm

Sterile Lancet

10
pcs

Rev. A1
Rel.: 2025/06/06

Beijing Ruicheng Medical Supplies Co.,Ltd.
Building 5, No. 8, Yanqi West Road, Yanqi
Economic Development Zone, Hualou
District 101407 Beijing P.R. China

Lotus NL B.V.
Koningin Julianaplein 10, le Verd, 2595AA,
The Hague, Netherlands

LOT

CE

0197

STERILE | R

Manufactured for Wondfo Biotech by Beijing Ruicheng Medical Supplies Co., Ltd.

Inverted cup

 Foshan Yuyang Medical Instrument Co.,Ltd.
Room 501,5th Floor, Building B, No. 1 Kesheng Road,Lishui
Town,Nanhai District, Foshan City,Guangdong, China

LOT

Manufactured for Wondfo Biotech by
Foshan Yuyang Medical Instrument Co.,Ltd.

Language: English
Rev:A1 Ref:2025/09/06

size:38*20mm

**Buffer for Wondfo Malaria
P.f (HRP2/pLDH) Test**

Lot: XXXXXXXX

Ref: W037P0028

Exp: YYYY-MM-DD

Vol:5ml

Store at:1-40°C

Version:A2

Rel.:2024/06/25

size:38*20mm

Instructions for use³

³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

Wondfo Malaria P.f (HRP2/pLDH) Test

CATALOG No.

W037P0022	W037P0028	W037P0030
W037P0024	W037P0033	W037P0035
W037P0037	W037P0041	


INTENDED USE

The Wondfo Malaria P.f (HRP2/pLDH) Test is an immunochromatographic, rapid diagnostic assay intended for the qualitative detection of *P. falciparum* specific histidine rich protein-2 (P.f HRP2) and *P. falciparum* specific pLDH (P.f pLDH) in human fingerstick whole blood and venous whole blood. It is intended to be used by laboratory professionals and at point-of-care by healthcare workers or trained lay providers. It is used as an aid for diagnosis of *P. falciparum* Malaria infection for symptomatic patients including paediatric testing. Results from the Wondfo Malaria Pf (HRP2/pLDH) Test should not be used as the sole basis for diagnosis.


For in vitro diagnostic use only. For professional use only.

BACKGROUND INFORMATION


Malaria remains one of the most serious tropical diseases in many regions of the world, mainly occurring in tropical and subtropical areas. Malaria is caused by a parasite that commonly infects the *Anopheles* mosquito which feeds on human blood. There are four kinds of malaria parasites that can infect humans: *P. falciparum*, *P. vivax*, *P. malariae*,

2. Test procedure for venous whole blood
A. Please gently mix the venous blood specimen by inverting the blood collection tube to make it homogenous. Take 5µL venous whole blood specimen with a pipette and drop it into the center of the sample well vertically.

B. Add 4 drops (about 100µL) of buffer into the buffer well vertically, and start timing.

3. Interpretation
A. Read the results after 15 minutes but no more than 30 minutes.

B. To avoid confusion, discard the test device after reading the result.

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WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* use only.
- Apply standard biosafety precautions for handling and disposal of potentially infective material.
 - Handle all specimens as potentially infectious.
 - Wear gloves while handling specimens and performing the test.
 - Avoid splashing and aerosol formation.
 - Clean up spills thoroughly using an appropriate disinfectant.
- Do not** reuse the test cassette, inverted cup, lancet, and alcohol prep pad.
- Do not** use the kit beyond the expiration date.
- Do not** use the test cassette if the pouch is punctured or improperly sealed.
- Do not** mix a buffer and a test cassette from different lots.
- Do not** open the sealed pouch until you are ready to conduct the assay. Once opened, use the test cassette within 1 hour.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials as biohazardous waste.
- Read the instructions for use completely before using the kit. Failure to follow directions may yield inaccurate test results.
- DISPOSAL OF DIAGNOSTIC DEVICE: The used device has risk of infection. The process for disposing the diagnostic device must follow local infectious waste disposal laws or laboratory regulations.
- Keep out of reach of children.
- The Buffer contains 2% Triton X-100, which presents no hazard to the user if normal laboratory safety precautions are followed. Wear protective gloves. In case of contact with skin


QUALITY CONTROL

- Internal control: This kit contains a built-in control feature as a procedural control, the control region (C). A colored line develops in the control region (C) after addition of the specimen and buffer. If a colored line does not develop in the C region, the test is invalid. Review the procedure and repeat the test with a new device.
- External control: The test kit does not provide external control.

INTERPRETATION OF RESULTS

Lines that you see	Picture	Read and Interpretation of Results
No line at 'C' (control)		Result: Invalid. Action: Take a new cassette and repeat the test. Please contact the local distributor or manufacturer if invalid results continuously appear.
Line at 'C' and No other line		Result: Negative. Action: No

Lines that you see	Picture	Read and Interpretation of Results
Line at 'C' and at 'T1' and/or 'T2'		Result: Positive. 1. Line at 'C', 'T1' and 'T2' indicates positive for P.f HRP2 and P.f pLDH. 2. Line at 'C' and 'T1' indicates positive for P.f HRP2. 3. Line at 'C' and 'T2' indicates positive for P.f pLDH. Note: Line intensity may vary from faint to strong intensity. Consider a faint test line as a positive result. Action: No

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

QUICK GUIDE


For fingerstick whole blood

Bring all specimens and kit components to room temperature (15°C- 40°C) before testing. When ready to test, remove the test cassette from the foil pouch by tearing the notch and place it on a flat, dry surface. Label the device with specimen ID number.



Material provided


- 1 Wash hands and put on protective gloves before testing. Wipe the fingertip to be lanced with an alcohol prep pad. Select the middle or ring finger for sampling and avoid calloused areas of the finger.

- 2 Prick the lateral side of fingertip with a lancet. Wipe off the first drop of blood with a sterile cotton swab/ball. Squeeze the fingertip again and allow a new drop of blood to form.


- 3 Dip the circular end of the inverted cup into the blood specimen (5µL).


! Fingerstick whole blood should be tested immediately.

4 Place the circular end of the inverted cup into the specimen well.


! Do not add into buffer well.

5 Break off the head of buffer along with the crease line and then twist it off.

6 Add 4 drops of buffer into the buffer well vertically and start timing.

! Do not add buffer into the specimen well.

7 Read the results after 15 minutes but no more than 30 minutes. To avoid confusion, discard the test device after reading the result.

15-30 mins

! Do not wait more than 30 minutes to read the result.

HOW TO INTERPRET THE TEST

Lines that you see	Picture	Read and Interpretation of Results
No line at 'C' (control)		<p>Result: Invalid.</p> <p>Action: Take a new cassette and repeat the test. Please contact the local distributor or manufacturer if invalid results continuously appear.</p>
Line at 'C' and No other line		<p>Result: Negative.</p> <p>Action: No</p>
Line at 'C' and at 'T1' and/or 'T2'		<p>Result: Positive.</p> <ol style="list-style-type: none"> Line at 'C', 'T1' and 'T2' indicates positive for P.f HRP2 and P.f pLDH. Line at 'C' and 'T1' indicates positive for P.f HRP2. Line at 'C' and 'T2' indicates positive for P.f pLDH. <p>Note: Line intensity may vary from faint to strong intensity. Consider a faint test line as a positive result.</p> <p>Action: No</p>

DISPOSE

! IMPORTANT

Put all used disposable components into the medical waste container. Dispose according to local infectious waste disposal laws and regulations.