WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Bioline Malaria Ag P.f/P.f/P.v
WHO reference number: PQDx 0297-012-00

Bioline Malaria Ag P.f/P.f/P.v¹ with product codes **05FK120** and **05FK123**, manufactured by **Abbott Diagnostics Korea Inc²**, **CE-marked regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 04 December 2018.

Summary of WHO Prequalification Assessment for the Bioline Malaria Ag P.f/P.f/P.v

	Date	Outcome
Prequalification listing	4 December 2018	listed
Dossier assessment	15 May 2018	MR
Site inspection(s) of the quality management system	10 May 2023	MR
Product performance evaluation	2018	MR

MR: Meets Requirements

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product, for which the WHO has been notified and has undertaken a review. The amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Updated performance characteristics in the IFU.	29 January 2019
3.0	The product name was changed from SD BIOLINE Malaria Ag P.f/P.f/P.v to Bioline Malaria Ag P.f/P.f/P.v. The manufacturer's name changed from Standard Diagnostics Inc. to Abbott Diagnostics Korea Inc. Manufacturing site inspection	26 November 2020

¹ Formerly called SD BIOLINE Malaria Ag P.f/P.f/P.v to Bioline Malaria Ag P.f/P.f/P.v.

² Formerly called Standard Diagnostics Inc.

	commitments for prequalification were fulfilled and removed from the public report.	
4.0	Updated performance evaluation results and addition of the limitations of the test on the detection of pLDH in line with the Global Malaria Programme recommendations.	13 January 2025.

Intended use:

According to the claim of intended use from Abbott Diagnostics Korea Inc, *“The Bioline Malaria Ag P.f/P.f/P.v test kit is an in vitro rapid, qualitative test for the detection of histidine rich protein II (HRP-II) antigen and plasmodium lactate dehydrogenase (pLDH) from Malaria Plasmodium falciparum and pLDH from Malaria Plasmodium vivax in human whole blood. Bioline Malaria Ag P.f/P.f/P.v test is intended for professional use, only for an initial screening test as an aid to diagnosis of clinical malaria disease and reactive specimens should be confirmed by a supplemental assay such as microscopic examination of blood smear.”*

Assay description:

According to the claim of assay description from Abbott Diagnostics Korea Inc, *“the Bioline Malaria Ag P.f/P.f/P.v test utilizes the principle of immunochromatography. The Bioline Malaria Ag P.f/P.f/P.v test kit contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to histidine-rich protein II (HRP-II) of P. falciparum in test line 1(T1), mouse monoclonal antibodies specific to plasmodium lactate dehydrogenase (pLDH) of P. falciparum in test line 2(T2) and mouse monoclonal antibodies specific to pLDH of P. vivax in test line 3(T3) region. This test device has a letter of T1, T2, T3 and C as “Test Line” and “Control Line” on the surface of the device. All the Test Lines and Control Line in result window are not visible before applying any specimens. The Control Line is used for procedural control. Control Line should always appear if the test procedure is performed properly and the test reagents of control line are working. The principle for reactive results is as follows. As the P.f (HRP-II, pLDH) and P.v reactive specimen flow through the membrane assembly of the cassette after addition of the assay diluent, they form a complex with mouse monoclonal antibodies specific to P.f HRP-II- gold conjugates and mouse monoclonal antibodies specific to pan pLDH -gold conjugates. This complex move further on the membrane to the test line region where it is immobilized by mouse monoclonal P.f HRP-II antibodies, P.f pLDH antibodies and P.v pLDH antibodies coated on the membrane leading to formation of a red-purple colored band in the respective regions which confirms a reactive test result.*

Absence of a colored band in the test region indicates a non-reactive test result for the corresponding Antigen”.

Test kit contents:

Test kit contents:			
Product code	Tests/kit	Component	Quantity
05FK120	25 T/kit	Test devise with desiccant in individual foil pouches	25
		Assay diluent vial	1 × 5 ml/vial
		Disposable inverted cup	25 (5 µL/cup)
		Sterile lancet	25
		Alcohol swabs	25
		Instruction for Use	1
05FK123	1 test/kit (bundle kit)	Each kit contains all necessary accessories (Lancet, alcohol swab, disposable applicator, and assay diluent) for a single test	25 each of 1 test/kit

Items required but not provided:

- Protective gloves
- Timer
- Biohazard container

Storage:

The test kit should be stored at 1°C to 40 °C.

Shelf-life upon manufacture:

24 months.

Warnings/limitations:

Refer to the current version of the manufacturer's instructions for use attached to this public report.

The WHO Global Malaria Programme recommends the use of tests detecting pLDH and a panel detection of at least 75% on the panel of nonHRP2 expressing specimens in areas with a high prevalence of pfhrp2/3 gene deletions. This test does not meet these criteria and is not recommended for use in such settings.³

³ <https://iris.who.int/bitstream/handle/10665/258972/WHO-HTM-GMP-2017.18-eng.pdf?sequence=1>

Prioritization for prequalification:

Based on the results of the WHO product testing of malaria RDTs for Round 8, Bioline Malaria Ag P.f/P.f/P.v was given priority for WHO prequalification assessment.

Dossier assessment

Abbott Diagnostics Korea Inc. submitted a product dossier for Bioline Malaria Ag P.f/P.f/P.v as per the *“Instructions for compilation of a product dossier”* (PQDx_018 v1). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 15 May 2018.

Based on the product dossier screening and assessment findings, the product dossier for Bioline Malaria Ag P.f/P.f/P.v meets WHO prequalification requirements.

Manufacturing site inspection

An inspection or a desk assessment of Abbott Diagnostics Korea Inc. located at 46, Hagal-ro 15 beon-gil, Giheung-gu, Yongin-si, Gyeonggi-do, 17099, Republic of Korea and 65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, 17099, Republic of Korea was conducted on 10-12 May 2023. At the time of inspection, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of a product of consistent quality. Routine inspections of the Manufacturing site will be conducted with copies of the WHO Public Inspection Report (WHOPIR) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current inspection or desk assessment performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection or desk assessment can be found at:
<https://extranet.who.int/pgweb/vitro-diagnostics/who-public-inspection-reports>

All published WHOPIRs are with the agreement of the manufacturer.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 12 October 2023.

Product performance evaluation

Bioline Malaria Ag P.f/P.f/P.v was evaluated in the eighth⁴ round of WHO product testing of RDTs for malaria antigen detection, completed in 2018.

The following results were obtained: *Plasmodium falciparum* panel detection score against HRP2-expressing parasites of 89.0% at 200 parasites/μl, *P. falciparum* detection score against HRP2-negative parasites of 42.5%⁵ at 200 parasites/μL, *P. vivax* panel detection score of 97.1% at 200 parasites/μl, false-positive rates of 0% for clean negatives, 0% for *P. falciparum* at 200 parasites/μl, 0% for *P. vivax* at 200 parasites/μl, 0% for *P. falciparum* at 2000 to 5000 parasites/μl, 0% for *P. vivax* at 2000 to 5000 parasites/μl and invalid rate of 0%.

Summary performance characteristics	Panel detection score (%)			False positive rate (%)			Invalid rate (%)
	200 parasites/μl			200 parasites/μl		Clean negatives	
	<i>Pf</i> (wild-type)	<i>Pv</i>	<i>Pf</i> (hrp2/3 deletion)	<i>Pf</i>	<i>Pv</i>		
Bioline Malaria Ag P.f/P.f/P.v	89 (89 for HRP2, 62 for Pf-pLDH)	97.1	42.5	0	0	0	0

Based on these results, the performance evaluation for Bioline Malaria Ag P.f/P.f/P.v meets the current laboratory evaluation requirements for prequalification. See Warning/limitations sections regarding use in areas with pfhrp2/3 gene deletions.

⁴ <https://www.who.int/publications/i/item/9789241514965/>

⁵ Results presented in the report “Malaria Rapid Diagnostic Test Performance – Results of WHO product testing of malaria RDTs: round 8 (2016-2018)” by WHO, FIND and CDC indicate a PDS of 20.0% on the HRP2-negative *P. falciparum* panel. This is due to an error in data analysis, which was later corrected, leading to a corrected PDS of 42.5%.


Labelling

- 1. Labels**
- 2. Instructions for use**

**Note: Labelling has been changed per site name change and product rebranding. But temporary labelling of a legacy brand (SD BIOLINE) will be used in the market according to registration status in each country.*

1. Labels

1.1 Package box for 05FK120

1.2 Package box for 05FK123

EN
Rapid, Immunochromatographic test for the qualitative detection of Histidine-rich protein 2 (HRP2) and Lactate dehydrogenase (pLDH) from *Malaria Plasmodium falciparum* and Lactate dehydrogenase (pLDH) from *Malaria Plasmodium vivax* in human whole blood

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CONTENTS :

- 25 Test devices with desiccant in individual foil pouch
- 25 Assay diluents
- 25 Disposable inverted cups (5 µl),
25 Sterile lancets, 25 Alcohol swabs
- 25 Summarized instructions for use
- 1 Instructions for use

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

PRINT AREA


MT Promed Consulting GmbH
Altenhofstrasse 80,
66386 St. Ingbert, Germany

© 2019 Abbott. All rights reserved. All trademarks
referenced are trademarks of either the Abbott group of
companies or their respective owners.


1.3 Device pouch for 05FK120, 05FK123


1.4 Outer Pouch for 05FK123

1.5 Assay diluent label

1.6 Test device image

2. Instructions for use⁶

⁶ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

2.1 IFU for 05FK120

Bioline™
Malaria Ag P.f/P.f/P.v

Malaria HRP-II (P.f), pLDH (P.f) and pLDH (P.v) Antigen Rapid Test
Test de diagnostic rapide en une etape de l'antigène la protéine HRP-II, de la pLDH de P.f et de la pLDH de P.v
Prueba en un solo paso de deteccion rápida de antígenos de HRP-II (P.f), pLDH (P.f) y pLDH (P.v) de la malaria
Teste rápido de uma etapa do antígeno Malaria HRP-II (P.f), pLDH (P.f) e pLDH (P.v)

PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

1 **EN** Open the package and look for the following:

- Test device with desiccant in individual foil pouch
- Assay diluent
- Instructions for use
- Disposable inverted cup (5 µl)
- Sterile lancet
- Alcohol swab

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1 **EN** Clean the area to be lanced with an alcohol swab.

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2 **EN** Prick the lateral side of the finger with the sterile lancet provided. Then, safely dispose of the lancet immediately after.

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3 **EN** Using a disposable inverted cup (5 µl) provided, dip the circular end of a inverted cup into the blood specimen.

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

4 **EN** Dispense 5 µl of drawn blood into round specimen well touching pad.

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EN Use Inverted cup : Let the circular end of the inverted cup touch the pad, then press down lightly.

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5 **EN** Hold assay diluent bottle vertically and dispense 4 drops of assay diluent into the assay diluent well. Exactly, 4 drops should be added. Do not let bottle tip touch device in order to avoid cross-contamination.

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
6 **EN** Interpret test results 15 minutes (up to 30 minutes) after adding assay diluent. Reading outside of this time frame (before 15 min or after 30 min) may provide false results.

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
REACTIVE / RÉACTIF / REACTIVO / REATIVO

EN ⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

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EN 1. "C", "T1" (Two lines) : HRP-II reactive

2. "C", "T2" (Two lines) : pLDH reactive

3. "C", "T1" and "T2" (Three lines) : pLDH and HRP-II reactive

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

EN "C" and "T3" (Two lines) : pLDH of P.v reactive

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
3. "C", "T1", "T2" e "T3" (quatro linhas)

INVALID / NON VALIDE / NO VALIDO / INVÁLIDO

EN If the Control Line "C" is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

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ENGLISH

Intended use

The Bioline™ Malaria Ag P.f/P.f/P.v test kit is a rapid, qualitative test for the detection of **histidine-rich protein II (HRP-II) antigen and plasmodium lactate dehydrogenase (pLDH) from Malaria Plasmodium falciparum and pLDH from Malaria Plasmodium vivax** in human whole blood.

REF 05FK123

Materials provided

- Test device
- Assay diluent
- Inverted cup (5 µl)
- Sterile lancet
- Alchol swab

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the kit or its components.
Note : When stored at refrigerator, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to the test. Do not open the pouch whilst components come to room temperature.
2. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
3. Do not use the test kit beyond its expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.
6. When transporting or storing the kit, avoid exposure to high temperature (above 45 °C) for a period longer than 1 week.

Specimen collection using a lancet

1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the end of the fingertip and pierce with a sterile lancet provided.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. With a disposable inverted cup (5 µl) provided, dip the circular end of an inverted cup into the blood specimen and carefully place the circular end of the inverted cup into the round specimen well.

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© 2019 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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
- 1

EN

Dispense 5 µl of drawn blood into round specimen well touching pad.

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
EN

Twist and pull cap to open assay diluent.

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- 3

EN

Dispense all of the assay diluent from the diluent tube into the square well of test device.

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
EN

Interpret test results 15 minutes (up to 30 minutes) after adding assay diluent. Reading outside of this time frame (before 15 min or after 30 min) may provide false results.

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⚠

Caution

: Even if the presence of faint test lines occur, the result is considered positive.

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
EN

1. “C”, “T1” (Two bands) : HRP2 Reactive
2. “C”, “T2” (Two bands) : pLDH Reactive
3. “C”, “T1” and “T2” (Three bands) : pLDH and HRP2 Reactive

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

- EN

1. “C”, “T1” and “T3” (Three bands)
2. “C”, “T2” and “T3” (Three bands)
3. “C”, “T1”, “T2” and “T3” (Four bands)
- FR

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

- Invalid / Non valide

EN

If the control line “C” is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

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- 1

EN

Open the package and look for the following:

1. Test device with desiccant in individual foil pouch

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- 2

EN

First, read carefully the instruction on how to use the Bioline™ Malaria Ag P.f/P.f/P.v test kit.
- ES

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- 3

EN

Next, look at the expiration date on the back of the foil pouch. If the expiration date has passed, use another lot. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.
- FR

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- 4

EN

Open the foil pouch and look for the following:

1. Result window

2. Specimen well

3. Assay diluent well

Then, label the device with the patient identifi er.
- FR

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
- 1

EN

Clean the area to be lanced with an alcohol swab.
- FR

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2

EN

Prick the lateral side of the finger with the sterile lancet provided. Then, safely dispose of the lancet immediately after.

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- 3


EN

Using a disposable inverted cup (5 µl) provided, dip the circular end of a inverted cup into the blood specimen.
- FR

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- EN

Use Inverted cup : Let the circular end of the inverted cup touch the pad, then press down lightly.
- FR

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- 5

EN

1. Twist and pull cap to open assay diluent.

2. Dispense all of the assay diluent from the diluent tube into the square well of test device.
- FR

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REACTIVE / RÉACTIF / REACTIVO / REATIVO

- EN

⚠ Caution: The presence of any test line, no matter how faint, the result is considered reactive.
- FR

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1. "C", "T1" (Two lines) : HRP-II reactive

2. "C", "T2" (Two lines) : pLDH reactive

3. "C", "T1" and "T2" (Three lines) : pLDH and HRP-II reactive
- FR

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- EN

"C" and "T3" (Two lines) : pLDH of P.v reactive
- FR

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INVALID / NON VALIDE / NO VALIDO / INVÁLIDO

- EN

If the Control Line "C" is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.
- FR

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