**WHO Prequalification of In Vitro Diagnostics  
PUBLIC REPORT**

**Product: STANDARD Q Malaria P.f/Pan Ag Test  
WHO reference number: PQDx 0347-117-00**

**STANDARD Q Malaria P.f/Pan Ag Test** with product code **09MAL30D**, manufactured by **SD Biosensor, Inc, CE-Mark regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 6 March 2020.

**Summary of WHO prequalification assessment for STANDARD Q Malaria P.f/Pan Ag Test**

|  | Date                              | Outcome |
|--|-----------------------------------|---------|
| <b>Prequalification listing</b>                        | 6 March 2020                      | listed  |
| <b>Dossier assessment</b>                              | 31 January 2019                   | MR      |
| <b>Site inspection(s) of quality management system</b> | 24-26 May 2023<br>14-16 June 2023 | MR      |
| <b>Product performance evaluation</b>                  | 2020                              | MR      |

MR: Meets Requirements

**Report amendments and/or product changes**

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

| Public report amendment | Summary of amendment   | Date of report amendment |
|-------------------------|--|--------------------------|
| 2.0                     | 1. Change of IFU from IFU in English to multilanguage IFU (EN, ES, FR, PT).<br>2. Change of outer package.   | 18 May 2021.             |
| 3.0                     | 1. Addition of manufacturer site for Cutting, Assembly, Buffer Preparation & dispensing, sealing, printing, packing and shipping process.<br>2. Addition of suppliers for raw materials.<br>(1) Upper device, (2) Lower device, (3) Aluminum Pouch, (4) Inverted cup, (5) Lancet, (6) Alcohol swab”. | 19 February 2024.        |

|     |  |               |
|-----|--|---------------|
| 4.0 | The addition of the manufacturing site is located in the public report of Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India. | 23 April 2024 |
|-----|--|---------------|

### Intended use

According to the claim of intended use from SD Biosensor, Inc, “*STANDARD Q Malaria P.f/Pan Ag Test is a rapid and membrane based immunochromatography for the qualitative detection of Plasmodium falciparum (P. falciparum) specific Histidine Rich Protein 2 (HRP-2) and Plasmodium species (P. falciparum, vivax, ovale and malariae) specific Plasmodium lactate dehydrogenase (pLDH) in human capillary and venous whole blood specimens of patients suspected of having malaria. STANDARD Q Malaria P.f/Pan Ag Test is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating in point-of-care settings in resource-limited lower- and middle-income countries. This product is not intended for self-testing.*”

### Assay description

According to the claim of assay description from SD Biosensor, Inc., “*STANDARD Q Malaria P.f/Pan Ag Test contains two pre-coated lines, “P.f” (P. falciparum), “Pan” (Plasmodium species: P. falciparum, vivax, ovale and malariae) as test lines and “C” as control line on the surface of the nitrocellulose membrane. The test lines and control line in the result window of the test device are not visible before applying any specimens. Monoclonal anti-P. falciparum HRP-2 is coated on the P.f test line, monoclonal anti-Malaria pLDH is coated on the Pan test line, and monoclonal anti-chicken IgY is coated on the control line region. During the test, the P. falciparum specific HRP-2 antigen and/or Plasmodium species specific pLDH in the specimen react to the gold-conjugated monoclonal anti-Malaria HRP-2 and/or gold-conjugated monoclonal anti-Malaria pLDH, and then bind to them respectively. Any P. falciparum specific HRP-2 antigen-antibody gold particle complex and/or Plasmodium species specific pLDH antigen-antibody gold particle complex also migrate with the buffer and are immobilized by monoclonal anti-P. falciparum HRP-2 and/ or monoclonal anti-Malaria pLDH at the two individual test lines to formation of violet test colored band(s) which confirms a positive result. Absence of this violet colored band indicates a negative result. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working. ”*

### Test kit contents

| Component   | 25 tests<br>(product code 09MAL30D) |
|---|-------------------------------------|
| Test device (individually in a foil pouch with desiccant) | 25                                  |
| Buffer Bottle   | 1 x 4 mL                            |
| Inverted cup (5µL)  | Pack of 25                          |

|                      |            |
|----------------------|------------|
| Sterile lancet       | Pack of 25 |
| Alcohol swab         | Pack of 25 |
| Instructions for use | 1          |

**Items required but not provided**

- Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood
- Micropipette and tip
- Timer
- PPE (Personal Protective Equipment)
- Pen/pencil
- Extra lancets and alcohol swabs
- Biosafety sharps container
- Biohazard container
- Sterile gauze

**Storage**

The test kit should be stored at 2-40°C.

**Shelf-life upon manufacture**

24 months.

**Warnings/limitations**

Refer to current version of the manufacturer's instructions for use.

**Prioritization for Prequalification**

Based on the established eligibility criteria, the STANDARD Q Malaria P.f/Pan Ag Test was given priority for the WHO prequalification assessment.

**Dossier assessment**

SD Biosensor, Inc. submitted a product dossier for STANDARD Q Malaria P.f/Pan Ag Test as per the "*Instructions for compilation of a product dossier*" (PQDx\_018 version 3). The information (data and documentation) submitted in the product dossier was reviewed by WHO staff and external technical experts (assessors) appointed by WHO.

The manufacturer's responses to the nonconformities found during dossier screening and assessment findings were accepted on 31 January 2019.

## Commitment for Prequalification

The manufacturer committed to submit the in-use stability study by the end of April 2021. WHO will follow up on implementing this commitment at the next re-inspection or date assigned in the commitment review letter.

Based on the product dossier screening and assessment findings, the product dossier for the STANDARD Q Malaria P.f/Pan Ag Test meets WHO prequalification requirements.

## Manufacturing site inspection

An inspection of **SD Biosensor Inc.**, located at *74, Osongsaeengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea* was conducted between 24-26 May 2023 and Plot no. 38, Sector 4, IMT Manesar, Gurugram, Haryana 122052, India between 14-16 June 2023. At the time of considering the product application for Prequalification, the Manufacturer of the product had a well-established quality management system and manufacturing practices in place that would support the manufacture of the product of consistent quality. Routine inspections of the Manufacturer will be conducted with copies of these WHO Public Inspection Reports (WHOPIRs) published on the WHO Prequalification web page as per Resolution WHA57.14 of the World Health Assembly. Note that a WHOPIR reflects the information on the most current inspection performed at a manufacturing site for in vitro diagnostic products and gives a summary of the inspection findings.

Information on the most current inspection can be found at:

<https://extranet.who.int/prequal/inspection-services/prequalification-reports/whopirs-vitro-diagnostics>

All published WHOPIRs are with the manufacturer's agreement.

## Product performance evaluation

**STANDARD Q Malaria P.f/Pan Ag Test<sup>1</sup>** was included in the eighth round of WHO product testing of RDTs for malaria antigen detection, which was completed in 2018. The product was evaluated against a *Plasmodium falciparum* cultured line panel, *P. falciparum* wild-type parasite panel, *P. vivax* wild-type parasite panel and a negative panel. Thermal stability was assessed after 2 months of storage at elevated temperature and humidity, and a descriptive ease of use assessment was recorded.

<sup>1</sup> The product code of the product included in the WHO product testing of malaria RDTs was 09MAL30B, which corresponded to the same product, with a different kit configuration (i.e. number of tests) <https://www.who.int/malaria/publications/atoz/9789241514965/en/>

Based on the demonstrated *P. falciparum* panel detection score (88.0% at 200 parasites/ $\mu$ l), *P. vivax* panel detection score (100% at 200 parasites/ $\mu$ l), false-positive rates (0% for clean negatives, 0% for *P. falciparum* at 200 parasites/ $\mu$ l, 0% for *P. vivax* at 200 parasites/ $\mu$ l, 0.5% for *P. falciparum* at 2000 parasites/ $\mu$ l, 0% for *P. vivax* at 2000 / $\mu$ l) and invalid rate (0%), STANDARD Q Malaria P.f/Pan Ag Test meets the current laboratory evaluation requirements for prequalification.

| Summary performance characteristics | Panel detection score  |     | False positive rate (%) |    |                 | Invalid rate (%) |
|-------------------------------------|------------------------|-----|-------------------------|----|-----------------|------------------|
|                                     | 200 parasites/ $\mu$ l |     | 200 parasites/ $\mu$ l  |    | Clean negatives |                  |
|                                     | Pf                     | Pv  | Pf                      | Pv |                 |                  |
| STANDARD Q Malaria P.f/Pan Ag Test  | 88.0                   | 100 | 0                       | 0  | 0               | 0                |

Based on these results, the STANDARD Q Malaria P.f/Pan Ag Test meets the current performance evaluation requirements for prequalification.

## **Labelling**



- 1. Labels**
- 2. Instructions for use**

## **1. Labels**

### **1.1 Kit packaging artwork**

# STANDARD Q Malaria P.f/Pan Ag Test 25T

# STANDARD Q Malaria P.f/Pan Ag

STANDARD™ Q Malaria P.f/Pan Ag Test

[www.sdbiosensor.com](http://www.sdbiosensor.com)

25 TESTS  
KIT



SD BIOSENSOR


# STANDARD Q Malaria P.f/Pan Ag Test 25T


## **1.2 Foil pouch labels**



15mm Notch

15mm Notch

## 1.3 Buffer labels

**1.4 Inverted cup(5  $\mu\text{l}$ ) label**

# Inverted cup (5 $\mu\text{l}$ )

LOT No. : XXXXXX

EXP : YYYY.MM.DD.

Quantity : 25PCS

L15IC1ENR2  
Issue date : 2020.02

**1.5 Sterile Lancet label**

# Disposable Sterile Lancets

LOT No. : XXXXXX

Product code : 01GL25

MFG Date : YYYY.MM.DD.

25PCS

EXP Date : YYYY.MM.DD.

28G

**INTENDED USE**

To obtain a capillary blood specimen from the fingertip.

B06LCT2ENR2

**INSTRUCTIONS FOR USE**

Issue Date : 2020.02

To use, twist-off the protective cap.

**CAUTION**

The lancet is guaranteed sterile while protective cap is sealed to the base.

Do not use if the seal has been damaged or broken.

Manufactured by

Beijing Ruicheng Medical Supplies Co., Ltd.  
No. 558 Zhangzikou, Yangsong Town,  
Huairou District, 101400 Beijing, China

EC   REP

Authorized Representative

Lotus NL B.V.  
Koningin Julianaplein 10, 1e Verd,  
2595AA, The Hague, Netherlands.  
Tel: +31645171879, +31626669008

STERILE R

# Disposable Sterile Lancets

LOT No. : X00000X

Product code : 01GL27

MFG Date : YYYY.MM.DD.

25PCS

EXP Date : YYYY.MM.DD.

28G

## INTENDED USE

To obtain a capillary blood specimen from the fingertip.

B05SL2ENRO

## INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

Issue Date : 2020.07

## CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base.

Do not use if the seal has been damaged or broken.

### Manufactured by

Tianjin Huahong Technology Co., Ltd.  
A01, Plant B No. 278, Hangkong Road,  
Tianjin Pilot Free Trade Zone  
(Air Port Industrial Park),  
300308 Tianjin, China

### Authorized Representative

Shanghai International Holding Corp.  
GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany  
Tel:+49-40-2513175

STERILE R

**SURFACE JOB**

220MM

110MM

110MM

100MM

## **Disposable Sterile Lancets**

 LOT No. : XXXXXXXX

Product code : SM.01

 MFG Date : YYYY-MM

**25PCS**

 EXP Date : YYYY-MM

**28G**

### INTENDED USE

To obtain a capillary blood specimen from the fingertip.

### INSTRUCTIONS FOR USE

To use, twist-off the protective cap.

Issue Date : 2023.09

### CAUTION

The lancet is guaranteed sterile while protective cap is sealed to the base.  
Do not use if the seal has been damaged or broken.

 Manufactured by

 EC REP Authorized Representative

Phoenix Innovative Healthcare Manufacturers Pvt. Ltd.  
EL-209, Shil Mahape Road, Electronic Zone, MIDC,  
TTC Industrial Area, Mahape, Navi Mumbai - 400 710,  
Maharashtra, India

Med Path Gm6H  
Miles-van-der-Rohe-Strasse 8,80807  
Munich, Germany

 CE 0197

 STERILE

25 Pcs.


### **1.6 Alcohol swab labels**

## Item Code

IM2006

## Description

FA Alcohol Swab-S (8Line)

Front

PANTONE Reflex Blue U


PANTONE Orange 021 U

2@ 3@ 4@ 5@ 6@ 7@ 8@

Back

2@ 3@ 4@ 5@ 6@ 7@ 8@

**1.7 Cassette image**

## **2. Instructions for use<sup>2</sup>**

---

<sup>2</sup> English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

EN

Cat. No.: 09MAL300 REF: QMAL03G

**STANDARD Q Malaria P.f/Pan Ag**  
 STANDARD™ Q Malaria P.f/Pan Ag Test

DE, IT

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

SD BIOSENSOR

**KIT CONTENTS****MATERIALS REQUIRED BUT NOT PROVIDED**

- ① Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood.
- ② Micropipette and tip ③ Timer
- ④ PPE (Personal Protective Equipment) ⑤ Pen/pencil
- ⑥ Extra lancets and alcohol swabs ⑦ BioSafety sharps container
- ⑧ Biohazard container ⑨ Sterile gauze

**SPECIMEN COLLECTION AND PREPARATION****■ Capillary whole blood**

1. Capillary whole blood should be collected aseptically by fingerstick.
2. Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
3. Squeeze the end of the fingertip and pierce with a sterile lancet.
4. Wipe the first drop with a sterile gauze.
5. Collect capillary whole blood and it must be tested immediately after collection.


**■ Venous whole blood**

1. Collect venous whole blood into a commercially available anti-coagulant tube such as heparin, EDTA, or sodium citrate by venipuncture.
2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 72 hours after collection.
3. Do not use hemolyzed blood specimen.

**CAUTION**  
Anticoagulants such as heparin, EDTA, or sodium citrate do not affect the test result.  
As known relevant interference, hemolytic specimen, rheumatoid factors-contained specimen and lipemic, icteric specimen may lead to inaccurate results.  
Use separate inverted cup for each specimen in order to avoid cross-contamination of either specimens which can cause erroneous results.

**PREPARATION AND TEST PROCEDURE****■ Preparation**

1. Carefully read the instructions for using the STANDARD Q Malaria P.f/Pan Ag Test.
2. Look at the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.
3. Open the foil pouch, and check the test device and the desiccant inside the foil pouch.


**CAUTION**  
Bring the kit contents and the specimens to room temperature before testing.

**■ Test Procedure**

Capillary and venous whole blood can be used for the STANDARD Q Malaria P.f/Pan Ag Test.

**1. Test with a sterile lancet**

1. Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.
2. Dry and pierce the wiped fingertip with a sterile lancet to bleed.

If the control band ("C" Control line) is not visible within the window, the result is considered invalid. The directions may not have been followed correctly or the test device may have deteriorated. Re-test with a new specimen and a new test device.

\*NOTE • Pf.: Plasmodium falciparum, • Pv.: Plasmodium vivax, • Po.: Plasmodium ovale, • Pm.: Plasmodium malariae


**PERFORMANCE CHARACTERISTICS****■ Limit of Detection (LOD)**

|                                   | P.f Line              |                       | Pan Line         |                  |                     |
|-----------------------------------|-----------------------|-----------------------|------------------|------------------|---------------------|
|                                   | Plasmodium falciparum | Plasmodium falciparum | Plasmodium vivax | Plasmodium ovale | Plasmodium malariae |
| Limit of detection (parasites/µL) | 163                   | 327                   | 152              | 35               | 224                 |


**■ Potential Interfering Substances**

STANDARD Q Malaria P.f/Pan Ag Test was evaluated with the following interfering substances present in specimen in order to assess their potential effect on the assay performance as per CLSI guideline EP7-A2. It was not affected by interfering materials such as artemether-lumefantrine (malaria drug), doxycycline hydate (malaria drug), lamivudine (retroviral medication), acetaminophen, erythromycin, ibuprofen, acetylsalicylic acid, bilirubin (unconjugate), caffeine, ethanol. There was no interference with Human anti-mouse antibody, whole blood of pregnant women, infant and neonate specimen and whole blood having elevated levels of hemoglobin, elevated levels of C-reactive protein and lipidemic specimen.


3. Wipe away the first drop of blood with sterile gauze.


4. Take the inverted cup (5µl) provided to collect 5µl of the capillary whole blood specimen by dipping the circular end of the inverted cup into the specimen.


5. To apply the collected capillary whole blood specimen to the specimen well(circle) of the test device, place the circular end of the inverted cup on the specimen well pad and then press down lightly.

**1.1 Test by venipuncture**


1. Gently invert 4-5 times for proper mixing.

2. Using a micropipette, draw 5µl of venous whole blood specimen.

3. Apply the collected venous whole blood specimen to the specimen well(circle) of the test device.

**2. Dropping of Buffer**


1. Hold the buffer bottle at 90° (perpendicular) to the test device and not touching the specimen well to avoid contamination. Add 3 drops (90µl) of the buffer into the buffer well(square) of the test device.

**3. Reading Time**

1. Read the test result in 15 minutes. The test can be read up to 30 minutes.

2. Do not read test results after 30 minutes. It may give false results.

3. If the T line does not appear until 15 minutes, make sure to observe the appearance of the T line until 30minutes.

**4. KIT STORAGE AND STABILITY**

1. Store the sealed pouch and the buffer provided in the kit at 2-40°C/36-104°F out of the direct sunlight for the duration of its shelf life.

2. Do not open the foil pouch until you are ready to perform a test.

3. Close the buffer cap tightly after using, and then store it at 2-40°C/36-104°F out of the direct sunlight. It is stable until the expiry date of the kit and the buffer label after opening its cap, if it is tightly closed.

**5. WARNINGS AND PRECAUTIONS**

1. Do not freeze.

2. Do not use beyond the expiration date.

3. Do not use the test kit.

4. Do not use the test kit if the pouch is damaged or the seal is broken.

5. Do not use the test kit if the buffer is damaged or the seal is broken.

6. Do not smoke, drink or eat while handling specimen.

7. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly afterwards.

8. Clean up spills thoroughly using an appropriate disinfectant.

9. Handle all specimens as if they contain infectious agents.

10. Observe established precautions against microbiological hazards throughout testing procedures.

11. Dispose of all specimens and materials used to perform the test as biohazard waste. Laboratory chemical and bio-hazard wastes must be handled and discarded in accordance with all local, state, and national regulations.

12. In the pouch, there is a desiccant containing a humidity indicator. The humidity indicator should be yellow. If the humidity indicator is green, throw away the test device and take another test device. Throw away the desiccant in the nonsharps (non-infectious) disposal.

13. Before testing, check that the lot number of buffer on the bottle matches the lot number of the package and the pouch you are using.

14. The buffer contains 0.09% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with large quantities of water.

15. For in vitro diagnostic use only.

**5. INTERPRETATION OF TEST RESULT****\* "C" Control Line | "Pan" (Pf., Pv., Po. and Pm.) Line | "P.f" P.f Line****Negative****Positive****Pf. Positive or mixed positive of Pv., Po. and/or Pm. with P.f.****Positive Pf., Pv., Po., or Pm. positive or mixed positive of these****Invalid****Caution****The presence of only one colored band ("C" Control line) within the result window indicates a negative result.****Positive****Pf. Positive or mixed positive of Pv., Po. and/or Pm. with P.f.****Pf., Pv., Po., or Pm. positive or mixed positive of these****Invalid****Caution****The presence of any line, no matter how faint the result is considered positive.****Invalid****Caution****If the control band ("C" Control line) is not visible within the window, the result is considered invalid. The directions may not have been followed correctly or the test device may have deteriorated. Re-test with a new specimen and a new test device.****Note****Pf.: Plasmodium falciparum, Pv.: Plasmodium vivax, Po.: Plasmodium ovale, Pm.: Plasmodium malariae**

1. A colored band will appear in the top section of the result window to show that the test is working properly. This band is the control line (C).
2. Colored bands will appear in the middle and lower section of the result window. These bands are the *Plasmodium falciparum*, *vivax*, *ovale* and *malariae* line (P) and the *Plasmodium falciparum* line (Pf).

**\*This test only indicates the presence of Histidine Rich Protein-2 (HRP2) of *Plasmodium falciparum* and *Plasmodium lactate dehydrogenase* (PLDH) of *Plasmodium species* (*Plasmodium falciparum*, *vivax*, *ovale*, and *malariae*) in human capillary and venous whole blood specimens. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.**

**QUALITY CONTROL**

1. The test procedure, precautions, and interpretation of results for this test must be followed strictly when testing.

2. STANDARD Q Malaria P.f/Pan Ag Test is designed for use only on human capillary and venous whole blood specimens.
3. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- LIMITATION OF TEST**
1. The test procedure, precautions, and interpretation of results for this test must be followed strictly when testing.
2. STANDARD Q Malaria P.f/Pan Ag Test is designed for use only on human capillary and venous whole blood specimens.

3. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

**CLINICAL SENSITIVITY**

1. Pf. confirmed specimen on Pf line: 99.58% (476/478, 95% CI: 98.50-99.95%)

2. Pf. confirmed specimen on Pan line: 56.69% (271/478, 95% CI: 52.12-61.19%)

3. Pv., Po. and Pm. confirmed specimen on Pan line: 100% (129/129, 95% CI: 97.18-100%)

4. Pv., Po. and Pm. co-infection confirmed specimen on Pf line: 100% (11/11, 95% CI: 71.51-100%)
<div data-bbox="40

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


**BIBLIOGRAPHY**

- World Malaria Report 2015, World Health Organization, p.8-12, 2015.
- Confronting Plasmodium vivax malaria, World Health Organization, p.1-3, 2015.
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition, U.S. Department of Health and Human Services, CDC, NIH, Washington DC, p.182-184, 2007.
- Guidelines for the treatment of malaria, World Health Organization, p.5-6, 2006.
- Recommended selection criteria for procurement of malaria rapid diagnostic tests, World Health Organization, p.1-4, 2015.
- C. Beadle, G.W. Long, Mc Elroy, P. Olliaro, C.M. Naing, S.M. Maret, A.O. Oloso, Diagnosis of malaria by detection of Plasmodium falciparum HRP-2 antigen with a rapid dipstick antigen-capture assay, *The Lancet*, Vol. 343, Issue 8897, 1994, p.564-568.
- A. Moody, Rapid diagnostic tests for malaria parasites, *Clin Microbiol Rev*, Vol. 15, Issue 1, 2002, p.66-78.
- B. LB Ochola, P. Vounatsou, T. Smith, M. Mabaso, C. Newton, The reliability of diagnostic techniques in the diagnosis and management of malaria in the absence of a gold standard, *Lancet Infect Dis*, Vol. 6, Issue 9, 2006, p.582-588.
- K. Abba, J.J. Deeks, P. Olliaro, C.M. Naing, S.M. Jackson et al, Rapid diagnostic tests for diagnosing uncomplicated Plasmodium malaria in endemic areas, *Cochrane Database Syst Rev*, Vol. 18, Issue 12, 2011.
- M. Ahmed, M. Hossain, A.S. Suzuman, M. Alam, A. Khan et al, Rapid Diagnosis of malaria by Antigen Detection, *Bangladesh Journal of Medical Microbiology* Vol. 15, Issue 1, 2009, p. 14-16.

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diagnostic techniques in the diagnosis and management of malaria in the absence of a gold standard, *Lancet Infect Dis*, Vol. 6, Issue 9, 2006, p.582-588.

6. C. Beadle, G.W. Long, Mc Elroy, P. Olliaro, C.M. Naing, S.M. Maret, A.O. Oloso, Diagnosis of malaria by detection of Plasmodium falciparum HRP-2 antigen with a rapid dipstick antigen-capture assay, *The Lancet*, Vol. 343, Issue 8897, 1994, p.564-568.

7. K. Abba, J.J. Deeks, P. Olliaro, C.M. Naing, S.M. Jackson et al, Rapid diagnostic tests for diagnosing uncomplicated Plasmodium malaria in endemic areas, *Cochrane Database Syst Rev*, Vol. 18, Issue 12, 2011.

8. A. Moody, Rapid diagnostic tests for malaria parasites, *Clin Microbiol Rev*, Vol. 15, Issue 1, 2002, p.66-78.

9. K. Abba, J.J. Deeks, P. Olliaro, C.M. Naing, S.M. Jackson et al, Rapid diagnostic tests for diagnosing uncomplicated Plasmodium malaria in endemic areas, *Cochrane Database Syst Rev*, Vol. 18, Issue 12, 2011.

10. M. Ahmed, M. Hossain, A.S. Suzuman, M. Alam, A. Khan et al, Rapid Diagnosis of malaria by Antigen Detection, *Bangladesh Journal of Medical Microbiology* Vol. 15, Issue 1, 2009, p. 14-16.

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
STANDARD™ Q Malaria P.f/Pan Ag Test

DE, IT

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST © SD BIOSENSOR

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
- ① Anti-coagulant tube containing heparin, EDTA or sodium citrate for collection of venous whole blood
- ② Micropipette and tip ③ Timer
- ④ PPE (Personal Protective Equipment) ⑤ Pen/pencil
- ⑥ Extra lancets and alcohol swabs ⑦ Biosafety sharps container
- ⑧ Biobranchard container ⑨ Sterile gauze

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
**■ Venous whole blood**

1. Collect venous whole blood into a commercially available anti-coagulant tube such as heparin, EDTA, or sodium citrate.
2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 72 hours after collection.
3. Do not use hemolyzed blood specimen.

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

**1. Test with a sterile lancet**

1. Select the finger that is not calloused and clean the fingertip by wiping with an alcohol swab.


2. Dry and pierce the wiped fingertip with a sterile lancet to bleed.

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

4. Take the inverted cup (5μl) provided to collect 5μl of the capillary whole blood specimen by dipping the circular end of the inverted cup into the specimen.


5. To apply the collected capillary whole blood specimen to the specimen well/circle of the test device, place the circular end of the inverted cup on the specimen well pad and then press down lightly.

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
2. Using a micropipette, draw 5μl of venous whole blood specimen.

3. Apply the collected venous whole blood specimen to the specimen well/circle of the test device.

**2. Dropping of Buffer**


1. Hold the buffer bottle at 90° (perpendicular) to the test device and not touching the specimen well to avoid contamination. Add 3 drops (90μl) of the buffer into the buffer well (square) of the test device.

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2. Do not read test results after 30 minutes. If the T line does not appear until 15 minutes, make sure to observe the appearance of the T line until 30 minutes.

**INTERPRETATION OF TEST RESULT**

\* "C" Control Line | "Pan" | "P.f., P.v., P.o. and P.m." Line | "P.f." P.f Line

Negative

Yellow : Valid  
Green : Invalid

<Desiccant>

The presence of only one colored band ("C" Control line) within the result window indicates a negative result.

**Positive**

P.f. Positive or  
mixed positive of  
P.v., P.o. and/or P.m.  
with P.f.

P.f., P.v., P.o., or P.m.  
positive or  
mixed positive of these

**Invalid**

C Pan P.f.

The presence of any line, no matter how faint the result is considered positive.

**Caution**

If the control band ("C" Control line) is not visible within the result window, the result is considered invalid. The directions may not have been followed correctly or the test device may have deteriorated. Re-test with a new specimen and a new test device.

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