WHO Prequalification of In Vitro Diagnostics PUBLIC REPORT

Product: Wondfo HIV Self-Test WHO reference number: PQDx 0357-004-01

Wondfo HIV Self-Test with product codes W006P0058, W006P0059, and W006P0060, manufactured by Guangzhou Wondfo Biotech Co., Ltd, Rest-of-World regulatory version, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 13 July 2022.

Summary of WHO prequalification assessment for Wondfo HIV Self-Test¹

	Date	Outcome
Prequalification listing for the Wondfo HIV Self-Test	13 July 2022	Listed
Prequalification listing for the One Step HIV1/2 Whole Blood/Serum/Plasma Test	29 November 2018	Listed
Dossier assessment	26 April 2018	MR
Site inspection(s) of the quality management system	12 October 2018	MR
Product performance evaluation	First-quarter of 2018	MR

MR: Meets Requirements

*Change notification

In 2021, Guangzhou Wondfo Biotech Co., Ltd submitted a change notification for their prequalified product One Step HIV1/2 Whole Blood/Serum/Plasma Test to introduce a new configuration with an intended use specific for HIV self-testing (Wondfo HIV Self-Test). The new format was adapted from the corresponding professional use product (One Step HIV1/2 Whole Blood/Serum/Plasma Test), for which a WHO prequalification assessment has already occurred. Guangzhou Wondfo Biotech Co., Ltd generated additional data to meet requirements for self-testing as set out in the WHO Technical Specifications Series document TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and self-testing².

¹ Dossier assessment and product performance evaluation for the Wondfo HIV Self-Test were considered from the previous assessment of professional use product, One Step HIV1/2 Whole Blood/Serum/Plasma Test, which was prequalified in 2018. Based on the product dossier assessment and product performance evaluation, Wondfo HIV Self-Test meets WHO prequalification requirements. Please refer to <https://extranet.who.int/prequal/WHOPR/public-report-one-step-hiv12-whole-bloodserumplasma-test-pqdx-0357-004-00>

² <https://apps.who.int/iris/bitstream/handle/10665/251857/9789241511742-eng.pdf;jsessionid=153ABC9D88E7623A1AD1DF946A22B4C8?sequence=1>

Report amendments and product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which the WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Public report amendment	Summary of amendment	Date of report amendment
2.0	<p>Added local languages to IFU. The changes included:</p> <ol style="list-style-type: none"> 1. The length of the box for W006P0059(20T/kit) increased by 2 cm, from 230X185X105mm to 250X185X105mm; another artwork of the box has not changed. 2. Two optional sizes of IFU, the prequalified 420*285 mm and the additional 570*420 mm. 	19 July 2024

Intended use:

According to the intended use claim from Guangzhou Wondfo Biotech Co., Ltd, *“the Wondfo HIV Self-Test is a single-use in vitro diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as a self-test and/or by medical professionals. For in vitro diagnostic use only.”*

Assay Description:

According to the claim of assay description from Guangzhou Wondfo Biotech Co., Ltd, *“Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into their respective wells, they migrate along the device by capillary action. When the levels of HIV-1 or HIV-2 antibodies are at or above the target cut-off level, HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen(gp36/41). The complex is captured by HIV recombinant antigen (gp41and gp36) immobilized in the test region (T). This produces a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV 1/2 antibody are zero or below the target cut-off level, there is no visible colored band in the test region (T). This indicates a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C), if the test has been performed properly.”*

Test kit contents

Catalog No.		W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)	1	1x20	1x100
	Desiccant (pcs)	1	1x20	1x100
Accessories	Dropper (pcs)	1	1x20	1x100
	Buffer (vial)	1	1x20	1x100
	IFU (pcs)	1	1x20	1x100
	Blood Lancet for Single Use (pcs)	1	1x20	1x100
	Alcohol prep pad (pcs)	1	1x20	1x100
	Cotton swab (pcs)	1	1x20	1x100
	Disposal bag (pcs)	1	/	/

Materials required but not provided:

- Timer
- Tissue
- Bandage.

Storage

The test kit must be stored at 2 -30 °C.

Shelf-life upon manufacture

24 months (test kit and buffer).

Warnings/limitations

Refer to manufacturer's instructions for use.

Labelling

- 1. Labels**
- 2. Instructions for use**

1. Labels

1.1 Labels for product code W006P0058

Rev. A4 Rel.: 2023/03/31

Wondfo

Leading POCT Manufacturer

Contents:

- | | |
|--------------------|-------------------------------|
| 1 Test cassette | 1 Desiccant pouch |
| 1 Dropper | 1 Instruction for use (IFU) |
| 1 Buffer | 1 Cotton swab |
| 1 Alcohol Prep Pad | 1 Blood Lancet for Single Use |
| 1 Disposal bag | |


REF W006P0058
LOT WXXXXXXXXXX
 YYYY-MM-DD
 YYYY-MM-DD

Wondfo

Leading POCT Manufacturer

Wondfo HIV Self-Test

**FOR IN VITRO DIAGNOSTIC USE ONLY
FOR SELF-TESTING USE**


Wondfo
Leading POCT Manufacturer


Guangzhou Wondfo Biotech Co., Ltd.
Add: No. 8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063 Website: en.wondfo.com

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXXXX which is nine characters.

Front

Back

Foil Pouch for W006P0058

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Test Cassette Foil Pouch for W006P0058


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front

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
Back


110mm

70mm

Disposal Bag for W006P0058


Round Sticker

1.2 Labels for product code W006P0059

230 x185x 105 mm

Note: The lot number have three formats. Format 1 is WXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Wondfo
S'aaäq * ÄUUÖVÄ' äj' -e&C' !Ä!

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishihe Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: en.wondfo.com

Wondfo
S'aaäq * ÄUUÖVÄ' äj' -e&C' !Ä!

IVD

Wondfo
S'aaäq * ÄUUÖVÄ' äj' -e&C' !Ä!

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishihe Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: en.wondfo.com

IVD

Wondfo
S'aaäq * ÄUUÖVÄ' äj' -e&C' !Ä!


REF W006P0059
LOT WXXXXXXXX
EXP YYYY-MM-DD
EXP YYYY-MM-DD

IVD

Health
We're Working For You

250 x 185 x 105mm

Note: The lot number have three formats. Format 1 is WXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishiwan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: en.wondfo.com

IVD

Wondfo
Leading POCT Manufacturer

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lishiwan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: en.wondfo.com

Version: A1
Rel.: 2023/03/31

Wondfo
Leading POCT Manufacturer

REF	W006P0059
LOT	WXXXXXXXX
	YYYY-MM-DD
	YYYY-MM-DD

IVD

Wondfo
Leading POCT Manufacturer

We Are Working For Your Health

IVD

Wondfo HIV Self-Test

CONTENTS:


- 20 Individual pouches, each containing:
 - 1 Test cassette
 - 1 Desiccant pouch
 - 1 Dropper
 - 1 Instruction for use (IFU)
 - 1 Blood Lancet for Single Use
 - 1 Alcohol Prep Pad
 - 1 Buffer
 - 1 Cotton Swab


Version: A4
Rel.: 2023/03/31

Secondary Packaging Label of W006P0059

Front

Back

Foil Pouch for W006P0059

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Test Cassette Foil Pouch for W006P0059


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front

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
Back

110mm

70mm

Round Sticker

1.3 Labels for product code W006P0060

Wondfo HIV Self-Test

Wondfo
Leading POCT Manufacturer

IVD

Version: A3
Rel.: 2023/03/31

CONTENTS:

100 Individual pouches, each containing:

1 Test cassette	1 Desiccant pouch
1 Cotton swab	1 Alcohol prep pad
1 Blood Lancet for Single Use	1 Buffer
1 Instruction for use (IFU)	1 Dropper

REF

W006P0060

LOT

WXXXXXXXXX

YYYY-MM-DD

YYYY-MM-DD

Guangzhou Wondfo Biotech Co., Ltd.


No. 8 Lizhishan Road, Science City, Huangpu District,
510663, Guangzhou, P.R. China
Tel: 0086-20-32296083 Fax: 0086-20-32296063
Website: en.wondfo.com


Secondary Packaging Label for W006P0060

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front


Back

Foil Pouch for W006P0060

Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Test Cassette Foil Pouch for W006P0060


Note: The lot number have three formats. Format 1 is WXXXXXXXXXX which is eleven characters, format 2 is WXXXXXXXXX which is ten characters, and format 3 is WXXXXXXXX which is nine characters.

Front

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Back

110mm

70mm

2.0 Instructions for use³


³ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

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Contents

Watch the operation video

Materials may be required but not provided.

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Components		Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)		1	1x20	1x100
	Desiccant (pcs)		1	1x20	1x100
	Dropper (pcs)		1	1x20	1x100
	Buffer (vial)		1	1x20	1x100
Accessories	IFU (pcs)		1	1x20	1x100
	Blood Lancet for Single Use (pcs)		1	1x20	1x100
	Alcohol prep pad (pcs)		1	1x20	1x100
	Cotton swab (pcs)		1	1x20	1x100
	Disposal bag (pcs)		1	/	/

One test strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid) and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

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SYMBOLS KEY

IVD	In vitro diagnostic medical device	Consult instructions for use	Use-by date
Contains sufficient for <n> tests	Date of manufacture	Keep dry	Keep away from sunlight
LOT	Batch code	Temperature limit	Do not re-use
Manufacturer	Caution	Do not use if package is damaged and consult instructions for use	REF

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RÉFÉRENCES

1. WHO. TGS-5 Designing instruction for use for in vitro diagnostic medical devices, Geneva: World Health Organization ; 2019.
2. Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
3. WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or selftesting, Geneva: World Health Organization; 2016.

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

Comment lire le test

One strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid) and rabbit IgG polyclonal antibody-gold colloid. Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

STORAGE AND STABILITY

- The test kit can be stored at 2-30 °C for 24 months.
- Use the test cassette within 1 hour after opening the pouch.
- Keep away from sunlight, moisture and heat.
- Use the kit at 10-30 °C.

LIMITATIONS OF THE PROCEDURE

- The test is designed for detecting human fingerstick whole blood.
- The test is limited to the qualitative detection of HIV-1 and HIV-2 antibodies.
- The assay procedure and result interpretation must be followed closely when testing. For optimal test performance, proper specimen collection is critical. Failure to follow the procedure may lead to inaccurate test results.
- False-negative results can occur in the following conditions:
 - Patients exposed to HIV less than 3 months.
 - Patients under HIV treatment (Antiretroviral therapy).
 - If the quantity of antibodies for HIV present in the specimen is below the detection limit of the assay.
- False-positive results can occur in the following conditions:
 - Patients have participated in a HIV vaccine clinical trial.
 - The presence of the control line only means that migration of added liquid occurred. It does not guarantee that:
 - The correct specimen has been used.
 - The specimen has been applied correctly.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but instead should be determined by a healthcare provider in conjunction with clinical findings and the results from other laboratory tests and evaluations. Results from the Wondfo HIV Self-Test should not be used as the sole basis for diagnosis.

PERFORMANCE CHARACTERISTICS

In the clinical study, 900 participants whose HIV status were unknown were given the Wondfo HIV Self-Test to test. The results were compared to the 4th generation laboratory test. The laboratory results showed that a total of 77 participants were HIV positive, 822 participants were HIV negative and 1 undetermined. A total of 43 participants (5 HIV positive, 37 HIV negative and 1 undetermined) were excluded from the performance analysis. The comparison of results was as follows:

- 95.8% of participants (69/72) correctly reported the result as positive. This means that 3 participants infected with HIV reported negative result. This is called a false negative.
- 99.6% of participants (782/785) correctly reported the result as negative. This means that 3 participants not infected with HIV reported positive result. This is called a false positive.
- 4.7% of participants (42/900) failed to obtain a result. 1 participant's HIV infection status was not confirmed during the clinical study so it was excluded from the analysis.

REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices. Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or selftesting. Geneva: World Health Organization; 2016.

SYMBOLS KEY

	In vitro diagnostic medical device		Consult instructions for use		Use-by date
	Contains sufficient for <n> tests		Date of manufacture		Keep dry
	Batch code		Temperature limit		Keep away from sunlight
	Manufacturer		Do not re-use		Catalogue number
	Caution		Do not use if package is damaged and consult instructions for use		

Manufacturer information

Guangzhou Wondfo Biotech Co., Ltd.
 Add: No.8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China
 Tel: +86-20-32296083 400-888-5268 (Toll Free)
 Fax: +86-20-32296063
 E-mail: global@wondfo.com.cn
 Website: en.wondfo.com

Please contact the manufacturer or your local distributor if you have any questions related to the product.

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RÉFÉRENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices. Geneva: World Health Organization ; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or selftesting. Geneva: World Health Organization; 2016.

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
HIV Self-Test

How to use am

- You must to follow how dem dey do de test je-je to get beta result.
- Siddon for clean place wey get light and make sure say yu get tins wey u need first for de test.
- Na only one-time dem dey use de test. No open am until yu ready do de test.


Things way dey inside

Some things wey yu go need but wey no dey inside.

Preparation

Questions and answers

1. Way thing be HIV?

The human immunodeficiency virus (HIV) dey target cells of the immune system, wey dem dey call CD4 cells, wey dey help body fight infection. Inside the CD4 cell, HIV dey double come dey damage and destroys the cell. If better treatment of a combination of antiretroviral (ARV) drugs no dey, the immune system no go strong again, and e no go fit fight infection and disease.

2. How dem they transmit HIV?

HIV dey stay for body of people wey get HIV, for blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV dey spread tru:

- When yu collect blood from person wey get HIV.
- When yu use needles, syringes, other injecting equipment, surgical equipment, or things wey sharp wey person don use before.
- From mama wey get belle to pikin wey dey her belle when she dey born or wen she dey give am breast milk.

If person wey get HIV dey take ART (Antiretroviral therapy) medicine well, e go dey kill HIV for body, and e go hard to give another person HIV.

3. Why I suppose test my self?

Person wey don expose to HIV suppose test. de type of test-kit wey dem dey use pass dey see antibodies wey the person produce as part of their immune response to take fight HIV. Plenty times, person body go produce antibodies to HIV within 28 days of infection. For this time person fit see negative wey no be true negative result if dem use antibody test. Person wey dey at risk suppose test sharp sharp, if na negative test result inside 28 days of exposure person suppose still do another test after 3 month to sure of de result.

4. Wetin I suppose do if my result na positive?

Yu go need go hospital make professional person confam am for yu. If dem don confam am, Person for hospital go tell you wetin yu suppose do next.

Product information

CATALOG NO.

W006P0058 W006P0059 W006P0060

INTENDED USE

The Wondfo HIV self-test na to test yourself for only one time use for invitro diagnostics for fingerstick whole blood detection of HIV-1/2. Dem suppose use am as self-test and/or by medical professionals.

For *in vitro* diagnostic use only.

PRINCIPLES OF THE PROCEDURE

Wondfo HIV Self-Test dey use double antigen sandwich immunochromatography method. Once yu put blood specimen and the buffer for their respective wells, dem go move along the device by capillary action. The HIV-1 and/or HIV-2 antibodies go join to the colloidal gold-HIV antigen (gp36/41), and the complex is then captured by the HIV recombinant antigen (gp41 and gp36) immobilized in the test region (T). When the levels of HIV-1 or HIV-2 antibodies are at or above de limit of detection (LOD) of de assay, e produce a visible colored line in the test region (T) and show a positive result. When the levels of HIV-1 or HIV-2 antibodies are zero or below de LOD, no visible colored line go show for the test region (T) showing a negative result. To show say de test kit dey work as an internal procedure control, a colored line go show for de control region (C).

WARNINGS AND CAUTIONS

No use de test kit if expiration date don pass.
No use am if e don burst, tear or e dey open. No test your self if you neva reach 12years.

No use am if you get bleeding problem.
No use am if you already know say u get HIV.
No open the test things if yu neva ready for de test.

Wash your hand and clean am dry before you start the test. The place wey you go do the test must get beta light if you want better result.

Things way they inside

E get 3 configurations of the test kits, 1 test/kit, 20 tests/kit and 100 tests/kit. The things way they inside they listed for down.

Components	Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)	1	1x20	1x100
	Desiccant (pcs)	1	1x20	1x100
Accessories	Dropper (pcs)	1	1x20	1x100
	Buffer (vial)	1	1x20	1x100
	IFU (pcs)	1	1x20	1x100
	Blood Lancet for Single Use (pcs)	1	1x20	1x100
	Alcohol prep pad (pcs)	1	1x20	1x100
	Cotton swab (pcs)	1	1x20	1x100
	Disposal bag (pcs)	1	/	/

One test strip get: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

How to keep am make e no spoil

- Keep de test kit for store at 2-30 °C for 24 months.
- Use de test cassette for inside 1 hour after yu open de packet.
- No keep am for where sunlight, moist and heat dey.
- Use de kit at 10-30 °C.

LIMITATIONS OF THE PROCEDURE

- de test na for detecting human fingerstick whole blood.
- de test na for qualitative detection of HIV-1 and HIV-2 antibodies.
- de test assay procedure and result interpretation must be followed well when testing. For beta test performance, yu must collect blood specimen well well. if yu no follow how to do am well well, yu no go get beta result.
- Lie lie results fit happen if:
 - Patients exposed to HIV less than 3 months.
 - Patients dey collect HIV medicine (Antiretroviral therapy).
 - If de number of antibodies of HIV wey dey de blood specimen no reach de amount wey de test go fit see.
- Lie lie-positive results fit happen if:
 - Patients don do HIV vaccine clinical trial.
- de colored line wey show for "C" only means say dem add another thing inside the test. It does not guarantee that:
 - de correct blood specimen has been used.
 - de blood specimen has been applied correctly.
- As all diagnostic tests be, a correct test no be only for result wey yu do one time, but professional healthcare provider also with wetin happen for hospital and de results from other laboratory tests and checkings. No be only results from de Wondfo HIV Self-Test you go use finalize say you don get HIV.

PERFORMANCE CHARACTERISTICS

For de clinical study, dem been give 900 participants wey no know dem HIV status Wondfo HIV Self-Test to test. Dem come compare de results to de 4th generation laboratory test. The laboratory results show say 77 participants be HIV positive, 822 participants be HIV negative and 1 undetermined. 43 of de participants (5 HIV positive, 37 HIV negative and 1 undetermined) been no follow for de performance analysis, de comparison of results be:

- 95.8% of participants (69/72) been correctly report de result as positive, which mean say 3 participants wey get HIV reported negative result. This is called a false negative.
- 99.6% of participants (782/785) been correctly report de result as negative, which mean say 3 participants wey no get HIV reported positive result. This is called a false positive.
- 4.7% of participants (42/900) no get result, 1 participant's HIV infection status no dey confam during de clinical study so e no dey inside de analysis.

REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVds for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or selftesting, Geneva: World Health Organization; 2016.

SYMBOLS KEY

IVD	Consult instructions for use	Use-by date
Contains sufficient for <n> tests	Date of manufacture	Keep dry
LOT	Temperature limit	Keep away from sunlight
Manufacturer	Do not re-use	Catalogue number
Caution	No use am if e package dey spoil, and check e instructions for how to use am.	

Manufacturer information

Guangzhou Wondfo Biotech Co., Ltd.

Add: No.8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R. China

Tel: +86-20-32296083 400-888-5268 (Toll Free)


Fax: +86-20-32296063

E-mail: global@wondfo.com.cn

Website: en.wondfo.com

Please contact the manufacturer or your local distributor if you have any questions related to the product.

How to use the test kit (for fingerstick whole blood use)

How to read the test

Positive (+)

if two color lines show for ("T" and "C")

Go clinic make dem confam de test **sharp sharp**.

Negative (-)

If color line no show for ("T")

Test again in **3 months**, if yu don dey expose to HIV or if yu get high risk of infection. Ask dem for hospital how yu fit protect your self.

Invalid

If color bandline no show for ("C")

Run another test or go clinic make dem confam de test.

DISPOSE



Make yu put all de things wey yu use do test inside dustbin or troway dem how dem dey troway local infectious waste or dustbin wey get blood and other things how government talk am.

Wondfo HIV Self-Test


Instruction for use

- You must follow the test procedure carefully to get an accurate result.
- Sit in a clean, well-lit area and ensure you have all contents before beginning the test.
- For single-use only. Do not open foil pouch containing cassette until ready to test.

Contents


Preparation


Watch the operation video

How to use the test kit (for fingerstick whole blood use)

How to read the test

Questions and answers

1. What is HIV?

The human immunodeficiency virus (HIV) targets cells of the immune system, called CD4 cells, which help the body respond to infection. Within the CD4 cell, HIV replicates and in turn, damages and destroys the cell. Without effective treatment of a combination of antiretroviral (ARV) drugs, the immune system will become weakened to the point that it can no longer fight infection and disease.

2. How is HIV transmitted?

HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:

- Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
- Blood transfusion of contaminated blood.
- Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments.
- From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.

If a person living with HIV is on ART (Antiretroviral therapy), which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

3. When do I need to test myself?

People who have been exposed to HIV or are at risk of HIV should seek testing. Most widely-used HIV diagnostic tests detect antibodies produced by the person as part of their immune response to fight HIV. In most cases, people develop antibodies to HIV within 28 days of infection. During this time, people may receive a false-negative result using an antibody test. While people at risk of HIV should test as soon as possible, negative test results within 28 days of exposure should be confirmed with additional testing no more than three months later.

4. What should I do if I get a positive result?

You need to follow up with a health care worker to get additional testing to confirm the result. At that time your local clinic and the provider will suggest the next steps that need to be taken.

Product information

CATALOG NO.

W006P0058 W006P0059 W006P0060

INTENDED USE

The Wondfo HIV Self-Test is a single-use *in vitro* diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as self-test and/or by medical professionals.

For *in vitro* diagnostic use only.

PRINCIPLES OF THE PROCEDURE

Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into the respective wells, they will migrate along the device by capillary action. The HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen (gp36/41), and the complex is then captured by the HIV recombinant antigen (gp41 and gp36) immobilized in the test region (T). When the levels of HIV-1 or HIV-2 antibodies are at or above the limit of detection (LOD) of the assay, it will produce a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV-1 or HIV-2 antibodies are zero or below the LOD, there is no visible colored band in the test region (T) indicating a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C).

WARNINGS AND CAUTIONS

- Do not** use if the test kit beyond expiration date.
 - Do not** use if the pouch is punctured or improperly sealed.
 - Do not** use for self-testing if you are under 12 years old.
 - Do not** use for self-testing if you have a bleeding disorder.
 - Do not** use for self-testing if you are already diagnosed as HIV positive.
 - Do not** open the pouch until you are ready to perform the test.
- Wash your hands and ensure that they are clean and dry before starting test. Adequate lighting is required to read the test results.

KIT CONTENT

There are 3 configurations of the test kits, 1 test/kit, 20 tests/kit and 100 tests/kit. The kit components are provided as below:

Components		Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)		1	1x20	1x100
	Desiccant (pcs)		1	1x20	1x100
Accessories	Dropper (pcs)		1	1x20	1x100
	Buffer (vial)		1	1x20	1x100
	IFU (pcs)		1	1x20	1x100
	Blood Lancet for Single Use (pcs)		1	1x20	1x100
	Alcohol prep pad (pcs)		1	1x20	1x100
	Cotton swab (pcs)		1	1x20	1x100
	Disposal bag (pcs)		1	/	/

One test strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

STORAGE AND STABILITY

- The test kit can be stored at 2-30 °C for 24 months.
- Use the test cassette within 1 hour after opening the pouch.
- Keep away from sunlight, moisture and heat.
- Use the kit at 10-30 °C.

LIMITATIONS OF THE PROCEDURE

- The test is designed for detecting human fingerstick whole blood.
- The test is limited to the qualitative detection of HIV-1 and HIV-2 antibodies.
- The assay procedure and result interpretation must be followed closely when testing. For optimal test performance, proper specimen collection is critical. Failure to follow the procedure may lead to inaccurate test results.
- False-negative results can occur in the following conditions:
 - Patients exposed to HIV less than 3 months.
 - Patients under HIV treatment (Antiretroviral therapy).
 - If the quantity of antibodies for HIV present in the specimen is below the detection limit of the assay.
- False-positive results can occur in the following conditions:
 - Patients have participated in a HIV vaccine clinical trial.
- The presence of the control line only means that migration of added liquid occurred. It does not guarantee that:
 - The correct specimen has been used.
 - The specimen has been applied correctly.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but instead should be determined by a healthcare provider in conjunction with clinical findings and the results from other laboratory tests and evaluations. Results from the Wondfo HIV Self-Test should not be used as the sole basis for diagnosis.

PERFORMANCE CHARACTERISTICS

In the clinical study, 900 participants whose HIV status were unknown were given the Wondfo HIV Self-Test to test. The results were compared to the 4th generation laboratory test. The laboratory results shown that a total of 77 participants were HIV positive, 822 participants were HIV negative and 1 undetermined. A total of 43 participants (5 HIV positive, 37 HIV negative and 1 undetermined) were excluded from the performance analysis. The comparison of results was as follows:

- 95.8% of participants (69/72) correctly reported the result as positive. This means that 3 participants infected with HIV reported negative result. This is called a false negative.
- 99.6% of participants (782/785) correctly reported the result as negative. This means that 3 participants not infected with HIV reported positive result. This is called a false positive.
- 4.7% of participants (42/900) failed to obtain a result, 1 participant's HIV infection status was not confirmed during the clinical study so it was excluded from the analysis.

REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or selftesting, Geneva: World Health Organization; 2016.

SYMBOLS KEY

	<i>In vitro</i> diagnostic medical device		Consult instructions for use		Use-by date
	Contains sufficient for <n> tests		Date of manufacture		Keep dry
	Batch code		Temperature limit		Keep away from sunlight
	Manufacturer		Do not re-use		Catalogue number
	Caution		Do not use if package is damaged and consult instructions for use		

Manufacturer information

Guangzhou Wondfo Biotech Co., Ltd.

Add: No.8 Lizhishan Road, Science City, Huangpu District, 510663,Guangzhou, P.R. China

Tel: +86-20-32296083 400-888-5268 (Toll Free)

Fax: +86-20-32296063

E-mail: global@wondfo.com.cn

Website: en.wondfo.com

Please contact the manufacturer or your local distributor if you have any questions related to the product.


Wondfo

HIV Self-Test


Malangizo akugwilitsira ntchito

- Muyeñera kutsatira ndondomeko yoyeserera mosamala kuti mupeze zotsatira zolondola
- Khalani pamalo oyera, owala bwino ndipo onetsetsani kuti muli ndi zones musanayambe mayeso.
- Kugwilitsa ntchito kamodzi kokha. Osatsegula thumba la zojambuluzo lomwe lili ndi makaseti mpaka mutakonzeka kuyesa.


Contenu

Momwe mungawilitsira ntchito zida zoyesera(kugwilitsira ntchito magari pachala)

Momwe mungawerengere zoyesera

Mafunso ndi mayankho

1. Kodi HIV ndichiyani?

Kachilombo ka HIV (Human immunodeficiency virus) ndi kachirombo ka HIV komwe kamalimbana ndi maselo a chitetezo cha mthupi, otchedwa CD4 cell, omwe amathandiza thupi kuyankha matenda. Mkati mwa CD4 cell, kachilombo ka HIV kamabwerezabwerezera, kenaka, kumawononga maselo. Popanda mankhwala ophatikizika a ma ARV, chitetezo chamthupi chimafooka kwambiri moti sichingathenso kulimbana ndi matenda.

2. Kodi HIV imafalikira bwani?

Kachilombo ka HIV kamapezeka m'madzi ena am'thupi mwa anthu omwe ali ndi kachilombo ka HIV, kuphatikiza magari, umuna, ukazi, madzi am'mimba ndi mkaka wa m'mawere. HIV imatha kufalikira ndi:

- Kugonana kosadziteza kumaliseche, kapena kumatako, komanso nthawi zina kudzera muku-gonana m' kamwa ndi munthu yemwe ali ndi kachilombo ka HIV.
- Kuthilidwa mwazi wamagazi opitsidwa.
- Kugawana singano, majakisoni, zida zina zobaya, zida zopangira opaleshoni kapena zakuthwa
- Kuchokera kwa amayi emwe ali ndi HIV kupita kwa mwana wake wakhandu panthawi yoyembekezera, yobereka kapena yoyamwitsa.

Ngati munthu yemwe ali ndi kachilombo ka HIV akumwa mankhwala a ART (Antiretroviral therapy), omwe amapondereza kwambiri kachilombo ka HIV m' thupi, mwayi wake wopatsira wina munthu umachepea kwambiri.

3. Ndiyenera kudziyesa lithi?

Anthu omwe ali pachiwopsezo chotenga kachilombo ka HIV kapena omwe ali pachiwopsezo chotenga kachilombo ka HIV ayenera kukayezetsa. Kachilombo ka HIV komwe amagwiritsidwa ntchito kwambiri amapeza ma antibodies opangidwa ndi munthu ngati gawo la chitetezo chawo cholimbana ndi HIV. Nthawi zambiri, anthu amapanga ma antibodies ku HIV pasanathe masiku 28 atadwala. Panthawi imeneyi, anthu amatha kulandira zotsatira zabadza pogwiritsa ntchito antibody test. Ngakhale kuti anthu omwe ali pachiwopsezo chotenga kachilombo ka HIV akuyenera kuyezetsa msanga, zotsatira zake zosonyeza kuti alibe HIV pasanathe masiku 28 atadzizikira ziyenera kutsimikizidwa ndi kuyezetsa kwina pasanathe miyezi itatu.

4. N' tchite chiyani nditapeza zotsatira zovomeleza?

Muyeñera kutsata wazachipatala kuti mukayezetse zina kuti mutsimikizire zotsatira zake. Panthawiyi chipatala cha kwanuko ndi wopereka chithandizo adzakuuzani njira zina zomwe muyenera kuchita.

Mau okhudza za chinthu

CATALOG NO.

W006P0058 W006P0059 W006P0060

ZOFUNIKA KUGWILITSA NTCHITO

Wondfo HIV Self-Test ndi njira imodzi yokha yodzoyesera yokha mu vitro yodzoyesera kuti muone magari athunthu a HIV-1/2 ndi chala. Amapangidwa kuti azigwiritsidwa ntchito ngati kudziyesa komanso/kapena ndi akatswiri zachipatala Kuyang' anira matenda ndi njira ya vitro kokha

MFUNDO YANJIRA YAMACHITIDWE

Wondfo HIV Self-Test amatengera njira ya ma antigen sandwich immunochromatography. Zitsanzo ndi buffer zikawonjezeredwa m'zitsime zomwezo, zimasantha motsatira chipangizocho ndi capillary action. Ma antibodies a HIV-1 ndi/kapena HIV-2 amamanga ku antigen ya colloidal gold-HIV (gp36/41), ndipo zovutazo zimagwidwa ndi kachilombo ka HIV (gp41 ndi gp36) osasunthika m'chigawo choyesera (T). Pamene milingo ya ma antibodies a HIV-1 kapena HIV-2 ali pamwamba kapena pamwamba pa malire a kudziwika (LOD) ya kuyesa, idzatulutsa gulu lowoneka lachikuda m'chigawo choyesera (T) ndikuwonetsa zotsatira zabwino. Pamene milingo ya ma antibodies a HIV-1 kapena HIV-2 ndi ziro kapena pansu pa LOD, palibe gulu lowoneka bwino m'dera loyetzetsa (T) lomwe likuwonetsa zotsatira zoyipa. Kutu ikhale yoyang'anira ndondomeko yamkati, mzere wachikuda udzawonekera pagawo lolamulira (C).

CHENJEZO NDI MALANGIZO

- **Osagwiritsa ntchito** ngati zida zoyeserera zitadutsa tsiku lotha ntchito.
 - **Osagwiritsa ntchito** ngati thumba laboola kapena losindikizidwa molakwika.
 - **Osagwiritsa ntchito** kudziyesa ngati muli ndi zaka zosakwana 12.
 - **Osagwiritsa ntchito** kudziyesa nokha ngati muli ndi vuto lotaya magari.
 - **Osadziyenza ngati** mwapezeka kale ndi kachilombo ka HIV.
 - **Osatsegula thumbalo** mpaka mutakonzeka kuyesa.
- Sambani m'manja ndipo onetsetsani kuti ali aukhondo ndi owuma musanayambe kuyezetsa. Kuunikira kokwanira kumafunika kuti muwerenge zotsatira za mayeso

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Pali masanjidwe atatu a zida zoyesera, mayeso amodzi / zida, mayeso 20 / zida ndi 100 mayeso / zida. Zigawo za kit zimaperkedwa monga pansipa:

Zigawo	No. Ya Ndandanda	W006P0058	W006P0059	W006P0060
Thumba la makaseti oyosera	Kaseti woyesera (pcs)	1	1x20	1x100
	Desiccant (pcs)	1	1x20	1x100
Zida	Chogwetsera (pcs)	1	1x20	1x100
	Mankhwara (vial)	1	1x20	1x100
	IFU (pcs)	1	1x20	1x100
	Chotengera magari chogwira ntchito kamodzi (pcs)	1	1x20	1x100
	Pad yoyesera yajang' ala (pcs)	1	1x20	1x100
	Thonje lopukutira (pcs)	1	1x20	1x100
	Disposal bag (pcs)	1	/	/

Mzere umodzi woyesera umaphatikizapo: Mankhwara a Gold conjugate (HIV gp41/gp36 ophatikizidwa ndi antigen ndiponso -gold colloid), Mayeso a(HIV gp41 ndi HIV gp36 yophatikizidwa) ndi Control line (Goat anti rabbit IgG polyclonal antibody).

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MAUMBONI

1. WHO. TGS-5 Designing instruction for use for in vitro diagnostic medical devices, Geneva: World Health Organization ; 2019.
2. Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
3. WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or self-testing, Geneva: World Health Organization; 2016.

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








Questions and answers

1. What is HIV?

The human immunodeficiency virus (HIV) targets cells of the immune system, called CD4 cells, which help the body respond to infection. Within the CD4 cell, HIV replicates and in turn, damages and destroys the cell. Without effective treatment of a combination of antiretroviral (ARV) drugs, the immune system will become weakened to the point that it can no longer fight infection and disease.

2. How is HIV transmitted?

- HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:
 - Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
 - Blood transfusion of contaminated blood.
 - Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments with HIV.
 - From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.
- If a person living with HIV is on ART (Antiretroviral therapy), which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

3. When do I need to test myself?

People who have been exposed to HIV or are at risk of HIV should seek testing. Most widely-used HIV diagnostic tests detect antibodies produced by the person as part of their immune response to fight HIV. In most cases, people develop antibodies to HIV within 28 days of infection. During this time, people may receive a false-negative result using an antibody test. While people at risk of HIV should test as soon as possible, negative test results within 28 days of exposure should be confirmed with additional testing no more than three months later.

4. What should I do if I get a positive result?

You need to follow up with a health care worker to get additional testing to confirm the result. At that time your local clinic and the provider will suggest the next steps that need to be taken.

Product information

CATALOG NO. W006P0058 W006P0059 W006P0060

INTENDED USE

The Wondfo HIV Self-Test is a single-use *in vitro* diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used as self-test and/or by medical professionals. For *in vitro* diagnostic use only.

PRINCIPLES OF THE PROCEDURE

Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into the respective wells, they will migrate along the device by capillary action. The HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen (gp36/41), and the complex is then captured by the HIV recombinant antigen (gp41 and gp36) immobilized in the test region (T). When the levels of HIV-1 or HIV-2 antibodies are at or above the limit of detection (LOD) of the assay, it will produce a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV-1 or HIV-2 antibodies are zero or below the LOD, there is no visible colored band in the test region (T) indicating a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C).

WARNINGS AND CAUTIONS

- Do not use if the test kit beyond expiration date.
 - Do not use if the pouch is punctured or improperly sealed.
 - Do not use for self-testing if you are under 12 years old.
 - Do not use for self-testing if you have a bleeding disorder.
 - Do not use for self-testing if you are already diagnosed as HIV positive.
 - Do not open the pouch until you are ready to perform the test.
- Wash your hands and ensure that they are clean and dry before starting test. Adequate lighting is required to read the test results.

KIT CONTENT

There are 3 configurations of the test kits, 1 test/kit, 20 tests/kit and 100 tests/kit. The kit components are provided as below:

Components	Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)	1	1x20	1x100
	Desiccant (pcs)	1	1x20	1x100
	Dropper (pcs)	1	1x20	1x100
	Buffer (vial)	1	1x20	1x100
	IFU (pcs)	1	1x20	1x100
Accessories	Blood Lancet for Single Use (pcs)	1	1x20	1x100
	Alcohol prep pad (pcs)	1	1x20	1x100
	Cotton swab (pcs)	1	1x20	1x100
	Disposal bag (pcs)	1	/	/

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One strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

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
REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices. Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDS for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or self-testing. Geneva: World Health Organization; 2016.

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
RÉFÉRENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices. Geneva: World Health Organization ; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDS for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or self-testing. Geneva: World Health Organization; 2016.

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


2. How is HIV transmitted?

- HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:
 - Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
 - Blood transfusion of contaminated blood.
 - Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments.
 - From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.
- If a person living with HIV is on ART (Antiretroviral therapy), which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

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

How to read the test

DISPOSE

Throw away all contents in the normal trash or dispose according to local infectious waste disposal laws and regulations.

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









BIWANDIKO EBYAKOZESIBWA

- WHO. TGS-5 Designing instruction for use for in vitro diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVds for self-testing. London, UK: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or self-testing. Geneva: World Health Organization; 2016.

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EBIHANDIKO EBIKORESIBWE


- WHO. TGS-5 Designing instruction for use for in vitro diagnostic medical devices. Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVds for self-testing. London, UK: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or self-testing. Geneva: World Health Organization; 2016.

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
2. How is HIV transmitted?

- HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:
 - Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
 - Blood transfusion of contaminated blood.
 - Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments.
 - From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.
- If a person living with HIV is on ART (Antiretroviral therapy), which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

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WARNINGS AND CAUTIONS


- Do not use if the test kit beyond expiration date.
- Do not use if the pouch is punctured or improperly sealed.
- Do not use for self-testing if you are under 12 years old.
- Do not use for self-testing if you have a bleeding disorder.
- Do not use for self-testing if you are already diagnosed as HIV positive.
- Do not open the pouch until you are ready to perform the test.
- Wash your hands and ensure that they are clean and dry before starting test.
- Adequate lighting is required to read the test results.

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
Components	Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)	1	1x20	1x100
	Desiccant (pcs)	1	1x20	1x100
	Dropper (pcs)	1	1x20	1x100
	Buffer (vial)	1	1x20	1x100
	IFU (pcs)	1	1x20	1x100
	Blood Lancet for Single Use (pcs)	1	1x20	1x100
Accessories	Alcohol prep pad (pcs)	1	1x20	1x100
	Cotton swab (pcs)	1	1x20	1x100
	Disposal bag (pcs)	1	/	/

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



REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing, London, UK: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or self-testing, Geneva: World Health Organization; 2016.

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JERERO

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing, London, UK: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or self-testing, Geneva: World Health Organization; 2016.

MIFANANIDZO NEZVAINOREVA

Ruzivo maererano nevagadziri vemuziyo uyu

Guangzhou Wondfo Biotech Co., Ltd.
 WEDZERA: 8 rue Lizhishan, Cité des Sciences, Arrondissement Huangpu, 510663, Guangzhou, République populaire de Chine
 Nhaba dzenharo: +86-20-32296083 400-888-5268 (appel gratuit)
 Fakisai: +86-20-32296063
 Imeri: global@wondfo.com.cn
 Webusaiti: en.wondfo.com

Kana mune mundo chero iyi zvayo maererano nemuziyo uyu, zvisai vagadziri vavo kana kuti vanoutengesha munyika menyika.