

**E008: Auto-disable syringe for fixed dose immunization****PQS code:** E008/050**Description:** single-use auto-disable needle-free jet injector**Manufacturer's reference:** Pharmajet Stratis®**Manufactured in:** United States**Company:** PharmaJet, Inc.**Address:** 400 Corporate Circle,  
Suite N  
Golden, Colorado 80401  
USA**Telephone:** +1 303 526-4278**Email:** [Chris.Cappello@pharmajet.com](mailto:Chris.Cappello@pharmajet.com)**Web address:** <http://www.pharmajet.com/>**Specifications**

| Reset station  | Injector                                 | Needle free syringe                                     | Filling adapter   |
|--|--|---|---|
| Reference manufacturer C501-20   | Reference manufacturer: C501-10          | Reference manufacturer: C501-34                         | Reference manufacturer: C501-43 & C501-46               |
| Activation: manual   | Administration: IM or SC                 | Nominal capacity and graduation: 0.5ml                  | Material of adaptor: Polycarbonate                      |
| Force applied:   | Energy source: spring                    | Main materials of syringe: Polycarbonate                | Compatibility only with the Stratis needle-free syringe |
| Life time cycles: 20,000 min   | Main material: ABS, stainless steel      | Prevented from reuse by: Plunger break before injection | Intended for one vial use for field-filling of syringes |
| Carrying case IP: 55   | Inadvertent actuation prevention feature | Primary packaging: paper blister packs                  | Primary packaging: paper blister packs                  |
| Warranty 2 years   | Warranty 2 years                         | Number of syringes per box: 100 units                   | Number of adapters per box: 75 units                    |
| The PharmaJet Stratis™ needle-free injection system is intended to deliver various medications and vaccines either intramuscularly or subcutaneously by means of narrow precise fluid stream, which penetrates body tissues. |  |   |   |
| Quality system Standard  | ISO 13485,                               | Product standard:                                       | ISO 7886-3; FDA 510(k); CE 556625                       |

**Comments:**

Single-use auto-disable needle-free jet injector

Note: Jet Injectors (single-use auto-disable needle-free syringe injectors) as medical devices (Class II) are regulated by National Regulatory Authorities. A Jet injector prequalified by the WHO PQS programme to deliver medications and vaccines is to be used according to label to deliver ONLY those medications and vaccines that have been approved by the relevant Authorities.

**Current PQS status:** *pre-qualified ; 31 Jan 2013***Valid until :** *May 2024***Note:** If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased