

WHO Immunisation Devices - Performance, Quality & Safety



E013: Injection devices for therapeutic purposes

Product name: VanishPoint II® Retractable 5ml

IMD-PQS Code:E013-008Status:PrequalifiedDate of acceptance:16 Jun 2005Appliance type:E013

Product description: RUP Retractable 5ml

Manufacturer: Retractable Technologies, Inc.

Manufacturers Reference: VanishPoint II® Retractable 5ml

Country of Manufacture: China

Address: 511 Lobo Lane, Little Elm, TX 75068

United States of America

Telephone: +1 972 294 1010

Email: kduesman@retractable.com

Website address: www.retractable.com

Valid until: 31 May 2026

Specifications

PRODUCT SPECIFICATIONS - MAIN

Intended purpose Medication variable dose Restrictive use for None

administration

Vaccine capacity (mL) 5.00 Graduations (ml) 0.20

Needle Size Other Needle Size (Other) 21Gx1.5" (38mm)

Needle FixationFixedNominal capacity (ml)5.00Packaging & SterilityYesOther needle size optionsYes

Assurance

Plunger aspirationYesNumber of piecesRUP mechanismRetraction of needle into barrelRUP Type

by spring

Location RUP activationEnd of injectionRUP Needle FixationFixedRUP Needle Sizes21G x 1 1/2" (38 mm)RUP Nominal Capacity (ml)5.00

RUP plunger Aspiration Multiple RUP Primary Packaging Paper blister pack

RUP Sub-category Retractable RUP

Comments The retraction of the needle is considered as a reuse prevention feature (RUP) and a Sharp Injury

Protection (SIP)

PRICE AND SHIPPING

Shipping volume (m3)0.045Shipping weight (kg)8.20Pieces per package600Minimum order600.00Price 1On requestYear of base price2005

Incoterms FCA

QUALITY STANDARD

Quality Standard ISO 13485 **Specification Reference** Applicable PQS specification:

ISO 7886-4 Applicable PQS

VP(s): ISO 7886-4

3.00

Other

PRODUCT SITES

Organization name Retractable Technologies, Inc.

Address 511 Lobo Lane, Little Elm, TX 75068

United States of America

Current PQS status:PrequalifiedValid until:31 May 2026

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.