Notes on the Design of Bioequivalence Study: Cabotegravir

Notes on the design of bioequivalence studies with products invited for submission to the WHO Prequalification Team – Medicines (PQT/MED) are issued to aid manufacturers with the development of their product dossier. Deviations from the approach suggested below can be considered acceptable if justified by sound scientific evidence.

The current notes should be read and followed in line with the general guidelines of submission of documentation for WHO prequalification. In particular, please consult the "Multisource (generic) pharmaceutical products: guidelines on registration requirements to establish interchangeability" in *Fifty-first Report of the WHO Expert Committee on Specifications for Pharmaceutical Preparations*. Geneva, World Health Organization, 2017. WHO Technical Report Series, No. 1003, Annex 6.

Below, additional specific guidance is provided on the invited immediate release products containing cabotegravir.

Pharmacokinetics of cabotegravir

Cabotegravir is rapidly absorbed following oral administration, with median Tmax at 3 hours post dose for the tablet formulation. Cabotegravir may be administered with or without food. Food increases the extent of absorption of cabotegravir. Bioavailability of cabotegravir is independent of meal content: high fat meals increased cabotegravir AUC_{0-inf} by 14% and increased C_{max} by 14% relative to fasting conditions. Cabotegravir has a mean terminal half-life of 41 h.

Cabotegravir IM injection exhibits absorption-limited (flip-flop) kinetics resulting from slow absorption from the gluteal muscle into the systemic circulation resulting in sustained plasma concentrations. Following a single intramuscular dose, plasma cabotegravir concentrations are detectable on the first day and gradually rise to reach maximum plasma concentration with a median Tmax of 7 days. Cabotegravir has been detected in plasma up to 52 weeks or longer after administration of a single injection. Plasma cabotegravir exposure increases in proportion or slightly less than in proportion to dose following single doses ranging from 100 to 800 mg. Cabotegravir mean apparent terminal phase half-life is absorption-rate limited and is estimated to be 5.6 to 11.5 weeks after a single dose IM injection. The significantly longer apparent half-life compared to oral reflects elimination from the injection site into the systemic circulation.

Guidance for the design of bioequivalence studies

Taking into account the pharmacokinetic properties of cabotegravir, the following guidance with regard to the study design should be taken into account.

<u>Study design</u>: A single-dose cross-over design is recommended for the oral tablet. A single dose with cross-over or parallel design may be employed for the intramuscular prolonged-release suspension for injection.

Dose: As the Eol only includes the strength 30 mg as tablet and 600 mg / 3 ml intramuscular prolonged-release suspension for injection, the bioequivalence study should be conducted with these strengths.

Fasted/fed: The bioequivalence study should be conducted in the fasting state as cabotegravir tablets can be taken irrespective of meals.



<u>Subjects</u>: Healthy adult subjects should be recruited. It is not necessary to include patients in the bioequivalence study.

Parent or metabolite data for assessment of bioequivalence: The parent drug is considered to best reflect the biopharmaceutical quality of the product. Therefore, bioequivalence should be based on the determination of the parent compound.

<u>Sample size</u>: Cabotegravir C_{max} seems to be moderately variable (intra-subject CV of 22% approx.) when administered orally. After intramuscular administration the expected inter-subject CV is around 55-60%. These data may facilitate the calculation of a sufficient sample size for a cross-over bioequivalence study.

Washout: Taking into account the elimination half-life of cabotegravir of 41 hours in healthy volunteers, a washout period of 2 -3 weeks is considered sufficient to prevent carry over in the bioequivalence study for the oral product. A wash-out period of 60 weeks would be necessary for the intramuscular injection; therefore, a parallel design is an acceptable alternative.

Blood sampling: The blood sampling in the bioequivalence study for the oral tablet should be intensive for the first hours after administration to properly characterize the C_{max} of cabotegravir. It is not necessary to take blood samples beyond 72 hours for the characterization of cabotegravir pharmacokinetics. For example, samples can be taken pre-dose and at 0.50, 1.00, 1.50, 2.00, 2.50, 3.00, 3.50, 4.00, 4.50, 5.00, 6.00, 8.00, 12.00, 16.00, 24.00, 36.00, 48.00 and 72 hours after dosing.

The blood sampling in the bioequivalence study for the intramuscular injection should be taken frequently around the 7th day and for at least 42 weeks. For example, samples can be taken pre-dose and at 4, 8, 16, 24, 48, 96, 120, 144, 168, 192 hours, 2, 3, 4, 6, 8, 12, 20, 28, 36, and 42 weeks after administration.

<u>Analytical considerations</u>: Information currently available to PQT/MED indicates that it is possible to measure cabotegravir in human plasma using LC-MS/MS analytical methodology. The bioanalytical method should be sufficiently sensitive to detect concentrations that are 5% of the C_{max} in most profiles of each formulation (test or comparator).

<u>Statistical considerations</u>: The data for cabotegravir should meet the following bioequivalence standards in a single-dose cross-over or parallel design study:

For the oral tablet:

- The 90% confidence interval of the relative mean AUC $_{0-72\,h}$ of the test to comparator product should be within 80.00 125.00%
- The 90% confidence interval of the relative mean C_{max} of the test to comparator product should be within 80.00 125.00%.

For the intramuscular prolonged-release suspension for injection:

- The 90% confidence interval of the relative mean AUC_{0-t} of the test to comparator product should be within 80.00 – 125.00%
- The 90% confidence interval of the relative mean AUC_{0-inf} of the test to comparator product should be within 80.00 125.00%
- The 90% confidence interval of the relative mean C_{max} of the test to comparator product should be within 80.00 125.00%.
- The 90% confidence intervals of partial AUCs (e.g., AUC_{0-8 weeks} and AUC_{8 weeks-t}) should also be submitted as supportive information.

Notes on the Design of Bioequivalence Study: Cabotegravir