

E003: Refrigerators and freezers

PQS code: E003/126

Vaccine Freezer Type of appliance:

Manufacturer's reference: HBD-86

China; People's Republic of Manufactured in: Qingdao Haier Biomedical Co., Ltd Company:

Address: No.280 Fengyuan Road

High-tech Zone Qingdao(266111)

China, People's Republic of

+86-532-88935593 Telephone:

Email: wenming.wang@haierbiomedical.com

Web address: http://www.haiermedical.com

Specifications

Hot 0°C Climate zone: Min rated ambient temp:

R600a Refrigerant: Electric Mains 220-240V 50/60Hz **Energy source:**

+43°C Appliance tested at: Ext dimensions (HxLxD) 788mm*717mm*872mm Performance at: +43°C Fuel and cycle type: Electric - compression

REFRIGERATOR **FREEZER** Freeze protection: Not tested

Vaccine storage capacity: 61 L

(Litres)

86 L Gross volume (Litres): Gross volume:

Holdover time (hours): 7 hrs 19 minutes Waterpack freezing capacity:

Energy consumption, stable 0.95

running (kWh/24 hours):

Waterpack storage capacity:

Energy consumption, cool **Energy consumption during**

down test (kWh/24 hours): freezing:

- "1. Display: Solar LED Temperature Display
- Low power consumption"
 Display: Solar LED Temperature Display
- 2. Low power consumption"

Comments:

Accessories:

1. Instructions 2. Key 3. Handle component 4. Basket 5. Alcohol thermometer

Spare parts (ref, price):

"1 No. Dry filter Part-0074180001 2 No. overload protector 0274300008B 1 No. Switch 0075050003

1 No. Thermo sensor-0274000357A
1 No. Thermostat with sensor 0274000352A
1 No. Compressor-0060705994A
1 No. Fan 0074000172

2 No. Indicator 0075050018

2 No. Displayer 0274000357 1 No. Sticker-0270501015B"

Shipping volume: 0.71cbm Shipping weight: 90kg Price year 2022 Incoterms EXW

Price \$1065 inclduing 30 DTR

Quality standard: - ISO 9001:2008 -Specification ref.

Verification report: WB-21-0193 Verification laboratory: **CHEARI**

Current PQS status: pre-qualified: 07 Jul 2022 Validity until: May, 2024

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.