

**E008: Auto-disable syringe for fixed dose immunization**

| | |
|----------------------------------|---|
| PQS code: | E008/093 |
| Description: | Auto-Disable Syringe 0.25ml |
| Manufacturer's reference: | SOL-GUARD Auto-Disable Syringe 0.25ml |
| Manufactured in: | China; People's Republic of |
| Company: | Sol-Millennium Medical Inc. |
| Address: | 315 Shawnee North Drive Suite 100 Suwanee GA USA |
| Telephone: | +41 76 675 65 86 |
| Email: | mbelluzzi@sol-m.com , ppasia@sol-m.com , yzhao@sol-m.com , yhong@sol-m.com |
| Web address: | www.sol-m.com |

Specifications

| | | | |
|------------------------------|--------------------|--------------------------------|-----------------------------------|
| Graduations: | 0.25ml | Quality standard: | ISO 13485, |
| Syringe material(s): | Polypropylene | Markings: | Device license (Canada) |
| Fixed needle size: | | Pieces per carton: | 2400 |
| AD mechanism: | metal clip | Volume per carton (m3): | 9.6 |
| AD location: | Start of injection | Weight per carton (kg): | 9.6 |
| Number of components: | 3 pieces | Minimum order: | 100,800 |
| Other needle options: | Yes | Incoterms: | FOB |
| Primary packaging: | Blister pack | Prices per unit: | <i>price available on request</i> |
| Year base price: | 2022 | | |

Comments: Graduation 0.25mL. Average deadspace 0.003mL
0.25ml G25x5/8"

Current PQS status: *pre-qualified ; 17 May 2022*

Valid until : *May 2024*

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased