



## E006: Temperature monitoring devices

<b>PQS code:</b>	<b>E006/042</b>
<b>Description:</b>	30-day electronic temperature logger
<b>Manufacturer's ref:</b>	HETL-01
<b>Manufactured in:</b>	China; People's Republic of
<b>Company:</b>	Qingdao Haier Biomedical Co., Ltd
<b>Address:</b>	No.280 Fengyuan Road High-tech Zone Qingdao(266111), China; People's Republic of

<b>Telephone:</b>	+86-532-88935593
<b>Email:</b>	<a href="mailto:wenming.wang@haierbiomedical.com">wenming.wang@haierbiomedical.com</a>
<b>Web address:</b>	<a href="http://www.haiermedical.com">http://www.haiermedical.com</a>

## Specifications

<b>Product type:</b>	30-day electronic refrigerator temperature logger	<b>Display medium:</b>	LCD
<b>Mode of operation:</b>		<b>Weight</b>	0.14 kg
<b>Minimum temperature:</b>	-20	<b>Main material:</b>	ABS
<b>Maximum temperature:</b>	+50	<b>Dimensions:</b>	13.1 x 2.4

### MECHANICAL AND ELECTRONIC DEVICES ONLY

<b>Accuracy:</b>	± 0.5°C	<b>Min logging interval:</b>	6 min
<b>Power source:</b>	Non-replaceable battery	<b>Alarm type:</b>	Visual
<b>Battery type:</b>	Non-rechargeable lithium battery	<b>User interface:</b>	device itself
<b>Battery shelf life:</b>	3 years	<b>Programmability:</b>	Factory programmed
<b>Activated life:</b>	2 years	<b>IP rating:</b>	IP65
<b>Memory capacity:</b>	8192 data points		

**Comments:** 30-day electronic refrigerator temperature logger, integrated sensor devices, reusable. It is used to monitor and log the vaccine and pharmacy storage environment, applied in vaccine refrigerators and cold rooms. This device can log the daily temperature data, and alarm indications displayed in the LCD display when the temperature exceeds. Data with 30 days can be viewed and downloaded by USB port to PC, and generate temperature charts

<b>Minimum order:</b>	5	<b>Year base price:</b>	2016
<b>Quality standard:</b>	- ISO 9001:2008 -	<b>Pieces per carton:</b>	5
<b>Price per unit:</b>	1-100 units: USD 23.00 101-200 units: USD 22.00; 201-500 units: USD 21.00 >501 units: USD 20.00	<b>Volume per carton:</b>	0.001m3
		<b>Weight per carton:</b>	0.14 kg
		<b>Incoterms:</b>	FCA
<b>Verification report:</b>	WB-16-243	<b>Verification laboratory:</b>	CHEARI
<b>Current PQS status:</b>	<i>pre-qualified ; 14 Jul 2016</i>	<b>Valid until:</b>	<i>May, 2024</i>

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.