



## E010: Waste management products

<b>PQS code:</b>	<b>E010/010</b>
<b>Description:</b>	Safety box 5L
<b>Manufacturer's reference:</b>	Sol-Millennium Safety Box 5L
<b>Manufactured in:</b>	China; People's Republic of
<b>Company:</b>	Shanghai Sol-Millennium Medical Products Co. Ltd.
<b>Address:</b>	Sol-Millennium Medical Products Co., Ltd. Unit 1801-1803, No. 2 Grand Gateway No. 3 Hongqiao Road Xuhui District Shanghai 200030 China
<b>Telephone:</b>	+86-21-68622586-17
<b>Email:</b>	<a href="mailto:ayu@sol-m.com">ayu@sol-m.com</a>
<b>Web address:</b>	<a href="http://www.sol-m.com">www.sol-m.com</a>

## Specifications

<b>Volume (Liter):</b>	5	<b>Syringe hole diameter (mm):</b>	38
<b>Nominal syringe capacity:</b>	150	<b>Material:</b>	cardboard
<b>Dimensions before assembly (mm):</b>	600 x 280 x 5.0	<b>Colour:</b>	Brown
<b>Dimensions after assembly (mm):</b>	288 x 165 x 130		
<b>Weight empty (g):</b>	270		
<b>Minimum order:</b>	5000	<b>Pieces per carton:</b>	25
<b>Year base price:</b>	2010	<b>Volume per carton:</b>	0.018 m3
<b>Price per unit:</b>	On request	<b>Weight per carton:</b>	8 kg
		<b>Incoterms:</b>	EXW
		<b>Quality standard:</b>	ISO 9001
<b>Verification report:</b>	110-29507 (August 2010)	<b>Verification laboratory:</b>	ForceTechnology
<b>Comments:</b>	Carton of 25 units		
<b>Current PQS status:</b>	<i>pre-qualified; 01 Nov 2010</i>	<b>Valid until :</b>	<i>May, 2024</i>

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased