



## E006: Temperature monitoring devices



<b>Product name:</b>	<b>VaxTag +</b>
<b>IMD-PQS Code:</b>	E006-099
<b>Status:</b>	Prequalified
<b>Date of acceptance:</b>	07 Oct 2024
<b>Product description:</b>	30 day temperature logger
<b>Manufacturer:</b>	LogTag North America Inc.
<b>Manufacturers Reference:</b>	VaxTag +
<b>Country of Manufacture:</b>	Vietnam
<b>Address:</b>	Unit 10, 10 Millpond Drive Lafayette New Jersey United States of America
<b>Telephone:</b>	+64 9 448 2311
<b>Email:</b>	<a href="mailto:quality@logtagrecorders.com">quality@logtagrecorders.com</a>
<b>Website address:</b>	<a href="http://www.logtagrecorders.com">www.logtagrecorders.com</a>
<b>Valid until:</b>	31 May 2025

## Specifications

### PRODUCT SPECIFICATIONS - MAIN

<b>Mode of operation</b>	Electronic	<b>Maximum logging intervals (mins)</b>	1.00
<b>External dimensions height (mm)</b>	10.70	<b>External dimensions length (mm)</b>	8.60
<b>External dimensions width (mm)</b>	54.50	<b>Material</b>	Polycarbonate
<b>Minimum working temperature (°C)</b>	-30.00	<b>Maximum working temperature (°C)</b>	65.00

### MECHANICAL AND ELECTRONIC DEVICES ONLY

<b>Accuracy (+/- °C)</b>	0.50	<b>Power source</b>	Non-replaceable battery
<b>Battery type</b>	Lithium Button Cell	<b>Battery shelf life (year)</b>	4
<b>Activated battery life</b>	24	<b>Alarm type</b>	Visual
<b>Memory capacity (readings)</b>	16129.00	<b>Programmability</b>	Factory Programmed
<b>User interface</b>	Push button operation	<b>Display medium</b>	Built-in LCD
<b>Ingress protection rating IP</b>	IP65		

<b>Comments</b>	Minimum logging interval: 30 seconds (fixed at 6 minutes for type 1) LogTag® vaxtag+ is an advanced temperature monitoring device specifically designed for the safe storage of vaccines in healthcare facilities and intermediate storage areas. Fully compliant with World Health Organization (WHO) standards under PQS specification E006/TR06.4, the vaxtag+ monitors storage conditions to safeguard vaccine integrity. The vaxtag+ stores up to 16,129 temperature readings to ensure the continuous monitoring for up to 67 days using the 6 minute logging interval. Default configuration is Type 1 (as per TR06.4); other configurations are also possible. Users can easily inspect violation history directly on the display or view charts and statistics, downloading the on-board PDF via the built-in USB connector. The data is compatible with the free LogTag Analyzer software, which allows users to visualize data in various formats, facilitating effective data management and archiving. The LCD display provides a clear view of current and historical temperature readings, including a 30-day summary of the highest and lowest temperatures. The vaxtag+ is also supported by LogTag Online, LogTag's 21CFRPart11 compliant data management system. Trust the reliability of the LogTag® vaxtag+ to maintain the safety and efficacy of your vaccination programs.
-----------------	--

### PRICE AND SHIPPING

<b>Shipping volume (m3)</b>	0.0051	<b>Units per package</b>	50.00
<b>Minimum order</b>	20.00	<b>Incoterms</b>	EXW
<b>Price 1</b>	20- 99 unit:40.80 USD	<b>Price 2</b>	> 100 unit:on request to manufacturer
<b>Year of base price</b>	2024		

**QUALITY STANDARD**

<b>Quality Standard</b>	Other	<b>Quality standard (Other)</b>	ISO 9001:2015
-------------------------	-------	---------------------------------	---------------

**VERIFICATION**

<b>Verification Laboratory</b>	0013X00004MI8OaQAJ	<b>Verification Report Reference</b>	A66604, Issue1
<b>Current PQS status:</b>	<i>Prequalified</i>	<b>Valid until:</b>	<i>31 May 2025</i>

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.