



E008: Injection devices for immunization



| | |
|---------------------------------|---|
| Product name: | Sol-Guard 0.5mL Auto Disable syringe |
| IMD-PQS Code: | E008-090 |
| Status: | Prequalified |
| Date of acceptance: | 20 Apr 2022 |
| Appliance type: | E008 |
| Product description: | 0.5mL Auto Disable syringe |
| Manufacturer: | Sol-Millennium Medical Inc. |
| Manufacturers Reference: | Sol-Guard 0.5mL Auto Disable syringe |
| Country of Manufacture: | China |
| Address: | 311 South Wacker Drive, Suite 4100 Chicago 60606 Illinois United States of America |
| Telephone: | +41 76 675 65 86 |
| Email: | jwang@sol-m.com |
| Website address: | www.sol-m.com |
| Valid until: | 31 May 2026 |

Specifications

PRODUCT SPECIFICATIONS - MAIN

| | | | |
|----------------------------------|---|--|-----------------------|
| Vaccine capacity (mL) | 0.50 | Graduations (ml) | 0.50 |
| Markings | Device license (Canada) | Syringe material(s) | Polypropylene |
| RUP mechanism | locking clip, locked trapped plunger | Packaging & Sterility Assurance | Yes |
| Other needle size options | Yes | AD Number pieces | 3 pieces |
| Auto-disabling at: | Start of injection | Needle Size | 23Gx1" (0.60 x 25 mm) |
| AD Primary packaging | Blister pack | | |

PRICE AND SHIPPING

| | | | |
|-----------------------------|----------------------------|--------------------------------|-----------|
| Shipping volume (m3) | 0.088 | Weight per package (Kg) | 9.6000 |
| Pieces per package | 2400 | Minimum order (unit) | 100800.00 |
| Price 1 | price available on request | Year of base price | 2022 |
| Incoterms | FOB | | |

QUALITY STANDARD

| | | | |
|-------------------------|-----------|--------------------------------|---|
| Quality Standard | ISO 13485 | Specification Reference | Applicable PQS specification: 7886-3 Applicable PQS VP(s): |
|-------------------------|-----------|--------------------------------|---|

PRODUCT SITES

| | |
|--------------------------|---|
| Organization name | Sol-Millennium Medical Inc. |
| Address | 311 South Wacker Drive, Suite 4100 Chicago 60606 Illinois United States of America |

| | | | |
|----------------------------|---------------------|---------------------|-------------|
| Current PQS status: | <i>Prequalified</i> | Valid until: | 31 May 2026 |
|----------------------------|---------------------|---------------------|-------------|

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.