



E010: Waste management equipment



Product name:	Sol-Millennium Safety Box 5L
IMD-PQS Code:	E010-010
Status:	Prequalified
Date of acceptance:	11 Jan 2010
Appliance type:	E010
Product description:	Safety box 5L
Manufacturer:	Shanghai Sol-Millennium Medical Products Co. Ltd.
Manufacturers Reference:	Sol-Millennium Safety Box 5L
Country of Manufacture:	China
Address:	Unit 2501-2506, No. 2 Grand Gateway No. 3 Hongqiao Road Xuhui District 20030 Shanghai China
Telephone:	+86-21-68622586-17
Email:	ichen@solm.com
Website address:	www.sol-m.com
Valid until:	31 May 2026

Specifications

PRODUCT SPECIFICATIONS - MAIN

Volume (Litre)	5.00	Capacity (# units)	150.00
Diameter syringe hole (mm)	38.00	Material	cardboard
Colour	Brown	Walls thickness (mm)	1.30
Weight fully assembled (grams)	270.00	Ext.Dimension After Assembly Height(mm)	130.00
Ext.Dimension After Assembly Length(mm)	288.00	Ext.Dimension After Assembly Width(mm)	165.00
Ext.Dimension Before Assembly Height(mm)	5.00	Ext.Dimension Before Assembly Length(mm)	600.00
Ext.Dimension Before Assembly Width(mm)	280.00		

PRICE AND SHIPPING

Shipping volume (m3)	0.018	Weight per Unit (Kg)	8.0000
Pieces per package	25	Minimum order	5000.00
Price 1	On request	Year of base price	2010
Incoterms	EXW		

QUALITY STANDARD

Quality Standard	Other	Quality standard (Other)	ISO 9001
Specification Reference	Applicable PQS specification: E010/SB01.1 Applicable PQS VP(s): E010/SB01-VP.1		

VERIFICATION

Verification Laboratory	ForceTechnology	Address	FORCE Technology Brøndby Park Allé 345 2605 Brøndby DenmarkDenmark
Telephone	+45 43250591	Email	njs@force.dk
Verification Report Reference	110-29507 (August 2010)		

PRODUCT SITES

Organization name	Shanghai Sol-Millennium Medical Products Co. Ltd.
Address	Unit 2501-2506, No. 2 Grand Gateway No. 3 Hongqiao Road Xuhui District 20030 Shanghai China

Current PQS status:

Prequalified

Valid until:

31 May 2026

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.