

WHO Immunisation Devices - Performance, Quality & Safety

E013: Injection devices for therapeutic purposes



Product name: Haiou RUP retractable safety syringe 3ml

IMD-PQS Code:E013-089Status:PrequalifiedDate of acceptance:06 May 2013

Appliance type: E013

Product description: RUP retractable syringe 3ml

Manufacturer:Guangdong Haiou Medical Apparatus Co., LtdManufacturers Reference:Haiou RUP retractable safety syringe 3ml

Address: Nanyuan Industrial Zone, North Liusha, Puning City,

Guangdong Province, China

VP(s): ISO 7886-4

31 May 2026

China

Telephone: +86-663-2758666

Email: market01@haiou.net.cn

Website address: http://www.haiou.net.cn

Valid until: 31 May 2026

Specifications

PRODUCT SITES
Organization name

Current PQS status:

Address

PRODUCT SPECIFICATIONS	- MAIN		
Intended purpose	Medication variable dose administration; Medication variable dose reconstitution	Restrictive use for	None
Graduations (ml)	0.10	Needle Size	Other
Needle Size (Other)	23G x 25mm	Needle Fixation	Fixed
Nominal capacity (ml)	3.00	Packaging & Sterility Assurance	Yes
Other needle size options	Yes	Plunger aspiration	Yes
Number of pieces	3.00	RUP mechanism	Plunger lock and needle retraction
RUP Type	Plunger lock	Location RUP activation	End of injection
RUP Needle Fixation	Fixed	RUP Needle Sizes	23G x 1" (25 mm)
RUP Nominal Capacity (ml)	3.00	RUP plunger Aspiration	Multiple
RUP Primary Packaging	Blister pack	RUP Sub-category	Retractable RUP
CE marking	Yes		
PRICE AND SHIPPING			
Shipping volume (m3)	0.118	Shipping weight (kg)	17.00
Pieces per package	2400	Minimum order	100000.00
Price 1	Upon request to manufacturer	Year of base price	2013
Incoterms	FOB		
QUALITY STANDARD			
Quality Standard	ISO 13485	Specification Reference	Applicable PQS specification: ISO 7886-4 Applicable PQS

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.

China

Prequalified

Guangdong Haiou Medical Apparatus Co., Ltd

Nanyuan Industrial Zone, North Liusha, Puning City, Guangdong Province, China

Valid until: