

WHO Immunisation Devices - Performance, Quality & Safety

E013: Injection devices for therapeutic purposes



Product name: Haiou RUP retractable syringe 10ml

IMD-PQS Code:E013-091Status:PrequalifiedDate of acceptance:06 May 2013

Appliance type: E013

Product description: RUP retractable syringe 10ml

Manufacturer: Guangdong Haiou Medical Apparatus Co., Ltd

Manufacturers Reference: Haiou RUP retractable syringe 10ml
Address: Nanyuan Industrial Zone, North Liusha,

Puning City Guangdong China

Telephone: +86-663-2758666

Email: market01@haiou.net.cn

Website address: http://www.haiou.net.cn

Valid until: 31 May 2026

Specifications

PRODUCT	SPECIFICATIONS	3 - MAIN
---------	----------------	----------

Intended purpose Medication variable dose Restrictive use for None administration: Medication

variable dose reconstitution

Graduations (ml)0.50Needle SizeOtherNeedle Size (Other)21G x 38mmNeedle FixationFixed

Nominal capacity (ml) 10.00 Packaging & Sterility Yes

Assurance

Other needle size options Yes Plunger aspiration Yes

Number of pieces 3.00 RUP mechanism Plunger lock and needle

retraction

RUP TypePlunger lockLocation RUP activationEnd of injectionRUP Needle FixationFixedRUP Needle Sizes21G x 1 1/2" (38 mm)

RUP Nominal Capacity (ml) 10.00 RUP plunger Aspiration Multiple

RUP Primary Packaging Blister pack RUP Sub-category Retractable RUP

CE marking Yes

PRICE AND SHIPPING

Shipping volume (m3)0.11Shipping weight (kg)17.00Pieces per package1200Minimum order100000.00Price 1Upon request to manufacturerYear of base price2013

Incoterms FOB

QUALITY STANDARD

Quality Standard ISO 13485 **Specification Reference** Applicable PQS specification:

ISO 7886-4 Applicable PQS VP(s): ISO 7886-4

VI (3). I

PRODUCT SITES

Organization name Guangdong Haiou Medical Apparatus Co., Ltd

Address Nanyuan Industrial Zone, North Liusha,

Puning City Guangdong China

Current PQS status: Prequalified Valid until: 31 May 2026

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.