

WHO Immunisation Devices - Performance, Quality & Safety



E013: Injection devices for therapeutic purposes

Product name: SHIFA syringe 10ml

IMD-PQS Code:E013-131Status:PrequalifiedDate of acceptance:20 Oct 2020Appliance type:E013

Product description: RUP syringe 10ml

Manufacturer: Changzhou Tongda Medical Appliance Co., Ltd

Manufacturers Reference: SHIFA syringe 10ml

Address: Sanhekou Street, Zhenglu Town, Tianning District,

Changzhou, Jiangsu, China. No. 208, Sanhe Road

Changzhou Jiangsu China

Telephone: +86 51988670522

Email: info@tongdamedical.com

Website address: http://www.tongdachina.com.cn

Valid until: 31 May 2026

Specifications

	TIONS - MAIN

Intended purpose Medication fixed dose Restrictive use for None

administration; Medication variable dose administration

Graduations (ml)0.50Needle SizeOtherNeedle Size (Other)22Gx1 1/4"Needle FixationDetached

Nominal capacity (ml) 10.00 Packaging & Sterility Yes

Assurance

Other needle size options Yes Plunger aspiration Yes

Number of pieces3.00RUP mechanismlocked and broken plunger

RUP TypePlunger lockLocation RUP activationEnd of injectionRUP Needle FixationNot FixedRUP Needle Sizes22G x 1 1/4" (30 mm)

RUP Nominal Capacity (ml) 10.00 RUP plunger Aspiration Multiple

RUP Primary Packaging Blister pack RUP Sub-category Simple RUP

CE marking Yes

Comments Needle at side available Pieces per carton:1600 or 1200

PRICE AND SHIPPING

Shipping volume (m3)0.139776Shipping weight (kg)18.94Pieces per package1600Minimum order200000.00Price 1price on requestYear of base price2020

Incoterms FOB

QUALITY STANDARD

Quality Standard ISO 13485 **Specification Reference** Applicable PQS specification:

ISO 7886-4 Applicable PQS VP(s): ISO 7886-4

PRODUCT SITES

Organization name Changzhou Tongda Medical Appliance Co., Ltd

Address Group 4, Ninghe village, Zhenglu Town, Tianning District

Changzhou Jiangsu China Current PQS status: Prequalified Valid until: 31 May 2026

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.