

Packaging information

BIMERVAX® emulsion single dose will be available for UN purchase in a secondary package for 10 and 20 vials.

Packaging configuration is presented in the table below:

	No. of SI	OV/carton	Dimensions (LxHxW)	Cold Chain Volume (cm³/dose)		Gross
	Box pallet va-	Box pallet va-	Box pallet va-Q-pal	Box pallet va-Q-one	Box pallet va-Q-pal	Box pallet va-Q-one	weight (kg)
	Q-pal	Q-one*					
Secondary Packaging	20		9.1 x 7.05 x 4.5 cm		14.43		n/a
Secondary Lackaging	1	.0	9.1 x 3.75 x	4.75 cm	16.2	20	n/a
Tertiary Packaging (20 SDV)	1 440	1 040	44.4 x 29.7 x 19.7 cm	47 x 46 x 38 cm	n/a	n/a	19.25
Tertiary Packaging (10 SDV)	1 320	920	45.4 x 30.0 x 21.7 cm	4 / X 40 X 38 CIII	n/a	n/a	18.59

Pallet shipper(s)	No. of shipping boxes/pallet	No. of doses/pallet	Pallet dimensions (LxHxW)	Gross weight (kg)
Box pallet va-Q-pal (20 SDV)	48	69 120	144.4 x 156.0 x 122.2 cm	692.08
Box pallet va-Q-pal (10 SDV)	40	52 800	144.4 X 130.0 X 122.2 CIII	584.96

Box pallet va-Q-pal: Single-use container with three layers. From the outside to the inside, there is the outer protective shell (honeycomb board), an insulation layer made of 40 mm thick vacuum insulation panels (VIPs) and an inner protective carboard leading to an outstanding temperature stability.

Box pallet va-Q-one: Single-use container with vacuum insulation panels (VIPs) and protective carboard leading to an outstanding temperature stability.

*va-Q-one: Is generally used for shipment of less units and may be shipped individually or over EUR pallets (1580 mm (Height) x 1200 mm (Width) x 1600 mm (Length) of 22 kg.